



## ***POWER TILT ACCESSORY***

**FUZE T50 | STELLAR GL | STELLAR LEAP | BENTLEY**

# **OWNER'S OPERATION AND MAINTENANCE MANUAL**

### **DEALER**

This manual must be given  
to the user of the  
wheelchair

### **USER**

Before using this  
wheelchair, read this  
entire manual and save for  
future reference

For more information regarding PDG products,  
parts & service, please visit  
[www.pdgmobility.com](http://www.pdgmobility.com)



Class 1 Medical Device

## 1. DEVICE DESCRIPTION

Thank you for purchasing a PDG mobility wheelchair Power Tilt Accessory.

Please do not install, maintain or operate this equipment without first reading and understanding this manual. Do not modify any parts on the unit in any way not specified by the manual. If you are unable to understand the warnings and instructions, contact a health care professional such as a doctor or therapist who is familiar with this type of product before attempting to use this equipment. Otherwise, injury or damage may result. Refer to this manual regularly for maximum safety and performance.

The Power Tilt is a wheelchair accessory that allows the attendant or user to tilt the wheelchair using a remote control switch. The accessory consists of an electric linear actuator, a battery assembly, a controller box, a tilt limiting assembly, and the switches used to control the system.

If you have any questions or comments about this manual, the safety and reliability of your wheelchair and the service you receive by us or your PDG supplier, please write or call us using the contact information below:

PDG PRODUCT DESIGN GROUP INC.  
#103- 318 East Kent Avenue South, Vancouver, BC V5X 4N6  
Ph: (604) 323-9220 Fax: (604) 323-9097  
E-mail: [info@pdgmobility.com](mailto:info@pdgmobility.com)

**SAVE THIS MANUAL FOR FUTURE REFERENCE**

## 2. TABLE OF CONTENTS

1.	Device Description.....	2
2.	Table of Contents .....	3
3.	Special Notes .....	4
4.	Safety Summary .....	4
5.	Power Tilt Schematic .....	5
6.	Power Tilt Installation .....	5
7.	Power Tilt Operation .....	6
8.	Switch Operation .....	6
9.	Installing and Removing Control Switches .....	7
10.	Adjusting Tilt Limit.....	8
11.	Using and Charging the Battery.....	9
12.	Installing and Removing the Battery .....	10
13.	Battery Charger .....	11
14.	Battery Charger Wiring Schematic .....	11
15.	Troubleshooting .....	12
16.	Servicing the Battery .....	14
17.	Disabling the Tilt Limit Switch .....	15
18.	Checking the Control Switches – 3 Position Switch .....	16
19.	Checking the Control Switches – Mono (Egg Switch).....	17
20.	Troubleshooting Flowchart – Power Tilt System .....	18
21.	Troubleshooting Flowchart –Battery.....	19
22.	Warranty Information .....	20
23.	Technical Specifications.....	21
24.	APPENDIX .....	21
25.	NOTES: .....	38

### 3. SPECIAL NOTES

**WARNING/CAUTION** — Notices as used in this manual apply to hazards or unsafe practices which could result in personal injury or property damage.

**NOTE** — THE INFORMATION CONTAINED IN THIS DOCUMENT IS SUBJECT TO CHANGE WITHOUT NOTICE. An updated version of this owner's manual may be available at [www.PDGMobility.com](http://www.PDGMobility.com)

**INDICATIONS FOR USE** — This product is intended to provide the benefits of tilt-in-space to a wheelchair user, wherever an attendant to manually tilt the chair is not available, or the user is unable to self-tilt.

### 4. SAFETY SUMMARY

**CAUTION:** Keep clear of lower frame and underside of upper frame when operating tilt mechanism


**WARNING:** Disconnect battery before adjusting or performing maintenance on wheelchair.

**CAUTION:** Do not open battery enclosure unless you are a competent technician

**CAUTION:** Do not step on battery enclosure

**CAUTION:** Do not immerse the battery in any liquid

## 5. POWER TILT SCHEMATIC



## 6. POWER TILT INSTALLATION

- See assembly drawings in the Appendix.
- Installation and service to the Power Tilt should be performed by a competent technician.


## 7. POWER TILT OPERATION

- **Do not** use the Power Tilt if it is behaving abnormally or erratically.
- **Do not** use the Power Tilt in wet conditions or immerse any of the components in liquid.
- **Do not** use the system with an occupant exceeding the maximum rated weight.
- **Inspect** the battery enclosure, connector, and cable before each use. Do not use the Power Tilt system if any of the components show wear or damage.
- **Disconnect** the power to the Power Tilt system prior to performing maintenance or adjustments on the wheelchair.
- **Caution** Keep clear of lower frame and underside of upper frame when operating tilt mechanism
- The Power Tilt system can be operated and adjusted by the wheelchair user or attendant.
- The Power Tilt operates an electric actuator to tilt the chair to any point up to the frame's tilt limit.
- The supplied actuator stops automatically at the maximum and minimum tilt. Maximum tilt can be further limited, see below.


## 8. SWITCH OPERATION

- The Power Tilt is controlled by digital input to the controller from the supplied switches – a toggle switch, a 2-button switch or an egg switch. To tilt the chair, operate the switches as shown. The Power Tilt moves while the switch is held and stops when the switch is released.
- The egg switch operates in Toggle Mode, where the tilt direction reverses each time the switch is activated.

Toggle switch (viewed from the side)




2 button



## 9. INSTALLING AND REMOVING CONTROL SWITCHES

- **Use only** PDG provided or PDG approved switches with the Power Tilt system.
- Up to two switches can be connected to the controller simultaneously.
- To install or remove a switch, pry off the end cap on the front of the controller enclosure to expose the switch jacks. Route the switch cable through the hole with a grommet in the side of the enclosure.
- **Do not** route the cable out the front of the enclosure.
- 3-position switches should be plugged directly into the controller. Mono switches (egg switch) must be used with a PDG-provided mono-to-stereo adapter cable (only supplied when required).




- Ensure the switch connector is fully seated into the jack. Replace the end cap. Please note the cap is designed to contact the switch connector.
- **Caution:** Partial insertion of the switch connector may result in an abnormal or a non-functioning system.

## 10. ADJUSTING TILT LIMIT

**Note:** Not available on Stellar LEAP or Bentley

- The Power Tilt system has an adjustable tilt limit setting and a tilt limit indicator scale.




### Tilt Limit Slider and Indicator (Stellar GL & Fuze T50)

- To adjust tilt limit, move the tilt limit slider to the desired degree of tilt, shown by the label.
  - Note: the tilt limit is adjusted relative to the lower frame, not relative to the ground.
  - Note: regardless of tilt limit position the power tilt will return to a level seat unless there is a system fault.

## 11. USING AND CHARGING THE BATTERY



- **Do not** use the Power Tilt with depleted batteries.
- **Do not** leave the batteries in a discharged state if not using the unit for some time. The batteries should be checked monthly and recharged as needed.
- **Protect** the battery from sources of ignition, such as flames and sparks.
- **Do not** charge the battery in direct sunlight.
- **Use only** PDG provided battery assemblies, and battery chargers with the Power Tilt system.
- **Unplug** the battery from the controller when not in use.
- **Do not** pull or yank on the battery cord to disconnect the cable. Do not carry the battery by the cord.
- **Do not** step on the battery.



- The Power Tilt is equipped with two modular lead-acid batteries intended to be alternated on and off the chair.
- Keep one battery on the chair and one on the battery charger, to ensure that a fully charged battery is always available.
- The Power Tilt controller consumes a small amount of power when plugged into the battery, even when the chair is not being used. If the chair is going to remain unused for several days, please unplug the power tilt to avoid draining the battery.
- Don't wait to change batteries until the battery runs down to zero and the Power Tilt stops working. To prolong battery life, change out the Power Tilt batteries at regular intervals, preferably daily.

## 12. INSTALLING AND REMOVING THE BATTERY

- To connect the battery module to the Power Tilt system, drop the latch receiver onto the mounting plate, located at the rear of the wheelchair, and plug in the cord.
- To disconnect the battery, first unplug the power connector, then simultaneously depress the release lever and pull the battery module directly upwards. Immediately plug the battery into the charger in order to prolong battery life.



### 13. BATTERY CHARGER

- **Read** the below instructions and the provided battery charger manufacturer's manual before using the battery charger.
- This battery charger is an off-board charger intended for the PDG 12V lead-acid battery only. The charging connector interfaces with the Power Tilt battery – see below schematic.
- To charge the battery, plug the power supply plug into a 120V wall socket and then plug the 3 pin XLR connector into the battery connector cable.
- Do not carry the battery charger on the wheelchair.
- Occasional use of the battery prior to the charging complete indication is acceptable if the need is urgent.

### 14. BATTERY CHARGER WIRING SCHEMATIC




## 15. TROUBLESHOOTING

Problem 	Possible Cause 	Solution 
Actuator whines but does not move	Low battery charge	Charge or swap out battery with a freshly charged battery.
	Actuator failed	Replace actuator.
	Frame movement obstructed	Ensure nothing is obstructing the movement of the upper frame.
Controller “clicks” and status LED turns on but actuator does not move, no whining sound heard from actuator	Actuator power cable unplugged	Plug in actuator power cable.
	Actuator cable damaged or actuator failed	Replace actuator.
Operate control switch, no sounds heard, actuator does not move	Control switch unplugged	Plug in control switch.
	Power cord unplugged	Plug in power cord.
	Tilt limit switch incorrectly adjusted	Adjust tilt limit switch.
	Battery dead	Charge or swap out battery with a fully charged battery.
	Tilt limit switch failed	See “Disabling the Tilt Limit Switch”.
	Control switch failed	See “Checking The Control Switches”.
	Controller failed	Replace controller
No LEDs lit on battery charger when charger is connected to the wall	Charger is not properly connected to the power outlet	Check all connections and cables.
	No power is being delivered to the wall outlet	Plug in another appliance to verify the power is working.
	Charger plugged into 240V power	Use 120V power only.
	Charger failed	Replace charger.

Green LED lit on battery charger when battery is initially connected, other LEDs not lighting up and then green LED starts blinking continuously	Charger fuse removed or blown	Open charger 5A fuse compartment and check that the fuse is not discolored. If the fuse is discolored, replace the charger.
	Battery fuse removed or blown, or battery failed.	Try charging the other battery, if the other battery charges then service the faulty battery.
	Battery not properly connected inside battery compartment	See "Servicing the Battery".
Battery runs out before 1 day of use	Battery is not being fully charged before use	Charge battery until the #5 LED on the charger is lit.
	Battery has failed	Replace battery, see "Battery Service".
Chair goes into tilt by itself	Control switch unplugged	Ensure control switch plug is fully seated and that the control switch cable is correctly routed.

## 16. SERVICING THE BATTERY

- **Do not** open the battery enclosure or service the battery without the manufacturer's written approval.
- **Always** wear eye protection when working with batteries. Lead acid batteries may explode.
- **Never** touch the contacts of the battery or short the contacts together. Shorting the contacts may cause the battery to overheat and explode.
- **Properly** dispose of used lead-acid batteries.
- **Do not** use third-party batteries with the Power Tilt system.
- Lead acid batteries have a finite service life and degrade over time. The life of the battery depends on the frequency of use, how deeply they are discharged, and the length of time spent in a discharged state before the batteries are recharged.
- To install a new lead-acid battery into the battery module:
  1. Disconnect the battery from the controller and remove the battery from the chair.
  2. Remove the 4 Phillips screws on the outside of the battery enclosure.
  3. Lift off the lid and pull the battery and wiring harness out of the enclosure.
  4. Disconnect the wiring harness from the battery terminals.
  5. Connect the new battery to the wiring harness (note polarity, see below schematic).
  6. Check the fuse before replacing the lid. To check the fuse, open the fuse compartment inside the battery enclosure and check that the fuse is present and not discolored or blackened. If the fuse is discolored, the fuse has failed and needs to be replaced.



## 17. DISABLING THE TILT LIMIT SWITCH

**WARNING:** Disconnect battery before adjusting or performing maintenance on wheelchair.

1. Identify the wires connecting the tilt limiter to the controller
2. If necessary, Gently Pull the wires for the tilt limiter clear of the controller enclosure.
3. Disconnect the quick-disconnect connections on the wires. Some force may be required. Hold the connectors, and not the wires. Pulling on the wires will likely result in damage to the wiring harness.
4. Remove the tilt limiter assembly from the wheelchair.
5. Reconnect the male and female connectors that are attached to the controller

The system will now travel through the full tilt range



Quick-disconnect connectors –  
Tilt limit sensor active





Quick-disconnect connectors –  
Tilt limit sensor disabled

## 18. CHECKING THE CONTROL SWITCHES – 3 POSITION SWITCH

**CAUTION:** Keep clear of lower frame and underside of upper frame when operating tilt mechanism


Verify that the switch is working correctly by using a digital multimeter to measure continuity between the three contacts on the switch plug. Contact 1 is closest to the cable end, and Button 1 is farthest from the cable end. See images below for reference.



1. Without depressing any of the buttons, check for continuity between all of the 3 contacts of the mini-plug connector (1-2, 1-3, and 2-3). If there is continuity then the switch has failed.
2. Depress Button 1 and check for continuity between (1-3). If there is no continuity then the switch has failed.
3. Depress Button 2 and check for continuity between (1-2). If there is no continuity then the switch has failed.

## 19. CHECKING THE CONTROL SWITCHES – MONO (EGG SWITCH)


1. Check for continuity between the two switch contacts with the switch not pushed. If there is continuity, then the switch has failed.
2. Check for continuity between the two switch contacts with the switch pushed. If there is no continuity, then the switch has failed.




Mono switch contacts

## 20. TROUBLESHOOTING FLOWCHART – POWER TILT SYSTEM

**WARNING:** Disconnect battery before adjusting or performing maintenance on wheelchair.



## 21. TROUBLESHOOTING FLOWCHART –BATTERY



## 22. WARRANTY INFORMATION

**PLEASE NOTE** – THE WARRANTY BELOW HAS BEEN DRAFTED TO COMPLY WITH FEDERAL LAW APPLICABLE TO PRODUCTS MANUFACTURED AFTER JULY 4<sup>th</sup>, 1975.

This warranty is extended only to the original purchaser/user of our products.

This warranty gives you specific legal rights and you may also have other legal rights, which vary from state to state.

PDG warrants its product, except for the seat cushion (which is not warranted), to be free from defects in materials and workmanship for a period of one (1) year from date of purchase. The side frames and cross members are warranted for the lifetime of the original purchaser/user. If within such warranty period any such product shall be proven to be defective, such product shall be repaired or replaced, at PDG's option.

This Warranty does not include any labor or shipping charges incurred in replacement part installation or repair of any such product. PDG's sole obligation and your exclusive remedy under this warranty shall be limited to such repair and/or replacement.

For warranty service, please contact the dealer from whom you purchased your PDG product. In the event you do not receive satisfactory warranty service, please write directly to PDG at the address on the back cover page, provide dealer's name, address, and date of purchase, indicate nature of the defect and, if the product is serialized, indicate the serial number. Do not return products to our factory without our prior consent. Limitations and exclusions: the foregoing warranty shall not apply to serial numbered products if the serial number has been removed or defaced, products subjected to negligence, accident, improper operation, maintenance or storage, products modified without PDG's express written consent including, but not limited to, modification through the use of unauthorized parts or attachments; products damaged by reason of repairs made to any component without the specific consent of PDG, or to a product damaged by circumstances beyond PDG's control, and such evaluation will be solely determined by PDG. The warranty shall not apply to problems arising from normal wear or failure to adhere to these instructions. The foregoing express warranty is exclusive and in lieu of any other warranties whatsoever, whether express or implied, including the implied warranties of merchantability and fitness for a particular purpose, and the sole remedy for violations of any warranty whatsoever, shall be limited to repair or replacement of the defective product pursuant to the terms contained herein. The application of any implied warranty whatsoever shall not extend beyond the duration of the express warranty provided herein. The manufacturer shall not be liable for any consequential or incidental damages whatsoever. This warranty shall be extended to comply with state/provincial laws and requirements.

## 23. TECHNICAL SPECIFICATIONS

Maximum rated occupant weight (Stellar GL, Stellar LEAP, and Fuze T50, Bentley)	250 lbs
Maximum rated occupant weight (Fuze T50 with 350 lbs Power Tilt upgrade)	350 lbs
Battery operation temperature range	-40° C to 60° C
Battery charging temperature range	-20° C to 50° C
Battery fuse and battery charger fuse	5A fast-acting blade type fuses
Battery voltage and capacity	12V 5.0 Ah
Battery cyclic endurance (IEC 60254-1)	200 charge/discharge cycles
Battery charging time	6 hours to 80%, 10 hours to 100% of capacity

## 24. APPENDIX

#	Description	Dwg #	Pages
A	Power Tilt – Fuze T50 – Chair Assembly	68950	22-25
B	Stellar GL - Power Tilt - UF/SF/LF Assembly	55661	26-29
C	Stellar LEAP - Power Tilt - UF/SF/LF Assembly	55651	30-33
D	Bentley - Power Tilt - UF/SF/LF Assembly	55685	34-37

1		2		3	
Rev	Revision description				
0	Initial Release				
Item	Part Name or Description			Number	250 QT
1	F2 - Power Tilt - Lower Frame Actuator Mount			68954	1
2	Power Tilt - Controller Enclosure - Asm			68468	1
3	5/16-18(3A) x 3" with Shank CBHCS - Black Magni 575			12899	4
4	5/16-18 (2B) Nylon Locknut - Thick - Black Zinc			12010	4
5	Power Tilt Actuator - LD20 - 200mm Stroke			36276	-
6	Power Tilt Actuator - LD3 - 200mm Stroke			12093	1
7	Power Tilt - Tilt Limit - Assembly			36279	1
8	5/16" ID x 1/2" OD x 1/16" Black Nylon Spacer			12512	6
9	5/16" ID x 3/4" OD x 5/16" L black nylon spacer			12012	1

### Actuator Mounting Block Assembly

Torque to 18 ft-lbs  
2 x TYP

### Actuator and Tilt Limit Assembly


Torque to 18 ft-lbs  
2 x TYP


PDG part number NA	Author Jim Forby
	Title: F2 - T50 - Power Tilt
© 2013 The information contained in this drawing is the sole property of PDG Mobility	Third angle projection 

Dimensioning and Tolerancing Standard ASME Y14.5-2009

4		5		6	
					ECO#
					1758
lbs	350 lbs				
QTY	QTY				
	1				
	1				
	4				
	4				
	1				
	-				
	1				
	6				
	1				
<p><b>Actuator Position Setup</b> Upper and lower frame set-up position must match</p> <p><b>0° - 40° Tilt Range</b> (Cannot be use with swing away front rigging)</p> <p><b>5° - 45° Tilt Range</b></p> <p><b>5° - 45° Tilt Range</b></p> <p><b>0° - 40° Tilt Range</b> (Cannot be used with swing away front rigging)</p>					
<p><b>Tilt Limiter Assembly</b></p> <p><b>Zip tie, then ensure that there is slack at max. tilt</b></p>					
Technical reference Torr Brown		Revised by Jim Formby		Approved by Torr Brown	
Lower Tilt - Assembly			File name (Drawing #) 68950	Status NA	
Tolerances unless otherwise stated: Angular: +/- 1° Other x +/-0.8, x.x +/-0.4, x.xx +/-0.15		Units mm	Sheet 1 of 2	Date of Issue 2020 Feb 05	Revision 0




4		5		6			
					ECO#		
					1758		
<b>Tilt BOM</b>							
Description	Number	250 lbs QTY	350 lbs QTY	A			
Lower Frame Actuator Mount	68954	1	1				
Controller Enclosure - Asm	68468	1	1				
Controller with Shank CBHCS - Black Magni 575	12899	4	4				
Locknut - Thick - Black Zinc	12010	4	4				
Actuator - LD20 - 200mm Stroke	36276	-	1				
Actuator - LD3 - 200mm Stroke	12093	1	-				
Limit - Assembly	36279	1	1				
Limit D x 1/16" Black Nylon Spacer	12512	6	6				
Limit D x 5/16" L black nylon spacer	12012	1	1				
Time Extension - Asm	68741	1	1				
Time Extension - Bridge	68353	1	1				
Driver Plate	68469	1	1				
Clamp	68451	2	2				
Limit Offset Spacer	36114	1	1				
Time - Masking Cap	68206	1	1	B			
Time - Rear Cap	68732	1	1				
Limit CBHCS - Black Magni 575	12898	4	4				
Head cap screw - black zinc	12121	4	4				
<b>Power supply, charger, connector and switch components</b>							
Description	Number	Toggle Switch QTY	Button Switch QTY			Egg Switch QTY	C
Power Supply Assembly	36248	2	2			2	
Power Supply	16056	1	1			1	
Power Supply 0.75A; 6 & 12 Volt	36693	1	1			1	
Power Supply XLR Connector	36694	1	1			1	
Power Supply Switch	36695	1	-	-			
Power Supply Switch	36696	-	1	-			
Power Supply Switch	36697	-	-	1			
Prepared by	Technical reference Torr Brown	Revised by Jim Formby	Approved by Torr Brown				
Lower Tilt - Assembly		File name (Drawing #) 68950		Status NA	D		
Notes	Tolerances unless otherwise stated: Angular: +/- 1° Other x +/-0.8, x.x +/-0.4, x.xx +/-0.15	Units mm	Sheet 2 of 2	Date of Issue 2020 Feb 05		Revision 0	




4	5	6
		ECO#
		1730

### Actuator, Secondary Strut and Tilt Limit Assembly



### Actuator, Secondary Strut and Tilt Limit Setup




Technical reference Torr Brown	Revised by Jim Formby	Approved by Torr Brown			
Power Tilt - UF-SF-LF Assembly		File name (Drawing #) 55661	Status NA		
Tolerances unless otherwise stated: Angular: +/- 1° Other x +/-0.8, x.x +/-0.4, x.xx +/-0.15		Units mm	Sheet 1 of 2	Date of Issue 2019 Jul 22	Revision 4

1	2	3
Rev	Revision description	
4	Changed battery mounting and updated BOM	

### Battery Mount Assembly (Left hand side of chair shown)

See power tilt manual (document 16056) for complete wiring schematic of the Power Tilt Accessory



### Other Parts: Battery, charger, connector and switch components

Item	Part Name or Description	Number	Toggle Switch QTY	BU
28	Power Tilt - Battery Assembly	36248	2	
29	Power Tilt - Manual	16056	1	
30	Power Tilt - Charger 0.75A; 6 & 12 Volt	36693	1	
31	Power Tilt -X-Connect XLR Connector	36694	1	
32	Power Tilt - Toggle Switch	36695	1	
33	Power Tilt - 2 Button Switch	36696	-	
34	Power Tilt - Egg Switch	36697	-	

PDG part number  
NA


Author  
Torr Brown

**pdg**<sup>®</sup>

Title:  
Stellar GL - Po

© 2013 The information contained in this drawing is the sole property of PDG Mobility

Third angle project




Dimensioning and Tolerancing Standard ASME Y14.5-2009



4		5		6	
					ECO#
					1730
<b>Stellar GL Power Tilt BOM</b>					
Part Name or Description		Number	QTY		A
Stellar GL - Power Tilt - SF Assembly		55662	1		
Stellar GL - Power Tilt - Controller Assembly		55664	1		
3/8-16 Nylon Locknut - Thick - Black Zinc		12021	6		
Oil Lite Bushing - Flanged - 1/2" ID x 5/8" OD x 1/2" L		12398	2		
Spring - Wave Disc - 0.650" ID x 0.855" OD x 0.060" Height		12480	2		
1/2" x 1/2" x 3/8-16 Socket Head Shoulder Bolt - Black Zinc		12241	2		
3/8-16 x 1" Hex Head Bolt		12408	4		
# 8 x 3/8" self-tapping phillips pan head screw		12067	3		
Power Tilt Actuator - LD20 - 100mm Stroke		36273	1		
G2 Secondary Gas Strut Assembly		55666	2		
Power Tilt - Tilt Limit - Assembly		36279	1		
Stellar Power Tilt - Operation Caution and Warning Label		16053	1		
5/16" ID x 3/4" OD x 5/16" L black nylon spacer		12012	8		
5/16" ID x 1/2" OD x 1/8" Black Nylon Spacer		12227	5		
5/16-18 x 5" Socket Head Cap Screw		12611	2		
5/16" ID x 1/2" OD x 1/16" Black Nylon Spacer		12512	2		
5/16-18(3A) x 2-3/4" with Shank CBHCS - Black Magni 575		12899	1		
5/16-18 (2B) Nylon Locknut - Thick - Black Zinc		12010	7		
5/16-18(3A) x 4" with Shank CBHCS - Black Magni 575		12899	2		
Accessory Extension Connector for Lower Frame		55218	1		
G2 - Power Tilt - Frame Extension for Battery Mount		55674	1		
Power Tilt - Receiver Plate		68469	1		
Stellar Power Tilt - Extension Cable - Assembly		36113	1		
Loop Clamp 3/4" ID x 2" LG		36690	1		
Lower Frame Tube Cap		55159	1		
5/16-18(3A) x 1-3/4" CBHCS - Black Magni 575		12898	2		
1/4-20 x 5/8" Socket Button Head - Black Zinc		12135	6		
Button Switch		Egg Switch	B		
QTY	QTY				
2	2				
1	1				
1	1				
1	1				
-	-				
1	-				
-	1				
Technical reference Torr Brown		Revised by Jim Formby		Approved by Torr Brown	
Power Tilt - UF-SF-LF Assembly			File name (Drawing #) 55661	Status NA	
Tolerances unless otherwise stated: Angular: +/- 1° Other x +/-0.8, x.x+/-0.4, x.xx+/-0.15		Units mm	Sheet 2 of 2	Date of Issue 2019 Jul 22	Revision 4

C


D





4		5		6	
					ECO#
					1730
Technical reference Torr Brown		Revised by Jim Formby		Approved by Torr Brown	
Power Tilt - UF-SF-LF Assembly			File name (Drawing #) 55651		Status NA
Tolerances unless otherwise stated: Angular: +/- 1° Other x +/-0.8, x.x +/-0.4, x.xx +/-0.15		Units mm	Sheet 1 of 2	Date of Issue 2019 Jul 22	Revision 3

Rev	1	2	3
3	Revision description Changed battery mounting, updated BOM		
			
Dimensioning and Tolerancing Standard ASME Y14.5-2009		PDG part number NA	Author Torr Brown
© 2013 The information contained in this drawing is the sole property of PDG Mobility		 Title: Stellar LEAP - Third angle project	

4		5		6	
				ECO#	
				1730	
Item	Part Name or Description	Number	QTY		
1	G2 - Lower Frame Assembly	55760	1		
2	Stellar LEAP - Power Tilt - Upper Frame Assembly	55655	1		
3	Stellar LEAP - Power Tilt - SF Assembly	55652	1		
4	G2 Secondary Gas Strut Assembly	55666	2		
5	Power Tilt Actuator - LD20 - 100mm Stroke	36273	1		
6	Stellar Power Tilt - Extension Cable - Assembly	36113	1		
7	Loop Clamp 3/4" ID x 2" LG	36690	1		
8	Stellar Power Tilt - Operation Caution and Warning Label	16053	1		
9	5/16" ID x 11/16" OD x 1/16" Washer - Black Zinc	12013	6		
10	5/16-18 (2B) Nylon Locknut - Thick - Black Zinc	12010	9		
11	G2 - Stellar GL - Upper Frame - Pivot Plate (RH/LH)	55734	1		
12	G2 - Stellar GL - Upper Frame - Pivot Plate (RH/LH)	55734	1		
13	5/16-18(3A) x 1" x 0.4" Shank FHCS - Black Magni 575	12897	4		
14	3/8-16 Nylon Locknut - Thick - Black Zinc	12021	6		
15	Oil Lite Bushing - Flanged - 1/2" ID x 5/8" OD x 1/2" L	12398	2		
16	Spring - Wave Disc - 0.650" ID x 0.855" OD x 0.060" Height	12480	2		
17	1/2" x 1/2" x 3/8-16 Socket Head Shoulder Bolt - Black Zinc	12241	2		
18	3/8-16 x 1" Hex Head Bolt	12408	4		
19	5/16" ID x 1/2" OD x 1/8" Black Nylon Spacer	12227	3		
20	5/16" ID x 3/4" OD x 5/16" L black nylon spacer	12012	8		
21	5/16-18 x 5" Socket Head Cap Screw	12611	2		
22	5/16-18(3A) x 3" with Shank CBHCS - Black Magni 575	12899	1		
23	1/4"ID x 1/2"OD x 0.06T Nylon Washer	12226	2		
24	Accessory Extension Connector for Lower Frame	55218	1		
25	G2 - Power Tilt - Frame Extension for Battery Mount	55674	1		
26	Power Tilt - Receiver Plate	68469	1		
27	5/16-18(3A) x 1-3/4" CBHCS - Black Magni 575	12898	2		
28	1/4-20 x 5/8" Socket Button Head - Black Zinc	12135	6		
29	Lower Frame Tube Cap	55159	1		
30	Power Tilt - Battery Assembly	36248	1		
31	Power Tilt - Charger 0.75A; 6 & 12 Volt	36693	1		
32	Power Tilt - X-Connect XLR Connector	36694	1		
33	Power Tilt Manual	16056	1		
Technical reference Torr Brown		Revised by Jim Formby		Approved by Torr Brown	
- Power Tilt - UF-SF-LF Assembly			File name (Drawing #) 55651	Status NA	
tion	Tolerances unless otherwise stated: Angular: +/- 1° Other x +/-0.8, x.x +/-0.4, x.xx +/-0.15	Units mm	Sheet 2 of 2	Date of Issue 2019 Jul 22	Revision 3



4		5		6	
					ECO#
					1751
Chang		Technical reference Andrew Chang		Revised by Jim Formby	
				Approved by Torr Brown	
er Tilt - UF-SF-LF - Asm				File name (Drawing #) 55685	
				Status NA	
Tolerances unless otherwise stated: Angular: +/- 1° Other: x +/-0.8, x.x +/-0.4, x.xx +/-0.15		Units mm		Date of Issue 2020 Jan 10	
		Sheet 1 of 2		Revision 0	

1	2	3
Rev	Revision description	
0	Initial Release	
		
<p>Notes:</p> <ol style="list-style-type: none"> <li>1. Material: As per detailed part</li> <li>2. Remove burrs &amp; sharp edges.</li> <li>3. Finish : As per detailed part</li> <li>4. Component weight [grams]: 23741.79</li> </ol> <p>Dimensioning and Tolerancing Standard ASME Y14.5-2009</p>		
PDG part number NA		Author Andrew Ch
		Title: Bentley - Pow
© 2013 The information contained in this drawing is the sole property of PDG Mobility		Third angle projection 

4		5		6	
					ECO#
					1751
Item	Part Name or Description	Number	QTY		
1	G2 - Lower Frame Assembly	55760	1		
2	Stellar GL - Power Tilt - Upper Frame Assembly	55668	1		
3	G2 Bentley - Power Tilt - Sub Frame - Asm	55686	1		
4	G2.1 - Bentley & GL - Upper Frame - Pivot Plate (LH/RH)	55916	1		
5	G2.1 - Bentley & GL - Upper Frame - Pivot Plate (RH/LH)	55916	1		
6	Power Tilt Actuator - LD20 - 100mm Stroke	36273	1		
7	G2 Secondary Gas Strut Assembly	55666	1		
8	Accessory Extension Connector for Lower Frame	55218	1		
9	G2 - Power Tilt - Frame Extension for Battery Mount	55674	1		
10	Power Tilt - Receiver Plate	68469	1		
11	Lower Frame Tube Cap	55159	1		
12	Loop Clamp 3/4" ID x 2" LG	36690	1		
13	3/8-16 x 1" Hex Head Bolt	12408	4		
14	1/2" x 1/2" x 3/8-16 Socket Head Shoulder Bolt - Black Zinc	12241	2		
15	3/8-16 Nylon Locknut - Thick - Black Zinc	12021	6		
16	Spring - Wave Disc - 0.650" ID x 0.855" OD x 0.060" Height	12480	2		
17	5/16-18(3A) x 1" x 0.4" Shank FHCS - Black Magni 575	12897	4		
18	5/16-18 (2B) Nylon Locknut - Thick - Black Zinc	12010	7		
19	5/16" ID x 11/16" OD x 1/16" Washer - Black Zinc	12013	4		
20	Oil Lite Bushing - Flanged - 1/2" ID x 5/8" OD x 1/2" L	12398	2		
21	5/16-18(3A) x 2-3/4" with Shank CBHCS - Black Magni 575	12899	1		
22	5/16-18(3A) x 4" with Shank CBHCS - Black Magni 575	12899	1		
23	5/16-18(3A) x 3-1/2" with Shank CBHCS - Black Magni 575	12899	1		
24	5/16-18 Nylon Locknut - Thin - Black Zinc	12066	2		
25	5/16-18(3A) x 1-3/4" CBHCS - Black Magni 575	12898	2		
26	1/4-20 x 5/8" Socket Button Head - Black Zinc	12135	7		
27	5/16" ID x 1/2" OD x 1/8" Black Nylon Spacer	12227	1		
28	5/16" ID x 3/4" OD x 5/16" L black nylon spacer	12012	2		
29	Power Tilt - Battery Assembly	36248	1		
30	Power Tilt - Charger 0.75A; 6 & 12 Volt	36693	1		
31	Power Tilt - X-Connect XLR Connector	36694	1		
ang	Technical reference Andrew Chang	Revised by Jim Formby	Approved by Torr Brown		
er Tilt - UF-SF-LF - Asm	File name (Drawing #) 55685		Status NA		
on	Tolerances unless otherwise stated: Angular: +/- 1° Other: x +/-0.8, x.x +/-0.4, x.xx +/-0.15	Units mm	Sheet 2 of 2	Date of Issue 2020 Jan 10	Revision 0

**25. NOTES:**

**NOTES:**



PDG PRODUCT DESIGN GROUP INC

103 – 318 East Kent Avenue South Vancouver, BC Canada V5X4N6

Phone: 604-323-9220 Fax: 604-323-9097 info@pdgmobility.com

[www.pdgmobility.com](http://www.pdgmobility.com)