



Janusz Krawczynski
PDG Product Design Group
#103-318 East Kent Ave, South,
Vancouver, BC, V5X 4N6,
Canada.

15th December 2010

Dear Janusz,

Reference: testing the Fuze T50 Manual wheelchair, our Laboratory reference is 491430.

This letter is to advise you that this Laboratory has finished partial testing your product to the following Standard method: AS3696-1:2008, AS3696-3:2008, AS3696-5:1989, AS3696-8:1998, ISO7176-7:1996 and the requirements of AS3695:1992 (except ignitibility, materials, finishes, effects of rain & storage conditions).

The official six (6) page report, dated 15th December 2010 and referencing job number 491430, identifies the sample, the test documents and the testing facility. All conformances with the Australian standards are indicated in the report.

The Fuze T50 Manual wheelchair as represented by the test sample, complies with the relevant requirements of the Australian Standard AS 3695 for fatigue strength tests conducted by the methods of two drum endurance and drop test only, for user mass 204 kg.

Please do not hesitate to contact us if you have any enquiries or comments regarding the test.

Kind regards,

Digitally signed by Wayne Wurfel
DN: CN = Wayne Wurfel, C = US,
O = NovitaTech Test Laboratory,
OU = Novita Children's Services
Reason: I attest to the accuracy
and integrity of this document
Date: 2010.12.15 13:54:29 +10'30'

Wayne Wurfel
Senior Test Technician, Test Laboratory
(NATA signatory)



NovitaTech (a division of Novita Children's Services)
171 Days Road, Regency Park, South Australia 5010
PO Box 2438, Regency Park, South Australia 5942
Telephone 1300 855 585 • Fax (08) 8243 8208
enquiries@novita.org.au • www.novitatech.org.au
ABN 53 137 537 636

