


STELLAR IMPACT



STELLAR IMPACT



The Stellar Impact was designed with the most demanding users in mind. It offers up to 45° of low knee rise tilt and a variety of available on-chair accessories. As such, the Stellar Impact is almost identical in function to the standard Stellar, but its ultra-tough steel construction makes it perfect for situations that require the most durable tilt chair available.


A Ultra-Tough Frame

Rugged steel construction results in a highly durable frame for the most demanding of users.



B Easily Adjustable Axle Position

Allows for fine adjustment of the chair's centre of gravity to allow for more stability and comfort.



C Force Adjustable Gas Struts

The impact absorbing gas strut and optional damper provide for superior energy dissipation.



D Dynamic Spring Back

Additional dampening is available with the PDG spring back.



Features

Ultra-tough construction	Designed for high fatigue applications
45° of manual tilt-in-space	Short wheelbase for maneuverability
Knees remain low in tilt	Default with spring back and compatible with Dynamic Rocker Back (DRB™)
Low seat-to-floor height	Transit tie-down available



SEAT WIDTH
14"–22"



DYNAMIC TILT
Up to 45°



SEAT DEPTH
17"–22"



BACK HEIGHT
25" & 30"



FRONT SEAT HEIGHT
14"–21"



MAX CAPACITY
250 lbs/113 kg



BACK ANGLE
90°–122°



CHAIR WEIGHT
64 lbs/29 kg