

Customer and Order Information

Purchase Order *	Tag	v.0601
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in *required

Frame Type

<input checked="" type="checkbox"/> 2401000	Bentley LT 20° Lightweight Manual Tilt Wheelchair Aluminum Frame. Not recommended for agitated users.	\$ 4,295
ⓘ For agitated or heavier users, choose the standard Bentley with steel frame. See separate order form.		

Frame Options

<input checked="" type="checkbox"/> 2402000	Lightweight Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/> 2402300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths					
<input type="checkbox"/> 1003014	14" Seat Width	<input type="checkbox"/> 1003015	15" Seat Width	<input type="checkbox"/> 1003016	16" Seat Width
<input type="checkbox"/> 1003017	17" Seat Width	<input type="checkbox"/> 1003018	18" Seat Width	<input type="checkbox"/> 1003019	19" Seat Width
<input type="checkbox"/> 1003020	20" Seat Width	<input type="checkbox"/> 1003021	21" Seat Width	<input type="checkbox"/> 1003022	22" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)					
<input type="checkbox"/> 1004016	16" Seat Depth	<input type="checkbox"/> 1004017	17" Seat Depth	<input type="checkbox"/> 1004018	18" Seat Depth
<input type="checkbox"/> 1004019	19" Seat Depth	<input type="checkbox"/> 1004020	20" Seat Depth		
Custom Seat Depths					
<input type="checkbox"/> 1004021	21" Seat Depth	<input type="checkbox"/> 1004022	22" Seat Depth	\$ 350	

Frame Color

PDG Premium Colors					
<input type="checkbox"/> 2405000	Tracer Red	<input type="checkbox"/> 1005001	Gloss Black	<input type="checkbox"/> 1005002	Silver
<input type="checkbox"/> 2405003	Black Shimmer	<input type="checkbox"/> 1005004	Orange Fusion	<input type="checkbox"/> 1005005	Karl's Imagination
<input type="checkbox"/> 2405006	Blue Shimmer	<input type="checkbox"/> 1005007	Gloss White	<input type="checkbox"/> 1005008	Dark Titanium
<input type="checkbox"/> 2405009	Geyser Green	<input type="checkbox"/> 1005010	Lollypop Blue	<input type="checkbox"/> 1005011	Lollypop Red
<input type="checkbox"/> 2405012	Bubble Gum Pink	<input type="checkbox"/> 1005013	Berry Purple	<input type="checkbox"/> 1005014	Safety Yellow
<input type="checkbox"/> 2405200	Custom Color Contact PDG for specifications.				\$ 300

Arm Type

<input type="checkbox"/> 2406000	No Armrest	
<input type="checkbox"/> 2406100	Single Post Fixed Height (12") Full Length Armpad only.	
Single Post Adjustable Height \$ 275		
<input type="checkbox"/> 2406111	8"-11"	<input type="checkbox"/> 1006114 10"-14"

Armrests

Regular Armrests

- 2407100 Desk Length - 2" x 10" 1007101 Full Length - 2" x 14"

Wide Armrests

May interfere with wheels in some tilt configurations.

\$ 200

- 2407110 Desk Length - 3" x 10" 1007111 Full Length - 3" x 12"

Gel Armrests

May interfere with wheels in some tilt configurations.

\$ 135

- 2407120 Full Length - 2" x 12" 1007121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad

May interfere with wheels in some tilt configurations.

\$ 125

- 2407130 Right (4.5" x 18") 1007131 Left (4.5" x 18")

- 2407310 **Clothing Guards** Available for 10"-14" Height range only.

\$ 110

Seat Pan

- 2408000 **No Seat Pan**

- 2408001 **Solid Aluminum Seat Pan**

\$ 120

Cushion

Refer to the enclosed Seating & Positioning order form for PDG's on-chair Cushion options. Standard PDG discounts apply.

Seat Height

➤ **Select desired Seat-to-Floor Height** STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High					
	12"	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	15"	16"	17"	18"	19"	20"
8" Casters	RW	Low					High					Extra-High					
	16"	15"	16"	17"	18"	19"	17"	18"	19"	17"	18"	19"	20"	17"	18"	19"	20"

Wheel size _____

Caster size _____

STF height _____

Range _____

Casters

Urethane

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2411061 6" x 1.25" Caster size
 2411063 6" x 1.50" Caster size
 2411081 8" x 1.25" Caster size
 2411084 8" x 2.00" Caster size

\$ 85

\$ 135

Pneumatic

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2411161 6" x 1.25" Caster size
 2411181 8" x 1.25" Caster size
 2411184 8" x 2.00" Caster size

\$ 110

\$ 120

\$ 150

- 2411260 **6" Soft Roll**

\$ 95

- 2411400 **Ultra Sport Fork** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.

\$ 395

Rear Wheels

Composite

- 2412012 12" Wheel size 1012016 16" Wheel size 1012020 20" Wheel size
 2412022 22" Wheel size 1012024 24" Wheel size

Spoke

\$ 75

- 2412120 20" Wheel size 1012122 22" Wheel size 1012124 24" Wheel size

Axle Type	
<input type="checkbox"/> 2413000	Quick Release Rear Axles
<input type="checkbox"/> 2413010	Fixed Rear Axles

Rear Tires	
<input type="checkbox"/> 2414000	Urethane
<input type="checkbox"/> 2414010	Pneumatic May result in a finished Seat Height 1/2" greater than stated. \$ 60
<input type="checkbox"/> 2414020	Pneumatic with Flat-Free Inserts May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95

Hand Rims	
<input type="checkbox"/> 2415000	No Hand Rims
<input type="checkbox"/> 2415010	Standard 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 2415020	Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 135

Wheel Locks	
<input type="checkbox"/> 2416010	Push to Lock
<input type="checkbox"/> 2416020	Pull to Lock
<input type="checkbox"/> 2416030	Attendant Wheel Lock (Push to Lock) Not compatible with TTDS. \$ 85
<input type="checkbox"/> 2416040	User & Attendant Wheel Lock Available only on 20", 22", or 24" Rear Wheels. \$ 125
<input type="checkbox"/> 2416100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock. \$ 60
<input type="checkbox"/> 2416200	Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts. \$ 450

Anti-Tippers	
<input checked="" type="checkbox"/> 2417010	Rear Anti-Tippers Required Safety Option. \$ 60

Back Post Style	
<input checked="" type="checkbox"/> 2418010	Push Bar Back
Back Post Height Seat Pan to Handle	
<input type="checkbox"/> 2418425	25" Height <input type="checkbox"/> 1018431 31" Height

Back Options	
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate. \$ 150	
<input type="checkbox"/> 2419010	Button on Outside <input type="checkbox"/> 1019011 Button on Inside
Seat to Back Angle	
<input type="checkbox"/> 2419090	Set at 90° <input type="checkbox"/> 1019098 Set at 98° <input type="checkbox"/> 1019106 Set at 106°
<input type="checkbox"/> 2419114	Set at 114° <input type="checkbox"/> 1019122 Set at 122°
<input type="checkbox"/> 2419210	10" Backrest Extension Requires the Angle/Height Adjustable Push Bar Handle. \$ 150
<input checked="" type="checkbox"/> 2419220	Fold Down Back Not compatible with Reclining Back.

Back Upholstery	
<input type="checkbox"/> 2420000	No Back Upholstery
Standard Back height measured from Seat Pan to top of Upholstery. \$ 175	
<input type="checkbox"/> 2420021	Set at 21" <input type="checkbox"/> 1020022 Set at 22" <input type="checkbox"/> 1020023 Set at 23" Not available with 25" Back.
Short 25" Back only. Back height measured from Seat Pan to top of Upholstery. \$ 175	
<input type="checkbox"/> 2420118	Set at 18" <input type="checkbox"/> 1020119 Set at 19" <input type="checkbox"/> 1020120 Set at 20"

See next page for more options >

<input type="checkbox"/>	Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.				\$ 175
<input type="checkbox"/>	2420224 Set at 24"	<input type="checkbox"/>	1020225 Set at 25"	<input type="checkbox"/>	1020226 Set at 26"
<input type="checkbox"/>	Adjustable Tension Back height measured from Seat Pan to top of Upholstery.				\$ 295
<input type="checkbox"/>	2420321 Set at 21"	<input type="checkbox"/>	1020322 Set at 22"	<input type="checkbox"/>	1020323 Set at 23"

Positioning Backrest

i Refer to the enclosed Seating & Positioning order form for PDG's on-chair Backrest options. Standard PDG discounts apply.

Tilt/Recline Options

<input type="checkbox"/>	Attendant Tilt only				
<input type="checkbox"/>	2423011 Trigger on Right Back Post	<input type="checkbox"/>	1023012 Trigger on Left Back Post		
<input type="checkbox"/>	User Tilt only Full Length Arm Pads only.				\$ 100
<input type="checkbox"/>	2423021 Trigger on Right Armrest	<input type="checkbox"/>	1023022 Trigger on Left Armrest		
<input type="checkbox"/>	Attendant & User Tilt Full Length Arm Pads only.				\$ 195
<input type="checkbox"/>	2423031 Trigger on Right Back Post	<input type="checkbox"/>	1023032 Trigger on Left Back Post		
<input type="checkbox"/>	2423041 Trigger on Right Armrest	<input type="checkbox"/>	1023042 Trigger on Left Armrest		
<input type="checkbox"/>	2423110 Trigger Extension				\$ 50
<input type="checkbox"/>	2423310 Reclining Backrest (Manual) 30° maximum recline. Not compatible with Backrest Extension, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery. Requires longer Frame. Maximum Seat Depth is 19".				\$ 895

Headrest

i Refer to the enclosed Seating & Positioning order form for PDG's on-chair Headrest options. Standard PDG discounts apply.

<input type="checkbox"/>	2424010 Headrest Mounting Plate Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back.	\$ 95
--------------------------	---	-------

Front Rigging

<input type="checkbox"/>	2425000 No Front Rigging	<input type="checkbox"/>	2425001 Sector Blocks only		
<input type="checkbox"/>	PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out				
<input type="checkbox"/>	2425410 Fixed 70°				\$ 95
<input type="checkbox"/>	2425420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) Adjust angle by loosening one bolt. Tool-free leg-length adjustment.				\$ 295
<input type="checkbox"/>	Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.				
<input type="checkbox"/>	2425130 Swing Away 60°	<input type="checkbox"/>	2425140 Swing Away 70°		
<input type="checkbox"/>	2425200 Swing Away Elevating Legrests Includes Calf Pads.				\$ 375
<input type="checkbox"/>	2425300 1" Offset Hangers for Swing Away				\$ 75
<input type="checkbox"/>	2425110 Standard Plug In 70° Legacy Design. Not compatible with Easy-Out.				

Foot Plates

<input type="checkbox"/>	2426010 Composite				
<input type="checkbox"/>	2426020 Aluminum Angle Adjustable				\$ 155
<input type="checkbox"/>	2426030 One Piece Fixed Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	2426040 One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	2426050 Interlocking Angle Adjustable Seat Widths of 14"-21" only.				\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	Calf Pad Not compatible with Gel Padding.				
<input type="checkbox"/>	2427012 For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.				\$ 200
<input type="checkbox"/>	2427010 For Swing Away Front Rigging				\$ 135

See next page for more options >

<input type="checkbox"/>	2427020	Gel Padding Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pad			
<input type="checkbox"/>	2427032	For Easy-Out™ Front Rigging	\$ 200
<input type="checkbox"/>	2427030	For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/>	2427040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	2427050	Padded Calf Strap	\$ 195

Extension Tubes

Regular Extension Tubes			
<input type="checkbox"/>	2428010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	2428020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75

Positioning Belts & Straps

Refer to the enclosed Seating & Positioning order form for PDG's on-chair Belts/Straps options. Standard PDG discounts apply.

Miscellaneous Options

<input type="checkbox"/>		One-Arm Drive Requires 24" Composite Wheels with Urethane Tires. 18" STF minimum. Limits tilt. Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.	\$ 3,500
<input type="checkbox"/>	2430010	Installed on Right side	<input type="checkbox"/>
	1030011	Installed on Left side	
<input type="checkbox"/>	2430050	Protective Lower Padding	\$ 825
<input type="checkbox"/>	2430080	Gel Hip Guides Reduces available Seat Width by 2".	\$ 375
<input type="checkbox"/>	2430170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/bentley/ to learn more about the **Bentley 20° Manual Tilt Wheelchair**



Cushion

1. Select the type ↓	2. Select the model ↓
--------------------------------	---------------------------------

General Use	Stealth Simplicity 350 lbs weight capacity. <input type="checkbox"/> 4009141 16"-20" width, 16"-20" depth \$ 310 <input type="checkbox"/> 4009142 22"-24" width, 18"-20" depth \$ 370
--------------------	--

Positioning	ROHO AirLITE 300 lbs weight capacity. <input type="checkbox"/> 4009120 14"-20" width, 14"-18" depth \$ 375
--------------------	--

Skin Protection & Positioning	Stealth Spectrum Foam 650 lbs weight capacity. <input type="checkbox"/> 4009143 16"-20" width, 16"-24" depth \$ 580 <input type="checkbox"/> 4009144 22"-32" width, 18"-28" depth \$ 630
	Stealth Tru-Comfort 2 450 lbs weight capacity. <input type="checkbox"/> 4009146 16"-20" width, 16"-24" depth \$ 515 <input type="checkbox"/> 4009147 22" width, 18"-24" depth \$ 575
	nxt BioFit 300 lbs weight capacity. <input type="checkbox"/> 4009101 14"-20" width, 16"-20" depth \$ 595 <input type="checkbox"/> 4009102 22"-24" width, 16"-20" depth \$ 625
	Matrx VI 300 lbs weight capacity. <input type="checkbox"/> 4009111 14"-20" width, 15"-20" depth \$ 680 <input type="checkbox"/> 4009112 22" width, 18"-20" depth \$ 745
	Matrx VI Heavy Duty 600 lbs weight capacity. <input type="checkbox"/> 4009115 18"-21" width, 18"-26" depth \$ 1,050 <input type="checkbox"/> 4009116 22"-30" width, 18"-26" depth \$ 1,250

Skin Protection & Adjustable	ROHO Single Valve High Profile No weight limit if cushion is properly sized. <input type="checkbox"/> 4009125 14"-20" width, 16"-22" depth \$ 840 <input type="checkbox"/> 4009126 22"-24" width, 16"-22" depth \$ 960
---	---

Skin Protection, Positioning & Adjustable	ROHO Quadro Select High Profile No weight limit if cushion is properly sized. <input type="checkbox"/> 4009121 14"-20" width, 16"-22" depth \$ 840 <input type="checkbox"/> 4009122 22"-24" width, 16"-22" depth \$ 960
--	--

Extra-Wide No weight limit if cushion is properly sized. Custom sizes available. Custom-made foam cushion with fluid-proof yet breathable fabric cover.	PDG Solar Basic High resiliency polyfoam with 2" medium-density base and 1" lower-density top.	<input type="checkbox"/> 4009200 20"-28" width, 16"-24" depth \$ 395 <input type="checkbox"/> 4009201 29"-36" width, 16"-24" depth \$ 495 <input type="checkbox"/> 4009202 37"-48" width, 16"-24" depth \$ 595
	PDG Solar Contour Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups.	<input type="checkbox"/> 4009210 20"-28" width, 16"-24" depth \$ 695 <input type="checkbox"/> 4009211 29"-36" width, 16"-24" depth \$ 795 <input type="checkbox"/> 4009212 37"-48" width, 16"-24" depth \$ 895
	PDG Solar Pressure Reduction All features of the Contour, plus added pincore latex rubber and containment foam outline.	<input type="checkbox"/> 4009220 20"-28" width, 16"-24" depth \$ 1,095 <input type="checkbox"/> 4009221 29"-36" width, 16"-24" depth \$ 1,195 <input type="checkbox"/> 4009222 37"-48" width, 16"-24" depth \$ 1,295

3. Specify the size →	Cushion Size Sizes based on availability from the cushion manufacturer. Specify Size (Width" x Depth") _____
------------------------------	--

Backrest

1. Select the type ↓	2. Select the contour depth ↓	3. Select the model ↓
-------------------------	----------------------------------	--------------------------

Positioning, posterior-lateral	2" depth	<input type="checkbox"/> 4021141 Stealth ADI AL Active 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 945
	4" depth	<input type="checkbox"/> 4021142 Stealth ADI AL Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,045
	6" depth	<input type="checkbox"/> 4021143 Stealth ADI AL Extra Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,145

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> 4021110 Matrx Elite 300 lbs weight capacity. 12"-20" width, 10"-20" height \$ 955
		<input type="checkbox"/> 4021111 Matrx Elite Heavy Duty 500 lbs weight capacity. 21" width, 16"/20" height \$ 1,195
		<input type="checkbox"/> 4021112 22"-24" width, 16"/20" height \$ 1,295
		<input type="checkbox"/> 4021114 25"-30" width, 16"/20" height \$ 1,995
	3.5" to 4" depth	<input type="checkbox"/> 4021120 ROHO Agility Min 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 900
		<input type="checkbox"/> 4021100 nxt Xtend 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021105 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021122 ROHO Agility Mid 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 930
	6" to 7" depth	<input type="checkbox"/> 4021101 nxt Xtend Deep 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021106 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021116 Matrx Elite Deep 300 lbs weight capacity. 14"-20" width, 12"-20" height \$ 1,145
		<input type="checkbox"/> 4021117 22"-24" width, 12"-20" height \$ 1,295
<input type="checkbox"/> 4021124 ROHO Agility Max 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 1,085		

4. Specify the size →	Backrest Size Sizes based on availability from the backrest manufacturer. Specify Size (Width" x Height") _____
-----------------------	---

Pediatric Seating System

4030400 **Diego Pediatric Seating System**

Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10 (25.40)	12 (30.48)	14 (35.56)	16 (40.64)
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓

A Head System Head Pad & Hardware

B Contoured Back System Base, Cushion & Cover

FlipLock Configurable as a between canes or full width set up.

Width Growth Adapted from our TFB System, the back wings allow for a pediatric growth adjustment.

Back Hardware Cover Easy to remove back cover for quick on-demand adjustments and a more aesthetic look.

Coolcore Covers Our Diego System uses the Coolcore cover technology which adapt all-year-round for a more comfortable experience.

Laterals Modular Laterals either Fixed or Swing Away

Contoured Seat System Base, Cushion, Cover & Hardware

diego a Stealth Product

WC-20 certified

Headrest

Stealth Comfort Plus Headrest	Stealth Hardware		
	<input type="checkbox"/> 4024110	Multi-Axis HMO475	\$ 225
	<input type="checkbox"/> 4024111	Link TWB480	\$ 350
	Comfort Plus Pad		\$ 275
	<input type="checkbox"/> 4024115	Large (10" x 6")	<input type="checkbox"/> 4024116 X-Large (14" x 6")
nxt MultiFit Head Support	<input type="checkbox"/> 4024100 nxt Hardware		\$ 285
	Multi-Adjustable Pad		\$ 275
	<input type="checkbox"/> 4024105	Contoured (9" x 5")	<input type="checkbox"/> 4024107 X-Large (14" x 5")
	<input type="checkbox"/> 4024106	Large (10" x 5")	

Positioning Belts & Straps

Bodypoint	Hip Belts		
	<input type="checkbox"/> 4029010	Center-Pull Auto Style	\$ 115
	<input type="checkbox"/> 4029015	Center-Pull Auto Style, Padded	\$ 215
	<input type="checkbox"/> 4029020	4-Point, Padded	\$ 235
	<input type="checkbox"/> 4029030	Chest Strap	\$ 165

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.

PDG proudly supplies on-chair seating and positioning products from these manufacturers:

