

Customer and Order Information

Purchase Order *	Tag	v.072318
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 1001000	Bentley 20° Manual Tilt Wheelchair	\$ 3,995
---	---	----------

Frame Options

<input type="checkbox"/> 1002000	Standard Base and Seat Frame 250 lbs Weight Capacity.	
<input type="checkbox"/> 1002050	Upgraded Base and Seat Frame 300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	\$ 200
<input type="checkbox"/> 1002100	Heavy Duty Base and Seat Frame 450 lbs Weight Capacity. Available for all sizes, but required for Seat Width of 23" or wider. Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	\$ 300
<input type="checkbox"/> 1002200	Extended Base with Rear Floor Jacks Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Locks, HD Anti-Tipper, IV Pole, O2 Holder.	\$ 600
<input type="checkbox"/> 1002300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths			
<input type="checkbox"/> 1003014	14" Seat Width	<input type="checkbox"/> 1003015	15" Seat Width
<input type="checkbox"/> 1003017	17" Seat Width	<input type="checkbox"/> 1003018	18" Seat Width
<input type="checkbox"/> 1003020	20" Seat Width	<input type="checkbox"/> 1003016	16" Seat Width
<input type="checkbox"/> 1003021	21" Seat Width	<input type="checkbox"/> 1003019	19" Seat Width
<input type="checkbox"/> 1003022	22" Seat Width		\$ 215
Custom Seat Widths Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.			
<input type="checkbox"/> 1003023	23" Seat Width	<input type="checkbox"/> 1003024	24" Seat Width
<input type="checkbox"/> 1003026	26" Seat Width	<input type="checkbox"/> 1003025	25" Seat Width
<input type="checkbox"/> 1003027	27" Seat Width	<input type="checkbox"/> 1003028	28" Seat Width
<input type="checkbox"/> 1003030	30" Seat Width	<input type="checkbox"/> 1003029	29" Seat Width
<input type="checkbox"/> 1003031	31" Seat Width	<input type="checkbox"/> 1003032	32" Seat Width
			\$ 625

Seat Depth

Standard Seat Depths (Adjustable)			
<input type="checkbox"/> 1004016	16" Seat Depth	<input type="checkbox"/> 1004017	17" Seat Depth
<input type="checkbox"/> 1004019	19" Seat Depth	<input type="checkbox"/> 1004018	18" Seat Depth
<input type="checkbox"/> 1004020	20" Seat Depth		
Custom Seat Depths			
<input type="checkbox"/> 1004021	21" Seat Depth	<input type="checkbox"/> 1004022	22" Seat Depth
			\$ 325


Frame Colour

PDG Premium Colours			
<input type="checkbox"/> 1005000	Tracer Red	<input type="checkbox"/> 1005001	Gloss Black
<input type="checkbox"/> 1005003	Black Shimmer	<input type="checkbox"/> 1005002	Silver
<input type="checkbox"/> 1005006	Blue Shimmer	<input type="checkbox"/> 1005004	Orange Fusion
<input type="checkbox"/> 1005009	Geyser Green	<input type="checkbox"/> 1005007	Gloss White
<input type="checkbox"/> 1005012	Bubble Gum Pink	<input type="checkbox"/> 1005008	Dark Titanium
<input type="checkbox"/> 1005010	Lollipop Blue	<input type="checkbox"/> 1005011	Lollipop Red
<input type="checkbox"/> 1005013	Berry Purple	<input type="checkbox"/> 1005014	Safety Yellow
<input type="checkbox"/> 1005200	Custom Colour Contact PDG for specifications.		\$ 300

Arm Type		
<input type="checkbox"/>	1006000	No Armrest
<input type="checkbox"/>	1006100	Single Post Fixed Height (12") Full Length Armpad only.
		Single Post Adjustable Height Not compatible with the Heavy Duty Base and/or Seat Width greater than 22". \$ 255
<input type="checkbox"/>	1006111	8"-11" <input type="checkbox"/> 1006114 10"-14"
<input type="checkbox"/>	1006200	Dual Post Fixed Height (12") Full Length Armpad only.
		Dual Post Adjustable Height \$ 255
<input type="checkbox"/>	1006211	8"-11" <input type="checkbox"/> 1006214 10"-14" <input type="checkbox"/> 1006216 12"-16"
		Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. \$ 255
<input type="checkbox"/>	1006310	Set at 10" <input type="checkbox"/> 1006311 Set at 11" <input type="checkbox"/> 1006312 Set at 12"
<input type="checkbox"/>	1006313	Set at 13" <input type="checkbox"/> 1006314 Set at 14"

Armrests		
		Regular Armpads
<input type="checkbox"/>	1007100	Desk Length - 2" x 10" <input type="checkbox"/> 1007101 Full Length - 2" x 14"
<input type="checkbox"/>	1007102	Full Length Premium Grade - 2" x 14" \$ 35
		Wide Armpads May interfere with wheels in some tilt configurations. \$ 115
<input type="checkbox"/>	1007110	Desk Length - 3" x 10" <input type="checkbox"/> 1007111 Full Length - 3" x 12"
		Gel Armpads May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1007120	Full Length - 2" x 12" <input type="checkbox"/> 1007121 Full Length - 3.5" x 12"
		Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1007130	Right (4.5" x 18") <input type="checkbox"/> 1007131 Left (4.5" x 18")
<input type="checkbox"/>	1007220	Armrest Wideners for Dual Post Armrests Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards. \$ 375
		Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.
<input type="checkbox"/>	1007310	For Single Post Armrests Available for 10"-14" Height range only. \$ 110
<input type="checkbox"/>	1007320	For Dual Post Armrests \$ 165
		Armrest Extensions Not compatible with Flip-up or Single Post Armrests. \$ 60
<input type="checkbox"/>	1007400	3" Extension - Forward <input type="checkbox"/> 1007401 3" Extension - Rearward

Seat Pan		
<input type="checkbox"/>	1008000	No Seat Pan
<input type="checkbox"/>	1008001	Solid Aluminum Seat Pan \$ 120

Cushion		
<input type="checkbox"/>	1009000	No Cushion
		nxt BioFit
<input type="checkbox"/>	1009101	14"-20" width, 16"-20" depth  \$ 595
<input type="checkbox"/>	1009102	22"-24" width, 16"-20" depth \$ 625
		Cushion Size Available in even number sizes only.
	1009400	Specify Size (Width" x Depth") _____

Seat Height

Select desired Seat-to-Floor Height STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High				
	12"	13"	14"	15"			15"	16"	17"			15"	16"	17"	18"	
16"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
20"	13"	14"	15"			15"	16"	17"			15"	16"	17"	18"		
22"		14"	15"	16"			16"	17"	18"			16"	17"	18"	19"	
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"

8" Casters	RW	Low			High			Extra-High							
	16"			15"	16"	17"			17"	18"	19"		17"	18"	19"
24"			15"	16"	17"			17"	18"	19"		17"	18"	19"	20"

Wheel size _____

Caster size _____

STF height _____

Range _____

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1011061 6" x 1.25" Caster size
- 1011063 6" x 1.50" Caster size \$ 85
- 1011081 8" x 1.25" Caster size
- 1011084 8" x 2.00" Caster size \$ 135

Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1011161 6" x 1.25" Caster size \$ 110
- 1011181 8" x 1.25" Caster size \$ 120
- 1011184 8" x 2.00" Caster size \$ 150

- 1011260 **6" Soft Roll** \$ 95

- 1011400 **Ultra Sport Fork** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. \$ 395

Rear Wheels

Composite

- 1012012 12" Wheel size
- 1012016 16" Wheel size
- 1012020 20" Wheel size
- 1012022 22" Wheel size
- 1012024 24" Wheel size

Spoke \$ 75

- 1012120 20" Wheel size
- 1012122 22" Wheel size
- 1012124 24" Wheel size

- 1012224 **24" Heavy Duty Spoke** Required for Heavy Duty Base and Seat Frame. \$ 350

Axle Type

- 1013000 **Quick Release Rear Axles**

- 1013010 **Fixed Rear Axles**

Rear Tires

- 1014000 **Urethane**

- 1014010 **Pneumatic** Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. \$ 60

- 1014020 **Pneumatic with Flat-Free Inserts** \$ 95

Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.

Hand Rims

- 1015000 **No Hand Rims**

- 1015010 **Standard** 20", 22", and 24" Rear Wheels only.

- 1015020 **Plastic Coated** 20", 22", and 24" Rear Wheels only. \$ 135

Wheel Locks			
<input type="checkbox"/>	1016010	Push to Lock	
<input type="checkbox"/>	1016020	Pull to Lock	
<input type="checkbox"/>	1016030	Attendant Wheel Lock (Push to Lock)	\$ 85
Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.			
<input type="checkbox"/>	1016040	User & Attendant Wheel Lock	\$ 125
Available only on 20", 22", or 24" Rear Wheels.			
<input type="checkbox"/>	1016100	Wheel Lock Extensions	\$ 60
Not compatible with Attendant Wheel Lock.			
<input type="checkbox"/>	1016200	Hill Holders	\$ 350
Requires Pneumatic Rear Tires with Flat-Free Inserts.			

Anti-Tippers			
<input checked="" type="checkbox"/>	1017010	Rear Anti-Tippers	\$ 60
Required Safety Option.			
<input type="checkbox"/>	1017020	Front Anti-Tippers	\$ 300
<input type="checkbox"/>	1017030	Heavy Duty Rear Anti-Tipper	\$ 600
Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.			


Back Post Style			
<input type="checkbox"/>	1018010	Push Bar Back	
<input type="checkbox"/>	1018020	Cane Handle Back	
Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.			
Back Post Height Seat Pan to Handle			
<input type="checkbox"/>	1018425	25" Height	<input type="checkbox"/>
			1018431 31" Height

Back Options			
Angle/Height Adjustable Push Bar Handle			\$ 150
Not compatible with Headrest Mounting Plate or Cane Handle Back.			
<input type="checkbox"/>	1019010	Button on Outside	<input type="checkbox"/>
			1019011 Button on Inside
Seat to Back Angle			
<input type="checkbox"/>	1019090	Set at 90°	<input type="checkbox"/>
			1019098 Set at 98°
			<input type="checkbox"/>
			1019106 Set at 106°
<input type="checkbox"/>	1019114	Set at 114°	<input type="checkbox"/>
			1019122 Set at 122°
<input type="checkbox"/>	1019210	10" Backrest Extension	\$ 135
Not compatible with Angle/Height Adjustable Push Bar Handle.			
<input type="checkbox"/>	1019220	Fold Down Back	
Not compatible with Reclining Back.			
<input checked="" type="checkbox"/>	1019230	Spring Back Suspension	
Not compatible with Reclining Back.			
<input type="checkbox"/>	1019240	Dynamic Rocker Back Assembly	\$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height. Requires longer Frame.			


Back Upholstery			
<input type="checkbox"/>	1020000	No Back Upholstery	
Standard Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	1020021	Set at 21"	<input type="checkbox"/>
			1020022 Set at 22"
Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	1020118	Set at 18"	<input type="checkbox"/>
			1020119 Set at 19"
			<input type="checkbox"/>
			1020120 Set at 20"
Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	1020224	Set at 24"	<input type="checkbox"/>
			1020225 Set at 25"
			<input type="checkbox"/>
			1020226 Set at 26"
Adjustable Tension Back height measured from Seat Pan to top of Upholstery.			\$ 295
<input type="checkbox"/>	1020321	Set at 21"	<input type="checkbox"/>
			1020322 Set at 22"
			<input type="checkbox"/>
			1020323 Set at 23"

Backrest			
<input type="checkbox"/>	1021000	No Backrest	

See next page for more Backrest options >

next Xtend		
<input type="checkbox"/>	1021100 16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1021105 22"-24" width, 16"-20" height	\$ 720
		
next Xtend Deep		
<input type="checkbox"/>	1021101 16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1021106 22"-24" width, 16"-20" height	\$ 720
Backrest Size Available in even number sizes only.		
<input type="checkbox"/>	1021400 Specify Size (Width" x Height") _____	

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1023011 Trigger on Right Back Post	<input type="checkbox"/>	1023012 Trigger on Left Back Post
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only.			\$ 100
<input type="checkbox"/>	1023021 Trigger on Right Armrest	<input type="checkbox"/>	1023022 Trigger on Left Armrest
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only.			\$ 195
<input type="checkbox"/>	1023031 Trigger on Right Back Post	<input type="checkbox"/>	1023032 Trigger on Left Back Post
<input type="checkbox"/>	1023041 Trigger on Right Armrest	<input type="checkbox"/>	1023042 Trigger on Left Armrest
<input type="checkbox"/>	1023110 Trigger Extension		\$ 50
<input type="checkbox"/>	1023410 Reclining Backrest (Manual)		\$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19".			

Headrest			
<input type="checkbox"/>	1024000 No Headrest		
<input type="checkbox"/>	1024010 Headrest Mounting Plate	\$ 75	
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	1024100 next MultiFit Head Support Hardware	\$ 235	
			
<input type="checkbox"/>	1024105 next MultiFit Head Support Pad	\$ 245	
<input type="checkbox"/>	1024105 Regular (9" x 5")	<input type="checkbox"/>	1024106 Large (10" x 5")
		<input type="checkbox"/>	1024107 X-Large (14" x 5")

Front Rigging			
<input type="checkbox"/>	1025000 No Front Rigging	<input type="checkbox"/>	1025001 Sector Blocks only
PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.			
<input type="checkbox"/>	1025410 Fixed 70°		\$ 95
<input type="checkbox"/>	1025420 Swing In/Away Universal Angle Adjustable (0°-110°)	Set at (use increments of 10°): _____	\$ 195
<input type="checkbox"/>	1025430 Swing In/Away Elevating Legrests	High knee pivot. Includes reinforced premium Calf Pad mount.	\$ 495
Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.			
<input type="checkbox"/>	1025110 Plug In 70°		
<input type="checkbox"/>	1025130 Swing Away 60°	<input type="checkbox"/>	1025140 Swing Away 70°
<input type="checkbox"/>	1025200 Swing Away Elevating Legrests		\$ 375
<input type="checkbox"/>	1025300 1" Offset Hangers for Swing Away		\$ 75

Foot Plates			
<input type="checkbox"/>	1026010 Composite		
<input type="checkbox"/>	1026020 Aluminum Angle Adjustable	\$ 155	
<input type="checkbox"/>	1026030 One Piece Fixed	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1026040 One Piece Flip-up	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1026050 Interlocking Angle Adjustable	Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1027010	Swing Away 60° & 70° Calf Pad	Not compatible with 60° & 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1027020	Swing Away 60° & 70° Gel Padding	Not compatible with 60° & 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1027030	Swing Away ELR Gel Knee Pad		\$ 100
<input type="checkbox"/>	1027040	Adjustable Length Heel Loops		\$ 50
<input type="checkbox"/>	1027050	Bodypoint - Padded Calf Strap	Sizes greater than 22" are not padded.	\$ 105
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.				
<input type="checkbox"/>	1027060	For Plug In Front Rigging		\$ 190
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging	Not available on 13" STF Height.	\$ 380

Extension Tubes

Regular Extension Tubes				
<input type="checkbox"/>	1028010	Short	14"-18" for Easy-Out SA Angle Adjustable 11"-15" for Easy Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1028020	Long	18"-24" for Easy-Out SA Angle Adjustable 15"-21" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. \$ 100				
<input type="checkbox"/>	1028110	Standard (13"-17")	<input type="checkbox"/>	1028120 Long (15"-19")

Positioning Straps

Bodypoint - Hip Belts				
<input type="checkbox"/>	1029010	Center-Pull Auto Style		\$ 105
<input type="checkbox"/>	1029015	Center-Pull Auto Style, Padded		\$ 155
<input type="checkbox"/>	1029020	4-Point, Padded		\$ 175
<input type="checkbox"/>	1029030	Bodypoint - Chest Strap		\$ 85

Miscellaneous Options

One-Arm Drive Maximum Seat Width of 22". 24" Composite Wheels with Urethane Tires required. 18" STF minimum. Limits tilt. \$ 3,500 Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.				
<input type="checkbox"/>	1030010	Installed on Right side	<input type="checkbox"/>	1030011 Installed on Left side
<input type="checkbox"/>	1030050	Protective Lower Padding		\$ 425
<input type="checkbox"/>	1030080	Gel Hip Guides	Reduces available Seat Width by 2".	\$ 325
<input type="checkbox"/>	1030110	IV Pole	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$ 295
<input type="checkbox"/>	1030120	O2 Holder	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$ 395
<input type="checkbox"/>	1030170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/bentley/ to learn more about the **Bentley 20° Manual Tilt Wheelchair**

