

### Customer and Order Information

Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

### Frame Type

<input checked="" type="checkbox"/> 1001000	<b>Bentley 20° Manual Tilt Wheelchair</b>	\$ 3,995
---	---	----------

### Frame Options

<input type="checkbox"/> 1002000	<b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.	
<input type="checkbox"/> 1002050	<b>Upgraded Base and Seat Frame</b> 300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	\$ 200
<input type="checkbox"/> 1002100	<b>Heavy Duty Base and Seat Frame</b> 450 lbs Weight Capacity. Available for all sizes, but required for Seat Width of 23" or wider. Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	\$ 300
<input type="checkbox"/> 1002200	<b>Extended Base with Rear Floor Jacks</b> Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Locks, HD Anti-Tipper, IV Pole, O2 Holder.	\$ 600
<input type="checkbox"/> 1002300	<b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

### Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

<b>Standard Seat Widths</b>			
<input type="checkbox"/> 1003014	14" Seat Width	<input type="checkbox"/> 1003015	15" Seat Width
<input type="checkbox"/> 1003017	17" Seat Width	<input type="checkbox"/> 1003018	18" Seat Width
<input type="checkbox"/> 1003020	20" Seat Width	<input type="checkbox"/> 1003016	16" Seat Width
<input type="checkbox"/> 1003021	21" Seat Width	<input type="checkbox"/> 1003019	19" Seat Width
<input type="checkbox"/> 1003022	22" Seat Width		\$ 215
<b>Custom Seat Widths</b> Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.			
<input type="checkbox"/> 1003023	23" Seat Width	<input type="checkbox"/> 1003024	24" Seat Width
<input type="checkbox"/> 1003026	26" Seat Width	<input type="checkbox"/> 1003025	25" Seat Width
<input type="checkbox"/> 1003027	27" Seat Width	<input type="checkbox"/> 1003028	28" Seat Width
<input type="checkbox"/> 1003030	30" Seat Width	<input type="checkbox"/> 1003029	29" Seat Width
<input type="checkbox"/> 1003031	31" Seat Width	<input type="checkbox"/> 1003032	32" Seat Width
			\$ 325
			\$ 325
			\$ 625
			\$ 625

### Seat Depth

<b>Standard Seat Depths (Adjustable)</b>			
<input type="checkbox"/> 1004016	16" Seat Depth	<input type="checkbox"/> 1004017	17" Seat Depth
<input type="checkbox"/> 1004019	19" Seat Depth	<input type="checkbox"/> 1004020	20" Seat Depth
<input type="checkbox"/> 1004018	18" Seat Depth		
<b>Custom Seat Depths</b>			
<input type="checkbox"/> 1004021	21" Seat Depth	<input type="checkbox"/> 1004022	22" Seat Depth
			\$ 325


### Frame Colour

<b>PDG Premium Colours</b>			
<input type="checkbox"/> 1005000	Tracer Red	<input type="checkbox"/> 1005001	Gloss Black
<input type="checkbox"/> 1005003	Black Shimmer	<input type="checkbox"/> 1005002	Silver
<input type="checkbox"/> 1005006	Blue Shimmer	<input type="checkbox"/> 1005004	Orange Fusion
<input type="checkbox"/> 1005009	Geyser Green	<input type="checkbox"/> 1005007	Gloss White
<input type="checkbox"/> 1005012	Bubble Gum Pink	<input type="checkbox"/> 1005008	Dark Titanium
<input type="checkbox"/> 1005010	Lollipop Blue	<input type="checkbox"/> 1005011	Lollipop Red
<input type="checkbox"/> 1005013	Berry Purple	<input type="checkbox"/> 1005014	Safety Yellow
<input type="checkbox"/> 1005200	<b>Custom Colour</b> Contact PDG for specifications.		\$ 300

Arm Type		
<input type="checkbox"/>	1006000	<b>No Armrest</b>
<input type="checkbox"/>	1006100	<b>Single Post Fixed Height (12")</b> Full Length Armpad only.
		<b>Single Post Adjustable Height</b> Not compatible with the Heavy Duty Base and/or Seat Width greater than 22". <span style="float: right;">\$ 255</span>
<input type="checkbox"/>	1006111	8"-11" <input type="checkbox"/>
<input type="checkbox"/>	1006114	10"-14"
<input type="checkbox"/>	1006200	<b>Dual Post Fixed Height (12")</b> Full Length Armpad only.
		<b>Dual Post Adjustable Height</b> <span style="float: right;">\$ 255</span>
<input type="checkbox"/>	1006211	8"-11" <input type="checkbox"/>
<input type="checkbox"/>	1006214	10"-14" <input type="checkbox"/>
<input type="checkbox"/>	1006216	12"-16"
		<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. <span style="float: right;">\$ 255</span>
<input type="checkbox"/>	1006310	Set at 10" <input type="checkbox"/>
<input type="checkbox"/>	1006311	Set at 11" <input type="checkbox"/>
<input type="checkbox"/>	1006312	Set at 12" <input type="checkbox"/>
<input type="checkbox"/>	1006313	Set at 13" <input type="checkbox"/>
<input type="checkbox"/>	1006314	Set at 14" <input type="checkbox"/>

Armrests		
		<b>Regular Armpads</b>
<input type="checkbox"/>	1007100	Desk Length - 2" x 10" <input type="checkbox"/>
<input type="checkbox"/>	1007101	Full Length - 2" x 14" <span style="float: right;">\$ 35</span>
<input type="checkbox"/>	1007102	Full Length Premium Grade - 2" x 14" <span style="float: right;">\$ 115</span>
		<b>Wide Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 115</span>
<input type="checkbox"/>	1007110	Desk Length - 3" x 10" <input type="checkbox"/>
<input type="checkbox"/>	1007111	Full Length - 3" x 12" <span style="float: right;">\$ 125</span>
		<b>Gel Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>
<input type="checkbox"/>	1007120	Full Length - 2" x 12" <input type="checkbox"/>
<input type="checkbox"/>	1007121	Full Length - 3.5" x 12" <span style="float: right;">\$ 125</span>
		<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>
<input type="checkbox"/>	1007130	Right (4.5" x 18") <input type="checkbox"/>
<input type="checkbox"/>	1007131	Left (4.5" x 18") <span style="float: right;">\$ 375</span>
<input type="checkbox"/>	1007220	<b>Armrest Wideners for Dual Post Armrests</b> Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards. <span style="float: right;">\$ 375</span>
		<b>Clothing Guards</b> Not compatible with Flip-up Armrests, Armrest Wideners.
<input type="checkbox"/>	1007310	For Single Post Armrests Available for 10"-14" Height range only. <span style="float: right;">\$ 110</span>
<input type="checkbox"/>	1007320	For Dual Post Armrests <span style="float: right;">\$ 165</span>
		<b>Armrest Extensions</b> Not compatible with Flip-up or Single Post Armrests. <span style="float: right;">\$ 60</span>
<input type="checkbox"/>	1007400	3" Extension - Forward <input type="checkbox"/>
<input type="checkbox"/>	1007401	3" Extension - Rearward

Seat Pan		
<input type="checkbox"/>	1008000	<b>No Seat Pan</b>
<input type="checkbox"/>	1008001	<b>Solid Aluminum Seat Pan</b> <span style="float: right;">\$ 120</span>

Cushion		
<input type="checkbox"/>	1009000	<b>No Cushion</b>
		<b>nxt BioFit</b>
<input type="checkbox"/>	1009101	14"-20" width, 16"-20" depth  <span style="float: right;">\$ 595</span>
<input type="checkbox"/>	1009102	22"-24" width, 16"-20" depth <span style="float: right;">\$ 625</span>
		<b>Cushion Size</b> Available in even number sizes only.
1009400	Specify Size (Width" x Depth") _____	

### Seat Height

**Select desired Seat-to-Floor Height** STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High					
	12"	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	15"	16"	17"	18"	19"	20"
12"																	
16"				15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
20"	13"	14"	15"				15"	16"	17"			15"	16"	17"	18"		
22"		14"	15"	16"				16"	17"	18"			16"	17"	18"	19"	
24"			15"	16"	17"				17"	18"	19"			17"	18"	19"	20"

8" Casters	RW	Low			High			Extra-High			
	16"	15"	16"	17"	17"	18"	19"	17"	18"	19"	20"
16"			15"	16"	17"			17"	18"	19"	20"
24"			15"	16"	17"			17"	18"	19"	20"

Wheel size \_\_\_\_\_

Caster size \_\_\_\_\_

STF height \_\_\_\_\_

Range \_\_\_\_\_

### Casters

**Urethane** 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1011061 6" x 1.25" Caster size
- 1011063 6" x 1.50" Caster size \$ 85
- 1011081 8" x 1.25" Caster size
- 1011084 8" x 2.00" Caster size \$ 135

**Pneumatic** 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1011161 6" x 1.25" Caster size \$ 110
- 1011181 8" x 1.25" Caster size \$ 120
- 1011184 8" x 2.00" Caster size \$ 150

- 1011260 **6" Soft Roll** \$ 95

- 1011400 **Ultra Sport Fork** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. \$ 395

### Rear Wheels

**Composite**

- 1012012 12" Wheel size
- 1012016 16" Wheel size
- 1012020 20" Wheel size
- 1012022 22" Wheel size
- 1012024 24" Wheel size

**Spoke** \$ 75

- 1012120 20" Wheel size
- 1012122 22" Wheel size
- 1012124 24" Wheel size

- 1012224 **24" Heavy Duty Spoke** Required for Heavy Duty Base and Seat Frame. \$ 350

### Axle Type

- 1013000 **Quick Release Rear Axles**

- 1013010 **Fixed Rear Axles**

### Rear Tires

- 1014000 **Urethane**

- 1014010 **Pneumatic** Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. \$ 60

- 1014020 **Pneumatic with Flat-Free Inserts** \$ 95

Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.

### Hand Rims

- 1015000 **No Hand Rims**

- 1015010 **Standard** 20", 22", and 24" Rear Wheels only.

- 1015020 **Plastic Coated** 20", 22", and 24" Rear Wheels only. \$ 135

Wheel Locks			
<input type="checkbox"/>	1016010	<b>Push to Lock</b>	
<input type="checkbox"/>	1016020	<b>Pull to Lock</b>	
<input type="checkbox"/>	1016030	<b>Attendant Wheel Lock (Push to Lock)</b>	\$ 85
Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.			
<input type="checkbox"/>	1016040	<b>User &amp; Attendant Wheel Lock</b>	\$ 125
Available only on 20", 22", or 24" Rear Wheels.			
<input type="checkbox"/>	1016100	<b>Wheel Lock Extensions</b>	\$ 60
Not compatible with Attendant Wheel Lock.			
<input type="checkbox"/>	1016200	<b>Hill Holders</b>	\$ 350
Requires Pneumatic Rear Tires with Flat-Free Inserts.			

Anti-Tippers			
<input checked="" type="checkbox"/>	1017010	<b>Rear Anti-Tippers</b>	\$ 60
Required Safety Option.			
<input type="checkbox"/>	1017020	<b>Front Anti-Tippers</b>	\$ 300
<input type="checkbox"/>	1017030	<b>Heavy Duty Rear Anti-Tipper</b>	\$ 600
Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.			

Back Post Style			
<input type="checkbox"/>	1018010	<b>Push Bar Back</b>	
<input type="checkbox"/>	1018020	<b>Cane Handle Back</b>	
Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.			
<b>Back Post Height</b> Seat Pan to Handle			
<input type="checkbox"/>	1018425	25" Height	<input type="checkbox"/>
	1018431	31" Height	

Back Options			
<b>Angle/Height Adjustable Push Bar Handle</b>			\$ 150
Not compatible with Headrest Mounting Plate or Cane Handle Back.			
<input type="checkbox"/>	1019010	Button on Outside	<input type="checkbox"/>
	1019011	Button on Inside	
<b>Seat to Back Angle</b>			
<input type="checkbox"/>	1019090	Set at 90°	<input type="checkbox"/>
	1019098	Set at 98°	<input type="checkbox"/>
	1019106	Set at 106°	
<input type="checkbox"/>	1019114	Set at 114°	<input type="checkbox"/>
	1019122	Set at 122°	
<input type="checkbox"/>	1019210	<b>10" Backrest Extension</b>	\$ 135
Not compatible with Angle/Height Adjustable Push Bar Handle.			
<input type="checkbox"/>	1019220	<b>Fold Down Back</b>	
Not compatible with Reclining Back.			
<input checked="" type="checkbox"/>	1019230	<b>Spring Back Suspension</b>	
Not compatible with Reclining Back.			
<input type="checkbox"/>	1019240	<b>Dynamic Rocker Back Assembly</b>	\$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height. Requires longer Frame.			

Back Upholstery			
<input type="checkbox"/>	1020000	<b>No Back Upholstery</b>	
<b>Standard</b> Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	1020021	Set at 21"	<input type="checkbox"/>
	1020022	Set at 22"	
Not available with 25" Back.			
<b>Short</b> 25" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	1020118	Set at 18"	<input type="checkbox"/>
	1020119	Set at 19"	<input type="checkbox"/>
	1020120	Set at 20"	
<b>Tall</b> 31" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	1020224	Set at 24"	<input type="checkbox"/>
	1020225	Set at 25"	<input type="checkbox"/>
	1020226	Set at 26"	
<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery.			\$ 295
<input type="checkbox"/>	1020321	Set at 21"	<input type="checkbox"/>
	1020322	Set at 22"	<input type="checkbox"/>
	1020323	Set at 23"	

Backrest			
<input type="checkbox"/>	1021000	<b>No Backrest</b>	

See next page for more Backrest options >

nxt Optima HA Active Contour			
<input type="checkbox"/>	1021100	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1021105	22"-24" width, 16"-20" height	\$ 720
nxt Optima HA Plus Contour			
<input type="checkbox"/>	1021101	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1021106	22"-24" width, 16"-20" height	\$ 720
<b>Backrest Size</b> Available in even number sizes only.			
<input type="checkbox"/>	1021400	Specify Size (Width" x Height") _____	



Tilt/Recline Options						
Attendant Tilt only						
<input type="checkbox"/>	1023011	Trigger on Right Back Post	<input type="checkbox"/>	1023012	Trigger on Left Back Post	
<b>User Tilt only</b> Not available on widths greater than 22". Full Length Arm Pads only.						\$ 100
<input type="checkbox"/>	1023021	Trigger on Right Armrest	<input type="checkbox"/>	1023022	Trigger on Left Armrest	
<b>Attendant &amp; User Tilt</b> Not available on widths greater than 22". Full Length Arm Pads only.						\$ 195
<input type="checkbox"/>	1023031	Trigger on Right Back Post	<input type="checkbox"/>	1023032	Trigger on Left Back Post	
<input type="checkbox"/>	1023041	Trigger on Right Armrest	<input type="checkbox"/>	1023042	Trigger on Left Armrest	
<input type="checkbox"/>	1023110	<b>Trigger Extension</b>				\$ 50
<input type="checkbox"/>	1023410	<b>Reclining Backrest (Manual)</b>				\$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19".						

Headrest								
<input type="checkbox"/>	1024000	<b>No Headrest</b>						
<input type="checkbox"/>	1024010	<b>Headrest Mounting Plate</b>	\$ 75					
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.								
<input type="checkbox"/>	1024100	<b>nxt MultiFit Head Support Hardware</b>	\$ 235					
<b>nxt MultiFit Head Support Pad</b>			\$ 245					
<input type="checkbox"/>	1024105	Regular (9" x 5")	<input type="checkbox"/>	1024106	Large (10" x 5")	<input type="checkbox"/>	1024107	X-Large (14" x 5")



Front Rigging						
<input type="checkbox"/>	1025000	<b>No Front Rigging</b>				
<b>PDG Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame.						
<input type="checkbox"/>	1025410	Fixed 70°	Available in March	\$ 95		
<input type="checkbox"/>	1025420	Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____	Available in April	\$ 195		
<input type="checkbox"/>	1025430	Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.	Available in April	\$ 495		
<b>Standard Front Rigging</b> Legacy Design. Not compatible with Easy-Out™.						
<input type="checkbox"/>	1025110	Plug In 70°				
<input type="checkbox"/>	1025130	Swing Away 60°	<input type="checkbox"/>	1025140	Swing Away 70°	
<input type="checkbox"/>	1025200	Swing Away Elevating Legrests				\$ 375
<input type="checkbox"/>	1025300	1" Offset Hangers for Swing Away				\$ 75

Foot Plates			
<input type="checkbox"/>	1026010	<b>Composite</b>	
<input type="checkbox"/>	1026020	<b>Aluminum Angle Adjustable</b>	\$ 155
<input type="checkbox"/>	1026030	<b>One Piece Fixed</b> Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1026040	<b>One Piece Flip-up</b> Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1026050	<b>Interlocking Angle Adjustable</b> Seat Widths of 14"-21" only.	\$ 195

### Footrest/Legrest Accessories

<input type="checkbox"/>	1027010	<b>Swing Away 60° &amp; 70° Calf Pad</b>	Not compatible with 60° & 70° Gel Padding.	\$	135
<input type="checkbox"/>	1027020	<b>Swing Away 60° &amp; 70° Gel Padding</b>	Not compatible with 60° & 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$	100
<input type="checkbox"/>	1027030	<b>Swing Away ELR Gel Knee Pad</b>		\$	100
<input type="checkbox"/>	1027040	<b>Adjustable Length Heel Loops</b>		\$	50
<input type="checkbox"/>	1027050	<b>Bodypoint - Padded Calf Strap</b>	Sizes greater than 22" are not padded.	\$	105
<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.					
<input type="checkbox"/>	1027060	For Plug In Front Rigging		\$	190
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging	Not available on 13" STF Height.	\$	380

### Extension Tubes

<b>Regular Extension Tubes</b>					
<input type="checkbox"/>	1028010	Short	13"-17" for Easy-Out SA Angle Adjustable 10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging		
<input type="checkbox"/>	1028020	Long	16"-22.5" for Easy-Out SA Angle Adjustable 12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$	75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only.					
<input type="checkbox"/>	1028110	Standard (13"-17")		<input type="checkbox"/>	1028120 Long (15"-19")

### Positioning Straps

<b>Bodypoint - Hip Belts</b>					
<input type="checkbox"/>	1029010	Center-Pull Auto Style		\$	105
<input type="checkbox"/>	1029015	Center-Pull Auto Style, Padded		\$	155
<input type="checkbox"/>	1029020	4-Point, Padded		\$	175
<input type="checkbox"/>	1029030	<b>Bodypoint - Chest Strap</b>		\$	85



### Miscellaneous Options

<b>One-Arm Drive</b> Maximum Seat Width of 22". 24" Composite Wheels with Urethane Tires required. 18" STF minimum. Limits tilt. Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.					
<input type="checkbox"/>	1030010	Installed on Right side	<input type="checkbox"/>	1030011	Installed on Left side
<input type="checkbox"/>	1030050	<b>Protective Lower Padding</b>		\$	425
<input type="checkbox"/>	1030080	<b>Gel Hip Guides</b>		\$	325
<input type="checkbox"/>	1030110	<b>IV Pole</b>	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	295
<input type="checkbox"/>	1030120	<b>O2 Holder</b>	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	395
<input type="checkbox"/>	1030170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.	\$	10

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/bentley/](http://pdgmobility.com/products/bentley/) to learn more about the **Bentley 20° Manual Tilt Wheelchair**

