

Customer and Order Information

| | | |
|--|-----|-----------|
| Purchase Order * | Tag | v.031518 |
| Dealer/Provider Name * | | |
| Purchaser Name and Contact Information * | | |
| Billing Address * | | |
| Shipping Address * | | *required |

Frame Type

Elevation Ultralight Manual Wheelchair Extended delivery lead times may apply. This wheelchair is not recommended for children under 12 years.

| | | |
|---|---|----------|
| <input checked="" type="checkbox"/> 2001000 | Ultralight Manual Rigid Wheelchair | \$ 3,900 |
| <input checked="" type="checkbox"/> 2001001 | Patient Controlled Anterior/Posterior Back Angle Recline System | \$ 850 |
| <input checked="" type="checkbox"/> 2001002 | Patient Controlled Anterior Tilt & Elevation System | \$ 1,800 |

Frame Options

2002000 **Standard Base and Seat Frame** 250 lbs Weight Capacity.

User Weight Range Required for selection of Gas Struts (Force measured in Newton - N)

| | | | |
|----------------------------------|----------------------------|----------------------------------|----------------------------|
| <input type="checkbox"/> 2002150 | 65 lbs to 90 lbs (150 N) | <input type="checkbox"/> 2002200 | 90 lbs to 140 lbs (200 N) |
| <input type="checkbox"/> 2002300 | 140 lbs to 200 lbs (300 N) | <input type="checkbox"/> 2002400 | 200 lbs to 250 lbs (400 N) |

Seat Width

Standard Seat Widths

| | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2003012 | 12" Seat Width | <input type="checkbox"/> 2003013 | 13" Seat Width | <input type="checkbox"/> 2003014 | 14" Seat Width |
| <input type="checkbox"/> 2003015 | 15" Seat Width | <input type="checkbox"/> 2003016 | 16" Seat Width | <input type="checkbox"/> 2003017 | 17" Seat Width |
| <input type="checkbox"/> 2003018 | 18" Seat Width | <input type="checkbox"/> 2003019 | 19" Seat Width | | |

Seat Depth

Standard Seat Depths

| | | | | | |
|----------------------------------|----------------|----------------------------------|----------------|----------------------------------|----------------|
| <input type="checkbox"/> 2004014 | 14" Seat Depth | <input type="checkbox"/> 2004015 | 15" Seat Depth | <input type="checkbox"/> 2004016 | 16" Seat Depth |
| <input type="checkbox"/> 2004017 | 17" Seat Depth | <input type="checkbox"/> 2004018 | 18" Seat Depth | <input type="checkbox"/> 2004019 | 19" Seat Depth |
| <input type="checkbox"/> 2004020 | 20" Seat Depth | | | | |

Frame Colour

PDG Premium Colours

| | | | | | |
|----------------------------------|----------------------|----------------------------------|---------------|----------------------------------|--------------------|
| <input type="checkbox"/> 2005000 | Tracer Red | <input type="checkbox"/> 2005001 | Gloss Black | <input type="checkbox"/> 2005002 | Silver |
| <input type="checkbox"/> 2005003 | Black Shimmer | <input type="checkbox"/> 2005004 | Orange Fusion | <input type="checkbox"/> 2005005 | Karl's Imagination |
| <input type="checkbox"/> 2005006 | Blue Shimmer | <input type="checkbox"/> 2005007 | Gloss White | <input type="checkbox"/> 2005008 | Dark Titanium |
| <input type="checkbox"/> 2005009 | Geyser Green | <input type="checkbox"/> 2005010 | Lollypop Blue | <input type="checkbox"/> 2005011 | Lollypop Red |
| <input type="checkbox"/> 2005012 | Bubble Gum Pink | <input type="checkbox"/> 2005013 | Berry Purple | <input type="checkbox"/> 2005014 | Safety Yellow |
| <input type="checkbox"/> 2005200 | Custom Colour | Contact PDG for specifications. | | | \$ 300 |

Arm Type

| | | |
|----------------------------------|----------------------|--|
| <input type="checkbox"/> 2006000 | No Armrest | |
| <input type="checkbox"/> 2006001 | Rigid Armrest | Not compatible with Push to Lock, Fold Down Back. \$ 150 |

Seat Option

| | | |
|----------------------------------|----------------------------|---|
| <input type="checkbox"/> 2008110 | Standard Seat Sling | |
| <input type="checkbox"/> 2008100 | No Seat Sling | Third-party Seat Sling must be installed prior to wheelchair use. |

Cushion

2009000 **No Cushion**

Seat Height

24" Wheel Size

2010018 18" STF 2010019 19" STF 2010020 20" STF

25" Wheel Size

Not compatible with 4" Caster with Frog Leg Uni-tine Fork.

2010118 18.5" STF 2010119 19.5" STF 2010120 20.5" STF

26" Wheel Size

Not compatible with Frog Leg Uni-tine Fork.

2010219 19" STF 1610220 20" STF 1610221 21" STF

Casters

Composite

2011040 4" x 1.0" Composite Sport 2011042 4" x 1.5" Composite Softroll 2011052 5" x 1.5" Composite Softroll

Aluminium

2011140 4" x 1.0" Aluminum Sport 2011142 4" x 1.5" Aluminum Softroll 2011152 5" x 1.5" Aluminum Softroll \$ 125

Caster Fork & Fork Suspension

2011410 Standard Aluminum Fork

2011420 Frog Legs Uni-tine Fork \$ 195

2011421 Frog Legs Ultra-Sport Suspension Forks Compatible with 4" Caster only. \$ 450

Rear Wheels

Spinergy Wheels

2012310 Spinergy Wire

2012311 Spinergy SPOX Everyday \$ 745

2012312 Spinergy SPOX Sport \$ 745

2012313 Spinergy LX Everyday \$ 895

2012314 Spinergy LX Sport \$ 895

2012315 **Spinergy Lite Extreme LX Flexirims** LX Wheels only. Not compatible with 26" Wheel size. \$ 1,245



Select Wheel Camber

2012400 0° Camber 2012402 2° Camber

2012404 4° Camber 2012406 6° Camber

Axle Type

2013000 **Quick Release Rear Axles**

Rear Tires

2014110 **Primo Pneumatic Tires**


2014120 **Schwalbe Marathon Plus Evolution**  \$ 125

Hand Rims

2015000 **No Hand Rims**

2015100 **Aluminum Clear Anodized**

2015110 **Plastic Coated** \$ 135

2015130 **BBraver Rough Rider**  \$ 750

1 → Black 24"-26" Wheel size White 24"-25" Wheel size Orange 24"-25" Wheel size

2 → Long Short

Wheel Locks

- 2016010 **Push to Lock** Not compatible with Rigid Armrest.
- 2016090 **Composite Scissor Locks**

Anti-Tippers

- 2017010 **Rear Anti-Tippers** \$ 145

Back Options

- 2019020 **Cane Style Back Post** \$ 195
- 2019221 **Rigid Back**
- 2019220 **Fold Down Back** Not compatible with Rigid Armrest. \$ 150
- Seat Back Height**
- 2019300 Short (8.5"-12") - Rigidizer Bar Height = 5.5"
- 2019301 Medium (10.5"-16") - Rigidizer Bar Height = 7.5"
- 2019302 Tall (16.5"-20") - Rigidizer Bar Height = 9.5"
- 2019310 **Custom Rigidizer Bar Height** _____ Must be at least 4.5' above Seat Sling. \$ 150
- 2019311 **Deep Rigidizer Bar** For installing third-party Backrest. \$ 150

Back Upholstery

- 2020010 **Standard Back Upholstery**
- 2020000 **No Back Upholstery** Third-party Back Upholstery must be installed prior to wheelchair use.

Footrest

- 2025010 **Standard Footrest**
- 2025011 **Angle Adjustable High Mount**
- 2026020 **Aluminum Angle Adjustable Foot Plates**
- 2026060 **Flip-up Foot Plates**
- 2026075 **Fixed Heel Retention Foot Plate** \$ 265

Calf Straps

- 2027050 **Bodypoint - Padded Calf Strap** \$ 105


Positioning Straps

- Bodypoint - Hip Belts**
Push button with flat-mount end fitting. Recommended for mounting through slots in side guards, to allow "anti-shear" feature to be used.
- 2029010 Center-Pull Auto Style \$ 105
- 2029015 Center-Pull Auto Style, Padded \$ 155
- 2029200 **Bodypoint - EvoFlex Pelvic Stabilizer** \$ 275
Push button recommended for when anti-shear feature is not used. Mounted directly to the seat tubes using the included seat tube mounting bracket.

Side Guards

- 2030400 **No Side Guards**
- Non-Lifting Side Guards**
- 2030410 Rigid - Small 2030411 Rigid - Regular \$ 295
- 2030415 Folding - Regular \$ 350
- Lifting (Anti-Shear) Side Guards**
- 2030420 Rigid - Regular \$ 600
- 2030425 Folding - Regular 2030426 Folding - Large Not compatible with Rigid Armrest. \$ 650

Miscellaneous Options

| | | | | |
|--------------------------|---------|----------------------------------|---|--------|
| <input type="checkbox"/> | 2030450 | Alber Wheel Mount Adapter |  | \$ 600 |
| Quad Triggers | | | | |
| <input type="checkbox"/> | 2030461 | Inboard Mounting | | \$ 115 |
| <input type="checkbox"/> | 2030462 | Outboard Mounting | | \$ 115 |
| <input type="checkbox"/> | 2030463 | Thumb Loops | | \$ 130 |

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to elevation.pdgmobility.com to learn more about the **Elevation Ultralight Manual Wheelchair**

