

Visit the new website: [ElevationbyPDG.com](http://ElevationbyPDG.com)

Customer and Order Information		
Purchase Order *	Tag	v.070118
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type		
<b>Elevation Ultralight Manual Wheelchair</b> Extended delivery lead times may apply. This wheelchair is not recommended for children under 12 years.		
<input checked="" type="checkbox"/> 2001000	Ultralight Manual Rigid Wheelchair	\$ 3,900
<input checked="" type="checkbox"/> 2001001	Patient Controlled Anterior/Posterior Back Angle Recline System	\$ 850
<input checked="" type="checkbox"/> 2001002	Patient Controlled Anterior Tilt & Elevation System	\$ 1,800

Frame Options		
<input checked="" type="checkbox"/> 2002000	<b>Standard Base and Seat Frame</b>	250 lbs Weight Capacity.
<b>User Weight Range</b> Required for selection of Gas Struts (Force measured in Newton - N)		
<input type="checkbox"/> 2002150	65 lbs to 90 lbs (150 N)	<input type="checkbox"/> 2002200 90 lbs to 140 lbs (200 N)
<input type="checkbox"/> 2002300	140 lbs to 200 lbs (300 N)	<input type="checkbox"/> 2002400 200 lbs to 250 lbs (400 N)

Seat Width		
<b>Standard Seat Widths</b>		
<input type="checkbox"/> 2003012	12" Seat Width	<input type="checkbox"/> 2003013 13" Seat Width
<input type="checkbox"/> 2003015	15" Seat Width	<input type="checkbox"/> 2003016 16" Seat Width
<input type="checkbox"/> 2003018	18" Seat Width	<input type="checkbox"/> 2003019 19" Seat Width
<input type="checkbox"/> 2003014	14" Seat Width	<input type="checkbox"/> 2003017 17" Seat Width

Seat Depth		
<b>Standard Seat Depths</b>		
<input type="checkbox"/> 2004014	14" Seat Depth	<input type="checkbox"/> 2004015 15" Seat Depth
<input type="checkbox"/> 2004017	17" Seat Depth	<input type="checkbox"/> 2004018 18" Seat Depth
<input type="checkbox"/> 2004020	20" Seat Depth	<input type="checkbox"/> 2004016 16" Seat Depth
		<input type="checkbox"/> 2004019 19" Seat Depth

Frame Colour		
<b>PDG Premium Colours</b>		
<input type="checkbox"/> 2005000	Tracer Red	<input type="checkbox"/> 2005001 Gloss Black
<input type="checkbox"/> 2005003	Black Shimmer	<input type="checkbox"/> 2005002 Silver
<input type="checkbox"/> 2005006	Blue Shimmer	<input type="checkbox"/> 2005004 Orange Fusion
<input type="checkbox"/> 2005009	Geyser Green	<input type="checkbox"/> 2005007 Gloss White
<input type="checkbox"/> 2005012	Bubble Gum Pink	<input type="checkbox"/> 2005010 Lollipop Blue
		<input type="checkbox"/> 2005008 Dark Titanium
		<input type="checkbox"/> 2005011 Lollipop Red
		<input type="checkbox"/> 2005013 Berry Purple
		<input type="checkbox"/> 2005014 Safety Yellow
<input type="checkbox"/> 2005200	<b>Custom Colour</b>	Contact PDG for specifications. <span style="float: right;">\$ 300</span>

Arm Type		
<input type="checkbox"/> 2006000	<b>No Armrest</b>	
<input type="checkbox"/> 2006001	<b>Rigid Armrest</b>	Not compatible with Push to Lock, Fold Down Back. <span style="float: right;">\$ 150</span>

Seat Option		
<input type="checkbox"/> 2008110	<b>Standard Seat Sling</b>	
<input type="checkbox"/> 2008100	<b>No Seat Sling</b>	Third-party Seat Sling must be installed prior to wheelchair use.

## Cushion

2009000 **No Cushion**

## Seat Height

### 24" Wheel Size

2010018 18" STF  2010019 19" STF  2010020 20" STF

### 25" Wheel Size

Not compatible with 4" Caster with Frog Leg Uni-tine Fork.

2010118 18.5" STF  2010119 19.5" STF  2010120 20.5" STF

### 26" Wheel Size

Not compatible with Frog Leg Uni-tine Fork.

2010219 19" STF  1610220 20" STF  1610221 21" STF

## Casters

### Composite

2011040 4" x 1.0" Composite Sport  2011042 4" x 1.5" Composite Softroll  2011052 5" x 1.5" Composite Softroll

### Aluminium

2011140 4" x 1.0" Aluminum Sport  2011142 4" x 1.5" Aluminum Softroll  2011152 5" x 1.5" Aluminum Softroll \$ 125

### Caster Fork & Fork Suspension

2011410 Standard Aluminum Fork

2011420 Frog Legs Uni-tine Fork \$ 195

2011421 Frog Legs Ultra-Sport Suspension Forks \$ 450  
Compatible with 4" Caster only.

## Rear Wheels

### Spinergey Wheels

2012310 Spinergey Wire

2012311 Spinergey SPOX Everyday \$ 745

2012312 Spinergey SPOX Sport \$ 745

2012313 Spinergey LX Everyday \$ 895

2012314 Spinergey LX Sport \$ 895

2012315 **Spinergey Lite Extreme LX Flexirims** \$ 1,245  
LX Wheels only. Not compatible with 26" Wheel size.



### Select Wheel Camber

2012400 0° Camber  2012402 2° Camber

2012404 4° Camber  2012406 6° Camber

## Axle Type

2013000 **Quick Release Rear Axles**

## Rear Tires

2014110 **Primo Pneumatic Tires**


2014120 **Schwalbe Marathon Plus Evolution**  \$ 125

## Hand Rims

2015000 **No Hand Rims**

2015100 **Aluminum Clear Anodized**

2015110 **Plastic Coated** \$ 135

2015130 **BBraver Rough Rider**  \$ 750

1 →  Black 24"-26" Wheel size  White 24"-25" Wheel size  Orange 24"-25" Wheel size

2 →  Long  Short

## Wheel Locks

- 2016010 **Push to Lock** Not compatible with Rigid Armrest.
- 2016090 **Composite Scissor Locks**

## Anti-Tippers

- 2017010 **Rear Anti-Tippers** \$ 145

## Back Options

- 2019020 **Cane Style Back Post** \$ 195
- 2019221 **Rigid Back**
- 2019220 **Fold Down Back** Not compatible with Rigid Armrest. \$ 150
- Seat Back Height**
- 2019300 Short (8.5"-12") - Rigidizer Bar Height = 5.5"
- 2019301 Medium (10.5"-16") - Rigidizer Bar Height = 7.5"
- 2019302 Tall (16.5"-20") - Rigidizer Bar Height = 9.5"
- 2019310 **Custom Rigidizer Bar Height** \_\_\_\_\_ Must be at least 4.5' above Seat Sling. \$ 150
- 2019311 **Deep Rigidizer Bar** For installing third-party Backrest. \$ 150

## Back Upholstery

- 2020010 **Standard Back Upholstery**
- 2020000 **No Back Upholstery** Third-party Back Upholstery must be installed prior to wheelchair use.

## Footrest

- 2025010 **Standard Footrest**
- 2025011 **Angle Adjustable High Mount**
- 2026020 **Aluminum Angle Adjustable Foot Plate**
- 2026060 **Flip-up Foot Plate**
- 2026075 **Fixed Heel Retention Foot Plate** \$ 265

## Calf Straps

- 2027050 **Bodypoint - Padded Calf Strap** \$ 105

## Positioning Straps

- Bodypoint - Hip Belts**  
Push button with flat-mount end fitting. Recommended for mounting through slots in side guards, to allow "anti-shear" feature to be used.
- 2029010 Center-Pull Auto Style \$ 105
- 2029015 Center-Pull Auto Style, Padded \$ 155
- 2029200 **Bodypoint - EvoFlex Pelvic Stabilizer** \$ 275  
Push button recommended for when anti-shear feature is not used. Mounted directly to the seat tubes using the included seat tube mounting bracket.

## Side Guards

- 2030400 **No Side Guards**
- Non-Lifting Side Guards**
- 2030410 Rigid - Small  2030411 Rigid - Regular \$ 295
- 2030415 Folding - Regular \$ 350
- Lifting (Anti-Shear) Side Guards**
- 2030420 Rigid - Regular \$ 600
- 2030425 Folding - Regular  2030426 Folding - Large Not compatible with Rigid Armrest. \$ 650

### Miscellaneous Options

<input type="checkbox"/>	2030450	<b>Alber Wheel Mount Adapter</b>		\$ 600
<b>Quad Triggers</b> Not compatible with Push to Lock Wheel Locks				
<input type="checkbox"/>	2030461	Inboard Mounting		\$ 115
<input type="checkbox"/>	2030462	Outboard Mounting		\$ 115
<input type="checkbox"/>	2030463	Thumb Loops		\$ 130

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [ElevationbyPDG.com](http://ElevationbyPDG.com) to learn more about the **Elevation Ultralight Manual Wheelchair**



  PDG Mobility

