

Customer and Order Information		
Purchase Order *	Tag	v.070118
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type	
<input checked="" type="checkbox"/> 1701000 Fuze JR 50° Pediatric Manual Tilt Wheelchair	\$ 3,595

Frame Options	
<input checked="" type="checkbox"/> 1702025 Adjustable/Growing Base and Seat Frame 180 lbs Weight Capacity.	
<input type="checkbox"/> 1702300 Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

Seat Width		Overall chair width is +8.5"
Adjustable Seat Widths		
<input type="checkbox"/> 1703112 12"-15" Range - Set at 12"	<input type="checkbox"/> 1703113 12"-15" Range - Set at 13"	<input type="checkbox"/> 1703114 12"-15" Range - Set at 14"
<input type="checkbox"/> 1703114 12"-15" Range - Set at 15"		
<input type="checkbox"/> 1703215 15"-18" Range - Set at 15"	<input type="checkbox"/> 1703216 15"-18" Range - Set at 16"	<input type="checkbox"/> 1703217 15"-18" Range - Set at 17"
<input type="checkbox"/> 1703218 15"-18" Range - Set at 18"		

Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 1704013 13" Seat Depth	<input type="checkbox"/> 1704014 14" Seat Depth
<input type="checkbox"/> 1704015 15" Seat Depth	<input type="checkbox"/> 1704016 16" Seat Depth
<input type="checkbox"/> 1704017 17" Seat Depth	<input type="checkbox"/> 1704018 18" Seat Depth

Frame Colour	
PDG Premium Colours	
<input type="checkbox"/> 1705000 Tracer Red	<input type="checkbox"/> 1705001 Gloss Black
<input type="checkbox"/> 1705002 Silver	<input type="checkbox"/> 1705003 Black Shimmer
<input type="checkbox"/> 1705004 Orange Fusion	<input type="checkbox"/> 1705005 Karl's Imagination
<input type="checkbox"/> 1705006 Blue Shimmer	<input type="checkbox"/> 1705007 Gloss White
<input type="checkbox"/> 1705008 Dark Titanium	<input type="checkbox"/> 1705009 Geyser Green
<input type="checkbox"/> 1705010 Lollypop Blue	<input type="checkbox"/> 1705011 Lollypop Red
<input type="checkbox"/> 1705012 Bubble Gum Pink	<input type="checkbox"/> 1705013 Berry Purple
<input type="checkbox"/> 1705014 Safety Yellow	
<input type="checkbox"/> 1705200 Custom Colour	Contact PDG for specifications.

Arm Type	
<input type="checkbox"/> 1706000 No Armrest	
<input type="checkbox"/> 1706100 Single Post Fixed Height (12")	
<input type="checkbox"/> 1706111 Single Post Adjustable Height (8"-11")	\$ 255
<input type="checkbox"/> 1506200 Dual Post Fixed Height (12")	Full Length Armpad only.
Dual Post Adjustable Height	
<input type="checkbox"/> 1706211 8"-11"	<input type="checkbox"/> 1706214 10"-14"
Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.	
<input type="checkbox"/> 1706308 Set at 8"	<input type="checkbox"/> 1706309 Set at 9"
<input type="checkbox"/> 1706310 Set at 10"	<input type="checkbox"/> 1706311 Set at 11"
<input type="checkbox"/> 1706312 Set at 12"	<input type="checkbox"/> 1706313 Set at 13"

Armrests

<input type="checkbox"/>	1707095	Curved Armpads	Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.		
Regular Armpads					
<input type="checkbox"/>	1707100	Desk Length - 2" x 10"	<input type="checkbox"/>	1707101	Full Length - 2" x 14"
<input type="checkbox"/>	1707102	Full Length Premium Grade - 2" x 14"			\$ 35
Wide Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707110	Desk Length - 3" x 10"	<input type="checkbox"/>	1707111	Full Length - 3" x 12"
					\$ 115
Gel Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707120	Full Length - 2" x 12"	<input type="checkbox"/>	1707121	Full Length - 3.5" x 12"
					\$ 125
Contour Arm Trough with Hand Pad					
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707130	Right (4.5" x 18")	<input type="checkbox"/>	1707131	Left (4.5" x 18")
					\$ 125

Seat Pan

<input type="checkbox"/>	1708000	No Seat Pan	
<input type="checkbox"/>	1708001	Solid Aluminum Seat Pan	\$ 120

Cushion

<input type="checkbox"/>	1709000	No Cushion	
nxt BioFit			
<input type="checkbox"/>	1709101	12"-18" width, 12"-18" depth	\$ 595
Cushion Size Available in even number sizes only.			
	1709400	Specify Size (Width" x Depth")	_____



Seat Height

Select desired Seat-to-Floor Height

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters							
12"	13"	14"						
16"		14"	15"	16"				
20"				16"	17"			
22"					17"			

RW	6" Casters & 6" Forks							
12"		14"				18"		
16"		14"	15"	16"				
20"				16"	17"	18"		
22"					17"	18"		
24"						18"		

RW	6" Casters & 8" Forks							
12"					18"	19"		
16"			15"	16"				
20"				16"	17"	18"		
22"					17"	18"	19"	
24"						18"	19"	

RW	8" Casters & 8" Forks							
12"						18"	19"	20"
16"				16"				20"
20"				16"	17"	18"		
22"					17"	18"	19"	
24"						18"	19"	20"

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

<input type="checkbox"/>	1711051	5" x 1.25" Caster size	
<input type="checkbox"/>	1711061	6" x 1.25" Caster size	
<input type="checkbox"/>	1711063	6" x 1.50" Caster size	\$ 85
<input type="checkbox"/>	1711081	8" x 1.25" Caster size	
<input type="checkbox"/>	1711084	8" x 2.00" Caster size	\$ 135

Pneumatic

<input type="checkbox"/>	1711161	6" x 1.25" Caster size	\$ 110
<input type="checkbox"/>	1711181	8" x 1.25" Caster size	\$ 120
<input type="checkbox"/>	1711184	8" x 2.00" Caster size	\$ 150

Rear Wheels

Composite								
<input type="checkbox"/>	1712012	12" Wheel size	<input type="checkbox"/>	1712016	16" Wheel size	<input type="checkbox"/>	1712020	20" Wheel size
<input type="checkbox"/>	1712022	22" Wheel size	<input type="checkbox"/>	1712024	24" Wheel size			
Spoke								\$ 75
<input type="checkbox"/>	1712120	20" Wheel size	<input type="checkbox"/>	1712122	22" Wheel size	<input type="checkbox"/>	1712124	24" Wheel size

Axle Type

<input type="checkbox"/>	1713000	Quick Release Rear Axles
<input type="checkbox"/>	1713010	Fixed Rear Axles

Rear Tires

<input type="checkbox"/>	1714000	Urethane	
<input type="checkbox"/>	1714010	Pneumatic	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. \$ 60
<input type="checkbox"/>	1714020	Pneumatic with Flat-Free Inserts	\$ 95
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.			

Hand Rims

<input type="checkbox"/>	1715000	No Hand Rims	
<input type="checkbox"/>	1715010	Standard	20", 22", and 24" Rear Wheels only.
<input type="checkbox"/>	1715020	Plastic Coated	20", 22", and 24" Rear Wheels only. \$ 135

Wheel Locks

<input type="checkbox"/>	1716005	Foot Operated Wheel Lock	\$ 200
<input type="checkbox"/>	1716010	Push to Lock	
<input type="checkbox"/>	1716020	Pull to Lock	
<input type="checkbox"/>	1716030	Attendant Wheel Lock (Push to Lock)	Not compatible with Wheel Lock Extensions. \$ 85
<input type="checkbox"/>	1716100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock. \$ 60

Anti-Tippers

<input checked="" type="checkbox"/>	1717010	Rear Anti-Tippers	Required Safety Option. \$ 60
-------------------------------------	---------	--------------------------	-------------------------------

Back Post Style

Push Bar Back					
<input type="checkbox"/>	1718010	Straight			
<input type="checkbox"/>	1718030	8° Bend	\$ 75		
Cane Handle Back Not compatible with Headrest Mounting Plate.					
<input type="checkbox"/>	1718020	Straight			
<input type="checkbox"/>	1718040	8° Bend	\$ 75		
Back Post Height Seat Pan to Handle.					
<input type="checkbox"/>	1718420	20" Height	<input type="checkbox"/>	1718425	25" Height

Back Options

Angle/Height Adjustable Push Bar Handle			Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150		
<input type="checkbox"/>	1719010	Button on Outside	<input type="checkbox"/>	1719011	Button on Inside
<input type="checkbox"/>	1719050	Removable Cane Handle Push Bar	Not compatible with Push Bar Back. \$ 140		

See next page for more Back Options >

Seat to Back Angle			
<input type="checkbox"/> 1719090	Set at 90°	<input type="checkbox"/> 1719096	Set at 96°
<input type="checkbox"/> 1719108	Set at 108°	<input type="checkbox"/> 1719114	Set at 114°
<input type="checkbox"/> 1719220	Fold Down Back	<input type="checkbox"/> 1719102	Set at 102°
<input type="checkbox"/> 1719221	Non Fold Down Back	<input type="checkbox"/> 1719120	Set at 120°
			\$ 75

Back Upholstery	
<input checked="" type="checkbox"/> 1720000	No Back Upholstery

Backrest	
<input type="checkbox"/> 1521000	No Backrest
nxt Xtend	
<input type="checkbox"/> 1521100	16"-18" width, 16"-20" height \$ 690
nxt Xtend Deep	
<input type="checkbox"/> 1521101	16"-16" width, 16"-20" height \$ 690
Backrest Size Available in even number sizes only.	
1521400	Specify Size (Width" x Height") _____

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/> 1723011	Trigger on Right Back Post	<input type="checkbox"/> 1723012	Trigger on Left Back Post
Tilt Limit			
<input type="checkbox"/> 1723220	Set at 20°	<input type="checkbox"/> 1723225	Set at 25°
<input type="checkbox"/> 1723235	Set at 35°	<input type="checkbox"/> 1723240	Set at 40°
<input type="checkbox"/> 1723250	Set at 50°	<input type="checkbox"/> 1723230	Set at 30°
<input type="checkbox"/> 1723245	Set at 45°		
<input type="checkbox"/> 1723410	Reclining Backrest (Manual)	\$ 875	
30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Requires 25" Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15" or 16" Seat Depth.			

Headrest					
<input type="checkbox"/> 1724000	No Headrest				
<input type="checkbox"/> 1724010	Headrest Mounting Plate \$ 75				
Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.					
<input type="checkbox"/> 1724100	nxt MultiFit Head Support Hardware \$ 235				
nxt MultiFit Head Support Pad \$ 245					
<input type="checkbox"/> 1724105	Regular (9" x 5")	<input type="checkbox"/> 1724106	Large (10" x 5")	<input type="checkbox"/> 1724107	X-Large (14" x 5")

Front Rigging			
<input type="checkbox"/> 1725000	No Front Rigging	<input type="checkbox"/> 1725001	Sector Blocks only
PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.			
<input type="checkbox"/> 1725410	Fixed 70°	\$ 95	
<input type="checkbox"/> 1725420	Swing In/Away Universal Angle Adjustable (0°-110°)	Set at (use increments of 10°): _____	\$ 195
<input type="checkbox"/> 1725430	Swing In/Away Elevating Legrests	High knee pivot. Includes reinforced premium Calf Pad mount.	\$ 495
Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.			
<input type="checkbox"/> 1725140	Swing Away 70°		
<input type="checkbox"/> 1525160	Swing Away 90°	Not compatible with 8" Casters.	\$ 55
<input type="checkbox"/> 1725200	Swing Away Elevating Legrest	\$ 375	
<input type="checkbox"/> 1725300	1" Offset Hangers for Swing Away	\$ 75	

Foot Plates

<input type="checkbox"/>	1726010	Composite Not available on 12" Seat Width.	
<input type="checkbox"/>	1726020	Aluminum Angle Adjustable	\$ 155
<input type="checkbox"/>	1726040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1726050	Interlocking Angle Adjustable Seat Widths of 14"-18" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1727015	Swing Away 70° Calf Pad Not compatible with 70° Gel Padding, Seat Widths of 12" and 13".	\$ 135
<input type="checkbox"/>	1727020	Swing Away 70° Gel Padding Not compatible with 70° Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1727030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1727040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	1727050	Bodypoint - Padded Calf Strap	\$ 105
<input type="checkbox"/>	1727065	Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.	\$ 380

Extension Tubes

<input type="checkbox"/>	1728010	Short 14"-18" for Easy-Out Fixed & SA Angle Adjustable 11"-15" for Easy-Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1728020	Long 18"-24" for Easy-Out Fixed & SA Angle Adjustable 15"-21" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
<input type="checkbox"/>	1728200	For 90° FR (11"-16") Not compatible with Easy-Out, Standard 70° or ELR	\$ 75

Positioning Straps

		Bodypoint - Pediatric Hip Belts	
<input type="checkbox"/>	1729017	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1729022	4-Point, Padded	\$ 175
<input type="checkbox"/>	1729030	Bodypoint - Chest Strap	\$ 85



Miscellaneous Options

<input type="checkbox"/>	1730160	Clear Lap Tray (Slide On) Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	\$ 490
<input type="checkbox"/>	1730170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/fuze-T50Jr to learn more about the **Fuze JR 50° Manual Tilt Wheelchair**

