

Customer and Order Information	
Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address * <span style="float: right;">*required</span>	

Frame Type	
<input checked="" type="checkbox"/> 1701000 <b>Fuze JR 50° Pediatric Manual Tilt Wheelchair</b>	\$ 3,595

Frame Options	
<input checked="" type="checkbox"/> 1702025 <b>Adjustable/Growing Base and Seat Frame</b> 180 lbs Weight Capacity.	
<input type="checkbox"/> 1702300 <b>Transit Tie Down System (TTDS)</b>	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.	

Seat Width		Overall chair width is +8.5"
<b>Adjustable Seat Widths</b>		
<input type="checkbox"/> 1703112 12"-15" Range - Set at 12"	<input type="checkbox"/> 1703113 12"-15" Range - Set at 13"	<input type="checkbox"/> 1703114 12"-15" Range - Set at 14"
<input type="checkbox"/> 1703114 12"-15" Range - Set at 15"		
<input type="checkbox"/> 1703215 15"-18" Range - Set at 15"	<input type="checkbox"/> 1703216 15"-18" Range - Set at 16"	<input type="checkbox"/> 1703217 15"-18" Range - Set at 17"
<input type="checkbox"/> 1703218 15"-18" Range - Set at 18"		

Seat Depth	
<b>Standard Seat Depths (Adjustable)</b>	
<input type="checkbox"/> 1704013 13" Seat Depth	<input type="checkbox"/> 1704014 14" Seat Depth
<input type="checkbox"/> 1704015 15" Seat Depth	<input type="checkbox"/> 1704016 16" Seat Depth
<input type="checkbox"/> 1704017 17" Seat Depth	<input type="checkbox"/> 1704018 18" Seat Depth

Frame Colour	
<b>PDG Premium Colours</b>	
<input type="checkbox"/> 1705000 Tracer Red	<input type="checkbox"/> 1705001 Gloss Black
<input type="checkbox"/> 1705002 Silver	<input type="checkbox"/> 1705003 Black Shimmer
<input type="checkbox"/> 1705004 Orange Fusion	<input type="checkbox"/> 1705005 Karl's Imagination
<input type="checkbox"/> 1705006 Blue Shimmer	<input type="checkbox"/> 1705007 Gloss White
<input type="checkbox"/> 1705008 Dark Titanium	<input type="checkbox"/> 1705009 Geyser Green
<input type="checkbox"/> 1705010 Lollypop Blue	<input type="checkbox"/> 1705011 Lollypop Red
<input type="checkbox"/> 1705012 Bubble Gum Pink	<input type="checkbox"/> 1705013 Berry Purple
<input type="checkbox"/> 1705014 Safety Yellow	
<input type="checkbox"/> 1705200 <b>Custom Colour</b>	Contact PDG for specifications.

Arm Type	
<input type="checkbox"/> 1706000 <b>No Armrest</b>	
<input type="checkbox"/> 1706100 <b>Single Post Fixed Height (12")</b>	
<input type="checkbox"/> 1706111 <b>Single Post Adjustable Height (8"-11")</b>	\$ 255
<input type="checkbox"/> 1506200 <b>Dual Post Fixed Height (12")</b>	Full Length Armpad only.
<b>Dual Post Adjustable Height</b>	
<input type="checkbox"/> 1706211 8"-11"	\$ 255
<input type="checkbox"/> 1706214 10"-14"	
<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.	
<input type="checkbox"/> 1706308 Set at 8"	\$ 255
<input type="checkbox"/> 1706309 Set at 9"	
<input type="checkbox"/> 1706310 Set at 10"	
<input type="checkbox"/> 1706311 Set at 11"	
<input type="checkbox"/> 1706312 Set at 12"	
<input type="checkbox"/> 1706313 Set at 13"	

## Armrests

<input type="checkbox"/>	1707095	<b>Curved Armpads</b>	Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.		
<b>Regular Armpads</b>					
<input type="checkbox"/>	1707100	Desk Length - 2" x 10"	<input type="checkbox"/>	1707101	Full Length - 2" x 14"
<input type="checkbox"/>	1707102	Full Length Premium Grade - 2" x 14"			\$ 35
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707110	Desk Length - 3" x 10"	<input type="checkbox"/>	1707111	Full Length - 3" x 12"
					\$ 115
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707120	Full Length - 2" x 12"	<input type="checkbox"/>	1707121	Full Length - 3.5" x 12"
					\$ 125
<b>Contour Arm Trough with Hand Pad</b>					
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707130	Right (4.5" x 18")	<input type="checkbox"/>	1707131	Left (4.5" x 18")
					\$ 125

## Seat Pan

<input type="checkbox"/>	1708000	<b>No Seat Pan</b>	
<input type="checkbox"/>	1708001	<b>Solid Seat Pan</b>	\$ 120

## Cushion

<input type="checkbox"/>	1709000	<b>No Cushion</b>	
<b>nxt BioFit</b>			
<input type="checkbox"/>	1709101	12"-18" width, 12"-18" depth	\$ 595
<b>Cushion Size</b> Available in even number sizes only.			
1709400	Specify Size (Width" x Depth") _____		



## Seat Height

### Select desired Seat-to-Floor Height

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 8" Forks					
12"					18"	19"
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	8" Casters & 8" Forks					
12"					18"	19"
16"			16"			20"
20"			16"	17"	18"	
22"				17"	18"	19"
24"					18"	19"

Wheel size \_\_\_\_\_

Caster size \_\_\_\_\_

Fork size \_\_\_\_\_

STF height \_\_\_\_\_

## Casters

### Urethane

<input type="checkbox"/>	1711051	5" x 1.25" Caster size	
<input type="checkbox"/>	1711061	6" x 1.25" Caster size	
<input type="checkbox"/>	1711063	6" x 1.50" Caster size	\$ 85
<input type="checkbox"/>	1711081	8" x 1.25" Caster size	
<input type="checkbox"/>	1711084	8" x 2.00" Caster size	\$ 135

### Pneumatic

<input type="checkbox"/>	1711161	6" x 1.25" Caster size	\$ 110
<input type="checkbox"/>	1711181	8" x 1.25" Caster size	\$ 120
<input type="checkbox"/>	1711184	8" x 2.00" Caster size	\$ 150

## Rear Wheels

### Composite

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 1712012 12" Wheel size | <input type="checkbox"/> 1712016 16" Wheel size | <input type="checkbox"/> 1712020 20" Wheel size |
| <input type="checkbox"/> 1712022 22" Wheel size | <input type="checkbox"/> 1712024 24" Wheel size |   |

### Spoke

- |   |   |   |       |
|---|---|---|-------|
| <input type="checkbox"/> 1712120 20" Wheel size | <input type="checkbox"/> 1712122 22" Wheel size | <input type="checkbox"/> 1712124 24" Wheel size | \$ 75 |
|---|---|---|-------|

## Axle Type

- |  |
|--|
| <input type="checkbox"/> 1713000 <b>Quick Release Rear Axles</b> |
| <input type="checkbox"/> 1713010 <b>Fixed Rear Axles</b>         |

## Rear Tires

- |  |  |       |
|--|--|-------|
| <input type="checkbox"/> 1714000 <b>Urethane</b>                         |  |       |
| <input type="checkbox"/> 1714010 <b>Pneumatic</b>                        | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.                                   | \$ 60 |
| <input type="checkbox"/> 1714020 <b>Pneumatic with Flat-Free Inserts</b> | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. | \$ 95 |

## Hand Rims

- |  |                                     |       |
|--|-------------------------------------|-------|
| <input type="checkbox"/> 1715000 <b>No Hand Rims</b>   |                                     |       |
| <input type="checkbox"/> 1715010 <b>Standard</b>       | 20", 22", and 24" Rear Wheels only. |       |
| <input type="checkbox"/> 1715020 <b>Plastic Coated</b> | 20", 22", and 24" Rear Wheels only. | \$ 95 |

## Wheel Locks

- |   |  |        |
|---|--|--------|
| <input type="checkbox"/> 1716005 <b>Foot Operated Wheel Lock</b>            |  | \$ 200 |
| <input type="checkbox"/> 1716010 <b>Push to Lock</b>                        |  |        |
| <input type="checkbox"/> 1716020 <b>Pull to Lock</b>                        |  |        |
| <input type="checkbox"/> 1716030 <b>Attendant Wheel Lock (Push to Lock)</b> | Not compatible with Wheel Lock Extensions. | \$ 85  |
| <input type="checkbox"/> 1716100 <b>Wheel Lock Extensions</b>               | Not compatible with Attendant Wheel Lock.  | \$ 60  |

## Anti-Tippers

- |  |                         |       |
|--|-------------------------|-------|
| <input checked="" type="checkbox"/> 1717010 <b>Rear Anti-Tippers</b> | Required Safety Option. | \$ 60 |
|--|-------------------------|-------|

## Back Post Style

### Push Bar Back

- |   |  |       |
|---|--|-------|
| <input type="checkbox"/> 1718010 Straight |  |       |
| <input type="checkbox"/> 1718030 8° Bend  |  | \$ 75 |

### Cane Handle Back

Not compatible with Headrest Mounting Plate.

- |   |  |       |
|---|--|-------|
| <input type="checkbox"/> 1718020 Straight |  |       |
| <input type="checkbox"/> 1718040 8° Bend  |  | \$ 75 |

### Back Post Height

Seat Pan to Handle.

- |   |   |
|---|---|
| <input type="checkbox"/> 1718420 20" Height | <input type="checkbox"/> 1718425 25" Height |
|---|---|


## Back Options

<input type="checkbox"/> 1719010 <b>Angle/Height Adjustable Push Bar Handle</b> Button on Outside		<input type="checkbox"/> 1719011 Button on Inside		Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150	
<input type="checkbox"/> 1719050 <b>Removable Cane Handle Push Bar</b>		Not compatible with Push Bar Back.				\$ 140	
<b>Seat to Back Angle</b>							
<input type="checkbox"/> 1719090	Set at 90°	<input type="checkbox"/> 1719096	Set at 96°	<input type="checkbox"/> 1719102	Set at 102°		
<input type="checkbox"/> 1719108	Set at 108°	<input type="checkbox"/> 1719114	Set at 114°	<input type="checkbox"/> 1719120	Set at 120°		
<input type="checkbox"/> 1719220 <b>Fold Down Back</b>							
<input type="checkbox"/> 1719221 <b>Non Fold Down Back</b>							\$ 75

## Back Upholstery

<input checked="" type="checkbox"/> 1720000 <b>No Back Upholstery</b>
---


## Backrest

<input type="checkbox"/> 1521000 <b>No Backrest</b>			
<b>nxt Optima HA Active Contour</b>			
<input type="checkbox"/> 1521100	16"-18" width, 16"-20" height		\$ 690
<b>nxt Optima HA Plus Contour</b>			
<input type="checkbox"/> 1521101	16"-16" width, 16"-20" height		\$ 690
<b>Backrest Size</b> Available in even number sizes only.			
1521400	Specify Size (Width" x Height")	_____	

## Tilt/Recline Options

<b>Attendant Tilt only</b>				
<input type="checkbox"/> 1723011	Trigger on Right Back Post	<input type="checkbox"/> 1723012	Trigger on Left Back Post	
<b>Tilt Limit</b>				
<input type="checkbox"/> 1723220	Set at 20°	<input type="checkbox"/> 1723225	Set at 25°	
<input type="checkbox"/> 1723235	Set at 35°	<input type="checkbox"/> 1723240	Set at 40°	
<input type="checkbox"/> 1723250	Set at 50°	<input type="checkbox"/> 1723230	Set at 30°	
<input type="checkbox"/> 1723245	Set at 45°			
<input type="checkbox"/> 1723410 <b>Reclining Backrest (Manual)</b>				\$ 875
30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Requires 25" Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15" or 16" Seat Depth.				

## Headrest

<input type="checkbox"/> 1724000 <b>No Headrest</b>			
<input type="checkbox"/> 1724010 <b>Headrest Mounting Plate</b>	Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.		
<input type="checkbox"/> 1724100 <b>nxt MultiFit Head Support Hardware</b>		\$ 155	
<b>nxt MultiFit Head Support Pad</b>			
<input type="checkbox"/> 1724105	Regular (9" x 5")	<input type="checkbox"/> 1724106	
	Large (10" x 5")	<input type="checkbox"/> 1724107	
		X-Large (14" x 5")	\$ 195

## Front Rigging

<input type="checkbox"/> 1725000 <b>No Front Rigging</b>			
<input type="checkbox"/> 1725001 <b>Sector Blocks only</b>			
<input type="checkbox"/> 1725140 <b>Swing Away 70°</b>			
<input type="checkbox"/> 1725160 <b>Swing Away 90°</b>	Not compatible with 8" Casters.		
		\$ 55	
<input type="checkbox"/> 1725200 <b>Swing Away Elevating Legrest</b>			\$ 375

### Foot Plates

<input type="checkbox"/>	1726010	<b>Composite</b>	Not available on 12" Seat Width.	
<input type="checkbox"/>	1726020	<b>Aluminum Angle Adjustable</b>		\$ 155
<input type="checkbox"/>	1726040	<b>One Piece Flip-up</b>	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1726050	<b>Interlocking Angle Adjustable</b>	Seat Widths of 14"-18" only.	\$ 195

### Footrest/Legrest Accessories

<input type="checkbox"/>	1727015	<b>Swing Away 70° Calf Pad</b>	Not compatible with 70° Gel Padding, Seat Widths of 12" and 13".	\$ 135
<input type="checkbox"/>	1727020	<b>Swing Away 70° Gel Padding</b>	Not compatible with 70° Calf Pad.	\$ 100
<input type="checkbox"/>	1727030	<b>Swing Away ELR Gel Knee Pad</b>		\$ 100
<input type="checkbox"/>	1727040	<b>Adjustable Length Heel Loops</b>		\$ 50
<input type="checkbox"/>	1727050	<b>Bodypoint - Padded Calf Strap</b>		\$ 105
<input type="checkbox"/>	1727065	<b>Miller's Dynamic Footrest Mounts</b>		\$ 380

### Extension Tubes

		<b>Regular Extension Tubes</b>	Not compatible with Swing Away 90°.	
<input type="checkbox"/>	1728010	Standard (13"-17")		
<input type="checkbox"/>	1728020	Long (15"-19")		\$ 75

### Positioning Straps

		<b>Bodypoint - Pediatric Hip Belts</b>		
<input type="checkbox"/>	1729017	Center-Pull Auto Style, Padded		\$ 155
<input type="checkbox"/>	1729022	4-Point, Padded		\$ 175
<input type="checkbox"/>	1729030	<b>Bodypoint - Chest Strap</b>		\$ 85

### Miscellaneous Options

<input type="checkbox"/>	1730160	<b>Clear Lap Tray (Slide On)</b>	Not compatible with 8"- 11" Dual Post Arms, Desk Length Pads, Gel Arm Pads. Seat Widths of 15"-18" only.	\$ 490
<input type="checkbox"/>	1730170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.	\$ 10

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/fuze-T50Jr](http://pdgmobility.com/products/fuze-T50Jr) to learn more about the **Fuze JR 50° Manual Tilt Wheelchair**

