

Customer and Order Information		
Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type	
<input checked="" type="checkbox"/> 1701000 Fuze JR 50° Pediatric Manual Tilt Wheelchair	\$ 3,595

Frame Options	
<input checked="" type="checkbox"/> 1702025 Adjustable/Growing Base and Seat Frame 180 lbs Weight Capacity.	
<input type="checkbox"/> 1702300 Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.

Seat Width		Overall chair width is +8.5"
Adjustable Seat Widths		
<input type="checkbox"/> 1703112 12"-15" Range - Set at 12"	<input type="checkbox"/> 1703113 12"-15" Range - Set at 13"	<input type="checkbox"/> 1703114 12"-15" Range - Set at 14"
<input type="checkbox"/> 1703114 12"-15" Range - Set at 15"		
<input type="checkbox"/> 1703215 15"-18" Range - Set at 15"	<input type="checkbox"/> 1703216 15"-18" Range - Set at 16"	<input type="checkbox"/> 1703217 15"-18" Range - Set at 17"
<input type="checkbox"/> 1703218 15"-18" Range - Set at 18"		

Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 1704013 13" Seat Depth	<input type="checkbox"/> 1704014 14" Seat Depth
<input type="checkbox"/> 1704015 15" Seat Depth	<input type="checkbox"/> 1704016 16" Seat Depth
<input type="checkbox"/> 1704017 17" Seat Depth	<input type="checkbox"/> 1704018 18" Seat Depth

Frame Colour	
PDG Premium Colours	
<input type="checkbox"/> 1705000 Tracer Red	<input type="checkbox"/> 1705001 Gloss Black
<input type="checkbox"/> 1705002 Silver	<input type="checkbox"/> 1705003 Black Shimmer
<input type="checkbox"/> 1705004 Orange Fusion	<input type="checkbox"/> 1705005 Karl's Imagination
<input type="checkbox"/> 1705006 Blue Shimmer	<input type="checkbox"/> 1705007 Gloss White
<input type="checkbox"/> 1705008 Dark Titanium	<input type="checkbox"/> 1705009 Geyser Green
<input type="checkbox"/> 1705010 Lollypop Blue	<input type="checkbox"/> 1705011 Lollypop Red
<input type="checkbox"/> 1705012 Bubble Gum Pink	<input type="checkbox"/> 1705013 Berry Purple
<input type="checkbox"/> 1705014 Safety Yellow	
<input type="checkbox"/> 1705200 Custom Colour Contact PDG for specifications.	

Arm Type	
<input type="checkbox"/> 1706000 No Armrest	
<input type="checkbox"/> 1706100 Single Post Fixed Height (12")	
<input type="checkbox"/> 1706111 Single Post Adjustable Height (8"-11")	\$ 255
<input type="checkbox"/> 1506200 Dual Post Fixed Height (12") Full Length Armpad only.	
<input type="checkbox"/> 1706211 Dual Post Adjustable Height	\$ 255
<input type="checkbox"/> 1706211 8"-11"	<input type="checkbox"/> 1706214 10"-14"
<input type="checkbox"/> 1706308 Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.	\$ 255
<input type="checkbox"/> 1706309 Set at 8"	<input type="checkbox"/> 1706310 Set at 9"
<input type="checkbox"/> 1706311 Set at 11"	<input type="checkbox"/> 1706312 Set at 12"
<input type="checkbox"/> 1706313 Set at 10"	<input type="checkbox"/> 1706314 Set at 13"

Armrests

<input type="checkbox"/>	1707095	Curved Armpads	Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.		
Regular Armpads					
<input type="checkbox"/>	1707100	Desk Length - 2" x 10"	<input type="checkbox"/>	1707101	Full Length - 2" x 14"
<input type="checkbox"/>	1707102	Full Length Premium Grade - 2" x 14"			\$ 35
Wide Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707110	Desk Length - 3" x 10"	<input type="checkbox"/>	1707111	Full Length - 3" x 12"
					\$ 115
Gel Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707120	Full Length - 2" x 12"	<input type="checkbox"/>	1707121	Full Length - 3.5" x 12"
					\$ 125
Contour Arm Trough with Hand Pad					
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707130	Right (4.5" x 18")	<input type="checkbox"/>	1707131	Left (4.5" x 18")
					\$ 125

Seat Pan

<input type="checkbox"/>	1708000	No Seat Pan			
<input type="checkbox"/>	1708001	Solid Aluminum Seat Pan			\$ 120

Cushion

<input type="checkbox"/>	1709000	No Cushion			
nxt BioFit					
<input type="checkbox"/>	1709101	12"-18" width, 12"-18" depth			\$ 595
Cushion Size Available in even number sizes only.					
	1709400	Specify Size (Width" x Depth")	_____		

Seat Height

Select desired Seat-to-Floor Height

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 8" Forks					
12"					18"	19"
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	8" Casters & 8" Forks					
12"					18"	19"
16"			16"			20"
20"			16"	17"	18"	
22"				17"	18"	19"
24"					18"	19"

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

<input type="checkbox"/>	1711051	5" x 1.25" Caster size			
<input type="checkbox"/>	1711061	6" x 1.25" Caster size			
<input type="checkbox"/>	1711063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/>	1711081	8" x 1.25" Caster size			
<input type="checkbox"/>	1711084	8" x 2.00" Caster size			\$ 135

Pneumatic

<input type="checkbox"/>	1711161	6" x 1.25" Caster size			\$ 110
<input type="checkbox"/>	1711181	8" x 1.25" Caster size			\$ 120
<input type="checkbox"/>	1711184	8" x 2.00" Caster size			\$ 150

Rear Wheels

Composite

- | | | |
|---|---|---|
| <input type="checkbox"/> 1712012 12" Wheel size | <input type="checkbox"/> 1712016 16" Wheel size | <input type="checkbox"/> 1712020 20" Wheel size |
| <input type="checkbox"/> 1712022 22" Wheel size | <input type="checkbox"/> 1712024 24" Wheel size | |

Spoke

- | | | | |
|---|---|---|-------|
| <input type="checkbox"/> 1712120 20" Wheel size | <input type="checkbox"/> 1712122 22" Wheel size | <input type="checkbox"/> 1712124 24" Wheel size | \$ 75 |
|---|---|---|-------|

Axle Type

- | |
|--|
| <input type="checkbox"/> 1713000 Quick Release Rear Axles |
| <input type="checkbox"/> 1713010 Fixed Rear Axles |

Rear Tires

- | | | |
|--|--|-------|
| <input type="checkbox"/> 1714000 Urethane | | |
| <input type="checkbox"/> 1714010 Pneumatic | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 1714020 Pneumatic with Flat-Free Inserts | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. | \$ 95 |

Hand Rims

- | | | |
|--|-------------------------------------|--------|
| <input type="checkbox"/> 1715000 No Hand Rims | | |
| <input type="checkbox"/> 1715010 Standard | 20", 22", and 24" Rear Wheels only. | |
| <input type="checkbox"/> 1715020 Plastic Coated | 20", 22", and 24" Rear Wheels only. | \$ 135 |

Wheel Locks

- | | | |
|---|--|--------|
| <input type="checkbox"/> 1716005 Foot Operated Wheel Lock | | \$ 200 |
| <input type="checkbox"/> 1716010 Push to Lock | | |
| <input type="checkbox"/> 1716020 Pull to Lock | | |
| <input type="checkbox"/> 1716030 Attendant Wheel Lock (Push to Lock) | Not compatible with Wheel Lock Extensions. | \$ 85 |
| <input type="checkbox"/> 1716100 Wheel Lock Extensions | Not compatible with Attendant Wheel Lock. | \$ 60 |

Anti-Tippers

- | | | |
|--|-------------------------|-------|
| <input checked="" type="checkbox"/> 1717010 Rear Anti-Tippers | Required Safety Option. | \$ 60 |
|--|-------------------------|-------|

Back Post Style

Push Bar Back

- | | | |
|---|--|-------|
| <input type="checkbox"/> 1718010 Straight | | |
| <input type="checkbox"/> 1718030 8° Bend | | \$ 75 |

Cane Handle Back

Not compatible with Headrest Mounting Plate.

- | | | |
|---|--|-------|
| <input type="checkbox"/> 1718020 Straight | | |
| <input type="checkbox"/> 1718040 8° Bend | | \$ 75 |

Back Post Height

Seat Pan to Handle.

- | | |
|---|---|
| <input type="checkbox"/> 1718420 20" Height | <input type="checkbox"/> 1718425 25" Height |
|---|---|

Back Options

Angle/Height Adjustable Push Bar Handle

Not compatible with Headrest Mounting Plate or Cane Handle Back.

- | | | |
|--|---|--------|
| <input type="checkbox"/> 1719010 Button on Outside | <input type="checkbox"/> 1719011 Button on Inside | \$ 150 |
| <input type="checkbox"/> 1719050 Removable Cane Handle Push Bar | Not compatible with Push Bar Back. | \$ 140 |

See next page for more Back Options >

Seat to Back Angle			
<input type="checkbox"/> 1719090	Set at 90°	<input type="checkbox"/> 1719096	Set at 96°
<input type="checkbox"/> 1719108	Set at 108°	<input type="checkbox"/> 1719114	Set at 114°
<input type="checkbox"/> 1719220	Fold Down Back	<input type="checkbox"/> 1719102	Set at 102°
<input type="checkbox"/> 1719221	Non Fold Down Back	<input type="checkbox"/> 1719120	Set at 120°
			\$ 75

Back Upholstery	
<input checked="" type="checkbox"/> 1720000	No Back Upholstery

Backrest	
<input type="checkbox"/> 1521000	No Backrest
<input type="checkbox"/> 1521100	nxt Optima HA Active Contour 16"-18" width, 16"-20" height
<input type="checkbox"/> 1521101	nxt Optima HA Plus Contour 16"-16" width, 16"-20" height
<input type="checkbox"/> 1521400	Backrest Size Available in even number sizes only. Specify Size (Width" x Height") _____
\$ 690	
\$ 690	

Tilt/Recline Options	
Attendant Tilt only	
<input type="checkbox"/> 1723011	Trigger on Right Back Post
<input type="checkbox"/> 1723012	Trigger on Left Back Post
Tilt Limit	
<input type="checkbox"/> 1723220	Set at 20°
<input type="checkbox"/> 1723225	Set at 25°
<input type="checkbox"/> 1723230	Set at 30°
<input type="checkbox"/> 1723235	Set at 35°
<input type="checkbox"/> 1723240	Set at 40°
<input type="checkbox"/> 1723245	Set at 45°
<input type="checkbox"/> 1723410	Reclining Backrest (Manual)
30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Requires 25" Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15" or 16" Seat Depth.	
\$ 875	

Headrest	
<input type="checkbox"/> 1724000	No Headrest
<input type="checkbox"/> 1724010	Headrest Mounting Plate
Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	
<input type="checkbox"/> 1724100	nxt MultiFit Head Support Hardware
<input type="checkbox"/> 1724105	nxt MultiFit Head Support Pad
<input type="checkbox"/> 1724106	Regular (9" x 5")
<input type="checkbox"/> 1724107	Large (10" x 5")
<input type="checkbox"/> 1724108	X-Large (14" x 5")
\$ 75	
\$ 235	
\$ 245	

Front Rigging	
<input type="checkbox"/> 1725000	No Front Rigging
PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.	
<input type="checkbox"/> 1725410	Fixed 70°
<input type="checkbox"/> 1725420	Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____
<input type="checkbox"/> 1725430	Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.
Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.	
<input type="checkbox"/> 1725140	Swing Away 70°
<input type="checkbox"/> 1525160	Swing Away 90° Not compatible with 8" Casters.
<input type="checkbox"/> 1725200	Swing Away Elevating Legrest
<input type="checkbox"/> 1725300	1" Offset Hangers for Swing Away
Available in March \$ 95	
Available in April \$ 195	
Available in April \$ 495	
\$ 55	
\$ 375	
\$ 75	

Foot Plates

<input type="checkbox"/>	1726010	Composite	Not available on 12" Seat Width.	
<input type="checkbox"/>	1726020	Aluminum Angle Adjustable		\$ 155
<input type="checkbox"/>	1726040	One Piece Flip-up	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1726050	Interlocking Angle Adjustable	Seat Widths of 14"-18" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1727015	Swing Away 70° Calf Pad	Not compatible with 70° Gel Padding, Seat Widths of 12" and 13".	\$ 135
<input type="checkbox"/>	1727020	Swing Away 70° Gel Padding	Not compatible with 70° Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1727030	Swing Away ELR Gel Knee Pad		\$ 100
<input type="checkbox"/>	1727040	Adjustable Length Heel Loops		\$ 50
<input type="checkbox"/>	1727050	Bodypoint - Padded Calf Strap		\$ 105
<input type="checkbox"/>	1727065	Miller's Dynamic Footrest Mounts	Not compatible with Easy-Out Front Rigging.	\$ 380

Extension Tubes

<input type="checkbox"/>	1728010	Short	13"-17" for Easy-Out SA Angle Adjustable 10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1728020	Long	16"-22.5" for Easy-Out SA Angle Adjustable 12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
<input type="checkbox"/>	1728200	For 90° FR (11"-16")	Not compatible with Easy-Out, Standard 70° or ELR	\$ 75

Positioning Straps

		Bodypoint - Pediatric Hip Belts		
<input type="checkbox"/>	1729017	Center-Pull Auto Style, Padded		\$ 155
<input type="checkbox"/>	1729022	4-Point, Padded		\$ 175
<input type="checkbox"/>	1729030	Bodypoint - Chest Strap		\$ 85

Miscellaneous Options

<input type="checkbox"/>	1730160	Clear Lap Tray (Slide On)	Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	\$ 490
<input type="checkbox"/>	1730170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/fuze-T50Jr to learn more about the **Fuze JR 50° Manual Tilt Wheelchair**

