

### Customer and Order Information

Purchase Order *	Tag	v.0301
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

### Frame Type

<input checked="" type="checkbox"/> 1601000	<b>Fuze T20 20° Manual Tilt Wheelchair</b>	\$ 3,995
---	--	----------

### Frame Options

<input type="checkbox"/> 1602010	<b>Fixed Base and Seat Frame</b>	250 lbs Weight Capacity. Select Fixed Seat Width below.	
<input type="checkbox"/> 1602025	<b>Adjustable Base and Seat Frame</b>	250 lbs Weight Capacity. Select Adjustable Seat Width below.	\$ 100
<input type="checkbox"/> 1602300	<b>Transit Tie Down System (TTDS)</b>	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.	
<b>User Weight Range</b> Required for selection of Gas Struts (Force measured in Newton - N)			
<input type="checkbox"/> 1602406	Up to 150 lbs (600 N)	<input type="checkbox"/> 1602410	150 lbs to 200 lbs (1000 N)
		<input type="checkbox"/> 1602412	200 lbs to 250 lbs (1200 N)

### Seat Width

Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

<b>Fixed Seat Widths</b> Order with Fixed Base and Seat Frame.			
<input type="checkbox"/> 1603015	15" Seat Width	<input type="checkbox"/> 1603016	16" Seat Width
<input type="checkbox"/> 1603017	17" Seat Width	<input type="checkbox"/> 1603018	18" Seat Width
<input type="checkbox"/> 1603019	19" Seat Width	<input type="checkbox"/> 1603020	20" Seat Width
<input type="checkbox"/> 1603021	21" Seat Width	<input type="checkbox"/> 1603022	22" Seat Width
<b>Adjustable Seat Widths</b> Order with Adjustable Base and Seat Frame.			
<input type="checkbox"/> 1603113	13"-15" Range - Set at 13"	<input type="checkbox"/> 1603114	13"-15" Range - Set at 14"
<input type="checkbox"/> 1603115	13"-15" Range - Set at 15"	<input type="checkbox"/> 1603215	15"-20" Range - Set at 15"
<input type="checkbox"/> 1603216	15"-20" Range - Set at 16"	<input type="checkbox"/> 1603217	15"-20" Range - Set at 17"
<input type="checkbox"/> 1603218	15"-20" Range - Set at 18"	<input type="checkbox"/> 1603219	15"-20" Range - Set at 19"
<input type="checkbox"/> 1603220	15"-20" Range - Set at 20"	<input type="checkbox"/> 1603321	21"-22" Range - Set at 21"
<input type="checkbox"/> 1603322	21"-22" Range - Set at 22"		

### Seat Depth

<b>Standard Seat Depths (Adjustable)</b>			
<input type="checkbox"/> 1604015	15" Seat Depth	<input type="checkbox"/> 1604016	16" Seat Depth
<input type="checkbox"/> 1604017	17" Seat Depth	<input type="checkbox"/> 1604018	18" Seat Depth
<input type="checkbox"/> 1604019	19" Seat Depth	<input type="checkbox"/> 1604020	20" Seat Depth

### Frame Color

<b>PDG Premium Colors</b>			
<input type="checkbox"/> 1605000	Tracer Red	<input type="checkbox"/> 1605001	Gloss Black
<input type="checkbox"/> 1605002	Silver	<input type="checkbox"/> 1605003	Black Shimmer
<input type="checkbox"/> 1605004	Orange Fusion	<input type="checkbox"/> 1605005	Karl's Imagination
<input type="checkbox"/> 1605006	Blue Shimmer	<input type="checkbox"/> 1605007	Gloss White
<input type="checkbox"/> 1605008	Dark Titanium	<input type="checkbox"/> 1605009	Geyser Green
<input type="checkbox"/> 1605010	Lollipop Blue	<input type="checkbox"/> 1605011	Lollipop Red
<input type="checkbox"/> 1605012	Bubble Gum Pink	<input type="checkbox"/> 1605013	Berry Purple
<input type="checkbox"/> 1605014	Safety Yellow		
<input type="checkbox"/> 1605200	<b>Custom Color</b>	Contact PDG for specifications.	
			\$ 300

Arm Type		
<input type="checkbox"/>	1606000	<b>No Armrest</b>
<input type="checkbox"/>	1606100	<b>Single Post Fixed Height (12")</b>
	<b>Single Post Adjustable Height</b>	8"-11" not compatible with 22" and 24" Rear Wheels. <span style="float: right;">\$ 275</span>
<input type="checkbox"/>	1606111	8"-11" <input type="checkbox"/>
<input type="checkbox"/>	1606114	10"-14"
<input type="checkbox"/>	1606200	<b>Dual Post Fixed Height (12")</b> Full Length Armpad only.
	<b>Dual Post Adjustable Height</b>	8"-11" not compatible with 22" and 24" Rear Wheels. <span style="float: right;">\$ 275</span>
<input type="checkbox"/>	1606211	8"-11" <input type="checkbox"/>
<input type="checkbox"/>	1606214	10"-14" <input type="checkbox"/>
<input type="checkbox"/>	1606216	12"-16"
	<b>Locking Flip-up</b>	Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery. <span style="float: right;">\$ 295</span>
<input type="checkbox"/>	1606310	Set at 10" <input type="checkbox"/>
<input type="checkbox"/>	1606311	Set at 11" <input type="checkbox"/>
<input type="checkbox"/>	1606312	Set at 12" <input type="checkbox"/>
<input type="checkbox"/>	1606313	Set at 13" <input type="checkbox"/>
<input type="checkbox"/>	1606314	Set at 14" <input type="checkbox"/>

Armrests		
<input type="checkbox"/>	1607095	<b>Curved Armpads</b> Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests, 22" and 24" Rear Wheels.
	<b>Regular Armpads</b>	
<input type="checkbox"/>	1607100	Desk Length - 2" x 10" <input type="checkbox"/>
<input type="checkbox"/>	1607101	Full Length - 2" x 14"
	<b>Wide Armpads</b>	May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 200</span>
<input type="checkbox"/>	1607110	Desk Length - 3" x 10" <input type="checkbox"/>
<input type="checkbox"/>	1607111	Full Length - 3" x 12"
	<b>Gel Armpads</b>	May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 135</span>
<input type="checkbox"/>	1607120	Full Length - 2" x 12" <input type="checkbox"/>
<input type="checkbox"/>	1607121	Full Length - 3.5" x 12"
	<b>Contour Arm Trough with Hand Pad</b>	May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>
<input type="checkbox"/>	1607130	Right (4.5" x 18") <input type="checkbox"/>
<input type="checkbox"/>	1607131	Left (4.5" x 18")

Seat Pan		
<input type="checkbox"/>	1608000	<b>No Seat Pan</b>
<input type="checkbox"/>	1608001	<b>Solid Aluminum Seat Pan</b> <span style="float: right;">\$ 120</span>

Cushion		
	Refer to the enclosed Seating & Positioning order form for PDG's on-chair Cushion options. Standard PDG discounts apply.	

### Seat Height

#### Select desired Seat-to-Floor Height

##### Standard T20 Options

RW	5" Casters					
12"	13"					
16"	13"	14"	15"			
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	

RW	6" Casters & 6" Forks					
16"		14"	15"			
20"		14"	15"	16"	17"	
22"		14"	15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	6" Casters & 8" Forks					
16"			15"			
20"			15"	16"	17"	
22"			15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	8" Casters & 8" Forks					
20"				16"	17"	
22"				16"	17"	18"
24"				16"	17"	18"

##### Narrow T20 Options

Reduces the overall width of the Chair by 2". Recommended only if the width of the chair is a concern. Some STF combinations offer less than 20° of tilt, but all offer a minimum of 12°.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	8" Casters & 8" Forks					
16"				16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

- Standard  
 Narrow

Wheel size \_\_\_\_\_

Caster size \_\_\_\_\_

Fork size \_\_\_\_\_

STF height \_\_\_\_\_

### Casters

#### Urethane

- 1611051 5" x 1.25" Caster size
- 1611061 6" x 1.25" Caster size
- 1611063 6" x 1.50" Caster size \$ 85
- 1611081 8" x 1.25" Caster size
- 1611084 8" x 2.00" Caster size \$ 135

#### Pneumatic

- 1611161 6" x 1.25" Caster size \$ 110
- 1611181 8" x 1.25" Caster size \$ 120
- 1611184 8" x 2.00" Caster size \$ 150
- 1611260 **6" Soft Roll** \$ 95

### Rear Wheels

#### Composite

- 1612012 12" Wheel size
- 1612016 16" Wheel size
- 1612020 20" Wheel size
- 1612022 22" Wheel size
- 1612024 24" Wheel size

#### Spoke

- 1612120 20" Wheel size \$ 75
- 1612122 22" Wheel size
- 1612124 24" Wheel size
- 1612224 **24" Heavy Duty Spoke** \$ 350

### Axle Type

- 1613000 **Quick Release Rear Axles**
- 1613010 **Fixed Rear Axles**

Rear Tires			
<input type="checkbox"/>	1614000	<b>Urethane</b>	
<input type="checkbox"/>	1614010	<b>Pneumatic</b> May result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/>	1614020	<b>Pneumatic with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95

Hand Rims			
<input type="checkbox"/>	1615000	<b>No Hand Rims</b>	
<input type="checkbox"/>	1615010	<b>Standard</b> 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1615020	<b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only.	\$ 135

Wheel Locks			
<input type="checkbox"/>	1616010	<b>Push to Lock</b>	
<input type="checkbox"/>	1616020	<b>Pull to Lock</b>	
<input type="checkbox"/>	1616030	<b>Attendant Wheel Lock (Push to Lock)</b> Not compatible with Wheel Lock Extensions.	\$ 85
<input type="checkbox"/>	1616040	<b>User &amp; Attendant Wheel Lock</b>	\$ 125
<input type="checkbox"/>	1616100	<b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	\$ 60

Anti-Tippers			
<input checked="" type="checkbox"/>	1617010	<b>Rear Anti-Tippers</b> Required Safety Option.	\$ 60

Back Post Style			
<b>Push Bar Back</b>			
<input type="checkbox"/>	1618010	Straight	
<input type="checkbox"/>	1618030	8° Bend	\$ 75
<b>Cane Handle Back</b> Not compatible with Headrest Mounting Plate.			
<input type="checkbox"/>	1618020	Straight	
<input type="checkbox"/>	1618040	8° Bend	\$ 75
<b>Back Post Height</b> Seat Pan to Handle.			
<input type="checkbox"/>	1618420	20" Height	
<input type="checkbox"/>	1618425	25" Height	

Back Options			
<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.			\$ 150
<input type="checkbox"/>	1619010	Button on Outside	
<input type="checkbox"/>	1619011	Button on Inside	
<input type="checkbox"/>	1619050	<b>Removable Cane Handle Push Bar</b> Not compatible with Push Bar Back.	\$ 140
<b>Seat to Back Angle</b>			
<input type="checkbox"/>	1619090	Set at 90°	
<input type="checkbox"/>	1619096	Set at 96°	
<input type="checkbox"/>	1619102	Set at 102°	
<input type="checkbox"/>	1619108	Set at 108°	
<input type="checkbox"/>	1619114	Set at 114°	
<input type="checkbox"/>	1619120	Set at 120°	
<input type="checkbox"/>	1619220	<b>Fold Down Back</b>	

Back Upholstery			
<input type="checkbox"/>	1620000	<b>No Back Upholstery</b>	
<b>Standard</b> Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.			\$ 175
<input type="checkbox"/>	1620018	Set at 18"	
<input type="checkbox"/>	1620019	Set at 19"	
<input type="checkbox"/>	1620020	Set at 20"	
<input type="checkbox"/>	1620021	Set at 21"	
<input type="checkbox"/>	1620022	Set at 22"	
<input type="checkbox"/>	1620023	Set at 23"	
<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.			\$ 295
<input type="checkbox"/>	1620321	Set at 21"	
<input type="checkbox"/>	1620322	Set at 22"	
<input type="checkbox"/>	1620323	Set at 23"	

### Positioning Backrest

**Refer to the enclosed Seating & Positioning order form for PDG's on-chair Backrest options. Standard PDG discounts apply.**

### Tilt/Recline Options

#### Attendant Tilt only

- 1623011 Trigger on Right Back Post       1623012 Trigger on Left Back Post

### Headrest

**Refer to the enclosed Seating & Positioning order form for PDG's on-chair Headrest options. Standard PDG discounts apply.**

- 1624010 **Headrest Mounting Plate** \$ 95  
 Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.

### Front Rigging

- 1625000 **No Front Rigging**       1625001 **Sector Blocks only**
- PDG Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame. [pdgmobility.com/easy-out](http://pdgmobility.com/easy-out)
- 1625410 Fixed 70° \$ 95
- 1625420 Universal Angle Adjustable (0°-110°)      Set at angle: \_\_\_\_\_ (use increments of 10°) \$ 295  
 Adjust angle by loosening one bolt. Tool-free leg-length adjustment.
- Swing Away Front Rigging** Legacy Design. Not compatible with Easy-Out.
- 1625140 Swing Away 70°
- 1625160 Swing Away 90° Includes 11"-16" Extension Tubes. Not compatible with 8" Casters. \$ 75
- 1625200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375
- 1625300 1" Offset Hangers for Swing Away \$ 75
- 1625110 **Standard Plug In 70°** Legacy Design. Not compatible with Easy-Out.

### Foot Plates

- 1626010 **Composite**
- 1626020 **Aluminum Angle Adjustable** \$ 155
- 1626050 **Interlocking Angle Adjustable** Seat Widths of 14"-21" only. \$ 195

### Footrest/Legrest Accessories

- Calf Pad** Not compatible with Gel Padding.
- 1627012 For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount. \$ 200
- 1627010 For Swing Away Front Rigging \$ 135
- 1627020 **Gel Padding** Not compatible with Calf Pad or Easy-Out Front Rigging. \$ 125
- Gel Knee Pad**
- 1627032 For Easy-Out™ Front Rigging \$ 200
- 1627030 For Swing Away ELR Front Rigging \$ 100
- 1627040 **Adjustable Length Heel Loops** \$ 50
- 1627050 **Padded Calf Strap** \$ 195

### Extension Tubes

- Regular Extension Tubes** No selection required for 90° Front Rigging
- 1628010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging
- 1628020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging \$ 75

Positioning Straps		
	Refer to the enclosed Seating & Positioning order form for PDG's on-chair Belts/Straps options. Standard PDG discounts apply.	
<input type="checkbox"/>	1629005 <b>Lap Belt Mount</b>	\$ 55

Miscellaneous Options		
<input type="checkbox"/>	1630080 <b>Gel Hip Guides</b> Not compatible with Dual Post Arms. Reduces available Seat Width by 2".	\$ 875
<input type="checkbox"/>	1630110 <b>IV Pole</b> 15" STF and higher.	\$ 295
<input type="checkbox"/>	1630120 <b>O2 Holder</b> 15" STF and higher.	\$ 395
<input type="checkbox"/>	1630160 <b>Clear Lap Tray (Slide On)</b> Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	\$ 490
<input type="checkbox"/>	1630170 <b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10
<input type="checkbox"/>	1630301 <b>Fuze T20 to T50 Conversion Kit</b>	\$ 995

Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/fuze-T20/](http://pdgmobility.com/products/fuze-T20/) to learn more about the **Fuze T20 20° Manual Tilt Wheelchair**



## Cushion

1. Select the type ↓	2. Select the model ↓
-------------------------	--------------------------

<b>General Use</b>	<p><b>Stealth Simplicity</b> 350 lbs weight capacity.</p> <p><input type="checkbox"/> 4009141 16"-20" width, 16"-20" depth \$ 310</p> <p><input type="checkbox"/> 4009142 22"-24" width, 18"-20" depth \$ 370</p>
--------------------	---

<b>Positioning</b>	<p><b>ROHO AirLITE</b> 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009120 14"-20" width, 14"-18" depth \$ 375</p>
--------------------	---

<b>Skin Protection &amp; Positioning</b>	<p><b>Stealth Spectrum Foam</b> 650 lbs weight capacity.</p> <p><input type="checkbox"/> 4009143 16"-20" width, 16"-24" depth \$ 580</p> <p><input type="checkbox"/> 4009144 22"-32" width, 18"-28" depth \$ 630</p>
	<p><b>Stealth Tru-Comfort 2</b> 450 lbs weight capacity.</p> <p><input type="checkbox"/> 4009146 16"-20" width, 16"-24" depth \$ 515</p> <p><input type="checkbox"/> 4009147 22" width, 18"-24" depth \$ 575</p>
	<p><b>nxt BioFit</b> 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009101 14"-20" width, 16"-20" depth \$ 595</p> <p><input type="checkbox"/> 4009102 22"-24" width, 16"-20" depth \$ 625</p>
	<p><b>Matrx VI</b> 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009111 14"-20" width, 15"-20" depth \$ 680</p> <p><input type="checkbox"/> 4009112 22" width, 18"-20" depth \$ 745</p>
	<p><b>Matrx VI Heavy Duty</b> 600 lbs weight capacity.</p> <p><input type="checkbox"/> 4009115 18"-21" width, 18"-26" depth \$ 1,050</p> <p><input type="checkbox"/> 4009116 22"-30" width, 18"-26" depth \$ 1,250</p>

<b>Skin Protection &amp; Adjustable</b>	<p><b>ROHO Single Valve High Profile</b> No weight limit if cushion is properly sized.</p> <p><input type="checkbox"/> 4009125 14"-20" width, 16"-22" depth \$ 840</p> <p><input type="checkbox"/> 4009126 22"-24" width, 16"-22" depth \$ 960</p>
---	--

<b>Skin Protection, Positioning &amp; Adjustable</b>	<p><b>ROHO Quadro Select High Profile</b> No weight limit if cushion is properly sized.</p> <p><input type="checkbox"/> 4009121 14"-20" width, 16"-22" depth \$ 840</p> <p><input type="checkbox"/> 4009122 22"-24" width, 16"-22" depth \$ 960</p>
--	---

<b>Extra-Wide</b>  No weight limit if cushion is properly sized. Custom sizes available.  Custom-made foam cushion with fluid-proof yet breathable fabric cover.	<p><b>PDG Solar Basic</b> High resiliency polyfoam with 2" medium-density base and 1" lower-density top.</p> <p><input type="checkbox"/> 4009200 20"-28" width, 16"-24" depth \$ 395</p> <p><input type="checkbox"/> 4009201 29"-36" width, 16"-24" depth \$ 495</p> <p><input type="checkbox"/> 4009202 37"-48" width, 16"-24" depth \$ 595</p>
	<p><b>PDG Solar Contour</b> Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups.</p> <p><input type="checkbox"/> 4009210 20"-28" width, 16"-24" depth \$ 695</p> <p><input type="checkbox"/> 4009211 29"-36" width, 16"-24" depth \$ 795</p> <p><input type="checkbox"/> 4009212 37"-48" width, 16"-24" depth \$ 895</p>
	<p><b>PDG Solar Pressure Reduction</b> All features of the Contour, plus added pincore latex rubber and containment foam outline.</p> <p><input type="checkbox"/> 4009220 20"-28" width, 16"-24" depth \$ 1,095</p> <p><input type="checkbox"/> 4009221 29"-36" width, 16"-24" depth \$ 1,195</p> <p><input type="checkbox"/> 4009222 37"-48" width, 16"-24" depth \$ 1,295</p>

3. Specify the size →	<p><b>Cushion Size</b> Sizes based on availability from the cushion manufacturer.</p> <p>Specify Size (Width" x Depth") _____</p>
-----------------------	---

## Backrest

1. Select the type ↓	2. Select the contour depth ↓	3. Select the model ↓
-------------------------	----------------------------------	--------------------------

Positioning, posterior-lateral	2" depth	<input type="checkbox"/> 4021141 <b>Stealth ADI AL Active</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 945
	4" depth	<input type="checkbox"/> 4021142 <b>Stealth ADI AL Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,045
	6" depth	<input type="checkbox"/> 4021143 <b>Stealth ADI AL Extra Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,145

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> 4021110 <b>Matrx Elite</b> 300 lbs weight capacity. 12"-20" width, 10"-20" height \$ 955
		<input type="checkbox"/> 4021111 <b>Matrx Elite Heavy Duty</b> 500 lbs weight capacity. 21" width, 16"/20" height \$ 1,195
		<input type="checkbox"/> 4021112 22"-24" width, 16"/20" height \$ 1,295
		<input type="checkbox"/> 4021114 25"-30" width, 16"/20" height \$ 1,995
	3.5" to 4" depth	<input type="checkbox"/> 4021120 <b>ROHO Agility Min</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 900
		<input type="checkbox"/> 4021100 <b>nxt Xtend</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021105 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021122 <b>ROHO Agility Mid</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 930
	6" to 7" depth	<input type="checkbox"/> 4021101 <b>nxt Xtend Deep</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021106 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021116 <b>Matrx Elite Deep</b> 300 lbs weight capacity. 14"-20" width, 12"-20" height \$ 1,145
		<input type="checkbox"/> 4021117 22"-24" width, 12"-20" height \$ 1,295
<input type="checkbox"/> 4021124 <b>ROHO Agility Max</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 1,085		

4. Specify the size →	<p><b>Backrest Size</b> Sizes based on availability from the backrest manufacturer.</p> <p>Specify Size (Width" x Height") _____</p>
-----------------------	--

## Pediatric Seating System

4030400 **Diego Pediatric Seating System**

Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10 (25.40)	12 (30.48)	14 (35.56)	16 (40.64)
Width				
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓

**A Head System** Head Pad & Hardware

**B Contoured Back System** Base, Cushion & Cover

**FlipLock** Configurable as a between canes or full width set up.

**Width Growth** Adapted from our TFB System, the back wings allow for a pediatric growth adjustment.

**Back Hardware Cover** Easy to remove back cover for quick on-demand adjustments and a more aesthetic look.

**Coolcore Covers** Our Diego System uses the Coolcore cover technology which adapt all-year-round for a more comfortable experience.

**Laterals** Modular Laterals either Fixed or Swing Away

**Contoured Seat System** Base, Cushion, Cover & Hardware

**diego** a Stealth Product

**WC-20** Certified



## Headrest

Stealth Comfort Plus Headrest	<b>Stealth Hardware</b>		
	<input type="checkbox"/> 4024110	Multi-Axis HMO475	\$ 225
	<input type="checkbox"/> 4024111	Link TWB480	\$ 350
	<b>Comfort Plus Pad</b>		\$ 275
	<input type="checkbox"/> 4024115	Large (10" x 6")	<input type="checkbox"/> 4024116 X-Large (14" x 6")
nxt MultiFit Head Support	<input type="checkbox"/> 4024100 <b>nxt Hardware</b>		\$ 285
	<b>Multi-Adjustable Pad</b>		\$ 275
	<input type="checkbox"/> 4024105	Contoured (9" x 5")	<input type="checkbox"/> 4024107 X-Large (14" x 5")
	<input type="checkbox"/> 4024106	Large (10" x 5")	

## Positioning Belts & Straps

Bodypoint	<b>Hip Belts</b>		
	<input type="checkbox"/> 4029010	Center-Pull Auto Style	\$ 115
	<input type="checkbox"/> 4029015	Center-Pull Auto Style, Padded	\$ 215
	<input type="checkbox"/> 4029020	4-Point, Padded	\$ 235
	<input type="checkbox"/> 4029030	<b>Chest Strap</b>	\$ 165

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.

PDG proudly supplies on-chair seating and positioning products from these manufacturers:

