

Customer and Order Information		
Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type	
<input checked="" type="checkbox"/> 1601000 Fuze T20 20° Manual Tilt Wheelchair	\$ 3,995

Frame Options	
<input type="checkbox"/> 1602010 Fixed Base and Seat Frame 250 lbs Weight Capacity. Select Fixed Seat Width below.	
<input type="checkbox"/> 1602025 Adjustable Base and Seat Frame 250 lbs Weight Capacity. Select Adjustable Seat Width below.	\$ 100
<input type="checkbox"/> 1602300 Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.	
User Weight Range Required for selection of Gas Struts (Force measured in Newton - N)	
<input type="checkbox"/> 1602406 Up to 150 lbs (600 N)	<input type="checkbox"/> 1602410 150 lbs to 200 lbs (1000 N)
<input type="checkbox"/> 1602412 200 lbs to 250 lbs (1200 N)	

Seat Width		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW
Fixed Seat Widths Order with Fixed Base and Seat Frame.		
<input type="checkbox"/> 1603015 15" Seat Width	<input type="checkbox"/> 1603016 16" Seat Width	<input type="checkbox"/> 1603017 17" Seat Width
<input type="checkbox"/> 1603018 18" Seat Width	<input type="checkbox"/> 1603019 19" Seat Width	<input type="checkbox"/> 1603020 20" Seat Width
<input type="checkbox"/> 1603021 21" Seat Width	<input type="checkbox"/> 1603022 22" Seat Width	\$ 215
Adjustable Seat Widths Order with Adjustable Base and Seat Frame.		
<input type="checkbox"/> 1603112 12"-15" Range - Set at 12"	<input type="checkbox"/> 1603113 12"-15" Range - Set at 13"	<input type="checkbox"/> 1603114 12"-15" Range - Set at 14"
<input type="checkbox"/> 1603115 12"-15" Range - Set at 15"		
<input type="checkbox"/> 1603215 15"-20" Range - Set at 15"	<input type="checkbox"/> 1603216 15"-20" Range - Set at 16"	<input type="checkbox"/> 1603217 15"-20" Range - Set at 17"
<input type="checkbox"/> 1603218 15"-20" Range - Set at 18"	<input type="checkbox"/> 1603219 15"-20" Range - Set at 19"	<input type="checkbox"/> 1603220 15"-20" Range - Set at 20"
<input type="checkbox"/> 1603321 21"-22" Range - Set at 21"	<input type="checkbox"/> 1603322 21"-22" Range - Set at 22"	\$ 215


Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 1604015 15" Seat Depth	<input type="checkbox"/> 1604016 16" Seat Depth
<input type="checkbox"/> 1604017 17" Seat Depth	<input type="checkbox"/> 1604018 18" Seat Depth
<input type="checkbox"/> 1604019 19" Seat Depth	<input type="checkbox"/> 1604020 20" Seat Depth

Frame Colour	
PDG Premium Colours	
<input type="checkbox"/> 1605000 Tracer Red	<input type="checkbox"/> 1605001 Gloss Black
<input type="checkbox"/> 1605002 Silver	<input type="checkbox"/> 1605003 Black Shimmer
<input type="checkbox"/> 1605004 Orange Fusion	<input type="checkbox"/> 1605005 Karl's Imagination
<input type="checkbox"/> 1605006 Blue Shimmer	<input type="checkbox"/> 1605007 Gloss White
<input type="checkbox"/> 1605008 Dark Titanium	<input type="checkbox"/> 1605009 Geyser Green
<input type="checkbox"/> 1605010 Lollypop Blue	<input type="checkbox"/> 1605011 Lollypop Red
<input type="checkbox"/> 1605012 Bubble Gum Pink	<input type="checkbox"/> 1605013 Berry Purple
<input type="checkbox"/> 1605014 Safety Yellow	
<input type="checkbox"/> 1605200 Custom Colour Contact PDG for specifications.	\$ 300

Arm Type		
<input type="checkbox"/>	1606000	No Armrest
<input type="checkbox"/>	1606100	Single Post Fixed Height (12")
	Single Post Adjustable Height	8"-11" not compatible with 22" and 24" Rear Wheels. \$ 255
<input type="checkbox"/>	1606111	8"-11" <input type="checkbox"/>
<input type="checkbox"/>	1606114	10"-14"
<input type="checkbox"/>	1606200	Dual Post Fixed Height (12") Full Length Armpad only.
	Dual Post Adjustable Height	8"-11" not compatible with 22" and 24" Rear Wheels. \$ 255
<input type="checkbox"/>	1606211	8"-11" <input type="checkbox"/>
<input type="checkbox"/>	1606214	10"-14" <input type="checkbox"/>
<input type="checkbox"/>	1606216	12"-16"
	Locking Flip-up	Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery. \$ 255
<input type="checkbox"/>	1606310	Set at 10" <input type="checkbox"/>
<input type="checkbox"/>	1606311	Set at 11" <input type="checkbox"/>
<input type="checkbox"/>	1606312	Set at 12" <input type="checkbox"/>
<input type="checkbox"/>	1606313	Set at 13" <input type="checkbox"/>
<input type="checkbox"/>	1606314	Set at 14" <input type="checkbox"/>

Armrests		
<input type="checkbox"/>	1607095	Curved Armpads Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests, 22" and 24" Rear Wheels.
	Regular Armpads	
<input type="checkbox"/>	1607100	Desk Length - 2" x 10" <input type="checkbox"/>
<input type="checkbox"/>	1607101	Full Length - 2" x 14"
<input type="checkbox"/>	1607102	Full Length Premium Grade - 2" x 14" \$ 35
	Wide Armpads	May interfere with wheels in some tilt configurations. \$ 115
<input type="checkbox"/>	1607110	Desk Length - 3" x 10" <input type="checkbox"/>
<input type="checkbox"/>	1607111	Full Length - 3" x 12"
	Gel Armpads	May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1607120	Full Length - 2" x 12" <input type="checkbox"/>
<input type="checkbox"/>	1607121	Full Length - 3.5" x 12"
	Contour Arm Trough with Hand Pad	May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1607130	Right (4.5" x 18") <input type="checkbox"/>
<input type="checkbox"/>	1607131	Left (4.5" x 18")

Seat Pan		
<input type="checkbox"/>	1608000	No Seat Pan
<input type="checkbox"/>	1608001	Solid Aluminum Seat Pan \$ 120

Cushion		
<input type="checkbox"/>	1609000	No Cushion
	nxt BioFit	
<input type="checkbox"/>	1609101	14"-20" width, 16"-20" depth  \$ 595
<input type="checkbox"/>	1609102	22"-24" width, 16"-20" depth \$ 625
	Cushion Size	Available in even number sizes only.
1609400	Specify Size (Width" x Depth")	_____

Seat Height

Select desired Seat-to-Floor Height

Standard T20 Options

RW	5" Casters					
12"	13"					
16"	13"	14"	15"			
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	

RW	6" Casters & 6" Forks					
16"		14"	15"			
20"		14"	15"	16"	17"	
22"		14"	15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	6" Casters & 8" Forks					
16"			15"			
20"			15"	16"	17"	
22"			15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	8" Casters & 8" Forks					
20"				16"	17"	
22"				16"	17"	18"
24"				16"	17"	18"

Narrow T20 Options

Reduces the overall width of the Chair by 2".
Recommended only if the width of the chair is a concern. Some STF combinations offer less than 20° of tilt, but all offer a minimum of 12°.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	8" Casters & 8" Forks					
16"				16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

- Standard
 Narrow

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

- 1611051 5" x 1.25" Caster size
- 1611061 6" x 1.25" Caster size
- 1611063 6" x 1.50" Caster size \$ 85
- 1611081 8" x 1.25" Caster size
- 1611084 8" x 2.00" Caster size \$ 135

Pneumatic

- 1611161 6" x 1.25" Caster size \$ 110
- 1611181 8" x 1.25" Caster size \$ 120
- 1611184 8" x 2.00" Caster size \$ 150
- 1611260 **6" Soft Roll** \$ 95

Rear Wheels

Composite

- 1612012 12" Wheel size 1612016 16" Wheel size 1612020 20" Wheel size
- 1612022 22" Wheel size 1612024 24" Wheel size

Spoke

- 1612120 20" Wheel size 1612122 22" Wheel size 1612124 24" Wheel size \$ 75
- 1612224 **24" Heavy Duty Spoke** \$ 350

Axle Type

- 1613000 **Quick Release Rear Axles**
- 1613010 **Fixed Rear Axles**

Rear Tires			
<input type="checkbox"/>	1614000	Urethane	
<input type="checkbox"/>	1614010	Pneumatic Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/>	1614020	Pneumatic with Flat-Free Inserts Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.	\$ 95

Hand Rims			
<input type="checkbox"/>	1615000	No Hand Rims	
<input type="checkbox"/>	1615010	Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1615020	Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 135

Wheel Locks			
<input type="checkbox"/>	1616010	Push to Lock	
<input type="checkbox"/>	1616020	Pull to Lock	
<input type="checkbox"/>	1616030	Attendant Wheel Lock (Push to Lock) Not compatible with Wheel Lock Extensions.	\$ 85
<input type="checkbox"/>	1616040	User & Attendant Wheel Lock	\$ 125
<input type="checkbox"/>	1616100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	\$ 60

Anti-Tippers			
<input checked="" type="checkbox"/>	1617010	Rear Anti-Tippers Required Safety Option.	\$ 60

Back Post Style			
Push Bar Back			
<input type="checkbox"/>	1618010	Straight	
<input type="checkbox"/>	1618030	8° Bend	\$ 75
Cane Handle Back Not compatible with Headrest Mounting Plate.			
<input type="checkbox"/>	1618020	Straight	
<input type="checkbox"/>	1618040	8° Bend	\$ 75
Back Post Height Seat Pan to Handle.			
<input type="checkbox"/>	1618420	20" Height	
<input type="checkbox"/>	1618425	25" Height	

Back Options			
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.			\$ 150
<input type="checkbox"/>	1619010	Button on Outside	
<input type="checkbox"/>	1619011	Button on Inside	
<input type="checkbox"/>	1619050	Removable Cane Handle Push Bar Not compatible with Push Bar Back.	\$ 140
Seat to Back Angle			
<input type="checkbox"/>	1619090	Set at 90°	
<input type="checkbox"/>	1619096	Set at 96°	
<input type="checkbox"/>	1619102	Set at 102°	
<input type="checkbox"/>	1619108	Set at 108°	
<input type="checkbox"/>	1619114	Set at 114°	
<input type="checkbox"/>	1619120	Set at 120°	
<input type="checkbox"/>	1619220	Fold Down Back	
<input type="checkbox"/>	1619221	Non Fold Down Back	\$ 75

Back Upholstery			
<input type="checkbox"/>	1620000	No Back Upholstery	
Standard Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.			\$ 155
<input type="checkbox"/>	1620018	Set at 18"	
<input type="checkbox"/>	1620019	Set at 19"	
<input type="checkbox"/>	1620020	Set at 20"	
<input type="checkbox"/>	1620021	Set at 21"	
<input type="checkbox"/>	1620022	Set at 22"	
<input type="checkbox"/>	1620023	Set at 23"	
Adjustable Tension Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.			\$ 295
<input type="checkbox"/>	1620321	Set at 21"	
<input type="checkbox"/>	1620322	Set at 22"	
<input type="checkbox"/>	1620323	Set at 23"	

Backrest			
<input type="checkbox"/>	1621000	No Backrest	
nxt Optima HA Active Contour			
<input type="checkbox"/>	1621100	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1621105	22"-24" width, 16"-20" height	\$ 720
nxt Optima HA Plus Contour			
<input type="checkbox"/>	1621101	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1621106	22"-24" width, 16"-20" height	\$ 720
Backrest Size Available in even number sizes only.			
	1621400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1623011	Trigger on Right Back Post	
<input type="checkbox"/>	1623012	Trigger on Left Back Post	

Headrest			
<input type="checkbox"/>	1624000	No Headrest	
<input type="checkbox"/>	1624010	Headrest Mounting Plate Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 75
<input type="checkbox"/>	1624100	nxt MultiFit Head Support Hardware	\$ 235
nxt MultiFit Head Support Pad			
<input type="checkbox"/>	1624105	Regular (9" x 5")	
<input type="checkbox"/>	1624106	Large (10" x 5")	
<input type="checkbox"/>	1624107	X-Large (14" x 5")	\$ 245

Front Rigging			
<input type="checkbox"/>	1625000	No Front Rigging	
PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.			
<input type="checkbox"/>	1625410	Fixed 70°	Available in March \$ 95
<input type="checkbox"/>	1625420	Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____	Available in April \$ 195
<input type="checkbox"/>	1625430	Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.	Available in April \$ 495
Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.			
<input type="checkbox"/>	1625140	Swing Away 70°	
<input type="checkbox"/>	1625160	Swing Away 90° Not compatible with 8" Casters.	\$ 55
<input type="checkbox"/>	1625200	Swing Away Elevating Legrest	\$ 375
<input type="checkbox"/>	1625300	1" Offset Hangers for Swing Away	\$ 75

Foot Plates			
<input type="checkbox"/>	1626010	Composite	
<input type="checkbox"/>	1626020	Aluminum Angle Adjustable	\$ 155
<input type="checkbox"/>	1626050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1627010	Swing Away 70° Calf Pad Not compatible with 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1627020	Swing Away 70° Gel Padding Not compatible with 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1627030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1627040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	1627050	Bodypoint - Padded Calf Strap	\$ 105

Extension Tubes			
<input type="checkbox"/>	1528010	Short 13"-17" for Easy-Out SA Angle Adjustable 10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1528020	Long 16"-22.5" for Easy-Out SA Angle Adjustable 12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
<input type="checkbox"/>	1628200	For 90° FR (11"-16") Not compatible with Easy-Out, Standard 70° or ELR	\$ 75

Positioning Straps			
<input type="checkbox"/>	1629005	Lap Belt Mount	\$ 55
Bodypoint - Hip Belts			
<input type="checkbox"/>	1629010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1629015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1629020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1629030	Bodypoint - Chest Strap	\$ 85



Miscellaneous Options			
<input type="checkbox"/>	1630080	Gel Hip Guides	\$ 325
<input type="checkbox"/>	1630110	IV Pole 15" STF and higher.	\$ 295
<input type="checkbox"/>	1630120	O2 Holder 15" STF and higher.	\$ 395
<input type="checkbox"/>	1630160	Clear Lap Tray (Slide On) Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	\$ 490
<input type="checkbox"/>	1630170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10
<input type="checkbox"/>	1630301	Fuze T20 to T50 Conversion Kit	\$ 995

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/фуze-T20/ to learn more about the **Fuze T20 20° Manual Tilt Wheelchair**

