

Customer and Order Information

Purchase Order *	Tag	v.040518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 1501000	Fuze T50 50° Manual Tilt Wheelchair	\$ 3,995
---	--	----------

Frame Options

Fixed Base and Seat Frame Select Fixed Seat Width below.

<input type="checkbox"/> 1502010	300 lbs Weight Capacity	
<input type="checkbox"/> 1502011	350 lbs Weight Capacity	Not compatible with Transit Tie Down System. Requires Flat-free Rear Tires or Heavy Duty Spoke Rear Wheels.
<input type="checkbox"/> 1502025	Adjustable Base and Seat Frame	250 lbs Weight Capacity. Select Adjustable Seat Width below. \$ 100
<input type="checkbox"/> 1502300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.

Seat Width Overall chair width is +7.5" for 12/16" RW and +8.5" for 20/22/24" RW

Fixed Seat Widths Order with Fixed Base and Seat Frame.

<input type="checkbox"/> 1503015	15" Seat Width	<input type="checkbox"/> 1503016	16" Seat Width	<input type="checkbox"/> 1503017	17" Seat Width	
<input type="checkbox"/> 1503018	18" Seat Width	<input type="checkbox"/> 1503019	19" Seat Width	<input type="checkbox"/> 1503020	20" Seat Width	
<input type="checkbox"/> 1503021	21" Seat Width	<input type="checkbox"/> 1503022	22" Seat Width			\$ 215
<input type="checkbox"/> 1503023	23" Seat Width	<input type="checkbox"/> 1503024	24" Seat Width			\$ 325

Adjustable Seat Widths Order with Adjustable Base and Seat Frame.

<input type="checkbox"/> 1503112	12"-15" Range - Set at 12"	<input type="checkbox"/> 1503113	12"-15" Range - Set at 13"	<input type="checkbox"/> 1503114	12"-15" Range - Set at 14"	
<input type="checkbox"/> 1503115	12"-15" Range - Set at 15"					
<input type="checkbox"/> 1503215	15"-20" Range - Set at 15"	<input type="checkbox"/> 1503216	15"-20" Range - Set at 16"	<input type="checkbox"/> 1503217	15"-20" Range - Set at 17"	
<input type="checkbox"/> 1503218	15"-20" Range - Set at 18"	<input type="checkbox"/> 1503219	15"-20" Range - Set at 19"	<input type="checkbox"/> 1503220	15"-20" Range - Set at 20"	
<input type="checkbox"/> 1503321	21"-22" Range - Set at 21"	<input type="checkbox"/> 1503322	21"-22" Range - Set at 22"			\$ 215

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 1504015	15" Seat Depth	<input type="checkbox"/> 1504016	16" Seat Depth	<input type="checkbox"/> 1504017	17" Seat Depth	
<input type="checkbox"/> 1504018	18" Seat Depth	<input type="checkbox"/> 1504019	19" Seat Depth	<input type="checkbox"/> 1504020	20" Seat Depth	

Custom Seat Depths Requires 25" Push Bar Back Straight, Angle/Height Adjustable Push Bar. Not compatible with Reclining or Fold Down Back. \$ 325

<input type="checkbox"/> 1504021	21" Seat Depth	<input type="checkbox"/> 1504022	22" Seat Depth	
----------------------------------	----------------	----------------------------------	----------------	--

Frame Colour

PDG Premium Colours

<input type="checkbox"/> 1505000	Tracer Red	<input type="checkbox"/> 1505001	Gloss Black	<input type="checkbox"/> 1505002	Silver
<input type="checkbox"/> 1505003	Black Shimmer	<input type="checkbox"/> 1505004	Orange Fusion	<input type="checkbox"/> 1505005	Karl's Imagination
<input type="checkbox"/> 1505006	Blue Shimmer	<input type="checkbox"/> 1505007	Gloss White	<input type="checkbox"/> 1505008	Dark Titanium
<input type="checkbox"/> 1505009	Geyser Green	<input type="checkbox"/> 1505010	Lollypop Blue	<input type="checkbox"/> 1505011	Lollypop Red
<input type="checkbox"/> 1505012	Bubble Gum Pink	<input type="checkbox"/> 1505013	Berry Purple	<input type="checkbox"/> 1505014	Safety Yellow
<input type="checkbox"/> 1505200	Custom Colour	Contact PDG for specifications.			\$ 300

Arm Type		
<input type="checkbox"/>	1506000	No Armrest
<input type="checkbox"/>	1506100	Single Post Fixed Height (12")
<input type="checkbox"/>	1506114	Single Post Adjustable Height (10"-14") \$ 255
<input type="checkbox"/>	1506200	Dual Post Fixed Height (12") Full Length Armpad only.
		Dual Post Adjustable Height \$ 255
<input type="checkbox"/>	1506214	10"-14"
<input type="checkbox"/>	1506216	12"-16"
		Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery. \$ 255
<input type="checkbox"/>	1506310	Set at 10"
<input type="checkbox"/>	1506311	Set at 11"
<input type="checkbox"/>	1506312	Set at 12"
<input type="checkbox"/>	1506313	Set at 13"
<input type="checkbox"/>	1506314	Set at 14"

Armrests		
<input type="checkbox"/>	1507095	Curved Armpads Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.
		Regular Armpads
<input type="checkbox"/>	1507100	Desk Length - 2" x 10"
<input type="checkbox"/>	1507101	Full Length - 2" x 14"
<input type="checkbox"/>	1507102	Full Length Premium Grade - 2" x 14" \$ 35
		Wide Armpads May interfere with wheels in some tilt configurations. \$ 115
<input type="checkbox"/>	1507110	Desk Length - 3" x 10"
<input type="checkbox"/>	1507111	Full Length - 3" x 12"
		Gel Armpads May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1507120	Full Length - 2" x 12"
<input type="checkbox"/>	1507121	Full Length - 3.5" x 12"
		Contour Arm Trough with Hand Pad \$ 125
		Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.
<input type="checkbox"/>	1507130	Right (4.5" x 18")
<input type="checkbox"/>	1507131	Left (4.5" x 18")

Seat Pan		
<input type="checkbox"/>	1508000	No Seat Pan
<input type="checkbox"/>	1508001	Solid Aluminum Seat Pan \$ 120

Cushion		
<input type="checkbox"/>	1509000	No Cushion
		nxt BioFit
<input type="checkbox"/>	1509101	12"-20" width, 14"-20" depth \$ 595
<input type="checkbox"/>	1509102	22"-24" width, 14"-20" depth \$ 625
		Cushion Size Available in even number sizes only.
<input type="checkbox"/>	1509400	Specify Size (Width" x Depth") _____



Seat Height

Select desired Seat-to-Floor Height

Standard T50 Options

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
12"					18"	19"
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	8" Casters & 8" Forks					
12"					18"	19"
16"				16"		20"
20"				16"	17"	18"
22"					17"	18"
24"					18"	19"

RW	6" Soft Roll & Ultra Sport Forks					
16"				16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

Wide T50 Options

Increases the overall width of the Chair by 2".

RW	5" Casters					
12"	13"					
16"	13"	14"	15"			
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	

RW	6" Casters & 6" Forks					
16"		14"	15"			
20"		14"	15"	16"	17"	
22"		14"	15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	6" Casters & 8" Forks					
16"			15"			
20"			15"	16"	17"	
22"			15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	8" Casters & 8" Forks					
20"				16"	17"	
22"				16"	17"	18"
24"				16"	17"	18"

- Standard
 Wide

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

- 1511051 5" x 1.25" Caster size
- 1511061 6" x 1.25" Caster size
- 1511063 6" x 1.50" Caster size \$ 85
- 1511081 8" x 1.25" Caster size
- 1511084 8" x 2.00" Caster size \$ 135

Pneumatic

- 1511161 6" x 1.25" Caster size \$ 110
- 1511181 8" x 1.25" Caster size \$ 120
- 1511184 8" x 2.00" Caster size \$ 150

- 1511260 **6" Soft Roll** \$ 95

- 1511400 **Ultra Sport Fork** Not compatible with 12" Rear Wheels. Requires 6" Soft Roll Caster. \$ 395

Rear Wheels

Composite

- | | | |
|---|---|---|
| <input type="checkbox"/> 1512012 12" Wheel size | <input type="checkbox"/> 1512016 16" Wheel size | <input type="checkbox"/> 1512020 20" Wheel size |
| <input type="checkbox"/> 1512022 22" Wheel size | <input type="checkbox"/> 1512024 24" Wheel size | |

Spoke

- | | | | |
|--|---|---|--------|
| <input type="checkbox"/> 1512120 20" Wheel size | <input type="checkbox"/> 1512122 22" Wheel size | <input type="checkbox"/> 1512124 24" Wheel size | \$ 75 |
| <input type="checkbox"/> 1512224 24" Heavy Duty Spoke | | | \$ 350 |

Axle Type

- | |
|--|
| <input type="checkbox"/> 1513000 Quick Release Rear Axles |
| <input type="checkbox"/> 1513010 Fixed Rear Axles |

Rear Tires

- | | | |
|--|--|-------|
| <input type="checkbox"/> 1514000 Urethane | | |
| <input type="checkbox"/> 1514010 Pneumatic | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 1514020 Pneumatic with Flat-Free Inserts | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. Required for 350 lbs Weight Capacity. | \$ 95 |

Hand Rims

- | | | |
|--|-------------------------------------|--------|
| <input type="checkbox"/> 1515000 No Hand Rims | | |
| <input type="checkbox"/> 1515010 Standard | 20", 22", and 24" Rear Wheels only. | |
| <input type="checkbox"/> 1515020 Plastic Coated | 20", 22", and 24" Rear Wheels only. | \$ 135 |

Wheel Locks

- | | | |
|---|--|--------|
| <input type="checkbox"/> 1516005 Foot Operated Wheel Lock | Not compatible with One Arm Drive. | \$ 200 |
| <input type="checkbox"/> 1516010 Push to Lock | | |
| <input type="checkbox"/> 1516020 Pull to Lock | | |
| <input type="checkbox"/> 1516030 Attendant Wheel Lock (Push to Lock) | Not compatible with Wheel Lock Extensions. | \$ 85 |
| <input type="checkbox"/> 1516040 User & Attendant Wheel Lock | | \$ 125 |
| <input type="checkbox"/> 1516100 Wheel Lock Extensions | Not compatible with Attendant Wheel Lock. | \$ 60 |

Anti-Tippers

- | | | |
|--|-------------------------|-------|
| <input checked="" type="checkbox"/> 1517010 Rear Anti-Tippers | Required Safety Option. | \$ 60 |
|--|-------------------------|-------|

Back Post Style

Push Bar Back

- | | | |
|---|--|-------|
| <input type="checkbox"/> 1518010 Straight | | |
| <input type="checkbox"/> 1518030 8° Bend | | \$ 75 |

Cane Handle Back

Not compatible with Headrest Mounting Plate.

- | | | |
|---|--|-------|
| <input type="checkbox"/> 1518020 Straight | | |
| <input type="checkbox"/> 1518040 8° Bend | | \$ 75 |

Back Post Height

Seat Pan to Handle.

- | | |
|---|---|
| <input type="checkbox"/> 1518420 20" Height | <input type="checkbox"/> 1518425 25" Height |
|---|---|

Back Options

<input type="checkbox"/> 1519010 Angle/Height Adjustable Push Bar Handle Button on Outside		<input type="checkbox"/> 1519011 Angle/Height Adjustable Push Bar Handle Button on Inside		Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 1519050 Removable Cane Handle Push Bar		Not compatible with Push Bar Back.				\$ 140
Seat to Back Angle						
<input type="checkbox"/> 1519090	Set at 90°	<input type="checkbox"/> 1519096	Set at 96°	<input type="checkbox"/> 1519102	Set at 102°	
<input type="checkbox"/> 1519108	Set at 108°	<input type="checkbox"/> 1519114	Set at 114°	<input type="checkbox"/> 1519120	Set at 120°	
<input type="checkbox"/> 1519220	Fold Down Back					
<input type="checkbox"/> 1519221	Non Fold Down Back					\$ 75

Back Upholstery

<input type="checkbox"/> 1520000	No Back Upholstery					
Standard		Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.				\$ 155
<input type="checkbox"/> 1520018	Set at 18"	<input type="checkbox"/> 1520019	Set at 19"	<input type="checkbox"/> 1520020	Set at 20"	
<input type="checkbox"/> 1520021	Set at 21"	<input type="checkbox"/> 1520022	Set at 22"	<input type="checkbox"/> 1520023	Set at 23"	
Adjustable Tension		Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.				\$ 295
<input type="checkbox"/> 1520321	Set at 21"	<input type="checkbox"/> 1520322	Set at 22"	<input type="checkbox"/> 1520323	Set at 23"	

Backrest

<input type="checkbox"/> 1521000	No Backrest					
nxt Optima HA Active Contour						
<input type="checkbox"/> 1521100	16"-20" width, 16"-20" height					\$ 690
<input type="checkbox"/> 1521105	22"-24" width, 16"-20" height					\$ 720
nxt Optima HA Plus Contour						
<input type="checkbox"/> 1521101	16"-20" width, 16"-20" height					\$ 690
<input type="checkbox"/> 1521106	22"-24" width, 16"-20" height					\$ 720
Backrest Size Available in even number sizes only.						
1521400	Specify Size (Width" x Height") _____					



Tilt/Recline Options

Attendant Tilt only						
<input type="checkbox"/> 1523011	Trigger on Right Back Post	<input type="checkbox"/> 1523012	Trigger on Left Back Post			
Tilt Limit						
<input type="checkbox"/> 1523220	Set at 20°	<input type="checkbox"/> 1523225	Set at 25°	<input type="checkbox"/> 1523230	Set at 30°	
<input type="checkbox"/> 1523235	Set at 35°	<input type="checkbox"/> 1523240	Set at 40°	<input type="checkbox"/> 1523245	Set at 45°	
<input type="checkbox"/> 1523250	Set at 50°					
<input type="checkbox"/> 1523410	Reclining Backrest (Manual)					\$ 875
30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Tension Adjustable Upholstery. Requires 25" Angle/Height Adjustable Push Bar with 8" Bend. Not available with 15", 16", 21", or 22" Seat Depth.						
<input type="checkbox"/> 1523400	Power Tilt 250 lbs Weight Capacity. 4° to 45° tilt range. Includes 2 batteries for off-chair charging. Minimum Seat Width is 15".					\$ 1,695
	Toggle Switch	<input type="checkbox"/> 1523421	Installed on Right side	<input type="checkbox"/> 1523422	Installed on Left side	
	Button Switch	<input type="checkbox"/> 1523431	Installed on Right side	<input type="checkbox"/> 1523432	Installed on Left side	
	Egg Switch	<input type="checkbox"/> 1523440	Not installed			
<input type="checkbox"/> 1523445	Egg Switch Swing Away Mounting Hardware					\$ 625
Not compatible with all Headrests. Specify make and model to be used: _____						

Headrest			
<input type="checkbox"/>	1524000	No Headrest	
<input type="checkbox"/>	1524010	Headrest Mounting Plate Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 75
<input type="checkbox"/>	1524100	nxt MultiFit Head Support Hardware	\$ 235
<input type="checkbox"/>		nxt MultiFit Head Support Pad	\$ 245
<input type="checkbox"/>	1524105	Regular (9" x 5")	
<input type="checkbox"/>	1524106	Large (10" x 5")	
<input type="checkbox"/>	1524107	X-Large (14" x 5")	

Front Rigging			
<input type="checkbox"/>	1525000	No Front Rigging	
<input type="checkbox"/>		PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.	
<input type="checkbox"/>	1525410	Fixed 70°	Available in March \$ 95
<input type="checkbox"/>	1525420	Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____	Available in April \$ 195
<input type="checkbox"/>	1525430	Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.	Available in April \$ 495
<input type="checkbox"/>		Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.	
<input type="checkbox"/>	1525140	Swing Away 70°	
<input type="checkbox"/>	1525160	Swing Away 90° Not compatible with 8" Casters.	\$ 55
<input type="checkbox"/>	1525200	Swing Away Elevating Legrest	\$ 375
<input type="checkbox"/>	1525300	1" Offset Hangers for Swing Away	\$ 75

Foot Plates			
<input type="checkbox"/>	1526010	Composite	
<input type="checkbox"/>	1526020	Aluminum Angle Adjustable	\$ 155
<input type="checkbox"/>	1526040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1526050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1527015	Swing Away 70° Calf Pad Not compatible with 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1527020	Swing Away 70° Gel Padding Not compatible with 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1527030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1527040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	1527050	Bodypoint - Padded Calf Strap	\$ 105
<input type="checkbox"/>	1527060	Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.	\$ 190

Extension Tubes			
<input type="checkbox"/>	1528010	Short 13"-17" for Easy-Out SA Angle Adjustable 10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1528020	Long 16"-22.5" for Easy-Out SA Angle Adjustable 12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
<input type="checkbox"/>	1628200	For 90° FR (11"-16") Not compatible with Easy-Out, Standard 70° or ELR	\$ 75

Positioning Straps			
<input type="checkbox"/>	1529005	Lap Belt Mount	\$ 55
<input type="checkbox"/>		Bodypoint - Hip Belts	
<input type="checkbox"/>	1529010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1529015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1529020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1529030	Bodypoint - Chest Strap	\$ 85

Miscellaneous Options		
One-Arm Drive Requires 24" Composite Rear Wheels and forward wheel position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks.		
<input type="checkbox"/> 1530010	Installed on Right side	<input type="checkbox"/> 1530011 Installed on Left side
<input type="checkbox"/> 1530060	Stealth Flip Down Abductor	\$ 495
<input type="checkbox"/> 1530070	QR Bracket - Accessory Mount	\$ 200
<input type="checkbox"/> 1530080	Gel Hip Guides Not compatible with Dual Post Arms.	\$ 325
<input type="checkbox"/> 1530110	IV Pole	\$ 295
<input type="checkbox"/> 1530120	O2 Holder May add to overall width of the Chair.	\$ 395
<input type="checkbox"/> 1530150	Ventilator Tray Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 16". Not compatible with O2 Holder, Power Tilt.	\$ 295
<input type="checkbox"/> 1530152	Ventilator Bracket Compatible with LTV1 (laptop) ventilator.	\$ 295
<input type="checkbox"/> 1530154	Ventilator Battery Tray 16" minimum Seat Width. Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with STF restrictions.	\$ 95
<input type="checkbox"/> 1530160	Clear Lap Tray (Slide On) Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	\$ 490
<input type="checkbox"/> 1530170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10
<input type="checkbox"/> 1530300	Fuze T50 to T20 Conversion Kit	\$ 995

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/fuze-T50/ to learn more about the **Fuze T50 50° Manual Tilt Wheelchair**

