

Customer and Order Information	
Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address * *required	

Frame Type	
<input checked="" type="checkbox"/> 1501000 Fuze T50 50° Manual Tilt Wheelchair	\$ 3,995

Frame Options	
Fixed Base and Seat Frame Select Fixed Seat Width below.	
<input type="checkbox"/> 1502010 300 lbs Weight Capacity	
<input type="checkbox"/> 1502011 350 lbs Weight Capacity Not compatible with Transit Tie Down System. Requires Flat-free Rear Tires or Heavy Duty Spoke Rear Wheels.	
<input type="checkbox"/> 1502025 Adjustable Base and Seat Frame 250 lbs Weight Capacity. Select Adjustable Seat Width below.	\$ 100
<input type="checkbox"/> 1502300 Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.	

Seat Width		Overall chair width is +7.5" for 12/16" RW and +8.5" for 20/22/24" RW
Fixed Seat Widths Order with Fixed Base and Seat Frame.		
<input type="checkbox"/> 1503013 13" Seat Width	<input type="checkbox"/> 1503014 14" Seat Width	<input type="checkbox"/> 1503015 15" Seat Width
<input type="checkbox"/> 1503016 16" Seat Width	<input type="checkbox"/> 1503017 17" Seat Width	<input type="checkbox"/> 1503018 18" Seat Width
<input type="checkbox"/> 1503019 19" Seat Width	<input type="checkbox"/> 1503020 20" Seat Width	
<input type="checkbox"/> 1503021 21" Seat Width	<input type="checkbox"/> 1503022 22" Seat Width	\$ 215
<input type="checkbox"/> 1503023 23" Seat Width	<input type="checkbox"/> 1503024 24" Seat Width	\$ 325
Adjustable Seat Widths Order with Adjustable Base and Seat Frame.		
<input type="checkbox"/> 1503113 13"-15" Range - Set at 13"	<input type="checkbox"/> 1503114 13"-15" Range - Set at 14"	<input type="checkbox"/> 1503115 13"-15" Range - Set at 15"
<input type="checkbox"/> 1503215 15"-20" Range - Set at 15"	<input type="checkbox"/> 1503216 15"-20" Range - Set at 16"	<input type="checkbox"/> 1503217 15"-20" Range - Set at 17"
<input type="checkbox"/> 1503218 15"-20" Range - Set at 18"	<input type="checkbox"/> 1503219 15"-20" Range - Set at 19"	<input type="checkbox"/> 1503220 15"-20" Range - Set at 20"
<input type="checkbox"/> 1503321 21"-22" Range - Set at 21"	<input type="checkbox"/> 1503322 21"-22" Range - Set at 22"	\$ 215


Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 1504015 15" Seat Depth	<input type="checkbox"/> 1504016 16" Seat Depth
<input type="checkbox"/> 1504018 18" Seat Depth	<input type="checkbox"/> 1504019 19" Seat Depth
<input type="checkbox"/> 1504017 17" Seat Depth	<input type="checkbox"/> 1504020 20" Seat Depth
Custom Seat Depths Requires 25" Push Bar Back Straight, Angle/Height Adjustable Push Bar. Not compatible with Reclining or Fold Down Back.	
<input type="checkbox"/> 1504021 21" Seat Depth	<input type="checkbox"/> 1504022 22" Seat Depth
	\$ 325

Frame Colour	
PDG Premium Colours	
<input type="checkbox"/> 1505000 Tracer Red	<input type="checkbox"/> 1505001 Gloss Black
<input type="checkbox"/> 1505003 Black Shimmer	<input type="checkbox"/> 1505004 Orange Fusion
<input type="checkbox"/> 1505006 Blue Shimmer	<input type="checkbox"/> 1505007 Gloss White
<input type="checkbox"/> 1505009 Geyser Green	<input type="checkbox"/> 1505010 Lollypop Blue
<input type="checkbox"/> 1505012 Bubble Gum Pink	<input type="checkbox"/> 1505013 Berry Purple
<input type="checkbox"/> 1505002 Silver	<input type="checkbox"/> 1505005 Karl's Imagination
<input type="checkbox"/> 1505008 Dark Titanium	<input type="checkbox"/> 1505011 Lollypop Red
<input type="checkbox"/> 1505014 Safety Yellow	
<input type="checkbox"/> 1505200 Custom Colour Contact PDG for specifications.	\$ 300

Arm Type		
<input type="checkbox"/>	1506000	No Armrest
<input type="checkbox"/>	1506100	Single Post Fixed Height (12")
<input type="checkbox"/>	1506114	Single Post Adjustable Height (10"-14") \$ 255
<input type="checkbox"/>	1506200	Dual Post Fixed Height (12") Full Length Armpad only.
		Dual Post Adjustable Height \$ 255
<input type="checkbox"/>	1506214	10"-14"
<input type="checkbox"/>	1506216	12"-16"
		Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery. \$ 255
<input type="checkbox"/>	1506310	Set at 10"
<input type="checkbox"/>	1506311	Set at 11"
<input type="checkbox"/>	1506312	Set at 12"
<input type="checkbox"/>	1506313	Set at 13"
<input type="checkbox"/>	1506314	Set at 14"

Armrests		
<input type="checkbox"/>	1507095	Curved Armpads Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.
		Regular Armpads
<input type="checkbox"/>	1507100	Desk Length - 2" x 10"
<input type="checkbox"/>	1507101	Full Length - 2" x 14"
<input type="checkbox"/>	1507102	Full Length Premium Grade - 2" x 14" \$ 35
		Wide Armpads May interfere with wheels in some tilt configurations. \$ 115
<input type="checkbox"/>	1507110	Desk Length - 3" x 10"
<input type="checkbox"/>	1507111	Full Length - 3" x 12"
		Gel Armpads May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1507120	Full Length - 2" x 12"
<input type="checkbox"/>	1507121	Full Length - 3.5" x 12"
		Contour Arm Trough with Hand Pad \$ 125
		Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.
<input type="checkbox"/>	1507130	Right (4.5" x 18")
<input type="checkbox"/>	1507131	Left (4.5" x 18")

Seat Pan		
<input type="checkbox"/>	1508000	No Seat Pan
<input type="checkbox"/>	1508001	Solid Seat Pan \$ 120

Cushion		
<input type="checkbox"/>	1509000	No Cushion
		nxt BioFit
<input type="checkbox"/>	1509101	12"-20" width, 14"-20" depth  \$ 595
<input type="checkbox"/>	1509102	22"-24" width, 14"-20" depth \$ 625
		Cushion Size Available in even number sizes only.
1509400	Specify Size (Width" x Depth") _____	

Seat Height

Select desired Seat-to-Floor Height

Standard T50 Options

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
12"					18"	19"
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	8" Casters & 8" Forks					
12"					18"	19"
16"				16"		20"
20"				16"	17"	18"
22"					17"	18"
24"					18"	19"

RW	6" Soft Roll & Ultra Sport Forks					
16"				16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

Wide T50 Options

Increases the overall width of the Chair by 2".

RW	5" Casters					
12"	13"					
16"	13"	14"	15"			
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	

RW	6" Casters & 6" Forks					
16"		14"	15"			
20"		14"	15"	16"	17"	
22"		14"	15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	6" Casters & 8" Forks					
16"			15"			
20"			15"	16"	17"	
22"			15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	8" Casters & 8" Forks					
20"				16"	17"	
22"				16"	17"	18"
24"				16"	17"	18"

- Standard
 Wide

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

- | | | | | |
|--------------------------|---------|------------------------|----|-----|
| <input type="checkbox"/> | 1511051 | 5" x 1.25" Caster size | | |
| <input type="checkbox"/> | 1511061 | 6" x 1.25" Caster size | | |
| <input type="checkbox"/> | 1511063 | 6" x 1.50" Caster size | \$ | 85 |
| <input type="checkbox"/> | 1511081 | 8" x 1.25" Caster size | | |
| <input type="checkbox"/> | 1511084 | 8" x 2.00" Caster size | \$ | 135 |

Pneumatic

- | | | | | |
|--------------------------|---------|--|----|-----|
| <input type="checkbox"/> | 1511161 | 6" x 1.25" Caster size | \$ | 110 |
| <input type="checkbox"/> | 1511181 | 8" x 1.25" Caster size | \$ | 120 |
| <input type="checkbox"/> | 1511184 | 8" x 2.00" Caster size | \$ | 150 |
| <input type="checkbox"/> | 1511260 | 6" Soft Roll | \$ | 95 |
| <input type="checkbox"/> | 1511400 | Ultra Sport Fork Not compatible with 12" Rear Wheels. Requires 6" Soft Roll Caster. | \$ | 395 |

Rear Wheels

Composite			
<input type="checkbox"/> 1512012	12" Wheel size	<input type="checkbox"/> 1512016	16" Wheel size
<input type="checkbox"/> 1512022	22" Wheel size	<input type="checkbox"/> 1512024	24" Wheel size
<input type="checkbox"/> 1512020	20" Wheel size		
Spoke			\$ 75
<input type="checkbox"/> 1512120	20" Wheel size	<input type="checkbox"/> 1512122	22" Wheel size
<input type="checkbox"/> 1512224	24" Heavy Duty Spoke	<input type="checkbox"/> 1512124	24" Wheel size
			\$ 350

Axle Type

<input type="checkbox"/> 1513000	Quick Release Rear Axles
<input type="checkbox"/> 1513010	Fixed Rear Axles

Rear Tires

<input type="checkbox"/> 1514000	Urethane	
<input type="checkbox"/> 1514010	Pneumatic	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. \$ 60
<input type="checkbox"/> 1514020	Pneumatic with Flat-Free Inserts	\$ 95
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. Required for 350 lbs Weight Capacity.		

Hand Rims

<input type="checkbox"/> 1515000	No Hand Rims
<input type="checkbox"/> 1515010	Standard 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 1515020	Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks

<input type="checkbox"/> 1516005	Foot Operated Wheel Lock	Not compatible with One Arm Drive, O2 Holder. \$ 200
<input type="checkbox"/> 1516010	Push to Lock	
<input type="checkbox"/> 1516020	Pull to Lock	
<input type="checkbox"/> 1516030	Attendant Wheel Lock (Push to Lock)	Not compatible with Wheel Lock Extensions. \$ 85
<input type="checkbox"/> 1516040	User & Attendant Wheel Lock	\$ 125
<input type="checkbox"/> 1516100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock. \$ 60

Anti-Tippers

<input checked="" type="checkbox"/> 1517010	Rear Anti-Tippers	Required Safety Option. \$ 60
---	--------------------------	-------------------------------

Back Post Style

Push Bar Back		
<input type="checkbox"/> 1518010	Straight	
<input type="checkbox"/> 1518030	8° Bend	\$ 75
Cane Handle Back Not compatible with Headrest Mounting Plate.		
<input type="checkbox"/> 1518020	Straight	
<input type="checkbox"/> 1518040	8° Bend	\$ 75
Back Post Height Seat Pan to Handle.		
<input type="checkbox"/> 1518420	20" Height	<input type="checkbox"/> 1518425 25" Height

Back Options		
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 1519010 Button on Outside	<input type="checkbox"/> 1519011 Button on Inside	
Removable Cane Handle Push Bar Not compatible with Push Bar Back.		\$ 140
Seat to Back Angle		
<input type="checkbox"/> 1519090 Set at 90°	<input type="checkbox"/> 1519096 Set at 96°	<input type="checkbox"/> 1519102 Set at 102°
<input type="checkbox"/> 1519108 Set at 108°	<input type="checkbox"/> 1519114 Set at 114°	<input type="checkbox"/> 1519120 Set at 120°
<input type="checkbox"/> 1519220 Fold Down Back		
<input type="checkbox"/> 1519221 Non Fold Down Back		\$ 75

Back Upholstery		
<input type="checkbox"/> 1520000 No Back Upholstery		
Standard Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		\$ 155
<input type="checkbox"/> 1520018 Set at 18"	<input type="checkbox"/> 1520019 Set at 19"	<input type="checkbox"/> 1520020 Set at 20"
<input type="checkbox"/> 1520021 Set at 21"	<input type="checkbox"/> 1520022 Set at 22"	<input type="checkbox"/> 1520023 Set at 23"
Adjustable Tension Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		\$ 295
<input type="checkbox"/> 1520321 Set at 21"	<input type="checkbox"/> 1520322 Set at 22"	<input type="checkbox"/> 1520323 Set at 23"

Backrest		
<input type="checkbox"/> 1521000 No Backrest		
nxt Optima HA Active Contour		
<input type="checkbox"/> 1521100 16"-20" width, 16"-20" height		\$ 690
<input type="checkbox"/> 1521105 22"-24" width, 16"-20" height		\$ 720
nxt Optima HA Plus Contour		
<input type="checkbox"/> 1521101 16"-20" width, 16"-20" height		\$ 690
<input type="checkbox"/> 1521106 22"-24" width, 16"-20" height		\$ 720
Backrest Size Available in even number sizes only.		
1521400 Specify Size (Width" x Height")	_____	

Tilt/Recline Options		
Attendant Tilt only		
<input type="checkbox"/> 1523011 Trigger on Right Back Post	<input type="checkbox"/> 1523012 Trigger on Left Back Post	
Tilt Limit		
<input type="checkbox"/> 1523220 Set at 20°	<input type="checkbox"/> 1523225 Set at 25°	<input type="checkbox"/> 1523230 Set at 30°
<input type="checkbox"/> 1523235 Set at 35°	<input type="checkbox"/> 1523240 Set at 40°	<input type="checkbox"/> 1523245 Set at 45°
<input type="checkbox"/> 1523250 Set at 50°		
<input type="checkbox"/> 1523410 Reclining Backrest (Manual)		\$ 875
30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Tension Adjustable Upholstery. Requires 25" Angle/Height Adjustable Push Bar with 8" Bend. Not available with 15", 16", 21", or 22" Seat Depth.		
Power Tilt 4° to 45° Tilt range, 250 lbs Weight Capacity.		\$ 1,695
Toggle Switch	<input type="checkbox"/> 1523421 Switch installed on Right side	<input type="checkbox"/> 1523422 Switch installed on Left side
Button Switch	<input type="checkbox"/> 1523431 Switch installed on Right side	<input type="checkbox"/> 1523432 Switch installed on Left side

Headrest		
<input type="checkbox"/> 1524000 No Headrest		
<input type="checkbox"/> 1524010 Headrest Mounting Plate		\$ 75
Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.		
<input type="checkbox"/> 1524100 nxt MultiFit Head Support Hardware		\$ 155
nxt MultiFit Head Support Pad		\$ 195
<input type="checkbox"/> 1524105 Regular (9" x 5")	<input type="checkbox"/> 1524106 Large (10" x 5")	<input type="checkbox"/> 1524107 X-Large (14" x 5")

Front Rigging

<input type="checkbox"/>	1525000	No Front Rigging		
<input type="checkbox"/>	1525001	Sector Blocks only		
<input type="checkbox"/>	1525110	Plug in 70°	Includes required Frame upgrade.	\$ 300
<input type="checkbox"/>	1525120	Plug in Adjustable Knee Angle	Includes required Frame upgrade.	\$ 1,275
<input type="checkbox"/>	1525140	Swing Away 70°		
<input type="checkbox"/>	1525160	Swing Away 90°	Not compatible with 8" Casters.	\$ 55
<input type="checkbox"/>	1525200	Swing Away Elevating Legrest		\$ 375
<input type="checkbox"/>	1525300	1" Offset Hangers for Swing Away		\$ 75

Foot Plates

<input type="checkbox"/>	1526010	Composite		
<input type="checkbox"/>	1526020	Aluminum Angle Adjustable		\$ 155
<input type="checkbox"/>	1526040	One Piece Flip-up	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1526050	Interlocking Angle Adjustable	Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1527015	Swing Away 70° Calf Pad	Not compatible with 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1527020	Swing Away 70° Gel Padding	Not compatible with 70° Calf Pad.	\$ 100
<input type="checkbox"/>	1527030	Swing Away ELR Gel Knee Pad		\$ 100
<input type="checkbox"/>	1527040	Adjustable Length Heel Loops		\$ 50
<input type="checkbox"/>	1527050	Bodypoint - Padded Calf Strap		\$ 105
Miller's Dynamic Footrest Mounts				
<input type="checkbox"/>	1527060	For Plug In Front Rigging		\$ 190
<input type="checkbox"/>	1527065	For Swing Away 70° Front Rigging	Not available on 13" STF Height.	\$ 380

Extension Tubes

Regular Extension Tubes Not compatible with Swing Away 90°.				
<input type="checkbox"/>	1528010	Standard (13"-17")		
<input type="checkbox"/>	1528020	Long (15"-19")		\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only.				
<input type="checkbox"/>	1528110	Standard (13"-17")	<input type="checkbox"/>	1528120 Long (15"-19")
<input type="checkbox"/>	1528200	90° Extension Tubes (11"-16")	Not compatible with 70° or ELR.	\$ 75

Positioning Straps

<input type="checkbox"/>	1529005	Lap Belt Mount		\$ 55
Bodypoint - Hip Belts				
<input type="checkbox"/>	1529010	Center-Pull Auto Style		\$ 105
<input type="checkbox"/>	1529015	Center-Pull Auto Style, Padded		\$ 155
<input type="checkbox"/>	1529020	4-Point, Padded		\$ 175
<input type="checkbox"/>	1529030	Bodypoint - Chest Strap		\$ 85

Miscellaneous Options

<input type="checkbox"/>	1530010 One-Arm Drive		\$ 3,500
	Requires 24" Composite Rear Wheel. Forward Wheel Position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks.		
<input type="checkbox"/>	1530010 Installed on Right side	<input type="checkbox"/>	1530011 Installed on Left side
<input type="checkbox"/>	1530060 Stealth Flip Down Abductor		\$ 495
<input type="checkbox"/>	1530070 QR Bracket - Accessory Mount		\$ 200
<input type="checkbox"/>	1530080 Gel Hip Guides	Not compatible with Dual Post Arms.	\$ 325
<input type="checkbox"/>	1530110 IV Pole		\$ 295
<input type="checkbox"/>	1530120 O2 Holder	May add to overall width of the Chair. Not compatible with Foot Operated Wheel Lock.	\$ 395
<input type="checkbox"/>	1530150 Ventilator Tray	16" minimum Seat Width. Reduces ground clearance. Chair may be unstable when unoccupied. Not compatible with O2 Holder.	\$ 295
<input type="checkbox"/>	1530152 Ventilator Bracket	Compatible with LTV1 (laptop) ventilator.	\$ 295
<input type="checkbox"/>	1530154 Ventilator Battery Tray	16" minimum Seat Width. Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with STF restrictions.	\$ 95
<input type="checkbox"/>	1530160 Clear Lap Tray (Slide On)	Not compatible with 8"- 11" Dual Post Arms, Desk Length Pads, Gel Arm Pads. Seat Widths of 15"-20" only.	\$ 490
<input type="checkbox"/>	1530170 Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10
<input type="checkbox"/>	1530300 Fuze T50 to T20 Conversion Kit		\$ 995

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/fuze-T50 to learn more about the Fuze T50 50° Manual Tilt Wheelchair

