

Customer and Order Information

Purchase Order *	Tag	v.061518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 1501000	Fuze T50 50° Manual Tilt Wheelchair	\$ 3,995
---	--	----------

Frame Options

Fixed Base and Seat Frame Select Fixed Seat Width below.

<input type="checkbox"/> 1502010	300 lbs Weight Capacity	
<input type="checkbox"/> 1502011	350 lbs Weight Capacity	Not compatible with Transit Tie Down System. Requires Flat-free Rear Tires or Heavy Duty Spoke Rear Wheels.
<input type="checkbox"/> 1502025	Adjustable Base and Seat Frame	250 lbs Weight Capacity. Select Adjustable Seat Width below. \$ 100
<input type="checkbox"/> 1502300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.

Seat Width Overall chair width is +7.5" for 12/16" RW and +8.5" for 20/22/24" RW

Fixed Seat Widths Order with Fixed Base and Seat Frame.

<input type="checkbox"/> 1503015	15" Seat Width	<input type="checkbox"/> 1503016	16" Seat Width	<input type="checkbox"/> 1503017	17" Seat Width	
<input type="checkbox"/> 1503018	18" Seat Width	<input type="checkbox"/> 1503019	19" Seat Width	<input type="checkbox"/> 1503020	20" Seat Width	
<input type="checkbox"/> 1503021	21" Seat Width	<input type="checkbox"/> 1503022	22" Seat Width			\$ 215
<input type="checkbox"/> 1503023	23" Seat Width	<input type="checkbox"/> 1503024	24" Seat Width			\$ 325

Adjustable Seat Widths Order with Adjustable Base and Seat Frame.

<input type="checkbox"/> 1503112	12"-15" Range - Set at 12"	<input type="checkbox"/> 1503113	12"-15" Range - Set at 13"	<input type="checkbox"/> 1503114	12"-15" Range - Set at 14"	
<input type="checkbox"/> 1503115	12"-15" Range - Set at 15"					
<input type="checkbox"/> 1503215	15"-20" Range - Set at 15"	<input type="checkbox"/> 1503216	15"-20" Range - Set at 16"	<input type="checkbox"/> 1503217	15"-20" Range - Set at 17"	
<input type="checkbox"/> 1503218	15"-20" Range - Set at 18"	<input type="checkbox"/> 1503219	15"-20" Range - Set at 19"	<input type="checkbox"/> 1503220	15"-20" Range - Set at 20"	
<input type="checkbox"/> 1503321	21"-22" Range - Set at 21"	<input type="checkbox"/> 1503322	21"-22" Range - Set at 22"			\$ 215

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 1504015	15" Seat Depth	<input type="checkbox"/> 1504016	16" Seat Depth	<input type="checkbox"/> 1504017	17" Seat Depth	
<input type="checkbox"/> 1504018	18" Seat Depth	<input type="checkbox"/> 1504019	19" Seat Depth	<input type="checkbox"/> 1504020	20" Seat Depth	

Custom Seat Depths Requires 25" Push Bar Back Straight, Angle/Height Adjustable Push Bar. Not compatible with Reclining or Fold Down Back. \$ 325

<input type="checkbox"/> 1504021	21" Seat Depth	<input type="checkbox"/> 1504022	22" Seat Depth		
----------------------------------	----------------	----------------------------------	----------------	--	--

Frame Colour

PDG Premium Colours

<input type="checkbox"/> 1505000	Tracer Red	<input type="checkbox"/> 1505001	Gloss Black	<input type="checkbox"/> 1505002	Silver
<input type="checkbox"/> 1505003	Black Shimmer	<input type="checkbox"/> 1505004	Orange Fusion	<input type="checkbox"/> 1505005	Karl's Imagination
<input type="checkbox"/> 1505006	Blue Shimmer	<input type="checkbox"/> 1505007	Gloss White	<input type="checkbox"/> 1505008	Dark Titanium
<input type="checkbox"/> 1505009	Geyser Green	<input type="checkbox"/> 1505010	Lollypop Blue	<input type="checkbox"/> 1505011	Lollypop Red
<input type="checkbox"/> 1505012	Bubble Gum Pink	<input type="checkbox"/> 1505013	Berry Purple	<input type="checkbox"/> 1505014	Safety Yellow
<input type="checkbox"/> 1505200	Custom Colour	Contact PDG for specifications.			\$ 300

Arm Type			
<input type="checkbox"/>	1506000	No Armrest	
<input type="checkbox"/>	1506100	Single Post Fixed Height (12")	
<input type="checkbox"/>	1506114	Single Post Adjustable Height (10"-14")	\$ 255
<input type="checkbox"/>	1506200	Dual Post Fixed Height (12") Full Length Armpad only.	
		Dual Post Adjustable Height	\$ 255
<input type="checkbox"/>	1506214	10"-14"	<input type="checkbox"/>
			<input type="checkbox"/> 1506216 12"-16"
		Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.	\$ 255
<input type="checkbox"/>	1506310	Set at 10"	<input type="checkbox"/>
			<input type="checkbox"/> 1506311 Set at 11"
			<input type="checkbox"/> 1506312 Set at 12"
<input type="checkbox"/>	1506313	Set at 13"	<input type="checkbox"/>
			<input type="checkbox"/> 1506314 Set at 14"

Armrests			
<input type="checkbox"/>	1507095	Curved Armpads Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.	
		Regular Armpads	
<input type="checkbox"/>	1507100	Desk Length - 2" x 10"	<input type="checkbox"/>
			<input type="checkbox"/> 1507101 Full Length - 2" x 14"
<input type="checkbox"/>	1507102	Full Length Premium Grade - 2" x 14"	\$ 35
		Wide Armpads May interfere with wheels in some tilt configurations.	\$ 115
<input type="checkbox"/>	1507110	Desk Length - 3" x 10"	<input type="checkbox"/>
			<input type="checkbox"/> 1507111 Full Length - 3" x 12"
		Gel Armpads May interfere with wheels in some tilt configurations.	\$ 125
<input type="checkbox"/>	1507120	Full Length - 2" x 12"	<input type="checkbox"/>
			<input type="checkbox"/> 1507121 Full Length - 3.5" x 12"
		Contour Arm Trough with Hand Pad	\$ 125
		Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.	
<input type="checkbox"/>	1507130	Right (4.5" x 18")	<input type="checkbox"/>
			<input type="checkbox"/> 1507131 Left (4.5" x 18")

Seat Pan			
<input type="checkbox"/>	1508000	No Seat Pan	
<input type="checkbox"/>	1508001	Solid Aluminum Seat Pan	\$ 120

Cushion			
<input type="checkbox"/>	1509000	No Cushion	
		nxt BioFit	
<input type="checkbox"/>	1509101	12"-20" width, 14"-20" depth	\$ 595
<input type="checkbox"/>	1509102	22"-24" width, 14"-20" depth	\$ 625
		Cushion Size Available in even number sizes only.	
<input type="checkbox"/>	1509400	Specify Size (Width" x Depth") _____	



Seat Height

Select desired Seat-to-Floor Height

Standard T50 Options

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
12"					18"	19"
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	8" Casters & 8" Forks					
12"					18"	19"
16"				16"		20"
20"				16"	17"	18"
22"					17"	18"
24"					18"	19"

RW	6" Soft Roll & Ultra Sport Forks					
16"				16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

Wide T50 Options

Increases the overall width of the Chair by 2".

RW	5" Casters					
12"	13"					
16"	13"	14"	15"			
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	

RW	6" Casters & 6" Forks					
16"		14"	15"			
20"		14"	15"	16"	17"	
22"		14"	15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	6" Casters & 8" Forks					
16"			15"			
20"			15"	16"	17"	
22"			15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	8" Casters & 8" Forks					
20"				16"	17"	
22"				16"	17"	18"
24"				16"	17"	18"

- Standard
 Wide

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

- 1511051 5" x 1.25" Caster size
- 1511061 6" x 1.25" Caster size
- 1511063 6" x 1.50" Caster size \$ 85
- 1511081 8" x 1.25" Caster size
- 1511084 8" x 2.00" Caster size \$ 135

Pneumatic

- 1511161 6" x 1.25" Caster size \$ 110
- 1511181 8" x 1.25" Caster size \$ 120
- 1511184 8" x 2.00" Caster size \$ 150

- 1511260 **6" Soft Roll** \$ 95

- 1511400 **Ultra Sport Fork** Not compatible with 12" Rear Wheels. Requires 6" Soft Roll Caster. \$ 395

Rear Wheels

Composite

- | | | |
|---|---|---|
| <input type="checkbox"/> 1512012 12" Wheel size | <input type="checkbox"/> 1512016 16" Wheel size | <input type="checkbox"/> 1512020 20" Wheel size |
| <input type="checkbox"/> 1512022 22" Wheel size | <input type="checkbox"/> 1512024 24" Wheel size | |

Spoke

- | | | | |
|--|---|---|--------|
| <input type="checkbox"/> 1512120 20" Wheel size | <input type="checkbox"/> 1512122 22" Wheel size | <input type="checkbox"/> 1512124 24" Wheel size | \$ 75 |
| <input type="checkbox"/> 1512224 24" Heavy Duty Spoke | | | \$ 350 |

Axle Type

- | |
|--|
| <input type="checkbox"/> 1513000 Quick Release Rear Axles |
| <input type="checkbox"/> 1513010 Fixed Rear Axles |

Rear Tires

- | | | |
|--|--|-------|
| <input type="checkbox"/> 1514000 Urethane | | |
| <input type="checkbox"/> 1514010 Pneumatic | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 1514020 Pneumatic with Flat-Free Inserts | Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. Required for 350 lbs Weight Capacity. | \$ 95 |

Hand Rims

- | | | |
|--|-------------------------------------|--------|
| <input type="checkbox"/> 1515000 No Hand Rims | | |
| <input type="checkbox"/> 1515010 Standard | 20", 22", and 24" Rear Wheels only. | |
| <input type="checkbox"/> 1515020 Plastic Coated | 20", 22", and 24" Rear Wheels only. | \$ 135 |

Wheel Locks

- | | | |
|---|--|--------|
| <input type="checkbox"/> 1516005 Foot Operated Wheel Lock | Not compatible with One Arm Drive. | \$ 200 |
| <input type="checkbox"/> 1516010 Push to Lock | | |
| <input type="checkbox"/> 1516020 Pull to Lock | | |
| <input type="checkbox"/> 1516030 Attendant Wheel Lock (Push to Lock) | Not compatible with Wheel Lock Extensions. | \$ 85 |
| <input type="checkbox"/> 1516040 User & Attendant Wheel Lock | | \$ 125 |
| <input type="checkbox"/> 1516100 Wheel Lock Extensions | Not compatible with Attendant Wheel Lock. | \$ 60 |

Anti-Tippers

- | | | |
|--|-------------------------|-------|
| <input checked="" type="checkbox"/> 1517010 Rear Anti-Tippers | Required Safety Option. | \$ 60 |
|--|-------------------------|-------|

Back Post Style

Push Bar Back

- | | | |
|---|--|-------|
| <input type="checkbox"/> 1518010 Straight | | |
| <input type="checkbox"/> 1518030 8° Bend | | \$ 75 |

Cane Handle Back

Not compatible with Headrest Mounting Plate.

- | | | |
|---|--|-------|
| <input type="checkbox"/> 1518020 Straight | | |
| <input type="checkbox"/> 1518040 8° Bend | | \$ 75 |

Back Post Height

Seat Pan to Handle.

- | | |
|---|---|
| <input type="checkbox"/> 1518420 20" Height | <input type="checkbox"/> 1518425 25" Height |
|---|---|

Back Options		
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 1519010 Button on Outside	<input type="checkbox"/> 1519011 Button on Inside	
Removable Cane Handle Push Bar Not compatible with Push Bar Back.		\$ 140
Seat to Back Angle		
<input type="checkbox"/> 1519090 Set at 90°	<input type="checkbox"/> 1519096 Set at 96°	<input type="checkbox"/> 1519102 Set at 102°
<input type="checkbox"/> 1519108 Set at 108°	<input type="checkbox"/> 1519114 Set at 114°	<input type="checkbox"/> 1519120 Set at 120°
<input type="checkbox"/> 1519220 Fold Down Back		
<input type="checkbox"/> 1519221 Non Fold Down Back		\$ 75

Back Upholstery		
<input type="checkbox"/> 1520000 No Back Upholstery		
Standard Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		\$ 155
<input type="checkbox"/> 1520018 Set at 18"	<input type="checkbox"/> 1520019 Set at 19"	<input type="checkbox"/> 1520020 Set at 20"
<input type="checkbox"/> 1520021 Set at 21"	<input type="checkbox"/> 1520022 Set at 22"	<input type="checkbox"/> 1520023 Set at 23"
Adjustable Tension Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		\$ 295
<input type="checkbox"/> 1520321 Set at 21"	<input type="checkbox"/> 1520322 Set at 22"	<input type="checkbox"/> 1520323 Set at 23"

Backrest		
<input type="checkbox"/> 1521000 No Backrest		
nxt Optima HA Active Contour		
<input type="checkbox"/> 1521100 16"-20" width, 16"-20" height		\$ 690
<input type="checkbox"/> 1521105 22"-24" width, 16"-20" height		\$ 720
nxt Optima HA Plus Contour		
<input type="checkbox"/> 1521101 16"-20" width, 16"-20" height		\$ 690
<input type="checkbox"/> 1521106 22"-24" width, 16"-20" height		\$ 720
Backrest Size Available in even number sizes only.		
1521400 Specify Size (Width" x Height")	_____	

Tilt/Recline Options		
Attendant Tilt only		
<input type="checkbox"/> 1523011 Trigger on Right Back Post	<input type="checkbox"/> 1523012 Trigger on Left Back Post	
Tilt Limit		
<input type="checkbox"/> 1523220 Set at 20°	<input type="checkbox"/> 1523225 Set at 25°	<input type="checkbox"/> 1523230 Set at 30°
<input type="checkbox"/> 1523235 Set at 35°	<input type="checkbox"/> 1523240 Set at 40°	<input type="checkbox"/> 1523245 Set at 45°
<input type="checkbox"/> 1523250 Set at 50°		
<input type="checkbox"/> 1523410 Reclining Backrest (Manual)		\$ 875
30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Tension Adjustable Upholstery. Requires 25" Angle/Height Adjustable Push Bar with 8" Bend. Not available with 15", 16", 21", or 22" Seat Depth.		
Power Tilt 0° to 40° tilt range. Includes 2 batteries for off-chair charging.		
<input type="checkbox"/> 1523400 250 lbs Weight Capacity		\$ 2,195
<input type="checkbox"/> 1523401 350 lbs Weight Capacity Requires Fixed Base and Seat Frame		\$ 2,495
Toggle Switch	<input type="checkbox"/> 1523421 Installed on Right side	<input type="checkbox"/> 1523422 Installed on Left side
Button Switch	<input type="checkbox"/> 1523431 Installed on Right side	<input type="checkbox"/> 1523432 Installed on Left side
Egg Switch	<input type="checkbox"/> 1523440 Not installed	
Optional with Egg Switch:	<input type="checkbox"/> 1523445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models. \$ 625
Specify Headrest to be used: _____		

Headrest			
<input type="checkbox"/>	1524000	No Headrest	
<input type="checkbox"/>	1524010	Headrest Mounting Plate Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 75
<input type="checkbox"/>	1524100	nxt MultiFit Head Support Hardware	\$ 235
<input type="checkbox"/>		nxt MultiFit Head Support Pad	\$ 245
<input type="checkbox"/>	1524105	Regular (9" x 5")	
<input type="checkbox"/>	1524106	Large (10" x 5")	
<input type="checkbox"/>	1524107	X-Large (14" x 5")	

Front Rigging			
<input type="checkbox"/>	1525000	No Front Rigging	
<input type="checkbox"/>	1525001	Sector Blocks only	
<input type="checkbox"/>		PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.	
<input type="checkbox"/>	1525410	Fixed 70°	\$ 95
<input type="checkbox"/>	1525420	Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____	\$ 195
<input type="checkbox"/>	1525430	Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.	\$ 495
<input type="checkbox"/>		Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.	
<input type="checkbox"/>	1525140	Swing Away 70°	
<input type="checkbox"/>	1525160	Swing Away 90° Not compatible with 8" Casters.	\$ 55
<input type="checkbox"/>	1525200	Swing Away Elevating Legrest	\$ 375
<input type="checkbox"/>	1525300	1" Offset Hangers for Swing Away	\$ 75

Foot Plates			
<input type="checkbox"/>	1526010	Composite	
<input type="checkbox"/>	1526020	Aluminum Angle Adjustable	\$ 155
<input type="checkbox"/>	1526040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1526050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1527015	Swing Away 70° Calf Pad Not compatible with 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1527020	Swing Away 70° Gel Padding Not compatible with 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1527030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1527040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	1527050	Bodypoint - Padded Calf Strap	\$ 105
<input type="checkbox"/>	1527060	Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.	\$ 190

Extension Tubes			
<input type="checkbox"/>	1528010	Short 14"-18" for Easy-Out Fixed & SA Angle Adjustable 11"-15" for Easy-Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1528020	Long 18"-24" for Easy-Out Fixed & SA Angle Adjustable 15"-21" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
<input type="checkbox"/>	1628200	For 90° FR (11"-16") Not compatible with Easy-Out, Standard 70° or ELR	\$ 75

Positioning Straps			
<input type="checkbox"/>	1529005	Lap Belt Mount	\$ 55
<input type="checkbox"/>		Bodypoint - Hip Belts	
<input type="checkbox"/>	1529010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1529015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1529020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1529030	Bodypoint - Chest Strap	\$ 85

Miscellaneous Options		
One-Arm Drive Requires 24" Composite Rear Wheels and forward wheel position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks.		
<input type="checkbox"/> 1530010	Installed on Right side	<input type="checkbox"/> 1530011 Installed on Left side
<input type="checkbox"/> 1530060	Stealth Flip Down Abductor	
<input type="checkbox"/> 1530070	QR Bracket - Accessory Mount	
<input type="checkbox"/> 1530080	Gel Hip Guides Not compatible with Dual Post Arms. Reduces available Seat Width by 2".	
<input type="checkbox"/> 1530110	IV Pole	
<input type="checkbox"/> 1530120	O2 Holder May add to overall width of the Chair.	
<input type="checkbox"/> 1530150	Ventilator Tray Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 16". Not compatible with O2 Holder, Power Tilt.	
<input type="checkbox"/> 1530152	Ventilator Bracket Compatible with LTV1 (laptop) ventilator.	
<input type="checkbox"/> 1530154	Ventilator Battery Tray	
16" minimum Seat Width. Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with STF restrictions.		
<input type="checkbox"/> 1530160	Clear Lap Tray (Slide On) Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	
<input type="checkbox"/> 1530170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	
<input type="checkbox"/> 1530300	Fuze T50 to T20 Conversion Kit	

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/fuze-T50/ to learn more about the **Fuze T50 50° Manual Tilt Wheelchair**

