

Customer and Order Information

Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address *	<i>*required</i>

Frame Type

<input checked="" type="checkbox"/> 1901000 Stellar GL 45° Manual Tilt Wheelchair	\$ 3,995
--	----------

Frame Options

<input checked="" type="checkbox"/> 1902000 Standard Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/> 1902300 Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths		
<input type="checkbox"/> 1903014 14" Seat Width	<input type="checkbox"/> 1903015 15" Seat Width	<input type="checkbox"/> 1903016 16" Seat Width
<input type="checkbox"/> 1903017 17" Seat Width	<input type="checkbox"/> 1903018 18" Seat Width	<input type="checkbox"/> 1903019 19" Seat Width
<input type="checkbox"/> 1903020 20" Seat Width		
<input type="checkbox"/> 1903021 21" Seat Width	<input type="checkbox"/> 1903022 22" Seat Width	\$ 215

Seat Depth

Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 1904016 16" Seat Depth	<input type="checkbox"/> 1904017 17" Seat Depth	<input type="checkbox"/> 1904018 18" Seat Depth
<input type="checkbox"/> 1904019 19" Seat Depth	<input type="checkbox"/> 1904020 20" Seat Depth	
Custom Seat Depths		
<input type="checkbox"/> 1904021 21" Seat Depth	<input type="checkbox"/> 1904022 22" Seat Depth	\$ 325

Frame Colour

PDG Premium Colours		
<input type="checkbox"/> 1905000 Tracer Red	<input type="checkbox"/> 1905001 Gloss Black	<input type="checkbox"/> 1905002 Silver
<input type="checkbox"/> 1905003 Black Shimmer	<input type="checkbox"/> 1905004 Orange Fusion	<input type="checkbox"/> 1905005 Karl's Imagination
<input type="checkbox"/> 1905006 Blue Shimmer	<input type="checkbox"/> 1905007 Gloss White	<input type="checkbox"/> 1905008 Dark Titanium
<input type="checkbox"/> 1905009 Geyser Green	<input type="checkbox"/> 1905010 Lollipop Blue	<input type="checkbox"/> 1905011 Lollipop Red
<input type="checkbox"/> 1905012 Bubble Gum Pink	<input type="checkbox"/> 1905013 Berry Purple	<input type="checkbox"/> 1905014 Safety Yellow
<input type="checkbox"/> 1905200 Custom Colour Contact PDG for specifications.		\$ 300

Arm Type

<input type="checkbox"/> 1906000 No Armrest		
<input type="checkbox"/> 1906100 Single Post Fixed Height (12") Full Length Armpad only.		
Single Post Adjustable Height		\$ 255
<input type="checkbox"/> 1906111 8"-11"	<input type="checkbox"/> 1906114 10"-14"	
<input type="checkbox"/> 1906200 Dual Post Fixed Height (12") Full Length Armpad only.		
Dual Post Adjustable Height		\$ 255
<input type="checkbox"/> 1906211 8"-11"	<input type="checkbox"/> 1906214 10"-14"	<input type="checkbox"/> 1906216 12"-16"

See next page for more Arm Type options >

Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.			\$ 255
<input type="checkbox"/> 1906310	Set at 10"	<input type="checkbox"/> 1906311	Set at 11"
<input type="checkbox"/> 1906312	Set at 12"	<input type="checkbox"/> 1906313	Set at 13"
<input type="checkbox"/> 1906314	Set at 14"		

Armrests

Regular Armrests			
<input type="checkbox"/> 1907100	Desk Length - 2" x 10"	<input type="checkbox"/> 1907101	Full Length - 2" x 14"
<input type="checkbox"/> 1907102	Full Length Premium Grade - 2" x 14"		\$ 35
Wide Armrests May interfere with wheels in some tilt configurations.			\$ 115
<input type="checkbox"/> 1907110	Desk Length - 3" x 10"	<input type="checkbox"/> 1907111	Full Length - 3" x 12"
Gel Armrests May interfere with wheels in some tilt configurations.			\$ 125
<input type="checkbox"/> 1907120	Full Length - 2" x 12"	<input type="checkbox"/> 1907121	Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.			\$ 125
<input type="checkbox"/> 1907130	Right (4.5" x 18")	<input type="checkbox"/> 1907131	Left (4.5" x 18")
<input type="checkbox"/> 1907220	Armrest Wideners for Dual Post Armrests Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.		\$ 375
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.			
<input type="checkbox"/> 1907310	For Single Post Armrests Available for 10"-14" Height range only.		\$ 110
<input type="checkbox"/> 1907320	For Dual Post Armrests		\$ 165
Armrest Extensions Not compatible with Flip-up or Single Post Armrests.			\$ 60
<input type="checkbox"/> 1907400	3" Extension - Forward	<input type="checkbox"/> 1907401	3" Extension - Rearward

Seat Pan

<input type="checkbox"/> 1908000	No Seat Pan	
<input type="checkbox"/> 1908001	Solid Seat Pan	\$ 120

Cushion

<input type="checkbox"/> 1909000	No Cushion		
nxt BioFit			
<input type="checkbox"/> 1909101	14"-20" width, 16"-20" depth		\$ 595
<input type="checkbox"/> 1909102	22"-24" width, 16"-20" depth		\$ 625
Cushion Size Available in even number sizes only.			
<input type="checkbox"/> 1909400	Specify Size (Width" x Depth") _____		

Seat Height

Select desired Seat-to-Floor Height STF Height is adjustable within the range of the selected Rear Wheel size.																																																													
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>20"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td></td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table>	RW	6" Casters						12"	13"	14"	15"	16"			16"			15"	16"	17"	18"	20"	13"	14"	15"	16"			22"		14"	15"	16"	17"		24"			15"	16"	17"	18"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="4">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>24"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table>	RW	8" Casters				16"		15"	16"	17"	18"	24"		15"	16"	17"	18"	Wheel size _____ Caster size _____ STF height _____
RW	6" Casters																																																												
12"	13"	14"	15"	16"																																																									
16"			15"	16"	17"	18"																																																							
20"	13"	14"	15"	16"																																																									
22"		14"	15"	16"	17"																																																								
24"			15"	16"	17"	18"																																																							
RW	8" Casters																																																												
16"		15"	16"	17"	18"																																																								
24"		15"	16"	17"	18"																																																								

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/> 1911061	6" x 1.25" Caster size	
<input type="checkbox"/> 1911063	6" x 1.50" Caster size	\$ 85
<input type="checkbox"/> 1911081	8" x 1.25" Caster size	
<input type="checkbox"/> 1911084	8" x 2.00" Caster size	\$ 135

See next page for more Casters options >

Back Options

<input type="checkbox"/> 1919010 Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150	
<input type="checkbox"/> 1919010 Button on Outside	<input type="checkbox"/> 1919011 Button on Inside
Seat to Back Angle	
<input type="checkbox"/> 1919090 Set at 90°	<input type="checkbox"/> 1919098 Set at 98°
<input type="checkbox"/> 1919114 Set at 114°	<input type="checkbox"/> 1919122 Set at 122°
<input type="checkbox"/> 1919210 10" Backrest Extension Not compatible with Angle/Height Adjustable Push Bar Handle. \$ 135	
<input type="checkbox"/> 1919220 Fold Down Back Not compatible with Reclining Back.	
<input type="checkbox"/> 1919230 Spring Back Suspension Not compatible with Reclining Back. \$ 105	
<input type="checkbox"/> 1919240 Dynamic Rocker Back Assembly \$ 2,995	
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.	

Back Upholstery

<input type="checkbox"/> 1920000 No Back Upholstery	
<input type="checkbox"/> 1920000 Standard Back height measured from Seat Pan to top of Upholstery. \$ 155	
<input type="checkbox"/> 1920021 Set at 21"	<input type="checkbox"/> 1920022 Set at 22"
Short 25" Back only. Back height measured from Seat Pan to top of Upholstery. \$ 155	
<input type="checkbox"/> 1920118 Set at 18"	<input type="checkbox"/> 1920119 Set at 19"
<input type="checkbox"/> 1920120 Set at 20"	
Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery. \$ 155	
<input type="checkbox"/> 1920224 Set at 24"	<input type="checkbox"/> 1920225 Set at 25"
<input type="checkbox"/> 1920226 Set at 26"	
Adjustable Tension Back height measured from Seat Pan to top of Upholstery. \$ 295	
<input type="checkbox"/> 1920321 Set at 21"	<input type="checkbox"/> 1920322 Set at 22"
<input type="checkbox"/> 1920323 Set at 23"	

Backrest

<input type="checkbox"/> 1921000 No Backrest	
nxt Optima HA Active Contour	
<input type="checkbox"/> 1921100 16"-20" width, 16"-20" height \$ 690	
<input type="checkbox"/> 1921105 22"-24" width, 16"-20" height \$ 720	
nxt Optima HA Plus Contour	
<input type="checkbox"/> 1921101 16"-20" width, 16"-20" height \$ 690	
<input type="checkbox"/> 1921106 22"-24" width, 16"-20" height \$ 720	
Backrest Size Available in even number sizes only.	
1921400 Specify Size (Width" x Height") _____	

Tilt/Recline Options

Attendant Tilt only	
<input type="checkbox"/> 1923011 Trigger on Right Back Post	<input type="checkbox"/> 1923012 Trigger on Left Back Post
Tilt Limit	
<input type="checkbox"/> 1923220 Set at 20°	<input type="checkbox"/> 1923230 Set at 30°
<input type="checkbox"/> 1923240 Set at 40°	
<input type="checkbox"/> 1923410 Reclining Backrest (Manual) \$ 875	
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.	

Headrest			
<input type="checkbox"/>	1924000	No Headrest	
<input type="checkbox"/>	1924010	Headrest Mounting Plate Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 75
<input type="checkbox"/>	1924100	nxt MultiFit Head Support Hardware	\$ 155
		nxt MultiFit Head Support Pad	\$ 195
<input type="checkbox"/>	1924105	Regular (9" x 5")	
<input type="checkbox"/>	1924106	Large (10" x 5")	
<input type="checkbox"/>	1924107	X-Large (14" x 5")	

Front Rigging			
<input type="checkbox"/>	1925000	No Front Rigging	
<input type="checkbox"/>	1925001	Sector Blocks only	
<input type="checkbox"/>	1925110	Plug In 70°	
<input type="checkbox"/>	1925120	Plug In Adjustable Knee Angle	\$ 975
<input type="checkbox"/>	1925130	Swing Away 60°	
<input type="checkbox"/>	1925140	Swing Away 70°	
<input type="checkbox"/>	1925200	Swing Away Elevating Legrest	\$ 375
<input type="checkbox"/>	1925300	1" Offset Hangers for Swing Away	\$ 75

Foot Plates			
<input type="checkbox"/>	1926010	Composite	
<input type="checkbox"/>	1926020	Aluminum Angle Adjustable	\$ 155
<input type="checkbox"/>	1926030	One Piece Fixed Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1926040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1926050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1927010	Swing Away 60° & 70° Calf Pad Not compatible with 60° & 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1927020	Swing Away 60° & 70° Gel Padding Not compatible with 60° & 70° Swing Away Calf Pad.	\$ 100
<input type="checkbox"/>	1927030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1927040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	1927050	Bodypoint - Padded Calf Strap	\$ 105
		Miller's Dynamic Footrest Mounts	
<input type="checkbox"/>	1927060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1927065	For Swing Away 70° Front Rigging Not available on 13" STF Height.	\$ 380

Extension Tubes			
		Regular Extension Tubes	
<input type="checkbox"/>	1928010	Standard (13"-17")	
<input type="checkbox"/>	1928020	Long (15"-19")	\$ 75
		Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only.	\$ 100
<input type="checkbox"/>	1928110	Standard (13"-17")	
<input type="checkbox"/>	1928120	Long (15"-19")	

Positioning Straps

Bodypoint - Hip Belts

<input type="checkbox"/>	1929010	Center-Pull Auto Style		\$	105
<input type="checkbox"/>	1929015	Center-Pull Auto Style, Padded		\$	155
<input type="checkbox"/>	1929020	4-Point, Padded		\$	175
<input type="checkbox"/>	1929030	Bodypoint - Chest Strap		\$	85

Miscellaneous Options

<input type="checkbox"/>	1930080	Gel Hip Guides		\$	325
<input type="checkbox"/>	1930110	IV Pole	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$	295
<input type="checkbox"/>	1930120	O2 Holder	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$	395
<input type="checkbox"/>	1930170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$	10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar-GL to learn more about the Stellar GL 45° Manual Tilt Wheelchair

