

Customer and Order Information

Purchase Order *	Tag	v.0601
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in *required

Frame Type

<input checked="" type="checkbox"/> 2501000	Stellar GLT 45° Lightweight Manual Tilt Wheelchair Aluminum Frame. Not recommended for agitated users.	\$ 4,295
<input type="checkbox"/>	For agitated or heavier users, choose the standard Stellar GL with steel frame. See separate order form.	
<input type="checkbox"/>	Consider the Stellar Impact for highly agitated users and the Stellar HD for bariatric users. See separate order forms.	

Frame Options

<input checked="" type="checkbox"/> 2502000	Lightweight Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/> 2502300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths		
<input type="checkbox"/> 2503014 14" Seat Width	<input type="checkbox"/> 1903015 15" Seat Width	<input type="checkbox"/> 1903016 16" Seat Width
<input type="checkbox"/> 2503017 17" Seat Width	<input type="checkbox"/> 1903018 18" Seat Width	<input type="checkbox"/> 1903019 19" Seat Width
<input type="checkbox"/> 2503020 20" Seat Width	<input type="checkbox"/> 1903021 21" Seat Width	<input type="checkbox"/> 1903022 22" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 2504016 16" Seat Depth	<input type="checkbox"/> 1904017 17" Seat Depth	<input type="checkbox"/> 1904018 18" Seat Depth
<input type="checkbox"/> 2504019 19" Seat Depth	<input type="checkbox"/> 1904020 20" Seat Depth	
Custom Seat Depths		
<input type="checkbox"/> 2504021 21" Seat Depth	<input type="checkbox"/> 1904022 22" Seat Depth	\$ 350

Frame Color

PDG Premium Colors		
<input type="checkbox"/> 2505000 Tracer Red	<input type="checkbox"/> 1905001 Gloss Black	<input type="checkbox"/> 1905002 Silver
<input type="checkbox"/> 2505003 Black Shimmer	<input type="checkbox"/> 1905004 Orange Fusion	<input type="checkbox"/> 1905005 Karl's Imagination
<input type="checkbox"/> 2505006 Blue Shimmer	<input type="checkbox"/> 1905007 Gloss White	<input type="checkbox"/> 1905008 Dark Titanium
<input type="checkbox"/> 2505009 Geyser Green	<input type="checkbox"/> 1905010 Lollypop Blue	<input type="checkbox"/> 1905011 Lollypop Red
<input type="checkbox"/> 2505012 Bubble Gum Pink	<input type="checkbox"/> 1905013 Berry Purple	<input type="checkbox"/> 1905014 Safety Yellow
<input type="checkbox"/> 2505200 Custom Color Contact PDG for specifications.		\$ 300

Arm Type

<input type="checkbox"/> 2506000 No Armrest		
<input type="checkbox"/> 2506100 Single Post Fixed Height (12") Full Length Armpad only.		
Single Post Adjustable Height		\$ 275
<input type="checkbox"/> 2506111 8"-11"	<input type="checkbox"/> 1906114 10"-14"	

Armrests

Regular Armrests		
<input type="checkbox"/> 2507100	Desk Length - 2" x 10"	<input type="checkbox"/> 1907101 Full Length - 2" x 14"
Wide Armrests May interfere with wheels in some tilt configurations.		\$ 200
<input type="checkbox"/> 2507110	Desk Length - 3" x 10"	<input type="checkbox"/> 1907111 Full Length - 3" x 12"
Gel Armrests May interfere with wheels in some tilt configurations.		\$ 135
<input type="checkbox"/> 2507120	Full Length - 2" x 12"	<input type="checkbox"/> 1907121 Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.		\$ 125
<input type="checkbox"/> 2507130	Right (4.5" x 18")	<input type="checkbox"/> 1907131 Left (4.5" x 18")
<input type="checkbox"/> 2507310	Clothing Guards Available for 10"-14" Height range only.	\$ 110

Seat Pan

<input type="checkbox"/> 2508000	No Seat Pan	
<input type="checkbox"/> 2508001	Solid Aluminum Seat Pan	\$ 120

Cushion

Refer to the enclosed Seating & Positioning order form for PDG's on-chair Cushion options. Standard PDG discounts apply.

Seat Height

Select desired Seat-to-Floor Height STF Height is adjustable within the range of the selected Rear Wheel size.

RW	6" Casters					
12"	13"	14"	15"	16"		
16"			15"	16"	17"	18"
20"	13"	14"	15"	16"		
22"		14"	15"	16"	17"	
24"			15"	16"	17"	18"

RW	8" Casters				
16"		15"	16"	17"	18"
24"		15"	16"	17"	18"

Wheel size _____

Caster size _____

STF height _____

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/> 2511061	6" x 1.25" Caster size	
<input type="checkbox"/> 2511063	6" x 1.50" Caster size	\$ 85
<input type="checkbox"/> 2511081	8" x 1.25" Caster size	
<input type="checkbox"/> 2511084	8" x 2.00" Caster size	\$ 135
Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/> 2511161	6" x 1.25" Caster size	\$ 110
<input type="checkbox"/> 2511181	8" x 1.25" Caster size	\$ 120
<input type="checkbox"/> 2511184	8" x 2.00" Caster size	\$ 150
<input type="checkbox"/> 2511260	6" Soft Roll	\$ 95
<input type="checkbox"/> 2511400	Ultra Sport Fork Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	\$ 395

Rear Wheels

Composite		
<input type="checkbox"/> 2512012	12" Wheel size	<input type="checkbox"/> 1912016 16" Wheel size
<input type="checkbox"/> 2512022	22" Wheel size	<input type="checkbox"/> 1912024 24" Wheel size
<input type="checkbox"/> 1912020	20" Wheel size	
Spoke		\$ 75
<input type="checkbox"/> 2512120	20" Wheel size	<input type="checkbox"/> 1912122 22" Wheel size
<input type="checkbox"/> 1912124	24" Wheel size	

Axle Type

<input type="checkbox"/> 2513000	Quick Release Rear Axles
<input type="checkbox"/> 2513010	Fixed Rear Axles

Rear Tires			
<input type="checkbox"/>	2514000	Urethane	
<input type="checkbox"/>	2514010	Pneumatic May result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/>	2514020	Pneumatic with Flat-Free Inserts May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95

Hand Rims			
<input type="checkbox"/>	2515000	No Hand Rims	
<input type="checkbox"/>	2515010	Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	2515020	Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 135

Wheel Locks			
<input type="checkbox"/>	2516010	Push to Lock	
<input type="checkbox"/>	2516020	Pull to Lock	
<input type="checkbox"/>	2516030	Attendant Wheel Lock (Push to Lock)	\$ 85
<input type="checkbox"/>	2516040	User & Attendant Wheel Lock Available only on 20", 22", or 24" Rear Wheels.	\$ 125
<input type="checkbox"/>	2516100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/>	2516200	Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	2517010	Rear Anti-Tippers Required Safety Option.	\$ 60

Back Post Style			
<input checked="" type="checkbox"/>	2518010	Push Bar Back	
Back Post Height Seat Pan to Handle			
<input type="checkbox"/>	2518425	25" Height	<input type="checkbox"/>
	1918431	31" Height	

Back Options			
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.			\$ 150
<input type="checkbox"/>	2519010	Button on Outside	<input type="checkbox"/>
	1919011	Button on Inside	
Seat to Back Angle			
<input type="checkbox"/>	2519090	Set at 90°	<input type="checkbox"/>
	1919098	Set at 98°	<input type="checkbox"/>
	1919106	Set at 106°	
<input type="checkbox"/>	2519114	Set at 114°	<input type="checkbox"/>
	1919122	Set at 122°	
<input type="checkbox"/>	2519210	10" Backrest Extension Requires the Angle/Height Adjustable Push Bar Handle.	\$ 150
<input checked="" type="checkbox"/>	2519220	Fold Down Back Not compatible with Reclining Back.	

Back Upholstery			
<input type="checkbox"/>	2520000	No Back Upholstery	
Standard Back height measured from Seat Pan to top of Upholstery.			\$ 175
<input type="checkbox"/>	2520021	Set at 21"	<input type="checkbox"/>
	1920022	Set at 22"	<input type="checkbox"/>
	1920023	Set at 23"	Not available with 25" Back.
Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 175
<input type="checkbox"/>	2520118	Set at 18"	<input type="checkbox"/>
	1920119	Set at 19"	<input type="checkbox"/>
	1920120	Set at 20"	
Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 175
<input type="checkbox"/>	2520224	Set at 24"	<input type="checkbox"/>
	1920225	Set at 25"	<input type="checkbox"/>
	1920226	Set at 26"	
Adjustable Tension Back height measured from Seat Pan to top of Upholstery.			\$ 295
<input type="checkbox"/>	2520321	Set at 21"	<input type="checkbox"/>
	1920322	Set at 22"	<input type="checkbox"/>
	1920323	Set at 23"	

Positioning Backrest	
	Refer to the enclosed Seating & Positioning order form for PDG's on-chair Backrest options. Standard PDG discounts apply.

Tilt/Recline Options

Attendant Tilt only

- 2523011 Trigger on Right Back Post 1923012 Trigger on Left Back Post

Tilt Limit

- 2523220 Set at 20° 1923230 Set at 30° 1923240 Set at 40°

- 2523310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.

- 2523400 **Power Tilt** 250 lbs Weight Capacity. 0° to 45° tilt range. Includes 2 batteries for off-chair charging. Minimum Seat Width is 16". \$ 1,695

Toggle Switch 1923421 Installed on Right side 1923422 Installed on Left side

Button Switch 1923431 Installed on Right side 1923432 Installed on Left side

Egg Switch 1923440 Not installed

Optional for Egg Switch: 1923445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625

Specify Headrest to be used: _____

Headrest

-  Refer to the enclosed Seating & Positioning order form for PDG's on-chair Headrest options. Standard PDG discounts apply.

- 2524010 **Headrest Mounting Plate** \$ 95

Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.

Front Rigging

- 2525000 **No Front Rigging** 2525001 **Sector Blocks only**

PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

- 2525410 Fixed 70° \$ 95

- 2525420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) \$ 295

Adjust angle by loosening one bolt. Tool-free leg-length adjustment.

Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.

- 2525130 Swing Away 60° 2525140 Swing Away 70°

- 2525200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375

- 2525300 1" Offset Hangers for Swing Away \$ 75

- 2525110 **Standard Plug In 70°** Legacy Design. Not compatible with Easy-Out.

Foot Plates

- 2526010 **Composite**

- 2526020 **Aluminum Angle Adjustable** \$ 155

- 2526030 **One Piece Fixed** Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350

- 2526040 **One Piece Flip-up** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395

- 2526050 **Interlocking Angle Adjustable** Seat Widths of 14"-21" only. \$ 195

Footrest/Legrest Accessories

Calf Pad Not compatible with Gel Padding.

- 2527012 For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount. \$ 200

- 2527010 For Swing Away Front Rigging \$ 135

- 2527020 **Gel Padding** Not compatible with Calf Pad or Easy-Out Front Rigging. \$ 125

Gel Knee Pad

- 2527032 For Easy-Out™ Front Rigging \$ 200

- 2527030 For Swing Away ELR Front Rigging \$ 100

- 2527040 **Adjustable Length Heel Loops** \$ 50

- 2527050 **Padded Calf Strap** \$ 195

Extension Tubes

Regular Extension Tubes

<input type="checkbox"/>	2528010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	2528020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75

Positioning Straps

Refer to the enclosed Seating & Positioning order form for PDG's on-chair Belts/Straps options. Standard PDG discounts apply.

Miscellaneous Options

<input type="checkbox"/>	2530080	Gel Hip Guides	Reduces available Seat Width by 2".		\$ 375
<input type="checkbox"/>	2530170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.		\$ 10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar-GL/ to learn more about the **Stellar GL 45° Manual Tilt Wheelchair**



Cushion

1. Select the type ↓	2. Select the model ↓
--------------------------------	---------------------------------

General Use	Stealth Simplicity 350 lbs weight capacity. <input type="checkbox"/> 4009141 16"-20" width, 16"-20" depth \$ 310 <input type="checkbox"/> 4009142 22"-24" width, 18"-20" depth \$ 370
--------------------	--

Positioning	ROHO AirLITE 300 lbs weight capacity. <input type="checkbox"/> 4009120 14"-20" width, 14"-18" depth \$ 375
--------------------	--

Skin Protection & Positioning	Stealth Spectrum Foam 650 lbs weight capacity. <input type="checkbox"/> 4009143 16"-20" width, 16"-24" depth \$ 580 <input type="checkbox"/> 4009144 22"-32" width, 18"-28" depth \$ 630
	Stealth Tru-Comfort 2 450 lbs weight capacity. <input type="checkbox"/> 4009146 16"-20" width, 16"-24" depth \$ 515 <input type="checkbox"/> 4009147 22" width, 18"-24" depth \$ 575
	nxt BioFit 300 lbs weight capacity. <input type="checkbox"/> 4009101 14"-20" width, 16"-20" depth \$ 595 <input type="checkbox"/> 4009102 22"-24" width, 16"-20" depth \$ 625
	Matrx VI 300 lbs weight capacity. <input type="checkbox"/> 4009111 14"-20" width, 15"-20" depth \$ 680 <input type="checkbox"/> 4009112 22" width, 18"-20" depth \$ 745
	Matrx VI Heavy Duty 600 lbs weight capacity. <input type="checkbox"/> 4009115 18"-21" width, 18"-26" depth \$ 1,050 <input type="checkbox"/> 4009116 22"-30" width, 18"-26" depth \$ 1,250

Skin Protection & Adjustable	ROHO Single Valve High Profile No weight limit if cushion is properly sized. <input type="checkbox"/> 4009125 14"-20" width, 16"-22" depth \$ 840 <input type="checkbox"/> 4009126 22"-24" width, 16"-22" depth \$ 960
---	---

Skin Protection, Positioning & Adjustable	ROHO Quadro Select High Profile No weight limit if cushion is properly sized. <input type="checkbox"/> 4009121 14"-20" width, 16"-22" depth \$ 840 <input type="checkbox"/> 4009122 22"-24" width, 16"-22" depth \$ 960
--	--

Extra-Wide No weight limit if cushion is properly sized. Custom sizes available. Custom-made foam cushion with fluid-proof yet breathable fabric cover.	PDG Solar Basic High resiliency polyfoam with 2" medium-density base and 1" lower-density top.	<input type="checkbox"/> 4009200 20"-28" width, 16"-24" depth \$ 395 <input type="checkbox"/> 4009201 29"-36" width, 16"-24" depth \$ 495 <input type="checkbox"/> 4009202 37"-48" width, 16"-24" depth \$ 595
	PDG Solar Contour Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups.	<input type="checkbox"/> 4009210 20"-28" width, 16"-24" depth \$ 695 <input type="checkbox"/> 4009211 29"-36" width, 16"-24" depth \$ 795 <input type="checkbox"/> 4009212 37"-48" width, 16"-24" depth \$ 895
	PDG Solar Pressure Reduction All features of the Contour, plus added pincore latex rubber and containment foam outline.	<input type="checkbox"/> 4009220 20"-28" width, 16"-24" depth \$ 1,095 <input type="checkbox"/> 4009221 29"-36" width, 16"-24" depth \$ 1,195 <input type="checkbox"/> 4009222 37"-48" width, 16"-24" depth \$ 1,295

3. Specify the size →	Cushion Size Sizes based on availability from the cushion manufacturer. Specify Size (Width" x Depth") _____
------------------------------	--

Backrest

1. Select the type ↓	2. Select the contour depth ↓	3. Select the model ↓
-------------------------	----------------------------------	--------------------------

Positioning, posterior-lateral	2" depth	<input type="checkbox"/> 4021141 Stealth ADI AL Active 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 945
	4" depth	<input type="checkbox"/> 4021142 Stealth ADI AL Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,045
	6" depth	<input type="checkbox"/> 4021143 Stealth ADI AL Extra Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,145

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> 4021110 Matrx Elite 300 lbs weight capacity. 12"-20" width, 10"-20" height \$ 955
		<input type="checkbox"/> 4021111 Matrx Elite Heavy Duty 500 lbs weight capacity. 21" width, 16"/20" height \$ 1,195
		<input type="checkbox"/> 4021112 22"-24" width, 16"/20" height \$ 1,295
		<input type="checkbox"/> 4021114 25"-30" width, 16"/20" height \$ 1,995
	3.5" to 4" depth	<input type="checkbox"/> 4021120 ROHO Agility Min 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 900
		<input type="checkbox"/> 4021100 nxt Xtend 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021105 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021122 ROHO Agility Mid 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 930
	6" to 7" depth	<input type="checkbox"/> 4021101 nxt Xtend Deep 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021106 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021116 Matrx Elite Deep 300 lbs weight capacity. 14"-20" width, 12"-20" height \$ 1,145
		<input type="checkbox"/> 4021117 22"-24" width, 12"-20" height \$ 1,295
<input type="checkbox"/> 4021124 ROHO Agility Max 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 1,085		

4. Specify the size →	Backrest Size Sizes based on availability from the backrest manufacturer. Specify Size (Width" x Height") _____
-----------------------	---

Pediatric Seating System

4030400 **Diego Pediatric Seating System**

Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10	12	14	16
Width	(25.40)	(30.48)	(35.56)	(40.64)
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓

A Head System Head Pad & Hardware

B Contoured Back System Base, Cushion & Cover

FlipLock Configurable as a between canes or full width set up.

CoolCore Covers Our Diego System uses the CoolCore cover technology which adapt all-year-round for a more comfortable experience.

Laterals Modular Laterals either Fixed or Swing Away

Contoured Seat System Base, Cushion, Cover & Hardware

Back Hardware Cover Easy to remove back cover for quick on-demand adjustments and a more aesthetic look.

Width Growth Adapted from our TFB System, the back wings allow for a pediatric growth adjustment.

Multi-Axis Hardware The HMO501, included as part of a package for the Diego Head System, allows growth adjustability and on-demand configurations.

WC-20[®] Certified Some requirements need to be met for safety purposes. Please verify!

Headrest

Stealth Comfort Plus Headrest	Stealth Hardware		
	<input type="checkbox"/> 4024110	Multi-Axis HMO475	\$ 225
	<input type="checkbox"/> 4024111	Link TWB480	\$ 350
	Comfort Plus Pad		\$ 275
	<input type="checkbox"/> 4024115	Large (10" x 6")	<input type="checkbox"/> 4024116 X-Large (14" x 6")
nxt MultiFit Head Support	<input type="checkbox"/> 4024100 nxt Hardware		\$ 285
	Multi-Adjustable Pad		\$ 275
	<input type="checkbox"/> 4024105	Contoured (9" x 5")	<input type="checkbox"/> 4024107 X-Large (14" x 5")
	<input type="checkbox"/> 4024106	Large (10" x 5")	

Positioning Belts & Straps

Bodypoint	Hip Belts		
	<input type="checkbox"/> 4029010	Center-Pull Auto Style	\$ 115
	<input type="checkbox"/> 4029015	Center-Pull Auto Style, Padded	\$ 215
	<input type="checkbox"/> 4029020	4-Point, Padded	\$ 235
	<input type="checkbox"/> 4029030	Chest Strap	\$ 165

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.

PDG proudly supplies on-chair seating and positioning products from these manufacturers:

