

Customer and Order Information

Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 2301000	Stellar IMPACT 45° Manual Tilt Wheelchair Extended delivery lead times may apply.	\$ 3,995
<input checked="" type="checkbox"/> 2301201	Fatigue-resistant Reinforced Steel Frame Package Ultra-durable steel reinforcement of frame and other critical components.	\$ 1,095

Frame Options

<input checked="" type="checkbox"/> 2302000	Standard Base and Seat Frame 250 lbs Weight Capacity.	
<input type="checkbox"/> 2302200	Extended Base with Rear Floor Jacks Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Locks, HD Anti-Tipper, IV Pole, O2 Holder.	\$ 600
<input type="checkbox"/> 2302300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base.	

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths		
<input type="checkbox"/> 2303014 14" Seat Width	<input type="checkbox"/> 2303015 15" Seat Width	<input type="checkbox"/> 2303016 16" Seat Width
<input type="checkbox"/> 2303017 17" Seat Width	<input type="checkbox"/> 2303018 18" Seat Width	<input type="checkbox"/> 2303019 19" Seat Width
<input type="checkbox"/> 2303020 20" Seat Width		
<input type="checkbox"/> 2303021 21" Seat Width	<input type="checkbox"/> 2303022 22" Seat Width	\$ 215

Seat Depth

Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 2304017 17" Seat Depth	<input type="checkbox"/> 2304018 18" Seat Depth	<input type="checkbox"/> 2304019 19" Seat Depth
<input type="checkbox"/> 2304020 20" Seat Depth		
Custom Seat Depths		
<input type="checkbox"/> 2304021 21" Seat Depth	<input type="checkbox"/> 2304022 22" Seat Depth	\$ 325

Frame Colour

PDG Premium Colours		
<input type="checkbox"/> 2305000 Tracer Red	<input type="checkbox"/> 2305001 Gloss Black	<input type="checkbox"/> 2305002 Silver
<input type="checkbox"/> 2305003 Black Shimmer	<input type="checkbox"/> 2305004 Orange Fusion	<input type="checkbox"/> 2305005 Karl's Imagination
<input type="checkbox"/> 2305006 Blue Shimmer	<input type="checkbox"/> 2305007 Gloss White	<input type="checkbox"/> 2305008 Dark Titanium
<input type="checkbox"/> 2305009 Geyser Green	<input type="checkbox"/> 2305010 Lollypop Blue	<input type="checkbox"/> 2305011 Lollypop Red
<input type="checkbox"/> 2305012 Bubble Gum Pink	<input type="checkbox"/> 2305013 Berry Purple	<input type="checkbox"/> 2305014 Safety Yellow
<input type="checkbox"/> 2305200 Custom Colour Contact PDG for specifications.		\$ 300

Arm Type

<input type="checkbox"/> 2306000 No Armrest		
<input type="checkbox"/> 2306100 Single Post Fixed Height (12") Full Length Armpad only.		
Single Post Adjustable Height		\$ 255
<input type="checkbox"/> 2306111 8"-11"	<input type="checkbox"/> 2306114 10"-14"	

See next page for more Arm Type options >

<input type="checkbox"/> 2306200	Dual Post Fixed Height (12")	Full Length Armpad only.			
	Dual Post Adjustable Height				\$ 255
<input type="checkbox"/> 2306211	8"-11"	<input type="checkbox"/> 2306214	10"-14"	<input type="checkbox"/> 2306216	12"-16"
	Locking Flip-up	Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.			\$ 255
<input type="checkbox"/> 2306310	Set at 10"	<input type="checkbox"/> 2306311	Set at 11"	<input type="checkbox"/> 2306312	Set at 12"
<input type="checkbox"/> 2306313	Set at 13"	<input type="checkbox"/> 2306314	Set at 14"		

Armrests					
Regular Armpads					
<input type="checkbox"/> 2307100	Desk Length - 2" x 10"	<input type="checkbox"/> 2307101	Full Length - 2" x 14"		
<input type="checkbox"/> 2307102	Full Length Premium Grade - 2" x 14"				\$ 35
Wide Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2307110	Desk Length - 3" x 10"	<input type="checkbox"/> 2307111	Full Length - 3" x 12"		\$ 115
Gel Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2307120	Full Length - 2" x 12"	<input type="checkbox"/> 2307121	Full Length - 3.5" x 12"		\$ 125
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 2307130	Right (4.5" x 18")	<input type="checkbox"/> 2307131	Left (4.5" x 18")		\$ 125
<input type="checkbox"/> 2307220	Armrest Wideners for Dual Post Armrests	Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.			\$ 375
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.					
<input type="checkbox"/> 2307310	For Single Post Armrests	Available for 10"-14" Height range only.			\$ 110
<input type="checkbox"/> 2307320	For Dual Post Armrests				\$ 165
Armrest Extensions Not compatible with Flip-up or Single Post Armrests.					
<input type="checkbox"/> 2307400	3" Extension - Forward	<input type="checkbox"/> 2307401	3" Extension - Rearward		\$ 60

Seat Pan	
<input type="checkbox"/> 2308000	No Seat Pan
<input type="checkbox"/> 2308001	Solid Aluminum Seat Pan
	\$ 120

Cushion	
<input type="checkbox"/> 2309000	No Cushion
next BioFit	
<input type="checkbox"/> 2309101	14"-20" width, 16"-20" depth
<input type="checkbox"/> 2309102	22" width, 16"-20" depth
	\$ 595
	\$ 625
Cushion Size Available in even number sizes only.	
2309400	Specify Size (Width" x Depth") _____

Seat Height																																																																																		
Select desired Seat-to-Floor Height																																																																																		
STF Height is adjustable within the range of the selected Rear Wheel size.																																																																																		
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="8">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td></td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> <td>21"</td> </tr> <tr> <td>20"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td></td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> <td></td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> <td>21"</td> </tr> </tbody> </table>	RW	6" Casters								12"	14"	15"	16"	17"	18"	19"			16"			16"	17"	18"	19"	20"	21"	20"	14"	15"	16"	17"	18"	19"			22"		15"	16"	17"	18"	19"	20"		24"			16"	17"	18"	19"	20"	21"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="8">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> <td>21"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> <td>21"</td> </tr> </tbody> </table>	RW	8" Casters								16"			16"	17"	18"	19"	20"	21"	24"			16"	17"	18"	19"	20"	21"
RW	6" Casters																																																																																	
12"	14"	15"	16"	17"	18"	19"																																																																												
16"			16"	17"	18"	19"	20"	21"																																																																										
20"	14"	15"	16"	17"	18"	19"																																																																												
22"		15"	16"	17"	18"	19"	20"																																																																											
24"			16"	17"	18"	19"	20"	21"																																																																										
RW	8" Casters																																																																																	
16"			16"	17"	18"	19"	20"	21"																																																																										
24"			16"	17"	18"	19"	20"	21"																																																																										
	Wheel size _____																																																																																	
	Caster size _____																																																																																	
	STF height _____																																																																																	

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2311061	6" x 1.25" Caster size			
<input type="checkbox"/>	2311063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/>	2311081	8" x 1.25" Caster size			
<input type="checkbox"/>	2311084	8" x 2.00" Caster size			\$ 135

Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2311161	6" x 1.25" Caster size			\$ 110
<input type="checkbox"/>	2311181	8" x 1.25" Caster size			\$ 120
<input type="checkbox"/>	2311184	8" x 2.00" Caster size			\$ 150

<input type="checkbox"/>	2311260	6" Soft Roll			\$ 95
--------------------------	---------	---------------------	--	--	-------

<input type="checkbox"/>	2311400	Ultra Sport Fork	Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.		\$ 395
--------------------------	---------	-------------------------	---	--	--------

Rear Wheels

Composite

<input type="checkbox"/>	2312012	12" Wheel size	<input type="checkbox"/>	2312016	16" Wheel size	<input type="checkbox"/>	2312020	20" Wheel size
<input type="checkbox"/>	2312022	22" Wheel size	<input type="checkbox"/>	2312024	24" Wheel size			

Spoke

<input type="checkbox"/>	2312120	20" Wheel size	<input type="checkbox"/>	2312122	22" Wheel size	<input type="checkbox"/>	2312124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------

<input type="checkbox"/>	2312224	24" Heavy Duty Spoke							\$ 350
--------------------------	---------	-----------------------------	--	--	--	--	--	--	--------

Axle Type

<input type="checkbox"/>	2313000	Quick Release Rear Axles			
--------------------------	---------	---------------------------------	--	--	--

<input type="checkbox"/>	2313010	Fixed Rear Axles			
--------------------------	---------	-------------------------	--	--	--

Rear Tires

<input type="checkbox"/>	2314000	Urethane			
--------------------------	---------	-----------------	--	--	--

<input type="checkbox"/>	2314010	Pneumatic	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.		\$ 60
--------------------------	---------	------------------	--	--	-------

<input type="checkbox"/>	2314020	Pneumatic with Flat-Free Inserts	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.		\$ 95
--------------------------	---------	---	--	--	-------

Hand Rims

<input type="checkbox"/>	2315000	No Hand Rims			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	2315010	Standard	20", 22", and 24" Rear Wheels only.		
--------------------------	---------	-----------------	-------------------------------------	--	--

<input type="checkbox"/>	2315020	Plastic Coated	20", 22", and 24" Rear Wheels only.		\$ 135
--------------------------	---------	-----------------------	-------------------------------------	--	--------

Wheel Locks

<input type="checkbox"/>	2316010	Push to Lock			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	2316020	Pull to Lock			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	2316030	Attendant Wheel Lock (Push to Lock)	Not compatible with TTDS.		\$ 85
--------------------------	---------	--	---------------------------	--	-------

<input type="checkbox"/>	2316040	User & Attendant Wheel Lock	Available only on 20", 22", or 24" Rear Wheels.		\$ 125
--------------------------	---------	--	---	--	--------

<input type="checkbox"/>	2316100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.		\$ 60
--------------------------	---------	------------------------------	---	--	-------

<input type="checkbox"/>	2316200	Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts.		\$ 350
--------------------------	---------	---------------------	---	--	--------

Anti-Tippers

<input checked="" type="checkbox"/>	2317010	Rear Anti-Tippers Required Safety Option.	\$ 60
<input type="checkbox"/>	2317020	Front Anti-Tippers	\$ 300
<input type="checkbox"/>	2317030	Heavy Duty Rear Anti-Tipper Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	\$ 600

Back Post Style

<input type="checkbox"/>	2318010	Push Bar Back	
<input type="checkbox"/>	2318020	Cane Handle Back Not compatible with Headrest Mounting Plate.	
Back Post Height Seat Pan to Handle			
<input type="checkbox"/>	2318425	25" Height	<input type="checkbox"/> 2318431 31" Height

Back Options

Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.			\$ 150
<input type="checkbox"/>	2319010	Button on Outside	<input type="checkbox"/> 2319011 Button on Inside
Seat to Back Angle			
<input type="checkbox"/>	2319090	Set at 90°	<input type="checkbox"/> 2319098 Set at 98° <input type="checkbox"/> 2319106 Set at 106°
<input type="checkbox"/>	2319114	Set at 114°	<input type="checkbox"/> 2319122 Set at 122°
<input type="checkbox"/>	2319210	10" Backrest Extension Not compatible with Angle/Height Adjustable Push Bar Handle.	\$ 135
<input type="checkbox"/>	2319220	Fold Down Back Not compatible with Reclining Back.	
<input checked="" type="checkbox"/>	2319230	Spring Back Suspension Not compatible with Reclining Back.	
<input type="checkbox"/>	2319240	Dynamic Rocker Back Assembly Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.	\$ 2,995

Back Upholstery

<input type="checkbox"/>	2320000	No Back Upholstery	
Standard Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	2320021	Set at 21"	<input type="checkbox"/> 2320022 Set at 22" Not available with 25" Back.
Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	2320118	Set at 18"	<input type="checkbox"/> 2320119 Set at 19" <input type="checkbox"/> 2320120 Set at 20"
Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.			\$ 155
<input type="checkbox"/>	2320224	Set at 24"	<input type="checkbox"/> 2320225 Set at 25" <input type="checkbox"/> 2320226 Set at 26"
Adjustable Tension Back height measured from Seat Pan to top of Upholstery.			\$ 295
<input type="checkbox"/>	2320321	Set at 21"	<input type="checkbox"/> 2320322 Set at 22" <input type="checkbox"/> 2320323 Set at 23"

Backrest

<input type="checkbox"/>	2321000	No Backrest	
nxt Optima HA Active Contour			
<input type="checkbox"/>	2321100	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	2321105	22" width, 16"-20" height	\$ 720
nxt Optima HA Plus Contour			
<input type="checkbox"/>	2321101	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	2321106	22" width, 16"-20" height	\$ 720
Backrest Size Available in even number sizes only.			
2321400	Specify Size (Width" x Height") _____		



Tilt/Recline Options

Attendant Tilt only				
<input type="checkbox"/>	2323011	Trigger on Right Back Post	<input type="checkbox"/>	2323012 Trigger on Left Back Post
User Tilt only		Full Length Arm Pads only.		\$ 100
<input type="checkbox"/>	2323021	Trigger on Right Armrest	<input type="checkbox"/>	2323022 Trigger on Left Armrest
Attendant & User Tilt		Full Length Arm Pads only.		\$ 195
<input type="checkbox"/>	2323031	Trigger on Right Back Post	<input type="checkbox"/>	2323032 Trigger on Left Back Post
<input type="checkbox"/>	2323041	Trigger on Right Armrest	<input type="checkbox"/>	2323042 Trigger on Left Armrest
<input type="checkbox"/>	2323110	Trigger Extension		\$ 50
Tilt Limit				
<input type="checkbox"/>	2323220	Set at 20°	<input type="checkbox"/>	2323225 Set at 25°
<input type="checkbox"/>	2323235	Set at 35°	<input type="checkbox"/>	2323240 Set at 40°
<input type="checkbox"/>	2323230	Set at 30°	<input type="checkbox"/>	2323245 Set at 45°
<input type="checkbox"/>	2323410	Reclining Backrest (Manual)		\$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.				

Headrest

<input type="checkbox"/>	2324000	No Headrest		
<input type="checkbox"/>	2324010	Headrest Mounting Plate		\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.				
<input type="checkbox"/>	2324100	next MultiFit Head Support Hardware		\$ 235
		next MultiFit Head Support Pad		\$ 245
<input type="checkbox"/>	2324105	Regular (9" x 5")	<input type="checkbox"/>	2324106 Large (10" x 5")
			<input type="checkbox"/>	2324107 X-Large (14" x 5")

Front Rigging

<input type="checkbox"/>	2325000	No Front Rigging		
		PDG Easy-Out™ Front Rigging	Effortless insertion and removal. No Hanger Bracket remains on frame.	
<input type="checkbox"/>	2325410	Fixed 70°		Available in March \$ 95
<input type="checkbox"/>	2325420	Swing In/Away Universal Angle Adjustable (0°-110°)	Set at (use increments of 10°): _____	Available in April \$ 195
<input type="checkbox"/>	2325430	Swing In/Away Elevating Legrests	High knee pivot. Includes reinforced premium Calf Pad mount.	Available in April \$ 495
		Standard Front Rigging	Legacy Design. Not compatible with Easy-Out™.	
<input type="checkbox"/>	2325110	Plug In 70°		
<input type="checkbox"/>	2325130	Swing Away 60°	<input type="checkbox"/>	2325140 Swing Away 70°
<input type="checkbox"/>	2325200	Swing Away Elevating Legrest		\$ 375
<input type="checkbox"/>	2325300	1" Offset Hangers for Swing Away		\$ 75

Foot Plates

<input type="checkbox"/>	2326010	Composite		
<input type="checkbox"/>	2326020	Aluminum Angle Adjustable		\$ 155
<input type="checkbox"/>	2326030	One Piece Fixed	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	2326040	One Piece Flip-up	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	2326050	Interlocking Angle Adjustable	Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	2327010	Swing Away 60° & 70° Calf Pad	Not compatible with 60° & 70° Gel Padding.	\$ 135
<input type="checkbox"/>	2327020	Swing Away 60° & 70° Gel Padding	Not compatible with 60° & 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	2327030	Swing Away ELR Gel Knee Pad		\$ 100
<input type="checkbox"/>	2327040	Adjustable Length Heel Loops		\$ 50

See next page for more Footrest/Legrest Accessories options >

<input type="checkbox"/>	2327050	Bodypoint - Padded Calf Strap	Sizes greater than 22" are not padded.	\$	105
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.					
<input type="checkbox"/>	2327060	For Plug In Front Rigging		\$	190
<input type="checkbox"/>	2327065	For Swing Away 70° Front Rigging	Not available on 13" STF Height.	\$	380

Extension Tubes

Regular Extension Tubes

<input type="checkbox"/>	2328010	Short	13"-17" for Easy-Out SA Angle Adjustable	10"-13.5" for Easy Out SA ELR	13"-17" for Standard Front Rigging	
<input type="checkbox"/>	2328020	Long	16"-22.5" for Easy-Out SA Angle Adjustable	12.5"-19.5" for Easy-Out SA ELR	15"-19" for Standard Front Rigging	\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only.						\$ 100
<input type="checkbox"/>	2328110	Standard (13"-17")		<input type="checkbox"/>	2328120	Long (15"-19")

Positioning Straps

Bodypoint - Hip Belts

<input type="checkbox"/>	2329010	Center-Pull Auto Style		\$	105
<input type="checkbox"/>	2329015	Center-Pull Auto Style, Padded		\$	155
<input type="checkbox"/>	2329020	4-Point, Padded		\$	175
<input type="checkbox"/>	2329030	Bodypoint - Chest Strap		\$	85

Miscellaneous Options

<input type="checkbox"/>	2330410	Strut Extension Damper	Optional Damper for superior energy dissipation.	\$	255
<input type="checkbox"/>	2330050	Protective Lower Padding		\$	425
<input type="checkbox"/>	2330080	Gel Hip Guides		\$	325
<input type="checkbox"/>	2330110	IV Pole	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	295
<input type="checkbox"/>	2330120	O2 Holder	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	395
<input type="checkbox"/>	2330170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$	10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar-impact/ to learn more about the **Stellar IMPACT 45° Manual Tilt Wheelchair**

