

## Customer and Order Information

Purchase Order *	Tag	v.0301
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <span style="float: right;">*required</span>

## Frame Type

<input checked="" type="checkbox"/> 2301000	<b>Stellar IMPACT 45° Manual Tilt Wheelchair</b> Extended delivery lead times may apply.	\$ 3,995
<input checked="" type="checkbox"/> 2301201	<b>Fatigue-resistant Reinforced Steel Frame Package</b> Ultra-durable steel reinforcement of frame and other critical components.	\$ 1,095
<b>i</b> Consider the <b>Stellar</b> or the <b>Stellar HD</b> for users over 250 lbs. Custom reinforcement available. See separate order forms.		

## Frame Options

<input checked="" type="checkbox"/> 2302000	<b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.	
<input type="checkbox"/> 2302200	<b>Extended Base with Rear Floor Jacks</b> Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.	\$ 900
<input type="checkbox"/> 2302300	<b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory-installed transport brackets. See Owners Manual for complete details. Not compatible with Extended Base.	\$ 250

## Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

<b>Standard Seat Widths</b>					
<input type="checkbox"/> 2303014	14" Seat Width	<input type="checkbox"/> 2303015	15" Seat Width	<input type="checkbox"/> 2303016	16" Seat Width
<input type="checkbox"/> 2303017	17" Seat Width	<input type="checkbox"/> 2303018	18" Seat Width	<input type="checkbox"/> 2303019	19" Seat Width
<input type="checkbox"/> 2303020	20" Seat Width	<input type="checkbox"/> 2303021	21" Seat Width	<input type="checkbox"/> 2303022	22" Seat Width

## Seat Depth

<b>Standard Seat Depths (Adjustable)</b>					
<input type="checkbox"/> 2304017	17" Seat Depth	<input type="checkbox"/> 2304018	18" Seat Depth	<input type="checkbox"/> 2304019	19" Seat Depth
<input type="checkbox"/> 2304020	20" Seat Depth				
<b>Custom Seat Depths</b>					
<input type="checkbox"/> 2304021	21" Seat Depth	<input type="checkbox"/> 2304022	22" Seat Depth	\$ 350	

## Frame Color

<b>PDG Premium Colors</b>					
<input type="checkbox"/> 2305000	Tracer Red	<input type="checkbox"/> 2305001	Gloss Black	<input type="checkbox"/> 2305002	Silver
<input type="checkbox"/> 2305003	Black Shimmer	<input type="checkbox"/> 2305004	Orange Fusion	<input type="checkbox"/> 2305005	Karl's Imagination
<input type="checkbox"/> 2305006	Blue Shimmer	<input type="checkbox"/> 2305007	Gloss White	<input type="checkbox"/> 2305008	Dark Titanium
<input type="checkbox"/> 2305009	Geyser Green	<input type="checkbox"/> 2305010	Lollypop Blue	<input type="checkbox"/> 2305011	Lollypop Red
<input type="checkbox"/> 2305012	Bubble Gum Pink	<input type="checkbox"/> 2305013	Berry Purple	<input type="checkbox"/> 2305014	Safety Yellow
<input type="checkbox"/> 2305200	<b>Custom Color</b> Contact PDG for specifications.				\$ 300

## Arm Type

<input type="checkbox"/> 2306000	<b>No Armrest</b>	
<input type="checkbox"/> 2306100	<b>Single Post Fixed Height (12")</b> Full Length Armpad only.	
<b>Single Post Adjustable Height</b>		\$ 275
<input type="checkbox"/> 2306111	8"-11"	<input type="checkbox"/> 2306114 10"-14"

See next page for more options >

<input type="checkbox"/>	2306200	<b>Dual Post Fixed Height (12")</b>	Full Length Armpad only.				
		<b>Dual Post Adjustable Height</b>				\$	275
<input type="checkbox"/>	2306211	8"-11"		<input type="checkbox"/>	2306214	10"-14"	
				<input type="checkbox"/>	2306216	12"-16"	
		<b>Locking Flip-up</b>	Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.			\$	295
<input type="checkbox"/>	2306310	Set at 10"		<input type="checkbox"/>	2306311	Set at 11"	
<input type="checkbox"/>	2306313	Set at 13"		<input type="checkbox"/>	2306314	Set at 14"	

Armrests							
		<b>Regular Armpads</b>					
<input type="checkbox"/>	2307100	Desk Length - 2" x 10"		<input type="checkbox"/>	2307101	Full Length - 2" x 14"	
		<b>Wide Armpads</b>	May interfere with wheels in some tilt configurations.			\$	200
<input type="checkbox"/>	2307110	Desk Length - 3" x 10"		<input type="checkbox"/>	2307111	Full Length - 3" x 12"	
		<b>Gel Armpads</b>	May interfere with wheels in some tilt configurations.			\$	135
<input type="checkbox"/>	2307120	Full Length - 2" x 12"		<input type="checkbox"/>	2307121	Full Length - 3.5" x 12"	
		<b>Contour Arm Trough with Hand Pad</b>	May interfere with wheels in some tilt configurations.			\$	125
<input type="checkbox"/>	2307130	Right (4.5" x 18")		<input type="checkbox"/>	2307131	Left (4.5" x 18")	
<input type="checkbox"/>	2307220	<b>Armrest Wideners for Dual Post Armrests</b>	Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.			\$	375
		<b>Clothing Guards</b>	Not compatible with Flip-up Armrests, Armrest Wideners.				
<input type="checkbox"/>	2307310	For Single Post Armrests	Available for 10"-14" Height range only.			\$	110
<input type="checkbox"/>	2307320	For Dual Post Armrests				\$	165
		<b>Armrest Extensions</b>	Not compatible with Flip-up or Single Post Armrests.			\$	75
<input type="checkbox"/>	2307400	3" Extension - Forward		<input type="checkbox"/>	2307401	3" Extension - Rearward	

Seat Pan							
<input type="checkbox"/>	2308000	<b>No Seat Pan</b>					
<input type="checkbox"/>	2308001	<b>Solid Aluminum Seat Pan</b>				\$	120

**Cushion**

Refer to the enclosed Seating & Positioning order form for PDG's on-chair Cushion options. Standard PDG discounts apply.

Seat Height								
		<b>Select desired Seat-to-Floor Height</b>	STF Height is adjustable within the range of the selected Rear Wheel size.					
	<b>RW</b>	<b>6" Casters</b>						
	12"	14"	15"	16"	17"	18"	19"	
	16"			16"	17"	18"	19"	20" 21"
	20"	14"	15"	16"	17"	18"	19"	
	22"		15"	16"	17"	18"	19"	20"
	24"			16"	17"	18"	19"	20" 21"
	<b>RW</b>	<b>8" Casters</b>						
	16"			16"	17"	18"	19"	20" 21"
	24"			16"	17"	18"	19"	20" 21"
								Wheel size _____
								Caster size _____
								STF height _____

Casters								
		<b>Urethane</b>	8" Casters are compatible only with 16" and 24" Rear Wheels.					
<input type="checkbox"/>	2311061	6" x 1.25" Caster size						
<input type="checkbox"/>	2311063	6" x 1.50" Caster size					\$ 85	
<input type="checkbox"/>	2311081	8" x 1.25" Caster size						
<input type="checkbox"/>	2311084	8" x 2.00" Caster size					\$ 135	
		<b>Pneumatic</b>	8" Casters are compatible only with 16" and 24" Rear Wheels.					
<input type="checkbox"/>	2311161	6" x 1.25" Caster size					\$ 110	
<input type="checkbox"/>	2311181	8" x 1.25" Caster size					\$ 120	
<input type="checkbox"/>	2311184	8" x 2.00" Caster size					\$ 150	

See next page for more options >

<input type="checkbox"/>	2311260	<b>6" Soft Roll</b>		\$	95
<input type="checkbox"/>	2311400	<b>Ultra Sport Fork</b>	Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	\$	395

## Rear Wheels

<b>Composite</b>							
<input type="checkbox"/>	2312012	12" Wheel size	<input type="checkbox"/>	2312016	16" Wheel size		
<input type="checkbox"/>	2312022	22" Wheel size	<input type="checkbox"/>	2312024	24" Wheel size		
<input type="checkbox"/>	2312020	20" Wheel size					
<b>Spoke</b>					\$	75	
<input type="checkbox"/>	2312120	20" Wheel size	<input type="checkbox"/>	2312122	22" Wheel size		
<input type="checkbox"/>	2312124	24" Wheel size					
<input type="checkbox"/>	2312224	<b>24" Heavy Duty Spoke</b>				\$	350

## Axle Type

<input type="checkbox"/>	2313000	<b>Quick Release Rear Axles</b>			
<input type="checkbox"/>	2313010	<b>Fixed Rear Axles</b>			

## Rear Tires

<input type="checkbox"/>	2314000	<b>Urethane</b>			
<input type="checkbox"/>	2314010	<b>Pneumatic</b>	May result in a finished Seat Height 1/2" greater than stated.	\$	60
<input type="checkbox"/>	2314020	<b>Pneumatic with Flat-Free Inserts</b>	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$	95

## Hand Rims

<input type="checkbox"/>	2315000	<b>No Hand Rims</b>			
<input type="checkbox"/>	2315010	<b>Standard</b>	20", 22", and 24" Rear Wheels only.		
<input type="checkbox"/>	2315020	<b>Plastic Coated</b>	20", 22", and 24" Rear Wheels only.	\$	135

## Wheel Locks

<input type="checkbox"/>	2316010	<b>Push to Lock</b>			
<input type="checkbox"/>	2316020	<b>Pull to Lock</b>			
<input type="checkbox"/>	2316030	<b>Attendant Wheel Lock (Push to Lock)</b>		\$	85
<input type="checkbox"/>	2316040	<b>User &amp; Attendant Wheel Lock</b>	Available only on 20", 22", or 24" Rear Wheels.	\$	125
<input type="checkbox"/>	2316100	<b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.	\$	60
<input type="checkbox"/>	2316200	<b>Hill Holders</b>	Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$	450

## Anti-Tippers

<input checked="" type="checkbox"/>	2317010	<b>Rear Anti-Tippers</b>	Required Safety Option.	\$	60
<input type="checkbox"/>	2317020	<b>Front Anti-Tippers</b>		\$	795
<input type="checkbox"/>	2317030	<b>Heavy Duty Rear Anti-Tipper</b>	Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	\$	895

## Back Post Style

<input type="checkbox"/>	2318010	<b>Push Bar Back</b>			
<input type="checkbox"/>	2318020	<b>Cane Handle Back</b>	Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.		
<b>Back Post Height</b> Seat Pan to Handle					
<input type="checkbox"/>	2318425	25" Height	<input type="checkbox"/>	2318431	31" Height

## Back Options

<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 2319010	Button on Outside	<input type="checkbox"/> 2319011 Button on Inside
<b>Seat to Back Angle</b>		
<input type="checkbox"/> 2319090	Set at 90°	<input type="checkbox"/> 2319098 Set at 98°
<input type="checkbox"/> 2319114	Set at 114°	<input type="checkbox"/> 2319106 Set at 106°
<input type="checkbox"/> 2319122	Set at 122°	
<input type="checkbox"/> 2319210	<b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Push Bar Handle.	\$ 150
<input type="checkbox"/> 2319220	<b>Fold Down Back</b> Not compatible with Reclining Back.	
<input checked="" type="checkbox"/> 2319230	<b>Spring Back Suspension</b> Not compatible with Reclining Back.	
<input type="checkbox"/> 2319240	<b>Dynamic Rocker Back Assembly</b> Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.	\$ 2,995

## Back Upholstery

<input type="checkbox"/> 2320000	<b>No Back Upholstery</b>	
<b>Standard</b> Back height measured from Seat Pan to top of Upholstery.		\$ 175
<input type="checkbox"/> 2320021	Set at 21"	<input type="checkbox"/> 2320022 Set at 22"
<input type="checkbox"/> 2320023	Set at 23"	Not available with 25" Back.
<b>Short</b> 25" Back only. Back height measured from Seat Pan to top of Upholstery.		\$ 175
<input type="checkbox"/> 2320118	Set at 18"	<input type="checkbox"/> 2320119 Set at 19"
<input type="checkbox"/> 2320120	Set at 20"	
<b>Tall</b> 31" Back only. Back height measured from Seat Pan to top of Upholstery.		\$ 175
<input type="checkbox"/> 2320224	Set at 24"	<input type="checkbox"/> 2320225 Set at 25"
<input type="checkbox"/> 2320226	Set at 26"	
<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery.		\$ 295
<input type="checkbox"/> 2320321	Set at 21"	<input type="checkbox"/> 2320322 Set at 22"
<input type="checkbox"/> 2320323	Set at 23"	

## Positioning Backrest

**Refer to the enclosed Seating & Positioning order form for PDG's on-chair Backrest options. Standard PDG discounts apply.**

## Tilt/Recline Options

<b>Attendant Tilt only</b>		
<input type="checkbox"/> 2323011	Trigger on Right Back Post	<input type="checkbox"/> 2323012 Trigger on Left Back Post
<b>User Tilt only</b> Not available on widths greater than 22". Full Length Arm Pads only.		\$ 100
<input type="checkbox"/> 2323021	Trigger on Right Armrest	<input type="checkbox"/> 2323022 Trigger on Left Armrest
<b>Attendant &amp; User Tilt</b> Not available on widths greater than 22". Full Length Arm Pads only.		\$ 195
<input type="checkbox"/> 2323031	Trigger on Right Back Post	<input type="checkbox"/> 2323032 Trigger on Left Back Post
<input type="checkbox"/> 2323041	Trigger on Right Armrest	<input type="checkbox"/> 2323042 Trigger on Left Armrest
<input type="checkbox"/> 2323110	<b>Trigger Extension</b>	\$ 50
<b>Tilt Limit</b>		
<input type="checkbox"/> 2323220	Set at 20°	<input type="checkbox"/> 2323225 Set at 25°
<input type="checkbox"/> 2323230	Set at 30°	
<input type="checkbox"/> 2323235	Set at 35°	<input type="checkbox"/> 2323240 Set at 40°
<input type="checkbox"/> 2323245	Set at 45°	
<input type="checkbox"/> 2323310	<b>Reclining Backrest (Manual)</b> 30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.	\$ 895

## Headrest

**Refer to the enclosed Seating & Positioning order form for PDG's on-chair Headrest options. Standard PDG discounts apply.**

<input type="checkbox"/> 2324010	<b>Headrest Mounting Plate</b> Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 95
----------------------------------	---	-------

## Front Rigging

<input type="checkbox"/>	2325000	<b>No Front Rigging</b>	<input type="checkbox"/>	2325001	<b>Sector Blocks only</b>		
<b>PDG Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>							
<input type="checkbox"/>	2325410	Fixed 70°					\$ 95
<input type="checkbox"/>	2325420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)				\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment.							
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.							
<input type="checkbox"/>	2325130	Swing Away 60°	<input type="checkbox"/>	2325140	Swing Away 70°		
<input type="checkbox"/>	2325200	Swing Away Elevating Legrests	Includes Calf Pads.				\$ 375
<input type="checkbox"/>	2325300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	2325110	<b>Standard Plug In 70°</b>	Legacy Design. Not compatible with Easy-Out.				

## Foot Plates

<input type="checkbox"/>	2326010	<b>Composite</b>					
<input type="checkbox"/>	2326020	<b>Aluminum Angle Adjustable</b>					\$ 155
<input type="checkbox"/>	2326030	<b>One Piece Fixed</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	2326040	<b>One Piece Flip-up</b>	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	2326050	<b>Interlocking Angle Adjustable</b>	Seat Widths of 14"-21" only.				\$ 195

## Footrest/Legrest Accessories

<b>Calf Pad</b> Not compatible with Gel Padding.							
<input type="checkbox"/>	2327012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.				\$ 200
<input type="checkbox"/>	2327010	For Swing Away Front Rigging					\$ 135
<input type="checkbox"/>	2327020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.				\$ 125
<b>Gel Knee Pad</b>							
<input type="checkbox"/>	2327032	For Easy-Out™ Front Rigging					\$ 200
<input type="checkbox"/>	2327030	For Swing Away ELR Front Rigging					\$ 100
<input type="checkbox"/>	2327040	<b>Adjustable Length Heel Loops</b>					\$ 50
<input type="checkbox"/>	2327050	<b>Padded Calf Strap</b>	Sizes greater than 22" are not padded.				\$ 195
<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.							
<input type="checkbox"/>	2327060	For Standard Plug In Front Rigging					\$ 345
<input type="checkbox"/>	2327065	For Swing Away 70° Front Rigging					\$ 475

## Extension Tubes

<b>Regular Extension Tubes</b>							
<input type="checkbox"/>	2328010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	2328020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging			\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.							
<input type="checkbox"/>	2328110	Standard (13"-17")	<input type="checkbox"/>	2328120	Long (15"-19")		\$ 100

## Positioning Straps

**Refer to the enclosed Seating & Positioning order form for PDG's on-chair Belts/Straps options. Standard PDG discounts apply.**

### Miscellaneous Options

<input type="checkbox"/>	2330410	<b>Strut Extension Damper</b>	Optional Damper for superior energy dissipation.	\$	255
<input type="checkbox"/>	2330050	<b>Protective Lower Padding</b>		\$	825
<input type="checkbox"/>	2330080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".	\$	375
<input type="checkbox"/>	2330110	<b>IV Pole</b>	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	295
<input type="checkbox"/>	2330120	<b>O2 Holder</b>	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	395
<input type="checkbox"/>	2330170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.	\$	10

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/stellar-impact/](http://pdgmobility.com/products/stellar-impact/) to learn more about the **Stellar IMPACT 45° Manual Tilt Wheelchair**



## Cushion

<b>1. Select the type</b> ↓	<b>2. Select the model</b> ↓
--------------------------------	---------------------------------

<b>General Use</b>	<b>Stealth Simplicity</b> 350 lbs weight capacity. <input type="checkbox"/> 4009141 16"-20" width, 16"-20" depth \$ 310 <input type="checkbox"/> 4009142 22"-24" width, 18"-20" depth \$ 370
--------------------	--

<b>Positioning</b>	<b>ROHO AirLITE</b> 300 lbs weight capacity. <input type="checkbox"/> 4009120 14"-20" width, 14"-18" depth \$ 375
--------------------	--

<b>Skin Protection &amp; Positioning</b>	<b>Stealth Spectrum Foam</b> 650 lbs weight capacity. <input type="checkbox"/> 4009143 16"-20" width, 16"-24" depth \$ 580 <input type="checkbox"/> 4009144 22"-32" width, 18"-28" depth \$ 630
	<b>Stealth Tru-Comfort 2</b> 450 lbs weight capacity. <input type="checkbox"/> 4009146 16"-20" width, 16"-24" depth \$ 515 <input type="checkbox"/> 4009147 22" width, 18"-24" depth \$ 575
	<b>nxt BioFit</b> 300 lbs weight capacity. <input type="checkbox"/> 4009101 14"-20" width, 16"-20" depth \$ 595 <input type="checkbox"/> 4009102 22"-24" width, 16"-20" depth \$ 625
	<b>Matrx VI</b> 300 lbs weight capacity. <input type="checkbox"/> 4009111 14"-20" width, 15"-20" depth \$ 680 <input type="checkbox"/> 4009112 22" width, 18"-20" depth \$ 745
	<b>Matrx VI Heavy Duty</b> 600 lbs weight capacity. <input type="checkbox"/> 4009115 18"-21" width, 18"-26" depth \$ 1,050 <input type="checkbox"/> 4009116 22"-30" width, 18"-26" depth \$ 1,250

<b>Skin Protection &amp; Adjustable</b>	<b>ROHO Single Valve High Profile</b> No weight limit if cushion is properly sized. <input type="checkbox"/> 4009125 14"-20" width, 16"-22" depth \$ 840 <input type="checkbox"/> 4009126 22"-24" width, 16"-22" depth \$ 960
---	---

<b>Skin Protection, Positioning &amp; Adjustable</b>	<b>ROHO Quadro Select High Profile</b> No weight limit if cushion is properly sized. <input type="checkbox"/> 4009121 14"-20" width, 16"-22" depth \$ 840 <input type="checkbox"/> 4009122 22"-24" width, 16"-22" depth \$ 960
--	--

<b>Extra-Wide</b>  No weight limit if cushion is properly sized. Custom sizes available.  Custom-made foam cushion with fluid-proof yet breathable fabric cover.	<b>PDG Solar Basic</b> High resiliency polyfoam with 2" medium-density base and 1" lower-density top. <input type="checkbox"/> 4009200 20"-28" width, 16"-24" depth \$ 395 <input type="checkbox"/> 4009201 29"-36" width, 16"-24" depth \$ 495 <input type="checkbox"/> 4009202 37"-48" width, 16"-24" depth \$ 595
	<b>PDG Solar Contour</b> Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups. <input type="checkbox"/> 4009210 20"-28" width, 16"-24" depth \$ 695 <input type="checkbox"/> 4009211 29"-36" width, 16"-24" depth \$ 795 <input type="checkbox"/> 4009212 37"-48" width, 16"-24" depth \$ 895
	<b>PDG Solar Pressure Reduction</b> All features of the Contour, plus added pincore latex rubber and containment foam outline. <input type="checkbox"/> 4009220 20"-28" width, 16"-24" depth \$ 1,095 <input type="checkbox"/> 4009221 29"-36" width, 16"-24" depth \$ 1,195 <input type="checkbox"/> 4009222 37"-48" width, 16"-24" depth \$ 1,295

<b>3. Specify the size</b> →	<b>Cushion Size</b> Sizes based on availability from the cushion manufacturer. Specify Size (Width" x Depth") _____
------------------------------	--

### Backrest

1. Select the type ↓	2. Select the contour depth ↓	3. Select the model ↓
-------------------------	----------------------------------	--------------------------

Positioning, posterior-lateral	2" depth	<input type="checkbox"/> 4021141 <b>Stealth ADI AL Active</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 945
	4" depth	<input type="checkbox"/> 4021142 <b>Stealth ADI AL Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,045
	6" depth	<input type="checkbox"/> 4021143 <b>Stealth ADI AL Extra Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,145

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> 4021110 <b>Matrx Elite</b> 300 lbs weight capacity. 12"-20" width, 10"-20" height \$ 955
		<input type="checkbox"/> 4021111 <b>Matrx Elite Heavy Duty</b> 500 lbs weight capacity. 21" width, 16"/20" height \$ 1,195
		<input type="checkbox"/> 4021112 22"-24" width, 16"/20" height \$ 1,295
		<input type="checkbox"/> 4021114 25"-30" width, 16"/20" height \$ 1,995
	3.5" to 4" depth	<input type="checkbox"/> 4021120 <b>ROHO Agility Min</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 900
		<input type="checkbox"/> 4021100 <b>nxt Xtend</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021105 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021122 <b>ROHO Agility Mid</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 930
	6" to 7" depth	<input type="checkbox"/> 4021101 <b>nxt Xtend Deep</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 690
		<input type="checkbox"/> 4021106 22"-24" width, 16"-20" height \$ 720
		<input type="checkbox"/> 4021116 <b>Matrx Elite Deep</b> 300 lbs weight capacity. 14"-20" width, 12"-20" height \$ 1,145
		<input type="checkbox"/> 4021117 22"-24" width, 12"-20" height \$ 1,295
<input type="checkbox"/> 4021124 <b>ROHO Agility Max</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 1,085		

4. Specify the size →	<b>Backrest Size</b> Sizes based on availability from the backrest manufacturer. Specify Size (Width" x Height") _____
-----------------------	---

### Pediatric Seating System

4030400 **Diego Pediatric Seating System**  
 Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10	12	14	16
Width	(25.40)	(30.48)	(35.56)	(40.64)
Height / Depth	10 (25.40)	12 (30.48)	14 (35.56)	16 (40.64)
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓

**A Head System** Head Pad & Hardware

**B Contoured Back System** Base, Cushion & Cover

**FlipLock** Configurable as a between canes or full width set up.

**CoolCore Covers** Our Diego System uses the CoolCore cover technology which adapt all-year-round for a more comfortable experience.

**Laterals** Modular Laterals either Fixed or Swing Away

**Contoured Seat System** Base, Cushion, Cover & Hardware

**Back Hardware Cover** Easy to remove back cover for quick on-demand adjustments and a more aesthetic look.

**Width Growth** Adapted from our TFB System, the back wings allow for a pediatric growth adjustment.

**Multi-Axis Hardware** The HMO501, included as part of a package for the Diego Head System, allows growth adjustability and on-demand configurations.

**diego** a Stealth Product

**WC-20** Certified



## Headrest

Stealth Comfort Plus Headrest	<b>Stealth Hardware</b>		
	<input type="checkbox"/> 4024110	Multi-Axis HMO475	\$ 225
	<input type="checkbox"/> 4024111	Link TWB480	\$ 350
	<b>Comfort Plus Pad</b>		\$ 275
	<input type="checkbox"/> 4024115	Large (10" x 6")	<input type="checkbox"/> 4024116 X-Large (14" x 6")
nxt MultiFit Head Support	<input type="checkbox"/> 4024100 <b>nxt Hardware</b>		\$ 285
	<b>Multi-Adjustable Pad</b>		\$ 275
	<input type="checkbox"/> 4024105	Contoured (9" x 5")	<input type="checkbox"/> 4024107 X-Large (14" x 5")
	<input type="checkbox"/> 4024106	Large (10" x 5")	

## Positioning Belts & Straps

Bodypoint	<b>Hip Belts</b>		
	<input type="checkbox"/> 4029010	Center-Pull Auto Style	\$ 115
	<input type="checkbox"/> 4029015	Center-Pull Auto Style, Padded	\$ 215
	<input type="checkbox"/> 4029020	4-Point, Padded	\$ 235
	<input type="checkbox"/> 4029030	<b>Chest Strap</b>	\$ 165

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.

PDG proudly supplies on-chair seating and positioning products from these manufacturers:

