

Customer and Order Information

Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

Stellar LEAP Manual Posterior & Anterior Tilt Wheelchair Extended delivery lead times may apply.

<input checked="" type="checkbox"/> 2201000	Manual Wheelchair with 20° Posterior Tilt	\$ 3,995
<input checked="" type="checkbox"/> 2201101	25° Anterior Tilt with Linked Back 30° Anterior Tilt available with reduced Posterior Tilt range of 15°.	\$ 1,050
<input type="checkbox"/> 2201102	30° Dynamic Back Recline	\$ 500

Frame Options

<input checked="" type="checkbox"/> 2202000	Standard Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/> 2202300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths

<input type="checkbox"/> 2203014	14" Seat Width	<input type="checkbox"/> 2203015	15" Seat Width	<input type="checkbox"/> 2203016	16" Seat Width
<input type="checkbox"/> 2203017	17" Seat Width	<input type="checkbox"/> 2203018	18" Seat Width	<input type="checkbox"/> 2203019	19" Seat Width
<input type="checkbox"/> 2203020	20" Seat Width				
<input type="checkbox"/> 2203021	21" Seat Width	<input type="checkbox"/> 2203022	22" Seat Width		\$ 215

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 2204016	16" Seat Depth	<input type="checkbox"/> 2204017	17" Seat Depth	<input type="checkbox"/> 2204018	18" Seat Depth
<input type="checkbox"/> 2204019	19" Seat Depth	<input type="checkbox"/> 2204020	20" Seat Depth		
Custom Seat Depths					
<input type="checkbox"/> 2204021	21" Seat Depth	<input type="checkbox"/> 2204022	22" Seat Depth		\$ 325

Frame Colour

PDG Premium Colours

<input type="checkbox"/> 2205000	Tracer Red	<input type="checkbox"/> 2205001	Gloss Black	<input type="checkbox"/> 2205002	Silver
<input type="checkbox"/> 2205003	Black Shimmer	<input type="checkbox"/> 2205004	Orange Fusion	<input type="checkbox"/> 2205005	Karl's Imagination
<input type="checkbox"/> 2205006	Blue Shimmer	<input type="checkbox"/> 2205007	Gloss White	<input type="checkbox"/> 2205008	Dark Titanium
<input type="checkbox"/> 2205009	Geyser Green	<input type="checkbox"/> 2205010	Lollypop Blue	<input type="checkbox"/> 2205011	Lollypop Red
<input type="checkbox"/> 2205012	Bubble Gum Pink	<input type="checkbox"/> 2205013	Berry Purple	<input type="checkbox"/> 2205014	Safety Yellow
<input type="checkbox"/> 2205200	Custom Colour Contact PDG for specifications.				\$ 300

Arm Type		
<input type="checkbox"/> 2206000	No Armrest	
<input type="checkbox"/> 2206100	Single Post Fixed Height (12") Full Length Armpad only.	
Single Post Adjustable Height \$ 255		
<input type="checkbox"/> 2206111 8"-11"	<input type="checkbox"/> 2206114 10"-14"	
<input type="checkbox"/> 2206200	Dual Post Fixed Height (12") Full Length Armpad only.	
Dual Post Adjustable Height \$ 255		
<input type="checkbox"/> 2206211 8"-11"	<input type="checkbox"/> 2206214 10"-14"	<input type="checkbox"/> 2206216 12"-16"
Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery. \$ 255		
<input type="checkbox"/> 2206310 Set at 10"	<input type="checkbox"/> 2206311 Set at 11"	<input type="checkbox"/> 2206312 Set at 12"
<input type="checkbox"/> 2206313 Set at 13"	<input type="checkbox"/> 2206314 Set at 14"	

Armrests	
Regular Armpads	
<input type="checkbox"/> 2207100 Desk Length - 2" x 10"	<input type="checkbox"/> 2207101 Full Length - 2" x 14"
<input type="checkbox"/> 2207102 Full Length Premium Grade - 2" x 14"	\$ 35
Wide Armpads May interfere with wheels in some tilt configurations. \$ 115	
<input type="checkbox"/> 2207110 Desk Length - 3" x 10"	<input type="checkbox"/> 2207111 Full Length - 3" x 12"
Gel Armpads May interfere with wheels in some tilt configurations. \$ 125	
<input type="checkbox"/> 2207120 Full Length - 2" x 12"	<input type="checkbox"/> 2207121 Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125	
<input type="checkbox"/> 2207130 Right (4.5" x 18")	<input type="checkbox"/> 2207131 Left (4.5" x 18")
Armrest Extensions Not compatible with Flip-up or Single Post Armrests. \$ 60	
<input type="checkbox"/> 2207400 3" Extension - Forward	<input type="checkbox"/> 2207401 3" Extension - Rearward

Seat Pan	
<input type="checkbox"/> 2208000	No Seat Pan
<input type="checkbox"/> 2208001	Solid Aluminum Seat Pan \$ 120

Cushion	
<input type="checkbox"/> 2209000	No Cushion
nxt BioFit	
<input type="checkbox"/> 2209101 14"-20" width, 16"-20" depth	
<input type="checkbox"/> 2209102 22"-24" width, 16"-20" depth	
Cushion Size Available in even number sizes only.	
2209400	Specify Size (Width" x Depth") _____

Seat Height																																																				
Select desired Seat-to-Floor Height STF Height is adjustable within the range of the selected Rear Wheel size.																																																				
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17" 18"</td> </tr> <tr> <td>20"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17" 18"</td> </tr> </tbody> </table>	RW	6" Casters					12"	13"	14"	15"	16"		16"			15"	16"	17" 18"	20"	13"	14"	15"	16"		22"		14"	15"	16"	17"	24"			15"	16"	17" 18"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="4">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td>15"</td> <td>16"</td> <td>17" 18"</td> </tr> <tr> <td>24"</td> <td></td> <td>15"</td> <td>16"</td> <td>17" 18"</td> </tr> </tbody> </table>	RW	8" Casters				16"		15"	16"	17" 18"	24"		15"	16"	17" 18"
RW	6" Casters																																																			
12"	13"	14"	15"	16"																																																
16"			15"	16"	17" 18"																																															
20"	13"	14"	15"	16"																																																
22"		14"	15"	16"	17"																																															
24"			15"	16"	17" 18"																																															
RW	8" Casters																																																			
16"		15"	16"	17" 18"																																																
24"		15"	16"	17" 18"																																																
	Wheel size _____																																																			
	Caster size _____																																																			
	STF height _____																																																			

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2211061	6" x 1.25" Caster size			
<input type="checkbox"/>	2211063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/>	2211081	8" x 1.25" Caster size			
<input type="checkbox"/>	2211084	8" x 2.00" Caster size			\$ 135

Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2211161	6" x 1.25" Caster size			\$ 110
<input type="checkbox"/>	2211181	8" x 1.25" Caster size			\$ 120
<input type="checkbox"/>	2211184	8" x 2.00" Caster size			\$ 150

<input type="checkbox"/>	2211260	6" Soft Roll			\$ 95
--------------------------	---------	---------------------	--	--	-------

<input type="checkbox"/>	2211400	Ultra Sport Fork	Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.		\$ 395
--------------------------	---------	-------------------------	---	--	--------

Rear Wheels

Composite

<input type="checkbox"/>	2212012	12" Wheel size	<input type="checkbox"/>	2212016	16" Wheel size	<input type="checkbox"/>	2212020	20" Wheel size
<input type="checkbox"/>	2212022	22" Wheel size	<input type="checkbox"/>	2212024	24" Wheel size			

Spoke

<input type="checkbox"/>	2212120	20" Wheel size	<input type="checkbox"/>	2212122	22" Wheel size	<input type="checkbox"/>	2212124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------

<input type="checkbox"/>	2212224	24" Heavy Duty Spoke							\$ 350
--------------------------	---------	-----------------------------	--	--	--	--	--	--	--------

Axle Type

<input type="checkbox"/>	2213000	Quick Release Rear Axles			
--------------------------	---------	---------------------------------	--	--	--

<input type="checkbox"/>	2213010	Fixed Rear Axles			
--------------------------	---------	-------------------------	--	--	--

Rear Tires

<input type="checkbox"/>	2214000	Urethane			
--------------------------	---------	-----------------	--	--	--

<input type="checkbox"/>	2214010	Pneumatic	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.		\$ 60
--------------------------	---------	------------------	--	--	-------

<input type="checkbox"/>	2214020	Pneumatic with Flat-Free Inserts	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.		\$ 95
--------------------------	---------	---	--	--	-------

Hand Rims

<input type="checkbox"/>	2215000	No Hand Rims			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	2215010	Standard	20", 22", and 24" Rear Wheels only.		
--------------------------	---------	-----------------	-------------------------------------	--	--

<input type="checkbox"/>	2215020	Plastic Coated	20", 22", and 24" Rear Wheels only.		\$ 135
--------------------------	---------	-----------------------	-------------------------------------	--	--------

Wheel Locks

<input type="checkbox"/>	2216010	Push to Lock			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	2216020	Pull to Lock			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	2216030	Attendant Wheel Lock (Push to Lock)			\$ 85
--------------------------	---------	--	--	--	-------

<input type="checkbox"/>	2216040	User & Attendant Wheel Lock	Available only on 20", 22", or 24" Rear Wheels.		\$ 125
--------------------------	---------	--	---	--	--------

<input type="checkbox"/>	2216100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.		\$ 60
--------------------------	---------	------------------------------	---	--	-------

<input type="checkbox"/>	2216200	Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts.		\$ 350
--------------------------	---------	---------------------	---	--	--------

Anti-Tippers

<input checked="" type="checkbox"/>	2217010	Rear Anti-Tippers	Required Safety Option.	\$	60
<input type="checkbox"/>	2217020	Front Anti-Tippers		\$	300
<input type="checkbox"/>	2217030	Heavy Duty Rear Anti-Tipper		\$	600
Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.					

Back Post Style

<input type="checkbox"/>	2218010	Push Bar Back			
<input type="checkbox"/>	2218020	Cane Handle Back	Not compatible with Headrest Mounting Plate.		
Back Post Height Seat Pan to Handle					
<input type="checkbox"/>	2218420	20" Height	<input type="checkbox"/>	2218425	25" Height
			<input type="checkbox"/>	2218431	31" Height

Back Options

Angle/Height Adjustable Push Bar Handle				Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$	150
<input type="checkbox"/>	2219010	Button on Outside	<input type="checkbox"/>	2219011	Button on Inside	
Seat to Back Angle						
<input type="checkbox"/>	2219090	Set at 90°	<input type="checkbox"/>	2219098	Set at 98°	<input type="checkbox"/>
			<input type="checkbox"/>	2219106	Set at 106°	
<input type="checkbox"/>	2219114	Set at 114°	<input type="checkbox"/>	2219122	Set at 122°	
<input type="checkbox"/>	2219210	10" Backrest Extension	Not compatible with Angle/Height Adjustable Push Bar Handle.			\$
<input checked="" type="checkbox"/>	2219220	Fold Down Back				\$
						135

Back Upholstery

<input type="checkbox"/>	2220000	No Back Upholstery				
Standard				Back height measured from Seat Pan to top of Upholstery.	\$	155
<input type="checkbox"/>	2220018	Set at 18"	<input type="checkbox"/>	2220019	Set at 19"	<input type="checkbox"/>
			<input type="checkbox"/>	2220020	Set at 20"	
<input type="checkbox"/>	2220021	Set at 21"	<input type="checkbox"/>	2220022	Set at 22"	<input type="checkbox"/>
			<input type="checkbox"/>	2220023	Set at 23"	
Adjustable Tension				Back height measured from Seat Pan to top of Upholstery.	\$	295
<input type="checkbox"/>	2220321	Set at 21"	<input type="checkbox"/>	2220322	Set at 22"	<input type="checkbox"/>
			<input type="checkbox"/>	2220323	Set at 23"	

Backrest


<input type="checkbox"/>	2221000	No Backrest					
nxt Optima HA Active Contour							
<input type="checkbox"/>	2221100	16"-20" width, 16"-20" height				\$	
						690	
<input type="checkbox"/>	2221105	22"-24" width, 16"-20" height				\$	
						720	
nxt Optima HA Plus Contour							
<input type="checkbox"/>	2221101	16"-20" width, 16"-20" height				\$	
						690	
<input type="checkbox"/>	2221106	22"-24" width, 16"-20" height				\$	
						720	
Backrest Size Available in even number sizes only.							
	2221400	Specify Size (Width" x Height")	_____				



Tilt/Recline Options

Attendant Tilt only			
<input type="checkbox"/> 2223011	Trigger on Right Back Post	<input type="checkbox"/> 2223012	Trigger on Left Back Post
User Tilt only Full Length Arm Pads only.			\$ 100
<input type="checkbox"/> 2223021	Trigger on Right Armrest	<input type="checkbox"/> 2223022	Trigger on Left Armrest
Attendant & User Tilt Full Length Arm Pads only.			\$ 195
<input type="checkbox"/> 2223031	Trigger on Right Back Post	<input type="checkbox"/> 2223032	Trigger on Left Back Post
<input type="checkbox"/> 2223041	Trigger on Right Armrest	<input type="checkbox"/> 2223042	Trigger on Left Armrest
Anterior Tilt Limit			
<input type="checkbox"/> 2223220	Set at 20°	<input type="checkbox"/> 2223225	Set at 25°
<input type="checkbox"/> 2223230			Set at 30°
<input type="checkbox"/> 2223400	Power Tilt 250 lbs Weight Capacity. 30° Anterior to 20° Posterior tilt range. Includes 2 batteries for off-chair charging. Minimum Seat Width is 16".		\$ 1,695
	Toggle Switch	<input type="checkbox"/> 2223421	Installed on Right side
	Button Switch	<input type="checkbox"/> 2223431	Installed on Right side
	Egg Switch	<input type="checkbox"/> 2223440	Not installed
<input type="checkbox"/> 2223445	Egg Switch Swing Away Mounting Hardware		\$ 625
Not compatible with all Headrests. Specify make and model to be used: _____			

Headrest

<input type="checkbox"/> 2224000	No Headrest		
<input type="checkbox"/> 2224010	Headrest Mounting Plate		\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/> 2224100	next MultiFit Head Support Hardware		\$ 235
	next MultiFit Head Support Pad		\$ 245
<input type="checkbox"/> 2224105	Regular (9" x 5")	<input type="checkbox"/> 2224106	Large (10" x 5")
<input type="checkbox"/> 2224107			X-Large (14" x 5")

Front Rigging

<input type="checkbox"/> 1025000	No Front Rigging		
	PDG Easy-Out™ Front Rigging	Effortless insertion and removal. No Hanger Bracket remains on frame.	
<input type="checkbox"/> 2225410	Fixed 70°		Available in March
<input type="checkbox"/> 2225420	Swing In/Away Universal Angle Adjustable (0°-110°)	Set at (use increments of 10°): _____	Available in April \$ 195
<input type="checkbox"/> 2225430	Swing In/Away Elevating Legrests	High knee pivot. Includes reinforced premium Calf Pad mount.	Available in April \$ 495

Foot Plates

<input type="checkbox"/> 2226010	Composite		
<input type="checkbox"/> 2226020	Aluminum Angle Adjustable		\$ 155
<input type="checkbox"/> 2226030	One Piece Fixed	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/> 2226040	One Piece Flip-up	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/> 2226050	Interlocking Angle Adjustable	Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/> 2227040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/> 2227050	Bodypoint - Padded Calf Strap	\$ 105

Extension Tubes

Regular Extension Tubes			
<input type="checkbox"/> 2228010	Short	13"-17" for Easy-Out SA Angle Adjustable	10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging
<input type="checkbox"/> 2228020	Long	16"-22.5" for Easy-Out SA Angle Adjustable	12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging
			\$ 75

Positioning Straps

Bodypoint - Hip Belt

<input type="checkbox"/>	2229010	Center-Pull Auto Style		\$ 105
<input type="checkbox"/>	2229015	Center-Pull Auto Style, Padded		\$ 155
<input type="checkbox"/>	2229020	4-Point, Padded		\$ 175
<input type="checkbox"/>	2229030	Bodypoint - Chest Strap		\$ 85

Miscellaneous Options

<input type="checkbox"/>	2230110	IV Pole Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$ 295
<input type="checkbox"/>	2230120	O2 Holder Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$ 395
<input type="checkbox"/>	2230170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar-leap/ to learn more about the **Stellar LEAP Manual Tilt Wheelchair**

