

Customer and Order Information		
Purchase Order *	Tag	v.072318
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type	
<input checked="" type="checkbox"/> 1201000 <b>Stellar 45° Manual Tilt Wheelchair</b>	\$ 3,995

Frame Options	
<input type="checkbox"/> 1202000 <b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.	
<input type="checkbox"/> 1202050 <b>Upgraded Base and Seat Frame</b> 300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	\$ 200
<input type="checkbox"/> 1202100 <b>Heavy Duty Base and Seat Frame</b> 450 lbs Weight Capacity. Required for Seat Width of 23" or 24" Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	\$ 300
<input type="checkbox"/> 1202200 <b>Extended Base with Rear Floor Jacks</b> Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Locks, HD Anti-Tipper, IV Pole, O2 Holder.	\$ 600
<input type="checkbox"/> 1202201 <b>Shortened Strut Bracket</b> Increases heel clearance for foot propelling. Available only for Seat Width of 23" & 24".	\$ 125
<input type="checkbox"/> 1202300 <b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW
<b>Standard Seat Widths</b>		
<input type="checkbox"/> 1203014 14" Seat Width	<input type="checkbox"/> 1203015 15" Seat Width	<input type="checkbox"/> 1203016 16" Seat Width
<input type="checkbox"/> 1203017 17" Seat Width	<input type="checkbox"/> 1203018 18" Seat Width	<input type="checkbox"/> 1203019 19" Seat Width
<input type="checkbox"/> 1203020 20" Seat Width		
<input type="checkbox"/> 1203021 21" Seat Width	<input type="checkbox"/> 1203022 22" Seat Width	\$ 215
<b>Custom Seat Widths</b> Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1203023 23" Seat Width	<input type="checkbox"/> 1203024 24" Seat Width	\$ 325
<b>Use Stellar HD Order Form for Seat Widths greater than 24"</b>		

Seat Depth	
<b>Standard Seat Depths (Adjustable)</b>	
<input type="checkbox"/> 1204016 16" Seat Depth	<input type="checkbox"/> 1204017 17" Seat Depth
<input type="checkbox"/> 1204019 19" Seat Depth	<input type="checkbox"/> 1204020 20" Seat Depth
<input type="checkbox"/> 1204021 21" Seat Depth	<input type="checkbox"/> 1204022 22" Seat Depth
<b>Custom Seat Depths</b>	
<input type="checkbox"/> 1204021 21" Seat Depth	<input type="checkbox"/> 1204022 22" Seat Depth
	\$ 325

Frame Colour	
<b>PDG Premium Colours</b>	
<input type="checkbox"/> 1205000 Tracer Red	<input type="checkbox"/> 1205001 Gloss Black
<input type="checkbox"/> 1205003 Black Shimmer	<input type="checkbox"/> 1205004 Orange Fusion
<input type="checkbox"/> 1205006 Blue Shimmer	<input type="checkbox"/> 1205007 Gloss White
<input type="checkbox"/> 1205009 Geyser Green	<input type="checkbox"/> 1205010 Lollypop Blue
<input type="checkbox"/> 1205012 Bubble Gum Pink	<input type="checkbox"/> 1205013 Berry Purple
<input type="checkbox"/> 1205002 Silver	<input type="checkbox"/> 1205005 Karl's Imagination
<input type="checkbox"/> 1205008 Dark Titanium	<input type="checkbox"/> 1205011 Lollypop Red
<input type="checkbox"/> 1205014 Safety Yellow	
<input type="checkbox"/> 1205200 <b>Custom Colour</b> Contact PDG for specifications.	\$ 300

Arm Type		
<input type="checkbox"/>	1206000	<b>No Armrest</b>
<input type="checkbox"/>	1206100	<b>Single Post Fixed Height (12")</b> Full Length Armpad only.
		<b>Single Post Adjustable Height</b> Not compatible with the Heavy Duty Base and/or Seat Width greater than 22". <span style="float: right;">\$ 255</span>
<input type="checkbox"/>	1206111	8"-11" <input type="checkbox"/> 1206114 10"-14"
<input type="checkbox"/>	1206200	<b>Dual Post Fixed Height (12")</b> Full Length Armpad only.
		<b>Dual Post Adjustable Height</b> <span style="float: right;">\$ 255</span>
<input type="checkbox"/>	1206211	8"-11" <input type="checkbox"/> 1206214 10"-14" <input type="checkbox"/> 1206216 12"-16"
		<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. <span style="float: right;">\$ 255</span>
<input type="checkbox"/>	1206310	Set at 10" <input type="checkbox"/> 1206311 Set at 11" <input type="checkbox"/> 1206312 Set at 12"
<input type="checkbox"/>	1206313	Set at 13" <input type="checkbox"/> 1206314 Set at 14"

Armrests		
<b>Regular Armpads</b>		
<input type="checkbox"/>	1207100	Desk Length - 2" x 10" <input type="checkbox"/> 1207101 Full Length - 2" x 14"
<input type="checkbox"/>	1207102	Full Length Premium Grade - 2" x 14" <span style="float: right;">\$ 35</span>
		<b>Wide Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 115</span>
<input type="checkbox"/>	1207110	Desk Length - 3" x 10" <input type="checkbox"/> 1207111 Full Length - 3" x 12"
		<b>Gel Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>
<input type="checkbox"/>	1207120	Full Length - 2" x 12" <input type="checkbox"/> 1207121 Full Length - 3.5" x 12"
		<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>
<input type="checkbox"/>	1207130	Right (4.5" x 18") <input type="checkbox"/> 1207131 Left (4.5" x 18")
<input type="checkbox"/>	1207220	<b>Armrest Wideners for Dual Post Armrests</b> Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards. <span style="float: right;">\$ 375</span>
		<b>Clothing Guards</b> Not compatible with Flip-up Armrests, Armrest Wideners.
<input type="checkbox"/>	1207310	For Single Post Armrests Available for 10"-14" Height range only. <span style="float: right;">\$ 110</span>
<input type="checkbox"/>	1207320	For Dual Post Armrests <span style="float: right;">\$ 165</span>
		<b>Armrest Extensions</b> Not compatible with Flip-up or Single Post Armrests. <span style="float: right;">\$ 60</span>
<input type="checkbox"/>	1207400	3" Extension - Forward <input type="checkbox"/> 1207401 3" Extension - Rearward

Seat Pan		
<input type="checkbox"/>	1208000	<b>No Seat Pan</b>
<input type="checkbox"/>	1208001	<b>Solid Aluminum Seat Pan</b> <span style="float: right;">\$ 120</span>

Cushion		
<input type="checkbox"/>	1209000	<b>No Cushion</b>
<b>nxt BioFit</b>		
<input type="checkbox"/>	1209101	14"-20" width, 16"-20" depth <span style="float: right;">\$ 595</span>
<input type="checkbox"/>	1209102	22"-24" width, 16"-20" depth <span style="float: right;">\$ 625</span>
		<b>Cushion Size</b> Available in even number sizes only.
1209400	Specify Size (Width" x Depth") _____	

Seat Height																																																																			
<b>Select desired Seat-to-Floor Height</b> STF Height is adjustable within the range of the selected Rear Wheel size.																																																																			
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>20"</td> </tr> <tr> <td>20"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>20"</td> </tr> </tbody> </table>	RW	6" Casters						12"	14"	15"	16"	17"	18"		16"			16"	17"	18"	20"	20"	14"	15"	16"	17"	18"		22"		15"	16"	17"	18"	19"	24"			16"	17"	18"	20"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> </tbody> </table>	RW	8" Casters						16"			16"	17"	18"	19"	20"	24"			16"	17"	18"	19"	20"	Wheel size _____ Caster size _____ STF height _____
RW	6" Casters																																																																		
12"	14"	15"	16"	17"	18"																																																														
16"			16"	17"	18"	20"																																																													
20"	14"	15"	16"	17"	18"																																																														
22"		15"	16"	17"	18"	19"																																																													
24"			16"	17"	18"	20"																																																													
RW	8" Casters																																																																		
16"			16"	17"	18"	19"	20"																																																												
24"			16"	17"	18"	19"	20"																																																												

## Casters

**Urethane** 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1211061	6" x 1.25" Caster size			
<input type="checkbox"/>	1211063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/>	1211081	8" x 1.25" Caster size			
<input type="checkbox"/>	1211084	8" x 2.00" Caster size			\$ 135

**Pneumatic** 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1211161	6" x 1.25" Caster size			\$ 110
<input type="checkbox"/>	1211181	8" x 1.25" Caster size			\$ 120
<input type="checkbox"/>	1211184	8" x 2.00" Caster size			\$ 150

<input type="checkbox"/>	1211260	<b>6" Soft Roll</b>			\$ 95
--------------------------	---------	---------------------	--	--	-------

<input type="checkbox"/>	1211400	<b>Ultra Sport Fork</b>	Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.		\$ 395
--------------------------	---------	-------------------------	---	--	--------

## Rear Wheels

**Composite**

<input type="checkbox"/>	1212012	12" Wheel size	<input type="checkbox"/>	1212016	16" Wheel size	<input type="checkbox"/>	1212020	20" Wheel size
<input type="checkbox"/>	1212022	22" Wheel size	<input type="checkbox"/>	1212024	24" Wheel size			

**Spoke**

<input type="checkbox"/>	1212120	20" Wheel size	<input type="checkbox"/>	1212122	22" Wheel size	<input type="checkbox"/>	1212124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------

<input type="checkbox"/>	1212224	<b>24" Heavy Duty Spoke</b>	Required for Heavy Duty Base and Seat Frame.						\$ 350
--------------------------	---------	-----------------------------	--	--	--	--	--	--	--------

## Axle Type

<input type="checkbox"/>	1213000	<b>Quick Release Rear Axles</b>			
--------------------------	---------	---------------------------------	--	--	--

<input type="checkbox"/>	1213010	<b>Fixed Rear Axles</b>			
--------------------------	---------	-------------------------	--	--	--

## Rear Tires

<input type="checkbox"/>	1214000	<b>Urethane</b>			
--------------------------	---------	-----------------	--	--	--

<input type="checkbox"/>	1214010	<b>Pneumatic</b>	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.		\$ 60
--------------------------	---------	------------------	--	--	-------

<input type="checkbox"/>	1214020	<b>Pneumatic with Flat-Free Inserts</b>	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.		\$ 95
--------------------------	---------	---	--	--	-------

## Hand Rims

<input type="checkbox"/>	1215000	<b>No Hand Rims</b>			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	1215010	<b>Standard</b>	20", 22", and 24" Rear Wheels only.		
--------------------------	---------	-----------------	-------------------------------------	--	--

<input type="checkbox"/>	1215020	<b>Plastic Coated</b>	20", 22", and 24" Rear Wheels only.		\$ 135
--------------------------	---------	-----------------------	-------------------------------------	--	--------

## Wheel Locks

<input type="checkbox"/>	1216010	<b>Push to Lock</b>			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	1216020	<b>Pull to Lock</b>			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	1216030	<b>Attendant Wheel Lock (Push to Lock)</b>	Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.		\$ 85
--------------------------	---------	--	---	--	-------

<input type="checkbox"/>	1216040	<b>User &amp; Attendant Wheel Lock</b>	Available only on 20", 22", or 24" Rear Wheels.		\$ 125
--------------------------	---------	--	---	--	--------

<input type="checkbox"/>	1216100	<b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.		\$ 60
--------------------------	---------	------------------------------	---	--	-------

<input type="checkbox"/>	1216200	<b>Hill Holders</b>	Requires Pneumatic Rear Tires with Flat-Free Inserts.		\$ 350
--------------------------	---------	---------------------	---	--	--------

Anti-Tippers		
<input checked="" type="checkbox"/>	1217010	<b>Rear Anti-Tippers</b> Required Safety Option. \$ 60
<input type="checkbox"/>	1217020	<b>Front Anti-Tippers</b> \$ 300
<input type="checkbox"/>	1217030	<b>Heavy Duty Rear Anti-Tipper</b> \$ 600 Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.

Back Post Style		
<input type="checkbox"/>	1218010	<b>Push Bar Back</b>
<input type="checkbox"/>	1218020	<b>Cane Handle Back</b> Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.
<b>Back Post Height</b> Seat Pan to Handle		
<input type="checkbox"/>	1218425	25" Height
<input type="checkbox"/>	1218431	31" Height

Back Options		
<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/>	1219010	Button on Outside
<input type="checkbox"/>	1219011	Button on Inside
<b>Seat to Back Angle</b>		
<input type="checkbox"/>	1219090	Set at 90°
<input type="checkbox"/>	1219098	Set at 98°
<input type="checkbox"/>	1219106	Set at 106°
<input type="checkbox"/>	1219114	Set at 114°
<input type="checkbox"/>	1219122	Set at 122°
<input type="checkbox"/>	1219210	<b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Push Bar Handle. \$ 135
<input type="checkbox"/>	1219220	<b>Fold Down Back</b> Not compatible with Reclining Back.
<input type="checkbox"/>	1219230	<b>Spring Back Suspension</b> Not compatible with Reclining Back. \$ 105
<input type="checkbox"/>	1219240	<b>Dynamic Rocker Back Assembly</b> \$ 2,995 Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.

Back Upholstery		
<input type="checkbox"/>	1220000	<b>No Back Upholstery</b>
<b>Standard</b> Back height measured from Seat Pan to top of Upholstery.		\$ 155
<input type="checkbox"/>	1220021	Set at 21"
<input type="checkbox"/>	1220022	Set at 22" Not available with 25" Back.
<b>Short</b> 25" Back only. Back height measured from Seat Pan to top of Upholstery.		\$ 155
<input type="checkbox"/>	1220118	Set at 18"
<input type="checkbox"/>	1220119	Set at 19"
<input type="checkbox"/>	1220120	Set at 20"
<b>Tall</b> 31" Back only. Back height measured from Seat Pan to top of Upholstery.		\$ 155
<input type="checkbox"/>	1220224	Set at 24"
<input type="checkbox"/>	1220225	Set at 25"
<input type="checkbox"/>	1220226	Set at 26"
<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery.		\$ 295
<input type="checkbox"/>	1220321	Set at 21"
<input type="checkbox"/>	1220322	Set at 22"
<input type="checkbox"/>	1220323	Set at 23"


Backrest		
<input type="checkbox"/>	1221000	<b>No Backrest</b>
<b>nxt Xtend</b>		
<input type="checkbox"/>	1221100	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/>	1221105	22"-24" width, 16"-20" height \$ 720
<b>nxt Xtend Deep</b>		
<input type="checkbox"/>	1221101	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/>	1221106	22"-24" width, 16"-20" height \$ 720
<b>Backrest Size</b> Available in even number sizes only.		
1221400	Specify Size (Width" x Height") _____	



## Tilt/Recline Options

<b>Attendant Tilt only</b>			
<input type="checkbox"/> 1223011	Trigger on Right Back Post	<input type="checkbox"/> 1223012	Trigger on Left Back Post
<b>User Tilt only</b> Not available on widths greater than 22". Full Length Arm Pads only.			\$ 100
<input type="checkbox"/> 1223021	Trigger on Right Armrest	<input type="checkbox"/> 1223022	Trigger on Left Armrest
<b>Attendant &amp; User Tilt</b> Not available on widths greater than 22". Full Length Arm Pads only.			\$ 195
<input type="checkbox"/> 1223031	Trigger on Right Back Post	<input type="checkbox"/> 1223032	Trigger on Left Back Post
<input type="checkbox"/> 1223041	Trigger on Right Armrest	<input type="checkbox"/> 1223042	Trigger on Left Armrest
<input type="checkbox"/> 1223110	<b>Trigger Extension</b>		\$ 50
<b>Tilt Limit</b>			
<input type="checkbox"/> 1223220	Set at 20°	<input type="checkbox"/> 1223225	Set at 25°
<input type="checkbox"/> 1223235	Set at 35°	<input type="checkbox"/> 1223240	Set at 40°
<input type="checkbox"/> 1223230	Set at 30°	<input type="checkbox"/> 1223245	Set at 45°
<input type="checkbox"/> 1223410	<b>Reclining Backrest (Manual)</b>		\$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.			

## Headrest

<input type="checkbox"/> 1224000	<b>No Headrest</b>		
<input type="checkbox"/> 1224010	<b>Headrest Mounting Plate</b>		\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/> 1224100	<b>next MultiFit Head Support Hardware</b>		\$ 235
<b>next MultiFit Head Support Pad</b>			\$ 245
<input type="checkbox"/> 1224105	Regular (9" x 5")	<input type="checkbox"/> 1224106	Large (10" x 5")
		<input type="checkbox"/> 1224107	X-Large (14" x 5")

## Front Rigging

<input type="checkbox"/> 1225000	<b>No Front Rigging</b>	<input type="checkbox"/> 1225001	<b>Sector Blocks only</b>
<b>PDG Easy-Out™ Front Rigging</b>		Effortless insertion and removal. No Hanger Bracket remains on frame.	
<input type="checkbox"/> 1225410	Fixed 70°		\$ 95
<input type="checkbox"/> 1225420	Swing In/Away Universal Angle Adjustable (0°-110°)	Set at (use increments of 10°): _____	\$ 195
<input type="checkbox"/> 1225430	Swing In/Away Elevating Legrests	High knee pivot. Includes reinforced premium Calf Pad mount.	\$ 495
<b>Standard Front Rigging</b>		Legacy Design. Not compatible with Easy-Out™.	
<input type="checkbox"/> 1225110	Plug In 70°		
<input type="checkbox"/> 1225130	Swing Away 60°	<input type="checkbox"/> 1225140	Swing Away 70°
<input type="checkbox"/> 1225200	Swing Away Elevating Legrest		\$ 375
<input type="checkbox"/> 1225300	1" Offset Hangers for Swing Away		\$ 75

## Foot Plates

<input type="checkbox"/> 1226010	<b>Composite</b>		
<input type="checkbox"/> 1226020	<b>Aluminum Angle Adjustable</b>		\$ 155
<input type="checkbox"/> 1226030	<b>One Piece Fixed</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/> 1226040	<b>One Piece Flip-up</b>	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/> 1226050	<b>Interlocking Angle Adjustable</b>	Seat Widths of 14"-21" only.	\$ 195

## Footrest/Legrest Accessories

<input type="checkbox"/>	1227010	<b>Swing Away 60° &amp; 70° Calf Pad</b>	Not compatible with 60° & 70° Gel Padding.	\$	135
<input type="checkbox"/>	1227020	<b>Swing Away 60° &amp; 70° Gel Padding</b>	Not compatible with 60° & 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$	100
<input type="checkbox"/>	1227030	<b>Swing Away ELR Gel Knee Pad</b>		\$	100
<input type="checkbox"/>	1227040	<b>Adjustable Length Heel Loops</b>		\$	50
<input type="checkbox"/>	1227050	<b>Bodypoint - Padded Calf Strap</b>	Sizes greater than 22" are not padded.	\$	105
<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.					
<input type="checkbox"/>	1227060	For Plug In Front Rigging		\$	190
<input type="checkbox"/>	1227065	For Swing Away 70° Front Rigging	Not available on 13" STF Height.	\$	380

## Extension Tubes

<b>Regular Extension Tubes</b>					
<input type="checkbox"/>	1228010	Short	14"-18" for Easy-Out SA Angle Adjustable	11"-15" for Easy Out SA ELR	13"-17" for Standard Front Rigging
<input type="checkbox"/>	1228020	Long	18"-24" for Easy-Out SA Angle Adjustable	15"-21" for Easy-Out SA ELR	15"-19" for Standard Front Rigging
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only.					
<input type="checkbox"/>	1228110	Standard (13"-17")		<input type="checkbox"/>	1228120 Long (15"-19")

## Positioning Straps

<b>Bodypoint - Hip Belts</b>						
<input type="checkbox"/>	1229010	Center-Pull Auto Style			\$	105
<input type="checkbox"/>	1229015	Center-Pull Auto Style, Padded			\$	155
<input type="checkbox"/>	1229020	4-Point, Padded			\$	175
<input type="checkbox"/>	1229030	<b>Bodypoint - Chest Strap</b>			\$	85

## Miscellaneous Options

<input type="checkbox"/>	1230050	<b>Protective Lower Padding</b>		\$	425
<input type="checkbox"/>	1230080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".	\$	325
<input type="checkbox"/>	1230110	<b>IV Pole</b>	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	295
<input type="checkbox"/>	1230120	<b>O2 Holder</b>	Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$	395
<input type="checkbox"/>	1230170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.	\$	10

## Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/stellar/](http://pdgmobility.com/products/stellar/) to learn more about the **Stellar 45° Manual Tilt Wheelchair**

