

Customer and Order Information		
Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type	
<input checked="" type="checkbox"/> 1201000 Stellar 45° Manual Tilt Wheelchair	\$ 3,995

Frame Options	
<input type="checkbox"/> 1202000 Standard Base and Seat Frame 250 lbs Weight Capacity.	
<input type="checkbox"/> 1202050 Upgraded Base and Seat Frame 300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	\$ 200
<input type="checkbox"/> 1202100 Heavy Duty Base and Seat Frame 450 lbs Weight Capacity. Required for Seat Width of 23" or 24" Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	\$ 300
<input type="checkbox"/> 1202200 Extended Base with Rear Floor Jacks Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Locks, HD Anti-Tipper, IV Pole, O2 Holder.	\$ 600
<input type="checkbox"/> 1202201 Shortened Strut Bracket Increases heel clearance for foot propelling. Available only for Seat Width of 23" & 24".	\$ 125
<input type="checkbox"/> 1202300 Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW
Standard Seat Widths		
<input type="checkbox"/> 1203014 14" Seat Width	<input type="checkbox"/> 1203015 15" Seat Width	<input type="checkbox"/> 1203016 16" Seat Width
<input type="checkbox"/> 1203017 17" Seat Width	<input type="checkbox"/> 1203018 18" Seat Width	<input type="checkbox"/> 1203019 19" Seat Width
<input type="checkbox"/> 1203020 20" Seat Width		
<input type="checkbox"/> 1203021 21" Seat Width	<input type="checkbox"/> 1203022 22" Seat Width	\$ 215
Custom Seat Widths Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1203023 23" Seat Width	<input type="checkbox"/> 1203024 24" Seat Width	\$ 325
Use Stellar HD Order Form for Seat Widths greater than 24"		

Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 1204016 16" Seat Depth	<input type="checkbox"/> 1204017 17" Seat Depth
<input type="checkbox"/> 1204019 19" Seat Depth	<input type="checkbox"/> 1204020 20" Seat Depth
<input type="checkbox"/> 1204021 21" Seat Depth	<input type="checkbox"/> 1204022 22" Seat Depth
	\$ 325

Frame Colour	
PDG Premium Colours	
<input type="checkbox"/> 1205000 Tracer Red	<input type="checkbox"/> 1205001 Gloss Black
<input type="checkbox"/> 1205003 Black Shimmer	<input type="checkbox"/> 1205004 Orange Fusion
<input type="checkbox"/> 1205006 Blue Shimmer	<input type="checkbox"/> 1205007 Gloss White
<input type="checkbox"/> 1205009 Geyser Green	<input type="checkbox"/> 1205010 Lollypop Blue
<input type="checkbox"/> 1205012 Bubble Gum Pink	<input type="checkbox"/> 1205013 Berry Purple
<input type="checkbox"/> 1205002 Silver	<input type="checkbox"/> 1205005 Karl's Imagination
	<input type="checkbox"/> 1205008 Dark Titanium
	<input type="checkbox"/> 1205011 Lollypop Red
	<input type="checkbox"/> 1205014 Safety Yellow
<input type="checkbox"/> 1205200 Custom Colour Contact PDG for specifications.	\$ 300

Arm Type		
<input type="checkbox"/>	1206000	No Armrest
<input type="checkbox"/>	1206100	Single Post Fixed Height (12") Full Length Armpad only.
		Single Post Adjustable Height Not compatible with the Heavy Duty Base and/or Seat Width greater than 22". \$ 255
<input type="checkbox"/>	1206111	8"-11" <input type="checkbox"/> 1206114 10"-14"
<input type="checkbox"/>	1206200	Dual Post Fixed Height (12") Full Length Armpad only.
		Dual Post Adjustable Height \$ 255
<input type="checkbox"/>	1206211	8"-11" <input type="checkbox"/> 1206214 10"-14" <input type="checkbox"/> 1206216 12"-16"
		Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. \$ 255
<input type="checkbox"/>	1206310	Set at 10" <input type="checkbox"/> 1206311 Set at 11" <input type="checkbox"/> 1206312 Set at 12"
<input type="checkbox"/>	1206313	Set at 13" <input type="checkbox"/> 1206314 Set at 14"

Armrests		
Regular Armpads		
<input type="checkbox"/>	1207100	Desk Length - 2" x 10" <input type="checkbox"/> 1207101 Full Length - 2" x 14"
<input type="checkbox"/>	1207102	Full Length Premium Grade - 2" x 14" \$ 35
		Wide Armpads May interfere with wheels in some tilt configurations. \$ 115
<input type="checkbox"/>	1207110	Desk Length - 3" x 10" <input type="checkbox"/> 1207111 Full Length - 3" x 12"
		Gel Armpads May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1207120	Full Length - 2" x 12" <input type="checkbox"/> 1207121 Full Length - 3.5" x 12"
		Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125
<input type="checkbox"/>	1207130	Right (4.5" x 18") <input type="checkbox"/> 1207131 Left (4.5" x 18")
<input type="checkbox"/>	1207220	Armrest Wideners for Dual Post Armrests Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards. \$ 375
		Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.
<input type="checkbox"/>	1207310	For Single Post Armrests Available for 10"-14" Height range only. \$ 110
<input type="checkbox"/>	1207320	For Dual Post Armrests \$ 165
		Armrest Extensions Not compatible with Flip-up or Single Post Armrests. \$ 60
<input type="checkbox"/>	1207400	3" Extension - Forward <input type="checkbox"/> 1207401 3" Extension - Rearward

Seat Pan		
<input type="checkbox"/>	1208000	No Seat Pan
<input type="checkbox"/>	1208001	Solid Aluminum Seat Pan \$ 120

Cushion		
<input type="checkbox"/>	1209000	No Cushion
nxt BioFit		
<input type="checkbox"/>	1209101	14"-20" width, 16"-20" depth \$ 595
<input type="checkbox"/>	1209102	22"-24" width, 16"-20" depth \$ 625
		Cushion Size Available in even number sizes only.
1209400	Specify Size (Width" x Depth") _____	

Seat Height																																																																			
Select desired Seat-to-Floor Height STF Height is adjustable within the range of the selected Rear Wheel size.																																																																			
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>20"</td> </tr> <tr> <td>20"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>20"</td> </tr> </tbody> </table>	RW	6" Casters						12"	14"	15"	16"	17"	18"		16"			16"	17"	18"	20"	20"	14"	15"	16"	17"	18"		22"		15"	16"	17"	18"	19"	24"			16"	17"	18"	20"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> </tbody> </table>	RW	8" Casters						16"			16"	17"	18"	19"	20"	24"			16"	17"	18"	19"	20"	Wheel size _____ Caster size _____ STF height _____
RW	6" Casters																																																																		
12"	14"	15"	16"	17"	18"																																																														
16"			16"	17"	18"	20"																																																													
20"	14"	15"	16"	17"	18"																																																														
22"		15"	16"	17"	18"	19"																																																													
24"			16"	17"	18"	20"																																																													
RW	8" Casters																																																																		
16"			16"	17"	18"	19"	20"																																																												
24"			16"	17"	18"	19"	20"																																																												

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1211061	6" x 1.25" Caster size			
<input type="checkbox"/>	1211063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/>	1211081	8" x 1.25" Caster size			
<input type="checkbox"/>	1211084	8" x 2.00" Caster size			\$ 135

Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1211161	6" x 1.25" Caster size			\$ 110
<input type="checkbox"/>	1211181	8" x 1.25" Caster size			\$ 120
<input type="checkbox"/>	1211184	8" x 2.00" Caster size			\$ 150

<input type="checkbox"/>	1211260	6" Soft Roll			\$ 95
--------------------------	---------	---------------------	--	--	-------

<input type="checkbox"/>	1211400	Ultra Sport Fork	Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.		\$ 395
--------------------------	---------	-------------------------	---	--	--------

Rear Wheels

Composite

<input type="checkbox"/>	1212012	12" Wheel size	<input type="checkbox"/>	1212016	16" Wheel size	<input type="checkbox"/>	1212020	20" Wheel size
<input type="checkbox"/>	1212022	22" Wheel size	<input type="checkbox"/>	1212024	24" Wheel size			

Spoke

<input type="checkbox"/>	1212120	20" Wheel size	<input type="checkbox"/>	1212122	22" Wheel size	<input type="checkbox"/>	1212124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------

<input type="checkbox"/>	1212224	24" Heavy Duty Spoke	Required for Heavy Duty Base and Seat Frame.						\$ 350
--------------------------	---------	-----------------------------	--	--	--	--	--	--	--------

Axle Type

<input type="checkbox"/>	1213000	Quick Release Rear Axles			
--------------------------	---------	---------------------------------	--	--	--

<input type="checkbox"/>	1213010	Fixed Rear Axles			
--------------------------	---------	-------------------------	--	--	--

Rear Tires

<input type="checkbox"/>	1214000	Urethane			
--------------------------	---------	-----------------	--	--	--

<input type="checkbox"/>	1214010	Pneumatic	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.		\$ 60
--------------------------	---------	------------------	--	--	-------

<input type="checkbox"/>	1214020	Pneumatic with Flat-Free Inserts	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.		\$ 95
--------------------------	---------	---	--	--	-------

Hand Rims

<input type="checkbox"/>	1215000	No Hand Rims			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	1215010	Standard	20", 22", and 24" Rear Wheels only.		
--------------------------	---------	-----------------	-------------------------------------	--	--

<input type="checkbox"/>	1215020	Plastic Coated	20", 22", and 24" Rear Wheels only.		\$ 135
--------------------------	---------	-----------------------	-------------------------------------	--	--------

Wheel Locks

<input type="checkbox"/>	1216010	Push to Lock			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	1216020	Pull to Lock			
--------------------------	---------	---------------------	--	--	--

<input type="checkbox"/>	1216030	Attendant Wheel Lock (Push to Lock)	Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.		\$ 85
--------------------------	---------	--	---	--	-------

<input type="checkbox"/>	1216040	User & Attendant Wheel Lock	Available only on 20", 22", or 24" Rear Wheels.		\$ 125
--------------------------	---------	--	---	--	--------

<input type="checkbox"/>	1216100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.		\$ 60
--------------------------	---------	------------------------------	---	--	-------


<input type="checkbox"/>	1216200	Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts.		\$ 350
--------------------------	---------	---------------------	---	--	--------

Anti-Tippers		
<input checked="" type="checkbox"/>	1217010	Rear Anti-Tippers Required Safety Option. \$ 60
<input type="checkbox"/>	1217020	Front Anti-Tippers \$ 300
<input type="checkbox"/>	1217030	Heavy Duty Rear Anti-Tipper \$ 600 Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.

Back Post Style		
<input type="checkbox"/>	1218010	Push Bar Back
<input type="checkbox"/>	1218020	Cane Handle Back Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.
Back Post Height Seat Pan to Handle		
<input type="checkbox"/>	1218425	25" Height <input type="checkbox"/> 1218431 31" Height

Back Options		
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150		
<input type="checkbox"/>	1219010	Button on Outside <input type="checkbox"/> 1219011 Button on Inside
Seat to Back Angle		
<input type="checkbox"/>	1219090	Set at 90° <input type="checkbox"/> 1219098 Set at 98° <input type="checkbox"/> 1219106 Set at 106°
<input type="checkbox"/>	1219114	Set at 114° <input type="checkbox"/> 1219122 Set at 122°
<input type="checkbox"/>	1219210	10" Backrest Extension Not compatible with Angle/Height Adjustable Push Bar Handle. \$ 135
<input type="checkbox"/>	1219220	Fold Down Back Not compatible with Reclining Back.
<input type="checkbox"/>	1219230	Spring Back Suspension Not compatible with Reclining Back. \$ 105
<input type="checkbox"/>	1219240	Dynamic Rocker Back Assembly \$ 2,995 Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.

Back Upholstery		
<input type="checkbox"/>	1220000	No Back Upholstery \$ 155
Standard Back height measured from Seat Pan to top of Upholstery.		
<input type="checkbox"/>	1220021	Set at 21" <input type="checkbox"/> 1220022 Set at 22" Not available with 25" Back. \$ 155
Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.		
<input type="checkbox"/>	1220118	Set at 18" <input type="checkbox"/> 1220119 Set at 19" <input type="checkbox"/> 1220120 Set at 20" \$ 155
Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.		
<input type="checkbox"/>	1220224	Set at 24" <input type="checkbox"/> 1220225 Set at 25" <input type="checkbox"/> 1220226 Set at 26" \$ 155
Adjustable Tension Back height measured from Seat Pan to top of Upholstery.		
<input type="checkbox"/>	1220321	Set at 21" <input type="checkbox"/> 1220322 Set at 22" <input type="checkbox"/> 1220323 Set at 23" \$ 295

Backrest		
<input type="checkbox"/>	1221000	No Backrest
nxt Optima HA Active Contour		
<input type="checkbox"/>	1221100	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/>	1221105	22"-24" width, 16"-20" height \$ 720
nxt Optima HA Plus Contour 		
<input type="checkbox"/>	1221101	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/>	1221106	22"-24" width, 16"-20" height \$ 720
Backrest Size Available in even number sizes only.		
1221400	Specify Size (Width" x Height") _____	

Tilt/Recline Options

Attendant Tilt only			
<input type="checkbox"/> 1223011	Trigger on Right Back Post	<input type="checkbox"/> 1223012	Trigger on Left Back Post
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only.			\$ 100
<input type="checkbox"/> 1223021	Trigger on Right Armrest	<input type="checkbox"/> 1223022	Trigger on Left Armrest
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only.			\$ 195
<input type="checkbox"/> 1223031	Trigger on Right Back Post	<input type="checkbox"/> 1223032	Trigger on Left Back Post
<input type="checkbox"/> 1223041	Trigger on Right Armrest	<input type="checkbox"/> 1223042	Trigger on Left Armrest
<input type="checkbox"/> 1223110	Trigger Extension		\$ 50
Tilt Limit			
<input type="checkbox"/> 1223220	Set at 20°	<input type="checkbox"/> 1223225	Set at 25°
<input type="checkbox"/> 1223230		<input type="checkbox"/> 1223230	Set at 30°
<input type="checkbox"/> 1223235	Set at 35°	<input type="checkbox"/> 1223240	Set at 40°
<input type="checkbox"/> 1223245		<input type="checkbox"/> 1223245	Set at 45°
<input type="checkbox"/> 1223410	Reclining Backrest (Manual)		\$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.			

Headrest

<input type="checkbox"/> 1224000	No Headrest		
<input type="checkbox"/> 1224010	Headrest Mounting Plate		\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/> 1224100	next MultiFit Head Support Hardware		\$ 235
next MultiFit Head Support Pad			\$ 245
<input type="checkbox"/> 1224105	Regular (9" x 5")	<input type="checkbox"/> 1224106	Large (10" x 5")
<input type="checkbox"/> 1224107	X-Large (14" x 5")		

Front Rigging

<input type="checkbox"/> 1225000	No Front Rigging		
PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.			
<input type="checkbox"/> 1225410	Fixed 70°	Available in March	\$ 95
<input type="checkbox"/> 1225420	Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____	Available in April	\$ 195
<input type="checkbox"/> 1225430	Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.	Available in April	\$ 495
Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.			
<input type="checkbox"/> 1225110	Plug In 70°		
<input type="checkbox"/> 1225130	Swing Away 60°	<input type="checkbox"/> 1225140	Swing Away 70°
<input type="checkbox"/> 1225200	Swing Away Elevating Legrest		\$ 375
<input type="checkbox"/> 1225300	1" Offset Hangers for Swing Away		\$ 75

Foot Plates

<input type="checkbox"/> 1226010	Composite		
<input type="checkbox"/> 1226020	Aluminum Angle Adjustable		\$ 155
<input type="checkbox"/> 1226030	One Piece Fixed Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.		\$ 350
<input type="checkbox"/> 1226040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.		\$ 395
<input type="checkbox"/> 1226050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.		\$ 195

Footrest/Legrest Accessories		
<input type="checkbox"/>	1227010	Swing Away 60° & 70° Calf Pad Not compatible with 60° & 70° Gel Padding. \$ 135
<input type="checkbox"/>	1227020	Swing Away 60° & 70° Gel Padding Not compatible with 60° & 70° Swing Away Calf Pad, Easy-Out Front Rigging. \$ 100
<input type="checkbox"/>	1227030	Swing Away ELR Gel Knee Pad \$ 100
<input type="checkbox"/>	1227040	Adjustable Length Heel Loops \$ 50
<input type="checkbox"/>	1227050	Bodypoint - Padded Calf Strap Sizes greater than 22" are not padded. \$ 105
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.		
<input type="checkbox"/>	1227060	For Plug In Front Rigging \$ 190
<input type="checkbox"/>	1227065	For Swing Away 70° Front Rigging Not available on 13" STF Height. \$ 380

Extension Tubes		
Regular Extension Tubes		
<input type="checkbox"/>	1228010	Short <u>13"-17"</u> for Easy-Out SA Angle Adjustable <u>10"-13.5"</u> for Easy Out SA ELR <u>13"-17"</u> for Standard Front Rigging
<input type="checkbox"/>	1228020	Long <u>16"-22.5"</u> for Easy-Out SA Angle Adjustable <u>12.5"-19.5"</u> for Easy-Out SA ELR <u>15"-19"</u> for Standard Front Rigging \$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. \$ 100		
<input type="checkbox"/>	1228110	Standard (13"-17")
<input type="checkbox"/>	1228120	Long (15"-19")

Positioning Straps		
Bodypoint - Hip Belts		
<input type="checkbox"/>	1229010	Center-Pull Auto Style \$ 105
<input type="checkbox"/>	1229015	Center-Pull Auto Style, Padded \$ 155
<input type="checkbox"/>	1229020	4-Point, Padded \$ 175
<input type="checkbox"/>	1229030	Bodypoint - Chest Strap \$ 85



Miscellaneous Options		
<input type="checkbox"/>	1230050	Protective Lower Padding \$ 425
<input type="checkbox"/>	1230080	Gel Hip Guides \$ 325
<input type="checkbox"/>	1230110	IV Pole Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks. \$ 295
<input type="checkbox"/>	1230120	O2 Holder Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks. \$ 395
<input type="checkbox"/>	1230170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting. \$ 10

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar/ to learn more about the **Stellar 45° Manual Tilt Wheelchair**

