

## Customer and Order Information

Purchase Order *	Tag	v.0601
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <span style="float: right;">*required</span>

## Frame Type

<input checked="" type="checkbox"/> 2401000	<b>Bentley LT 20° Lightweight Manual Tilt Wheelchair</b> Aluminum Frame. Not recommended for agitated users.	E1161	\$ 3,295
<span style="color: orange;">ⓘ</span> For agitated or heavier users, choose the standard <b>Bentley</b> with steel frame. See separate order form.			

## Frame Options

<input checked="" type="checkbox"/> 2402000	<b>Lightweight Base and Seat Frame</b> 250 lbs Weight Capacity.
<input type="checkbox"/> 2402300	<b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

## Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

<b>Standard Seat Widths</b>			
<input type="checkbox"/> 1003014	14" Seat Width	<input type="checkbox"/> 1003015	15" Seat Width
<input type="checkbox"/> 1003017	17" Seat Width	<input type="checkbox"/> 1003018	18" Seat Width
<input type="checkbox"/> 1003020	20" Seat Width	<input type="checkbox"/> 1003021	21" Seat Width
<input type="checkbox"/> 1003016	16" Seat Width	<input type="checkbox"/> 1003019	19" Seat Width
<input type="checkbox"/> 1003022	22" Seat Width	E2201	\$ 350

## Seat Depth

<b>Standard Seat Depths (Adjustable)</b>			
<input type="checkbox"/> 1004016	16" Seat Depth	<input type="checkbox"/> 1004017	17" Seat Depth
<input type="checkbox"/> 1004019	19" Seat Depth	<input type="checkbox"/> 1004018	18" Seat Depth
<input type="checkbox"/> 1004020	20" Seat Depth	E2203	\$ 500
<b>Custom Seat Depths</b>			
<input type="checkbox"/> 1004021	21" Seat Depth	E2203	\$ 500
<input type="checkbox"/> 1004022	22" Seat Depth	E2204	\$ 750

## Frame Color

<b>PDG Premium Colors</b>			
<input type="checkbox"/> 2405000	Tracer Red	<input type="checkbox"/> 1005001	Gloss Black
<input type="checkbox"/> 2405003	Black Shimmer	<input type="checkbox"/> 1005002	Silver
<input type="checkbox"/> 2405006	Blue Shimmer	<input type="checkbox"/> 1005004	Orange Fusion
<input type="checkbox"/> 2405009	Geyser Green	<input type="checkbox"/> 1005007	Gloss White
<input type="checkbox"/> 2405012	Bubble Gum Pink	<input type="checkbox"/> 1005008	Dark Titanium
<input type="checkbox"/> 1005010	Lollypop Blue	<input type="checkbox"/> 1005011	Lollypop Red
<input type="checkbox"/> 1005013	Berry Purple	<input type="checkbox"/> 1005014	Safety Yellow
<input type="checkbox"/> 2405200	<b>Custom Color</b> Contact PDG for specifications.		\$ 300

## Arm Type

<input type="checkbox"/> 2406000	<b>No Armrest</b>
<input type="checkbox"/> 2406100	<b>Single Post Fixed Height (12")</b> Full Length Armpad only.
<b>Single Post Adjustable Height</b> <span style="float: right;">E0973 \$ 275</span>	
<input type="checkbox"/> 2406111	8"-11"
<input type="checkbox"/> 1006114	10"-14"

## Armrests

### Regular Armrests

- 2407100 Desk Length - 2" x 10"       1007101 Full Length - 2" x 14"

### Wide Armrests

May interfere with wheels in some tilt configurations.

\$ 200

- 2407110 Desk Length - 3" x 10"       1007111 Full Length - 3" x 12"

### Gel Armrests

May interfere with wheels in some tilt configurations.

\$ 135

- 2407120 Full Length - 2" x 12"       1007121 Full Length - 3.5" x 12"

### Contour Arm Trough with Hand Pad

May interfere with wheels in some tilt configurations.

E2209 \$ 125

- 2407130 Right (4.5" x 18")       1007131 Left (4.5" x 18")

- 2407310 **Clothing Guards** Available for 10"-14" Height range only.

\$ 110

## Seat Pan

- 2408000 **No Seat Pan**

- 2408001 **Solid Aluminum Seat Pan**

E2231 \$ 120

## Cushion

Refer to the enclosed Seating & Positioning order form for PDG's on-chair Cushion options. Standard PDG discounts apply.

## Seat Height

➤ **Select desired Seat-to-Floor Height** STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High					
	12"	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	15"	16"	17"	18"	19"	20"
8" Casters	RW	Low					High					Extra-High					
	16"	15"	16"	17"	18"	19"	17"	18"	19"	17"	18"	19"	20"	17"	18"	19"	20"

Wheel size \_\_\_\_\_

Caster size \_\_\_\_\_

STF height \_\_\_\_\_

Range \_\_\_\_\_

## Casters

### Urethane

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2411061 6" x 1.25" Caster size  
 2411063 6" x 1.50" Caster size  
 2411081 8" x 1.25" Caster size  
 2411084 8" x 2.00" Caster size

\$ 85

\$ 135

### Pneumatic

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2411161 6" x 1.25" Caster size  
 2411181 8" x 1.25" Caster size  
 2411184 8" x 2.00" Caster size

E2214 & E2215 \$ 110

E2214 & E2215 \$ 120

E2214 & E2215 \$ 150

- 2411260 **6" Soft Roll**

\$ 95

- 2411400 **Ultra Sport Fork** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.

E1015 \$ 395

## Rear Wheels

### Composite

- 2412012 12" Wheel size       1012016 16" Wheel size       1012020 20" Wheel size  
 2412022 22" Wheel size       1012024 24" Wheel size

### Spoke

\$ 75

- 2412120 20" Wheel size       1012122 22" Wheel size       1012124 24" Wheel size

Axle Type	
<input type="checkbox"/> 2413000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/> 2413010	<b>Fixed Rear Axles</b>

Rear Tires	
<input type="checkbox"/> 2414000	<b>Urethane</b>
<input type="checkbox"/> 2414010	<b>Pneumatic</b> May result in a finished Seat Height 1/2" greater than stated. <span style="float: right;">E2211 &amp; E2212 \$ 60</span>
<input type="checkbox"/> 2414020	<b>Pneumatic with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. <span style="float: right;">E2211 &amp; E2213 \$ 95</span>

Hand Rims	
<input type="checkbox"/> 2415000	<b>No Hand Rims</b>
<input type="checkbox"/> 2415010	<b>Standard</b> 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 2415020	<b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only. <span style="float: right;">\$ 135</span>

Wheel Locks	
<input type="checkbox"/> 2416010	<b>Push to Lock</b>
<input type="checkbox"/> 2416020	<b>Pull to Lock</b>
<input type="checkbox"/> 2416030	<b>Attendant Wheel Lock (Push to Lock)</b> Not compatible with TTDS. <span style="float: right;">\$ 85</span>
<input type="checkbox"/> 2416040	<b>User &amp; Attendant Wheel Lock</b> Available only on 20", 22", or 24" Rear Wheels. <span style="float: right;">\$ 125</span>
<input type="checkbox"/> 2416100	<b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock. <span style="float: right;">E0961 \$ 60</span>
<input type="checkbox"/> 2416200	<b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts. <span style="float: right;">E0974 \$ 450</span>

Anti-Tippers	
<input checked="" type="checkbox"/> 2417010	<b>Rear Anti-Tippers</b> Required Safety Option. <span style="float: right;">E0971 \$ 60</span>

Back Post Style	
<input checked="" type="checkbox"/> 2418010	<b>Push Bar Back</b>
<b>Back Post Height</b> Seat Pan to Handle	
<input type="checkbox"/> 2418425	25" Height <input type="checkbox"/> 1018431 31" Height

Back Options	
<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate. <span style="float: right;">\$ 150</span>	
<input type="checkbox"/> 2419010	Button on Outside <input type="checkbox"/> 1019011 Button on Inside
<b>Seat to Back Angle</b>	
<input type="checkbox"/> 2419090	Set at 90° <input type="checkbox"/> 1019098 Set at 98° <input type="checkbox"/> 1019106 Set at 106°
<input type="checkbox"/> 2419114	Set at 114° <input type="checkbox"/> 1019122 Set at 122°
<input type="checkbox"/> 2419210	<b>10" Backrest Extension</b> Requires the Angle/Height Adjustable Push Bar Handle. <span style="float: right;">E0966 \$ 150</span>
<input checked="" type="checkbox"/> 2419220	<b>Fold Down Back</b> Not compatible with Reclining Back.

Back Upholstery	
<input type="checkbox"/> 2420000	<b>No Back Upholstery</b>
<b>Standard</b> Back height measured from Seat Pan to top of Upholstery. <span style="float: right;">E2611 / E2612 \$ 175</span>	
<input type="checkbox"/> 2420021	Set at 21" <input type="checkbox"/> 1020022 Set at 22" <input type="checkbox"/> 1020023 Set at 23" Not available with 25" Back.
<b>Short</b> 25" Back only. Back height measured from Seat Pan to top of Upholstery. <span style="float: right;">E2611 / E2612 \$ 175</span>	
<input type="checkbox"/> 2420118	Set at 18" <input type="checkbox"/> 1020119 Set at 19" <input type="checkbox"/> 1020120 Set at 20"

See next page for more options >

<input type="checkbox"/> 2420224	<b>Tall</b> 31" Back only. Back height measured from Seat Pan to top of Upholstery. Set at 24"	<input type="checkbox"/> 1020225	Set at 25"	<input type="checkbox"/> 1020226	Set at 26"	E2611 / E2612	\$ 175
<input type="checkbox"/> 2420321	<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery. Set at 21"	<input type="checkbox"/> 1020322	Set at 22"	<input type="checkbox"/> 1020323	Set at 23"	E2611 / E2612	\$ 295

### Positioning Backrest

**i** Refer to the enclosed Seating & Positioning order form for PDG's on-chair Backrest options. Standard PDG discounts apply.

### Tilt/Recline Options

<input type="checkbox"/> 2423011	<b>Attendant Tilt only</b> Trigger on Right Back Post	<input type="checkbox"/> 1023012	Trigger on Left Back Post				
<input type="checkbox"/> 2423021	<b>User Tilt only</b> Full Length Arm Pads only. Trigger on Right Armrest	<input type="checkbox"/> 1023022	Trigger on Left Armrest				\$ 100
<input type="checkbox"/> 2423031	<b>Attendant &amp; User Tilt</b> Full Length Arm Pads only. Trigger on Right Back Post	<input type="checkbox"/> 1023032	Trigger on Left Back Post				\$ 195
<input type="checkbox"/> 2423041	Trigger on Right Armrest	<input type="checkbox"/> 1023042	Trigger on Left Armrest				
<input type="checkbox"/> 2423110	<b>Trigger Extension</b>						\$ 50
<input type="checkbox"/> 2423310	<b>Reclining Backrest (Manual)</b> 30° maximum recline. Not compatible with Backrest Extension, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery. Requires longer Frame. Maximum Seat Depth is 19".					E1225	\$ 895

### Headrest

**i** Refer to the enclosed Seating & Positioning order form for PDG's on-chair Headrest options. Standard PDG discounts apply.

<input type="checkbox"/> 2424010	<b>Headrest Mounting Plate</b> Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back.						\$ 95
----------------------------------	---	--	--	--	--	--	-------

### Front Rigging

<input type="checkbox"/> 2425000	<b>No Front Rigging</b>	<input type="checkbox"/> 2425001	<b>Sector Blocks only</b>				
<input type="checkbox"/> 2425410	<b>PDG Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a> Fixed 70°						\$ 95
<input type="checkbox"/> 2425420	Universal Angle Adjustable (0°-110°) Adjust angle by loosening one bolt. Tool-free leg-length adjustment.		Set at angle: _____ (use increments of 10°)				\$ 295
<input type="checkbox"/> 2425130	<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out. Swing Away 60°	<input type="checkbox"/> 2425140	Swing Away 70°				
<input type="checkbox"/> 2425200	Swing Away Elevating Legrests Includes Calf Pads.					E0990	\$ 375
<input type="checkbox"/> 2425300	1" Offset Hangers for Swing Away						\$ 75
<input type="checkbox"/> 2425110	<b>Standard Plug In 70°</b> Legacy Design. Not compatible with Easy-Out.						

### Foot Plates

<input type="checkbox"/> 2426010	<b>Composite</b>						
<input type="checkbox"/> 2426020	<b>Aluminum Angle Adjustable</b>					K0040	\$ 155
<input type="checkbox"/> 2426030	<b>One Piece Fixed</b> Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.						\$ 350
<input type="checkbox"/> 2426040	<b>One Piece Flip-up</b> Requires selection of Composite or Aluminum Angle Adjustable.						\$ 395
<input type="checkbox"/> 2426050	<b>Interlocking Angle Adjustable</b> Seat Widths of 14"-21" only.					K0040	\$ 195

## Footrest/Legrest Accessories

<b>Calf Pad</b> Not compatible with Gel Padding.		
<input type="checkbox"/> 2427012	For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/> 2427010	For Swing Away Front Rigging	\$ 135
<input type="checkbox"/> 2427020	<b>Gel Padding</b> Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
<b>Gel Knee Pad</b>		
<input type="checkbox"/> 2427032	For Easy-Out™ Front Rigging	\$ 200
<input type="checkbox"/> 2427030	For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/> 2427040	<b>Adjustable Length Heel Loops</b>	E0951 \$ 50
<input type="checkbox"/> 2427050	<b>Padded Calf Strap</b>	K0038 \$ 195

## Extension Tubes

<b>Regular Extension Tubes</b>		
<input type="checkbox"/> 2428010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/> 2428020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75

## Positioning Belts & Straps

Refer to the enclosed Seating & Positioning order form for PDG's on-chair Belts/Straps options. Standard PDG discounts apply.

## Miscellaneous Options

<b>One-Arm Drive</b> Requires 24" Composite Wheels with Urethane Tires. 18" STF minimum. Limits tilt. Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.		E0958	\$ 3,500
<input type="checkbox"/> 2430010	Installed on Right side	<input type="checkbox"/> 1030011	Installed on Left side
<input type="checkbox"/> 2430050	<b>Protective Lower Padding</b>		\$ 825
<input type="checkbox"/> 2430080	<b>Gel Hip Guides</b> Reduces available Seat Width by 2".	E0956 & E1028	\$ 375
<input type="checkbox"/> 2430170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.		\$ 10

## Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/bentley/](http://pdgmobility.com/products/bentley/) to learn more about the **Bentley 20° Manual Tilt Wheelchair**



### Cushion

1. Select the type ↓	2. Select the model ↓
-------------------------	--------------------------

<b>General Use</b>	<p><b>Stealth Simplicity</b> 350 lbs weight capacity.</p> <p><input type="checkbox"/> 4009141 16"-20" width, 16"-20" depth <span style="float: right;">E2601 \$ 230</span></p> <p><input type="checkbox"/> 4009142 22"-24" width, 18"-20" depth <span style="float: right;">E2602 \$ 290</span></p>
--------------------	---

<b>Positioning</b>	<p><b>ROHO AirLITE</b> 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009120 14"-20" width, 14"-18" depth <span style="float: right;">E2605 \$ 375</span></p>
--------------------	--

<b>Skin Protection &amp; Positioning</b>	<p><b>Stealth Spectrum Foam</b> 650 lbs weight capacity.</p> <p><input type="checkbox"/> 4009143 16"-20" width, 16"-24" depth <span style="float: right;">E2607 \$ 550</span></p> <p><input type="checkbox"/> 4009144 22"-32" width, 18"-28" depth <span style="float: right;">E2608 \$ 625</span></p>
	<p><b>Stealth Tru-Comfort 2</b> 450 lbs weight capacity.</p> <p><input type="checkbox"/> 4009146 16"-20" width, 16"-24" depth <span style="float: right;">E2607 \$ 490</span></p> <p><input type="checkbox"/> 4009147 22" width, 18"-24" depth <span style="float: right;">E2608 \$ 590</span></p>
	<p><b>nxt BioFit</b> 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009101 14"-20" width, 16"-20" depth <span style="float: right;">E2607 \$ 555</span></p> <p><input type="checkbox"/> 4009102 22"-24" width, 16"-20" depth <span style="float: right;">E2608 \$ 585</span></p>
	<p><b>Matrx VI</b> 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009111 14"-20" width, 15"-20" depth <span style="float: right;">E2607 \$ 680</span></p> <p><input type="checkbox"/> 4009112 22" width, 18"-20" depth <span style="float: right;">E2608 \$ 745</span></p>
	<p><b>Matrx VI Heavy Duty</b> 600 lbs weight capacity.</p> <p><input type="checkbox"/> 4009115 18"-21" width, 18"-26" depth <span style="float: right;">E2607 \$ 1,050</span></p> <p><input type="checkbox"/> 4009116 22"-30" width, 18"-26" depth <span style="float: right;">E2608 \$ 1,250</span></p>

<b>Skin Protection &amp; Adjustable</b>	<p><b>ROHO Single Valve High Profile</b> No weight limit if cushion is properly sized.</p> <p><input type="checkbox"/> 4009125 14"-20" width, 16"-22" depth <span style="float: right;">E2622 \$ 840</span></p> <p><input type="checkbox"/> 4009126 22"-24" width, 16"-22" depth <span style="float: right;">E2623 \$ 960</span></p>
---	--

<b>Skin Protection, Positioning &amp; Adjustable</b>	<p><b>ROHO Quadro Select High Profile</b> No weight limit if cushion is properly sized.</p> <p><input type="checkbox"/> 4009121 14"-20" width, 16"-22" depth <span style="float: right;">E2624 \$ 840</span></p> <p><input type="checkbox"/> 4009122 22"-24" width, 16"-22" depth <span style="float: right;">E2625 \$ 960</span></p>
--	---

<b>Extra-Wide</b>	<p><b>PDG Solar Basic</b></p> <p>High resiliency polyfoam with 2" medium-density base and 1" lower-density top.</p>	<p><input type="checkbox"/> 4009200 20"-28" width, 16"-24" depth <span style="float: right;">\$ 395</span></p> <p><input type="checkbox"/> 4009201 29"-36" width, 16"-24" depth <span style="float: right;">\$ 495</span></p> <p><input type="checkbox"/> 4009202 37"-48" width, 16"-24" depth <span style="float: right;">\$ 595</span></p>
	<p><b>PDG Solar Contour</b></p> <p>Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups.</p>	<p><input type="checkbox"/> 4009210 20"-28" width, 16"-24" depth <span style="float: right;">\$ 695</span></p> <p><input type="checkbox"/> 4009211 29"-36" width, 16"-24" depth <span style="float: right;">\$ 795</span></p> <p><input type="checkbox"/> 4009212 37"-48" width, 16"-24" depth <span style="float: right;">\$ 895</span></p>
	<p><b>PDG Solar Pressure Reduction</b></p> <p>All features of the Contour, plus added pincore latex rubber and containment foam outline.</p>	<p><input type="checkbox"/> 4009220 20"-28" width, 16"-24" depth <span style="float: right;">\$ 1,095</span></p> <p><input type="checkbox"/> 4009221 29"-36" width, 16"-24" depth <span style="float: right;">\$ 1,195</span></p> <p><input type="checkbox"/> 4009222 37"-48" width, 16"-24" depth <span style="float: right;">\$ 1,295</span></p>

No weight limit if cushion is properly sized. Custom sizes available.

Custom-made foam cushion with fluid-proof yet breathable fabric cover.

3. Specify the size →	<p><b>Cushion Size</b> Sizes based on availability from the cushion manufacturer.</p> <p>Specify Size (Width" x Depth") _____</p>
-----------------------	---

### Backrest

1. Select the type ↓	2. Select the contour depth ↓	3. Select the model ↓
-------------------------	----------------------------------	--------------------------

Positioning, posterior-lateral	2" depth	<input type="checkbox"/> 4021141 <b>Stealth ADI AL Active</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height <span style="float:right">E2615 \$ 790</span>
	4" depth	<input type="checkbox"/> 4021142 <b>Stealth ADI AL Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height <span style="float:right">E2615 \$ 875</span>
	6" depth	<input type="checkbox"/> 4021143 <b>Stealth ADI AL Extra Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height <span style="float:right">E2615 \$ 960</span>

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> 4021110 <b>Matrx Elite</b> 300 lbs weight capacity. 12"-20" width, 10"-20" height <span style="float:right">E2620 \$ 955</span>
		<input type="checkbox"/> 4021111 <b>Matrx Elite Heavy Duty</b> 500 lbs weight capacity. 21" width, 16"/20" height <span style="float:right">E2620 \$ 1,195</span>
		<input type="checkbox"/> 4021112 22"-24" width, 16"/20" height <span style="float:right">E2621 \$ 1,295</span>
		<input type="checkbox"/> 4021114 25"-30" width, 16"/20" height <span style="float:right">E2621 \$ 1,995</span>
	3.5" to 4" depth	<input type="checkbox"/> 4021120 <b>ROHO Agility Min</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height <span style="float:right">E2620 \$ 900</span>
		<input type="checkbox"/> 4021100 <b>nxt Xtend</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height <span style="float:right">E2620 \$ 635</span>
		<input type="checkbox"/> 4021105 22"-24" width, 16"-20" height <span style="float:right">E2621 \$ 660</span>
		<input type="checkbox"/> 4021122 <b>ROHO Agility Mid</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height <span style="float:right">E2620 \$ 930</span>
	6" to 7" depth	<input type="checkbox"/> 4021101 <b>nxt Xtend Deep</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height <span style="float:right">E2620 \$ 635</span>
		<input type="checkbox"/> 4021106 22"-24" width, 16"-20" height <span style="float:right">E2621 \$ 660</span>
		<input type="checkbox"/> 4021116 <b>Matrx Elite Deep</b> 300 lbs weight capacity. 14"-20" width, 12"-20" height <span style="float:right">E2620 \$ 1,145</span>
		<input type="checkbox"/> 4021117 22"-24" width, 12"-20" height <span style="float:right">E2621 \$ 1,295</span>
<input type="checkbox"/> 4021124 <b>ROHO Agility Max</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height <span style="float:right">E2620 \$ 1,085</span>		

4. Specify the size →	<b>Backrest Size</b> Sizes based on availability from the backrest manufacturer. Specify Size (Width" x Height") _____
-----------------------	---

### Pediatric Seating System

4030400 **Diego Pediatric Seating System**  
 Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10	12	14	16
Width	(25.40)	(30.48)	(35.56)	(40.64)
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓

**A Head System** Head Pad & Hardware

**B Contoured Back System** Base, Cushion & Cover

**C Contoured Seat System** Base, Cushion, Cover & Hardware

**FlipLock** Configurable as a between cases or full width set up.

**Width Growth** Adapted from our TEB System, the back wings allow for a pediatric growth adjustment.

**Back Hardware Cover** Easy to remove back cover for quick on-demand adjustments and a more aesthetic look.

**Multi-Axis Hardware** The IMOSOL included as part of a package for the Diego Head System, allows growth adjustability and on-demand configurations.

**Coolcore Covers** Our Diego System uses the Coolcore cover technology which adapt all-year-round for a more comfortable experience.

**Laterals** Modular Laterals either Fixed or Swing Away

**diego** a Stealth Product

**WC-20** Certified

### Headrest

Stealth Comfort Plus Headrest	<b>Stealth Hardware</b>		
	<input type="checkbox"/> 4024110	Multi-Axis HMO475	E1028 \$ 225
	<input type="checkbox"/> 4024111	Link TWB480	E1028 \$ 350
	<b>Comfort Plus Pad</b>		
<input type="checkbox"/> 4024116	Large (10" x 6")	<input type="checkbox"/> 4024117	X-Large (14" x 6")
nxt MultiFit Head Support	<b>nxt Hardware</b>		
	<input type="checkbox"/> 4024100		E1028 \$ 285
	<b>Multi-Adjustable Pad</b>		
<input type="checkbox"/> 4024105	Contoured (9" x 5")	E0955 \$ 275	
<input type="checkbox"/> 4024106	Large (10" x 5")	<input type="checkbox"/> 4024107	X-Large (14" x 5")

### Positioning Belts & Straps

Bodypoint	<b>Hip Belts</b>		
	<input type="checkbox"/> 4029010	Center-Pull Auto Style	E0978 \$ 115
	<input type="checkbox"/> 4029015	Center-Pull Auto Style, Padded	E0978 \$ 215
	<input type="checkbox"/> 4029020	4-Point, Padded	E0978 \$ 235
<input type="checkbox"/> 4029030	<b>Chest Strap</b>	E0960 \$ 165	

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.

PDG proudly supplies on-chair seating and positioning products from these manufacturers:

