

Customer and Order Information	
Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address * *required	

Frame Type	
<input checked="" type="checkbox"/> 1001000 Bentley 20° Manual Tilt Wheelchair	E1161 \$ 2,995

Frame Options	
<input type="checkbox"/> 1002000 Standard Base and Seat Frame 250 lbs Weight Capacity.	
<input type="checkbox"/> 1002050 Upgraded Base and Seat Frame 300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	K0108 \$ 200
<input type="checkbox"/> 1002100 Heavy Duty Base and Seat Frame 450 lbs Weight Capacity. Available for all sizes, but required for Seat Width of 23" or wider. Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	K0108 \$ 300
<input type="checkbox"/> 1002200 Extended Base with Rear Floor Jacks Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant wheel lock, HD Anti-Tipper, IV Pole, O2 Holder.	K0108 \$ 600
<input type="checkbox"/> 1002300 Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW
Standard Seat Widths		
<input type="checkbox"/> 1003014 14" Seat Width	<input type="checkbox"/> 1003015 15" Seat Width	<input type="checkbox"/> 1003016 16" Seat Width
<input type="checkbox"/> 1003017 17" Seat Width	<input type="checkbox"/> 1003018 18" Seat Width	<input type="checkbox"/> 1003019 19" Seat Width
<input type="checkbox"/> 1003020 20" Seat Width	<input type="checkbox"/> 1003021 21" Seat Width	<input type="checkbox"/> 1003022 22" Seat Width E2201 \$ 325
Custom Seat Widths Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1003023 23" Seat Width E2201 \$ 325		
<input type="checkbox"/> 1003024 24" Seat Width	<input type="checkbox"/> 1003025 25" Seat Width	<input type="checkbox"/> 1003026 26" Seat Width E2202 \$ 525
<input type="checkbox"/> 1003027 27" Seat Width E2202 \$ 525		
<input type="checkbox"/> 1003028 28" Seat Width	<input type="checkbox"/> 1003029 29" Seat Width	<input type="checkbox"/> 1003030 30" Seat Width K0108 \$ 925
<input type="checkbox"/> 1003031 31" Seat Width K0108 \$ 925	<input type="checkbox"/> 1003032 32" Seat Width	

Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 1004016 16" Seat Depth	<input type="checkbox"/> 1004017 17" Seat Depth
<input type="checkbox"/> 1004019 19" Seat Depth	<input type="checkbox"/> 1004018 18" Seat Depth
<input type="checkbox"/> 1004020 20" Seat Depth E2203 \$ 475	
Custom Seat Depths	
<input type="checkbox"/> 1004021 21" Seat Depth E2203 \$ 475	
<input type="checkbox"/> 1004022 22" Seat Depth E2204 \$ 750	

Frame Color	
PDG Premium Colors	
<input type="checkbox"/> 1005000 Tracer Red	<input type="checkbox"/> 1005001 Gloss Black
<input type="checkbox"/> 1005003 Black Shimmer	<input type="checkbox"/> 1005004 Orange Fusion
	<input type="checkbox"/> 1005002 Silver
	<input type="checkbox"/> 1005005 Karl's Imagination

See next page for more Frame Color options >

<input type="checkbox"/> 1005006 Blue Shimmer	<input type="checkbox"/> 1005007 Gloss White	<input type="checkbox"/> 1005008 Dark Titanium
<input type="checkbox"/> 1005009 Geyser Green	<input type="checkbox"/> 1005010 Lollypop Blue	<input type="checkbox"/> 1005011 Lollypop Red
<input type="checkbox"/> 1005012 Bubble Gum Pink	<input type="checkbox"/> 1005013 Berry Purple	<input type="checkbox"/> 1005014 Safety Yellow
<input type="checkbox"/> 1005200 Custom Color Contact PDG for specifications.		\$ 300

Arm Type		
<input type="checkbox"/> 1006000 No Armrest		
<input type="checkbox"/> 1006100 Single Post Fixed Height (12")	Full Length Armpad only.	
Single Post Adjustable Height Not compatible with the Heavy Duty Base and/or Seat Width greater than 22".		<i>E0973</i> \$ 255
<input type="checkbox"/> 1006111 8"-11"	<input type="checkbox"/> 1006114 10"-14"	
<input type="checkbox"/> 1006200 Dual Post Fixed Height (12")	Full Length Armpad only.	
Dual Post Adjustable Height		<i>E0973</i> \$ 255
<input type="checkbox"/> 1006211 8"-11"	<input type="checkbox"/> 1006214 10"-14"	<input type="checkbox"/> 1006216 12"-16"
Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.		<i>E0973</i> \$ 255
<input type="checkbox"/> 1006310 Set at 10"	<input type="checkbox"/> 1006311 Set at 11"	<input type="checkbox"/> 1006312 Set at 12"
<input type="checkbox"/> 1006313 Set at 13"	<input type="checkbox"/> 1006314 Set at 14"	

Armrests		
Regular Armpads		
<input type="checkbox"/> 1007100 Desk Length - 2" x 10"	<input type="checkbox"/> 1007101 Full Length - 2" x 14"	
<input type="checkbox"/> 1007102 Full Length Premium Grade - 2" x 14"		\$ 35
Wide Armpads May interfere with wheels in some tilt configurations.		\$ 115
<input type="checkbox"/> 1007110 Desk Length - 3" x 10"	<input type="checkbox"/> 1007111 Full Length - 3" x 12"	
Gel Armpads May interfere with wheels in some tilt configurations.		\$ 125
<input type="checkbox"/> 1007120 Full Length - 2" x 12"	<input type="checkbox"/> 1007121 Full Length - 3.5" x 12"	
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.		<i>E2209</i> \$ 125
<input type="checkbox"/> 1007130 Right (4.5" x 18")	<input type="checkbox"/> 1007131 Left (4.5" x 18")	
<input type="checkbox"/> 1007220 Armrest Wideners for Dual Post Armrests	Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.	\$ 375
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.		
<input type="checkbox"/> 1007310 For Single Post Armrests Available for 10"-14" Height range only.		<i>K0108</i> \$ 110
<input type="checkbox"/> 1007320 For Dual Post Armrests		<i>K0108</i> \$ 165
Armrest Extensions Not compatible with Flip-up or Single Post Armrests.		\$ 60
<input type="checkbox"/> 1007400 3" Extension - Forward	<input type="checkbox"/> 1007401 3" Extension - Rearward	

Seat Pan		
<input type="checkbox"/> 1008000 No Seat Pan		
<input type="checkbox"/> 1008001 Solid Seat Pan		<i>E2231</i> \$ 120

Cushion		
<input type="checkbox"/> 1009000 No Cushion		
nxt BioFit		
<input type="checkbox"/> 1009101 14"-20" width, 16"-20" depth		<i>E2607</i> \$ 550
<input type="checkbox"/> 1009102 22"-24" width, 16"-20" depth		<i>E2608</i> \$ 600
Matrix VI		
<input type="checkbox"/> 1009111 14"-20" width, 16"-22" depth		<i>E2607</i> \$ 550
<input type="checkbox"/> 1009112 22" width, 18"-22" depth		<i>E2608</i> \$ 600
<input type="checkbox"/> 1009115 HD - 18"-30" width, 18"-22" depth		<i>E2608</i> \$ 1,150

See next page for more Cushion options >

ROHO Quadro Select High Profile				
<input type="checkbox"/>	1009121	14"-20" width, 16"-22" depth	<i>E2624</i>	\$ 650
<input type="checkbox"/>	1009122	22"-24" width, 16"-22" depth	<i>E2625</i>	\$ 750
ROHO High Profile				
<input type="checkbox"/>	1009125	14"-20" width, 16"-22" depth	<i>E2622</i>	\$ 650
<input type="checkbox"/>	1009126	22"-24" width, 16"-22" depth	<i>E2623</i>	\$ 750
Varilite Reflex				
<input type="checkbox"/>	1009131	14"-20" width, 16"-20" depth	<i>E2603</i>	\$ 375
Cushion Size				
1009400	Specify Size (Width" x Depth") _____			

Seat Height

Select desired Seat-to-Floor Height STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High				
	12"	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	15"	16"	17"	18"	19"
16"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
20"	13"	14"	15"			15"	16"	17"			15"	16"	17"	18"		
22"		14"	15"	16"			16"	17"	18"			16"	17"	18"	19"	
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"

8" Casters	RW	Low			High			Extra-High								
	16"			15"	16"	17"			17"	18"	19"			17"	18"	19"
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"

Wheel size _____
Caster size _____
STF height _____
Range _____

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1011061	6" x 1.25" Caster size		
<input type="checkbox"/>	1011063	6" x 1.50" Caster size		\$ 85
<input type="checkbox"/>	1011081	8" x 1.25" Caster size		
<input type="checkbox"/>	1011084	8" x 2.00" Caster size		\$ 135

Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1011161	6" x 1.25" Caster size	<i>E2214 & E2215</i>	\$ 110
<input type="checkbox"/>	1011181	8" x 1.25" Caster size	<i>E2214 & E2215</i>	\$ 120
<input type="checkbox"/>	1011184	8" x 2.00" Caster size	<i>E2214 & E2215</i>	\$ 150
<input type="checkbox"/>	1011260	6" Soft Roll		\$ 95
<input type="checkbox"/>	1011400	Ultra Sport Fork Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	<i>E1015</i>	\$ 395

Rear Wheels

Composite

<input type="checkbox"/>	1012012	12" Wheel size	<input type="checkbox"/>	1012016	16" Wheel size	<input type="checkbox"/>	1012020	20" Wheel size
<input type="checkbox"/>	1012022	22" Wheel size	<input type="checkbox"/>	1012024	24" Wheel size			

Spoke

<input type="checkbox"/>	1012120	20" Wheel size	<input type="checkbox"/>	1012122	22" Wheel size	<input type="checkbox"/>	1012124	24" Wheel size	\$ 75
<input type="checkbox"/>	1012224	24" Heavy Duty Spoke Required for Heavy Duty Base and Seat Frame							\$ 350

Axle Type

<input type="checkbox"/>	1013000	Quick Release Rear Axles
<input type="checkbox"/>	1013010	Fixed Rear Axles

Rear Tires			
<input type="checkbox"/>	1014000	Urethane	
<input type="checkbox"/>	1014010	Pneumatic Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/>	1014020	Pneumatic with Flat-Free Inserts Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.	E2211 & E2213 \$ 95

Hand Rims			
<input type="checkbox"/>	1015000	No Hand Rims	
<input type="checkbox"/>	1015010	Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1015020	Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 95

Wheel Locks			
<input type="checkbox"/>	1016010	Push to Lock	
<input type="checkbox"/>	1016020	Pull to Lock	
<input type="checkbox"/>	1016030	Attendant Wheel Lock (Push to Lock) Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.	\$ 85
<input type="checkbox"/>	1016040	User & Attendant Wheel Lock Available only on 20", 22", or 24" Rear Wheels.	\$ 125
<input type="checkbox"/>	1016100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 60
<input type="checkbox"/>	1016200	Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 350

Anti-Tippers			
<input checked="" type="checkbox"/>	1017010	Rear Anti-Tippers Required Safety Option.	E0971 \$ 60
<input type="checkbox"/>	1017020	Front Anti-Tippers	E0971 \$ 300
<input type="checkbox"/>	1017030	Heavy Duty Rear Anti-Tipper Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	E0971 \$ 600

Back Post Style			
<input type="checkbox"/>	1018010	Push Bar Back	
<input type="checkbox"/>	1018020	Cane Handle Back Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
Back Post Height Seat Pan to Handle			
<input type="checkbox"/>	1018425	25" Height	<input type="checkbox"/>
			<input type="checkbox"/> 1018431 31" Height

Back Options			
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150			
<input type="checkbox"/>	1019010	Button on Outside	<input type="checkbox"/>
			<input type="checkbox"/> 1019011 Button on Inside
Seat to Back Angle			
<input type="checkbox"/>	1019090	Set at 90°	<input type="checkbox"/>
			<input type="checkbox"/> 1019098 Set at 98° <input type="checkbox"/> 1019106 Set at 106°
<input type="checkbox"/>	1019114	Set at 114°	<input type="checkbox"/>
			<input type="checkbox"/> 1019122 Set at 122°
<input type="checkbox"/>	1019210	10" Backrest Extension Not compatible with Angle/Height Adjustable Push Bar Handle.	E0966 \$ 135
<input type="checkbox"/>	1019220	Fold Down Back Not compatible with Reclining Back.	
<input checked="" type="checkbox"/>	1019230	Spring Back Suspension Not compatible with Reclining Back.	
<input type="checkbox"/>	1019240	Dynamic Rocker Back Assembly Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height. Requires longer Frame.	K0108 \$ 2,995

Back Upholstery			
<input type="checkbox"/>	1020000	No Back Upholstery	
		Standard Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1020021	Set at 21"	<input type="checkbox"/>
<input type="checkbox"/>	1020022	Set at 22"	
		Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1020118	Set at 18"	<input type="checkbox"/>
<input type="checkbox"/>	1020119	Set at 19"	<input type="checkbox"/>
<input type="checkbox"/>	1020120	Set at 20"	
		Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1020224	Set at 24"	<input type="checkbox"/>
<input type="checkbox"/>	1020225	Set at 25"	<input type="checkbox"/>
<input type="checkbox"/>	1020226	Set at 26"	
		Adjustable Tension Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 295
<input type="checkbox"/>	1020321	Set at 21"	<input type="checkbox"/>
<input type="checkbox"/>	1020322	Set at 22"	<input type="checkbox"/>
<input type="checkbox"/>	1020323	Set at 23"	

Backrest			
<input type="checkbox"/>	1021000	No Backrest	
		nxt Optima HA Active Contour	
<input type="checkbox"/>	1021100	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1021105	22"-24" width, 16"-20" height	E2621 \$ 720
		nxt Optima HA Plus Contour	
<input type="checkbox"/>	1021101	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1021106	22"-24" width, 16"-20" height	E2621 \$ 720
		Matrx Elite	
<input type="checkbox"/>	1021110	14"-20" width, 10"-20" height	E2620 \$ 825
<input type="checkbox"/>	1021112	HD - 22"-26" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/>	1021114	HD - 26"-30" width, 16"-20" height	E2621 \$ 1,595
<input type="checkbox"/>	1021116	Deep - 14"-20" width, 12"-20" height	E2620 \$ 950
		ROHO Agility	
<input type="checkbox"/>	1021120	Min Contour 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1021122	Mid Contour 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1021124	Max Contour 14"-20" width, 10"-18" height	E2620 \$ 950
		Varilite Icon	
<input type="checkbox"/>	1021130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1021132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
		Backrest Size	
<input type="checkbox"/>	1021400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
		Attendant Tilt only	
<input type="checkbox"/>	1023011	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1023012	Trigger on Left Back Post	
		User Tilt only Not available on widths greater than 22". Full Length Arm Pads only.	\$ 100
<input type="checkbox"/>	1023021	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1023022	Trigger on Left Armrest	
		Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only.	\$ 195
<input type="checkbox"/>	1023031	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1023032	Trigger on Left Back Post	
<input type="checkbox"/>	1023041	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1023042	Trigger on Left Armrest	
<input type="checkbox"/>	1023110	Trigger Extension	\$ 50
<input type="checkbox"/>	1023410	Reclining Backrest (Manual)	E1225 \$ 875
		30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery. Requires longer Frame.	

Headrest			
<input type="checkbox"/>	1024000	No Headrest	
<input type="checkbox"/>	1024010	Headrest Mounting Plate Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 75
<input type="checkbox"/>	1024100	nxt MultiFit Head Support Hardware	E1028 \$ 155
		nxt MultiFit Head Support Pad	E0955 \$ 195
<input type="checkbox"/>	1024105	Regular (9" x 5")	
<input type="checkbox"/>	1024106	Large (10" x 5")	
<input type="checkbox"/>	1024107	X-Large (14" x 5")	
		Stealth Headrest Hardware	E1028 \$ 195
<input type="checkbox"/>	1024105	Multi-Axis HMO475	
<input type="checkbox"/>	1024105	Link TWB480	
		Stealth Headrest Pad	E0955 \$ 200
<input type="checkbox"/>	1024105	10" width	
<input type="checkbox"/>	1024106	14" width	

Front Rigging			
<input type="checkbox"/>	1025000	No Front Rigging	
<input type="checkbox"/>	1025001	Sector Blocks only	
<input type="checkbox"/>	1025110	Plug In 70°	
<input type="checkbox"/>	1025120	Plug In Adjustable Knee Angle	\$ 975
<input type="checkbox"/>	1025130	Swing Away 60°	
<input type="checkbox"/>	1025140	Swing Away 70°	
<input type="checkbox"/>	1025200	Swing Away Elevating Legrest	E0990 \$ 375
<input type="checkbox"/>	1025300	1" Offset Hangers for Swing Away	\$ 75

Foot Plates			
<input type="checkbox"/>	1026010	Composite	
<input type="checkbox"/>	1026020	Aluminum Angle Adjustable	K0040 \$ 155
<input type="checkbox"/>	1026030	One Piece Fixed Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1026040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1026050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.	K0040 \$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1027010	Swing Away 60° & 70° Calf Pad Not compatible with 60° & 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1027020	Swing Away 60° & 70° Gel Padding Not compatible with 60° & 70° Swing Away Calf Pad.	\$ 100
<input type="checkbox"/>	1027030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1027040	Adjustable Length Heel Loops	E0951 \$ 50
<input type="checkbox"/>	1027050	Bodypoint - Padded Calf Strap Sizes greater than 22" are not padded.	K0038 \$ 105
		Miller's Dynamic Footrest Mounts	
<input type="checkbox"/>	1027060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging Not available on 13" STF Height.	\$ 380

Extension Tubes			
		Regular Extension Tubes	
<input type="checkbox"/>	1028010	Standard (13"-17")	
<input type="checkbox"/>	1028020	Long (15"-19")	\$ 75
		Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only.	\$ 100
<input type="checkbox"/>	1028110	Standard (13"-17")	
<input type="checkbox"/>	1028120	Long (15"-19")	

Positioning Straps

Bodypoint - Hip Belts



<input type="checkbox"/>	1029010	Center-Pull Auto Style	<i>E0978</i>	\$ 105
<input type="checkbox"/>	1029015	Center-Pull Auto Style, Padded	<i>E0978</i>	\$ 155
<input type="checkbox"/>	1029020	4-Point, Padded	<i>E0978</i>	\$ 175
<input type="checkbox"/>	1029030	Bodypoint - Chest Strap	<i>E0960</i>	\$ 85

Miscellaneous Options

One-Arm Drive Maximum Seat Width of 22". 24" Composite Wheels required. 18" STF minimum. Limits tilt.

E0958 \$ 3,500

<input type="checkbox"/>	1030010	Installed on Right side	<input type="checkbox"/>	1030011	Installed on Left side	
<input type="checkbox"/>	1030050	Protective Lower Padding				\$ 425
<input type="checkbox"/>	1030080	Gel Hip Guides			<i>E0956 & E1028</i>	\$ 325
<input type="checkbox"/>	1030110	IV Pole <small>Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.</small>			<i>E1029</i>	\$ 295
<input type="checkbox"/>	1030120	O2 Holder <small>Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.</small>			<i>E2208</i>	\$ 395
<input type="checkbox"/>	1030170	Adhesive Back Hook Strip, 2" x 30" <small>For Cushion mounting.</small>				\$ 10

Order Notes

HCPSC codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPSC codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/bentley to learn more about the Bentley 20° Manual Tilt Wheelchair



PDG Mobility