

Customer and Order Information		
Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type		
<input checked="" type="checkbox"/> 1001000	<b>Bentley 20° Manual Tilt Wheelchair</b>	E1161 \$ 2,995

Frame Options		
<input type="checkbox"/> 1002000	<b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.	
<input type="checkbox"/> 1002050	<b>Upgraded Base and Seat Frame</b> 300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	K0108 \$ 200
<input type="checkbox"/> 1002100	<b>Heavy Duty Base and Seat Frame</b> 450 lbs Weight Capacity. Available for all sizes, but required for Seat Width of 23" or wider. Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	K0108 \$ 300
<input type="checkbox"/> 1002200	<b>Extended Base with Rear Floor Jacks</b> Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Locks, HD Anti-Tipper, IV Pole, O2 Holder.	K0108 \$ 600
<input type="checkbox"/> 1002300	<b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW
<b>Standard Seat Widths</b>		
<input type="checkbox"/> 1003014 14" Seat Width	<input type="checkbox"/> 1003015 15" Seat Width	<input type="checkbox"/> 1003016 16" Seat Width
<input type="checkbox"/> 1003017 17" Seat Width	<input type="checkbox"/> 1003018 18" Seat Width	<input type="checkbox"/> 1003019 19" Seat Width
<input type="checkbox"/> 1003020 20" Seat Width	<input type="checkbox"/> 1003021 21" Seat Width	<input type="checkbox"/> 1003022 22" Seat Width
		E2201 \$ 325
<b>Custom Seat Widths</b> Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1003023 23" Seat Width		E2201 \$ 325
<input type="checkbox"/> 1003024 24" Seat Width	<input type="checkbox"/> 1003025 25" Seat Width	<input type="checkbox"/> 1003026 26" Seat Width
<input type="checkbox"/> 1003027 27" Seat Width		E2202 \$ 525
<input type="checkbox"/> 1003028 28" Seat Width	<input type="checkbox"/> 1003029 29" Seat Width	<input type="checkbox"/> 1003030 30" Seat Width
<input type="checkbox"/> 1003031 31" Seat Width	<input type="checkbox"/> 1003032 32" Seat Width	K0108 \$ 925
		K0108 \$ 925

Seat Depth		
<b>Standard Seat Depths (Adjustable)</b>		
<input type="checkbox"/> 1004016 16" Seat Depth	<input type="checkbox"/> 1004017 17" Seat Depth	<input type="checkbox"/> 1004018 18" Seat Depth
<input type="checkbox"/> 1004019 19" Seat Depth		
<input type="checkbox"/> 1004020 20" Seat Depth		E2203 \$ 475
<b>Custom Seat Depths</b>		
<input type="checkbox"/> 1004021 21" Seat Depth		E2203 \$ 475
<input type="checkbox"/> 1004022 22" Seat Depth		E2204 \$ 750

Frame Color		
<b>PDG Premium Colors</b>		
<input type="checkbox"/> 1005000 Tracer Red	<input type="checkbox"/> 1005001 Gloss Black	<input type="checkbox"/> 1005002 Silver
<input type="checkbox"/> 1005003 Black Shimmer	<input type="checkbox"/> 1005004 Orange Fusion	<input type="checkbox"/> 1005005 Karl's Imagination

See next page for more Frame Color options >

<input type="checkbox"/> 1005006 Blue Shimmer	<input type="checkbox"/> 1005007 Gloss White	<input type="checkbox"/> 1005008 Dark Titanium
<input type="checkbox"/> 1005009 Geyser Green	<input type="checkbox"/> 1005010 Lollypop Blue	<input type="checkbox"/> 1005011 Lollypop Red
<input type="checkbox"/> 1005012 Bubble Gum Pink	<input type="checkbox"/> 1005013 Berry Purple	<input type="checkbox"/> 1005014 Safety Yellow
<input type="checkbox"/> 1005200 <b>Custom Color</b> Contact PDG for specifications.		\$ 300

### Arm Type

<input type="checkbox"/> 1006000 <b>No Armrest</b>			
<input type="checkbox"/> 1006100 <b>Single Post Fixed Height (12")</b>	Full Length Armpad only.		
<b>Single Post Adjustable Height</b>		Not compatible with the Heavy Duty Base and/or Seat Width greater than 22".	E0973 \$ 255
<input type="checkbox"/> 1006111 8"-11"	<input type="checkbox"/> 1006114 10"-14"		
<input type="checkbox"/> 1006200 <b>Dual Post Fixed Height (12")</b>	Full Length Armpad only.		
<b>Dual Post Adjustable Height</b>			E0973 \$ 255
<input type="checkbox"/> 1006211 8"-11"	<input type="checkbox"/> 1006214 10"-14"	<input type="checkbox"/> 1006216 12"-16"	
<b>Locking Flip-up</b>		Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.	E0973 \$ 255
<input type="checkbox"/> 1006310 Set at 10"	<input type="checkbox"/> 1006311 Set at 11"	<input type="checkbox"/> 1006312 Set at 12"	
<input type="checkbox"/> 1006313 Set at 13"	<input type="checkbox"/> 1006314 Set at 14"		

### Armrests

<b>Regular Armpads</b>			
<input type="checkbox"/> 1007100 Desk Length - 2" x 10"	<input type="checkbox"/> 1007101 Full Length - 2" x 14"		
<input type="checkbox"/> 1007102 Full Length Premium Grade - 2" x 14"			\$ 35
<b>Wide Armpads</b>		May interfere with wheels in some tilt configurations.	\$ 115
<input type="checkbox"/> 1007110 Desk Length - 3" x 10"	<input type="checkbox"/> 1007111 Full Length - 3" x 12"		
<b>Gel Armpads</b>		May interfere with wheels in some tilt configurations.	\$ 125
<input type="checkbox"/> 1007120 Full Length - 2" x 12"	<input type="checkbox"/> 1007121 Full Length - 3.5" x 12"		
<b>Contour Arm Trough with Hand Pad</b>		May interfere with wheels in some tilt configurations.	E2209 \$ 125
<input type="checkbox"/> 1007130 Right (4.5" x 18")	<input type="checkbox"/> 1007131 Left (4.5" x 18")		
<input type="checkbox"/> 1007220 <b>Armrest Wideners for Dual Post Armrests</b>	Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.		\$ 375
<b>Clothing Guards</b>		Not compatible with Flip-up Armrests, Armrest Wideners.	
<input type="checkbox"/> 1007310 For Single Post Armrests	Available for 10"-14" Height range only.		K0108 \$ 110
<input type="checkbox"/> 1007320 For Dual Post Armrests			K0108 \$ 165
<b>Armrest Extensions</b>		Not compatible with Flip-up or Single Post Armrests.	\$ 60
<input type="checkbox"/> 1007400 3" Extension - Forward	<input type="checkbox"/> 1007401 3" Extension - Rearward		

### Seat Pan

<input type="checkbox"/> 1008000 <b>No Seat Pan</b>		
<input type="checkbox"/> 1008001 <b>Solid Aluminum Seat Pan</b>		E2231 \$ 120

### Cushion

<input type="checkbox"/> 1009000 <b>No Cushion</b>		
<b>next BioFit</b>		
<input type="checkbox"/> 1009101 14"-20" width, 16"-20" depth		E2607 \$ 550
<input type="checkbox"/> 1009102 22"-24" width, 16"-20" depth		E2608 \$ 600
<b>Matrix VI</b>		
<input type="checkbox"/> 1009111 14"-20" width, 16"-22" depth		E2607 \$ 550
<input type="checkbox"/> 1009112 22" width, 18"-22" depth		E2608 \$ 600
<input type="checkbox"/> 1009115 HD - 18"-30" width, 18"-22" depth		E2608 \$ 1,150

See next page for more Cushion options >

<b>ROHO Quadro Select High Profile</b>				
<input type="checkbox"/>	1009121	14"-20" width, 16"-22" depth	E2624	\$ 650
<input type="checkbox"/>	1009122	22"-24" width, 16"-22" depth	E2625	\$ 750
<b>ROHO High Profile</b>				
<input type="checkbox"/>	1009125	14"-20" width, 16"-22" depth	E2622	\$ 650
<input type="checkbox"/>	1009126	22"-24" width, 16"-22" depth	E2623	\$ 750
<b>Varilite Reflex</b>				
<input type="checkbox"/>	1009131	14"-20" width, 16"-20" depth	E2603	\$ 375
<b>Cushion Size</b>				
1009400	Specify Size (Width" x Depth") _____			



Seat Height																	
Select desired Seat-to-Floor Height																	
STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.																	
6" Casters	RW	Low					High					Extra-High					
	12"	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	15"	16"	17"	18"	19"	20"
	16"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
	20"	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	15"	16"	17"	18"	19"	20"
	22"		14"	15"	16"	17"			16"	17"	18"			16"	17"	18"	19"
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"	
8" Casters	RW	Low					High					Extra-High					
	16"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
	24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
Wheel size	_____																
Caster size	_____																
STF height	_____																
Range	_____																

Casters				
<b>Urethane</b> 8" Casters are compatible only with 16" and 24" Rear Wheels.				
<input type="checkbox"/>	1011061	6" x 1.25" Caster size		
<input type="checkbox"/>	1011063	6" x 1.50" Caster size		\$ 85
<input type="checkbox"/>	1011081	8" x 1.25" Caster size		
<input type="checkbox"/>	1011084	8" x 2.00" Caster size		\$ 135
<b>Pneumatic</b> 8" Casters are compatible only with 16" and 24" Rear Wheels.				
<input type="checkbox"/>	1011161	6" x 1.25" Caster size	E2214 & E2215	\$ 110
<input type="checkbox"/>	1011181	8" x 1.25" Caster size	E2214 & E2215	\$ 120
<input type="checkbox"/>	1011184	8" x 2.00" Caster size	E2214 & E2215	\$ 150
<input type="checkbox"/>	1011260	<b>6" Soft Roll</b>		\$ 95
<input type="checkbox"/>	1011400	<b>Ultra Sport Fork</b> Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	E1015	\$ 395

Rear Wheels				
<b>Composite</b>				
<input type="checkbox"/>	1012012	12" Wheel size	<input type="checkbox"/>	1012016 16" Wheel size
<input type="checkbox"/>	1012022	22" Wheel size	<input type="checkbox"/>	1012024 24" Wheel size
<input type="checkbox"/>			<input type="checkbox"/>	1012020 20" Wheel size
<b>Spoke</b> \$ 75				
<input type="checkbox"/>	1012120	20" Wheel size	<input type="checkbox"/>	1012122 22" Wheel size
<input type="checkbox"/>			<input type="checkbox"/>	1012124 24" Wheel size
<input type="checkbox"/>	1012224	<b>24" Heavy Duty Spoke</b> Required for Heavy Duty Base and Seat Frame.		\$ 350

Axle Type	
<input type="checkbox"/>	1013000 <b>Quick Release Rear Axles</b>
<input type="checkbox"/>	1013010 <b>Fixed Rear Axles</b>

Rear Tires			
<input type="checkbox"/>	1014000	<b>Urethane</b>	
<input type="checkbox"/>	1014010	<b>Pneumatic</b> Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/>	1014020	<b>Pneumatic with Flat-Free Inserts</b> Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.	E2211 & E2213 \$ 95

Hand Rims			
<input type="checkbox"/>	1015000	<b>No Hand Rims</b>	
<input type="checkbox"/>	1015010	<b>Standard</b> 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1015020	<b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only.	\$ 135

Wheel Locks			
<input type="checkbox"/>	1016010	<b>Push to Lock</b>	
<input type="checkbox"/>	1016020	<b>Pull to Lock</b>	
<input type="checkbox"/>	1016030	<b>Attendant Wheel Lock (Push to Lock)</b> Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.	\$ 85
<input type="checkbox"/>	1016040	<b>User &amp; Attendant Wheel Lock</b> Available only on 20", 22", or 24" Rear Wheels.	\$ 125
<input type="checkbox"/>	1016100	<b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	E0961 \$ 60
<input type="checkbox"/>	1016200	<b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 350

Anti-Tippers			
<input checked="" type="checkbox"/>	1017010	<b>Rear Anti-Tippers</b> Required Safety Option.	E0971 \$ 60
<input type="checkbox"/>	1017020	<b>Front Anti-Tippers</b>	E0971 \$ 300
<input type="checkbox"/>	1017030	<b>Heavy Duty Rear Anti-Tipper</b> Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	E0971 \$ 600

Back Post Style			
<input type="checkbox"/>	1018010	<b>Push Bar Back</b>	
<input type="checkbox"/>	1018020	<b>Cane Handle Back</b> Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
<b>Back Post Height</b> Seat Pan to Handle			
<input type="checkbox"/>	1018425	25" Height	<input type="checkbox"/> 1018431 31" Height

Back Options			
<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150			
<input type="checkbox"/>	1019010	Button on Outside	<input type="checkbox"/> 1019011 Button on Inside
<b>Seat to Back Angle</b>			
<input type="checkbox"/>	1019090	Set at 90°	<input type="checkbox"/> 1019098 Set at 98° <input type="checkbox"/> 1019106 Set at 106°
<input type="checkbox"/>	1019114	Set at 114°	<input type="checkbox"/> 1019122 Set at 122°
<input type="checkbox"/>	1019210	<b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Push Bar Handle.	E0966 \$ 135
<input type="checkbox"/>	1019220	<b>Fold Down Back</b> Not compatible with Reclining Back.	
<input checked="" type="checkbox"/>	1019230	<b>Spring Back Suspension</b> Not compatible with Reclining Back.	
<input type="checkbox"/>	1019240	<b>Dynamic Rocker Back Assembly</b> Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height. Requires longer Frame.	K0108 \$ 2,995

Back Upholstery			
<input type="checkbox"/>	1020000	<b>No Back Upholstery</b>	
		<b>Standard</b> Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1020021	Set at 21"	<input type="checkbox"/>
<input type="checkbox"/>	1020022	Set at 22"	Not available with 25" Back.
		<b>Short</b> 25" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1020118	Set at 18"	<input type="checkbox"/>
<input type="checkbox"/>	1020119	Set at 19"	<input type="checkbox"/>
<input type="checkbox"/>	1020120	Set at 20"	
		<b>Tall</b> 31" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1020224	Set at 24"	<input type="checkbox"/>
<input type="checkbox"/>	1020225	Set at 25"	<input type="checkbox"/>
<input type="checkbox"/>	1020226	Set at 26"	
		<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 295
<input type="checkbox"/>	1020321	Set at 21"	<input type="checkbox"/>
<input type="checkbox"/>	1020322	Set at 22"	<input type="checkbox"/>
<input type="checkbox"/>	1020323	Set at 23"	

Backrest			
<input type="checkbox"/>	1021000	<b>No Backrest</b>	
		<b>nxt Optima HA Active Contour</b>	
<input type="checkbox"/>	1021100	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1021105	22"-24" width, 16"-20" height	E2621 \$ 720
		<b>nxt Optima HA Plus Contour</b>	
<input type="checkbox"/>	1021101	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1021106	22"-24" width, 16"-20" height	E2621 \$ 720
		<b>Matrx Elite</b>	
<input type="checkbox"/>	1021110	14"-20" width, 10"-20" height	E2620 \$ 825
<input type="checkbox"/>	1021112	HD - 22"-26" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/>	1021114	HD - 26"-30" width, 16"-20" height	E2621 \$ 1,595
<input type="checkbox"/>	1021116	Deep - 14"-20" width, 12"-20" height	E2620 \$ 950
		<b>ROHO Agility</b>	
<input type="checkbox"/>	1021120	Min Contour 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1021122	Mid Contour 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1021124	Max Contour 14"-20" width, 10"-18" height	E2620 \$ 950
		<b>Varilite Icon</b>	
<input type="checkbox"/>	1021130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1021132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
		<b>Backrest Size</b>	
<input type="checkbox"/>	1021400	Specify Size (Width" x Height")	_____

Tilt/Recline Options			
		<b>Attendant Tilt only</b>	
<input type="checkbox"/>	1023011	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1023012	Trigger on Left Back Post	
		<b>User Tilt only</b> Not available on widths greater than 22". Full Length Arm Pads only.	\$ 100
<input type="checkbox"/>	1023021	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1023022	Trigger on Left Armrest	
		<b>Attendant &amp; User Tilt</b> Not available on widths greater than 22". Full Length Arm Pads only.	\$ 195
<input type="checkbox"/>	1023031	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1023032	Trigger on Left Back Post	
<input type="checkbox"/>	1023041	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1023042	Trigger on Left Armrest	
<input type="checkbox"/>	1023110	<b>Trigger Extension</b>	\$ 50
<input type="checkbox"/>	1023410	<b>Reclining Backrest (Manual)</b>	E1225 \$ 875
		30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19".	

Headrest			
<input type="checkbox"/>	1024000	<b>No Headrest</b>	
<input type="checkbox"/>	1024010	<b>Headrest Mounting Plate</b> Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 75
<input type="checkbox"/>	1024100	<b>nxt MultiFit Head Support Hardware</b>	E1028 \$ 235
		<b>nxt MultiFit Head Support Pad</b>	E0955 \$ 245
<input type="checkbox"/>	1024105	Regular (9" x 5")	
<input type="checkbox"/>	1024106	Large (10" x 5")	
<input type="checkbox"/>	1024107	X-Large (14" x 5")	
		<b>Stealth Headrest Hardware</b>	
<input type="checkbox"/>	1024110	Multi-Axis HMO475	E1028 \$ 150
<input type="checkbox"/>	1024111	Link TWB480	E1028 \$ 325
		<b>Stealth Headrest Pad</b>	E0955 \$ 245
<input type="checkbox"/>	1024115	10" width	
<input type="checkbox"/>	1024116	14" width	

Front Rigging			
<input type="checkbox"/>	1025000	<b>No Front Rigging</b>	
		<b>PDG Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame.	
<input type="checkbox"/>	1025410	Fixed 70°	Available in March \$ 95
<input type="checkbox"/>	1025420	Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____	Available in April \$ 195
<input type="checkbox"/>	1025430	Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.	Available in April E0990 \$ 495
		<b>Standard Front Rigging</b> Legacy Design. Not compatible with Easy-Out™.	
<input type="checkbox"/>	1025110	Plug In 70°	
<input type="checkbox"/>	1025130	Swing Away 60°	
<input type="checkbox"/>	1025140	Swing Away 70°	
<input type="checkbox"/>	1025200	Swing Away Elevating Legrests	E0990 \$ 375
<input type="checkbox"/>	1025300	1" Offset Hangers for Swing Away	\$ 75

Foot Plates			
<input type="checkbox"/>	1026010	<b>Composite</b>	
<input type="checkbox"/>	1026020	<b>Aluminum Angle Adjustable</b>	K0040 \$ 155
<input type="checkbox"/>	1026030	<b>One Piece Fixed</b> Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1026040	<b>One Piece Flip-up</b> Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1026050	<b>Interlocking Angle Adjustable</b> Seat Widths of 14"-21" only.	K0040 \$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1027010	<b>Swing Away 60° &amp; 70° Calf Pad</b> Not compatible with 60° & 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1027020	<b>Swing Away 60° &amp; 70° Gel Padding</b> Not compatible with 60° & 70° Swing Away Calf Pad, Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1027030	<b>Swing Away ELR Gel Knee Pad</b>	\$ 100
<input type="checkbox"/>	1027040	<b>Adjustable Length Heel Loops</b>	E0951 \$ 50
<input type="checkbox"/>	1027050	<b>Bodypoint - Padded Calf Strap</b> Sizes greater than 22" are not padded.	K0038 \$ 105
		<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.	
<input type="checkbox"/>	1027060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging Not available on 13" STF Height.	\$ 380

Extension Tubes			
		<b>Regular Extension Tubes</b>	
<input type="checkbox"/>	1028010	Short 13"-17" for Easy-Out SA Angle Adjustable 10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1028020	Long 16"-22.5" for Easy-Out SA Angle Adjustable 12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
		<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.	\$ 100
<input type="checkbox"/>	1028110	Standard (13"-17")	
<input type="checkbox"/>	1028120	Long (15"-19")	

### Positioning Straps

#### Bodypoint - Hip Belts



<input type="checkbox"/>	1029010	Center-Pull Auto Style	E0978	\$	105
<input type="checkbox"/>	1029015	Center-Pull Auto Style, Padded	E0978	\$	155
<input type="checkbox"/>	1029020	4-Point, Padded	E0978	\$	175
<input type="checkbox"/>	1029030	<b>Bodypoint - Chest Strap</b>	E0960	\$	85

### Miscellaneous Options

**One-Arm Drive** Maximum Seat Width of 22". 24" Composite Wheels with Urethane Tires required. 18" STF minimum. Limits tilt. E0958 \$ 3,500  
 Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.

<input type="checkbox"/>	1030010	Installed on Right side	<input type="checkbox"/>	1030011	Installed on Left side	
<input type="checkbox"/>	1030050	<b>Protective Lower Padding</b>				\$ 425
<input type="checkbox"/>	1030080	<b>Gel Hip Guides</b>			E0956 & E1028	\$ 325
<input type="checkbox"/>	1030110	<b>IV Pole</b> Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			E1029	\$ 295
<input type="checkbox"/>	1030120	<b>O2 Holder</b> Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			E2208	\$ 395
<input type="checkbox"/>	1030170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.				\$ 10

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.  
 Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/bentley/](http://pdgmobility.com/products/bentley/) to learn more about the **Bentley 20° Manual Tilt Wheelchair**

