

Customer and Order Information

Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address *	<i>*required</i>

Frame Type

<input checked="" type="checkbox"/> 1101000	Eclipse Bariatric Manual Wheelchair	Extended delivery lead times may apply.	<i>K0007</i>	\$ 2,795
---	--	---	--------------	----------

Frame Options

Heavy Duty Base and Seat Frame				
<input type="checkbox"/> 1102000	600 lbs Weight Capacity			
<input type="checkbox"/> 1102100	1000 lbs Weight Capacity		<i>K0108</i>	\$ 300

Seat Width Overall chair width is +8"

Seat Widths Seat Width is adjustable within the original range with the optional Adjustment Kit (sold separately below).				
<input type="checkbox"/> 1103001	20"-26" Range, Width set to _____		<i>E2201 / E2202</i>	\$ 625
<input type="checkbox"/> 1103002	26"-32" Range, Width set to _____		<i>E2202 / K0108</i>	\$ 925
<input type="checkbox"/> 1103003	30"-36" Range, Width set to _____		<i>K0108</i>	\$ 1,125
<input type="checkbox"/> 1103004	36"-42" Range, Width set to _____		<i>K0108</i>	\$ 1,125
<input type="checkbox"/> 1103005	42"-48" Range, Width set to _____		<i>K0108</i>	\$ 1,225
Seat Widths greater than 48" are available for Custom Chairs				

Seat Width Adjustment Kit

Optional Seat Width Adjustment Kit				<i>K0108</i>	\$ 750
This optional Kit is required to grow or shrink the Chair in the future (Width must remain within the original range). Kit includes Seat Pan, Back Posts, Back Upholstery, and Width Adjustment Spacers.					
<input type="checkbox"/> 1131001	20"-26" Range, Width set to _____				
<input type="checkbox"/> 1131002	26"-32" Range, Width set to _____				
<input type="checkbox"/> 1131003	30"-36" Range, Width set to _____				
<input type="checkbox"/> 1131004	36"-42" Range, Width set to _____				
<input type="checkbox"/> 1131005	42"-48" Range, Width set to _____				

Seat Depth

Seat Depths (Adjustable)					
<input type="checkbox"/> 1104016	16" Seat Depth	<input type="checkbox"/> 1104017	17" Seat Depth		\$ 325
<input type="checkbox"/> 1104018	18" Seat Depth	<input type="checkbox"/> 1104019	19" Seat Depth		\$ 325
<input type="checkbox"/> 1104020	20" Seat Depth	<input type="checkbox"/> 1104021	21" Seat Depth	<i>E2203</i>	\$ 475
<input type="checkbox"/> 1104022	22" Seat Depth			<i>E2204</i>	\$ 750
Custom Seat Depths					
<input type="checkbox"/> 1104023	23" Seat Depth	<input type="checkbox"/> 1104024	24" Seat Depth	<i>E2204</i>	\$ 750

Frame Color

PDG Premium Colors

<input type="checkbox"/> 1105000 Tracer Red	<input type="checkbox"/> 1105001 Gloss Black	<input type="checkbox"/> 1105002 Silver
<input type="checkbox"/> 1105003 Black Shimmer	<input type="checkbox"/> 1105004 Orange Fusion	<input type="checkbox"/> 1105005 Karl's Imagination
<input type="checkbox"/> 1105006 Blue Shimmer	<input type="checkbox"/> 1105007 Gloss White	<input type="checkbox"/> 1105008 Dark Titanium
<input type="checkbox"/> 1105009 Geyser Green	<input type="checkbox"/> 1105010 Lollypop Blue	<input type="checkbox"/> 1105011 Lollypop Red
<input type="checkbox"/> 1105012 Bubble Gum Pink	<input type="checkbox"/> 1105013 Berry Purple	<input type="checkbox"/> 1105014 Safety Yellow
<input type="checkbox"/> 1105200 Custom Color Contact PDG for specifications. \$ 300		

Arm Type

<input type="checkbox"/> 1106000 No Armrest		
<input type="checkbox"/> 1106200 Dual Post Fixed Height (12")		
Dual Post Adjustable Height		<i>E0973</i> \$ 255
<input type="checkbox"/> 1106214 10"-14"	<input type="checkbox"/> 1106216 12"-16"	

Armrests

Regular Armpads

<input type="checkbox"/> 1107100 Desk Length - 2" x 10"	<input type="checkbox"/> 1107101 Full Length - 2" x 14"	
<input type="checkbox"/> 1107102 Full Length Premium Grade - 2" x 14"		\$ 35

Wide Armpads

May interfere with wheels in some tilt configurations.

\$ 115

<input type="checkbox"/> 1107110 Desk Length - 3" x 10"	<input type="checkbox"/> 1107111 Full Length - 3" x 12"	
---	---	--

Gel Armpads

May interfere with wheels in some tilt configurations.

\$ 125

<input type="checkbox"/> 1107120 Full Length - 2" x 12"	<input type="checkbox"/> 1107121 Full Length - 3.5" x 12"	
---	---	--

<input type="checkbox"/> 1107250 Forward Guards		\$ 200
--	--	--------

<input type="checkbox"/> 1107400 3" Armrest Extension - Forward		\$ 60
--	--	-------

Seat Pan

<input checked="" type="checkbox"/> 1108001 Solid Seat Pan		<i>E2231</i> \$ 120
---	--	---------------------

Cushion

<input type="checkbox"/> 1109000 No Cushion		
--	--	--

Matrix VI Heavy Duty



<input type="checkbox"/> 1109115 22"-30" width, 18"-22" depth		<i>E2608</i> \$ 1,150
---	--	-----------------------

3" Basic Foam (Nylon Cover)

<input type="checkbox"/> 1109200 Up to 28" Width		<i>E2602</i> \$ 485
--	--	---------------------

<input type="checkbox"/> 1109201 Above 28" Width		<i>E2602</i> \$ 605
--	--	---------------------

3" Contour Foam (Nylon Cover)

<input type="checkbox"/> 1109210 Up to 28" Width		<i>E2602</i> \$ 835
--	--	---------------------

<input type="checkbox"/> 1109211 Above 28" Width		<i>E2602</i> \$ 1,090
--	--	-----------------------

3" Contour Pressure Reduction (Nylon Cover)

<input type="checkbox"/> 1109220 Up to 28" Width		<i>E2602</i> \$ 1,485
--	--	-----------------------

<input type="checkbox"/> 1109221 Above 28" Width		<i>E2602</i> \$ 1,900
--	--	-----------------------

Spare Cushion Cover

<input type="checkbox"/> 1109230 Basic		\$ 390
--	--	--------

<input type="checkbox"/> 1109231 Contour		\$ 490
--	--	--------

Cushion Size

2109400 Specify Size (Width" x Depth") _____

Seat Height

Front Seat Height

Angle Adjustable Foot Plates required for 14.5" and 15".

- | | | | | | | |
|----------------------------------|-----------------|----------------------------------|---------------|----------------------------------|---------------|--------|
| <input type="checkbox"/> 1110014 | 14.5" STF Front | <input type="checkbox"/> 1110015 | 15" STF Front | <input type="checkbox"/> 1110016 | 16" STF Front | |
| <input type="checkbox"/> 1110017 | 17" STF Front | <input type="checkbox"/> 1110018 | 18" STF Front | | | |
| <input type="checkbox"/> 1110019 | 19" STF Front | <input type="checkbox"/> 1110020 | 20" STF Front | | | \$ 200 |

Rear Seat Height

Angle Adjustable Foot Plates required for 14.5" and 15".

- | | | | | | | |
|----------------------------------|----------------|----------------------------------|--------------|----------------------------------|--------------|--------|
| <input type="checkbox"/> 1110114 | 14.5" STF Rear | <input type="checkbox"/> 1110115 | 15" STF Rear | <input type="checkbox"/> 1110116 | 16" STF Rear | |
| <input type="checkbox"/> 1110117 | 17" STF Rear | <input type="checkbox"/> 1110118 | 18" STF Rear | | | |
| <input type="checkbox"/> 1110119 | 19" STF Rear | <input type="checkbox"/> 1110120 | 20" STF Rear | | | \$ 200 |

Casters

- 1111050 **5" Caster With Heavy Duty Steel Fork**

Rear Wheels

- 1112224 **24" Heavy Duty Spoke**

Axle Type

- 1113000 **Quick Release Rear Axles**
- 1113010 **Fixed Rear Axles**

Rear Axle Position

- | | | | |
|----------------------------------|-----------------|----------------------------------|------------------|
| <input type="checkbox"/> 1113101 | Set at 1 (Rear) | <input type="checkbox"/> 1113102 | Set at 2 |
| <input type="checkbox"/> 1113103 | Set at 3 | <input type="checkbox"/> 1113104 | Set at 4 (Front) |

Rear Tires

- 1114000 **Urethane**

Hand Rims

- 1115000 **No Hand Rims**
- 1115010 **Standard**
- 1115020 **Plastic Coated** \$ 250

Wheel Locks

- 1116010 **Push to Lock**
- 1116020 **Pull to Lock**
- 1116100 **Wheel Lock Extensions** *E0961* \$ 60

Anti-Tippers

- 1117010 **Rear Anti-Tippers** Required Safety Option. *E0971* \$ 60

Back Post Style

- 1118030 **Push Bar Back, 8° Bend**
- 1118040 **Cane Handle Back, 8° Bend** Not compatible with Headrest Mounting Plate.

Back Post Height

Seat Pan to Handle.

- | | | | |
|----------------------------------|------------|----------------------------------|------------|
| <input type="checkbox"/> 1118425 | 25" Height | <input type="checkbox"/> 1118430 | 30" Height |
|----------------------------------|------------|----------------------------------|------------|

Back Options

Seat to Back Angle		
<input type="checkbox"/> 1119085	Set at 85°	<input type="checkbox"/> 1119090 Set at 90°
<input type="checkbox"/> 1119095	Set at 95°	<input type="checkbox"/> 1119100 Set at 100°
<input type="checkbox"/> 1119210	10" Backrest Extension	E0966 \$ 135

Back Upholstery

<input type="checkbox"/> 1120000	No Back Upholstery		
Standard Back height measured from Seat Pan to top of Upholstery		E2612	\$ 250
<input type="checkbox"/> 1120018	Set at 18"	<input type="checkbox"/> 1120019 Set at 19"	<input type="checkbox"/> 1120020 Set at 20"
Short (25" Back only) Back height measured from Seat Pan to top of Upholstery		E2612	\$ 250
<input type="checkbox"/> 1120116	Set at 16"	<input type="checkbox"/> 1120117 Set at 17"	
Tall (30" Back only) Back height measured from Seat Pan to top of Upholstery		E2612	\$ 250
<input type="checkbox"/> 1120221	Set at 21"	<input type="checkbox"/> 1120222 Set at 22"	<input type="checkbox"/> 1120223 Set at 23"

Backrest

<input type="checkbox"/> 1121000	No Backrest		
Matrix Elite			
<input type="checkbox"/> 1121112	HD - 22"-26" width, 16"-20" height		E2621 \$ 1,095
<input type="checkbox"/> 1121114	HD - 26"-30" width, 16"-20" height		E2621 \$ 1,595
Backrest Size			
1021400	Specify Size (Width" x Height")	_____	

Tilt/Recline Options

<input type="checkbox"/> 1123410	Reclining Backrest (Manual)	E1225	\$ 875
30° maximum recline. Not compatible with Backrest Extension. Reclining Backrest is 1" narrower than Seat Width. Back Frame is not removable.			

Headrest

<input type="checkbox"/> 1124000	No Headrest		
<input type="checkbox"/> 1124010	Headrest Mounting Plate Not required if mounting to a solid Back.		\$ 75
<input type="checkbox"/> 1124100	nxt MultiFit Head Support Hardware		E1028 \$ 155
nxt MultiFit Head Support Pad			E0955 \$ 195
<input type="checkbox"/> 1124105	Regular (9" x 5")	<input type="checkbox"/> 1124106 Large (10" x 5")	<input type="checkbox"/> 1124107 X-Large (14" x 5")
Stealth Headrest Hardware			E1028 \$ 195
<input type="checkbox"/> 1124105	Multi-Axis HMO475		
<input type="checkbox"/> 1124105	Link TWB480		
Stealth Headrest Pad			E0955 \$ 200
<input type="checkbox"/> 1124105	10" width	<input type="checkbox"/> 1124106 14" width	

Front Rigging

<input type="checkbox"/> 1125001	Sector Blocks only		
<input type="checkbox"/> 1125002	Center Mount Front Rigging Bracket Compatible with 60° only.		\$ 300
<input type="checkbox"/> 1125130	Swing Away 60°		
<input type="checkbox"/> 1125200	Swing Away Elevating Legrest	E0990	\$ 375

Foot Plates

<input type="checkbox"/>	1126010	Composite	Not compatible with 14.5" and 15" STF Height.		
<input type="checkbox"/>	1126020	Aluminum Angle Adjustable	Required for 14.5" and 15" STF Height.	K0040	\$ 165
<input type="checkbox"/>	1126040	One Piece Flip-up	Requires selection of Composite or Aluminum Angle Adjustable.		\$ 395

Footrest/Legrest Accessories

<input type="checkbox"/>	1127015	Oversized 9" x 10" Calf Pad	Not compatible with Seat Width of less than 25".		\$ 200
<input type="checkbox"/>	1127030	Swing Away ELR Gel Knee Pad			\$ 100
<input type="checkbox"/>	1127040	Adjustable Length Heel Loops		E0951	\$ 50
<input type="checkbox"/>	1127051	Bodypoint - Calf Strap		K0038	\$ 105

Extension Tubes

Regular Extension Tubes

<input type="checkbox"/>	1128010	Standard (13"-17")			
<input type="checkbox"/>	1128020	Long (15"-19")			\$ 75

Positioning Straps

Bodypoint - Hip Belts

<input type="checkbox"/>	1129010	Center-Pull Auto Style		E0978	\$ 105
<input type="checkbox"/>	1129015	Center-Pull Auto Style, Padded		E0978	\$ 155

Miscellaneous Options

<input type="checkbox"/>	1130170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.		\$ 10
--------------------------	---------	---	-----------------------	--	-------

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/eclipse to learn more about the Eclipse Bariatric Manual Wheelchair

