

Customer and Order Information		
Purchase Order *	Tag	v.060118
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type		
<input checked="" type="checkbox"/> 1701000	<b>Fuze JR 50° Pediatric Manual Tilt Wheelchair</b>	E1233 \$ 2,695

Frame Options		
<input checked="" type="checkbox"/> 1702025	<b>Adjustable/Growing Base and Seat Frame</b> 180 lbs Weight Capacity.	
<input type="checkbox"/> 1702300	<b>Transit Tie Down System (TTDS)</b>	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.		

Seat Width			Overall chair width is +8.5"
<b>Adjustable Seat Widths</b>			
<input type="checkbox"/> 1703112	12"-15" Range - Set at 12"	<input type="checkbox"/> 1703113	12"-15" Range - Set at 13"
<input type="checkbox"/> 1703114	12"-15" Range - Set at 15"	<input type="checkbox"/> 1703114	12"-15" Range - Set at 14"
<input type="checkbox"/> 1703215	15"-18" Range - Set at 15"	<input type="checkbox"/> 1703216	15"-18" Range - Set at 16"
<input type="checkbox"/> 1703217	15"-18" Range - Set at 17"	<input type="checkbox"/> 1703218	15"-18" Range - Set at 18"

Seat Depth			
<b>Standard Seat Depths (Adjustable)</b>			
<input type="checkbox"/> 1704013	13" Seat Depth	<input type="checkbox"/> 1704014	14" Seat Depth
<input type="checkbox"/> 1704015	15" Seat Depth	<input type="checkbox"/> 1704016	16" Seat Depth
<input type="checkbox"/> 1704017	17" Seat Depth	<input type="checkbox"/> 1704018	18" Seat Depth

Frame Color			
<b>PDG Premium Colors</b>			
<input type="checkbox"/> 1705000	Tracer Red	<input type="checkbox"/> 1705001	Gloss Black
<input type="checkbox"/> 1705002	Silver	<input type="checkbox"/> 1705003	Black Shimmer
<input type="checkbox"/> 1705004	Orange Fusion	<input type="checkbox"/> 1705005	Karl's Imagination
<input type="checkbox"/> 1705006	Blue Shimmer	<input type="checkbox"/> 1705007	Gloss White
<input type="checkbox"/> 1705008	Dark Titanium	<input type="checkbox"/> 1705009	Geyser Green
<input type="checkbox"/> 1705010	Lollypop Blue	<input type="checkbox"/> 1705011	Lollypop Red
<input type="checkbox"/> 1705012	Bubble Gum Pink	<input type="checkbox"/> 1705012	Bubble Gum Pink
<input type="checkbox"/> 1705013	Berry Purple	<input type="checkbox"/> 1705013	Berry Purple
<input type="checkbox"/> 1705014	Safety Yellow	<input type="checkbox"/> 1705014	Safety Yellow
<input type="checkbox"/> 1705200	<b>Custom Color</b> Contact PDG for specifications.		

Arm Type			
<input type="checkbox"/> 1706000	<b>No Armrest</b>		
<input type="checkbox"/> 1706100	<b>Single Post Fixed Height (12")</b>		
<input type="checkbox"/> 1706111	<b>Single Post Adjustable Height (8"-11")</b>	E0973 \$ 255	
<input type="checkbox"/> 1506200	<b>Dual Post Fixed Height (12")</b> Full Length Armpad only.		
<b>Dual Post Adjustable Height</b>		E0973 \$ 255	
<input type="checkbox"/> 1706211	8"-11"	<input type="checkbox"/> 1706214	10"-14"
<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.		E0973 \$ 255	
<input type="checkbox"/> 1706308	Set at 8"	<input type="checkbox"/> 1706309	Set at 9"
<input type="checkbox"/> 1706310	Set at 10"	<input type="checkbox"/> 1706311	Set at 11"
<input type="checkbox"/> 1706312	Set at 12"	<input type="checkbox"/> 1706313	Set at 13"

### Armrests

<input type="checkbox"/>	1707095	<b>Curved Armpads</b>	Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.		
<b>Regular Armpads</b>					
<input type="checkbox"/>	1707100	Desk Length - 2" x 10"	<input type="checkbox"/>	1707101	Full Length - 2" x 14"
<input type="checkbox"/>	1707102	Full Length Premium Grade - 2" x 14"			\$ 35
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707110	Desk Length - 3" x 10"	<input type="checkbox"/>	1707111	Full Length - 3" x 12"
					\$ 115
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707120	Full Length - 2" x 12"	<input type="checkbox"/>	1707121	Full Length - 3.5" x 12"
					\$ 125
<b>Contour Arm Trough with Hand Pad</b> <span style="float: right;">E2209</span>					
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707130	Right (4.5" x 18")	<input type="checkbox"/>	1707131	Left (4.5" x 18")
					\$ 125

### Seat Pan

<input type="checkbox"/>	1708000	<b>No Seat Pan</b>			
<input type="checkbox"/>	1708001	<b>Solid Aluminum Seat Pan</b>			E2231 \$ 120

### Cushion

<input type="checkbox"/>	1709000	<b>No Cushion</b>			
<b>nxt BioFit</b>					
<input type="checkbox"/>	1709101	12"-18" width, 12"-18" depth			E2607 \$ 595
<b>Matrix VI</b>					
<input type="checkbox"/>	1709111	14"-18" width, 16"-18" depth			E2607 \$ 550
<b>ROHO Quadro Select High Profile</b>					
<input type="checkbox"/>	1709121	14"-18" width, 16"-18" depth			E2624 \$ 650
<b>ROHO High Profile</b>					
<input type="checkbox"/>	1709125	14"-18" width, 16"-18" depth			E2622 \$ 650
<b>Varilite Reflex</b>					
<input type="checkbox"/>	1709131	14"-18" width, 16"-18" depth			E2603 \$ 375
<b>Cushion Size</b>					
<input type="checkbox"/>	1709400	Specify Size (Width" x Depth")	_____		

### Seat Height

#### Select desired Seat-to-Floor Height

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
12"					18"	19"
16"		15"	16"			
20"			16"	17"	18"	
22"					17"	18" 19"
24"						18" 19"

RW	8" Casters & 8" Forks					
12"					18"	19" 20"
16"			16"			20"
20"			16"	17"	18"	
22"				17"	18"	19"
24"					18"	19" 20"

Wheel size \_\_\_\_\_

Caster size \_\_\_\_\_

Fork size \_\_\_\_\_

STF height \_\_\_\_\_

### Casters

#### Urethane

<input type="checkbox"/>	1711051	5" x 1.25" Caster size		
<input type="checkbox"/>	1711061	6" x 1.25" Caster size		
<input type="checkbox"/>	1711063	6" x 1.50" Caster size		\$ 85
<input type="checkbox"/>	1711081	8" x 1.25" Caster size		
<input type="checkbox"/>	1711084	8" x 2.00" Caster size		\$ 135

#### Pneumatic

<input type="checkbox"/>	1711161	6" x 1.25" Caster size	E2214 & E2215	\$ 110
<input type="checkbox"/>	1711181	8" x 1.25" Caster size	E2214 & E2215	\$ 120
<input type="checkbox"/>	1711184	8" x 2.00" Caster size	E2214 & E2215	\$ 150

### Rear Wheels

#### Composite

<input type="checkbox"/>	1712012	12" Wheel size	<input type="checkbox"/>	1712016	16" Wheel size	<input type="checkbox"/>	1712020	20" Wheel size
<input type="checkbox"/>	1712022	22" Wheel size	<input type="checkbox"/>	1712024	24" Wheel size			

#### Spoke

<input type="checkbox"/>	1712120	20" Wheel size	<input type="checkbox"/>	1712122	22" Wheel size	<input type="checkbox"/>	1712124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------

### Axle Type

1713000 **Quick Release Rear Axles**

1713010 **Fixed Rear Axles**

### Rear Tires

<input type="checkbox"/>	1714000	<b>Urethane</b>		
<input type="checkbox"/>	1714010	<b>Pneumatic</b>	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/>	1714020	<b>Pneumatic with Flat-Free Inserts</b>	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.	E2211 & E2213 \$ 95

### Hand Rims

<input type="checkbox"/>	1715000	<b>No Hand Rims</b>		
<input type="checkbox"/>	1715010	<b>Standard</b>	20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1715020	<b>Plastic Coated</b>	20", 22", and 24" Rear Wheels only.	\$ 135

### Wheel Locks

<input type="checkbox"/>	1716005	<b>Foot Operated Wheel Lock</b>		\$ 200
<input type="checkbox"/>	1716010	<b>Push to Lock</b>		
<input type="checkbox"/>	1716020	<b>Pull to Lock</b>		
<input type="checkbox"/>	1716030	<b>Attendant Wheel Lock (Push to Lock)</b>	Not compatible with Wheel Lock Extensions.	\$ 85
<input type="checkbox"/>	1716100	<b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.	E0961 \$ 60

### Anti-Tippers

<input checked="" type="checkbox"/>	1717010	<b>Rear Anti-Tippers</b>	Required Safety Option.	E0971 \$ 60
-------------------------------------	---------	--------------------------	-------------------------	-------------

### Back Post Style

<b>Push Bar Back</b>		
<input type="checkbox"/> 1718010	Straight	
<input type="checkbox"/> 1718030	8° Bend	\$ 75
<b>Cane Handle Back</b> Not compatible with Headrest Mounting Plate.		
<input type="checkbox"/> 1718020	Straight	
<input type="checkbox"/> 1718040	8° Bend	\$ 75
<b>Back Post Height</b> Seat Pan to Handle.		
<input type="checkbox"/> 1718420	20" Height	<input type="checkbox"/> 1718425 25" Height

### Back Options

<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 1719010	Button on Outside	<input type="checkbox"/> 1719011 Button on Inside
<input type="checkbox"/> 1719050	<b>Removable Cane Handle Push Bar</b> Not compatible with Push Bar Back.	\$ 140
<b>Seat to Back Angle</b>		
<input type="checkbox"/> 1719090	Set at 90°	<input type="checkbox"/> 1719096 Set at 96°
<input type="checkbox"/> 1719108	Set at 108°	<input type="checkbox"/> 1719114 Set at 114°
<input type="checkbox"/> 1719102	Set at 102°	<input type="checkbox"/> 1719120 Set at 120°
<input type="checkbox"/> 1719220	<b>Fold Down Back</b>	
<input type="checkbox"/> 1719221	<b>Non Fold Down Back</b>	\$ 75

### Back Upholstery

<input checked="" type="checkbox"/> 1720000	<b>No Back Upholstery</b>
---	---------------------------

### Backrest

<input type="checkbox"/> 1721000	<b>No Backrest</b>		
<b>nxt Xtend</b>			
<input type="checkbox"/> 1721100	16"-18" width, 16"-20" height		E2620 \$ 690
<b>nxt Xtend Deep</b>			
<input type="checkbox"/> 1721101	16"-16" width, 16"-20" height		E2620 \$ 690
<b>Matrx Elite</b>			
<input type="checkbox"/> 1721110	14"-18" width, 10"-20" height		E2620 \$ 825
<input type="checkbox"/> 1721116	Deep - 14"-18" width, 12"-20" height		E2620 \$ 950
<b>ROHO Agility</b>			
<input type="checkbox"/> 1721120	Min Contour 14"-18" width, 10"-18" height		E2620 \$ 825
<input type="checkbox"/> 1721122	Mid Contour 14"-18" width, 10"-18" height		E2620 \$ 850
<input type="checkbox"/> 1721124	Max Contour 14"-18" width, 10"-18" height		E2620 \$ 950
<b>Varilite Icon</b>			
<input type="checkbox"/> 1721130	Tall - 14"-18" width, 14"-24" height		E2615 \$ 850
<input type="checkbox"/> 1721132	Deep - 14"-18" width, 17"-24" height		E2620 \$ 895
<b>Backrest Size</b>			
1521400	Specify Size (Width" x Height")	_____	

## Tilt/Recline Options

### Attendant Tilt only

- 1723011 Trigger on Right Back Post       1723012 Trigger on Left Back Post

### Tilt Limit

- 1723220 Set at 20°       1723225 Set at 25°       1723230 Set at 30°  
 1723235 Set at 35°       1723240 Set at 40°       1723245 Set at 45°  
 1723250 Set at 50°

- 1723410 **Reclining Backrest (Manual)** E1225 \$ 875  
 30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Requires 25" Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15" or 16" Seat Depth.

## Headrest

- 1724000 **No Headrest**
- 1724010 **Headrest Mounting Plate** \$ 75  
 Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.
- 1724100 **nxt MultiFit Head Support Hardware** E1028 \$ 235
- 1724100 **nxt MultiFit Head Support Pad**  E0955 \$ 245
- 1724105 Regular (9" x 5")       1724106 Large (10" x 5")       1724107 X-Large (14" x 5")
- Stealth Headrest Hardware**
- 1724110 Multi-Axis HMO475 E1028 \$ 150
- 1724111 Link TWB480 E1028 \$ 325
- Stealth Headrest Pad**  E0955 \$ 245
- 1724115 10" width       1724116 14" width

## Front Rigging

- 1725000 **No Front Rigging**       1725001 **Sector Blocks only**
- PDG Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame.
- 1725410 Fixed 70° \$ 95
- 1725420 Swing In/Away Universal Angle Adjustable (0°-110°)      Set at (use increments of 10°): \_\_\_\_\_ \$ 195
- 1725430 Swing In/Away Elevating Legrests      High knee pivot. Includes reinforced premium Calf Pad mount. E0990 \$ 495
- Standard Front Rigging** Legacy Design. Not compatible with Easy-Out™.
- 1725140 Swing Away 70°
- 1525160 Swing Away 90°      Not compatible with 8" Casters. \$ 55
- 1725200 Swing Away Elevating Legrests E0990 \$ 375
- 1725300 1" Offset Hangers for Swing Away \$ 75

## Foot Plates

- 1726010 **Composite**      Not available on 12" Seat Width.
- 1726020 **Aluminum Angle Adjustable** K0040 \$ 155
- 1726040 **One Piece Flip-up**      Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
- 1726050 **Interlocking Angle Adjustable**      Seat Widths of 14"-18" only. K0040 \$ 195

## Footrest/Legrest Accessories

- 1727015 **Swing Away 70° Calf Pad**      Not compatible with 70° Gel Padding, Seat Widths of 12" and 13". \$ 135
- 1727020 **Swing Away 70° Gel Padding**      Not compatible with 70° Calf Pad, Easy-Out Front Rigging. \$ 100
- 1727030 **Swing Away ELR Gel Knee Pad** \$ 100
- 1727040 **Adjustable Length Heel Loops** E0951 \$ 50
- 1727050 **Bodypoint - Padded Calf Strap** K0038 \$ 105
- 1727065 **Miller's Dynamic Footrest Mounts**      Not compatible with Easy-Out Front Rigging. \$ 380

Extension Tubes			
<input type="checkbox"/>	1728010	<b>Short</b> 14"-18" for Easy-Out Fixed & SA Angle Adjustable 11"-15" for Easy-Out SA ELR 13"-17" for Standard Front Rigging	
<input type="checkbox"/>	1728020	<b>Long</b> 18"-24" for Easy-Out Fixed & SA Angle Adjustable 15"-21" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$ 75
<input type="checkbox"/>	1728200	<b>For 90° FR (11"-16")</b> Not compatible with Easy-Out, Standard 70° or ELR	\$ 75

Positioning Straps			
<b>Bodypoint - Pediatric Hip Belts</b>			
<input type="checkbox"/>	1729017	Center-Pull Auto Style, Padded	E0978 \$ 155
<input type="checkbox"/>	1729022	4-Point, Padded	E0978 \$ 175
<input type="checkbox"/>	1729030	<b>Bodypoint - Chest Strap</b>	E0960 \$ 85

Miscellaneous Options			
<input type="checkbox"/>	1730160	<b>Clear Lap Tray (Slide On)</b> Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	E0950 \$ 490
<input type="checkbox"/>	1730170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/fuze-T50Jr/](http://pdgmobility.com/products/fuze-T50Jr/) to learn more about the **Fuze JR 50° Manual Tilt Wheelchair**

