

Customer and Order Information

Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 1701000	Fuze JR 50° Pediatric Manual Tilt Wheelchair	E1233	\$ 2,695
---	---	-------	----------

Frame Options

<input checked="" type="checkbox"/> 1702025	Adjustable/Growing Base and Seat Frame	180 lbs Weight Capacity.
<input type="checkbox"/> 1702300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.

Seat Width Overall chair width is +8.5"

Adjustable Seat Widths					
<input type="checkbox"/> 1703112	12"-15" Range - Set at 12"	<input type="checkbox"/> 1703113	12"-15" Range - Set at 13"	<input type="checkbox"/> 1703114	12"-15" Range - Set at 14"
<input type="checkbox"/> 1703114	12"-15" Range - Set at 15"				
<input type="checkbox"/> 1703215	15"-18" Range - Set at 15"	<input type="checkbox"/> 1703216	15"-18" Range - Set at 16"	<input type="checkbox"/> 1703217	15"-18" Range - Set at 17"
<input type="checkbox"/> 1703218	15"-18" Range - Set at 18"				

Seat Depth

Standard Seat Depths (Adjustable)					
<input type="checkbox"/> 1704013	13" Seat Depth	<input type="checkbox"/> 1704014	14" Seat Depth	<input type="checkbox"/> 1704015	15" Seat Depth
<input type="checkbox"/> 1704016	16" Seat Depth	<input type="checkbox"/> 1704017	17" Seat Depth	<input type="checkbox"/> 1704018	18" Seat Depth

Frame Color

PDG Premium Colors					
<input type="checkbox"/> 1705000	Tracer Red	<input type="checkbox"/> 1705001	Gloss Black	<input type="checkbox"/> 1705002	Silver
<input type="checkbox"/> 1705003	Black Shimmer	<input type="checkbox"/> 1705004	Orange Fusion	<input type="checkbox"/> 1705005	Karl's Imagination
<input type="checkbox"/> 1705006	Blue Shimmer	<input type="checkbox"/> 1705007	Gloss White	<input type="checkbox"/> 1705008	Dark Titanium
<input type="checkbox"/> 1705009	Geyser Green	<input type="checkbox"/> 1705010	Lollypop Blue	<input type="checkbox"/> 1705011	Lollypop Red
<input type="checkbox"/> 1705012	Bubble Gum Pink	<input type="checkbox"/> 1705013	Berry Purple	<input type="checkbox"/> 1705014	Safety Yellow
<input type="checkbox"/> 1705200	Custom Color	Contact PDG for specifications.			

Arm Type

<input type="checkbox"/> 1706000	No Armrest		
<input type="checkbox"/> 1706100	Single Post Fixed Height (12")		
<input type="checkbox"/> 1706111	Single Post Adjustable Height (8"-11")	E0973	\$ 255
<input type="checkbox"/> 1506200	Dual Post Fixed Height (12")	Full Length Armpad only.	
Dual Post Adjustable Height		E0973	\$ 255
<input type="checkbox"/> 1706211	8"-11"	<input type="checkbox"/> 1706214	10"-14"
Locking Flip-up		Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.	
<input type="checkbox"/> 1706308	Set at 8"	<input type="checkbox"/> 1706309	Set at 9"
<input type="checkbox"/> 1706310	Set at 10"	<input type="checkbox"/> 1706311	Set at 11"
<input type="checkbox"/> 1706312	Set at 12"	<input type="checkbox"/> 1706313	Set at 13"

Armrests

<input type="checkbox"/>	1707095	Curved Armpads	Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.		
Regular Armpads					
<input type="checkbox"/>	1707100	Desk Length - 2" x 10"	<input type="checkbox"/>	1707101	Full Length - 2" x 14"
<input type="checkbox"/>	1707102	Full Length Premium Grade - 2" x 14"			\$ 35
Wide Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707110	Desk Length - 3" x 10"	<input type="checkbox"/>	1707111	Full Length - 3" x 12"
					\$ 115
Gel Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707120	Full Length - 2" x 12"	<input type="checkbox"/>	1707121	Full Length - 3.5" x 12"
					\$ 125
Contour Arm Trough with Hand Pad					
					E2209 \$ 125
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.					
<input type="checkbox"/>	1707130	Right (4.5" x 18")	<input type="checkbox"/>	1707131	Left (4.5" x 18")

Seat Pan

<input type="checkbox"/>	1708000	No Seat Pan			
<input type="checkbox"/>	1708001	Solid Aluminum Seat Pan			E2231 \$ 120

Cushion

<input type="checkbox"/>	1709000	No Cushion			
nxt BioFit					
<input type="checkbox"/>	1709101	12"-18" width, 12"-18" depth			E2607 \$ 595
Matrix VI					
<input type="checkbox"/>	1709111	14"-18" width, 16"-18" depth			E2607 \$ 550
ROHO Quadro Select High Profile					
<input type="checkbox"/>	1709121	14"-18" width, 16"-18" depth			E2624 \$ 650
ROHO High Profile					
<input type="checkbox"/>	1709125	14"-18" width, 16"-18" depth			E2622 \$ 650
Varilite Reflex					
<input type="checkbox"/>	1709131	14"-18" width, 16"-18" depth			E2603 \$ 375
Cushion Size					
<input type="checkbox"/>	1709400	Specify Size (Width" x Depth")	_____		

Seat Height

Select desired Seat-to-Floor Height

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters				
12"	13"	14"			
16"		14"	15"	16"	
20"			16"	17"	
22"				17"	

RW	6" Casters & 6" Forks				
12"		14"			18"
16"		14"	15"	16"	
20"			16"	17"	18"
22"				17"	18"
24"					18"

RW	6" Casters & 8" Forks				
12"				18"	19"
16"		15"	16"		
20"			16"	17"	18"
22"				17"	18"
24"					18"

RW	8" Casters & 8" Forks				
12"				18"	19"
16"			16"		20"
20"			16"	17"	18"
22"				17"	18"
24"					18"

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

<input type="checkbox"/>	1711051	5" x 1.25" Caster size		
<input type="checkbox"/>	1711061	6" x 1.25" Caster size		
<input type="checkbox"/>	1711063	6" x 1.50" Caster size		\$ 85
<input type="checkbox"/>	1711081	8" x 1.25" Caster size		
<input type="checkbox"/>	1711084	8" x 2.00" Caster size		\$ 135

Pneumatic

<input type="checkbox"/>	1711161	6" x 1.25" Caster size	E2214 & E2215	\$ 110
<input type="checkbox"/>	1711181	8" x 1.25" Caster size	E2214 & E2215	\$ 120
<input type="checkbox"/>	1711184	8" x 2.00" Caster size	E2214 & E2215	\$ 150

Rear Wheels

Composite

<input type="checkbox"/>	1712012	12" Wheel size	<input type="checkbox"/>	1712016	16" Wheel size	<input type="checkbox"/>	1712020	20" Wheel size
<input type="checkbox"/>	1712022	22" Wheel size	<input type="checkbox"/>	1712024	24" Wheel size			

Spoke

<input type="checkbox"/>	1712120	20" Wheel size	<input type="checkbox"/>	1712122	22" Wheel size	<input type="checkbox"/>	1712124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------

Axle Type

1713000 **Quick Release Rear Axles**

1713010 **Fixed Rear Axles**

Rear Tires

<input type="checkbox"/>	1714000	Urethane		
<input type="checkbox"/>	1714010	Pneumatic	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/>	1714020	Pneumatic with Flat-Free Inserts	Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.	E2211 & E2213 \$ 95

Hand Rims

<input type="checkbox"/>	1715000	No Hand Rims		
<input type="checkbox"/>	1715010	Standard	20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1715020	Plastic Coated	20", 22", and 24" Rear Wheels only.	\$ 135

Wheel Locks

<input type="checkbox"/>	1716005	Foot Operated Wheel Lock		\$ 200
<input type="checkbox"/>	1716010	Push to Lock		
<input type="checkbox"/>	1716020	Pull to Lock		
<input type="checkbox"/>	1716030	Attendant Wheel Lock (Push to Lock)	Not compatible with Wheel Lock Extensions.	\$ 85
<input type="checkbox"/>	1716100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.	E0961 \$ 60

Anti-Tippers

<input checked="" type="checkbox"/>	1717010	Rear Anti-Tippers	Required Safety Option.	E0971 \$ 60
-------------------------------------	---------	--------------------------	-------------------------	-------------

Back Post Style

Push Bar Back		
<input type="checkbox"/> 1718010	Straight	
<input type="checkbox"/> 1718030	8° Bend	\$ 75
Cane Handle Back Not compatible with Headrest Mounting Plate.		
<input type="checkbox"/> 1718020	Straight	
<input type="checkbox"/> 1718040	8° Bend	\$ 75
Back Post Height Seat Pan to Handle.		
<input type="checkbox"/> 1718420	20" Height	<input type="checkbox"/> 1718425 25" Height


Back Options

Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 1719010	Button on Outside	<input type="checkbox"/> 1719011 Button on Inside
<input type="checkbox"/> 1719050	Removable Cane Handle Push Bar Not compatible with Push Bar Back.	\$ 140
Seat to Back Angle		
<input type="checkbox"/> 1719090	Set at 90°	<input type="checkbox"/> 1719096 Set at 96°
<input type="checkbox"/> 1719108	Set at 108°	<input type="checkbox"/> 1719114 Set at 114°
<input type="checkbox"/> 1719102	Set at 102°	<input type="checkbox"/> 1719120 Set at 120°
<input type="checkbox"/> 1719220	Fold Down Back	
<input type="checkbox"/> 1719221	Non Fold Down Back	\$ 75

Back Upholstery

<input checked="" type="checkbox"/> 1720000	No Back Upholstery
---	---------------------------

Backrest

<input type="checkbox"/> 1721000	No Backrest		
nxt Optima HA Active Contour			
<input type="checkbox"/> 1721100	16"-18" width, 16"-20" height		E2620 \$ 690
nxt Optima HA Plus Contour			
<input type="checkbox"/> 1721101	16"-16" width, 16"-20" height		E2620 \$ 690
Matrx Elite			
<input type="checkbox"/> 1721110	14"-18" width, 10"-20" height		E2620 \$ 825
<input type="checkbox"/> 1721116	Deep - 14"-18" width, 12"-20" height		E2620 \$ 950
ROHO Agility			
<input type="checkbox"/> 1721120	Min Contour 14"-18" width, 10"-18" height		E2620 \$ 825
<input type="checkbox"/> 1721122	Mid Contour 14"-18" width, 10"-18" height		E2620 \$ 850
<input type="checkbox"/> 1721124	Max Contour 14"-18" width, 10"-18" height		E2620 \$ 950
Varilite Icon			
<input type="checkbox"/> 1721130	Tall - 14"-18" width, 14"-24" height		E2615 \$ 850
<input type="checkbox"/> 1721132	Deep - 14"-18" width, 17"-24" height		E2620 \$ 895
Backrest Size			
1521400	Specify Size (Width" x Height")	_____	

Tilt/Recline Options

Attendant Tilt only

- 1723011 Trigger on Right Back Post 1723012 Trigger on Left Back Post

Tilt Limit

- 1723220 Set at 20° 1723225 Set at 25° 1723230 Set at 30°
 1723235 Set at 35° 1723240 Set at 40° 1723245 Set at 45°
 1723250 Set at 50°

- 1723410 **Reclining Backrest (Manual)** E1225 \$ 875
 30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Requires 25" Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15" or 16" Seat Depth.

Headrest

- 1724000 **No Headrest**
- 1724010 **Headrest Mounting Plate** \$ 75
 Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.
- 1724100 **nxt MultiFit Head Support Hardware** E1028 \$ 235
- 1724100 **nxt MultiFit Head Support Pad**  E0955 \$ 245
- 1724105 Regular (9" x 5") 1724106 Large (10" x 5") 1724107 X-Large (14" x 5")
- Stealth Headrest Hardware**
- 1724110 Multi-Axis HMO475 E1028 \$ 150
- 1724111 Link TWB480 E1028 \$ 325
- Stealth Headrest Pad**  E0955 \$ 245
- 1724115 10" width 1724116 14" width

Front Rigging

- 1725000 **No Front Rigging**
- PDG Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame.
- 1725410 Fixed 70° Available in March \$ 95
- 1725420 Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____ Available in April \$ 195
- 1725430 Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount. Available in April E0990 \$ 495
- Standard Front Rigging** Legacy Design. Not compatible with Easy-Out™.
- 1725140 Swing Away 70°
- 1525160 Swing Away 90° Not compatible with 8" Casters. \$ 55
- 1725200 Swing Away Elevating Legrests E0990 \$ 375
- 1725300 1" Offset Hangers for Swing Away \$ 75

Foot Plates

- 1726010 **Composite** Not available on 12" Seat Width.
- 1726020 **Aluminum Angle Adjustable** K0040 \$ 155
- 1726040 **One Piece Flip-up** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
- 1726050 **Interlocking Angle Adjustable** Seat Widths of 14"-18" only. K0040 \$ 195

Footrest/Legrest Accessories

- 1727015 **Swing Away 70° Calf Pad** Not compatible with 70° Gel Padding, Seat Widths of 12" and 13". \$ 135
- 1727020 **Swing Away 70° Gel Padding** Not compatible with 70° Calf Pad, Easy-Out Front Rigging. \$ 100
- 1727030 **Swing Away ELR Gel Knee Pad** \$ 100
- 1727040 **Adjustable Length Heel Loops** E0951 \$ 50
- 1727050 **Bodypoint - Padded Calf Strap** K0038 \$ 105
- 1727065 **Miller's Dynamic Footrest Mounts** Not compatible with Easy-Out Front Rigging. \$ 380

Extension Tubes

<input type="checkbox"/>	1728010	Short	13"-17" for Easy-Out SA Angle Adjustable 10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging		
<input type="checkbox"/>	1728020	Long	16"-22.5" for Easy-Out SA Angle Adjustable 12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging	\$	75
<input type="checkbox"/>	1728200	For 90° FR (11"-16")	Not compatible with Easy-Out, Standard 70° or ELR	\$	75

Positioning Straps

Bodypoint - Pediatric Hip Belts						
<input type="checkbox"/>	1729017	Center-Pull Auto Style, Padded		E0978	\$	155
<input type="checkbox"/>	1729022	4-Point, Padded		E0978	\$	175
<input type="checkbox"/>	1729030	Bodypoint - Chest Strap		E0960	\$	85

Miscellaneous Options

<input type="checkbox"/>	1730160	Clear Lap Tray (Slide On)	Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	E0950	\$	490
<input type="checkbox"/>	1730170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.		\$	10

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/fuze-T50Jr/ to learn more about the **Fuze JR 50° Manual Tilt Wheelchair**

