

Customer and Order Information

| | | |
|------------------------------------------|-----|-----------|
| Purchase Order * | Tag | v.031518 |
| Dealer/Provider Name * | | |
| Purchaser Name and Contact Information * | | |
| Billing Address * | | |
| Shipping Address * | | *required |

Frame Type

| | | | |
|---------------------------------------------|--------------------------------------------|-------|----------|
| <input checked="" type="checkbox"/> 1601000 | Fuze T20 20° Manual Tilt Wheelchair | E1161 | \$ 2,995 |
|---------------------------------------------|--------------------------------------------|-------|----------|

Frame Options

| | | | |
|----------------------------------------------------------------------------------------------|---------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|
| <input type="checkbox"/> 1602010 | Fixed Base and Seat Frame | 250 lbs Weight Capacity. Select Fixed Seat Width below. | |
| <input type="checkbox"/> 1602025 | Adjustable Base and Seat Frame | 250 lbs Weight Capacity. Select Adjustable Seat Width below. | \$ 100 |
| <input type="checkbox"/> 1602300 | Transit Tie Down System (TTDS) | TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. | |
| User Weight Range Required for selection of Gas Struts (Force measured in Newton - N) | | | |
| <input type="checkbox"/> 1602406 | Up to 150 lbs (600 N) | <input type="checkbox"/> 1602410 | 150 lbs to 200 lbs (1000 N) |
| | | <input type="checkbox"/> 1602412 | 200 lbs to 250 lbs (1200 N) |

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

| | | | |
|--------------------------------------------------------------------------|----------------------------|----------------------------------|----------------------------|
| Fixed Seat Widths Order with Fixed Base and Seat Frame. | | | |
| <input type="checkbox"/> 1603015 | 15" Seat Width | <input type="checkbox"/> 1603016 | 16" Seat Width |
| <input type="checkbox"/> 1603017 | 17" Seat Width | <input type="checkbox"/> 1603018 | 18" Seat Width |
| <input type="checkbox"/> 1603019 | 19" Seat Width | <input type="checkbox"/> 1603020 | 20" Seat Width |
| <input type="checkbox"/> 1603021 | 21" Seat Width | <input type="checkbox"/> 1603022 | 22" Seat Width |
| | | E2201 | \$ 325 |
| Adjustable Seat Widths Order with Adjustable Base and Seat Frame. | | | |
| <input type="checkbox"/> 1603112 | 12"-15" Range - Set at 12" | <input type="checkbox"/> 1603113 | 12"-15" Range - Set at 13" |
| <input type="checkbox"/> 1603114 | 12"-15" Range - Set at 14" | <input type="checkbox"/> 1603115 | 12"-15" Range - Set at 15" |
| <input type="checkbox"/> 1603215 | 15"-20" Range - Set at 15" | <input type="checkbox"/> 1603216 | 15"-20" Range - Set at 16" |
| <input type="checkbox"/> 1603217 | 15"-20" Range - Set at 17" | <input type="checkbox"/> 1603218 | 15"-20" Range - Set at 18" |
| <input type="checkbox"/> 1603219 | 15"-20" Range - Set at 19" | <input type="checkbox"/> 1603220 | 15"-20" Range - Set at 20" |
| | | E2201 | \$ 325 |
| <input type="checkbox"/> 1603321 | 21"-22" Range - Set at 21" | <input type="checkbox"/> 1603322 | 21"-22" Range - Set at 22" |
| | | E2201 | \$ 325 |

Seat Depth

| | | | |
|------------------------------------------|----------------|----------------------------------|----------------|
| Standard Seat Depths (Adjustable) | | | |
| <input type="checkbox"/> 1604015 | 15" Seat Depth | <input type="checkbox"/> 1604016 | 16" Seat Depth |
| <input type="checkbox"/> 1604017 | 17" Seat Depth | <input type="checkbox"/> 1604018 | 18" Seat Depth |
| <input type="checkbox"/> 1604019 | 19" Seat Depth | <input type="checkbox"/> 1604020 | 20" Seat Depth |
| | | E2203 | \$ 475 |





Frame Color

| | | | |
|----------------------------------|-----------------|----------------------------------|-----------------------------------------------------|
| PDG Premium Colors | | | |
| <input type="checkbox"/> 1605000 | Tracer Red | <input type="checkbox"/> 1605001 | Gloss Black |
| <input type="checkbox"/> 1605002 | Silver | <input type="checkbox"/> 1605003 | Black Shimmer |
| <input type="checkbox"/> 1605004 | Orange Fusion | <input type="checkbox"/> 1605005 | Karl's Imagination |
| <input type="checkbox"/> 1605006 | Blue Shimmer | <input type="checkbox"/> 1605007 | Gloss White |
| <input type="checkbox"/> 1605008 | Dark Titanium | <input type="checkbox"/> 1605009 | Geyser Green |
| <input type="checkbox"/> 1605010 | Lollypop Blue | <input type="checkbox"/> 1605011 | Lollypop Red |
| <input type="checkbox"/> 1605012 | Bubble Gum Pink | <input type="checkbox"/> 1605013 | Berry Purple |
| <input type="checkbox"/> 1605014 | Safety Yellow | <input type="checkbox"/> 1605200 | Custom Color Contact PDG for specifications. |
| | | | \$ 300 |

| Arm Type | | | |
|--------------------------|---------|---------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> | 1606000 | No Armrest | |
| <input type="checkbox"/> | 1606100 | Single Post Fixed Height (12") | |
| | | Single Post Adjustable Height | 8"-11" not compatible with 22" and 24" Rear Wheels. E0973 \$ 255 |
| <input type="checkbox"/> | 1606111 | 8"-11" | <input type="checkbox"/> 1606114 10"-14" |
| <input type="checkbox"/> | 1606200 | Dual Post Fixed Height (12") | Full Length Armpad only. |
| | | Dual Post Adjustable Height | 8"-11" not compatible with 22" and 24" Rear Wheels. E0973 \$ 255 |
| <input type="checkbox"/> | 1606211 | 8"-11" | <input type="checkbox"/> 1606214 10"-14" <input type="checkbox"/> 1606216 12"-16" |
| | | Locking Flip-up | Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery. E0973 \$ 255 |
| <input type="checkbox"/> | 1606310 | Set at 10" | <input type="checkbox"/> 1606311 Set at 11" <input type="checkbox"/> 1606312 Set at 12" |
| <input type="checkbox"/> | 1606313 | Set at 13" | <input type="checkbox"/> 1606314 Set at 14" |

| Armrests | | | |
|--------------------------|---------|-----------------------------------------|--------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> | 1607095 | Curved Armpads | Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests, 22" and 24" Rear Wheels. |
| | | Regular Armpads | |
| <input type="checkbox"/> | 1607100 | Desk Length - 2" x 10" | <input type="checkbox"/> 1607101 Full Length - 2" x 14" |
| <input type="checkbox"/> | 1607102 | Full Length Premium Grade - 2" x 14" | \$ 35 |
| | | Wide Armpads | May interfere with wheels in some tilt configurations. \$ 115 |
| <input type="checkbox"/> | 1607110 | Desk Length - 3" x 10" | <input type="checkbox"/> 1607111 Full Length - 3" x 12" |
| | | Gel Armpads | May interfere with wheels in some tilt configurations. \$ 125 |
| <input type="checkbox"/> | 1607120 | Full Length - 2" x 12" | <input type="checkbox"/> 1607121 Full Length - 3.5" x 12" |
| | | Contour Arm Trough with Hand Pad | May interfere with wheels in some tilt configurations. E2209 \$ 125 |
| <input type="checkbox"/> | 1607130 | Right (4.5" x 18") | <input type="checkbox"/> 1607131 Left (4.5" x 18") |

| Seat Pan | | | |
|--------------------------|---------|--------------------------------|--------------|
| <input type="checkbox"/> | 1608000 | No Seat Pan | |
| <input type="checkbox"/> | 1608001 | Solid Aluminum Seat Pan | E2231 \$ 120 |

| Cushion | | | |
|--------------------------|---------|----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> | 1609000 | No Cushion | |
| | | nxt BioFit | |
| <input type="checkbox"/> | 1609101 | 14"-20" width, 16"-20" depth |  E2607 \$ 595 |
| <input type="checkbox"/> | 1609102 | 22" width, 16"-20" depth | E2608 \$ 625 |
| | | Matrx VI | |
| <input type="checkbox"/> | 1609111 | 14"-20" width, 16"-20" depth |  E2607 \$ 550 |
| <input type="checkbox"/> | 1609112 | 22" width, 18"-20" depth | E2608 \$ 600 |
| <input type="checkbox"/> | 1609115 | HD - 18"-22" width, 18"-20" depth | E2608 \$ 1,150 |
| | | ROHO Quadro Select High Profile | |
| <input type="checkbox"/> | 1609121 | 14"-20" width, 16"-20" depth |  E2624 \$ 650 |
| <input type="checkbox"/> | 1609122 | 22" width, 16"-20" depth | E2625 \$ 750 |
| | | ROHO High Profile | |
| <input type="checkbox"/> | 1609125 | 14"-20" width, 16"-20" depth | E2622 \$ 650 |
| <input type="checkbox"/> | 1609126 | 22" width, 16"-20" depth | E2623 \$ 750 |
| | | Varilite Reflex | |
| <input type="checkbox"/> | 1609131 | 14"-20" width, 16"-20" depth |  E2603 \$ 375 |
| | | Cushion Size | |
| <input type="checkbox"/> | 1609400 | Specify Size (Width" x Depth") | |

Seat Height

Select desired Seat-to-Floor Height

Standard T20 Options

| RW | 5" Casters | | | | | | |
|-----|------------|-----|-----|-----|-----|--|--|
| 12" | 13" | | | | | | |
| 16" | 13" | 14" | 15" | | | | |
| 20" | 13" | 14" | 15" | 16" | 17" | | |
| 22" | | 14" | 15" | 16" | 17" | | |
| 24" | | | 15" | 16" | 17" | | |

| RW | 6" Casters & 6" Forks | | | | | | |
|-----|-----------------------|-----|-----|-----|-----|-----|--|
| 16" | | 14" | 15" | | | | |
| 20" | | 14" | 15" | 16" | 17" | | |
| 22" | | 14" | 15" | 16" | 17" | 18" | |
| 24" | | | 15" | 16" | 17" | 18" | |

| RW | 6" Casters & 8" Forks | | | | | | |
|-----|-----------------------|--|-----|-----|-----|-----|-----|
| 16" | | | 15" | | | | |
| 20" | | | 15" | 16" | 17" | | |
| 22" | | | 15" | 16" | 17" | 18" | |
| 24" | | | 15" | 16" | 17" | 18" | 19" |

| RW | 8" Casters & 8" Forks | | | | | | |
|-----|-----------------------|--|--|-----|-----|-----|-----|
| 20" | | | | 16" | 17" | | |
| 22" | | | | 16" | 17" | 18" | |
| 24" | | | | 16" | 17" | 18" | 19" |

Narrow T20 Options

Reduces the overall width of the Chair by 2". Recommended only if the width of the chair is a concern. Some STF combinations offer less than 20° of tilt, but all offer a minimum of 12°.

| RW | 5" Casters | | | | | | |
|-----|------------|-----|-----|-----|-----|--|--|
| 12" | 13" | 14" | | | | | |
| 16" | | 14" | 15" | 16" | | | |
| 20" | | | | 16" | 17" | | |
| 22" | | | | | 17" | | |

| RW | 6" Casters & 6" Forks | | | | | | |
|-----|-----------------------|-----|-----|-----|-----|-----|--|
| 12" | | 14" | | | | | |
| 16" | | 14" | 15" | 16" | | | |
| 20" | | | | 16" | 17" | 18" | |
| 22" | | | | | 17" | 18" | |
| 24" | | | | | | 18" | |

| RW | 6" Casters & 8" Forks | | | | | | |
|-----|-----------------------|--|-----|-----|-----|-----|-----|
| 16" | | | 15" | 16" | | | |
| 20" | | | | 16" | 17" | 18" | |
| 22" | | | | | 17" | 18" | 19" |
| 24" | | | | | | 18" | 19" |

| RW | 8" Casters & 8" Forks | | | | | | |
|-----|-----------------------|--|--|-----|-----|-----|-----|
| 16" | | | | 16" | | | |
| 20" | | | | 16" | 17" | 18" | |
| 22" | | | | | 17" | 18" | 19" |
| 24" | | | | | | 18" | 19" |

- Standard
 Narrow

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

- 1611051 5" x 1.25" Caster size
- 1611061 6" x 1.25" Caster size
- 1611063 6" x 1.50" Caster size \$ 85
- 1611081 8" x 1.25" Caster size
- 1611084 8" x 2.00" Caster size \$ 135

Pneumatic

- 1611161 6" x 1.25" Caster size E2214 & E2215 \$ 110
- 1611181 8" x 1.25" Caster size E2214 & E2215 \$ 120
- 1611184 8" x 2.00" Caster size E2214 & E2215 \$ 150

- 1611260 **6" Soft Roll** \$ 95

Rear Wheels

Composite

- 1612012 12" Wheel size
- 1612016 16" Wheel size
- 1612020 20" Wheel size
- 1612022 22" Wheel size
- 1612024 24" Wheel size

Spoke

- 1612120 20" Wheel size \$ 75
- 1612122 22" Wheel size
- 1612124 24" Wheel size

- 1612224 **24" Heavy Duty Spoke** \$ 350

Axle Type

- 1613000 **Quick Release Rear Axles**
- 1613010 **Fixed Rear Axles**

| Rear Tires | | | |
|--------------------------|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| <input type="checkbox"/> | 1614000 | Urethane | |
| <input type="checkbox"/> | 1614010 | Pneumatic Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. | E2211 & E2212 \$ 60 |
| <input type="checkbox"/> | 1614020 | Pneumatic with Flat-Free Inserts Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. | E2211 & E2213 \$ 95 |

| Hand Rims | | | |
|--------------------------|---------|-----------------------------------------------------------|--------|
| <input type="checkbox"/> | 1615000 | No Hand Rims | |
| <input type="checkbox"/> | 1615010 | Standard 20", 22", and 24" Rear Wheels only. | |
| <input type="checkbox"/> | 1615020 | Plastic Coated 20", 22", and 24" Rear Wheels only. | \$ 135 |

| Wheel Locks | | | |
|--------------------------|---------|---------------------------------------------------------------------------------------|-------------|
| <input type="checkbox"/> | 1616010 | Push to Lock | |
| <input type="checkbox"/> | 1616020 | Pull to Lock | |
| <input type="checkbox"/> | 1616030 | Attendant Wheel Lock (Push to Lock) Not compatible with Wheel Lock Extensions. | \$ 85 |
| <input type="checkbox"/> | 1616040 | User & Attendant Wheel Lock | \$ 125 |
| <input type="checkbox"/> | 1616100 | Wheel Lock Extensions Not compatible with Attendant Wheel Lock. | E0961 \$ 60 |

| Anti-Tippers | | | |
|-------------------------------------|---------|--------------------------------------------------|-------------|
| <input checked="" type="checkbox"/> | 1617010 | Rear Anti-Tippers Required Safety Option. | E0971 \$ 60 |

| Back Post Style | | | |
|----------------------------------------------------------------------|---------|------------|-------|
| Push Bar Back | | | |
| <input type="checkbox"/> | 1618010 | Straight | |
| <input type="checkbox"/> | 1618030 | 8° Bend | \$ 75 |
| Cane Handle Back Not compatible with Headrest Mounting Plate. | | | |
| <input type="checkbox"/> | 1618020 | Straight | |
| <input type="checkbox"/> | 1618040 | 8° Bend | \$ 75 |
| Back Post Height Seat Pan to Handle. | | | |
| <input type="checkbox"/> | 1618420 | 20" Height | |
| <input type="checkbox"/> | 1618425 | 25" Height | |

| Back Options | | | |
|-----------------------------------------------------------------------------------------------------------------|---------|--------------------------------------------------------------------------|--------|
| Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back. | | | |
| <input type="checkbox"/> | 1619010 | Button on Outside | |
| <input type="checkbox"/> | 1619011 | Button on Inside | \$ 150 |
| <input type="checkbox"/> | 1619050 | Removable Cane Handle Push Bar Not compatible with Push Bar Back. | \$ 140 |
| Seat to Back Angle | | | |
| <input type="checkbox"/> | 1619090 | Set at 90° | |
| <input type="checkbox"/> | 1619096 | Set at 96° | |
| <input type="checkbox"/> | 1619102 | Set at 102° | |
| <input type="checkbox"/> | 1619108 | Set at 108° | |
| <input type="checkbox"/> | 1619114 | Set at 114° | |
| <input type="checkbox"/> | 1619120 | Set at 120° | |
| <input type="checkbox"/> | 1619220 | Fold Down Back | |
| <input type="checkbox"/> | 1619221 | Non Fold Down Back | \$ 75 |

| Back Upholstery | | | |
|-------------------------------------------------------------------------------------------------------------------------|---------|---------------------------|--------------|
| <input type="checkbox"/> | 1620000 | No Back Upholstery | |
| Standard Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back. | | | |
| <input type="checkbox"/> | 1620018 | Set at 18" | |
| <input type="checkbox"/> | 1620019 | Set at 19" | |
| <input type="checkbox"/> | 1620020 | Set at 20" | |
| <input type="checkbox"/> | 1620021 | Set at 21" | |
| <input type="checkbox"/> | 1620022 | Set at 22" | |
| <input type="checkbox"/> | 1620023 | Set at 23" | |
| Adjustable Tension Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back. | | | |
| <input type="checkbox"/> | 1620321 | Set at 21" | E2611 \$ 295 |
| <input type="checkbox"/> | 1620322 | Set at 22" | |
| <input type="checkbox"/> | 1620323 | Set at 23" | |

Backrest

| | | | | |
|-------------------------------------|---------|-------------------------------------------|-------|----------------|
| <input type="checkbox"/> | 1621000 | No Backrest | | |
| nxt Optima HA Active Contour | | | | |
| <input type="checkbox"/> | 1621100 | 16"-20" width, 16"-20" height | | E2620 \$ 690 |
| <input type="checkbox"/> | 1621105 | 22" width, 16"-20" height | | E2621 \$ 720 |
| nxt Optima HA Plus Contour | | | | |
| <input type="checkbox"/> | 1621101 | 16"-20" width, 16"-20" height | | E2620 \$ 690 |
| <input type="checkbox"/> | 1621106 | 22" width, 16"-20" height | | E2621 \$ 720 |
| Matrx Elite | | | | |
| <input type="checkbox"/> | 1621110 | 14"-20" width, 10"-20" height | | E2620 \$ 825 |
| <input type="checkbox"/> | 1621112 | HD - 22" width, 16"-20" height | | E2621 \$ 1,095 |
| <input type="checkbox"/> | 1621116 | Deep - 14"-20" width, 12"-20" height | | E2620 \$ 950 |
| ROHO Agility | | | | |
| <input type="checkbox"/> | 1621120 | Min Contour 14"-20" width, 10"-18" height | | E2620 \$ 825 |
| <input type="checkbox"/> | 1621122 | Mid Contour 14"-20" width, 10"-18" height | | E2620 \$ 850 |
| <input type="checkbox"/> | 1621124 | Max Contour 14"-20" width, 10"-18" height | | E2620 \$ 950 |
| Varilite Icon | | | | |
| <input type="checkbox"/> | 1621130 | Tall - 14"-20" width, 14"-24" height | | E2615 \$ 850 |
| <input type="checkbox"/> | 1621132 | Deep - 14"-20" width, 17"-24" height | | E2620 \$ 895 |
| Backrest Size | | | | |
| | 1621400 | Specify Size (Width" x Height") | _____ | |

Tilt/Recline Options

| | | | | |
|----------------------------|---------|----------------------------|--------------------------|-----------------------------------|
| Attendant Tilt only | | | | |
| <input type="checkbox"/> | 1623011 | Trigger on Right Back Post | <input type="checkbox"/> | 1623012 Trigger on Left Back Post |

Headrest

| | | | | |
|-------------------------------------------------------------------------------------------------------------------|---------|-------------------------------------------|--------------------------|--------------------------|
| <input type="checkbox"/> | 1624000 | No Headrest | | |
| <input type="checkbox"/> | 1624010 | Headrest Mounting Plate | | \$ 75 |
| Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back. | | | | |
| <input type="checkbox"/> | 1624100 | nxt MultiFit Head Support Hardware | | E1028 \$ 235 |
| nxt MultiFit Head Support Pad | | | | |
| <input type="checkbox"/> | 1624105 | Regular (9" x 5") | <input type="checkbox"/> | 1624106 Large (10" x 5") |
| <input type="checkbox"/> | 1624107 | X-Large (14" x 5") | | E0955 \$ 245 |
| Stealth Headrest Hardware | | | | |
| <input type="checkbox"/> | 1624110 | Multi-Axis HMO475 | | E1028 \$ 150 |
| <input type="checkbox"/> | 1624111 | Link TWB480 | | E1028 \$ 325 |
| Stealth Headrest Pad | | | | |
| <input type="checkbox"/> | 1624115 | 10" width | <input type="checkbox"/> | 1624116 14" width |

Front Rigging

| | | | | | |
|--------------------------|---------|----------------------------------------------------|-----------------------------------------------------------------------|--------------------|--------------|
| <input type="checkbox"/> | 1625000 | No Front Rigging | | | |
| | | PDG Easy-Out™ Front Rigging | Effortless insertion and removal. No Hanger Bracket remains on frame. | | |
| <input type="checkbox"/> | 1625410 | Fixed 70° | | Available in March | \$ 95 |
| <input type="checkbox"/> | 1625420 | Swing In/Away Universal Angle Adjustable (0°-110°) | Set at (use increments of 10°): _____ | Available in April | \$ 195 |
| <input type="checkbox"/> | 1625430 | Swing In/Away Elevating Legrests | High knee pivot. Includes reinforced premium Calf Pad mount. | Available in April | E0990 \$ 495 |
| | | Standard Front Rigging | Legacy Design. Not compatible with Easy-Out™. | | |
| <input type="checkbox"/> | 1625140 | Swing Away 70° | | | |
| <input type="checkbox"/> | 1625160 | Swing Away 90° | Not compatible with 8" Casters. | | \$ 55 |
| <input type="checkbox"/> | 1625200 | Swing Away Elevating Legrests | | E0990 | \$ 375 |
| <input type="checkbox"/> | 1625300 | 1" Offset Hangers for Swing Away | | | \$ 75 |

Foot Plates

| | | | | | |
|--------------------------|---------|--------------------------------------|------------------------------|-------|--------|
| <input type="checkbox"/> | 1626010 | Composite | | | |
| <input type="checkbox"/> | 1626020 | Aluminum Angle Adjustable | | K0040 | \$ 155 |
| <input type="checkbox"/> | 1626050 | Interlocking Angle Adjustable | Seat Widths of 14"-21" only. | K0040 | \$ 195 |

Footrest/Legrest Accessories

| | | | | | |
|--------------------------|---------|--------------------------------------|----------------------------------------------------------------------|-------|--------|
| <input type="checkbox"/> | 1627010 | Swing Away 70° Calf Pad | Not compatible with 70° Gel Padding. | | \$ 135 |
| <input type="checkbox"/> | 1627020 | Swing Away 70° Gel Padding | Not compatible with 70° Swing Away Calf Pad, Easy-Out Front Rigging. | | \$ 100 |
| <input type="checkbox"/> | 1627030 | Swing Away ELR Gel Knee Pad | | | \$ 100 |
| <input type="checkbox"/> | 1627040 | Adjustable Length Heel Loops | | E0951 | \$ 50 |
| <input type="checkbox"/> | 1627050 | Bodypoint - Padded Calf Strap | | K0038 | \$ 105 |

Extension Tubes

| | | | | | |
|--------------------------|---------|-----------------------------|---------------------------------------------------------------------------------------------------------------|--|-------|
| <input type="checkbox"/> | 1528010 | Short | 13"-17" for Easy-Out SA Angle Adjustable 10"-13.5" for Easy Out SA ELR 13"-17" for Standard Front Rigging | | |
| <input type="checkbox"/> | 1528020 | Long | 16"-22.5" for Easy-Out SA Angle Adjustable 12.5"-19.5" for Easy-Out SA ELR 15"-19" for Standard Front Rigging | | \$ 75 |
| <input type="checkbox"/> | 1628200 | For 90° FR (11"-16") | Not compatible with Easy-Out, Standard 70° or ELR | | \$ 75 |

Positioning Straps

| | | | | | |
|--------------------------|---------|--------------------------------|--|-------|--------|
| <input type="checkbox"/> | 1629005 | Lap Belt Mount | | E1028 | \$ 55 |
| | | Bodypoint - Hip Belts | | | |
| <input type="checkbox"/> | 1629010 | Center-Pull Auto Style | | E0978 | \$ 105 |
| <input type="checkbox"/> | 1629015 | Center-Pull Auto Style, Padded | | E0978 | \$ 155 |
| <input type="checkbox"/> | 1629020 | 4-Point, Padded | | E0978 | \$ 175 |
| <input type="checkbox"/> | 1629030 | Bodypoint - Chest Strap | | E0960 | \$ 85 |



Miscellaneous Options

| | | | | | |
|--------------------------|---------|-------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------|
| <input type="checkbox"/> | 1630080 | Gel Hip Guides | Not compatible with Dual Post Arms. | E0956 & E1028 | \$ 325 |
| <input type="checkbox"/> | 1630110 | IV Pole | 15" STF and higher. | E1029 | \$ 295 |
| <input type="checkbox"/> | 1630120 | O2 Holder | 15" STF and higher. | E2208 | \$ 395 |
| <input type="checkbox"/> | 1630160 | Clear Lap Tray (Slide On) | Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray. | E0950 | \$ 490 |
| <input type="checkbox"/> | 1630170 | Adhesive Back Hook Strip, 2" x 30" | For Cushion mounting. | | \$ 10 |
| <input type="checkbox"/> | 1630301 | Fuze T20 to T50 Conversion Kit | | K0108 | \$ 995 |

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to
pdgmobility.com/products/fuze-T20/
to learn more about the
Fuze T20 20° Manual Tilt Wheelchair



  PDG Mobility