

## Customer and Order Information

Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address * <span style="float: right;">*required</span>	

## Frame Type

<input checked="" type="checkbox"/> 1501000 <b>Fuze T50 50° Manual Tilt Wheelchair</b>	E1161 \$ 2,995
--	----------------

## Frame Options

<b>Fixed Base and Seat Frame</b> Select Fixed Seat Width below.	
<input type="checkbox"/> 1502010 300 lbs Weight Capacity	
<input type="checkbox"/> 1502011 350 lbs Weight Capacity Not compatible with Transit Tie Down System. Requires Flat-free Rear Tires or Heavy Duty Spoke Rear Wheels.	
<input type="checkbox"/> 1502025 <b>Adjustable Base and Seat Frame</b> 250 lbs Weight Capacity. Select Adjustable Seat Width below.	\$ 100
<input type="checkbox"/> 1502300 <b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.	

## Seat Width Overall chair width is +7.5" for 12/16" RW and +8.5" for 20/22/24" RW

<b>Fixed Seat Widths</b> Order with Fixed Base and Seat Frame.		
<input type="checkbox"/> 1503013 13" Seat Width	<input type="checkbox"/> 1503014 14" Seat Width	<input type="checkbox"/> 1503015 15" Seat Width
<input type="checkbox"/> 1503016 16" Seat Width	<input type="checkbox"/> 1503017 17" Seat Width	<input type="checkbox"/> 1503018 18" Seat Width
<input type="checkbox"/> 1503019 19" Seat Width		
<input type="checkbox"/> 1503020 20" Seat Width	<input type="checkbox"/> 1503021 21" Seat Width	<input type="checkbox"/> 1503022 22" Seat Width <span style="float: right;">E2201 \$ 325</span>
<input type="checkbox"/> 1503023 23" Seat Width <span style="float: right;">E2201 \$ 325</span>		
<input type="checkbox"/> 1503024 24" Seat Width <span style="float: right;">E2202 \$ 525</span>		
<b>Adjustable Seat Widths</b> Order with Adjustable Base and Seat Frame.		
<input type="checkbox"/> 1503113 13"-15" Range - Set at 13"	<input type="checkbox"/> 1503114 13"-15" Range - Set at 14"	<input type="checkbox"/> 1503115 13"-15" Range - Set at 15"
<input type="checkbox"/> 1503215 15"-20" Range - Set at 15"	<input type="checkbox"/> 1503216 15"-20" Range - Set at 16"	<input type="checkbox"/> 1503217 15"-20" Range - Set at 17"
<input type="checkbox"/> 1503218 15"-20" Range - Set at 18"	<input type="checkbox"/> 1503219 15"-20" Range - Set at 19"	
<input type="checkbox"/> 1503220 15"-20" Range - Set at 20" <span style="float: right;">E2201 \$ 325</span>		
<input type="checkbox"/> 1503321 21"-22" Range - Set at 21" <span style="float: right;">E2201 \$ 325</span>	<input type="checkbox"/> 1503322 21"-22" Range - Set at 22"	

## Seat Depth

<b>Standard Seat Depths (Adjustable)</b>		
<input type="checkbox"/> 1504015 15" Seat Depth	<input type="checkbox"/> 1504016 16" Seat Depth	<input type="checkbox"/> 1504017 17" Seat Depth
<input type="checkbox"/> 1504018 18" Seat Depth	<input type="checkbox"/> 1504019 19" Seat Depth	
<input type="checkbox"/> 1504020 20" Seat Depth <span style="float: right;">E2203 \$ 475</span>		
<b>Custom Seat Depths</b> Requires 25" Push Bar Back Straight, Angle/Height Adjustable Push Bar. Not compatible with Reclining or Fold Down Back.		
<input type="checkbox"/> 1504021 21" Seat Depth <span style="float: right;">E2203 \$ 475</span>		
<input type="checkbox"/> 1504022 22" Seat Depth <span style="float: right;">E2204 \$ 750</span>		

## Frame Color

### PDG Premium Colors

<input type="checkbox"/> 1505000 Tracer Red	<input type="checkbox"/> 1505001 Gloss Black	<input type="checkbox"/> 1505002 Silver
<input type="checkbox"/> 1505003 Black Shimmer	<input type="checkbox"/> 1505004 Orange Fusion	<input type="checkbox"/> 1505005 Karl's Imagination
<input type="checkbox"/> 1505006 Blue Shimmer	<input type="checkbox"/> 1505007 Gloss White	<input type="checkbox"/> 1505008 Dark Titanium
<input type="checkbox"/> 1505009 Geyser Green	<input type="checkbox"/> 1505010 Lollypop Blue	<input type="checkbox"/> 1505011 Lollypop Red
<input type="checkbox"/> 1505012 Bubble Gum Pink	<input type="checkbox"/> 1505013 Berry Purple	<input type="checkbox"/> 1505014 Safety Yellow
<input type="checkbox"/> 1505200 <b>Custom Color</b> Contact PDG for specifications. <span style="float: right;">\$ 300</span>		

## Arm Type

<input type="checkbox"/> 1506000 <b>No Armrest</b>		
<input type="checkbox"/> 1506100 <b>Single Post Fixed Height (12")</b>		
<input type="checkbox"/> 1506114 <b>Single Post Adjustable Height (10"-14")</b>		E0973 \$ 255
<input type="checkbox"/> 1506200 <b>Dual Post Fixed Height (12")</b>	Full Length Armpad only.	
<b>Dual Post Adjustable Height</b>		E0973 \$ 255
<input type="checkbox"/> 1506214 10"-14"	<input type="checkbox"/> 1506216 12"-16"	
<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.		E0973 \$ 255
<input type="checkbox"/> 1506310 Set at 10"	<input type="checkbox"/> 1506311 Set at 11"	<input type="checkbox"/> 1506312 Set at 12"
<input type="checkbox"/> 1506313 Set at 13"	<input type="checkbox"/> 1506314 Set at 14"	

## Armrests

<input type="checkbox"/> 1507095 <b>Curved Armpads</b>	Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.	
<b>Regular Armpads</b>		
<input type="checkbox"/> 1507100 Desk Length - 2" x 10"	<input type="checkbox"/> 1507101 Full Length - 2" x 14"	
<input type="checkbox"/> 1507102 Full Length Premium Grade - 2" x 14"		\$ 35
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 115</span>		
<input type="checkbox"/> 1507110 Desk Length - 3" x 10"	<input type="checkbox"/> 1507111 Full Length - 3" x 12"	
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>		
<input type="checkbox"/> 1507120 Full Length - 2" x 12"	<input type="checkbox"/> 1507121 Full Length - 3.5" x 12"	
<b>Contour Arm Trough with Hand Pad</b>		E2209 \$ 125
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.		
<input type="checkbox"/> 1507130 Right (4.5" x 18")	<input type="checkbox"/> 1507131 Left (4.5" x 18")	

## Seat Pan

<input type="checkbox"/> 1508000 <b>No Seat Pan</b>	
<input type="checkbox"/> 1508001 <b>Solid Seat Pan</b>	E2231 \$ 120

## Cushion

<input type="checkbox"/> 1509000 <b>No Cushion</b>		
<b>nxt BioFit</b>		
<input type="checkbox"/> 1509101 14"-20" width, 16"-20" depth		E2607 \$ 595
<input type="checkbox"/> 1509102 22"-24" width, 16"-20" depth		E2608 \$ 625
<b>Matrix VI</b>		
<input type="checkbox"/> 1509111 14"-20" width, 16"-22" depth		E2607 \$ 550
<input type="checkbox"/> 1509112 22" width, 18"-22" depth		E2608 \$ 600
<input type="checkbox"/> 1509115 HD - 18"-24" width, 18"-22" depth		E2608 \$ 1,150

See next page for more Cushion options >

### ROHO Quadro Select High Profile

- 1509121 14"-20" width, 16"-22" depth
- 1509122 22"-24" width, 16"-22" depth

E2624 \$ 650  
E2625 \$ 750



### ROHO High Profile

- 1509125 14"-20" width, 16"-22" depth
- 1509126 22"-24" width, 16"-22" depth

E2622 \$ 650  
E2623 \$ 750

### Varilite Reflex

- 1509131 14"-20" width, 16"-20" depth



E2603 \$ 375

### Cushion Size

1509400 Specify Size (Width" x Depth") \_\_\_\_\_

## Seat Height

### Select desired Seat-to-Floor Height

#### Standard T50 Options

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
12"						18" 19"
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18" 19"
24"						18" 19"

RW	8" Casters & 8" Forks					
12"						18" 19" 20"
16"				16"		20"
20"				16"	17"	18"
22"					17"	18" 19"
24"						18" 19" 20"

RW	6" Soft Roll & Ultra Sport Forks					
16"				16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

#### Wide T50 Options

Increases the overall width of the Chair by 2".

- Standard
- Wide

Wheel size \_\_\_\_\_

Caster size \_\_\_\_\_

Fork size \_\_\_\_\_

STF height \_\_\_\_\_

RW	5" Casters					
12"	13"					
16"	13"	14"	15"			
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	

RW	6" Casters & 6" Forks					
16"		14"	15"			
20"		14"	15"	16"	17"	
22"		14"	15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	6" Casters & 8" Forks					
16"			15"			
20"			15"	16"	17"	
22"			15"	16"	17"	18"
24"			15"	16"	17"	18" 19"

RW	8" Casters & 8" Forks					
20"				16"	17"	
22"				16"	17"	18"
24"				16"	17"	18" 19"

## Casters

### Urethane

- 1511051 5" x 1.25" Caster size
- 1511061 6" x 1.25" Caster size
- 1511063 6" x 1.50" Caster size
- 1511081 8" x 1.25" Caster size
- 1511084 8" x 2.00" Caster size

\$ 85

\$ 135

See next page for more Casters options >

Pneumatic		
<input type="checkbox"/>	1511161 6" x 1.25" Caster size	<i>E2214 &amp; E2215</i> \$ 110
<input type="checkbox"/>	1511181 8" x 1.25" Caster size	<i>E2214 &amp; E2215</i> \$ 120
<input type="checkbox"/>	1511184 8" x 2.00" Caster size	<i>E2214 &amp; E2215</i> \$ 150
<input type="checkbox"/>	1511260 <b>6" Soft Roll</b>	\$ 95
<input type="checkbox"/>	1511400 <b>Ultra Sport Fork</b> Not compatible with 12" Rear Wheels. Requires 6" Soft Roll Caster.	<i>E1015</i> \$ 395

Rear Wheels		
Composite		
<input type="checkbox"/>	1512012 12" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1512022 22" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1512016 16" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1512024 24" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1512020 20" Wheel size	
Spoke		
<input type="checkbox"/>	1512120 20" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1512122 22" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1512124 24" Wheel size	\$ 75
<input type="checkbox"/>	1512224 <b>24" Heavy Duty Spoke</b>	\$ 350

Axle Type	
<input type="checkbox"/>	1513000 <b>Quick Release Rear Axles</b>
<input type="checkbox"/>	1513010 <b>Fixed Rear Axles</b>

Rear Tires		
<input type="checkbox"/>	1514000 <b>Urethane</b>	
<input type="checkbox"/>	1514010 <b>Pneumatic</b> Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	<i>E2211 &amp; E2212</i> \$ 60
<input type="checkbox"/>	1514020 <b>Pneumatic with Flat-Free Inserts</b> Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. Required for 350 lbs Weight Capacity.	<i>E2211 &amp; E2213</i> \$ 95

Hand Rims		
<input type="checkbox"/>	1515000 <b>No Hand Rims</b>	
<input type="checkbox"/>	1515010 <b>Standard</b> 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1515020 <b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only.	
		\$ 95

Wheel Locks		
<input type="checkbox"/>	1516005 <b>Foot Operated Wheel Lock</b> Not compatible with One Arm Drive, O2 Holder.	\$ 200
<input type="checkbox"/>	1516010 <b>Push to Lock</b>	
<input type="checkbox"/>	1516020 <b>Pull to Lock</b>	
<input type="checkbox"/>	1516030 <b>Attendant Wheel Lock (Push to Lock)</b> Not compatible with Wheel Lock Extensions.	\$ 85
<input type="checkbox"/>	1516040 <b>User &amp; Attendant Wheel Lock</b>	\$ 125
<input type="checkbox"/>	1516100 <b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	<i>E0961</i> \$ 60

Anti-Tippers		
<input checked="" type="checkbox"/>	1517010 <b>Rear Anti-Tippers</b> Required Safety Option.	<i>E0971</i> \$ 60

Back Post Style		
Push Bar Back		
<input type="checkbox"/>	1518010 Straight	
<input type="checkbox"/>	1518030 8° Bend	
		\$ 75

See next page for more Back Post Style options >

<b>Cane Handle Back</b> Not compatible with Headrest Mounting Plate.		
<input type="checkbox"/> 1518020	Straight	
<input type="checkbox"/> 1518040	8° Bend	\$ 75
<b>Back Post Height</b> Seat Pan to Handle.		
<input type="checkbox"/> 1518420	20" Height	<input type="checkbox"/> 1518425 25" Height

Back Options		
<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 1519010	Button on Outside	<input type="checkbox"/> 1519011 Button on Inside
<input type="checkbox"/> 1519050	<b>Removable Cane Handle Push Bar</b> Not compatible with Push Bar Back.	\$ 140
<b>Seat to Back Angle</b>		
<input type="checkbox"/> 1519090	Set at 90°	<input type="checkbox"/> 1519096 Set at 96°
<input type="checkbox"/> 1519108	Set at 108°	<input type="checkbox"/> 1519114 Set at 114°
<input type="checkbox"/> 1519102	Set at 102°	<input type="checkbox"/> 1519120 Set at 120°
<input type="checkbox"/> 1519220	<b>Fold Down Back</b>	
<input type="checkbox"/> 1519221	<b>Non Fold Down Back</b>	\$ 75

Back Upholstery		
<input type="checkbox"/> 1520000	<b>No Back Upholstery</b>	
<b>Standard</b> Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		<i>E2611 / E2612</i> \$ 155
<input type="checkbox"/> 1520018	Set at 18"	<input type="checkbox"/> 1520019 Set at 19"
<input type="checkbox"/> 1520021	Set at 21"	<input type="checkbox"/> 1520022 Set at 22"
<input type="checkbox"/> 1520020	Set at 20"	<input type="checkbox"/> 1520023 Set at 23"
<b>Adjustable Tension</b> Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		<i>E2611 / E2612</i> \$ 295
<input type="checkbox"/> 1520321	Set at 21"	<input type="checkbox"/> 1520322 Set at 22"
<input type="checkbox"/> 1520323	Set at 23"	

Backrest		
<input type="checkbox"/> 1521000	<b>No Backrest</b>	
<b>next Optima HA Active Contour</b>		
<input type="checkbox"/> 1521100	16"-20" width, 16"-20" height	<i>E2620</i> \$ 690
<input type="checkbox"/> 1521105	22"-24" width, 16"-20" height	<i>E2621</i> \$ 720
<b>next Optima HA Plus Contour</b>		
<input type="checkbox"/> 1521101	16"-20" width, 16"-20" height	<i>E2620</i> \$ 690
<input type="checkbox"/> 1521106	22"-24" width, 16"-20" height	<i>E2621</i> \$ 720
<b>Matrix Elite</b>		
<input type="checkbox"/> 1521110	14"-20" width, 10"-20" height	<i>E2620</i> \$ 825
<input type="checkbox"/> 1521112	HD - 22"-24" width, 16"-20" height	<i>E2621</i> \$ 1,095
<input type="checkbox"/> 1521116	Deep - 14"-20" width, 12"-20" height	<i>E2620</i> \$ 950
<b>ROHO Agility</b>		
<input type="checkbox"/> 1521120	Min Contour 14"-20" width, 10"-18" height	<i>E2620</i> \$ 825
<input type="checkbox"/> 1521122	Mid Contour 14"-20" width, 10"-18" height	<i>E2620</i> \$ 850
<input type="checkbox"/> 1521124	Max Contour 14"-20" width, 10"-18" height	<i>E2620</i> \$ 950
<b>Varilite Icon</b>		
<input type="checkbox"/> 1521130	Tall - 14"-20" width, 14"-24" height	<i>E2615</i> \$ 850
<input type="checkbox"/> 1521132	Deep - 14"-20" width, 17"-24" height	<i>E2620</i> \$ 895
<b>Backrest Size</b>		
<input type="checkbox"/> 1521400	Specify Size (Width" x Height")	

## Tilt/Recline Options

### Attendant Tilt only

- 1523011 Trigger on Right Back Post       1523012 Trigger on Left Back Post

### Tilt Limit

- 1523220 Set at 20°       1523225 Set at 25°       1523230 Set at 30°  
 1523235 Set at 35°       1523240 Set at 40°       1523245 Set at 45°  
 1523250 Set at 50°

- 1523410 **Reclining Backrest (Manual)** E1225 \$ 875  
 30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Tension Adjustable Upholstery. Requires 25" Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15", 16", 21", or 22" Seat Depth.

### Power Tilt 4° to 45° Tilt range, 250 lbs Weight Capacity. E1002 \$ 2,195

- Toggle Switch       1523421 Switch installed on Right side       1523422 Switch installed on Left side  
 Button Switch       1523431 Switch installed on Right side       1523432 Switch installed on Left side

## Headrest

- 1524000 **No Headrest**

- 1524010 **Headrest Mounting Plate** \$ 75  
 Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.

- 1524100 **nxt MultiFit Head Support Hardware** E1028 \$ 155

### nxt MultiFit Head Support Pad E0955 \$ 195

- 1524105 Regular (9" x 5")       1524106 Large (10" x 5")       1524107 X-Large (14" x 5")

### Stealth Headrest Hardware E1028 \$ 195

- 1524105 Multi-Axis HMO475  
 1524105 Link TWB480

### Stealth Headrest Pad E0955 \$ 200

- 1524105 10" width       1524106 14" width

## Front Rigging

- 1525000 **No Front Rigging**

- 1525001 **Sector Blocks only**

- 1525110 **Plug in 70°** Includes required Frame upgrade. \$ 300

- 1525120 **Plug in Adjustable Knee Angle** Includes required Frame upgrade. \$ 1,275

- 1525140 **Swing Away 70°**

- 1525160 **Swing Away 90°** Not compatible with 8" Casters. \$ 55

- 1525200 **Swing Away Elevating Legrest** E0990 \$ 375

- 1525300 **1" Offset Hangers for Swing Away** \$ 75

## Foot Plates

- 1526010 **Composite**

- 1526020 **Aluminum Angle Adjustable** K0040 \$ 155

- 1526040 **One Piece Flip-up** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395

- 1526050 **Interlocking Angle Adjustable** Seat Widths of 14"-21" only. K0040 \$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1527015	<b>Swing Away 70° Calf Pad</b> Not compatible with 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1527020	<b>Swing Away 70° Gel Padding</b> Not compatible with 70° Calf Pad.	\$ 100
<input type="checkbox"/>	1527030	<b>Swing Away ELR Gel Knee Pad</b>	\$ 100
<input type="checkbox"/>	1527040	<b>Adjustable Length Heel Loops</b>	E0951 \$ 50
<input type="checkbox"/>	1527050	<b>Bodypoint - Padded Calf Strap</b>	K0038 \$ 105
<b>Miller's Dynamic Footrest Mounts</b>			
<input type="checkbox"/>	1527060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1527065	For Swing Away 70° Front Rigging Not available on 13" STF Height.	\$ 380

Extension Tubes			
<b>Regular Extension Tubes</b> Not compatible with Swing Away 90°.			
<input type="checkbox"/>	1528010	Standard (13"-17")	
<input type="checkbox"/>	1528020	Long (15"-19")	\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1528110	Standard (13"-17")	
<input type="checkbox"/>	1528120	Long (15"-19")	\$ 100
<input type="checkbox"/>	1528200	<b>90° Extension Tubes (11"-16")</b> Not compatible with 70° or ELR.	\$ 75

Positioning Straps			
<input type="checkbox"/>	1529005	<b>Lap Belt Mount</b>	E1028 \$ 55
<b>Bodypoint - Hip Belts</b>			
<input type="checkbox"/>	1529010	Center-Pull Auto Style	E0978 \$ 105
<input type="checkbox"/>	1529015	Center-Pull Auto Style, Padded	E0978 \$ 155
<input type="checkbox"/>	1529020	4-Point, Padded	E0978 \$ 175
<input type="checkbox"/>	1529030	<b>Bodypoint - Chest Strap</b>	E0960 \$ 85



Miscellaneous Options			
<b>One-Arm Drive</b> E0958 \$ 3,500			
Requires 24" Composite Rear Wheel. Forward Wheel Position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks.			
<input type="checkbox"/>	1530010	Installed on Right side	<input type="checkbox"/>
<input type="checkbox"/>	1530011	Installed on Left side	
<input type="checkbox"/>	1530060	<b>Stealth Flip Down Abductor</b>	E0957 \$ 495
<input type="checkbox"/>	1530070	<b>QR Bracket - Accessory Mount</b>	E1028 \$ 200
<input type="checkbox"/>	1530080	<b>Gel Hip Guides</b> Not compatible with Dual Post Arms.	E0956 & E1028 \$ 325
<input type="checkbox"/>	1530110	<b>IV Pole</b>	E1029 \$ 295
<input type="checkbox"/>	1530120	<b>O2 Holder</b> May add to overall width of the Chair. Not compatible with Foot Operated Wheel Lock.	E2208 \$ 395
<input type="checkbox"/>	1530150	<b>Ventilator Tray</b> 16" minimum Seat Width. Reduces ground clearance. Chair may be unstable when unoccupied. Not compatible with O2 Holder.	E1029 \$ 295
<input type="checkbox"/>	1530152	<b>Ventilator Bracket</b> Compatible with LTV1 (laptop) ventilator.	E1029 \$ 295
<input type="checkbox"/>	1530154	<b>Ventilator Battery Tray</b>	E1029 \$ 95
16" minimum Seat Width. Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with STF restrictions.			
<input type="checkbox"/>	1530160	<b>Clear Lap Tray (Slide On)</b> Not compatible with 8"- 11" Dual Post Arms, Desk Length Pads, Gel Arm Pads. Seat Widths of 15"-20" only.	E0950 \$ 490
<input type="checkbox"/>	1530170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10
<input type="checkbox"/>	1530300	<b>Fuze T50 to T20 Conversion Kit</b>	K0108 \$ 995

## Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.  
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to  
[pdgmobility.com/products/fuze-T50](http://pdgmobility.com/products/fuze-T50)  
to learn more about the  
**Fuze T50 50° Manual Tilt Wheelchair**



  PDG Mobility