

Customer and Order Information

Purchase Order *	Tag	v.061518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 1501000	Fuze T50 50° Manual Tilt Wheelchair	E1161	\$ 2,995
---	--	-------	----------

Frame Options

Fixed Base and Seat Frame Select Fixed Seat Width below.			
<input type="checkbox"/> 1502010	300 lbs Weight Capacity		
<input type="checkbox"/> 1502011	350 lbs Weight Capacity	Not compatible with Transit Tie Down System. Requires Flat-free Rear Tires or Heavy Duty Spoke Rear Wheels.	
<input type="checkbox"/> 1502025	Adjustable Base and Seat Frame	250 lbs Weight Capacity. Select Adjustable Seat Width below.	\$ 100
<input type="checkbox"/> 1502300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.	

Seat Width Overall chair width is +7.5" for 12/16" RW and +8.5" for 20/22/24" RW

Fixed Seat Widths Order with Fixed Base and Seat Frame.			
<input type="checkbox"/> 1503015	15" Seat Width	<input type="checkbox"/> 1503016	16" Seat Width
<input type="checkbox"/> 1503017	17" Seat Width	<input type="checkbox"/> 1503018	18" Seat Width
<input type="checkbox"/> 1503019	19" Seat Width	<input type="checkbox"/> 1503020	20" Seat Width
<input type="checkbox"/> 1503021	21" Seat Width	<input type="checkbox"/> 1503022	22" Seat Width
<input type="checkbox"/> 1503023	23" Seat Width		E2201 \$ 325
<input type="checkbox"/> 1503024	24" Seat Width		E2202 \$ 525
Adjustable Seat Widths Order with Adjustable Base and Seat Frame.			
<input type="checkbox"/> 1503112	12"-15" Range - Set at 12"	<input type="checkbox"/> 1503113	12"-15" Range - Set at 13"
<input type="checkbox"/> 1503114	12"-15" Range - Set at 14"	<input type="checkbox"/> 1503115	12"-15" Range - Set at 15"
<input type="checkbox"/> 1503215	15"-20" Range - Set at 15"	<input type="checkbox"/> 1503216	15"-20" Range - Set at 16"
<input type="checkbox"/> 1503217	15"-20" Range - Set at 17"	<input type="checkbox"/> 1503218	15"-20" Range - Set at 18"
<input type="checkbox"/> 1503219	15"-20" Range - Set at 19"	<input type="checkbox"/> 1503220	15"-20" Range - Set at 20"
<input type="checkbox"/> 1503321	21"-22" Range - Set at 21"	<input type="checkbox"/> 1503322	21"-22" Range - Set at 22"
			E2201 \$ 325
			E2201 \$ 325

Seat Depth

Standard Seat Depths (Adjustable)			
<input type="checkbox"/> 1504015	15" Seat Depth	<input type="checkbox"/> 1504016	16" Seat Depth
<input type="checkbox"/> 1504017	17" Seat Depth	<input type="checkbox"/> 1504018	18" Seat Depth
<input type="checkbox"/> 1504019	19" Seat Depth	<input type="checkbox"/> 1504020	20" Seat Depth
			E2203 \$ 475
Custom Seat Depths Requires 25" Push Bar Back Straight, Angle/Height Adjustable Push Bar. Not compatible with Reclining or Fold Down Back.			
<input type="checkbox"/> 1504021	21" Seat Depth		E2203 \$ 475
<input type="checkbox"/> 1504022	22" Seat Depth		E2204 \$ 750

Frame Color

PDG Premium Colors

<input type="checkbox"/> 1505000 Tracer Red	<input type="checkbox"/> 1505001 Gloss Black	<input type="checkbox"/> 1505002 Silver
<input type="checkbox"/> 1505003 Black Shimmer	<input type="checkbox"/> 1505004 Orange Fusion	<input type="checkbox"/> 1505005 Karl's Imagination
<input type="checkbox"/> 1505006 Blue Shimmer	<input type="checkbox"/> 1505007 Gloss White	<input type="checkbox"/> 1505008 Dark Titanium
<input type="checkbox"/> 1505009 Geyser Green	<input type="checkbox"/> 1505010 Lollypop Blue	<input type="checkbox"/> 1505011 Lollypop Red
<input type="checkbox"/> 1505012 Bubble Gum Pink	<input type="checkbox"/> 1505013 Berry Purple	<input type="checkbox"/> 1505014 Safety Yellow
<input type="checkbox"/> 1505200 Custom Color Contact PDG for specifications. \$ 300		

Arm Type

<input type="checkbox"/> 1506000 No Armrest		
<input type="checkbox"/> 1506100 Single Post Fixed Height (12")		
<input type="checkbox"/> 1506114 Single Post Adjustable Height (10"-14")		E0973 \$ 255
<input type="checkbox"/> 1506200 Dual Post Fixed Height (12") Full Length Armpad only.		
Dual Post Adjustable Height		E0973 \$ 255
<input type="checkbox"/> 1506214 10"-14"	<input type="checkbox"/> 1506216 12"-16"	
Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.		E0973 \$ 255
<input type="checkbox"/> 1506310 Set at 10"	<input type="checkbox"/> 1506311 Set at 11"	<input type="checkbox"/> 1506312 Set at 12"
<input type="checkbox"/> 1506313 Set at 13"	<input type="checkbox"/> 1506314 Set at 14"	

Armrests

<input type="checkbox"/> 1507095 Curved Armpads Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.		
Regular Armpads		
<input type="checkbox"/> 1507100 Desk Length - 2" x 10"	<input type="checkbox"/> 1507101 Full Length - 2" x 14"	
<input type="checkbox"/> 1507102 Full Length Premium Grade - 2" x 14"		\$ 35
Wide Armpads May interfere with wheels in some tilt configurations.		
<input type="checkbox"/> 1507110 Desk Length - 3" x 10"	<input type="checkbox"/> 1507111 Full Length - 3" x 12"	\$ 115
Gel Armpads May interfere with wheels in some tilt configurations.		
<input type="checkbox"/> 1507120 Full Length - 2" x 12"	<input type="checkbox"/> 1507121 Full Length - 3.5" x 12"	\$ 125
Contour Arm Trough with Hand Pad		E2209 \$ 125
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.		
<input type="checkbox"/> 1507130 Right (4.5" x 18")	<input type="checkbox"/> 1507131 Left (4.5" x 18")	

Seat Pan

<input type="checkbox"/> 1508000 No Seat Pan		
<input type="checkbox"/> 1508001 Solid Aluminum Seat Pan		E2231 \$ 120

Cushion

<input type="checkbox"/> 1509000 No Cushion		
nxt BioFit		
<input type="checkbox"/> 1509101 14"-20" width, 16"-20" depth		E2607 \$ 595
<input type="checkbox"/> 1509102 22"-24" width, 16"-20" depth		E2608 \$ 625
Matrix VI		
<input type="checkbox"/> 1509111 14"-20" width, 16"-22" depth		E2607 \$ 550
<input type="checkbox"/> 1509112 22" width, 18"-22" depth		E2608 \$ 600
<input type="checkbox"/> 1509115 HD - 18"-24" width, 18"-22" depth		E2608 \$ 1,150

See next page for more Cushion options >

ROHO Quadro Select High Profile			
<input type="checkbox"/>	1509121 14"-20" width, 16"-22" depth	E2624	\$ 650
<input type="checkbox"/>	1509122 22"-24" width, 16"-22" depth	E2625	\$ 750
ROHO High Profile			
<input type="checkbox"/>	1509125 14"-20" width, 16"-22" depth	E2622	\$ 650
<input type="checkbox"/>	1509126 22"-24" width, 16"-22" depth	E2623	\$ 750
Varilite Reflex			
<input type="checkbox"/>	1509131 14"-20" width, 16"-20" depth	E2603	\$ 375
Cushion Size			
1509400	Specify Size (Width" x Depth")	_____	



Seat Height

Select desired Seat-to-Floor Height

Standard T50 Options

STF Height is approximate and can vary by as much as 1" depending on Rear Wheel type and pre-tilt setting. To resolve, STF can be adjusted down by removing the pre-tilt bumper.

RW	5" Casters					
12"	13"	14"				
16"		14"	15"	16"		
20"				16"	17"	
22"					17"	

RW	6" Casters & 6" Forks					
12"		14"				18"
16"		14"	15"	16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

RW	6" Casters & 8" Forks					
12"						18" 19"
16"			15"	16"		
20"				16"	17"	18"
22"					17"	18" 19"
24"						18" 19"

RW	8" Casters & 8" Forks					
12"					18"	19" 20"
16"				16"		20"
20"				16"	17"	18"
22"					17"	18" 19"
24"						18" 19" 20"

RW	6" Soft Roll & Ultra Sport Forks					
16"				16"		
20"				16"	17"	18"
22"					17"	18"
24"						18"

Wide T50 Options

Increases the overall width of the Chair by 2".

RW	5" Casters					
12"	13"					
16"	13"	14"	15"			
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	

RW	6" Casters & 6" Forks					
16"		14"	15"			
20"		14"	15"	16"	17"	
22"		14"	15"	16"	17"	18"
24"			15"	16"	17"	18"

RW	6" Casters & 8" Forks					
16"			15"			
20"			15"	16"	17"	
22"			15"	16"	17"	18"
24"			15"	16"	17"	18" 19"

RW	8" Casters & 8" Forks					
20"				16"	17"	
22"				16"	17"	18"
24"				16"	17"	18" 19"

- Standard
 Wide

Wheel size _____

Caster size _____

Fork size _____

STF height _____

Casters

Urethane

<input type="checkbox"/>	1511051 5" x 1.25" Caster size		
<input type="checkbox"/>	1511061 6" x 1.25" Caster size		
<input type="checkbox"/>	1511063 6" x 1.50" Caster size		\$ 85
<input type="checkbox"/>	1511081 8" x 1.25" Caster size		
<input type="checkbox"/>	1511084 8" x 2.00" Caster size		\$ 135

See next page for more Casters options >

Pneumatic				
<input type="checkbox"/>	1511161	6" x 1.25" Caster size	E2214 & E2215	\$ 110
<input type="checkbox"/>	1511181	8" x 1.25" Caster size	E2214 & E2215	\$ 120
<input type="checkbox"/>	1511184	8" x 2.00" Caster size	E2214 & E2215	\$ 150
<input type="checkbox"/>	1511260	6" Soft Roll		\$ 95
<input type="checkbox"/>	1511400	Ultra Sport Fork Not compatible with 12" Rear Wheels. Requires 6" Soft Roll Caster.	E1015	\$ 395

Rear Wheels										
Composite										
<input type="checkbox"/>	1512012	12" Wheel size	<input type="checkbox"/>	1512016	16" Wheel size	<input type="checkbox"/>	1512020	20" Wheel size		
<input type="checkbox"/>	1512022	22" Wheel size	<input type="checkbox"/>	1512024	24" Wheel size					
Spoke										\$ 75
<input type="checkbox"/>	1512120	20" Wheel size	<input type="checkbox"/>	1512122	22" Wheel size	<input type="checkbox"/>	1512124	24" Wheel size		
<input type="checkbox"/>	1512224	24" Heavy Duty Spoke								\$ 350

Axle Type	
<input type="checkbox"/>	1513000 Quick Release Rear Axles
<input type="checkbox"/>	1513010 Fixed Rear Axles

Rear Tires				
<input type="checkbox"/>	1514000	Urethane		
<input type="checkbox"/>	1514010	Pneumatic Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212	\$ 60
<input type="checkbox"/>	1514020	Pneumatic with Flat-Free Inserts Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels. Required for 350 lbs Weight Capacity.	E2211 & E2213	\$ 95

Hand Rims		
<input type="checkbox"/>	1515000 No Hand Rims	
<input type="checkbox"/>	1515010 Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1515020 Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 135

Wheel Locks			
<input type="checkbox"/>	1516005	Foot Operated Wheel Lock Not compatible with One Arm Drive.	\$ 200
<input type="checkbox"/>	1516010	Push to Lock	
<input type="checkbox"/>	1516020	Pull to Lock	
<input type="checkbox"/>	1516030	Attendant Wheel Lock (Push to Lock) Not compatible with Wheel Lock Extensions.	\$ 85
<input type="checkbox"/>	1516040	User & Attendant Wheel Lock	\$ 125
<input type="checkbox"/>	1516100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 60

Anti-Tippers			
<input checked="" type="checkbox"/>	1517010	Rear Anti-Tippers Required Safety Option.	E0971 \$ 60

Back Post Style		
Push Bar Back		
<input type="checkbox"/>	1518010 Straight	
<input type="checkbox"/>	1518030 8° Bend	\$ 75

See next page for more Back Post Style options >

Cane Handle Back Not compatible with Headrest Mounting Plate.		
<input type="checkbox"/> 1518020	Straight	
<input type="checkbox"/> 1518040	8° Bend	\$ 75
Back Post Height Seat Pan to Handle.		
<input type="checkbox"/> 1518420	20" Height	<input type="checkbox"/> 1518425 25" Height

Back Options		
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 1519010	Button on Outside	<input type="checkbox"/> 1519011 Button on Inside
<input type="checkbox"/> 1519050	Removable Cane Handle Push Bar Not compatible with Push Bar Back.	\$ 140
Seat to Back Angle		
<input type="checkbox"/> 1519090	Set at 90°	<input type="checkbox"/> 1519096 Set at 96°
<input type="checkbox"/> 1519108	Set at 108°	<input type="checkbox"/> 1519114 Set at 114°
<input type="checkbox"/> 1519102	Set at 102°	<input type="checkbox"/> 1519120 Set at 120°
<input type="checkbox"/> 1519220	Fold Down Back	
<input type="checkbox"/> 1519221	Non Fold Down Back	\$ 75

Back Upholstery		
<input type="checkbox"/> 1520000	No Back Upholstery	
Standard Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		E2611 / E2612 \$ 155
<input type="checkbox"/> 1520018	Set at 18"	<input type="checkbox"/> 1520019 Set at 19"
<input type="checkbox"/> 1520021	Set at 21"	<input type="checkbox"/> 1520022 Set at 22"
<input type="checkbox"/> 1520020	Set at 20"	<input type="checkbox"/> 1520023 Set at 23"
Adjustable Tension Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.		E2611 / E2612 \$ 295
<input type="checkbox"/> 1520321	Set at 21"	<input type="checkbox"/> 1520322 Set at 22"
<input type="checkbox"/> 1520323	Set at 23"	

Backrest		
<input type="checkbox"/> 1521000	No Backrest	
next Optima HA Active Contour		
<input type="checkbox"/> 1521100	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/> 1521105	22"-24" width, 16"-20" height	E2621 \$ 720
next Optima HA Plus Contour		
<input type="checkbox"/> 1521101	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/> 1521106	22"-24" width, 16"-20" height	E2621 \$ 720
Matrx Elite		
<input type="checkbox"/> 1521110	14"-20" width, 10"-20" height	E2620 \$ 825
<input type="checkbox"/> 1521112	HD - 22"-24" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/> 1521116	Deep - 14"-20" width, 12"-20" height	E2620 \$ 950
ROHO Agility		
<input type="checkbox"/> 1521120	Min Contour 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/> 1521122	Mid Contour 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/> 1521124	Max Contour 14"-20" width, 10"-18" height	E2620 \$ 950
Varilite Icon		
<input type="checkbox"/> 1521130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/> 1521132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
Backrest Size		
1521400	Specify Size (Width" x Height")	

Tilt/Recline Options

Attendant Tilt only

- 1523011 Trigger on Right Back Post 1523012 Trigger on Left Back Post

Tilt Limit



- 1523220 Set at 20° 1523225 Set at 25° 1523230 Set at 30°
 1523235 Set at 35° 1523240 Set at 40° 1523245 Set at 45°
 1523250 Set at 50°

- 1523410 **Reclining Backrest (Manual)** E1225 \$ 875
 30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back, Tension Adjustable Upholstery. Requires 25" Angle/Height Adjustable Push Bar with 8" Bend. Not available with 15", 16", 21", or 22" Seat Depth.

Power Tilt 0° to 40° tilt range. Includes 2 batteries for off-chair charging.

- 1523400 250 lbs Weight Capacity E1002 \$ 2,195
 1523401 350 lbs Weight Capacity Requires Fixed Base and Seat Frame E1002 \$ 2,495
- Toggle Switch 1523421 Installed on Right side 1523422 Installed on Left side
 Button Switch 1523431 Installed on Right side 1523432 Installed on Left side
 Egg Switch 1523440 Not installed
- Optional with Egg Switch: 1523445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. E1028 \$ 625
- Specify Headrest to be used: _____

Headrest

- 1524000 **No Headrest**
- 1524010 **Headrest Mounting Plate** \$ 75
 Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.
- 1524100 **nxt MultiFit Head Support Hardware** E1028 \$ 235
- nxt MultiFit Head Support Pad**  E0955 \$ 245
- 1524105 Regular (9" x 5") 1524106 Large (10" x 5") 1524107 X-Large (14" x 5")
- Stealth Headrest Hardware**
- 1524110 Multi-Axis HMO475 E1028 \$ 150
 1524111 Link TWB480 E1028 \$ 325
- Stealth Headrest Pad**  E0955 \$ 245
- 1524115 10" width 1524116 14" width

Front Rigging

- 1525000 **No Front Rigging** 1525001 **Sector Blocks only**
- PDG Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame.
- 1525410 Fixed 70° \$ 95
 1525420 Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____ \$ 195
 1525430 Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount. E0990 \$ 495
- Standard Front Rigging** Legacy Design. Not compatible with Easy-Out™.
- 1525140 Swing Away 70°
 1525160 Swing Away 90° Not compatible with 8" Casters. \$ 55
 1525200 Swing Away Elevating Legrests E0990 \$ 375
 1525300 1" Offset Hangers for Swing Away \$ 75

Foot Plates

- 1526010 **Composite**
- 1526020 **Aluminum Angle Adjustable** K0040 \$ 155
- 1526040 **One Piece Flip-up** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
- 1526050 **Interlocking Angle Adjustable** Seat Widths of 14"-21" only. K0040 \$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1527015	Swing Away 70° Calf Pad Not compatible with 70° Gel Padding.		\$	135
<input type="checkbox"/>	1527020	Swing Away 70° Gel Padding Not compatible with 70° Swing Away Calf Pad, Easy-Out Front Rigging.		\$	100
<input type="checkbox"/>	1527030	Swing Away ELR Gel Knee Pad			100
<input type="checkbox"/>	1527040	Adjustable Length Heel Loops	E0951	\$	50
<input type="checkbox"/>	1527050	Bodypoint - Padded Calf Strap	K0038	\$	105
<input type="checkbox"/>	1527060	Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.		\$	190

Extension Tubes

<input type="checkbox"/>	1528010	Short <u>14"-18"</u> for Easy-Out Fixed & SA Angle Adjustable <u>11"-15"</u> for Easy-Out SA ELR <u>13"-17"</u> for Standard Front Rigging			
<input type="checkbox"/>	1528020	Long <u>18"-24"</u> for Easy-Out Fixed & SA Angle Adjustable <u>15"-21"</u> for Easy-Out SA ELR <u>15"-19"</u> for Standard Front Rigging		\$	75
<input type="checkbox"/>	1628200	For 90° FR (11"-16") Not compatible with Easy-Out, Standard 70° or ELR		\$	75

Positioning Straps

<input type="checkbox"/>	1529005	Lap Belt Mount	E1028	\$	55
Bodypoint - Hip Belts					
<input type="checkbox"/>	1529010	Center-Pull Auto Style	E0978	\$	105
<input type="checkbox"/>	1529015	Center-Pull Auto Style, Padded	E0978	\$	155
<input type="checkbox"/>	1529020	4-Point, Padded	E0978	\$	175
<input type="checkbox"/>	1529030	Bodypoint - Chest Strap	E0960	\$	85



Miscellaneous Options

		One-Arm Drive Requires 24" Composite Rear Wheels and forward wheel position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks.	E0958	\$	3,500
<input type="checkbox"/>	1530010	Installed on Right side	<input type="checkbox"/>	1530011	Installed on Left side
<input type="checkbox"/>	1530060	Stealth Flip Down Abductor	E0957	\$	495
<input type="checkbox"/>	1530070	QR Bracket - Accessory Mount	E1028	\$	200
<input type="checkbox"/>	1530080	Gel Hip Guides Not compatible with Dual Post Arms. Reduces available Seat Width by 2".	E0956 & E1028	\$	325
<input type="checkbox"/>	1530110	IV Pole	E1029	\$	295
<input type="checkbox"/>	1530120	O2 Holder May add to overall width of the Chair.	E2208	\$	395
<input type="checkbox"/>	1530150	Ventilator Tray Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 16". Not compatible with O2 Holder, Power Tilt.	E1029	\$	295
<input type="checkbox"/>	1530152	Ventilator Bracket Compatible with LTV1 (laptop) ventilator.	E1029	\$	295
<input type="checkbox"/>	1530154	Ventilator Battery Tray 16" minimum Seat Width. Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with STF restrictions.	E1029	\$	95
<input type="checkbox"/>	1530160	Clear Lap Tray (Slide On) Not compatible with 8"- 11" Dual Post, Curved Armpads, Contour Arm Trough. Seat Widths of 15"-20" only. Premium Grade or Gel Armpads recommended to ensure snug fit of the Tray.	E0950	\$	490
<input type="checkbox"/>	1530170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.		\$	10
<input type="checkbox"/>	1530300	Fuze T50 to T20 Conversion Kit	K0108	\$	995

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to
pdgmobility.com/products/fuze-T50/
to learn more about the
Fuze T50 50° Manual Tilt Wheelchair

