

Cushion

1. Select the type ↓	2. Select the model ↓
-------------------------	--------------------------

General Use	<p>Stealth Simplicity 350 lbs weight capacity.</p> <p><input type="checkbox"/> 4009141 16"-20" width, 16"-20" depth E2601 \$ 230</p> <p><input type="checkbox"/> 4009142 22"-24" width, 18"-20" depth E2602 \$ 290</p>
--------------------	---

Positioning	<p>ROHO AirLITE 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009120 14"-20" width, 14"-18" depth E2605 \$ 375</p>
--------------------	--

Skin Protection & Positioning	<p>Stealth Spectrum Foam 650 lbs weight capacity.</p> <p><input type="checkbox"/> 4009143 16"-20" width, 16"-24" depth E2607 \$ 550</p> <p><input type="checkbox"/> 4009144 22"-32" width, 18"-28" depth E2608 \$ 625</p>
	<p>Stealth Tru-Comfort 2 450 lbs weight capacity.</p> <p><input type="checkbox"/> 4009146 16"-20" width, 16"-24" depth E2607 \$ 490</p> <p><input type="checkbox"/> 4009147 22" width, 18"-24" depth E2608 \$ 590</p>
	<p>nxt BioFit 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009101 14"-20" width, 16"-20" depth E2607 \$ 555</p> <p><input type="checkbox"/> 4009102 22"-24" width, 16"-20" depth E2608 \$ 585</p>
	<p>Matrx VI 300 lbs weight capacity.</p> <p><input type="checkbox"/> 4009111 14"-20" width, 15"-20" depth E2607 \$ 680</p> <p><input type="checkbox"/> 4009112 22" width, 18"-20" depth E2608 \$ 745</p>
	<p>Matrx VI Heavy Duty 600 lbs weight capacity.</p> <p><input type="checkbox"/> 4009115 18"-21" width, 18"-26" depth E2607 \$ 1,050</p> <p><input type="checkbox"/> 4009116 22"-30" width, 18"-26" depth E2608 \$ 1,250</p>

Skin Protection & Adjustable	<p>ROHO Single Valve High Profile No weight limit if cushion is properly sized.</p> <p><input type="checkbox"/> 4009125 14"-20" width, 16"-22" depth E2622 \$ 840</p> <p><input type="checkbox"/> 4009126 22"-24" width, 16"-22" depth E2623 \$ 960</p>
---	--

Skin Protection, Positioning & Adjustable	<p>ROHO Quadro Select High Profile No weight limit if cushion is properly sized.</p> <p><input type="checkbox"/> 4009121 14"-20" width, 16"-22" depth E2624 \$ 840</p> <p><input type="checkbox"/> 4009122 22"-24" width, 16"-22" depth E2625 \$ 960</p>
--	---

Extra-Wide	<p>PDG Solar Basic</p> <p>High resiliency polyfoam with 2" medium-density base and 1" lower-density top.</p>	<p><input type="checkbox"/> 4009200 20"-28" width, 16"-24" depth \$ 395</p> <p><input type="checkbox"/> 4009201 29"-36" width, 16"-24" depth \$ 495</p> <p><input type="checkbox"/> 4009202 37"-48" width, 16"-24" depth \$ 595</p>
	<p>PDG Solar Contour</p> <p>Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups.</p>	<p><input type="checkbox"/> 4009210 20"-28" width, 16"-24" depth \$ 695</p> <p><input type="checkbox"/> 4009211 29"-36" width, 16"-24" depth \$ 795</p> <p><input type="checkbox"/> 4009212 37"-48" width, 16"-24" depth \$ 895</p>
	<p>PDG Solar Pressure Reduction</p> <p>All features of the Contour, plus added pincore latex rubber and containment foam outline.</p>	<p><input type="checkbox"/> 4009220 20"-28" width, 16"-24" depth \$ 1,095</p> <p><input type="checkbox"/> 4009221 29"-36" width, 16"-24" depth \$ 1,195</p> <p><input type="checkbox"/> 4009222 37"-48" width, 16"-24" depth \$ 1,295</p>

No weight limit if cushion is properly sized. Custom sizes available.

Custom-made foam cushion with fluid-proof yet breathable fabric cover.

3. Specify the size →	<p>Cushion Size Sizes based on availability from the cushion manufacturer.</p> <p>Specify Size (Width" x Depth") _____</p>
-----------------------	---

Backrest

1. Select the type ↓	2. Select the contour depth ↓	3. Select the model ↓
-------------------------	----------------------------------	--------------------------

Positioning, posterior-lateral	2" depth	<input type="checkbox"/> 4021141 Stealth ADI AL Active 250 lbs weight capacity. 14"-20" width, 10"-20" height <i>E2615</i> \$ 790
	4" depth	<input type="checkbox"/> 4021142 Stealth ADI AL Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height <i>E2615</i> \$ 875
	6" depth	<input type="checkbox"/> 4021143 Stealth ADI AL Extra Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height <i>E2615</i> \$ 960

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> 4021110 Matrx Elite 300 lbs weight capacity. 12"-20" width, 10"-20" height <i>E2620</i> \$ 955
		<input type="checkbox"/> 4021111 Matrx Elite Heavy Duty 500 lbs weight capacity. 21" width, 16"/20" height <i>E2620</i> \$ 1,195
		<input type="checkbox"/> 4021112 22"-24" width, 16"/20" height <i>E2621</i> \$ 1,295
		<input type="checkbox"/> 4021114 25"-30" width, 16"/20" height <i>E2621</i> \$ 1,995
	3.5" to 4" depth	<input type="checkbox"/> 4021120 ROHO Agility Min 300 lbs weight capacity. 14"-20" width, 10"-18" height <i>E2620</i> \$ 900
		<input type="checkbox"/> 4021100 nxt Xtend 300 lbs weight capacity. 16"-20" width, 16"-20" height <i>E2620</i> \$ 635
		<input type="checkbox"/> 4021105 22"-24" width, 16"-20" height <i>E2621</i> \$ 660
		<input type="checkbox"/> 4021122 ROHO Agility Mid 300 lbs weight capacity. 14"-20" width, 10"-18" height <i>E2620</i> \$ 930
	6" to 7" depth	<input type="checkbox"/> 4021101 nxt Xtend Deep 300 lbs weight capacity. 16"-20" width, 16"-20" height <i>E2620</i> \$ 635
		<input type="checkbox"/> 4021106 22"-24" width, 16"-20" height <i>E2621</i> \$ 660
		<input type="checkbox"/> 4021116 Matrx Elite Deep 300 lbs weight capacity. 14"-20" width, 12"-20" height <i>E2620</i> \$ 1,145
		<input type="checkbox"/> 4021117 22"-24" width, 12"-20" height <i>E2621</i> \$ 1,295
<input type="checkbox"/> 4021124 ROHO Agility Max 300 lbs weight capacity. 14"-20" width, 10"-18" height <i>E2620</i> \$ 1,085		

4. Specify the size →	Backrest Size Sizes based on availability from the backrest manufacturer. Specify Size (Width" x Height") _____
-----------------------	---

Pediatric Seating System

4030400 **Diego Pediatric Seating System**
 Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10	12	14	16
Width	(25.40)	(30.48)	(35.56)	(40.64)
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓

A Head System Head Pad & Hardware

B FlipLock Configurable as a between cases or full width set up.

C Contoured Seat System Base, Cushion & Cover

Multi-Axis Hardware The IMOSOL included as part of a package for the Diego Head System, allows growth adjustability and on-demand configurations.

Contoured Back System Base, Cushion & Cover

Coolcore Covers Our Diego System uses the Coolcore cover technology which adapt all-year-round for a more comfortable experience.

Laterals Modular Laterals either Fixed or Swing Away

Contoured Seat System Base, Cushion, Cover & Hardware

Width Growth Adapted from our TEB System, the back wings allow for a pediatric growth adjustment.

Back Hardware Cover Easy to remove back cover for quick on-demand adjustments and a more aesthetic look.

diego a Stealth Product

WC-20 Certified Some requirements need to be met for safety purposes. Read warning 17!

Headrest

Stealth Comfort Plus Headrest	Stealth Hardware		
	<input type="checkbox"/> 4024110	Multi-Axis HMO475	E1028 \$ 225
	<input type="checkbox"/> 4024111	Link TWB480	E1028 \$ 350
	Comfort Plus Pad		
<input type="checkbox"/> 4024116	Large (10" x 6")	<input type="checkbox"/> 4024117 X-Large (14" x 6")	E0955 \$ 275
nxt MultiFit Head Support	nxt Hardware		
	<input type="checkbox"/> 4024100		E1028 \$ 285
	Multi-Adjustable Pad		
	<input type="checkbox"/> 4024105	Contoured (9" x 5")	E0955 \$ 275
<input type="checkbox"/> 4024106	Large (10" x 5")	<input type="checkbox"/> 4024107 X-Large (14" x 5")	

Positioning Belts & Straps

Bodypoint	Hip Belts		
	<input type="checkbox"/> 4029010	Center-Pull Auto Style	E0978 \$ 115
	<input type="checkbox"/> 4029015	Center-Pull Auto Style, Padded	E0978 \$ 215
	<input type="checkbox"/> 4029020	4-Point, Padded	E0978 \$ 235
	<input type="checkbox"/> 4029030	Chest Strap	E0960 \$ 165

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.

PDG proudly supplies on-chair seating and positioning products from these manufacturers:

