

Customer and Order Information

Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address *	*required

Frame Type

<input checked="" type="checkbox"/> 1901000 Stellar GL 45° Manual Tilt Wheelchair	E1161 \$ 2,995
--	----------------

Frame Options

<input checked="" type="checkbox"/> 1902000 Standard Base and Seat Frame 250 lbs Weight Capacity.	
<input type="checkbox"/> 1902300 Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths			
<input type="checkbox"/> 1903014 14" Seat Width	<input type="checkbox"/> 1903015 15" Seat Width	<input type="checkbox"/> 1903016 16" Seat Width	
<input type="checkbox"/> 1903017 17" Seat Width	<input type="checkbox"/> 1903018 18" Seat Width	<input type="checkbox"/> 1903019 19" Seat Width	
<input type="checkbox"/> 1903020 20" Seat Width	<input type="checkbox"/> 1903021 21" Seat Width	<input type="checkbox"/> 1903022 22" Seat Width	E2201 \$ 325

Seat Depth

Standard Seat Depths (Adjustable)			
<input type="checkbox"/> 1904016 16" Seat Depth	<input type="checkbox"/> 1904017 17" Seat Depth	<input type="checkbox"/> 1904018 18" Seat Depth	
<input type="checkbox"/> 1904019 19" Seat Depth			
<input type="checkbox"/> 1904020 20" Seat Depth			E2203 \$ 475
Custom Seat Depths			
<input type="checkbox"/> 1904021 21" Seat Depth			E2203 \$ 475
<input type="checkbox"/> 1904022 22" Seat Depth			E2204 \$ 750

Frame Color

PDG Premium Colors			
<input type="checkbox"/> 1905000 Tracer Red	<input type="checkbox"/> 1905001 Gloss Black	<input type="checkbox"/> 1905002 Silver	
<input type="checkbox"/> 1905003 Black Shimmer	<input type="checkbox"/> 1905004 Orange Fusion	<input type="checkbox"/> 1905005 Karl's Imagination	
<input type="checkbox"/> 1905006 Blue Shimmer	<input type="checkbox"/> 1905007 Gloss White	<input type="checkbox"/> 1905008 Dark Titanium	
<input type="checkbox"/> 1905009 Geyser Green	<input type="checkbox"/> 1905010 Lollypop Blue	<input type="checkbox"/> 1905011 Lollypop Red	
<input type="checkbox"/> 1905012 Bubble Gum Pink	<input type="checkbox"/> 1905013 Berry Purple	<input type="checkbox"/> 1905014 Safety Yellow	
<input type="checkbox"/> 1905200 Custom Color Contact PDG for specifications.			\$ 300

Arm Type

<input type="checkbox"/> 1906000 No Armrest	
<input type="checkbox"/> 1906100 Single Post Fixed Height (12") Full Length Armpad only.	
Single Post Adjustable Height	
<input type="checkbox"/> 1906111 8"-11"	<input type="checkbox"/> 1906114 10"-14"
<input type="checkbox"/> 1906200 Dual Post Fixed Height (12") Full Length Armpad only.	E0973 \$ 255

See next page for more Arm Type options >

Dual Post Adjustable Height			E0973	\$	255
<input type="checkbox"/> 1906211	8"-11"	<input type="checkbox"/> 1906214	10"-14"	<input type="checkbox"/> 1906216	12"-16"
Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.			E0973	\$	255
<input type="checkbox"/> 1906310	Set at 10"	<input type="checkbox"/> 1906311	Set at 11"	<input type="checkbox"/> 1906312	Set at 12"
<input type="checkbox"/> 1906313	Set at 13"	<input type="checkbox"/> 1906314	Set at 14"		

Armrests					
Regular Armpads					
<input type="checkbox"/> 1907100	Desk Length - 2" x 10"	<input type="checkbox"/> 1907101	Full Length - 2" x 14"		
<input type="checkbox"/> 1907102	Full Length Premium Grade - 2" x 14"			\$	35
Wide Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 1907110	Desk Length - 3" x 10"	<input type="checkbox"/> 1907111	Full Length - 3" x 12"		\$ 115
Gel Armpads May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 1907120	Full Length - 2" x 12"	<input type="checkbox"/> 1907121	Full Length - 3.5" x 12"		\$ 125
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.					
<input type="checkbox"/> 1907130	Right (4.5" x 18")	<input type="checkbox"/> 1907131	Left (4.5" x 18")	E2209	\$ 125
<input type="checkbox"/> 1907220	Armrest Wideners for Dual Post Armrests	Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.			\$ 375
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.					
<input type="checkbox"/> 1907310	For Single Post Armrests	Available for 10"-14" Height range only.			K0108 \$ 110
<input type="checkbox"/> 1907320	For Dual Post Armrests			K0108	\$ 165
Armrest Extensions Not compatible with Flip-up or Single Post Armrests.					
<input type="checkbox"/> 1907400	3" Extension - Forward	<input type="checkbox"/> 1907401	3" Extension - Rearward		\$ 60

Seat Pan					
<input type="checkbox"/> 1908000	No Seat Pan				
<input type="checkbox"/> 1908001	Solid Seat Pan			E2231	\$ 120

Cushion					
<input type="checkbox"/> 1909000	No Cushion				
nxt BioFit					
<input type="checkbox"/> 1909101	14"-20" width, 16"-20" depth			E2607	\$ 595
<input type="checkbox"/> 1909102	22" width, 16"-20" depth			E2608	\$ 625
Matrx VI					
<input type="checkbox"/> 1909111	14"-20" width, 16"-22" depth			E2607	\$ 550
<input type="checkbox"/> 1909112	22" width, 18"-22" depth			E2608	\$ 600
<input type="checkbox"/> 1909115	HD - 18"-22" width, 18"-22" depth			E2608	\$ 1,150
ROHO Quadro Select High Profile					
<input type="checkbox"/> 1909121	14"-20" width, 16"-22" depth			E2624	\$ 650
<input type="checkbox"/> 1909122	22" width, 16"-22" depth			E2625	\$ 750
ROHO High Profile					
<input type="checkbox"/> 1909125	14"-20" width, 16"-22" depth			E2622	\$ 650
<input type="checkbox"/> 1909126	22" width, 16"-22" depth			E2623	\$ 750
Varilite Reflex					
<input type="checkbox"/> 1909131	14"-20" width, 16"-20" depth			E2603	\$ 375
Cushion Size					
1909400	Specify Size (Width" x Depth")				

Seat Height

Select desired Seat-to-Floor Height STF Height is adjustable within the range of the selected Rear Wheel size.

RW	6" Casters					
12"	13"	14"	15"	16"		
16"			15"	16"	17"	18"
20"	13"	14"	15"	16"		
22"		14"	15"	16"	17"	
24"			15"	16"	17"	18"

RW	8" Casters					
16"			15"	16"	17"	18"
24"			15"	16"	17"	18"

Wheel size _____

Caster size _____

STF height _____

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1911061 6" x 1.25" Caster size
- 1911063 6" x 1.50" Caster size \$ 85
- 1911081 8" x 1.25" Caster size
- 1911084 8" x 2.00" Caster size \$ 135

Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1911161 6" x 1.25" Caster size *E2214 & E2215* \$ 110
- 1911181 8" x 1.25" Caster size *E2214 & E2215* \$ 120
- 1911184 8" x 2.00" Caster size *E2214 & E2215* \$ 150
- 1911260 **6" Soft Roll** \$ 95
- 1911400 **Ultra Sport Fork** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. *E1015* \$ 395

Rear Wheels

Composite

- 1912012 12" Wheel size
- 1912016 16" Wheel size
- 1912020 20" Wheel size
- 1912022 22" Wheel size
- 1912024 24" Wheel size

Spoke

- 1912120 20" Wheel size
- 1912122 22" Wheel size
- 1912124 24" Wheel size
- 1912224 **24" Heavy Duty Spoke** \$ 350

Axle Type

- 1913000 **Quick Release Rear Axles**
- 1913010 **Fixed Rear Axles**

Rear Tires

- 1914000 **Urethane**
- 1914010 **Pneumatic** Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. *E2211 & E2212* \$ 60
- 1914020 **Pneumatic with Flat-Free Inserts** *E2211 & E2213* \$ 95
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.

Hand Rims

- 1915000 **No Hand Rims**
- 1915010 **Standard** 20", 22", and 24" Rear Wheels only.
- 1915020 **Plastic Coated** 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks			
<input type="checkbox"/>	1916010	Push to Lock	
<input type="checkbox"/>	1916020	Pull to Lock	
<input type="checkbox"/>	1916030	Attendant Wheel Lock (Push to Lock)	\$ 85
<input type="checkbox"/>	1916040	User & Attendant Wheel Lock Available only on 20", 22", or 24" Rear Wheels.	\$ 125
<input type="checkbox"/>	1916100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 60
<input type="checkbox"/>	1916200	Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 350

Anti-Tippers			
<input checked="" type="checkbox"/>	1917010	Rear Anti-Tippers Required Safety Option.	E0971 \$ 60
<input type="checkbox"/>	1917020	Front Anti-Tippers	E0971 \$ 300
<input type="checkbox"/>	1917030	Heavy Duty Rear Anti-Tipper Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.	E0971 \$ 600

Back Post Style			
<input type="checkbox"/>	1918010	Push Bar Back	
<input type="checkbox"/>	1918020	Cane Handle Back Not compatible with Headrest Mounting Plate.	
Back Post Height Seat Pan to Handle			
<input type="checkbox"/>	1918425	25" Height	<input type="checkbox"/>
			1918431 31" Height

Back Options			
<input type="checkbox"/>		Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150
<input type="checkbox"/>	1919010	Button on Outside	<input type="checkbox"/>
			1919011 Button on Inside
Seat to Back Angle			
<input type="checkbox"/>	1919090	Set at 90°	<input type="checkbox"/>
			1919098 Set at 98°
			<input type="checkbox"/>
			1919106 Set at 106°
<input type="checkbox"/>	1919114	Set at 114°	<input type="checkbox"/>
			1919122 Set at 122°
<input type="checkbox"/>	1919210	10" Backrest Extension Not compatible with Angle/Height Adjustable Push Bar Handle.	E0966 \$ 135
<input type="checkbox"/>	1919220	Fold Down Back Not compatible with Reclining Back.	
<input type="checkbox"/>	1919230	Spring Back Suspension Not compatible with Reclining Back.	\$ 105
<input type="checkbox"/>	1919240	Dynamic Rocker Back Assembly Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.	K0108 \$ 2,995

Back Upholstery			
<input type="checkbox"/>	1920000	No Back Upholstery	
<input type="checkbox"/>		Standard Back height measured from Seat Pan to top of Upholstery.	E2611 \$ 155
<input type="checkbox"/>	1920021	Set at 21"	<input type="checkbox"/>
			1920022 Set at 22"
<input type="checkbox"/>		Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 \$ 155
<input type="checkbox"/>	1920118	Set at 18"	<input type="checkbox"/>
			1920119 Set at 19"
			<input type="checkbox"/>
			1920120 Set at 20"
<input type="checkbox"/>		Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 \$ 155
<input type="checkbox"/>	1920224	Set at 24"	<input type="checkbox"/>
			1920225 Set at 25"
			<input type="checkbox"/>
			1920226 Set at 26"
<input type="checkbox"/>		Adjustable Tension Back height measured from Seat Pan to top of Upholstery.	E2611 \$ 295
<input type="checkbox"/>	1920321	Set at 21"	<input type="checkbox"/>
			1920322 Set at 22"
			<input type="checkbox"/>
			1920323 Set at 23"

Backrest			
<input type="checkbox"/>	1921000	No Backrest	
nxt Optima HA Active Contour			
<input type="checkbox"/>	1921100	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1921105	22" width, 16"-20" height	E2621 \$ 720
nxt Optima HA Plus Contour			
<input type="checkbox"/>	1921101	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1921106	22" width, 16"-20" height	E2621 \$ 720
Matrx Elite			
<input type="checkbox"/>	1921110	14"-20" width, 10"-20" height	E2620 \$ 825
<input type="checkbox"/>	1921112	HD - 22" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/>	1921116	Deep - 14"-20" width, 12"-20" height	E2620 \$ 950
ROHO Agility			
<input type="checkbox"/>	1921120	Min Contour 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1921122	Mid Contour 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1921124	Max Contour 14"-20" width, 10"-18" height	E2620 \$ 950
Varilite Icon			
<input type="checkbox"/>	1921130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1921132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
Backrest Size			
	1921400	Specify Size (Width" x Height")	_____

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1923011	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1923012	Trigger on Left Back Post	
Tilt Limit			
<input type="checkbox"/>	1923220	Set at 20°	<input type="checkbox"/>
<input type="checkbox"/>	1923230	Set at 30°	<input type="checkbox"/>
<input type="checkbox"/>	1923240	Set at 40°	
<input type="checkbox"/>	1923410	Reclining Backrest (Manual)	\$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.			

Headrest			
<input type="checkbox"/>	1924000	No Headrest	
<input type="checkbox"/>	1924010	Headrest Mounting Plate	\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	1924100	nxt MultiFit Head Support Hardware	E1028 \$ 155
nxt MultiFit Head Support Pad			
<input type="checkbox"/>	1924105	Regular (9" x 5")	<input type="checkbox"/>
<input type="checkbox"/>	1924106	Large (10" x 5")	<input type="checkbox"/>
<input type="checkbox"/>	1924107	X-Large (14" x 5")	
<input type="checkbox"/>	1924105	Stealth Headrest Hardware	E1028 \$ 195
<input type="checkbox"/>	1924105	Multi-Axis HMO475	
<input type="checkbox"/>	1924105	Link TWB480	
Stealth Headrest Pad			
<input type="checkbox"/>	1924105	10" width	<input type="checkbox"/>
<input type="checkbox"/>	1924106	14" width	

Front Rigging			
<input type="checkbox"/>	1925000	No Front Rigging	
<input type="checkbox"/>	1925001	Sector Blocks only	
<input type="checkbox"/>	1925110	Plug In 70°	
<input type="checkbox"/>	1925120	Plug In Adjustable Knee Angle	\$ 975
<input type="checkbox"/>	1925130	Swing Away 60°	
<input type="checkbox"/>	1925140	Swing Away 70°	
<input type="checkbox"/>	1925200	Swing Away Elevating Legrest	E0990 \$ 375
<input type="checkbox"/>	1925300	1" Offset Hangers for Swing Away	\$ 75

Foot Plates			
<input type="checkbox"/>	1926010	Composite	
<input type="checkbox"/>	1926020	Aluminum Angle Adjustable	K0040 \$ 155
<input type="checkbox"/>	1926030	One Piece Fixed Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1926040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1926050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.	K0040 \$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1927010	Swing Away 60° & 70° Calf Pad Not compatible with 60° & 70° Gel Padding.	\$ 135
<input type="checkbox"/>	1927020	Swing Away 60° & 70° Gel Padding Not compatible with 60° & 70° Swing Away Calf Pad.	\$ 100
<input type="checkbox"/>	1927030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1927040	Adjustable Length Heel Loops	E0951 \$ 50
<input type="checkbox"/>	1927050	Bodypoint - Padded Calf Strap	K0038 \$ 105
Miller's Dynamic Footrest Mounts			
<input type="checkbox"/>	1927060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1927065	For Swing Away 70° Front Rigging Not available on 13" STF Height.	\$ 380

Extension Tubes			
Regular Extension Tubes			
<input type="checkbox"/>	1928010	Standard (13"-17")	
<input type="checkbox"/>	1928020	Long (15"-19")	\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1928110	Standard (13"-17")	
<input type="checkbox"/>	1928120	Long (15"-19")	\$ 100

Positioning Straps			
Bodypoint - Hip Belts			
<input type="checkbox"/>	1929010	Center-Pull Auto Style	E0978 \$ 105
<input type="checkbox"/>	1929015	Center-Pull Auto Style, Padded	E0978 \$ 155
<input type="checkbox"/>	1929020	4-Point, Padded	E0978 \$ 175
<input type="checkbox"/>	1929030	Bodypoint - Chest Strap	E0960 \$ 85



Miscellaneous Options

<input type="checkbox"/>	1930080	Gel Hip Guides		<i>E0956 & E1028</i>	\$ 325
<input type="checkbox"/>	1930110	IV Pole	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	<i>E1029</i>	\$ 295
<input type="checkbox"/>	1930120	O2 Holder	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	<i>E2208</i>	\$ 395
<input type="checkbox"/>	1930170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.		\$ 10

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar-GL to learn more about the **Stellar GL 45° Manual Tilt Wheelchair**



PDG Mobility