

Customer and Order Information

Purchase Order *	Tag	v.031518
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 2101000	Stellar HD 45° Bariatric Manual Tilt Wheelchair	Extended delivery lead times may apply.	E1161	\$ 2,995
---	--	---	-------	----------

Frame Options

Heavy Duty Base and Seat Frame				
<input type="checkbox"/> 2102045	450 lbs Weight Capacity. 45° of Tilt		K0108	\$ 600
<input type="checkbox"/> 2102060	600 lbs Weight Capacity. 25° of Tilt		K0108	\$ 900
User Weight Range Required for selection of Gas Struts (Force measured in Newton - N)				
<input type="checkbox"/> 2102212	Up to 250 lbs (2 x 1200 N)	<input type="checkbox"/> 2102308	250 lbs to 350 lbs (3 x 800 N)	
<input type="checkbox"/> 2102310	350 lbs to 425 lbs (3 x 1000 N)	<input type="checkbox"/> 2102312	425 lbs to 600 lbs (3 x 1200 N)	

Seat Width Overall chair width is +11"

Standard Seat Widths					
<input type="checkbox"/> 2103022	22" Seat Width		E2201	\$ 325	
<input type="checkbox"/> 2103024	24" Seat Width	<input type="checkbox"/> 2103026	26" Seat Width	E2202 \$ 525	
<input type="checkbox"/> 2103028	28" Seat Width	<input type="checkbox"/> 2103030	30" Seat Width	<input type="checkbox"/> 2103032	32" Seat Width K0108 \$ 925

Seat Depth

Standard Seat Depths (Adjustable)				
<input type="checkbox"/> 2104018	18" Seat Depth	<input type="checkbox"/> 2104019	19" Seat Depth	\$ 325
<input type="checkbox"/> 2104020	20" Seat Depth	<input type="checkbox"/> 2104021	21" Seat Depth	E2203 \$ 475
<input type="checkbox"/> 2104022	22" Seat Depth			E2204 \$ 750

Frame Color

PDG Premium Colors					
<input type="checkbox"/> 2105000	Tracer Red	<input type="checkbox"/> 2105001	Gloss Black	<input type="checkbox"/> 2105002	Silver
<input type="checkbox"/> 2105003	Black Shimmer	<input type="checkbox"/> 2105004	Orange Fusion	<input type="checkbox"/> 2105005	Karl's Imagination
<input type="checkbox"/> 2105006	Blue Shimmer	<input type="checkbox"/> 2105007	Gloss White	<input type="checkbox"/> 2105008	Dark Titanium
<input type="checkbox"/> 2105009	Geyser Green	<input type="checkbox"/> 2105010	Lollypop Blue	<input type="checkbox"/> 2105011	Lollypop Red
<input type="checkbox"/> 2105012	Bubble Gum Pink	<input type="checkbox"/> 2105013	Berry Purple	<input type="checkbox"/> 2105014	Safety Yellow
<input type="checkbox"/> 2105200	Custom Color	Contact PDG for specifications.		\$ 300	

Arm Type

<input type="checkbox"/> 2106000	No Armrest	
<input type="checkbox"/> 2106214	Dual Post Adjustable Height (10"-14")	E0973 \$ 255

Armrests

Regular Armpads				
<input type="checkbox"/> 2107100	Desk Length - 2" x 10"	<input type="checkbox"/> 2107101	Full Length - 2" x 14"	
<input type="checkbox"/> 2107102	Full Length Premium Grade - 2" x 14"			\$ 35

See next page for more Armrests options >

Wide Armpads May interfere with wheels in some tilt configurations.		\$ 115
<input type="checkbox"/> 2107110	Desk Length - 3" x 10"	<input type="checkbox"/> 2107111 Full Length - 3" x 12"
Gel Armpads May interfere with wheels in some tilt configurations.		\$ 125
<input type="checkbox"/> 2107120	Full Length - 2" x 12"	<input type="checkbox"/> 2107121 Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.		E2209 \$ 125
<input type="checkbox"/> 2107130	Right (4.5" x 18")	<input type="checkbox"/> 2107131 Left (4.5" x 18")
<input type="checkbox"/> 2107320	Clothing Guards	\$ 165

Seat Pan		
<input checked="" type="checkbox"/> 2108001	Solid Seat Pan	E2231 \$ 120

Cushion		
<input type="checkbox"/> 2109000	No Cushion	
<input type="checkbox"/> 2109115	Matrx VI Heavy Duty 22"-30" width, 18"-22" depth	E2608 \$ 1,150
		
3" Basic Foam (Nylon Cover)		
<input type="checkbox"/> 2109200	Up to 28" width	E2602 \$ 485
<input type="checkbox"/> 2109201	Above 28" width	E2602 \$ 605
3" Contour Foam (Nylon Cover)		
<input type="checkbox"/> 2109210	Up to 28" width	E2602 \$ 835
<input type="checkbox"/> 2109211	Above 28" width	E2602 \$ 1,090
3" Contour Pressure Reduction (Nylon Cover)		
<input type="checkbox"/> 2109220	Up to 28" width	E2602 \$ 1,485
<input type="checkbox"/> 2109221	Above 28" width	E2602 \$ 1,900
Spare Cushion Cover		
<input type="checkbox"/> 2109230	Basic	\$ 390
<input type="checkbox"/> 2109231	Contour	\$ 490
Cushion Size		
2109400	Specify Size (Width" x Depth") _____	

Seat Height		
Seat to Floor Height		
<input type="checkbox"/> 2110016	Set at 16" STF	<input type="checkbox"/> 2110017 Set at 17" STF
		<input type="checkbox"/> 2110018 Set at 18" STF

Casters		
<input checked="" type="checkbox"/> 2111050	5" Caster With Heavy Duty Steel Fork	

Rear Wheels		
<input checked="" type="checkbox"/> 2112224	24" Heavy Duty Spoke	\$ 350

Axle Type		
<input checked="" type="checkbox"/> 2113010	Fixed Rear Axles	

Rear Tires		
<input checked="" type="checkbox"/> 2114000	Urethane	

Hand Rims		
<input checked="" type="checkbox"/> 2115010	Standard	

Wheel Locks			
<input type="checkbox"/>	2116010	Push to Lock	
<input type="checkbox"/>	2116020	Pull to Lock	
<input type="checkbox"/>	2116100	Wheel Lock Extensions	E0961 \$ 60

Anti-Tippers			
<input checked="" type="checkbox"/>	2117010	Rear Anti-Tippers	Required Safety Option. E0971 \$ 60

Back Post Style			
<input checked="" type="checkbox"/>	2118010	Push Bar Back	
		Back Post Height	Seat Pan to Handle
<input type="checkbox"/>	2118425	25" Height	<input type="checkbox"/> 2118431 31" Height

Back Options			
<input checked="" type="checkbox"/>	2119090	Seat to Back Angle set at 90°	
<input type="checkbox"/>	2119240	Dynamic Rocker Back Assembly	K0108 \$ 2,995
Not compatible with 30° Recline. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.			

Back Upholstery			
<input type="checkbox"/>	2120000	No Back Upholstery	
		Adjustable Tension	Back height measured from Seat Pan to top of Upholstery. E2612 \$ 295
<input type="checkbox"/>	2120321	Set at 21"	<input type="checkbox"/> 2120322 Set at 22" <input type="checkbox"/> 2120323 Set at 23"

Backrest			
<input type="checkbox"/>	2121000	No Backrest	
		Matrx Elite	
<input type="checkbox"/>	2121112	HD - 22"-26" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/>	2121114	HD - 26"-30" width, 16"-20" height	E2621 \$ 1,595
		Backrest Size	
1021400	Specify Size (Width" x Height")		_____

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	2123011	Trigger on Right Back Post	<input type="checkbox"/> 2123012 Trigger on Left Back Post
		Tilt Limit	Maximum Tilt Limit for 600 lbs Weight Capacity is 25°
<input type="checkbox"/>	2123220	Set at 20°	<input type="checkbox"/> 2123225 Set at 25° <input type="checkbox"/> 2123230 Set at 30°
<input type="checkbox"/>	2123235	Set at 35°	<input type="checkbox"/> 2123240 Set at 40° <input type="checkbox"/> 2123245 Set at 45°
<input type="checkbox"/>	2123410	Reclining Backrest (Manual)	30° maximum recline. Not compatible with Tension Adjustable Upholstery. E1225 \$ 875

Headrest			
<input type="checkbox"/>	2124000	No Headrest	
<input type="checkbox"/>	2124010	Headrest Mounting Plate	Not required if mounting to a solid Back. \$ 75
<input type="checkbox"/>	2124100	nxt MultiFit Head Support Hardware	E1028 \$ 235
		nxt MultiFit Head Support Pad	E0955 \$ 245
<input type="checkbox"/>	2124105	Regular (9" x 5")	<input type="checkbox"/> 2124106 Large (10" x 5") <input type="checkbox"/> 2124107 X-Large (14" x 5")

See next page for more Headrest options >

Stealth Headrest Hardware			
<input type="checkbox"/>	2124110 Multi-Axis HMO475	<i>E1028</i>	\$ 150
<input type="checkbox"/>	2124111 Link TWB480	<i>E1028</i>	\$ 325
Stealth Headrest Pad			
<input type="checkbox"/>	2124115 10" width	<i>E0955</i>	\$ 200
<input type="checkbox"/>	2124116 14" width		



Front Rigging			
<input type="checkbox"/>	2125001 Sector Blocks only		
<input type="checkbox"/>	2125130 Swing Away 60°		
<input type="checkbox"/>	2125200 Swing Away Elevating Legrest	<i>E0990</i>	\$ 375

Foot Plates			
<input checked="" type="checkbox"/>	2126020 Aluminum Angle Adjustable	<i>K0040</i>	\$ 165

Footrest/Legrest Accessories			
<input type="checkbox"/>	2127030 Swing Away ELR Gel Knee Pad		\$ 100
<input type="checkbox"/>	2127040 Adjustable Length Heel Loops	<i>E0951</i>	\$ 50

Extension Tubes			
<input checked="" type="checkbox"/>	2128010 Regular Extension Tubes (13"-17")		

Positioning Straps			
Bodypoint - Hip Belts			
<input type="checkbox"/>	2129010 Center-Pull Auto Style	<i>E0978</i>	\$ 105
<input type="checkbox"/>	2129015 Center-Pull Auto Style, Padded	<i>E0978</i>	\$ 155



Miscellaneous Options			
<input type="checkbox"/>	2130120 O2 Holder	<i>E2208</i>	\$ 395
<input type="checkbox"/>	2130170 Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.		\$ 10

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar-HD/ to learn more about the Stellar HD 45° Manual Tilt Wheelchair

