

Customer and Order Information		
Purchase Order *	Tag	v.072318
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type		
<b>Stellar LEAP Manual Posterior &amp; Anterior Tilt Wheelchair</b> Extended delivery lead times may apply.		
<input checked="" type="checkbox"/>	2201000 Manual Wheelchair with 20° Posterior Tilt	E1161 \$ 2,995
<input type="checkbox"/>	2201101 25° Anterior Tilt with Linked Back 30° Anterior Tilt available with reduced Posterior Tilt range of 15°.	K0108 \$ 1,250
<input type="checkbox"/>	2201102 30° Dynamic Back Recline	E1225 \$ 500

Frame Options		
<input checked="" type="checkbox"/>	2202000 <b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.	
<input type="checkbox"/>	2202300 <b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.	

Seat Width		
Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW		
<b>Standard Seat Widths</b>		
<input type="checkbox"/>	2203014 14" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203015 15" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203016 16" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203017 17" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203018 18" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203019 19" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203020 20" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203021 21" Seat Width	<input type="checkbox"/>
<input type="checkbox"/>	2203022 22" Seat Width	E2201 \$ 325





Seat Depth		
<b>Standard Seat Depths (Adjustable)</b>		
<input type="checkbox"/>	2204016 16" Seat Depth	<input type="checkbox"/>
<input type="checkbox"/>	2204017 17" Seat Depth	<input type="checkbox"/>
<input type="checkbox"/>	2204018 18" Seat Depth	<input type="checkbox"/>
<input type="checkbox"/>	2204019 19" Seat Depth	<input type="checkbox"/>
<input type="checkbox"/>	2204020 20" Seat Depth	E2203 \$ 475
<b>Custom Seat Depths</b>		
<input type="checkbox"/>	2204021 21" Seat Depth	E2203 \$ 475
<input type="checkbox"/>	2204022 22" Seat Depth	E2204 \$ 750

Frame Color		
<b>PDG Premium Colors</b>		
<input type="checkbox"/>	2205000 Tracer Red	<input type="checkbox"/>
<input type="checkbox"/>	2205001 Gloss Black	<input type="checkbox"/>
<input type="checkbox"/>	2205002 Silver	<input type="checkbox"/>
<input type="checkbox"/>	2205003 Black Shimmer	<input type="checkbox"/>
<input type="checkbox"/>	2205004 Orange Fusion	<input type="checkbox"/>
<input type="checkbox"/>	2205005 Karl's Imagination	<input type="checkbox"/>
<input type="checkbox"/>	2205006 Blue Shimmer	<input type="checkbox"/>
<input type="checkbox"/>	2205007 Gloss White	<input type="checkbox"/>
<input type="checkbox"/>	2205008 Dark Titanium	<input type="checkbox"/>
<input type="checkbox"/>	2205009 Geyser Green	<input type="checkbox"/>
<input type="checkbox"/>	2205010 Lollypop Blue	<input type="checkbox"/>
<input type="checkbox"/>	2205011 Lollypop Red	<input type="checkbox"/>
<input type="checkbox"/>	2205012 Bubble Gum Pink	<input type="checkbox"/>
<input type="checkbox"/>	2205013 Berry Purple	<input type="checkbox"/>
<input type="checkbox"/>	2205014 Safety Yellow	<input type="checkbox"/>
<input type="checkbox"/>	2205200 <b>Custom Color</b> Contact PDG for specifications.	\$ 300

Arm Type			
<input type="checkbox"/>	2206000	<b>No Armrest</b>	
<input type="checkbox"/>	2206100	<b>Single Post Fixed Height (12")</b> Full Length Armpad only.	
		<b>Single Post Adjustable Height</b>	E0973 \$ 255
<input type="checkbox"/>	2206111	8"-11"	<input type="checkbox"/>
<input type="checkbox"/>	2206114	10"-14"	
<input type="checkbox"/>	2206200	<b>Dual Post Fixed Height (12")</b> Full Length Armpad only.	
		<b>Dual Post Adjustable Height</b>	E0973 \$ 255
<input type="checkbox"/>	2206211	8"-11"	<input type="checkbox"/>
<input type="checkbox"/>	2206214	10"-14"	<input type="checkbox"/>
<input type="checkbox"/>	2206216	12"-16"	
		<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Back Upholstery.	E0973 \$ 255
<input type="checkbox"/>	2206310	Set at 10"	<input type="checkbox"/>
<input type="checkbox"/>	2206311	Set at 11"	<input type="checkbox"/>
<input type="checkbox"/>	2206312	Set at 12"	
<input type="checkbox"/>	2206313	Set at 13"	<input type="checkbox"/>
<input type="checkbox"/>	2206314	Set at 14"	

Armrests			
		<b>Regular Armpads</b>	
<input type="checkbox"/>	2207100	Desk Length - 2" x 10"	<input type="checkbox"/>
<input type="checkbox"/>	2207101	Full Length - 2" x 14"	
<input type="checkbox"/>	2207102	Full Length Premium Grade - 2" x 14"	\$ 35
		<b>Wide Armpads</b> May interfere with wheels in some tilt configurations.	\$ 115
<input type="checkbox"/>	2207110	Desk Length - 3" x 10"	<input type="checkbox"/>
<input type="checkbox"/>	2207111	Full Length - 3" x 12"	
		<b>Gel Armpads</b> May interfere with wheels in some tilt configurations.	\$ 125
<input type="checkbox"/>	2207120	Full Length - 2" x 12"	<input type="checkbox"/>
<input type="checkbox"/>	2207121	Full Length - 3.5" x 12"	
		<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations.	E2209 \$ 125
<input type="checkbox"/>	2207130	Right (4.5" x 18")	<input type="checkbox"/>
<input type="checkbox"/>	2207131	Left (4.5" x 18")	
		<b>Armrest Extensions</b> Not compatible with Flip-up or Single Post Armrests.	\$ 60
<input type="checkbox"/>	2207400	3" Extension - Forward	<input type="checkbox"/>
<input type="checkbox"/>	2207401	3" Extension - Rearward	

Seat Pan			
<input type="checkbox"/>	2208000	<b>No Seat Pan</b>	
<input type="checkbox"/>	2208001	<b>Solid Aluminum Seat Pan</b>	E2231 \$ 120

Cushion			
<input type="checkbox"/>	2209000	<b>No Cushion</b>	
		<b>nxt BioFit</b>	
<input type="checkbox"/>	2209101	14"-20" width, 16"-20" depth	 E2607 \$ 595
<input type="checkbox"/>	2209102	22" width, 16"-20" depth	E2608 \$ 625
		<b>Matrx VI</b>	
<input type="checkbox"/>	2209111	14"-20" width, 16"-22" depth	 E2607 \$ 550
<input type="checkbox"/>	2209112	22" width, 18"-22" depth	E2608 \$ 600
<input type="checkbox"/>	2209115	HD - 18"-22" width, 18"-22" depth	E2608 \$ 1,150
		<b>ROHO Quadtro Select High Profile</b>	
<input type="checkbox"/>	2209121	14"-20" width, 16"-22" depth	 E2624 \$ 650
<input type="checkbox"/>	2209122	22" width, 16"-22" depth	E2625 \$ 750
		<b>ROHO High Profile</b>	
<input type="checkbox"/>	2209125	14"-20" width, 16"-22" depth	E2622 \$ 650
<input type="checkbox"/>	2209126	22" width, 16"-22" depth	E2623 \$ 750
		<b>Varilite Reflex</b>	
<input type="checkbox"/>	2209131	14"-20" width, 16"-20" depth	 E2603 \$ 375
		<b>Cushion Size</b>	
<input type="checkbox"/>	2209400	Specify Size (Width" x Depth")	_____

### Seat Height

#### Select desired Seat-to-Floor Height

STF Height is adjustable within the range of the selected Rear Wheel size.

RW	6" Casters					
12"	13"	14"	15"	16"		
16"			15"	16"	17"	18"
20"	13"	14"	15"	16"		
22"		14"	15"	16"	17"	
24"			15"	16"	17"	18"

RW	8" Casters					
16"			15"	16"	17"	18"
24"			15"	16"	17"	18"

Wheel size \_\_\_\_\_

Caster size \_\_\_\_\_

STF height \_\_\_\_\_

### Casters

#### Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211061 6" x 1.25" Caster size
- 2211063 6" x 1.50" Caster size \$ 85
- 2211081 8" x 1.25" Caster size
- 2211084 8" x 2.00" Caster size \$ 135

#### Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211161 6" x 1.25" Caster size *E2214 & E2215* \$ 110
- 2211181 8" x 1.25" Caster size *E2214 & E2215* \$ 120
- 2211184 8" x 2.00" Caster size *E2214 & E2215* \$ 150
- 2211260 **6" Soft Roll** \$ 95
- 2211400 **Ultra Sport Fork** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. *E1015* \$ 395

### Rear Wheels

#### Composite

- 2212012 12" Wheel size
- 2212016 16" Wheel size
- 2212020 20" Wheel size
- 2212022 22" Wheel size
- 2212024 24" Wheel size

#### Spoke

- 2212120 20" Wheel size \$ 75
- 2212122 22" Wheel size
- 2212124 24" Wheel size
- 2212224 **24" Heavy Duty Spoke** \$ 350

### Axle Type

- 2213000 **Quick Release Rear Axles**
- 2213010 **Fixed Rear Axles**

### Rear Tires

- 2214000 **Urethane**
  - 2214010 **Pneumatic** Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. *E2211 & E2212* \$ 60
  - 2214020 **Pneumatic with Flat-Free Inserts** *E2211 & E2213* \$ 95
- Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.

### Hand Rims

- 2215000 **No Hand Rims**
- 2215010 **Standard** 20", 22", and 24" Rear Wheels only.
- 2215020 **Plastic Coated** 20", 22", and 24" Rear Wheels only. \$ 135

### Wheel Locks

- 2216010 **Push to Lock**
- 2216020 **Pull to Lock**

See next page for more Wheel Locks options >

<input type="checkbox"/>	2216030	<b>Attendant Wheel Lock (Push to Lock)</b>		\$	85
<input type="checkbox"/>	2216040	<b>User &amp; Attendant Wheel Lock</b>	Available only on 20", 22", or 24" Rear Wheels.	\$	125
<input type="checkbox"/>	2216100	<b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.	E0961	\$ 60
<input type="checkbox"/>	2216200	<b>Hill Holders</b>	Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974	\$ 350

### Anti-Tippers

<input checked="" type="checkbox"/>	2217010	<b>Rear Anti-Tippers</b>	Required Safety Option.	E0971	\$ 60
<input type="checkbox"/>	2217020	<b>Front Anti-Tippers</b>		E0971	\$ 300
<input type="checkbox"/>	2217030	<b>Heavy Duty Rear Anti-Tipper</b>		E0971	\$ 600
Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.					

### Back Post Style

<input type="checkbox"/>	2218010	<b>Push Bar Back</b>			
<input type="checkbox"/>	2218020	<b>Cane Handle Back</b>	Not compatible with Headrest Mounting Plate.		
<b>Back Post Height</b> Seat Pan to Handle					
<input type="checkbox"/>	2218420	20" Height	<input type="checkbox"/>	2218425	25" Height
			<input type="checkbox"/>	2218431	31" Height

### Back Options

<b>Angle/Height Adjustable Push Bar Handle</b>				Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150
<input type="checkbox"/>	2219010	Button on Outside	<input type="checkbox"/>	2219011	Button on Inside
<b>Seat to Back Angle</b>					
<input type="checkbox"/>	2219090	Set at 90°	<input type="checkbox"/>	2219098	Set at 98°
			<input type="checkbox"/>	2219106	Set at 106°
<input type="checkbox"/>	2219114	Set at 114°	<input type="checkbox"/>	2219122	Set at 122°
<input type="checkbox"/>	2219210	<b>10" Backrest Extension</b>	Not compatible with Angle/Height Adjustable Push Bar Handle.	E0966	\$ 135
<input checked="" type="checkbox"/>	2219220	<b>Fold Down Back</b>			

### Back Upholstery

<input type="checkbox"/>	2220000	<b>No Back Upholstery</b>			
<b>Standard</b>				Back height measured from Seat Pan to top of Upholstery.	E2611 \$ 155
<input type="checkbox"/>	2220018	Set at 18"	<input type="checkbox"/>	2220019	Set at 19"
			<input type="checkbox"/>	2220020	Set at 20"
<input type="checkbox"/>	2220021	Set at 21"	<input type="checkbox"/>	2220022	Set at 22"
			<input type="checkbox"/>	2220023	Set at 23"
<b>Adjustable Tension</b>				Back height measured from Seat Pan to top of Upholstery.	E2611 \$ 295
<input type="checkbox"/>	2220321	Set at 21"	<input type="checkbox"/>	2220322	Set at 22"
			<input type="checkbox"/>	2220323	Set at 23"

### Backrest

<input type="checkbox"/>	2221000	<b>No Backrest</b>			
<b>next Xtend</b>					
<input type="checkbox"/>	2221100	16"-20" width, 16"-20" height		E2620	\$ 690
<input type="checkbox"/>	2221105	22" width, 16"-20" height		E2621	\$ 720
<b>next Xtend Deep</b>					
<input type="checkbox"/>	2221101	16"-20" width, 16"-20" height		E2620	\$ 690
<input type="checkbox"/>	2221106	22" width, 16"-20" height		E2621	\$ 720
<b>Matrix Elite</b>					
<input type="checkbox"/>	2221110	14"-20" width, 10"-20" height		E2620	\$ 825
<input type="checkbox"/>	2221112	HD - 22" width, 16"-20" height		E2621	\$ 1,095
<input type="checkbox"/>	2221116	Deep - 14"-20" width, 12"-20" height		E2620	\$ 950



See next page for more Backrest options >

<b>ROHO Agility</b>				
<input type="checkbox"/>	2221120 Min Contour 14"-20" width, 10"-18" height		E2620	\$ 825
<input type="checkbox"/>	2221122 Mid Contour 14"-20" width, 10"-18" height		E2620	\$ 850
<input type="checkbox"/>	2221124 Max Contour 14"-20" width, 10"-18" height		E2620	\$ 950
<b>Varilite Icon</b>				
<input type="checkbox"/>	2221130 Tall - 14"-20" width, 14"-24" height		E2615	\$ 850
<input type="checkbox"/>	2221132 Deep - 14"-20" width, 17"-24" height		E2620	\$ 895
<b>Backrest Size</b>				
<input type="checkbox"/>	2221400 Specify Size (Width" x Height") _____			

<b>Tilt/Recline Options</b>				
<b>Attendant Tilt only</b>				
<input type="checkbox"/>	2223011 Trigger on Right Back Post	<input type="checkbox"/>	2223012 Trigger on Left Back Post	
<b>User Tilt only</b> Full Length Arm Pads only.				\$ 100
<input type="checkbox"/>	2223021 Trigger on Right Armrest	<input type="checkbox"/>	2223022 Trigger on Left Armrest	
<b>Attendant &amp; User Tilt</b> Full Length Arm Pads only.				\$ 195
<input type="checkbox"/>	2223031 Trigger on Right Back Post	<input type="checkbox"/>	2223032 Trigger on Left Back Post	
<input type="checkbox"/>	2223041 Trigger on Right Armrest	<input type="checkbox"/>	2223042 Trigger on Left Armrest	
<b>Anterior Tilt Limit</b>				
<input type="checkbox"/>	2223220 Set at 20°	<input type="checkbox"/>	2223225 Set at 25°	<input type="checkbox"/>
<input type="checkbox"/>	2223230 Set at 30°			
<input type="checkbox"/>	2223400 <b>Power Tilt</b> 250 lbs Weight Capacity. 30° Anterior to 20° Posterior tilt range. Includes 2 batteries for off-chair charging. Minimum Seat Width is 16".			E1002 \$ 2,195
	Toggle Switch	<input type="checkbox"/>	2223421 Installed on Right side	<input type="checkbox"/>
	Button Switch	<input type="checkbox"/>	2223431 Installed on Right side	<input type="checkbox"/>
	Egg Switch	<input type="checkbox"/>	2223440 Not installed	
	Optional for Egg Switch:	<input type="checkbox"/>	2223445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models. E1028 \$ 625
<i>Specify Headrest to be used: _____</i>				

<b>Headrest</b>				
<input type="checkbox"/>	2224000 <b>No Headrest</b>			
<input type="checkbox"/>	2224010 <b>Headrest Mounting Plate</b>			\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.				
<input type="checkbox"/>	2224100 <b>next MultiFit Head Support Hardware</b>			E1028 \$ 235
<b>next MultiFit Head Support Pad</b>				
<input type="checkbox"/>	2224105 Regular (9" x 5")	<input type="checkbox"/>	2224106 Large (10" x 5")	<input type="checkbox"/>
	<input type="checkbox"/> 2224107 X-Large (14" x 5")			
<b>Stealth Headrest Hardware</b>				
<input type="checkbox"/>	2224110 Multi-Axis HMO475			E1028 \$ 150
<input type="checkbox"/>	2224111 Link TWB480			E1028 \$ 325
<b>Stealth Headrest Pad</b>				
<input type="checkbox"/>	2224115 10" width			E0955 \$ 245
<input type="checkbox"/>	2224116 14" width			

<b>Front Rigging</b>				
<input type="checkbox"/>	1025000 <b>No Front Rigging</b>			
<b>PDG Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame.				
<input type="checkbox"/>	2225410 Fixed 70°			\$ 95
<input type="checkbox"/>	2225420 Swing In/Away Universal Angle Adjustable (0°-110°) Set at (use increments of 10°): _____			\$ 195
<input type="checkbox"/>	2225430 Swing In/Away Elevating Legrests High knee pivot. Includes reinforced premium Calf Pad mount.			E0990 \$ 495

Foot Plates			
<input type="checkbox"/>	2226010	<b>Composite</b>	
<input type="checkbox"/>	2226020	<b>Aluminum Angle Adjustable</b>	K0040 \$ 155
<input type="checkbox"/>	2226030	<b>One Piece Fixed</b> Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	2226040	<b>One Piece Flip-up</b> Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	2226050	<b>Interlocking Angle Adjustable</b> Seat Widths of 14"-21" only.	K0040 \$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	2227040	<b>Adjustable Length Heel Loops</b>	E0951 \$ 50
<input type="checkbox"/>	2227050	<b>Bodypoint - Padded Calf Strap</b>	K0038 \$ 105

Extension Tubes			
<b>Regular Extension Tubes</b>			
<input type="checkbox"/>	2228010	Short 14"-18" for Easy-Out SA Angle Adjustable 11"-15" for Easy Out SA ELR	
<input type="checkbox"/>	2228020	Long 18"-24" for Easy-Out SA Angle Adjustable 15"-21" for Easy-Out SA ELR	\$ 75

Positioning Straps			
<b>Bodypoint - Hip Belt</b>			
<input type="checkbox"/>	2229010	Center-Pull Auto Style	E0978 \$ 105
<input type="checkbox"/>	2229015	Center-Pull Auto Style, Padded	E0978 \$ 155
<input type="checkbox"/>	2229020	4-Point, Padded	E0978 \$ 175
<input type="checkbox"/>	2229030	<b>Bodypoint - Chest Strap</b>	E0960 \$ 85



Miscellaneous Options			
<input type="checkbox"/>	2230110	<b>IV Pole</b> Not compatible with HD Anti Tipper, Attendant Wheel Locks.	E1029 \$ 295
<input type="checkbox"/>	2230120	<b>O2 Holder</b> Not compatible with HD Anti Tipper, Attendant Wheel Locks.	E2208 \$ 395
<input type="checkbox"/>	2230170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.



Scan this code or go to [pdgmobility.com/products/stellar-leap/](http://pdgmobility.com/products/stellar-leap/) to learn more about the **Stellar LEAP Manual Tilt Wheelchair**

