

Customer and Order Information

Purchase Order *	Tag	v.072318
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Billing Address *		
Shipping Address *		*required

Frame Type

<input checked="" type="checkbox"/> 1201000	Stellar 45° Manual Tilt Wheelchair	E1161	\$ 2,995
---	---	-------	----------

Frame Options

<input type="checkbox"/> 1202000	Standard Base and Seat Frame	250 lbs Weight Capacity.		
<input type="checkbox"/> 1202050	Upgraded Base and Seat Frame	300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	K0108	\$ 200
<input type="checkbox"/> 1202100	Heavy Duty Base and Seat Frame	450 lbs Weight Capacity. Required for Seat Width of 23" or 24" Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	K0108	\$ 300
<input type="checkbox"/> 1202200	Extended Base with Rear Floor Jacks	Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Locks, HD Anti-Tipper, IV Pole, O2 Holder.	K0108	\$ 600
<input type="checkbox"/> 1202201	Shortened Strut Bracket	Increases heel clearance for foot propelling. Available only for Seat Width of 23" & 24".		\$ 125
<input type="checkbox"/> 1202300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.		

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths					
<input type="checkbox"/> 1203014	14" Seat Width	<input type="checkbox"/> 1203015	15" Seat Width	<input type="checkbox"/> 1203016	16" Seat Width
<input type="checkbox"/> 1203017	17" Seat Width	<input type="checkbox"/> 1203018	18" Seat Width	<input type="checkbox"/> 1203019	19" Seat Width
<input type="checkbox"/> 1203020	20" Seat Width	<input type="checkbox"/> 1203021	21" Seat Width	<input type="checkbox"/> 1203022	22" Seat Width
				E2201	\$ 325
Custom Seat Widths Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.					
<input type="checkbox"/> 1203023	23" Seat Width			E2201	\$ 325
<input type="checkbox"/> 1203024	24" Seat Width			E2202	\$ 525
Use Stellar HD Order Form for Seat Widths greater than 24"					

Seat Depth

Standard Seat Depths (Adjustable)					
<input type="checkbox"/> 1204016	16" Seat Depth	<input type="checkbox"/> 1204017	17" Seat Depth	<input type="checkbox"/> 1204018	18" Seat Depth
<input type="checkbox"/> 1204019	19" Seat Depth				
<input type="checkbox"/> 1204020	20" Seat Depth			E2203	\$ 475
Custom Seat Depths					
<input type="checkbox"/> 1204021	21" Seat Depth			E2203	\$ 475
<input type="checkbox"/> 1204022	22" Seat Depth			E2204	\$ 750

Frame Color

PDG Premium Colors					
<input type="checkbox"/> 1205000	Tracer Red	<input type="checkbox"/> 1205001	Gloss Black	<input type="checkbox"/> 1205002	Silver
<input type="checkbox"/> 1205003	Black Shimmer	<input type="checkbox"/> 1205004	Orange Fusion	<input type="checkbox"/> 1205005	Karl's Imagination
<input type="checkbox"/> 1205006	Blue Shimmer	<input type="checkbox"/> 1205007	Gloss White	<input type="checkbox"/> 1205008	Dark Titanium

See next page for more Frame Color options >

<input type="checkbox"/> 1205009 Geyser Green	<input type="checkbox"/> 1205010 Lollypop Blue	<input type="checkbox"/> 1205011 Lollypop Red
<input type="checkbox"/> 1205012 Bubble Gum Pink	<input type="checkbox"/> 1205013 Berry Purple	<input type="checkbox"/> 1205014 Safety Yellow
<input type="checkbox"/> 1205200 Custom Color Contact PDG for specifications.		\$ 300

Arm Type

<input type="checkbox"/> 1206000 No Armrest			
<input type="checkbox"/> 1206100 Single Post Fixed Height (12")	Full Length Armpad only.		
Single Post Adjustable Height		Not compatible with the Heavy Duty Base and/or Seat Width greater than 22".	E0973 \$ 255
<input type="checkbox"/> 1206111 8"-11"	<input type="checkbox"/> 1206114 10"-14"		
<input type="checkbox"/> 1206200 Dual Post Fixed Height (12")	Full Length Armpad only.		
Dual Post Adjustable Height			E0973 \$ 255
<input type="checkbox"/> 1206211 8"-11"	<input type="checkbox"/> 1206214 10"-14"	<input type="checkbox"/> 1206216 12"-16"	
Locking Flip-up		Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.	E0973 \$ 255
<input type="checkbox"/> 1206310 Set at 10"	<input type="checkbox"/> 1206311 Set at 11"	<input type="checkbox"/> 1206312 Set at 12"	
<input type="checkbox"/> 1206313 Set at 13"	<input type="checkbox"/> 1206314 Set at 14"		

Armrests

Regular Armpads			
<input type="checkbox"/> 1207100 Desk Length - 2" x 10"	<input type="checkbox"/> 1207101 Full Length - 2" x 14"		
<input type="checkbox"/> 1207102 Full Length Premium Grade - 2" x 14"			\$ 35
Wide Armpads		May interfere with wheels in some tilt configurations.	\$ 115
<input type="checkbox"/> 1207110 Desk Length - 3" x 10"	<input type="checkbox"/> 1207111 Full Length - 3" x 12"		
Gel Armpads		May interfere with wheels in some tilt configurations.	\$ 125
<input type="checkbox"/> 1207120 Full Length - 2" x 12"	<input type="checkbox"/> 1207121 Full Length - 3.5" x 12"		
Contour Arm Trough with Hand Pad		May interfere with wheels in some tilt configurations.	E2209 \$ 125
<input type="checkbox"/> 1207130 Right (4.5" x 18")	<input type="checkbox"/> 1207131 Left (4.5" x 18")		
<input type="checkbox"/> 1207220 Armrest Wideners for Dual Post Armrests	Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.		\$ 375
Clothing Guards		Not compatible with Flip-up Armrests, Armrest Wideners.	
<input type="checkbox"/> 1207310 For Single Post Armrests	Available for 10"-14" Height range only.		K0108 \$ 110
<input type="checkbox"/> 1207320 For Dual Post Armrests			K0108 \$ 165
Armrest Extensions		Not compatible with Flip-up or Single Post Armrests.	\$ 60
<input type="checkbox"/> 1207400 3" Extension - Forward	<input type="checkbox"/> 1207401 3" Extension - Rearward		

Seat Pan

<input type="checkbox"/> 1208000 No Seat Pan		
<input type="checkbox"/> 1208001 Solid Aluminum Seat Pan		E2231 \$ 120

Cushion

<input type="checkbox"/> 1209000 No Cushion		
next BioFit		
<input type="checkbox"/> 1209101 14"-20" width, 16"-20" depth		E2607 \$ 595
<input type="checkbox"/> 1209102 22"-24" width, 16"-20" depth		E2608 \$ 625
Matrx VI		
<input type="checkbox"/> 1209111 14"-20" width, 16"-22" depth		E2607 \$ 550
<input type="checkbox"/> 1209112 22" width, 18"-22" depth		E2608 \$ 600
<input type="checkbox"/> 1209115 HD - 18"-24" width, 18"-22" depth		E2608 \$ 1,150

See next page for more Cushion options >

ROHO Quadro Select High Profile				
<input type="checkbox"/>	1209121	14"-20" width, 16"-22" depth	E2624	\$ 650
<input type="checkbox"/>	1209122	22"-24" width, 16"-22" depth	E2625	\$ 750
ROHO High Profile				
<input type="checkbox"/>	1209125	14"-20" width, 16"-22" depth	E2622	\$ 650
<input type="checkbox"/>	1209126	22"-24" width, 16"-22" depth	E2623	\$ 750
Varilite Reflex				
<input type="checkbox"/>	1209131	14"-20" width, 16"-20" depth	E2603	\$ 375
Cushion Size				
1209400	Specify Size (Width" x Depth") _____			



Seat Height					
Select desired Seat-to-Floor Height			STF Height is adjustable within the range of the selected Rear Wheel size.		
RW	6" Casters				Wheel size _____ Caster size _____ STF height _____
12"	14"	15"	16"	17"	
16"			16"	17"	
20"	14"	15"	16"	17"	
22"		15"	16"	17"	
24"			16"	17"	
RW	8" Casters				
16"			16"	17"	
24"			16"	17"	

Casters				
Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.				
<input type="checkbox"/>	1211061	6" x 1.25" Caster size		
<input type="checkbox"/>	1211063	6" x 1.50" Caster size		\$ 85
<input type="checkbox"/>	1211081	8" x 1.25" Caster size		
<input type="checkbox"/>	1211084	8" x 2.00" Caster size		\$ 135
Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.				
<input type="checkbox"/>	1211161	6" x 1.25" Caster size	E2214 & E2215	\$ 110
<input type="checkbox"/>	1211181	8" x 1.25" Caster size	E2214 & E2215	\$ 120
<input type="checkbox"/>	1211184	8" x 2.00" Caster size	E2214 & E2215	\$ 150
<input type="checkbox"/>	1211260	6" Soft Roll		\$ 95
<input type="checkbox"/>	1211400	Ultra Sport Fork Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	E1015	\$ 395

Rear Wheels				
Composite				
<input type="checkbox"/>	1212012	12" Wheel size	<input type="checkbox"/>	1212016 16" Wheel size
<input type="checkbox"/>	1212022	22" Wheel size	<input type="checkbox"/>	1212024 24" Wheel size
<input type="checkbox"/>			<input type="checkbox"/>	1212020 20" Wheel size
Spoke \$ 75				
<input type="checkbox"/>	1212120	20" Wheel size	<input type="checkbox"/>	1212122 22" Wheel size
<input type="checkbox"/>	1212224	24" Heavy Duty Spoke Required for Heavy Duty Base and Seat Frame.		\$ 350

Axle Type				
<input type="checkbox"/>	1213000	Quick Release Rear Axles		
<input type="checkbox"/>	1213010	Fixed Rear Axles		

Rear Tires				
<input type="checkbox"/>	1214000	Urethane		
<input type="checkbox"/>	1214010	Pneumatic Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212	\$ 60
<input type="checkbox"/>	1214020	Pneumatic with Flat-Free Inserts Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.	E2211 & E2213	\$ 95

Hand Rims			
<input type="checkbox"/>	1215000	No Hand Rims	
<input type="checkbox"/>	1215010	Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1215020	Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 135


Wheel Locks			
<input type="checkbox"/>	1216010	Push to Lock	
<input type="checkbox"/>	1216020	Pull to Lock	
<input type="checkbox"/>	1216030	Attendant Wheel Lock (Push to Lock) Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.	\$ 85
<input type="checkbox"/>	1216040	User & Attendant Wheel Lock Available only on 20", 22", or 24" Rear Wheels.	\$ 125
<input type="checkbox"/>	1216100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 60
<input type="checkbox"/>	1216200	Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 350

Anti-Tippers			
<input checked="" type="checkbox"/>	1217010	Rear Anti-Tippers Required Safety Option.	E0971 \$ 60
<input type="checkbox"/>	1217020	Front Anti-Tippers	E0971 \$ 300
<input type="checkbox"/>	1217030	Heavy Duty Rear Anti-Tipper Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	E0971 \$ 600

Back Post Style			
<input type="checkbox"/>	1218010	Push Bar Back	
<input type="checkbox"/>	1218020	Cane Handle Back Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
		Back Post Height Seat Pan to Handle	
<input type="checkbox"/>	1218425	25" Height	<input type="checkbox"/> 1218431 31" Height

Back Options			
		Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150
<input type="checkbox"/>	1219010	Button on Outside	<input type="checkbox"/> 1219011 Button on Inside
		Seat to Back Angle	
<input type="checkbox"/>	1219090	Set at 90°	<input type="checkbox"/> 1219098 Set at 98° <input type="checkbox"/> 1219106 Set at 106°
<input type="checkbox"/>	1219114	Set at 114°	<input type="checkbox"/> 1219122 Set at 122°
<input type="checkbox"/>	1219210	10" Backrest Extension Not compatible with Angle/Height Adjustable Push Bar Handle.	E0966 \$ 135
<input type="checkbox"/>	1219220	Fold Down Back Not compatible with Reclining Back.	
<input type="checkbox"/>	1219230	Spring Back Suspension Not compatible with Reclining Back.	\$ 105
<input type="checkbox"/>	1219240	Dynamic Rocker Back Assembly Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.	K0108 \$ 2,995

Back Upholstery			
<input type="checkbox"/>	1220000	No Back Upholstery	
		Standard Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1220021	Set at 21"	<input type="checkbox"/> 1220022 Set at 22" Not available with 25" Back.
		Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1220118	Set at 18"	<input type="checkbox"/> 1220119 Set at 19" <input type="checkbox"/> 1220120 Set at 20"
		Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 155
<input type="checkbox"/>	1220224	Set at 24"	<input type="checkbox"/> 1220225 Set at 25" <input type="checkbox"/> 1220226 Set at 26"
		Adjustable Tension Back height measured from Seat Pan to top of Upholstery.	E2611 / E2612 \$ 295
<input type="checkbox"/>	1220321	Set at 21"	<input type="checkbox"/> 1220322 Set at 22" <input type="checkbox"/> 1220323 Set at 23"

Backrest			
<input type="checkbox"/>	1221000	No Backrest	
nxt Xtend			
<input type="checkbox"/>	1221100	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1221105	22"-24" width, 16"-20" height	E2621 \$ 720
nxt Xtend Deep 			
<input type="checkbox"/>	1221101	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1221106	22"-24" width, 16"-20" height	E2621 \$ 720
Matrx Elite			
<input type="checkbox"/>	1221110	14"-20" width, 10"-20" height	E2620 \$ 825
<input type="checkbox"/>	1221112	HD - 22"-24" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/>	1221116	Deep - 14"-20" width, 12"-20" height	E2620 \$ 950
ROHO Agility			
<input type="checkbox"/>	1221120	Min Contour 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1221122	Mid Contour 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1221124	Max Contour 14"-20" width, 10"-18" height	E2620 \$ 950
Varilite Icon			
<input type="checkbox"/>	1221130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1221132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
Backrest Size			
	1221400	Specify Size (Width" x Height")	_____

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1223011	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1223012	Trigger on Left Back Post	
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. \$ 100			
<input type="checkbox"/>	1223021	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1223022	Trigger on Left Armrest	
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. \$ 195			
<input type="checkbox"/>	1223031	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1223032	Trigger on Left Back Post	
<input type="checkbox"/>	1223041	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1223042	Trigger on Left Armrest	
<input type="checkbox"/>	1223110	Trigger Extension	\$ 50
Tilt Limit			
<input type="checkbox"/>	1223220	Set at 20°	<input type="checkbox"/>
<input type="checkbox"/>	1223225	Set at 25°	<input type="checkbox"/>
<input type="checkbox"/>	1223230	Set at 30°	
<input type="checkbox"/>	1223235	Set at 35°	<input type="checkbox"/>
<input type="checkbox"/>	1223240	Set at 40°	<input type="checkbox"/>
<input type="checkbox"/>	1223245	Set at 45°	
<input type="checkbox"/>	1223410	Reclining Backrest (Manual)	E1225 \$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.			

Headrest			
<input type="checkbox"/>	1224000	No Headrest	
<input type="checkbox"/>	1224010	Headrest Mounting Plate	\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	1224100	nxt MultiFit Head Support Hardware	E1028 \$ 235
nxt MultiFit Head Support Pad 			
<input type="checkbox"/>	1224105	Regular (9" x 5")	<input type="checkbox"/>
<input type="checkbox"/>	1224106	Large (10" x 5")	<input type="checkbox"/>
<input type="checkbox"/>	1224107	X-Large (14" x 5")	E0955 \$ 245
Stealth Headrest Hardware			
<input type="checkbox"/>	1224110	Multi-Axis HMO475	E1028 \$ 150
<input type="checkbox"/>	1224111	Link TWB480	E1028 \$ 325
Stealth Headrest Pad 			
<input type="checkbox"/>	1224115	10" width	<input type="checkbox"/>
<input type="checkbox"/>	1224116	14" width	

Front Rigging

<input type="checkbox"/>	1225000	No Front Rigging	<input type="checkbox"/>	1225001	Sector Blocks only		
PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame.							
<input type="checkbox"/>	1225410	Fixed 70°					\$ 95
<input type="checkbox"/>	1225420	Swing In/Away Universal Angle Adjustable (0°-110°)	Set at (use increments of 10°): _____				\$ 195
<input type="checkbox"/>	1225430	Swing In/Away Elevating Legrests	High knee pivot. Includes reinforced premium Calf Pad mount.			E0990	\$ 495
Standard Front Rigging Legacy Design. Not compatible with Easy-Out™.							
<input type="checkbox"/>	1225110	Plug In 70°					
<input type="checkbox"/>	1225130	Swing Away 60°	<input type="checkbox"/>	1225140	Swing Away 70°		
<input type="checkbox"/>	1225200	Swing Away Elevating Legrests				E0990	\$ 375
<input type="checkbox"/>	1225300	1" Offset Hangers for Swing Away					\$ 75

Foot Plates

<input type="checkbox"/>	1226010	Composite					
<input type="checkbox"/>	1226020	Aluminum Angle Adjustable				K0040	\$ 155
<input type="checkbox"/>	1226030	One Piece Fixed	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	1226040	One Piece Flip-up	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	1226050	Interlocking Angle Adjustable	Seat Widths of 14"-21" only.			K0040	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1227010	Swing Away 60° & 70° Calf Pad	Not compatible with 60° & 70° Gel Padding.				\$ 135
<input type="checkbox"/>	1227020	Swing Away 60° & 70° Gel Padding	Not compatible with 60° & 70° Swing Away Calf Pad, Easy-Out Front Rigging.				\$ 100
<input type="checkbox"/>	1227030	Swing Away ELR Gel Knee Pad					\$ 100
<input type="checkbox"/>	1227040	Adjustable Length Heel Loops				E0951	\$ 50
<input type="checkbox"/>	1227050	Bodypoint - Padded Calf Strap	Sizes greater than 22" are not padded.			K0038	\$ 105
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.							
<input type="checkbox"/>	1227060	For Plug In Front Rigging					\$ 190
<input type="checkbox"/>	1227065	For Swing Away 70° Front Rigging	Not available on 13" STF Height.				\$ 380

Extension Tubes

Regular Extension Tubes							
<input type="checkbox"/>	1228010	Short	14"-18" for Easy-Out SA Angle Adjustable	11"-15" for Easy Out SA ELR	13"-17" for Standard Front Rigging		
<input type="checkbox"/>	1228020	Long	18"-24" for Easy-Out SA Angle Adjustable	15"-21" for Easy-Out SA ELR	15"-19" for Standard Front Rigging		\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.							
<input type="checkbox"/>	1228110	Standard (13"-17")	<input type="checkbox"/>	1228120	Long (15"-19")		\$ 100

Positioning Straps

Bodypoint - Hip Belts							
<input type="checkbox"/>	1229010	Center-Pull Auto Style				E0978	\$ 105
<input type="checkbox"/>	1229015	Center-Pull Auto Style, Padded				E0978	\$ 155
<input type="checkbox"/>	1229020	4-Point, Padded				E0978	\$ 175
<input type="checkbox"/>	1229030	Bodypoint - Chest Strap				E0960	\$ 85



Miscellaneous Options			
<input type="checkbox"/>	1230050	Protective Lower Padding	\$ 425
<input type="checkbox"/>	1230080	Gel Hip Guides Reduces available Seat Width by 2".	E0956 & E1028 \$ 325
<input type="checkbox"/>	1230110	IV Pole Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	E1029 \$ 295
<input type="checkbox"/>	1230120	O2 Holder Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	E2208 \$ 395
<input type="checkbox"/>	1230170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar/ to learn more about the **Stellar 45° Manual Tilt Wheelchair**

