

Customer and Order Information	
Purchase Order *	Tag
Dealer/Provider Name *	
Purchaser Name and Contact Information *	
Billing Address *	
Shipping Address *	*required

Frame Type	
<input checked="" type="checkbox"/> 1201000 Stellar 45° Manual Tilt Wheelchair	E1161 \$ 2,995

Frame Options	
<input type="checkbox"/> 1202000 Standard Base and Seat Frame 250 lbs Weight Capacity.	
<input type="checkbox"/> 1202050 Upgraded Base and Seat Frame 300 lbs Weight Capacity with addition of a secondary Gas Strut. Available for Seat Widths up to 22".	K0108 \$ 200
<input type="checkbox"/> 1202100 Heavy Duty Base and Seat Frame 450 lbs Weight Capacity. Required for Seat Width of 23" or 24" Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	K0108 \$ 300
<input type="checkbox"/> 1202200 Extended Base with Rear Floor Jacks Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant wheel lock, HD Anti-Tipper, IV Pole, O2 Holder.	K0108 \$ 600
<input type="checkbox"/> 1202201 Shortened Strut Bracket Increases heel clearance for foot propelling. Available only for Seat Width of 23" & 24".	\$ 125
<input type="checkbox"/> 1202300 Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	

Seat Width		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW
Standard Seat Widths		
<input type="checkbox"/> 1203014 14" Seat Width	<input type="checkbox"/> 1203015 15" Seat Width	<input type="checkbox"/> 1203016 16" Seat Width
<input type="checkbox"/> 1203017 17" Seat Width	<input type="checkbox"/> 1203018 18" Seat Width	<input type="checkbox"/> 1203019 19" Seat Width
<input type="checkbox"/> 1203020 20" Seat Width	<input type="checkbox"/> 1203021 21" Seat Width	<input type="checkbox"/> 1203022 22" Seat Width
		E2201 \$ 325
Custom Seat Widths Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1203023 23" Seat Width		E2201 \$ 325
<input type="checkbox"/> 1203024 24" Seat Width		E2202 \$ 525
Use Stellar HD Order Form for Seat Widths greater than 24"		

Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 1204016 16" Seat Depth	<input type="checkbox"/> 1204017 17" Seat Depth
<input type="checkbox"/> 1204019 19" Seat Depth	<input type="checkbox"/> 1204018 18" Seat Depth
<input type="checkbox"/> 1204020 20" Seat Depth	
	E2203 \$ 475
Custom Seat Depths	
<input type="checkbox"/> 1204021 21" Seat Depth	E2203 \$ 475
<input type="checkbox"/> 1204022 22" Seat Depth	E2204 \$ 750

Frame Color	
PDG Premium Colors	
<input type="checkbox"/> 1205000 Tracer Red	<input type="checkbox"/> 1205001 Gloss Black
<input type="checkbox"/> 1205003 Black Shimmer	<input type="checkbox"/> 1205002 Silver
<input type="checkbox"/> 1205006 Blue Shimmer	<input type="checkbox"/> 1205004 Orange Fusion
	<input type="checkbox"/> 1205005 Karl's Imagination
	<input type="checkbox"/> 1205007 Gloss White
	<input type="checkbox"/> 1205008 Dark Titanium

See next page for more Frame Color options >

<input type="checkbox"/> 1205009 Geyser Green	<input type="checkbox"/> 1205010 Lollypop Blue	<input type="checkbox"/> 1205011 Lollypop Red
<input type="checkbox"/> 1205012 Bubble Gum Pink	<input type="checkbox"/> 1205013 Berry Purple	<input type="checkbox"/> 1205014 Safety Yellow
<input type="checkbox"/> 1205200 Custom Color Contact PDG for specifications.		\$ 300

Arm Type		
<input type="checkbox"/> 1206000 No Armrest		
<input type="checkbox"/> 1206100 Single Post Fixed Height (12")	Full Length Armpad only.	
Single Post Adjustable Height Not compatible with the Heavy Duty Base and/or Seat Width greater than 22".		<i>E0973</i> \$ 255
<input type="checkbox"/> 1206111 8"-11"	<input type="checkbox"/> 1206114 10"-14"	
<input type="checkbox"/> 1206200 Dual Post Fixed Height (12")	Full Length Armpad only.	
Dual Post Adjustable Height		<i>E0973</i> \$ 255
<input type="checkbox"/> 1206211 8"-11"	<input type="checkbox"/> 1206214 10"-14"	<input type="checkbox"/> 1206216 12"-16"
Locking Flip-up Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.		<i>E0973</i> \$ 255
<input type="checkbox"/> 1206310 Set at 10"	<input type="checkbox"/> 1206311 Set at 11"	<input type="checkbox"/> 1206312 Set at 12"
<input type="checkbox"/> 1206313 Set at 13"	<input type="checkbox"/> 1206314 Set at 14"	

Armrests		
Regular Armpads		
<input type="checkbox"/> 1207100 Desk Length - 2" x 10"	<input type="checkbox"/> 1207101 Full Length - 2" x 14"	
<input type="checkbox"/> 1207102 Full Length Premium Grade - 2" x 14"		\$ 35
Wide Armpads May interfere with wheels in some tilt configurations.		\$ 115
<input type="checkbox"/> 1207110 Desk Length - 3" x 10"	<input type="checkbox"/> 1207111 Full Length - 3" x 12"	
Gel Armpads May interfere with wheels in some tilt configurations.		\$ 125
<input type="checkbox"/> 1207120 Full Length - 2" x 12"	<input type="checkbox"/> 1207121 Full Length - 3.5" x 12"	
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.		<i>E2209</i> \$ 125
<input type="checkbox"/> 1207130 Right (4.5" x 18")	<input type="checkbox"/> 1207131 Left (4.5" x 18")	
<input type="checkbox"/> 1207220 Armrest Wideners for Dual Post Armrests Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.		\$ 375
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.		
<input type="checkbox"/> 1207310 For Single Post Armrests Available for 10"-14" Height range only.		<i>K0108</i> \$ 110
<input type="checkbox"/> 1207320 For Dual Post Armrests		<i>K0108</i> \$ 165
Armrest Extensions Not compatible with Flip-up or Single Post Armrests.		\$ 60
<input type="checkbox"/> 1207400 3" Extension - Forward	<input type="checkbox"/> 1207401 3" Extension - Rearward	

Seat Pan		
<input type="checkbox"/> 1208000 No Seat Pan		
<input type="checkbox"/> 1208001 Solid Seat Pan		<i>E2231</i> \$ 120

Cushion		
<input type="checkbox"/> 1209000 No Cushion		
nxt BioFit		
<input type="checkbox"/> 1209101 14"-20" width, 16"-20" depth		<i>E2607</i> \$ 595
<input type="checkbox"/> 1209102 22"-24" width, 16"-20" depth		<i>E2608</i> \$ 625
Matrx VI		
<input type="checkbox"/> 1209111 14"-20" width, 16"-22" depth		<i>E2607</i> \$ 550
<input type="checkbox"/> 1209112 22" width, 18"-22" depth		<i>E2608</i> \$ 600
<input type="checkbox"/> 1209115 HD - 18"-24" width, 18"-22" depth		<i>E2608</i> \$ 1,150

See next page for more Cushion options >

ROHO Quadro Select High Profile				
<input type="checkbox"/>	1209121	14"-20" width, 16"-22" depth	<i>E2624</i>	\$ 650
<input type="checkbox"/>	1209122	22"-24" width, 16"-22" depth	<i>E2625</i>	\$ 750
ROHO High Profile				
<input type="checkbox"/>	1209125	14"-20" width, 16"-22" depth	<i>E2622</i>	\$ 650
<input type="checkbox"/>	1209126	22"-24" width, 16"-22" depth	<i>E2623</i>	\$ 750
Varilite Reflex				
<input type="checkbox"/>	1209131	14"-20" width, 16"-20" depth	<i>E2603</i>	\$ 375
Cushion Size				
1209400	Specify Size (Width" x Depth") _____			



Seat Height

Select desired Seat-to-Floor Height STF Height is adjustable within the range of the selected Rear Wheel size.

RW	6" Casters						
12"	14"	15"	16"	17"	18"	19"	20"
16"			16"	17"	18"	19"	20"
20"	14"	15"	16"	17"	18"		
22"		15"	16"	17"	18"	19"	
24"			16"	17"	18"	19"	20"

RW	8" Casters					
16"		16"	17"	18"	19"	20"
24"		16"	17"	18"	19"	20"

Wheel size _____

Caster size _____

STF height _____

Casters

Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1211061	6" x 1.25" Caster size		
<input type="checkbox"/>	1211063	6" x 1.50" Caster size		\$ 85
<input type="checkbox"/>	1211081	8" x 1.25" Caster size		
<input type="checkbox"/>	1211084	8" x 2.00" Caster size		\$ 135

Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	1211161	6" x 1.25" Caster size	<i>E2214 & E2215</i>	\$ 110
<input type="checkbox"/>	1211181	8" x 1.25" Caster size	<i>E2214 & E2215</i>	\$ 120
<input type="checkbox"/>	1211184	8" x 2.00" Caster size	<i>E2214 & E2215</i>	\$ 150

1211260 **6" Soft Roll** \$ 95

1211400 **Ultra Sport Fork** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. *E1015* \$ 395

Rear Wheels

Composite

<input type="checkbox"/>	1212012	12" Wheel size	<input type="checkbox"/>	1212016	16" Wheel size	<input type="checkbox"/>	1212020	20" Wheel size
<input type="checkbox"/>	1212022	22" Wheel size	<input type="checkbox"/>	1212024	24" Wheel size			

Spoke \$ 75

<input type="checkbox"/>	1212120	20" Wheel size	<input type="checkbox"/>	1212122	22" Wheel size	<input type="checkbox"/>	1212124	24" Wheel size
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------

1212224 **24" Heavy Duty Spoke** Required for Heavy Duty Base and Seat Frame \$ 350

Axle Type

1213000 **Quick Release Rear Axles**

1213010 **Fixed Rear Axles**

Rear Tires

1214000 **Urethane**

1214010 **Pneumatic** Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. *E2211 & E2212* \$ 60

1214020 **Pneumatic with Flat-Free Inserts** *E2211 & E2213* \$ 95

Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoke Wheels.

Hand Rims			
<input type="checkbox"/>	1215000	No Hand Rims	
<input type="checkbox"/>	1215010	Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1215020	Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 95

Wheel Locks			
<input type="checkbox"/>	1216010	Push to Lock	
<input type="checkbox"/>	1216020	Pull to Lock	
<input type="checkbox"/>	1216030	Attendant Wheel Lock (Push to Lock) Not compatible with TTDS or Width greater than 23" on chairs with Upgraded or Heavy Duty Base and Seat Frame.	\$ 85
<input type="checkbox"/>	1216040	User & Attendant Wheel Lock Available only on 20", 22", or 24" Rear Wheels.	\$ 125
<input type="checkbox"/>	1216100	Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 60
<input type="checkbox"/>	1216200	Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 350

Anti-Tippers			
<input checked="" type="checkbox"/>	1217010	Rear Anti-Tippers Required Safety Option.	E0971 \$ 60
<input type="checkbox"/>	1217020	Front Anti-Tippers	E0971 \$ 300
<input type="checkbox"/>	1217030	Heavy Duty Rear Anti-Tipper Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	E0971 \$ 600

Back Post Style			
<input type="checkbox"/>	1218010	Push Bar Back	
<input type="checkbox"/>	1218020	Cane Handle Back Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
Back Post Height Seat Pan to Handle			
<input type="checkbox"/>	1218425	25" Height	<input type="checkbox"/>
			<input type="checkbox"/> 1218431 31" Height

Back Options			
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.			\$ 150
<input type="checkbox"/>	1219010	Button on Outside	<input type="checkbox"/>
			<input type="checkbox"/> 1219011 Button on Inside
Seat to Back Angle			
<input type="checkbox"/>	1219090	Set at 90°	<input type="checkbox"/>
			<input type="checkbox"/> 1219098 Set at 98°
			<input type="checkbox"/> 1219106 Set at 106°
<input type="checkbox"/>	1219114	Set at 114°	<input type="checkbox"/>
			<input type="checkbox"/> 1219122 Set at 122°
<input type="checkbox"/>	1219210	10" Backrest Extension Not compatible with Angle/Height Adjustable Push Bar Handle.	E0966 \$ 135
<input type="checkbox"/>	1219220	Fold Down Back Not compatible with Reclining Back.	
<input type="checkbox"/>	1219230	Spring Back Suspension Not compatible with Reclining Back.	\$ 105
<input type="checkbox"/>	1219240	Dynamic Rocker Back Assembly Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller than selected Back Post Height.	K0108 \$ 2,995

Back Upholstery			
<input type="checkbox"/>	1220000	No Back Upholstery	
Standard Back height measured from Seat Pan to top of Upholstery.			E2611 / E2612 \$ 155
<input type="checkbox"/>	1220021	Set at 21"	<input type="checkbox"/>
			<input type="checkbox"/> 1220022 Set at 22"
Short 25" Back only. Back height measured from Seat Pan to top of Upholstery.			E2611 / E2612 \$ 155
<input type="checkbox"/>	1220118	Set at 18"	<input type="checkbox"/>
			<input type="checkbox"/> 1220119 Set at 19"
			<input type="checkbox"/> 1220120 Set at 20"
Tall 31" Back only. Back height measured from Seat Pan to top of Upholstery.			E2611 / E2612 \$ 155
<input type="checkbox"/>	1220224	Set at 24"	<input type="checkbox"/>
			<input type="checkbox"/> 1220225 Set at 25"
			<input type="checkbox"/> 1220226 Set at 26"
Adjustable Tension Back height measured from Seat Pan to top of Upholstery.			E2611 / E2612 \$ 295
<input type="checkbox"/>	1220321	Set at 21"	<input type="checkbox"/>
			<input type="checkbox"/> 1220322 Set at 22"
			<input type="checkbox"/> 1220323 Set at 23"

Backrest			
<input type="checkbox"/>	1221000	No Backrest	
nxt Optima HA Active Contour			
<input type="checkbox"/>	1221100	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1221105	22"-24" width, 16"-20" height	E2621 \$ 720
nxt Optima HA Plus Contour			
<input type="checkbox"/>	1221101	16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1221106	22"-24" width, 16"-20" height	E2621 \$ 720
Matrx Elite			
<input type="checkbox"/>	1221110	14"-20" width, 10"-20" height	E2620 \$ 825
<input type="checkbox"/>	1221112	HD - 22"-24" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/>	1221116	Deep - 14"-20" width, 12"-20" height	E2620 \$ 950
ROHO Agility			
<input type="checkbox"/>	1221120	Min Contour 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1221122	Mid Contour 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1221124	Max Contour 14"-20" width, 10"-18" height	E2620 \$ 950
Varilite Icon			
<input type="checkbox"/>	1221130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1221132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
Backrest Size			
	1221400	Specify Size (Width" x Height")	_____

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1223011	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1223012	Trigger on Left Back Post	
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. \$ 100			
<input type="checkbox"/>	1223021	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1223022	Trigger on Left Armrest	
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. \$ 195			
<input type="checkbox"/>	1223031	Trigger on Right Back Post	<input type="checkbox"/>
<input type="checkbox"/>	1223032	Trigger on Left Back Post	
<input type="checkbox"/>	1223041	Trigger on Right Armrest	<input type="checkbox"/>
<input type="checkbox"/>	1223042	Trigger on Left Armrest	
<input type="checkbox"/>	1223110	Trigger Extension	\$ 50
Tilt Limit			
<input type="checkbox"/>	1223220	Set at 20°	<input type="checkbox"/>
<input type="checkbox"/>	1223225	Set at 25°	<input type="checkbox"/>
<input type="checkbox"/>	1223230	Set at 30°	
<input type="checkbox"/>	1223235	Set at 35°	<input type="checkbox"/>
<input type="checkbox"/>	1223240	Set at 40°	<input type="checkbox"/>
<input type="checkbox"/>	1223245	Set at 45°	
<input type="checkbox"/>	1223410	Reclining Backrest (Manual)	E1225 \$ 875
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.			

Headrest			
<input type="checkbox"/>	1224000	No Headrest	
<input type="checkbox"/>	1224010	Headrest Mounting Plate	\$ 75
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	1224100	nxt MultiFit Head Support Hardware	E1028 \$ 155
nxt MultiFit Head Support Pad			
<input type="checkbox"/>	1224105	Regular (9" x 5")	<input type="checkbox"/>
<input type="checkbox"/>	1224106	Large (10" x 5")	<input type="checkbox"/>
<input type="checkbox"/>	1224107	X-Large (14" x 5")	
Stealth Headrest Hardware			
<input type="checkbox"/>	1224105	Multi-Axis HMO475	E1028 \$ 195
<input type="checkbox"/>	1224105	Link TWB480	
Stealth Headrest Pad			
<input type="checkbox"/>	1224105	10" width	<input type="checkbox"/>
<input type="checkbox"/>	1224106	14" width	

Front Rigging

<input type="checkbox"/>	1225000	No Front Rigging		
<input type="checkbox"/>	1225001	Sector Blocks only		
<input type="checkbox"/>	1225110	Plug In 70°		
<input type="checkbox"/>	1225120	Plug In Adjustable Knee Angle		\$ 975
<input type="checkbox"/>	1225130	Swing Away 60°		
<input type="checkbox"/>	1225140	Swing Away 70°		
<input type="checkbox"/>	1225200	Swing Away Elevating Legrest	E0990	\$ 375
<input type="checkbox"/>	1225300	1" Offset Hangers for Swing Away		\$ 75

Foot Plates

<input type="checkbox"/>	1226010	Composite		
<input type="checkbox"/>	1226020	Aluminum Angle Adjustable	K0040	\$ 155
<input type="checkbox"/>	1226030	One Piece Fixed Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.		\$ 350
<input type="checkbox"/>	1226040	One Piece Flip-up Requires selection of Composite or Aluminum Angle Adjustable.		\$ 395
<input type="checkbox"/>	1226050	Interlocking Angle Adjustable Seat Widths of 14"-21" only.	K0040	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	1227010	Swing Away 60° & 70° Calf Pad Not compatible with 60° & 70° Gel Padding.		\$ 135
<input type="checkbox"/>	1227020	Swing Away 60° & 70° Gel Padding Not compatible with 60° & 70° Swing Away Calf Pad.		\$ 100
<input type="checkbox"/>	1227030	Swing Away ELR Gel Knee Pad		\$ 100
<input type="checkbox"/>	1227040	Adjustable Length Heel Loops	E0951	\$ 50
<input type="checkbox"/>	1227050	Bodypoint - Padded Calf Strap Sizes greater than 22" are not padded.	K0038	\$ 105
Miller's Dynamic Footrest Mounts				
<input type="checkbox"/>	1227060	For Plug In Front Rigging		\$ 190
<input type="checkbox"/>	1227065	For Swing Away 70° Front Rigging Not available on 13" STF Height.		\$ 380

Extension Tubes

Regular Extension Tubes				
<input type="checkbox"/>	1228010	Standard (13"-17")		
<input type="checkbox"/>	1228020	Long (15"-19")		\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only.				
<input type="checkbox"/>	1228110	Standard (13"-17")	<input type="checkbox"/>	1228120 Long (15"-19")
				\$ 100

Positioning Straps

Bodypoint - Hip Belts				
<input type="checkbox"/>	1229010	Center-Pull Auto Style		E0978 \$ 105
<input type="checkbox"/>	1229015	Center-Pull Auto Style, Padded		E0978 \$ 155
<input type="checkbox"/>	1229020	4-Point, Padded		E0978 \$ 175
<input type="checkbox"/>	1229030	Bodypoint - Chest Strap		E0960 \$ 85

Miscellaneous Options			
<input type="checkbox"/>	1230050	Protective Lower Padding	\$ 425
<input type="checkbox"/>	1230080	Gel Hip Guides	<i>E0956 & E1028</i> \$ 325
<input type="checkbox"/>	1230110	IV Pole Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	<i>E1029</i> \$ 295
<input type="checkbox"/>	1230120	O2 Holder Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	<i>E2208</i> \$ 395
<input type="checkbox"/>	1230170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
 Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.



Scan this code or go to pdgmobility.com/products/stellar to learn more about the Stellar 45° Manual Tilt Wheelchair

