



USA Parts List

STELLAR HD

Effective January 1, 2018. Prices in U.S. Dollars.

Prices and descriptions contained in this manual are subject to change without prior notice.

Stellar HD Frame-Seat Pan

Frame Components

		<u>Retail Price</u>
16532	Stellar HD Lower Frame - specify width & color	\$1,200.00
16533	Stellar HD Upper Frame - specify width	\$1,200.00

Seat Pan

		<u>Retail Price</u>
18023	Stellar HD Seat Pan - 22"	\$150.00
18024	Stellar HD Seat Pan - 24"	\$150.00
18025	Stellar HD Seat Pan - 26"	\$150.00
18026	Stellar HD Seat Pan - 28"	\$150.00
18027	Stellar HD Seat Pan - 30"	\$150.00
18028	Stellar HD Seat Pan - 32"	\$150.00

Stellar HD Back Posts

Backpost Components

		<u>Retail Price</u>
5247	B/S Fold Down Bracket	\$63.00
15224	DRB Assy - PDG Bentley Stellar - Pair (New backpost req'd)	\$4,144.00

Pushbar Backpost 25"

		<u>Retail Price</u>
17890	B/S Push Bar 8° Bend Backpost 25" H x 22" W	\$176.00
17892	B/S Push Bar 8° Bend Backpost 25" H x 24" W	\$207.00
17894	B/S Push Bar 8° Bend Backpost 25" H x 26" W	\$207.00
17896	B/S Push Bar 8° Bend Backpost 25" H x 28" W	\$207.00
17898	B/S Push Bar 8° Bend Backpost 25" H x 30" W	\$178.00
17900	B/S Push Bar 8° Bend Backpost 25" H x 32" W	\$178.00

Pushbar Backpost 31"

		<u>Retail Price</u>
18297	B/S Push Bar 8° Bend Backpost 31" H x 22" W	\$207.00
18299	B/S Push Bar 8° Bend Backpost 31" H x 24" W	\$207.00
18301	B/S Push Bar 8° Bend Backpost 31" H x 26" W	\$207.00
18303	B/S Push Bar 8° Bend Backpost 31" H x 28" W	\$207.00
18305	B/S Push Bar 8° Bend Backpost 31" H x 30" W	\$207.00
18307	B/S Push Bar 8° Bend Backpost 31" H x 32" W	\$207.00

Reclining Back

		<u>Retail Price</u>
2722	B/S Recline back ASSEMBLY - Supply width	\$1,105.00
7301	Elastic Locking Gas Strut (60-228-250N)(add 5230 & 5231)	\$101.00
5230	Fuze T20 Strut Activator - Block	\$17.00

Stellar HD Back Upholstery

Back Upholstery

		<u>Retail Price</u>
1677	E600 upholstery 22"	\$487.00
1679	E600 upholstery 24"	\$487.00
1681	E600 upholstery 26"	\$487.00
1682	E600 upholstery 28"	\$557.00
1683	E600 upholstery 30"	\$557.00
1684	E600 upholstery 32"	\$557.00

Stellar HD Armrest

Arm Components

		<u>Retail Price</u>
3991	Armrest Latch Assy (Includes mounting hardware)	\$25.00
2880	Clothing Guard - 16" Deep - L	\$101.00
2881	Clothing Guard - 16" Deep - R	\$101.00
2882	Clothing Guard - 17" Deep - L	\$101.00
2883	Clothing Guard - 17" Deep - R	\$101.00
2884	Clothing Guard - 18" Deep - L	\$101.00
2885	Clothing Guard - 18" Deep - R	\$101.00
2886	Clothing Guard - 19" Deep - L	\$101.00
2887	Clothing Guard - 19" Deep - R	\$101.00
2888	Clothing Guard - 20" Deep - L	\$101.00
2889	Clothing Guard - 20" Deep - R	\$101.00
1502	Snap button (single button) 0.25"OD x 0.5"-1.25" Tube ID - Black Zinc	\$3.00

Arm Pads

		<u>Retail Price</u>
11406	Arm pad - 3" wide Desk Length 3" x 10" (EA)	\$140.00
11407	Arm pad - 3" wide Full Length 3" x 12" (EA)	\$140.00
14215	Arm pad - Contour Arm Trough w/Hand Pad (EA)	\$160.00
5102	Arm pad - Gel 2" x 12" (PR)	\$193.00
5101	Arm pad - Gel 3.5" x 12" (PR)	\$207.00
18001	Arm pad - Premium Grade Full Length	\$68.00
1505	Arm pad - Standard Desk Length 2" x 10" (EA)	\$21.00
1504	Arm pad - Standard Full Length 2" x 14" (EA)	\$20.00

Dual Post Arm Assembly

		<u>Retail Price</u>
3141	B/S Arm-Desk Length 10"-14" w/Pad & middle ASS'Y (PR)	\$373.00

3138	B/S Arm-Full Length 10"-14" w/Pad & middle ASS'Y (PR)	\$373.00
5248	B/S dual post middle arm rest (10"-16")	\$64.00
5250	B/S dual post upper arm - Desk Length (10-14")	\$70.00
5249	B/S dual post upper arm - Full Length (10-14")	\$70.00

Stellar HD Tilt

Tilt Components

		<u>Retail Price</u>
10246	Cable Clip	\$3.00
9840	Stellar gas strut ASSEMBLY complete (1000N)	\$173.00
1538	Stellar gas strut ASSEMBLY complete (1200N)	\$176.00
9839	Stellar gas strut ASSEMBLY complete (800N)	\$173.00
9398	Tilt Cable - 59" long	\$40.00
16654	Titl limiter Assembly	\$49.00
10906	Trigger Only (w/ mounting hardware)	\$19.00

Stellar HD Axle Plates & Anti-tippers

Anti-tipper

		<u>Retail Price</u>
16904	B/S Anti-Tipper (16"/22"/24" RW) - Assy - L	\$49.00
18226	B/S Anti-Tipper (16"/22"/24" RW) - Assy - R	\$49.00

Axle Plate & Axles

		<u>Retail Price</u>
14017	Stellar HD Axle Plate Assy - L	\$194.00
14018	Stellar HD Axle Plate Assy - R	\$194.00
13384	Stellar HD Fixed Axle 5/8-18(3A) x 4-5/16" GR5 C-HHCS - Clear Zinc	\$13.00

Stellar HD Rear Wheel

24" Rear Wheel

Retail Price

16090 24" HD Spoke Wheel, Urethane Tire, Standard Handrim

\$262.00

Wheel Components

Retail Price

2777 Bearings 5/8 ID x 1 3/8 OD (Magnesium Rim)

\$8.00

Stellar HD Casters

5" Caster

Retail Price

2139 5" Caster - Steel Fork Assy (Eclipse/Stellar HD Caster)

\$94.00

Stellar HD Wheel Lock

Wheel Lock Components

		<u>Retail Price</u>
1592	Wheel Lock Extention - Pair	\$81.00

Wheel Locks

		<u>Retail Price</u>
5245	B/S Wheel Lock Assembly - Pull - L	\$50.00
5246	B/S Wheel Lock Assembly - Pull - R	\$50.00
5243	B/S Wheel Lock Assembly - Push - L	\$50.00
5244	B/S Wheel Lock Assembly - Push - R	\$50.00

Stellar HD Front Rigging

Extension Tube and Footplate (Swing away)

		<u>Retail Price</u>
6921	Oversize AA Footplate Assy w/std lower ext. - L	\$147.00
6922	Oversize AA Footplate Assy w/std lower ext. - R	\$147.00

Extension Tubes

		<u>Retail Price</u>
6341	Extension Tube - SA (Channeled)(13"-17")	\$47.00

Footplate - Angle Adjustable

		<u>Retail Price</u>
17948	Footplate - Angle Adjustable - Oversized 20" & above (Not incl. clamp)	\$53.00

Front Rigging Components

		<u>Retail Price</u>
15173	5/16-18(3A) x 3/4" Soft Point Set Screw - Black Zinc	\$5.00
16330	Calf pad - bariatric (25" & wider)	\$682.00
16443	Calf Pad Assembly - Large (16"w & up)	\$29.00
17874	Calf Pad Mounting Bracket	\$8.00
7722	Calfpad Hanger Ass'y - L	\$35.00
7726	Calfpad Hanger Ass'y - R	\$35.00
16748	Clamp for Angle Adjustable Footplate	\$40.00
9325	Gel Knee Protector for ELR's (PR)	\$134.00
16761	Hanger Pivot Pin	\$33.00
10590	Heel Loop Adjustable Length (ea)	\$37.00
5971	Knuckle for Angle Adjustable Footplate	\$33.00
9274	PDG Swing Away Front Rigging Latch ASSY - L	\$50.00
9273	PDG Swing Away Front Rigging Latch ASSY - R	\$50.00
18034	Sector Block Assembly - L	\$16.00

Swing Away 60* Assembly AA Plate**Retail Price**

13718	PDG SA 60° Assy w/AA Footplate (20"w & up) - L	\$260.00
13719	PDG SA 60° Assy w/AA Footplate (20"w & up) - R	\$260.00

Swing Away ELR Assembly AA Plate**Retail Price**

13783	PDG SA ELR Assy w/AA Footplate (20"w & up) - L	\$420.00
13784	PDG SA ELR Assy w/AA Footplate (20"w & up) - R	\$420.00

Swing Away Upper Hangers**Retail Price**

10536	PDG SA 60° Assy, Upper Hanger - L	\$114.00
10537	PDG SA 60° Assy, Upper Hanger - R	\$114.00
9095	PDG SA ELR Assy, Upper Hanger - L	\$273.00
9094	PDG SA ELR Assy, Upper Hanger - R	\$273.00

Stellar HD Miscellaneous Options

Hardware

		<u>Retail Price</u>
13383	1/2" x 1/2" x 3/8-16(3A) SHSS - Black Zinc	\$6.00
1390	1/4-20(2B) Thick - Nylon Locknut - Black Zinc	\$2.00
1393	3/8-16(2B) Thick - Nylon Locknut - Black Zinc	\$2.00
17919	5/16-18(3A) x 2-1/2" with Shank CBHCS - Black Magni 575	\$2.00
9970	Oil Lite Bushing - Flanged - 1/2" ID x 5/8" OD x 1/2" L	\$7.00

Miscellaneous

		<u>Retail Price</u>
6584	Headrest Mounting Plate	\$107.00
3124	Nighthawk O2 Holder Assy (Specify Chair Type)	\$467.00

Stellar HD Seat Cushion

PDG Cushion

		<u>Retail Price</u>
1635	Cover only for PDG Basic Cushion	\$354.00
2594	Cover only for PDG Contour/Pressure reduction cushion	\$448.00
2280	Cushion - PDG basic, 14" - 28"	\$485.00
2281	Cushion - PDG basic, 29" and up	\$605.00
2282	Cushion - PDG contour, 14"-28"	\$835.00
2283	Cushion - PDG contour, 29" and up	\$1,090.00
2284	Cushion - PDG pressure relief, 14" - 28"	\$1,485.00
2285	Cushion - PDG pressure relief, 29" and up	\$1,900.00