

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type
<input checked="" type="checkbox"/> 1001000 Bentley 20° Manual Tilt Wheelchair \$ 3,995

Frame Options
<input type="checkbox"/> 1002000 Standard Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/> 1002050 Weight Capacity Upgrade \$ 200 Increases Weight Capacity to 300 lbs through addition of a secondary Gas Strut. Available for Seat Width of up to 22".
<input type="checkbox"/> 1002100 Heavy Duty Base and Seat Frame \$ 300 450 lbs Weight Capacity. Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Required for Seat Width of 23" or greater.
<input type="checkbox"/> 1002200 Extended Base with Rear Floor Jacks \$ 600 Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Incompatible with Attendant wheel lock, HD Anti-Tipper, IV Pole, O2 Holder.
<input type="checkbox"/> 1002300 Transit Tie Down System (TTDS) TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22".

Seat Width
Standard Seat Widths
<input type="checkbox"/> 1003014 14" Seat Width
<input type="checkbox"/> 1003015 15" Seat Width
<input type="checkbox"/> 1003016 16" Seat Width
<input type="checkbox"/> 1003017 17" Seat Width
<input type="checkbox"/> 1003018 18" Seat Width
<input type="checkbox"/> 1003019 19" Seat Width
<input type="checkbox"/> 1003020 20" Seat Width
<input type="checkbox"/> 1003021 21" Seat Width \$ 215
<input type="checkbox"/> 1003022 22" Seat Width \$ 215

Custom Seat Widths
Requires 1002100 Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.
<input type="checkbox"/> 1003023 23" Seat Width \$ 325
<input type="checkbox"/> 1003024 24" Seat Width \$ 325
<input type="checkbox"/> 1003025 25" Seat Width \$ 325
<input type="checkbox"/> 1003026 26" Seat Width \$ 325
<input type="checkbox"/> 1003027 27" Seat Width \$ 325
<input type="checkbox"/> 1003028 28" Seat Width \$ 625
<input type="checkbox"/> 1003029 29" Seat Width \$ 625
<input type="checkbox"/> 1003030 30" Seat Width \$ 625
<input type="checkbox"/> 1003031 31" Seat Width \$ 625
<input type="checkbox"/> 1003032 32" Seat Width \$ 625

Seat Depth
Standard Seat Depths (Adjustable)
<input type="checkbox"/> 1004016 16" Seat Depth
<input type="checkbox"/> 1004017 17" Seat Depth
<input type="checkbox"/> 1004018 18" Seat Depth
<input type="checkbox"/> 1004019 19" Seat Depth
<input type="checkbox"/> 1004020 20" Seat Depth
<input type="checkbox"/> 1004021 21" Seat Depth \$ 325
<input type="checkbox"/> 1004022 22" Seat Depth

Frame Colour
PDG Primary Colours
<input type="checkbox"/> 1005000 Tracer Red
<input type="checkbox"/> 1005001 Gloss Black
<input type="checkbox"/> 1005002 Silver
<input type="checkbox"/> 1005003 Black Shimmer
<input type="checkbox"/> 1005004 Orange Fusion
<input type="checkbox"/> 1005005 Karl's Imagination
<input type="checkbox"/> 1005006 Blue Shimmer
<input type="checkbox"/> 1005007 Gloss White
<input type="checkbox"/> 1005008 Dark Titanium
<input type="checkbox"/> 1005009 Geyser Green \$ 150
<input type="checkbox"/> 1005010 Lollipop Blue
<input type="checkbox"/> 1005011 Lollipop Red
<input type="checkbox"/> 1005200 Custom Colour \$ 300 Contact PDG for specifications.

Arm Type	
<input type="checkbox"/> 1006000	No Armrest
<input type="checkbox"/> 1006100	Single Post Fixed Height (12") Full Length Armpad only.
	Single Post Adjustable Height \$ 255 Not compatible with the Heavy Duty Base and/or Seat Width greater than 22".
<input type="checkbox"/> 1006111	8"-11"
<input type="checkbox"/> 1006114	10"-14"
<input type="checkbox"/> 1006200	Dual Post Fixed Height (12") Full Length Armpad only.
	Dual Post Adjustable Height \$ 255
<input type="checkbox"/> 1006211	8"-11"
<input type="checkbox"/> 1006214	10"-14"
<input type="checkbox"/> 1006216	12"-16"
	Locking Flip-up \$ 255 Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.
<input type="checkbox"/> 1006310	Set at 10"
<input type="checkbox"/> 1006311	Set at 11"
<input type="checkbox"/> 1006312	Set at 12"
<input type="checkbox"/> 1006313	Set at 13"
<input type="checkbox"/> 1006314	Set at 14"

Armrests	
	Regular Armpads
<input type="checkbox"/> 1007100	Desk Length - 2" x 10" (pair)
<input type="checkbox"/> 1007101	Full Length - 2" x 14" (pair)
	Wide Armpads \$ 115 May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1007110	Desk Length - 3" x 10" (pair)
<input type="checkbox"/> 1007111	Full Length - 3" x 12" (pair)
	Gel Armpads \$ 115 May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1007120	Full Length - 2" x 12" (pair)
<input type="checkbox"/> 1007121	Full Length - 3.5" x 12" (pair)
	Contour Arm Trough \$ 125 May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1007130	Right (4.5" x 18") with Hand Pad
<input type="checkbox"/> 1007131	Left (4.5" x 18") with Hand Pad
<input type="checkbox"/> 1007220	Armrest Wideners for Dual Post Armrests \$ 375 Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels or Clothing Guards.
	Clothing Guards Not compatible with Flip-up Armrests.
<input type="checkbox"/> 1007310	For Single Post Armrests \$ 110 10"-14" Height range only
<input type="checkbox"/> 1007320	For Dual Post Armrests \$ 165
	Armrest Extensions \$ 60 Not compatible with Flip-up or Single Post Armrests.
<input type="checkbox"/> 1007400	3" Extension - Forward
<input type="checkbox"/> 1007401	3" Extension - Rearward



Scan this code or visit www.pdgmobility.com to learn more about the Bentley Manual Tilt Wheelchair.

Seat Pan	
<input type="checkbox"/> 1008000	No Seat Pan
<input type="checkbox"/> 1008001	Solid Seat Pan \$ 120

Cushion	
<input type="checkbox"/> 1009000	No Cushion
	nxt BioFit
<input type="checkbox"/> 1009101	14"-20" width, 16"-20" depth \$ 595
<input type="checkbox"/> 1009102	22"-24" width, 16"-20" depth \$ 625
	Cushion Size Available in even number sizes only.
1009400	Specify Size (Width" x Depth") _____

Seat Height																																																																																		
Circle desired Seat-to-Floor Height or enter information below:																																																																																		
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="8">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> <tr> <td>20"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td></td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> </tbody> </table>	RW	6" Casters								12"	13"	14"	15"	16"	17"	18"			16"			15"	16"	17"	18"	19"	20"	20"	13"	14"	15"	16"	17"	18"			22"		14"	15"	16"	17"	18"	19"		24"			15"	16"	17"	18"	19"	20"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="8">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> </tbody> </table>	RW	8" Casters								16"			15"	16"	17"	18"	19"	20"	24"			15"	16"	17"	18"	19"	20"
RW	6" Casters																																																																																	
12"	13"	14"	15"	16"	17"	18"																																																																												
16"			15"	16"	17"	18"	19"	20"																																																																										
20"	13"	14"	15"	16"	17"	18"																																																																												
22"		14"	15"	16"	17"	18"	19"																																																																											
24"			15"	16"	17"	18"	19"	20"																																																																										
RW	8" Casters																																																																																	
16"			15"	16"	17"	18"	19"	20"																																																																										
24"			15"	16"	17"	18"	19"	20"																																																																										
Wheel size _____ Caster size _____ STF _____																																																																																		
STF Adjustment Range Select the desired range for STF adjustments. Acutal range depends on selected Wheel size.																																																																																		
<input type="checkbox"/> 1010005	Adjustment Range 13" to 17"																																																																																	
<input type="checkbox"/> 1010006	Adjustment Range 15" to 19"																																																																																	
<input type="checkbox"/> 1010007	Adjustment Range 15" to 20" Adds 1 lb to total weight of chair.																																																																																	

Casters	
	Urethane 8" Casters are compatible only with 16" and 24" Rear Wheels.
<input type="checkbox"/> 1011061	6" x 1.25" Caster size
<input type="checkbox"/> 1011063	6" x 1.50" Caster size \$ 85
<input type="checkbox"/> 1011081	8" x 1.25" Caster size
<input type="checkbox"/> 1011084	8" x 2.00" Caster size \$ 135
	Pneumatic 8" Casters are compatible only with 16" and 24" Rear Wheels.
<input type="checkbox"/> 1011161	6" x 1.25" Caster size \$ 110
<input type="checkbox"/> 1011181	8" x 1.25" Caster size \$ 120
<input type="checkbox"/> 1011184	8" x 2.00" Caster size \$ 150
<input type="checkbox"/> 1011260	6" Soft Roll \$ 95
<input type="checkbox"/> 1011400	Ultra Sport Fork \$ 395 Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.

Rear Wheels	
Composite	
<input type="checkbox"/> 1012012	12" Wheel size
<input type="checkbox"/> 1012016	16" Wheel size
<input type="checkbox"/> 1012020	20" Wheel size
<input type="checkbox"/> 1012022	22" Wheel size
<input type="checkbox"/> 1012024	24" Wheel size
Spoke \$ 75	
<input type="checkbox"/> 1012120	20" Wheel size
<input type="checkbox"/> 1012122	22" Wheel size
<input type="checkbox"/> 1012124	24" Wheel size
<input type="checkbox"/> 1012224	24" Heavy Duty Spoke \$ 350 Required for Heavy Duty Package.

Axle Type	
<input type="checkbox"/> 1013000	Quick Release Rear Axles
<input type="checkbox"/> 1013010	Fixed Rear Axles

Rear Tires	
<input type="checkbox"/> 1014000	Urethane
<input type="checkbox"/> 1014010	Pneumatic \$ 60 Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.
<input type="checkbox"/> 1014020	Pneumatic with Flat-Free Inserts \$ 95 Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.

Hand Rims	
<input type="checkbox"/> 1015000	No Hand Rims
<input type="checkbox"/> 1015010	Standard 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 1015020	Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks	
<input type="checkbox"/> 1016010	Push to Lock
<input type="checkbox"/> 1016020	Pull to Lock
<input type="checkbox"/> 1016030	Attendant Wheel Lock (Push to Lock) \$ 85 Not compatible with TTDS (23" & wider)
<input type="checkbox"/> 1016040	User & Attendant Wheel Lock \$ 125 Available only on 20", 22", or 24" Rear Wheels.
<input type="checkbox"/> 1016100	Wheel Lock Extensions \$ 60 Not compatible with Attendant Wheel Lock.
<input type="checkbox"/> 1016200	Hill Holders \$ 350 Requires Pneumatic Rear Tires with Flat-Free Inserts.

Anti-Tippers	
<input checked="" type="checkbox"/> 1017010	Rear Anti-Tippers \$ 60 Required Safety Option.
<input type="checkbox"/> 1017020	Front Anti-Tippers \$ 300
<input type="checkbox"/> 1017030	Heavy Duty Rear Anti-Tipper \$ 600 Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.

Back Post Style	
<input type="checkbox"/> 1018010	Push Bar Back
<input type="checkbox"/> 1018020	Cane Handle Back Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.
Back Post Height Seat Pan to Handle	
<input type="checkbox"/> 1018425	25" Height
<input type="checkbox"/> 1018431	31" Height

Back Options	
<input type="checkbox"/> 1019010	Angle/Height Adjustable Push Bar Handle \$ 150 Not compatible with Headrest Mounting Plate or Cane Handle Back.
<input type="checkbox"/> 1019010	Button on Outside
<input type="checkbox"/> 1019011	Button on Inside
Seat to Back Angle	
<input type="checkbox"/> 1019090	Set at 90°
<input type="checkbox"/> 1019098	Set at 98°
<input type="checkbox"/> 1019106	Set at 106°
<input type="checkbox"/> 1019114	Set at 114°
<input type="checkbox"/> 1019122	Set at 122°
<input type="checkbox"/> 1019210	Backrest Extension 10" high \$ 120 Not compatible with Angle/Height Adjustable Push Bar Handle.
<input type="checkbox"/> 1019220	Fold Down Back Not compatible with Reclining Back.
<input checked="" type="checkbox"/> 1019230	Spring Back Suspension Not compatible with Reclining Back.
<input type="checkbox"/> 1019240	Dynamic Rocker Back Assembly \$ 2,995 Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller. Requires longer Frame.

Back Upholstery	
<input type="checkbox"/> 1020000	No Back Upholstery
<input type="checkbox"/> 1020021	Standard \$ 125 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1020021	Set at 21"
<input type="checkbox"/> 1020022	Set at 22"
<input type="checkbox"/> 1020118	Short (25" Back only) \$ 125 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1020118	Set at 18"
<input type="checkbox"/> 1020119	Set at 19"
<input type="checkbox"/> 1020120	Set at 20"
<input type="checkbox"/> 1020224	Tall (31" Back only) \$ 125 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1020224	Set at 24"
<input type="checkbox"/> 1020225	Set at 25"
<input type="checkbox"/> 1020226	Set at 26"
<input type="checkbox"/> 1020321	Adjustable Tension \$ 250 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1020321	Set at 21"
<input type="checkbox"/> 1020322	Set at 22"
<input type="checkbox"/> 1020323	Set at 23"

Backrest			
<input type="checkbox"/>	1021000	No Backrest	
		nxt Optima HA Active Contour	
<input type="checkbox"/>	1021100	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1021105	22"-24" width, 16"-20" height	\$ 720
		nxt Optima HA Plus Contour	
<input type="checkbox"/>	1021101	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1021106	22"-24" width, 16"-20" height	\$ 720
		Backrest Size	
		Available in even number sizes only.	
<input type="checkbox"/>	1021400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
		Attendant Tilt only	
<input type="checkbox"/>	1023011	Installed on Right Back Post	
<input type="checkbox"/>	1023012	Installed on Left Back Post	
		User Tilt only	\$ 100
		Not available on widths greater than 22". Full Length Arm Pads only.	
<input type="checkbox"/>	1023021	Installed on Right Armrest	
<input type="checkbox"/>	1023022	Installed on Left Armrest	
		Attendant & User Tilt	\$ 195
		Not available on widths greater than 22". Full Length Arm Pads only.	
<input type="checkbox"/>	1023031	Back Post mount - Right side	
<input type="checkbox"/>	1023032	Back Post mount - Left side	
<input type="checkbox"/>	1023041	Armrest mount - Right side	
<input type="checkbox"/>	1023042	Armrest mount - Left side	
<input type="checkbox"/>	1023110	Trigger Extension	\$ 50
<input type="checkbox"/>	1023410	Reclining Backrest (Manual)	\$ 830
		Maximum recline = 30°. Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, 20"-22" seat depth. Requires longer Frame.	

Headrest			
<input type="checkbox"/>	1024000	No Headrest	
<input type="checkbox"/>	1024010	Headrest Mounting Plate	\$ 75
		Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	
<input type="checkbox"/>	1024100	nxt MultiFit Head Support Hardware	\$ 155
		nxt MultiFit Head Support Pad	\$ 195
<input type="checkbox"/>	1024105	Regular (9" x 5")	
<input type="checkbox"/>	1024106	Large (10" x 5")	
<input type="checkbox"/>	1024107	X-Large (14" x 5")	

Front Rigging			
<input type="checkbox"/>	1025000	No Front Rigging	
<input type="checkbox"/>	1025001	Sector Blocks only	
<input type="checkbox"/>	1025110	Plug In 70°	
<input type="checkbox"/>	1025120	Plug In Adjustable Knee Angle	\$ 975
<input type="checkbox"/>	1025130	Swing Away 60°	
<input type="checkbox"/>	1025140	Swing Away 70°	
<input type="checkbox"/>	1025200	Swing Away Elevating Legrest	\$ 355

Foot Plates			
<input type="checkbox"/>	1026010	Composite	
<input type="checkbox"/>	1026020	Angle Adjustable	\$ 140
<input type="checkbox"/>	1026030	One Piece Fixed	\$ 350
		Available on Plug In Front Rigging only.	
<input type="checkbox"/>	1026040	One Piece Flip-up	\$ 395
<input type="checkbox"/>	1026050	Interlocking Angle Adjustable	\$ 195
		Seat Widths of 14"-21" only.	

Footrest/Legrest Accessories			
<input type="checkbox"/>	1027010	Swing Away 60° & 70° Calf Pad	\$ 135
		Not compatible with 60° & 70° Gel Pad.	
<input type="checkbox"/>	1027020	Swing Away 60° & 70° Gel Pad	\$ 100
		Not compatible with 60° & 70° Swing Away Calf Pad.	
<input type="checkbox"/>	1027030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1027040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	1027050	Bodypoint - Padded Calf Strap	\$ 105
		Sizes greater than 22" are not padded.	
		Miller's Dynamic Footrest Mounts	
<input type="checkbox"/>	1027060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging	\$ 380
		Not available on 13" STF	

Extension Tubes			
		Regular Extension Tubes	
<input type="checkbox"/>	1028010	Standard (13"-17")	
<input type="checkbox"/>	1028020	Long (15"-19")	\$ 75
		Extra Strength Extension Tubes (Solid)	\$ 100
		Available on Plug In Front Rigging only.	
<input type="checkbox"/>	1028110	Standard (13"-17")	
<input type="checkbox"/>	1028120	Long (15"-19")	

Positioning Straps			
		Bodypoint - Hip Belts	
<input type="checkbox"/>	1029010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1029015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1029020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1029030	Bodypoint - Chest Strap	\$ 85

Miscellaneous Options			
		One-Arm Drive	\$ 3,500
		Maximum width of 22". 24" Composite Wheels required. 18" STF minimum. Limits tilt.	
<input type="checkbox"/>	1030010	Installed on Right side	
<input type="checkbox"/>	1030011	Installed on Left side	
<input type="checkbox"/>	1030050	Protective Lower Padding	\$ 425
<input type="checkbox"/>	1030110	IV Pole	\$ 295
		Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	
<input type="checkbox"/>	1030120	O2 Holder	\$ 295
		Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	
<input type="checkbox"/>	1030170	Adhesive Back Hook Strip, 2" x 30"	\$ 10
		For Cushion mounting.	