



Arm Type			
<input type="checkbox"/>	1006000	<b>No Armrest</b>	
<input type="checkbox"/>	1006100	<b>Single Post Fixed Height (12")</b>	
		Full Length Armpad only.	
		<b>Single Post Adjustable Height</b>	<i>E0973</i> \$ 255
		Not compatible with the Heavy Duty Package and Seat Widths greater than 22".	
<input type="checkbox"/>	1006111	8"-11"	
<input type="checkbox"/>	1006114	10"-14"	
<input type="checkbox"/>	1006200	<b>Dual Post Fixed Height (12")</b>	
		Full Length Armpad only.	
		<b>Dual Post Adjustable Height</b>	<i>E0973</i> \$ 255
<input type="checkbox"/>	1006211	8"-11"	
<input type="checkbox"/>	1006214	10"-14"	
<input type="checkbox"/>	1006216	12"-16"	
		<b>Locking Flip-up</b>	<i>E0973</i> \$ 255
		Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.	
<input type="checkbox"/>	1006310	Set at 10"	
<input type="checkbox"/>	1006311	Set at 11"	
<input type="checkbox"/>	1006312	Set at 12"	
<input type="checkbox"/>	1006313	Set at 13"	
<input type="checkbox"/>	1006314	Set at 14"	

Armrests			
		<b>Regular Armpads</b>	
<input type="checkbox"/>	1007100	Desk Length - 2" x 10" (pair)	
<input type="checkbox"/>	1007101	Full Length - 2" x 14" (pair)	
		<b>Wide Armpads</b>	\$ 115
		May interfere with wheels in some tilt configurations.	
<input type="checkbox"/>	1007110	Desk Length - 3" x 10" (pair)	
<input type="checkbox"/>	1007111	Full Length - 3" x 12" (pair)	
		<b>Gel Armpads</b>	\$ 115
		May interfere with wheels in some tilt configurations.	
<input type="checkbox"/>	1007120	Full Length - 2" x 12" (pair)	
<input type="checkbox"/>	1007121	Full Length - 3.5" x 12" (pair)	
		<b>Contour Arm Trough</b>	<i>E2209</i> \$ 125
		May interfere with wheels in some tilt configurations.	
<input type="checkbox"/>	1007130	Right (4.5" x 18") with Hand Pad	
<input type="checkbox"/>	1007131	Left (4.5" x 18") with Hand Pad	
<input type="checkbox"/>	1007220	<b>Armrest Wideners for Dual Post Armrests</b>	\$ 375
		Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels or Clothing Guards.	
		<b>Clothing Guards</b>	<i>K0108</i>
		Not compatible with Flip-up Armrests.	
<input type="checkbox"/>	1007310	For Single Post Armrests	\$ 110
		10"-14" Height range only	
<input type="checkbox"/>	1007320	For Dual Post Armrests	\$ 165
		<b>Armrest Extensions</b>	\$ 60
		Not compatible with Flip-up or Single Post Armrests.	
<input type="checkbox"/>	1007400	3" Extension - Forward	
<input type="checkbox"/>	1007401	3" Extension - Rearward	



Scan this code or visit [www.pdgmobility.com](http://www.pdgmobility.com) to learn more about the Bentley Manual Tilt Wheelchair.

Seat Pan			
<input type="checkbox"/>	1008000	<b>No Seat Pan</b>	
<input type="checkbox"/>	1008001	<b>Solid Seat Pan</b>	<i>E2231</i> \$ 120

Cushion			
<input type="checkbox"/>	1009000	<b>No Cushion</b>	
		<b>next BioFit</b>	
<input type="checkbox"/>	1009101	14"-20" width, 16"-20" depth	<i>E2607</i> \$ 550
<input type="checkbox"/>	1009102	22"-24" width, 16"-20" depth	<i>E2608</i> \$ 600
		<b>Matrx VI</b>	
<input type="checkbox"/>	1009111	14"-20" width, 16"-22" depth	<i>E2607</i> \$ 550
<input type="checkbox"/>	1009112	22" width, 18"-22" depth	<i>E2608</i> \$ 600
<input type="checkbox"/>	1009115	HD - 18"-30" width, 18"-22" depth	<i>E2608</i> \$ 1,150
		<b>ROHO Quadtro Select HP</b>	
<input type="checkbox"/>	1009121	14"-20" width, 16"-22" depth	<i>E2624</i> \$ 650
<input type="checkbox"/>	1009122	22"-24" width, 16"-22" depth	<i>E2625</i> \$ 750
		<b>ROHO High Profile</b>	
<input type="checkbox"/>	1009125	14"-20" width, 16"-22" depth	<i>E2622</i> \$ 650
<input type="checkbox"/>	1009126	22"-24" width, 16"-22" depth	<i>E2623</i> \$ 750
		<b>Varilite Reflex</b>	
<input type="checkbox"/>	1009131	14"-20" width, 16"-20" depth	<i>E2603</i> \$ 375
		<b>Cushion Size</b>	
	1009400	Specify Size (Width" x Depth") _____	

Seat Height																			
Circle desired Seat-to-Floor Height or enter information below:																			
<b>RW</b>	<b>6" Casters</b>							<b>RW</b>	<b>8" Casters</b>										
<b>12"</b>	13"	14"	15"	16"	17"	18"		<b>16"</b>			15"	16"	17"	18"	19"	20"			
<b>16"</b>			15"	16"	17"	18"	19"	20"	<b>20"</b>	13"	14"	15"	16"	17"	18"				
<b>20"</b>			14"	15"	16"	17"	18"	19"	<b>22"</b>			15"	16"	17"	18"	19"			
<b>24"</b>				15"	16"	17"	18"	19"	20"	<b>24"</b>				15"	16"	17"	18"	19"	20"
Wheel size _____ Caster size _____ STF _____																			
<b>STF Adjustment Range</b>																			
Select the desired range for STF adjustments. Acutal range depends on selected Wheel size.																			
<input type="checkbox"/>	1010005	Adjustment Range 13" to 17"																	
<input type="checkbox"/>	1010006	Adjustment Range 15" to 19"																	
<input type="checkbox"/>	1010007	Adjustment Range 15" to 20"	Adds 1 lb to total weight of chair.																

Casters	
<b>Urethane</b>	
8" Casters are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/> 1011061	6" x 1.25" Caster size
<input type="checkbox"/> 1011063	6" x 1.50" Caster size \$ 85
<input type="checkbox"/> 1011081	8" x 1.25" Caster size
<input type="checkbox"/> 1011084	8" x 2.00" Caster size \$ 135
<b>Pneumatic</b>	
8" Casters are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/> 1011161	6" x 1.25" Caster size \$ 110
<input type="checkbox"/> 1011181	8" x 1.25" Caster size \$ 120
<input type="checkbox"/> 1011184	8" x 2.00" Caster size \$ 150
<input type="checkbox"/> 1011260	<b>6" Soft Roll</b> \$ 95
<input type="checkbox"/> 1011400	<b>Ultra Sport Fork</b> E1015 \$ 395
Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	

Rear Wheels	
<b>Composite</b>	
<input type="checkbox"/> 1012012	12" Wheel size
<input type="checkbox"/> 1012016	16" Wheel size
<input type="checkbox"/> 1012020	20" Wheel size
<input type="checkbox"/> 1012022	22" Wheel size
<input type="checkbox"/> 1012024	24" Wheel size
<b>Spoke</b> \$ 75	
<input type="checkbox"/> 1012120	20" Wheel size
<input type="checkbox"/> 1012122	22" Wheel size
<input type="checkbox"/> 1012124	24" Wheel size
<input type="checkbox"/> 1012224	<b>24" Heavy Duty Spoke</b> \$ 350
Required for Heavy Duty Package.	

Axle Type	
<input type="checkbox"/> 1013000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/> 1013010	<b>Fixed Rear Axles</b>

Rear Tires	
<input type="checkbox"/> 1014000	<b>Urethane</b>
<input type="checkbox"/> 1014010	<b>Pneumatic</b> \$ 60
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.	
<input type="checkbox"/> 1014020	<b>Pneumatic with Flat-Free Inserts</b> E2213 \$ 95
Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.	

Hand Rims	
<input type="checkbox"/> 1015000	<b>No Hand Rims</b>
<input type="checkbox"/> 1015010	<b>Standard</b> 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 1015020	<b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks	
<input type="checkbox"/> 1016010	<b>Push to Lock</b>
<input type="checkbox"/> 1016020	<b>Pull to Lock</b>
<input type="checkbox"/> 1016030	<b>Attendant Wheel Lock (Push to Lock)</b> \$ 85
Not compatible with TTDS (23" & wider)	
<input type="checkbox"/> 1016040	<b>User &amp; Attendant Wheel Lock</b> \$ 125
Available only on 20", 22", or 24" Rear Wheels.	
<input type="checkbox"/> 1016100	<b>Wheel Lock Extensions</b> E0961 \$ 60
Not compatible with Attendant Wheel Lock.	
<input type="checkbox"/> 1016200	<b>Hill Holders</b> \$ 350
Requires Pneumatic Rear Tires with Flat-Free Inserts.	

Anti-Tippers	
<input checked="" type="checkbox"/> 1017010	<b>Rear Anti-Tippers</b> E0971 \$ 60
Required Safety Option.	
<input type="checkbox"/> 1017020	<b>Front Anti-Tippers</b> \$ 300
<input type="checkbox"/> 1017030	<b>Heavy Duty Rear Anti-Tipper</b> \$ 600
Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	

Back Post Style	
<input type="checkbox"/> 1018010	<b>Push Bar Back</b>
<input type="checkbox"/> 1018020	<b>Cane Handle Back</b>
Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
<b>Back Post Height</b>	
Seat Pan to Handle	
<input type="checkbox"/> 1018425	25" Height
<input type="checkbox"/> 1018431	31" Height

Back Options	
<b>Angle/Height Adjustable Push Bar Handle</b> \$ 150	
Not compatible with Headrest Mounting Plate or Cane Handle Back.	
<input type="checkbox"/> 1019010	Button on Outside
<input type="checkbox"/> 1019011	Button on Inside
<b>Seat to Back Angle</b>	
<input type="checkbox"/> 1019090	Set at 90°
<input type="checkbox"/> 1019098	Set at 98°
<input type="checkbox"/> 1019106	Set at 106°
<input type="checkbox"/> 1019114	Set at 114°
<input type="checkbox"/> 1019122	Set at 122°
<input type="checkbox"/> 1019210	<b>Backrest Extension 10" high</b> \$ 120
Not compatible with Angle/Height Adjustable Push Bar Handle.	
<input type="checkbox"/> 1019220	<b>Fold Down Back</b>
Not compatible with Reclining Back.	
<input checked="" type="checkbox"/> 1019230	<b>Spring Back Suspension</b>
Not compatible with Reclining Back.	
<input type="checkbox"/> 1019240	<b>Dynamic Rocker Back Assembly</b> K0108 \$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller. Requires longer Frame.	

Back Upholstery			
<input type="checkbox"/>	1020000	<b>No Back Upholstery</b>	
		<b>Standard</b>	\$ 125
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1020021	Set at 21"	
<input type="checkbox"/>	1020022	Set at 22"	
		<b>Short (25" Back only)</b>	\$ 125
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1020118	Set at 18"	
<input type="checkbox"/>	1020119	Set at 19"	
<input type="checkbox"/>	1020120	Set at 20"	
		<b>Tall (31" Back only)</b>	\$ 125
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1020224	Set at 24"	
<input type="checkbox"/>	1020225	Set at 25"	
<input type="checkbox"/>	1020226	Set at 26"	
		<b>Adjustable Tension</b>	\$ 250
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1020321	Set at 21"	
<input type="checkbox"/>	1020322	Set at 22"	
<input type="checkbox"/>	1020323	Set at 23"	

Backrest			
<input type="checkbox"/>	1021000	<b>No Backrest</b>	
		<b>nxt Optima HA</b>	
<input type="checkbox"/>	1021100	AC - 16"-20" width, 16"-20" height E2620	\$ 690
<input type="checkbox"/>	1021101	PC - 16"-20" width, 16"-20" height E2620	\$ 690
<input type="checkbox"/>	1021105	AC - 22"-24" width, 16"-20" height E2621	\$ 720
<input type="checkbox"/>	1021106	PC - 22"-24" width, 16"-20" height E2621	\$ 720
		<b>Matrx Elite</b>	
<input type="checkbox"/>	1021110	14"-20" width, 10"-20" height E2620	\$ 825
<input type="checkbox"/>	1021112	HD - 22"-26" width, 16"-20" height E2621	\$ 1,095
<input type="checkbox"/>	1021114	HD - 26"-30" width, 16"-20" height E2621	\$ 1,595
<input type="checkbox"/>	1021116	Deep - 14"-20" width, 12"-20" height E2620	\$ 950
		<b>ROHO Agility</b>	
<input type="checkbox"/>	1021120	Min - 14"-20" width, 10"-18" height E2620	\$ 825
<input type="checkbox"/>	1021122	Mid - 14"-20" width, 10"-18" height E2620	\$ 850
<input type="checkbox"/>	1021124	Max - 14"-20" width, 10"-18" height E2620	\$ 950
		<b>Varilite Icon</b>	
<input type="checkbox"/>	1021130	Tall - 14"-20" width, 14"-24" height E2615	\$ 850
<input type="checkbox"/>	1021132	Deep - 14"-20" width, 17"-24" height E2620	\$ 895
		<b>Backrest Size</b>	
	1021400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
		<b>Attendant Tilt only</b>	
<input type="checkbox"/>	1023011	Installed on Right Back Post	
<input type="checkbox"/>	1023012	Installed on Left Back Post	
		<b>User Tilt only</b>	\$ 100
		Not available on widths greater than 22". Full Length Arm Pads only.	
<input type="checkbox"/>	1023021	Installed on Right Armrest	
<input type="checkbox"/>	1023022	Installed on Left Armrest	
		<b>Attendant &amp; User Tilt</b>	\$ 195
		Not available on widths greater than 22". Full Length Arm Pads only.	
<input type="checkbox"/>	1023031	Back Post mount - Right side	
<input type="checkbox"/>	1023032	Back Post mount - Left side	
<input type="checkbox"/>	1023041	Armrest mount - Right side	
<input type="checkbox"/>	1023042	Armrest mount - Left side	
<input type="checkbox"/>	1023110	<b>Trigger Extension</b>	\$ 50
<input type="checkbox"/>	1023410	<b>Reclining Backrest (Manual)</b> E1225	\$ 830
		Maximum recline = 30°. Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, 20"-22" seat depth. Requires longer Frame.	

Headrest			
<input type="checkbox"/>	1024000	<b>No Headrest</b>	
<input type="checkbox"/>	1024010	<b>Headrest Mounting Plate</b>	\$ 75
		Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	
<input type="checkbox"/>	1024100	<b>nxt MultiFit Head Support Hardware</b> E1028	\$ 155
		<b>nxt MultiFit Head Support Pad</b> E0955	\$ 195
<input type="checkbox"/>	1024105	Regular (9" x 5")	
<input type="checkbox"/>	1024106	Large (10" x 5")	
<input type="checkbox"/>	1024107	X-Large (14" x 5")	
		<b>Stealth Headrest Hardware</b> E1028	
<input type="checkbox"/>	1024105	Multi-Axis HMO475	\$ 150
<input type="checkbox"/>	1024106	Link TWB480	\$ 325
		<b>Stealth Headrest Pad</b> E0955	\$ 200
<input type="checkbox"/>	1024105	10" width	
<input type="checkbox"/>	1024106	14" width	

Front Rigging			
<input type="checkbox"/>	1025000	<b>No Front Rigging</b>	
<input type="checkbox"/>	1025001	<b>Sector Blocks only</b>	
<input type="checkbox"/>	1025110	<b>Plug In 70°</b>	
<input type="checkbox"/>	1025120	<b>Plug In Adjustable Knee Angle</b>	\$ 975
<input type="checkbox"/>	1025130	<b>Swing Away 60°</b>	
<input type="checkbox"/>	1025140	<b>Swing Away 70°</b>	
<input type="checkbox"/>	1025200	<b>Swing Away Elevating Legrest</b> E0990	\$ 355

Foot Plates			
<input type="checkbox"/>	1026010	<b>Composite</b>	
<input type="checkbox"/>	1026020	<b>Angle Adjustable</b> <i>K0040</i>	\$ 140
<input type="checkbox"/>	1026030	<b>One Piece Fixed</b> Available on Plug In Front Rigging only.	\$ 350
<input type="checkbox"/>	1026040	<b>One Piece Flip-up</b>	\$ 395
<input type="checkbox"/>	1026050	<b>Interlocking Angle Adjustable</b> Seat Widths of 14"-21" only.	\$ 195

Footrest/Legrest Accessories			
<input type="checkbox"/>	1027010	<b>Swing Away 60° &amp; 70° Calf Pad</b> Not compatible with 60° & 70° Gel Pad.	\$ 135
<input type="checkbox"/>	1027020	<b>Swing Away 60° &amp; 70° Gel Pad</b> Not compatible with 60° & 70° Swing Away Calf Pad.	\$ 100
<input type="checkbox"/>	1027030	<b>Swing Away ELR Gel Knee Pad</b>	\$ 100
<input type="checkbox"/>	1027040	<b>Adjustable Length Heel Loops</b> <i>E0951</i>	\$ 50
<input type="checkbox"/>	1027050	<b>Bodypoint - Padded Calf Strap</b> <i>K0038</i> Sizes greater than 22" are not padded.	\$ 105
<b>Miller's Dynamic Footrest Mounts</b>			
<input type="checkbox"/>	1027060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging Not available on 13" STF	\$ 380

Extension Tubes			
<b>Regular Extension Tubes</b>			
<input type="checkbox"/>	1028010	Standard (13"-17")	
<input type="checkbox"/>	1028020	Long (15"-19")	\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> <span style="float: right;">\$ 100</span>			
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1028110	Standard (13"-17")	
<input type="checkbox"/>	1028120	Long (15"-19")	

Positioning Straps			
<b>Bodypoint - Hip Belt</b> <i>E0978</i>			
<input type="checkbox"/>	1029010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1029015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1029020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1029030	<b>Bodypoint - Chest Strap</b> <i>E0960</i>	\$ 85

Miscellaneous Options			
<b>One-Arm Drive</b> <i>E0958</i> <span style="float: right;">\$ 3,500</span>			
Maximum width of 22". 24" Composite Wheels required. 18" STF minimum. Limits tilt.			
<input type="checkbox"/>	1030010	Installed on Right side	
<input type="checkbox"/>	1030011	Installed on Left side	
<input type="checkbox"/>	1030050	<b>Protective Lower Padding</b>	\$ 425
<input type="checkbox"/>	1030110	<b>IV Pole</b> <i>E1029</i> <span style="float: right;">\$ 295</span>	
Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			
<input type="checkbox"/>	1030120	<b>O2 Holder</b> <i>E2208</i> <span style="float: right;">\$ 295</span>	
Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			
<input type="checkbox"/>	1030170	<b>Adhesive Back Hook Strip, 2" x 30"</b> <span style="float: right;">\$ 10</span>	
For Cushion mounting.			