

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type

Elevation Ultralight Manual Wheelchair

- | | | |
|-------------------------------------|---|----------|
| <input checked="" type="checkbox"/> | 2001000 Ultralight Manual Rigid Wheelchair | \$ 3,900 |
| <input checked="" type="checkbox"/> | 2001001 Patient Controlled Anterior/Posterior Back Angle Recline System | \$ 850 |
| <input checked="" type="checkbox"/> | 2001002 Patient Controlled Anterior Tilt & Elevation System | \$ 1,800 |

Frame Options

- | | | |
|---|---|--|
| <input checked="" type="checkbox"/> | 2002000 Standard Base and Seat Frame
250 lbs Weight Capacity. | |
| User Weight Range
Required for selection of Gas Struts. | | |
| <input type="checkbox"/> | 2002150 Up to 90 lbs (150 Newton Gas Struts) | |
| <input type="checkbox"/> | 2002200 90 lbs to 140 lbs (200 Newton Gas Struts) | |
| <input type="checkbox"/> | 2002300 140 lbs to 200 lbs (300 Newton Gas Struts) | |
| <input type="checkbox"/> | 2002400 200 lbs to 250 lbs (400 Newton Gas Struts) | |

Seat Width

Standard Seat Widths

- | | |
|--------------------------|------------------------|
| <input type="checkbox"/> | 2003012 12" Seat Width |
| <input type="checkbox"/> | 2003013 13" Seat Width |
| <input type="checkbox"/> | 2003014 14" Seat Width |
| <input type="checkbox"/> | 2003015 15" Seat Width |
| <input type="checkbox"/> | 2003016 16" Seat Width |
| <input type="checkbox"/> | 2003017 17" Seat Width |
| <input type="checkbox"/> | 2003018 18" Seat Width |
| <input type="checkbox"/> | 2003019 19" Seat Width |

Seat Depth

Standard Seat Depths

- | | |
|--------------------------|------------------------|
| <input type="checkbox"/> | 2004014 14" Seat Depth |
| <input type="checkbox"/> | 2004015 15" Seat Depth |
| <input type="checkbox"/> | 2004016 16" Seat Depth |
| <input type="checkbox"/> | 2004017 17" Seat Depth |
| <input type="checkbox"/> | 2004018 18" Seat Depth |
| <input type="checkbox"/> | 2004019 19" Seat Depth |
| <input type="checkbox"/> | 2004020 20" Seat Depth |

Frame Colour

PDG Primary Colours

- | | |
|--------------------------|----------------------------|
| <input type="checkbox"/> | 2005000 Tracer Red |
| <input type="checkbox"/> | 2005001 Gloss Black |
| <input type="checkbox"/> | 2005002 Silver |
| <input type="checkbox"/> | 2005003 Black Shimmer |
| <input type="checkbox"/> | 2005004 Orange Fusion |
| <input type="checkbox"/> | 2005005 Karl's Imagination |
| <input type="checkbox"/> | 2005006 Blue Shimmer |
| <input type="checkbox"/> | 2005007 Gloss White |
| <input type="checkbox"/> | 2005008 Dark Titanium |

PDG Premium Colours

- | | | |
|--------------------------|-----------------------|-----|
| <input type="checkbox"/> | 2005009 Geyser Green | N/C |
| <input type="checkbox"/> | 2005010 Lollipop Blue | |
| <input type="checkbox"/> | 2005011 Lollipop Red | |

- | | | |
|---------------------------------|------------------------------|--------|
| <input type="checkbox"/> | 2005200 Custom Colour | \$ 300 |
| Contact PDG for specifications. | | |

Arm Type

- | | |
|--------------------------|------------------------------|
| <input type="checkbox"/> | 2006000 No Armrest |
| <input type="checkbox"/> | 2006001 Rigid Armrest |

Seat Option

- | | |
|---|------------------------------------|
| <input type="checkbox"/> | 2008110 Standard Seat Sling |
| <input type="checkbox"/> | 2008100 No Seat Sling |
| If omitted, third-party Seat Sling must be installed by dealer. | |

Cushion

- | | |
|-------------------------------------|---------------------------|
| <input checked="" type="checkbox"/> | 2009000 No Cushion |
|-------------------------------------|---------------------------|

Seat Height

24" Wheel Size

- 2010018 18" STF
- 2010019 19" STF
- 2010020 20" STF

25" Wheel Size

Not compatible with 4" Caster with Frog Leg Uni-tine Fork

- 2010118 18.5" STF
- 2010119 19.5" STF
- 2010120 20.5" STF

26" Wheel Size

Not compatible with Frog Leg Uni-tine Fork

- 2010219 19" STF
- 2010220 20" STF
- 2010221 21" STF

Casters

Composite

- 2011040 4" x 1.0" Composite Sport
- 2011042 4" x 1.5" Composite Softroll
- 2011052 5" x 1.5" Composite Softroll

Aluminium

\$ 125

- 2011140 4" x 1.0" Aluminum Sport
- 2011142 4" x 1.5" Aluminum Softroll
- 2011152 5" x 1.5" Aluminum Softroll

Caster Fork & Fork Suspension

- 2011410 Standard Aluminum Fork
- 2011420 Frog Legs Uni-tine Fork \$ 195
- 2011421 Frog Legs Ultra-Sport Suspension Forks \$ 450
Compatible with 4" Caster only.

Rear Wheels

Spinergy Wheels

- 2012310 Spinergy Wire
- 2012311 Spinergy SPOX Everyday \$ 745
- 2012312 Spinergy SPOX Sport \$ 745
- 2012313 Spinergy LX Everyday \$ 895
- 2012314 Spinergy LX Sport \$ 895
- 2012315 **Spinergy Lite Extreme LX Flexrims** \$ 1,245
LX Wheels only. Not compatible with 26" Wheel size.

Select Wheel Camber

- 2012400 0° Camber
- 2012402 2° Camber
- 2012404 4° Camber
- 2012406 6° Camber

Axle Type

- 2013000 **Quick Release Rear Axles**



Scan this code or visit www.pdgmobility.com to learn more about the Elevation Ultralight Wheelchair.

Rear Tires

- 2014110 **Primo Pneumatic Tires**
- 2014120 **Schwalbe Marathon Plus Evolution** \$ 125

Hand Rims

- 2015000 **No Hand Rims**
- 2015100 **Aluminum Clear Anodized**
- 2015110 **Plastic Coated** \$ 135
- 2015120 **Natural Fit with Standard Grip** \$ 425

Wheel Locks

- 2016010 **Push to Lock**
- 2016090 **Composite Scissor Locks**

Anti-Tippers

- 2017010 **Rear Anti-Tippers** \$ 145

Back Options

- 2019221 **Rigid Back**
- 2019220 **Fold Down Back** \$ 150

Seat Back Height

- 2019300 Short (8.5"-12") - Rigidizer Bar Height = 5.5"
- 2019301 Medium (10.5"-16") - Rigidizer Bar Height = 7.5"
- 2019302 Tall (16.5"-20") - Rigidizer Bar Height = 9.5"
- 2019310 **Custom Rigidizer Bar Height:** _____ \$ 150
Must be at least 4.5" above Seat Sling
- 2019311 **Deep Rigidizer Bar** \$ 150
For installing third-party Backrest

Back Upholstery

- 2020010 **Standard Back Upholstery**
- 2020000 **No Back Upholstery**
If omitted, third-party Back Upholstery must be installed by dealer.

Footrest

- 2025010 **Standard Footrest**
- 2025011 **High Mount Footrest**

Foot Plates

- 2026020 **Angle Adjustable**
- 2026060 **Flip-up**

Calf Straps

- 2027050 **Bodypoint - Padded Calf Strap** \$ 105

Positioning Straps

- 2029005 **Standard Lap Belt**

Bodypoint - Hip Belts

Push button with flat-mount end fitting. Recommended for mounting through slots in side guards, to allow "anti-shear" feature to be used.

- 2029010 Center-Pull Auto Style \$ 105
- 2029015 Center-Pull Auto Style, Padded \$ 155
- 2029200 **Bodypoint - EvoFlex Pelvic Stabilizer** \$ 275

Push button recommended for when anti-shear feature is not used. Mounted directly to the seat tubes using the included seat tube mounting bracket.

Side Guards

- 2030400 **No Side Guards**

Non-Lifting Side Guards

- 2030410 Rigid - Small \$ 295
- 2030411 Rigid - Regular \$ 295
- 2030415 Folding - Regular \$ 350

Lifting (Anti-Shear) Side Guards

- 2030420 Rigid - Regular \$ 600
- 2030425 Folding - Regular \$ 650
- 2030426 Folding - Large \$ 650

Miscellaneous Options

- 2030450 **Alber Wheel Mount Adapter** \$ 600

Quad Triggers

- 2030461 Inboard Mounting \$ 115
- 2030462 Outboard Mounting \$ 115
- 2030463 Thumb Loops \$ 130