

| | |
|----------------------------------------|-----|
| Purchase Order | Tag |
| Dealer/Organization Name | |
| Purchaser Name and Contact Information | |
| Mailing Address | |
| Shipping Address | |

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type

Elevation Ultralight Manual Wheelchair

- | | | | |
|-------------------------------------|---------|------------------------------------------------------------------------------|----------|
| <input checked="" type="checkbox"/> | 2001000 | Ultralight Manual Rigid Wheelchair <i>K0005 / E1235</i> | \$ 3,900 |
| <input checked="" type="checkbox"/> | 2001001 | Patient Controlled Anterior/Posterior Back Angle Recline System <i>E1225</i> | \$ 850 |
| <input checked="" type="checkbox"/> | 2001002 | Patient Controlled Anterior Tilt & Elevation System | \$ 1,800 |

Frame Options

- | | | | |
|-------------------------------------------------------------------|---------|-----------------------------------------------------------------|--|
| <input checked="" type="checkbox"/> | 2002000 | Standard Base and Seat Frame 250 lbs Weight Capacity. | |
| User Weight Range Required for selection of Gas Struts. | | | |
| <input type="checkbox"/> | 2002150 | Up to 90 lbs (150 Newton Gas Struts) | |
| <input type="checkbox"/> | 2002200 | 90 lbs to 140 lbs (200 Newton Gas Struts) | |
| <input type="checkbox"/> | 2002300 | 140 lbs to 200 lbs (300 Newton Gas Struts) | |
| <input type="checkbox"/> | 2002400 | 200 lbs to 250 lbs (400 Newton Gas Struts) | |

Seat Width

Standard Seat Widths

- | | | |
|--------------------------|---------|----------------|
| <input type="checkbox"/> | 2003012 | 12" Seat Width |
| <input type="checkbox"/> | 2003013 | 13" Seat Width |
| <input type="checkbox"/> | 2003014 | 14" Seat Width |
| <input type="checkbox"/> | 2003015 | 15" Seat Width |
| <input type="checkbox"/> | 2003016 | 16" Seat Width |
| <input type="checkbox"/> | 2003017 | 17" Seat Width |
| <input type="checkbox"/> | 2003018 | 18" Seat Width |
| <input type="checkbox"/> | 2003019 | 19" Seat Width |

Seat Depth

Standard Seat Depths

- | | | | |
|--------------------------|---------|----------------|--------------|
| <input type="checkbox"/> | 2004014 | 14" Seat Depth | |
| <input type="checkbox"/> | 2004015 | 15" Seat Depth | |
| <input type="checkbox"/> | 2004016 | 16" Seat Depth | |
| <input type="checkbox"/> | 2004017 | 17" Seat Depth | |
| <input type="checkbox"/> | 2004018 | 18" Seat Depth | |
| <input type="checkbox"/> | 2004019 | 19" Seat Depth | |
| <input type="checkbox"/> | 2004020 | 20" Seat Depth | <i>E2201</i> |

Frame Color

PDG Primary Colors

- | | | |
|--------------------------|---------|--------------------|
| <input type="checkbox"/> | 2005000 | Tracer Red |
| <input type="checkbox"/> | 2005001 | Gloss Black |
| <input type="checkbox"/> | 2005002 | Silver |
| <input type="checkbox"/> | 2005003 | Black Shimmer |
| <input type="checkbox"/> | 2005004 | Orange Fusion |
| <input type="checkbox"/> | 2005005 | Karl's Imagination |
| <input type="checkbox"/> | 2005006 | Blue Shimmer |
| <input type="checkbox"/> | 2005007 | Gloss White |
| <input type="checkbox"/> | 2005008 | Dark Titanium |

PDG Premium Colors

- | | | | |
|--------------------------|---------|---------------|-----|
| <input type="checkbox"/> | 2005009 | Geyser Green | N/C |
| <input type="checkbox"/> | 2005010 | Lollipop Blue | |
| <input type="checkbox"/> | 2005011 | Lollipop Red | |

- | | | | |
|---------------------------------|---------|---------------------|--------|
| <input type="checkbox"/> | 2005200 | Custom Color | \$ 300 |
| Contact PDG for specifications. | | | |

Arm Type

- | | | |
|--------------------------|---------|----------------------|
| <input type="checkbox"/> | 2006000 | No Armrest |
| <input type="checkbox"/> | 2006001 | Rigid Armrest |

Seat Option

- | | | |
|-----------------------------------------------------------------|---------|----------------------------|
| <input type="checkbox"/> | 2008110 | Standard Seat Sling |
| <input type="checkbox"/> | 2008100 | No Seat Sling |
| If omitted, third-party Seat Sling must be installed by dealer. | | |

Cushion

- | | | |
|-------------------------------------|---------|-------------------|
| <input checked="" type="checkbox"/> | 2009000 | No Cushion |
|-------------------------------------|---------|-------------------|

Seat Height

24" Wheel Size

- 2010018 18" STF
- 2010019 19" STF
- 2010020 20" STF

25" Wheel Size

Not compatible with 4" Caster with Frog Leg Uni-tine Fork

- 2010118 18.5" STF
- 2010119 19.5" STF
- 2010120 20.5" STF

26" Wheel Size

Not compatible with Frog Leg Uni-tine Fork

- 2010219 19" STF
- 2010220 20" STF
- 2010221 21" STF

Casters

Composite

- 2011040 4" x 1.0" Composite Sport
- 2011042 4" x 1.5" Composite Softroll
- 2011052 5" x 1.5" Composite Softroll

Aluminium

\$ 125

- 2011140 4" x 1.0" Aluminum Sport
- 2011142 4" x 1.5" Aluminum Softroll
- 2011152 5" x 1.5" Aluminum Softroll

Caster Fork & Fork Suspension

- 2011410 Standard Aluminum Fork
- 2011420 Frog Legs Uni-tine Fork \$ 195
- 2011421 Frog Legs Ultra-Sport Suspension Forks *E1015* \$ 450
Compatible with 4" Caster only.

Rear Wheels

Spinergy Wheels

- 2012310 Spinergy Wire
- 2012311 Spinergy SPOX Everyday \$ 745
- 2012312 Spinergy SPOX Sport \$ 745
- 2012313 Spinergy LX Everyday \$ 895
- 2012314 Spinergy LX Sport \$ 895
- 2012315 **Spinergy Lite Extreme LX Flexirims** \$ 1,245
LX Wheels only. Not compatible with 26" Wheel size.

Wheel Camber

Select Wheel Camber

- 2012400 0° Camber
- 2012402 2° Camber
- 2012404 4° Camber
- 2012406 6° Camber

Axle Type

- 2013000 **Quick Release Rear Axles**



Scan this code or visit www.pdgmobility.com to learn more about the Elevation Ultralight Wheelchair.

Rear Tires

- 2014110 **Primo Pneumatic Tires**
- 2014120 **Schwalbe Marathon Plus Evolution** \$ 125

Hand Rims

- 2015000 **No Hand Rims**
- 2015100 **Aluminum Clear Anodized**
- 2015110 **Plastic Coated** \$ 135
- 2015120 **Natural Fit with Standard Grip** \$ 425

Wheel Locks

- 2016010 **Push to Lock**
- 2016090 **Composite Scissor Locks**

Anti-Tippers

- 2017010 **Rear Anti-Tippers** *E0971* \$ 145

Back Options

- 2019221 **Rigid Back**
- 2019220 **Fold Down Back** \$ 150
- Seat Back Height**
- 2019300 Short (8.5"-12") - Rigidizer Bar Height = 5.5"
- 2019301 Medium (10.5"-16") - Rigidizer Bar Height = 7.5"
- 2019302 Tall (16.5"-20") - Rigidizer Bar Height = 9.5"
- 2019310 **Custom Rigidizer Bar Height:** _____ \$ 150
Must be at least 4.5" above Seat Sling
- 2019311 **Deep Rigidizer Bar** \$ 150
For installing third-party Backrest

Back Upholstery

- 2020010 **Standard Back Upholstery**
- 2020000 **No Back Upholstery**
If omitted, third-party Back Upholstery must be installed by dealer.

Footrest

- 2025010 **Standard Footrest**
- 2025011 **High Mount Footrest**

Foot Plates

- 2026020 **Angle Adjustable** *K0040*
- 2026060 **Flip-up**

Calf Straps

- 2027050 **Bodypoint - Padded Calf Strap** *K0038* \$ 105

Positioning Straps

- 2029005 **Standard Lap Belt**
- Bodypoint - Hip Belts** *E0978*
Push button with flat-mount end fitting. Recommended for mounting through slots in side guards to allow "anti-shear" feature to be used.
- 2029010 Center-Pull Auto Style \$ 105
- 2029015 Center-Pull Auto Style, Padded \$ 155
- 2029200 **Bodypoint - EvoFlex Pelvic Stabilizer** \$ 275
Push button recommended for when anti-shear feature is not used. Mounted directly to the seat tubes using the included seat tube mounting bracket.

Side Guards

- 2030400 **No Side Guards**
- Non-Lifting Side Guards**
- 2030410 Rigid - Small \$ 295
- 2030411 Rigid - Regular \$ 295
- 2030415 Folding - Regular \$ 350
- Lifting (Anti-Shear) Side Guards**
- 2030420 Rigid - Regular \$ 600
- 2030425 Folding - Regular \$ 650
- 2030426 Folding - Large \$ 650

Miscellaneous Options

- 2030450 **Alber Wheel Mount Adapter** \$ 600
- Quad Triggers**
- 2030461 Inboard Mounting \$ 115
- 2030462 Outboard Mounting \$ 115
- 2030463 Thumb Loops \$ 130