

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type

- 1501000 **Fuze T50 50° Manual Tilt Wheelchair** \$ 3,995

Frame Options

- 1502010 **Fixed Base and Seat Frame**
350 lbs Weight Capacity. Select Fixed Seat Width below.
- 1502025 **Adjustable Base and Seat Frame** \$ 100
250 lbs Weight Capacity. Select Adjustable Seat Width below.
- 1502100 **Frame for Plug In Front Rigging** \$ 300
- 1502300 **Transit Tie Down System (TTDS)**
TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.

Seat Width

Fixed Seat Widths

Order with 1502010 Fixed Base and Seat Frame.

- 1503013 13" Seat Width
- 1503014 14" Seat Width
- 1503015 15" Seat Width
- 1503016 16" Seat Width
- 1503017 17" Seat Width
- 1503018 18" Seat Width
- 1503019 19" Seat Width
- 1503020 20" Seat Width
- 1503021 21" Seat Width \$ 215
- 1503022 22" Seat Width \$ 215

Adjustable Seat Widths

Order with 1502025 Adjustable Base and Seat Frame.

- 1503113 13"-15" Range - Set at 13"
- 1503114 13"-15" Range - Set at 14"
- 1503115 13"-15" Range - Set at 15"
- 1503116 15"-20" Range - Set at 15"
- 1503117 15"-20" Range - Set at 16"
- 1503118 15"-20" Range - Set at 17"
- 1503119 15"-20" Range - Set at 18"
- 1503120 15"-20" Range - Set at 19"
- 1503121 15"-20" Range - Set at 20"
- 1603321 21"-22" Range - Set at 21" \$ 215
- 1603322 21"-22" Range - Set at 22" \$ 215

Seat Depth

Standard Seat Depths (Adjustable)

- 1504015 15" Seat Depth
- 1504016 16" Seat Depth
- 1504017 17" Seat Depth
- 1504018 18" Seat Depth
- 1504019 19" Seat Depth
- 1504020 20" Seat Depth

Custom Seat Depths

\$ 325

Requires 25" Push bar Back Straight, Angle/Height Adjustable Push Bar Handle. Not compatible with Reclining or Fold Down Back.

- 1504021 21" Seat Depth
- 1504022 22" Seat Depth

Frame Colour

PDG Primary Colours

- 1505000 Tracer Red
- 1505001 Gloss Black
- 1505002 Silver
- 1505003 Black Shimmer
- 1505004 Orange Fusion
- 1505005 Karl's Imagination
- 1505006 Blue Shimmer
- 1505007 Gloss White
- 1505008 Dark Titanium

PDG Premium Colours

\$ 150

- 1505009 Geyser Green
- 1505010 Lollipop Blue
- 1505011 Lollipop Red

- 1505200 **Custom Colour** \$ 300

Contact PDG for specifications.

PDG Product Design Group Inc.
318 East Kent Avenue South
Vancouver BC V5X 4N6 Canada
Phone 604-323-9220 | 888-858-4422

Arm Type	
<input type="checkbox"/> 1506000	No Armrest
<input type="checkbox"/> 1506100	Single Post Fixed Height (12")
<input type="checkbox"/> 1506114	Single Post Adjustable Height (10"-14") \$ 255
<input type="checkbox"/> 1506200	Dual Post Fixed Height (12")
	Dual Post Adjustable Height \$ 255
<input type="checkbox"/> 1506214	10"-14"
<input type="checkbox"/> 1506216	12"-16"
	Locking Flip-up \$ 255
	Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.
<input type="checkbox"/> 1506310	Set at 10"
<input type="checkbox"/> 1506311	Set at 11"
<input type="checkbox"/> 1506312	Set at 12"
<input type="checkbox"/> 1506313	Set at 13"
<input type="checkbox"/> 1506314	Set at 14"

Armrests	
<input type="checkbox"/> 1507095	Curved Armpads (pair) Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.
	Regular Armpads
<input type="checkbox"/> 1507100	Desk Length - 2" x 10" (pair)
<input type="checkbox"/> 1507101	Full Length - 2" x 14" (pair)
	Wide Armpads \$ 115
	May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1507110	Desk Length - 3" x 10" (pair)
<input type="checkbox"/> 1507111	Full Length - 3" x 12" (pair)
	Gel Armpads \$ 115
	May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1507120	Full Length - 2" x 12" (pair)
<input type="checkbox"/> 1507121	Full Length - 3.5" x 12" (pair)
	Contour Arm Trough with Hand Pad \$ 125
	Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configurations.
<input type="checkbox"/> 1507130	Right (4.5" x 18")
<input type="checkbox"/> 1507131	Left (4.5" x 18")

Seat Pan	
<input type="checkbox"/> 1508000	No Seat Pan
<input type="checkbox"/> 1508001	Solid Seat Pan \$ 120

Cushion	
<input type="checkbox"/> 1509000	No Cushion
	nxt BioFit
<input type="checkbox"/> 1509101	12"-20" width, 14"-20" depth \$ 595
<input type="checkbox"/> 1509102	22"-24" width, 14"-20" depth \$ 625
	Cushion Size
	Available in even number sizes only.
1509400	Specify Size (Width" x Depth") _____



Scan this code or visit www.pdgmobility.com to learn more about the Fuze T50 Manual Tilt Wheelchair.

Seat Height																																																																			
Circle desired Seat-to-Floor Height or enter information below:																																																																			
Standard T50 Options	Wide T50 Options																																																																		
STF Height is approximate and can vary by as much as 1" depending on rear wheel type and pre-tilt setting (Standard T50 Chairs are shipped in 3° pre-tilt). To resolve, STF Height can be adjusted down by removing the pre-tilt collar/bumper.	Increases the overall width of the Chair by 2".																																																																		
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">5" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>13"</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>16"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td></td> <td></td> </tr> <tr> <td>20"</td> <td></td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> </tr> <tr> <td>22"</td> <td></td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td>15"</td> <td>16"</td> </tr> </tbody> </table>	RW	5" Casters					12"	13"					16"	13"	14"	15"			20"		13"	14"	15"	16"	22"			14"	15"	16"					15"	16"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">6" Casters & 6" Forks</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td>14"</td> <td>15"</td> <td></td> <td></td> </tr> <tr> <td>20"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> </tr> <tr> <td>22"</td> <td></td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td></td> <td>15"</td> <td>16"</td> </tr> </tbody> </table>	RW	6" Casters & 6" Forks					16"		14"	15"			20"		14"	15"	16"	17"	22"			14"	15"	16"	24"				15"	16"
RW	5" Casters																																																																		
12"	13"																																																																		
16"	13"	14"	15"																																																																
20"		13"	14"	15"	16"																																																														
22"			14"	15"	16"																																																														
				15"	16"																																																														
RW	6" Casters & 6" Forks																																																																		
16"		14"	15"																																																																
20"		14"	15"	16"	17"																																																														
22"			14"	15"	16"																																																														
24"				15"	16"																																																														
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">6" Casters & 8" Forks</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td></td> <td>14"</td> <td></td> <td></td> <td>18"</td> </tr> <tr> <td>16"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> </tr> <tr> <td>20"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>22"</td> <td></td> <td></td> <td></td> <td>17"</td> <td>18"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td></td> <td></td> <td>18"</td> </tr> </tbody> </table>	RW	6" Casters & 8" Forks					12"		14"			18"	16"		14"	15"	16"		20"			16"	17"	18"	22"				17"	18"	24"					18"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">6" Casters & 8" Forks</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td></td> <td></td> </tr> <tr> <td>20"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> </tr> <tr> <td>22"</td> <td></td> <td></td> <td></td> <td>15"</td> <td>16"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td></td> <td></td> <td>15"</td> </tr> </tbody> </table>	RW	6" Casters & 8" Forks					16"			15"			20"			15"	16"	17"	22"				15"	16"	24"					15"
RW	6" Casters & 8" Forks																																																																		
12"		14"			18"																																																														
16"		14"	15"	16"																																																															
20"			16"	17"	18"																																																														
22"				17"	18"																																																														
24"					18"																																																														
RW	6" Casters & 8" Forks																																																																		
16"			15"																																																																
20"			15"	16"	17"																																																														
22"				15"	16"																																																														
24"					15"																																																														
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">8" Casters & 8" Forks</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td></td> <td></td> <td></td> <td>18"</td> <td>19"</td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td></td> </tr> <tr> <td>20"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>22"</td> <td></td> <td></td> <td></td> <td>17"</td> <td>18"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td></td> <td></td> <td>18"</td> </tr> </tbody> </table>	RW	8" Casters & 8" Forks					12"				18"	19"	16"			15"	16"		20"			16"	17"	18"	22"				17"	18"	24"					18"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">8" Casters & 8" Forks</th> </tr> </thead> <tbody> <tr> <td>20"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td></td> <td></td> <td>16"</td> <td>17"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td></td> <td></td> <td>16"</td> </tr> </tbody> </table>	RW	8" Casters & 8" Forks					20"			16"	17"		22"				16"	17"	24"					16"						
RW	8" Casters & 8" Forks																																																																		
12"				18"	19"																																																														
16"			15"	16"																																																															
20"			16"	17"	18"																																																														
22"				17"	18"																																																														
24"					18"																																																														
RW	8" Casters & 8" Forks																																																																		
20"			16"	17"																																																															
22"				16"	17"																																																														
24"					16"																																																														
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="5">Ultra Sport Forks</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>20"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>22"</td> <td></td> <td></td> <td></td> <td>17"</td> <td>18"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td></td> <td></td> <td>18"</td> </tr> </tbody> </table>	RW	Ultra Sport Forks					16"			16"			20"			16"	17"	18"	22"				17"	18"	24"					18"																																					
RW	Ultra Sport Forks																																																																		
16"			16"																																																																
20"			16"	17"	18"																																																														
22"				17"	18"																																																														
24"					18"																																																														
Wheel size _____ Caster size _____ STF _____																																																																			

Casters		
Urethane		
<input type="checkbox"/>	1511051	5" x 1.25" Caster size
<input type="checkbox"/>	1511061	6" x 1.25" Caster size
<input type="checkbox"/>	1511063	6" x 1.50" Caster size \$ 85
<input type="checkbox"/>	1511081	8" x 1.25" Caster size
<input type="checkbox"/>	1511082	8" x 2.00" Caster size \$ 135
Pneumatic		
<input type="checkbox"/>	1511161	6" x 1.25" Caster size \$ 110
<input type="checkbox"/>	1511181	8" x 1.25" Caster size \$ 120
<input type="checkbox"/>	1511182	8" x 2.00" Caster size \$ 150
<input type="checkbox"/>	1511260	6" Soft Roll \$ 95
<input type="checkbox"/>	1511400	Ultra Sport Fork \$ 395 Not compatible with 12" Rear Wheels. Requires 6" Soft Roll Caster.

Rear Wheels		
Composite		
<input type="checkbox"/>	1512012	12" Wheel size
<input type="checkbox"/>	1512016	16" Wheel size
<input type="checkbox"/>	1512020	20" Wheel size
<input type="checkbox"/>	1512022	22" Wheel size
<input type="checkbox"/>	1512024	24" Wheel size
		Spoke \$ 75
<input type="checkbox"/>	1512120	20" Wheel size
<input type="checkbox"/>	1512122	22" Wheel size
<input type="checkbox"/>	1512124	24" Wheel size
<input type="checkbox"/>	1512224	24" Heavy Duty Spoke \$ 350

Axle Type	
<input type="checkbox"/>	1513000 Quick Release Rear Axles
<input type="checkbox"/>	1513010 Fixed Rear Axles

Rear Tires		
<input type="checkbox"/>	1514000	Urethane
<input type="checkbox"/>	1514010	Pneumatic \$ 60 Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.
<input type="checkbox"/>	1514020	Pneumatic with Flat-Free Inserts \$ 95 Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.

Hand Rims		
<input type="checkbox"/>	1515000	No Hand Rims
<input type="checkbox"/>	1515010	Standard 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/>	1515020	Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks		
<input type="checkbox"/>	1516005	Foot Operated Wheel Lock \$ 200 Not compatible with One Arm Drive
<input type="checkbox"/>	1516010	Push to Lock
<input type="checkbox"/>	1516020	Pull to Lock
<input type="checkbox"/>	1516030	Attendant Wheel Lock (Push to Lock) \$ 85 Not compatible with Wheel Lock Extensions.
<input type="checkbox"/>	1516040	User & Attendant Wheel Lock \$ 125
<input type="checkbox"/>	1516100	Wheel Lock Extensions \$ 60 Not compatible with Attendant Wheel Lock.

Anti-Tippers		
<input checked="" type="checkbox"/>	1517010	Rear Anti-Tippers \$ 60 Required Safety Option.

Back Post Style		
Push Bar Back		
<input type="checkbox"/>	1518010	Straight
<input type="checkbox"/>	1518030	8° Bend \$ 75
Cane Handle Back Not compatible with Headrest Mounting Plate.		
<input type="checkbox"/>	1518020	Straight
<input type="checkbox"/>	1518040	8° Bend \$ 75
Back Post Height Seat Pan to Handle		
<input type="checkbox"/>	1518420	20" Height
<input type="checkbox"/>	1518425	25" Height

Back Options		
Angle/Height Adjustable Push Bar Handle		\$ 150
Not compatible with Headrest Mounting Plate or Cane Handle Back.		
<input type="checkbox"/>	1519010	Button on Outside
<input type="checkbox"/>	1519011	Button on Inside
<input type="checkbox"/>	1519050	Removable Cane Handle Push Bar \$ 140 Not compatible with Push Bar Back.
Seat to Back Angle		
<input type="checkbox"/>	1519090	Set at 90°
<input type="checkbox"/>	1519096	Set at 96°
<input type="checkbox"/>	1519102	Set at 102°
<input type="checkbox"/>	1519108	Set at 108°
<input type="checkbox"/>	1519114	Set at 114°
<input type="checkbox"/>	1519120	Set at 120°
<input type="checkbox"/>	1519220	Fold Down Back
<input type="checkbox"/>	1519221	Non Fold Down Back \$ 75

Back Upholstery	
<input type="checkbox"/> 1520000	No Back Upholstery
	<p>Standard \$ 125</p> <p>Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.</p>
<input type="checkbox"/> 1520018	Set at 18"
<input type="checkbox"/> 1520019	Set at 19"
<input type="checkbox"/> 1520020	Set at 20"
<input type="checkbox"/> 1520021	Set at 21"
<input type="checkbox"/> 1520022	Set at 22"
<input type="checkbox"/> 1520023	Set at 23"
	<p>Adjustable Tension \$ 250</p> <p>Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.</p>
<input type="checkbox"/> 1520321	Set at 21"
<input type="checkbox"/> 1520322	Set at 22"
<input type="checkbox"/> 1520323	Set at 23"

Backrest	
<input type="checkbox"/> 1521000	No Backrest
	<p>nxt Optima HA Active Contour</p>
<input type="checkbox"/> 1521100	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/> 1521105	22"-24" width, 16"-20" height \$ 720
	<p>nxt Optima HA Plus Contour</p>
<input type="checkbox"/> 1521101	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/> 1521106	22"-24" width, 16"-20" height \$ 720
	<p>Backrest Size</p> <p>Available in even number sizes only.</p> <p>1521400 Specify Size (Width" x Height") _____</p>

Tilt/Recline Options	
	<p>Attendant Tilt only</p>
<input type="checkbox"/> 1523011	Installed on Right Back Post
<input type="checkbox"/> 1523012	Installed on Left Back Post
	<p>Tilt Limit</p>
<input type="checkbox"/> 1523220	Set at 20°
<input type="checkbox"/> 1523225	Set at 25°
<input type="checkbox"/> 1523230	Set at 30°
<input type="checkbox"/> 1523235	Set at 35°
<input type="checkbox"/> 1523240	Set at 40°
<input type="checkbox"/> 1523245	Set at 45°
<input type="checkbox"/> 1523250	Set at 50°
<input type="checkbox"/> 1523410	Reclining Backrest (Manual) \$ 830
	<p>Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back or Tension Adjustable Upholstery. Requires 25" Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15", 16", 21", or 22" Seat Depth.</p>
	<p>Power Tilt \$ 1,595</p> <p>Tilt range 4° to 45°, Max weight 250 lbs.</p>
<input type="checkbox"/> 1523421	Toggle Switch - Installed on Right side
<input type="checkbox"/> 1523422	Toggle Switch - Installed on Left side
<input type="checkbox"/> 1523431	Button Switch - Installed on Right side
<input type="checkbox"/> 1523432	Button Switch - Installed on Left side

Headrest	
<input type="checkbox"/> 1524000	No Headrest
<input type="checkbox"/> 1524010	Headrest Mounting Plate \$ 75
	<p>Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.</p>
<input type="checkbox"/> 1524100	nxt MultiFit Head Support Hardware \$ 155
	<p>nxt MultiFit Head Support Pad \$ 195</p>
<input type="checkbox"/> 1524105	Regular (9" x 5")
<input type="checkbox"/> 1524106	Large (10" x 5")
<input type="checkbox"/> 1524107	X-Large (14" x 5")

Front Rigging	
<input type="checkbox"/> 1525000	No Front Rigging
<input type="checkbox"/> 1525001	Sector Blocks only
<input type="checkbox"/> 1525110	Plug in 70°
	<p>Frame upgrade required. Select 1502100 above.</p>
<input type="checkbox"/> 1525120	Plug in Adjustable Knee Angle \$ 975
	<p>Frame upgrade required. Select 1502100 above.</p>
<input type="checkbox"/> 1525140	Swing Away 70°
<input type="checkbox"/> 1525160	Swing Away 90°
	<p>Not compatible with 8" Casters.</p>
<input type="checkbox"/> 1525200	Swing Away Elevating Legrest \$ 355

Foot Plates	
<input type="checkbox"/> 1526010	Composite
<input type="checkbox"/> 1526020	Angle Adjustable \$ 140
<input type="checkbox"/> 1526040	One Piece Flip-up \$ 395
<input type="checkbox"/> 1526050	Interlocking Angle Adjustable \$ 195
	<p>Seat Widths of 14"-21" only.</p>

Footrest/Legrest Accessories	
<input type="checkbox"/> 1527015	Swing Away 70° Calf Pad \$ 135
	<p>Not compatible with 70° Gel Pad.</p>
<input type="checkbox"/> 1527020	Swing Away 70° Gel Pad \$ 100
	<p>Not compatible with 70° Swing Away Calf Pad.</p>
<input type="checkbox"/> 1527030	Swing Away ELR Gel Knee Pad \$ 100
<input type="checkbox"/> 1527040	Adjustable Length Heel Loops \$ 50
<input type="checkbox"/> 1527050	Bodypoint - Padded Calf Strap \$ 105
	<p>Miller's Dynamic Footrest Mounts</p>
<input type="checkbox"/> 1527060	For Plug In Front Rigging \$ 190
<input type="checkbox"/> 1527065	For Swing Away 70° Front Rigging \$ 380
	<p>Not available on 13" STF</p>

Extension Tubes	
Regular Extension Tubes	
Not compatible with Swing Away 90°.	
<input type="checkbox"/> 1528010	Standard (13"-17")
<input type="checkbox"/> 1528020	Long (15"-19")
Extra Strength Extension Tubes (Solid)	
Available on Plug In Front Rigging only.	
<input type="checkbox"/> 1528110	Standard (13"-17")
<input type="checkbox"/> 1528120	Long (15"-19")
<input type="checkbox"/> 1528200	90° Extension Tubes (11"-16")
Not compatible with 70° or ELR.	

Positioning Straps	
Bodypoint - Hip Belts	
<input type="checkbox"/> 1529010	Center-Pull Auto Style
<input type="checkbox"/> 1529015	Center-Pull Auto Style, Padded
<input type="checkbox"/> 1529020	4-Point, Padded
<input type="checkbox"/> 1529030	Bodypoint - Chest Strap

Miscellaneous Options	
One-Arm Drive	
Requires 24" Composite Rear Wheel. Forward Wheel Position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks.	
<input type="checkbox"/> 1530010	Installed on Right side
<input type="checkbox"/> 1530011	Installed on Left side
<input type="checkbox"/> 1530060	Stealth Flip Down Abductor
<input type="checkbox"/> 1530070	QR Bracket - Accessory Mount
Not compatible with TTDS.	
<input type="checkbox"/> 1530080	Gel Hip Guides
<input type="checkbox"/> 1530110	IV Pole
<input type="checkbox"/> 1530120	O2 Holder
May add to overall width of the Chair.	
<input type="checkbox"/> 1530150	Ventilator Tray
Minimum Seat Width is 16". Reduces ground clearance. May become unstable when no occupant in the chair. Not compatible with O2 holder.	
<input type="checkbox"/> 1530152	Ventilator Bracket
Compatible with LTV1 (lap tops) ventilator.	
<input type="checkbox"/> 1530154	Ventilator Battery Tray
Minimum Seat Width is 16". Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with seat height restrictions.	
<input type="checkbox"/> 1530160	Clear Lap Tray (Slide On)
Not compatible with 8"- 11" Dual Post Arms, Desk Length Pads, or Gel Arm Pads. Seat Widths of 15"-20" only.	
<input type="checkbox"/> 1530170	Adhesive Back Hook Strip, 2" x 30"
For Cushion mounting.	
<input type="checkbox"/> 1530300	Fuze T50 to T20 Conversion Kit

Total MSRP:

Price will update when form is completed electronically.