

| | |
|----------------------------------------|-----|
| Purchase Order | Tag |
| Dealer/Organization Name | |
| Purchaser Name and Contact Information | |
| Mailing Address | |
| Shipping Address | |

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type

- 1501000 **Fuze T50 50° Manual Tilt Wheelchair** E1161 \$ 2,995

Frame Options

- 1502010 **Fixed Base and Seat Frame**
350 lbs Weight Capacity. Select Fixed Seat Width below.
- 1502025 **Adjustable Base and Seat Frame** \$ 100
250 lbs Weight Capacity. Select Adjustable Seat Width below.
- 1502100 **Frame for Plug In Front Rigging** \$ 300
- 1502300 **Transit Tie Down System (TTDS)**
TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.

Seat Width

Fixed Seat Widths

Order with 1502010 Fixed Base and Seat Frame.

- 1503013 13" Seat Width
- 1503014 14" Seat Width
- 1503015 15" Seat Width
- 1503016 16" Seat Width
- 1503017 17" Seat Width
- 1503018 18" Seat Width
- 1503019 19" Seat Width
- 1503020 20" Seat Width E2201 \$ 325
- 1503021 21" Seat Width E2201 \$ 325
- 1503022 22" Seat Width E2201 \$ 325

Adjustable Seat Widths

Order with 1502025 Adjustable Base and Seat Frame.

- 1503113 13"-15" Range - Set at 13"
- 1503114 13"-15" Range - Set at 14"
- 1503115 13"-15" Range - Set at 15"
- 1503116 15"-20" Range - Set at 15"
- 1503117 15"-20" Range - Set at 16"
- 1503118 15"-20" Range - Set at 17"
- 1503119 15"-20" Range - Set at 18"
- 1503120 15"-20" Range - Set at 19"
- 1503220 15"-20" Range - Set at 20" E2201 \$ 325
- 1503321 21"-22" Range - Set at 21" E2201 \$ 325
- 1503322 21"-22" Range - Set at 22" E2201 \$ 325

Seat Depth

Standard Seat Depths (Adjustable)

- 1504015 15" Seat Depth
- 1504016 16" Seat Depth
- 1504017 17" Seat Depth
- 1504018 18" Seat Depth
- 1504019 19" Seat Depth
- 1504020 20" Seat Depth E2203 \$ 475

Custom Seat Depths

Requires 25" Push bar Back Straight, Angle/Height Adjustable Push Bar Handle. Not compatible with Reclining or Fold Down Back.

- 1504021 21" Seat Depth E2203 \$ 475
- 1504022 22" Seat Depth E2204 \$ 750

Frame Color

PDG Primary Colors

- 1505000 Tracer Red
- 1505001 Gloss Black
- 1505002 Silver
- 1505003 Black Shimmer
- 1505004 Orange Fusion
- 1505005 Karl's Imagination
- 1505006 Blue Shimmer
- 1505007 Gloss White
- 1505008 Dark Titanium

PDG Premium Colors

\$ 150

- 1505009 Geyser Green
- 1505010 Lollipop Blue
- 1505011 Lollipop Red

- 1505200 **Custom Color** \$ 300

Contact PDG for specifications.

PDG Product Design Group Inc.
318 East Kent Avenue South
Vancouver BC V5X 4N6 Canada
Phone 604-323-9220 | 888-858-4422

| Arm Type | | | |
|------------------------------------------------------------------------------------------------------------------------------------|---------|------------------------------------------------|--------------|
| <input type="checkbox"/> | 1506000 | No Armrest | |
| <input type="checkbox"/> | 1506100 | Single Post Fixed Height (12") | |
| <input type="checkbox"/> | 1506114 | Single Post Adjustable Height (10"-14") | E0973 \$ 255 |
| <input type="checkbox"/> | 1506200 | Dual Post Fixed Height (12") | |
| | | Dual Post Adjustable Height | E0973 \$ 255 |
| <input type="checkbox"/> | 1506214 | 10"-14" | |
| <input type="checkbox"/> | 1506216 | 12"-16" | |
| | | Locking Flip-up | E0973 \$ 255 |
| Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. | | | |
| <input type="checkbox"/> | 1506310 | Set at 10" | |
| <input type="checkbox"/> | 1506311 | Set at 11" | |
| <input type="checkbox"/> | 1506312 | Set at 12" | |
| <input type="checkbox"/> | 1506313 | Set at 13" | |
| <input type="checkbox"/> | 1506314 | Set at 14" | |

| Armrests | | | |
|--------------------------------------------------------------------------------------------------------------------------------------|---------|----------------------------------|--------------|
| <input type="checkbox"/> | 1507095 | Curved Armpads (pair) | |
| Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. | | | |
| | | Regular Armpads | |
| <input type="checkbox"/> | 1507100 | Desk Length - 2" x 10" (pair) | |
| <input type="checkbox"/> | 1507101 | Full Length - 2" x 14" (pair) | |
| | | Wide Armpads | \$ 115 |
| May interfere with wheels in some tilt configurations. | | | |
| <input type="checkbox"/> | 1507110 | Desk Length - 3" x 10" (pair) | |
| <input type="checkbox"/> | 1507111 | Full Length - 3" x 12" (pair) | |
| | | Gel Armpads | \$ 115 |
| May interfere with wheels in some tilt configurations. | | | |
| <input type="checkbox"/> | 1507120 | Full Length - 2" x 12" (pair) | |
| <input type="checkbox"/> | 1507121 | Full Length - 3.5" x 12" (pair) | |
| | | Contour Arm Trough | E2209 \$ 125 |
| Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configurations. | | | |
| <input type="checkbox"/> | 1507130 | Right (4.5" x 18") with Hand Pad | |
| <input type="checkbox"/> | 1507131 | Left (4.5" x 18") with Hand Pad | |

| Seat Pan | | | |
|--------------------------|---------|-----------------------|--------------|
| <input type="checkbox"/> | 1508000 | No Seat Pan | |
| <input type="checkbox"/> | 1508001 | Solid Seat Pan | E2231 \$ 120 |

| Cushion | | | |
|--------------------------|---------|-----------------------------------|----------------|
| <input type="checkbox"/> | 1509000 | No Cushion | |
| | | nxt BioFit | |
| <input type="checkbox"/> | 1509101 | 14"-20" width, 16"-20" depth | E2607 \$ 550 |
| <input type="checkbox"/> | 1509102 | 22"-24" width, 16"-20" depth | E2608 \$ 600 |
| | | Matrx VI | |
| <input type="checkbox"/> | 1509111 | 14"-20" width, 16"-22" depth | E2607 \$ 550 |
| <input type="checkbox"/> | 1509112 | 22" width, 18"-22" depth | E2608 \$ 600 |
| <input type="checkbox"/> | 1509115 | HD - 18"-30" width, 18"-22" depth | E2608 \$ 1,150 |



Scan this code or visit www.pdgmobility.com to learn more about the Fuze T50 Manual Tilt Wheelchair.

| | | | |
|-------------------------------|--------------------------------------|------------------------------|--------------|
| ROHO Quadtro Select HP | | | |
| <input type="checkbox"/> | 1509121 | 14"-20" width, 16"-22" depth | E2624 \$ 650 |
| <input type="checkbox"/> | 1509122 | 22"-24" width, 16"-22" depth | E2625 \$ 750 |
| ROHO High Profile | | | |
| <input type="checkbox"/> | 1509125 | 14"-20" width, 16"-22" depth | E2622 \$ 650 |
| <input type="checkbox"/> | 1509126 | 22"-24" width, 16"-22" depth | E2623 \$ 750 |
| Varilite Reflex | | | |
| <input type="checkbox"/> | 1509131 | 14"-20" width, 16"-20" depth | E2603 \$ 375 |
| Cushion Size | | | |
| 1509400 | Specify Size (Width" x Depth") _____ | | |

| Seat Height | | | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|-------------------------|-----|
| Circle desired Seat-to-Floor Height or enter information below: | | | |
| Standard T50 Options | | Wide T50 Options | |
| STF Height is approximate and can vary by as much as 1" depending on rear wheel type and pre-tilt setting (Standard T50 Chairs are shipped in 3" pre-tilt). To resolve, STF Height can be adjusted down by removing the pre-tilt collar/bumper. | | | |
| Increases the overall width of the Chair by 2". | | | |
| RW | 5" Casters | | |
| 12" | 13" | 14" | 15" |
| 16" | 13" | 14" | 15" |
| 20" | 13" | 14" | 15" |
| 22" | 14" | 15" | 16" |
| 24" | 14" | 15" | 16" |
| RW | 6" Casters & 6" Forks | | |
| 12" | 14" | 15" | 16" |
| 16" | 14" | 15" | 16" |
| 20" | 14" | 15" | 16" |
| 22" | 14" | 15" | 16" |
| 24" | 14" | 15" | 16" |
| RW | 6" Casters & 8" Forks | | |
| 12" | 14" | 15" | 16" |
| 16" | 14" | 15" | 16" |
| 20" | 14" | 15" | 16" |
| 22" | 14" | 15" | 16" |
| 24" | 14" | 15" | 16" |
| RW | 8" Casters & 8" Forks | | |
| 12" | 14" | 15" | 16" |
| 16" | 14" | 15" | 16" |
| 20" | 14" | 15" | 16" |
| 22" | 14" | 15" | 16" |
| 24" | 14" | 15" | 16" |
| RW | Ultra Sport Forks | | |
| 16" | 16" | 17" | 18" |
| 20" | 16" | 17" | 18" |
| 22" | 16" | 17" | 18" |
| 24" | 16" | 17" | 18" |
| Wheel size _____ Caster size _____ STF _____ | | | |

| Casters | | |
|--------------------------|---------|------------------------------------------------------------------------------------------------------------|
| Urethane | | |
| <input type="checkbox"/> | 1511051 | 5" x 1.25" Caster size |
| <input type="checkbox"/> | 1511061 | 6" x 1.25" Caster size |
| <input type="checkbox"/> | 1511063 | 6" x 1.50" Caster size \$ 85 |
| <input type="checkbox"/> | 1511081 | 8" x 1.25" Caster size |
| <input type="checkbox"/> | 1511082 | 8" x 2.00" Caster size \$ 135 |
| Pneumatic | | |
| <input type="checkbox"/> | 1511161 | 6" x 1.25" Caster size \$ 110 |
| <input type="checkbox"/> | 1511181 | 8" x 1.25" Caster size \$ 120 |
| <input type="checkbox"/> | 1511182 | 8" x 2.00" Caster size \$ 150 |
| <input type="checkbox"/> | 1511260 | 6" Soft Roll \$ 95 |
| <input type="checkbox"/> | 1511400 | Ultra Sport Fork E1015 \$ 395 Not compatible with 12" Rear Wheels. Requires 6" Soft Roll Caster. |

| Rear Wheels | | |
|--------------------------|---------|------------------------------------|
| Composite | | |
| <input type="checkbox"/> | 1512012 | 12" Wheel size |
| <input type="checkbox"/> | 1512016 | 16" Wheel size |
| <input type="checkbox"/> | 1512020 | 20" Wheel size |
| <input type="checkbox"/> | 1512022 | 22" Wheel size |
| <input type="checkbox"/> | 1512024 | 24" Wheel size |
| | | \$ 75 |
| Spoke | | |
| <input type="checkbox"/> | 1512120 | 20" Wheel size |
| <input type="checkbox"/> | 1512122 | 22" Wheel size |
| <input type="checkbox"/> | 1512124 | 24" Wheel size |
| <input type="checkbox"/> | 1512224 | 24" Heavy Duty Spoke \$ 350 |

| Axle Type | |
|--------------------------|-----------------------------------------|
| <input type="checkbox"/> | 1513000 Quick Release Rear Axles |
| <input type="checkbox"/> | 1513010 Fixed Rear Axles |

| Rear Tires | | |
|--------------------------|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> | 1514000 | Urethane |
| <input type="checkbox"/> | 1514010 | Pneumatic \$ 60 Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated. |
| <input type="checkbox"/> | 1514020 | Pneumatic with Flat-Free Inserts E2213 \$ 95 Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels. |

| Hand Rims | | |
|--------------------------|---------|-----------------------------------------------------------------|
| <input type="checkbox"/> | 1515000 | No Hand Rims |
| <input type="checkbox"/> | 1515010 | Standard 20", 22", and 24" Rear Wheels only. |
| <input type="checkbox"/> | 1515020 | Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 95 |

| Wheel Locks | | |
|--------------------------|---------|------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> | 1516005 | Foot Operated Wheel Lock \$ 200 Not compatible with One Arm Drive |
| <input type="checkbox"/> | 1516010 | Push to Lock |
| <input type="checkbox"/> | 1516020 | Pull to Lock |
| <input type="checkbox"/> | 1516030 | Attendant Wheel Lock (Push to Lock) \$ 85 Not compatible with Wheel Lock Extensions. |
| <input type="checkbox"/> | 1516040 | User & Attendant Wheel Lock \$ 125 |
| <input type="checkbox"/> | 1516100 | Wheel Lock Extensions E0961 \$ 60 Not compatible with Attendant Wheel Lock. |

| Anti-Tippers | | |
|-------------------------------------|---------|-----------------------------------------------------------------|
| <input checked="" type="checkbox"/> | 1517010 | Rear Anti-Tippers E0971 \$ 60 Required Safety Option. |

| Back Post Style | | |
|-------------------------------------------------------------------------|---------|---------------|
| Push Bar Back | | |
| <input type="checkbox"/> | 1518010 | Straight |
| <input type="checkbox"/> | 1518030 | 8° Bend \$ 75 |
| Cane Handle Back Not compatible with Headrest Mounting Plate. | | |
| <input type="checkbox"/> | 1518020 | Straight |
| <input type="checkbox"/> | 1518040 | 8° Bend \$ 75 |
| Back Post Height Seat Pan to Handle | | |
| <input type="checkbox"/> | 1518420 | 20" Height |
| <input type="checkbox"/> | 1518425 | 25" Height |

| Back Options | | |
|------------------------------------------------------------------|---------|------------------------------------------------------------------------------------|
| Angle/Height Adjustable Push Bar Handle | | \$ 150 |
| Not compatible with Headrest Mounting Plate or Cane Handle Back. | | |
| <input type="checkbox"/> | 1519010 | Button on Outside |
| <input type="checkbox"/> | 1519011 | Button on Inside |
| <input type="checkbox"/> | 1519050 | Removable Cane Handle Push Bar \$ 140 Not compatible with Push Bar Back. |
| Seat to Back Angle | | |
| <input type="checkbox"/> | 1519090 | Set at 90° |
| <input type="checkbox"/> | 1519096 | Set at 96° |
| <input type="checkbox"/> | 1519102 | Set at 102° |
| <input type="checkbox"/> | 1519108 | Set at 108° |
| <input type="checkbox"/> | 1519114 | Set at 114° |
| <input type="checkbox"/> | 1519120 | Set at 120° |
| <input type="checkbox"/> | 1519220 | Fold Down Back |
| <input type="checkbox"/> | 1519221 | Non Fold Down Back \$ 75 |

| Back Upholstery | |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> 1520000 | No Back Upholstery |
| | <p>Standard \$ 125</p> <p>Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.</p> |
| <input type="checkbox"/> 1520018 | Set at 18" |
| <input type="checkbox"/> 1520019 | Set at 19" |
| <input type="checkbox"/> 1520020 | Set at 20" |
| <input type="checkbox"/> 1520021 | Set at 21" |
| <input type="checkbox"/> 1520022 | Set at 22" |
| <input type="checkbox"/> 1520023 | Set at 23" |
| | <p>Adjustable Tension \$ 250</p> <p>Back height measured from Seat Pan to top of Upholstery. Not available with Cane Handle Back.</p> |
| <input type="checkbox"/> 1520321 | Set at 21" |
| <input type="checkbox"/> 1520322 | Set at 22" |
| <input type="checkbox"/> 1520323 | Set at 23" |

| Backrest | |
|----------------------------------|----------------------------------------------------------|
| <input type="checkbox"/> 1521000 | No Backrest |
| | <p>nxt Optima HA</p> |
| <input type="checkbox"/> 1521100 | AC - 16"-20" width, 16"-20" height <i>E2620</i> \$ 690 |
| <input type="checkbox"/> 1521101 | PC - 16"-20" width, 16"-20" height <i>E2620</i> \$ 690 |
| <input type="checkbox"/> 1521105 | AC - 22"-24" width, 16"-20" height <i>E2621</i> \$ 720 |
| <input type="checkbox"/> 1521106 | PC - 22"-24" width, 16"-20" height <i>E2621</i> \$ 720 |
| | <p>Matrx Elite</p> |
| <input type="checkbox"/> 1521110 | 14"-20" width, 10"-20" height <i>E2620</i> \$ 825 |
| <input type="checkbox"/> 1521112 | HD - 22"-26" width, 16"-20" height <i>E2621</i> \$ 1,095 |
| <input type="checkbox"/> 1521116 | Deep - 14"-20" width, 12"-20" height <i>E2620</i> \$ 950 |
| | <p>ROHO Agility</p> |
| <input type="checkbox"/> 1521120 | Min - 14"-20" width, 10"-18" height <i>E2620</i> \$ 825 |
| <input type="checkbox"/> 1521122 | Mid - 14"-20" width, 10"-18" height <i>E2620</i> \$ 850 |
| <input type="checkbox"/> 1521124 | Max - 14"-20" width, 10"-18" height <i>E2620</i> \$ 950 |
| | <p>Varilite Icon</p> |
| <input type="checkbox"/> 1521130 | Tall - 14"-20" width, 14"-24" height <i>E2615</i> \$ 850 |
| <input type="checkbox"/> 1521132 | Deep - 14"-20" width, 17"-24" height <i>E2620</i> \$ 895 |
| | <p>Backrest Size</p> |
| 1521400 | Specify Size (Width" x Height") _____ |

| Tilt/Recline Options | |
|----------------------------------|-----------------------------------|
| | <p>Attendant Tilt only</p> |
| <input type="checkbox"/> 1523011 | Installed on Right Back Post |
| <input type="checkbox"/> 1523012 | Installed on Left Back Post |
| | <p>Tilt Limit</p> |
| <input type="checkbox"/> 1523220 | Set at 20° |
| <input type="checkbox"/> 1523225 | Set at 25° |
| <input type="checkbox"/> 1523230 | Set at 30° |
| <input type="checkbox"/> 1523235 | Set at 35° |
| <input type="checkbox"/> 1523240 | Set at 40° |
| <input type="checkbox"/> 1523245 | Set at 45° |
| <input type="checkbox"/> 1523250 | Set at 50° |

| | |
|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> 1523410 | Reclining Backrest (Manual) <i>E1225</i> \$ 830 |
| | Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back or Tension Adjustable Upholstery. Requires 25° Angle/Height Adjustable Push Bar with 8° Bend. Not available with 15", 16", 21", or 22" Seat Depth. |
| | <p>Power Tilt <i>E1002</i> \$ 2,095</p> <p>Tilt range 4° to 45°, Max weight 250 lbs.</p> |
| <input type="checkbox"/> 1523421 | Toggle Switch - Installed on Right side |
| <input type="checkbox"/> 1523422 | Toggle Switch - Installed on Left side |
| <input type="checkbox"/> 1523431 | Button Switch - Installed on Right side |
| <input type="checkbox"/> 1523432 | Button Switch - Installed on Left side |

| Headrest | |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> 1524000 | No Headrest |
| <input type="checkbox"/> 1524010 | Headrest Mounting Plate \$ 75 |
| | Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back. |
| <input type="checkbox"/> 1524100 | nxt MultiFit Head Support Hardware <i>E1028</i> \$ 155 |
| | <p>nxt MultiFit Head Support Pad <i>E0955</i> \$ 195</p> |
| <input type="checkbox"/> 1524105 | Regular (9" x 5") |
| <input type="checkbox"/> 1524106 | Large (10" x 5") |
| <input type="checkbox"/> 1524107 | X-Large (14" x 5") |
| | <p>Stealth Headrest Hardware <i>E1028</i></p> |
| <input type="checkbox"/> 1524110 | Multi-Axis HMO475 \$ 150 |
| <input type="checkbox"/> 1524111 | Link TWB480 \$ 325 |
| | <p>Stealth Headrest Pad <i>E0955</i> \$ 200</p> |
| <input type="checkbox"/> 1524115 | 10" width |
| <input type="checkbox"/> 1524116 | 14" width |

| Front Rigging | |
|----------------------------------|---------------------------------------------------------|
| <input type="checkbox"/> 1525000 | No Front Rigging |
| <input type="checkbox"/> 1525001 | Sector Blocks only |
| <input type="checkbox"/> 1525110 | Plug in 70° |
| | Frame upgrade required. Select 1502100 above. |
| <input type="checkbox"/> 1525120 | Plug in Adjustable Knee Angle \$ 975 |
| | Frame upgrade required. Select 1502100 above. |
| <input type="checkbox"/> 1525140 | Swing Away 70° |
| <input type="checkbox"/> 1525160 | Swing Away 90° |
| | Not compatible with 8" Casters. |
| <input type="checkbox"/> 1525200 | Swing Away Elevating Legrest <i>E0990</i> \$ 355 |

| Foot Plates | |
|----------------------------------|---------------------------------------------|
| <input type="checkbox"/> 1526010 | Composite |
| <input type="checkbox"/> 1526020 | Angle Adjustable <i>K0040</i> \$ 140 |
| <input type="checkbox"/> 1526040 | One Piece Flip-up \$ 395 |
| <input type="checkbox"/> 1526050 | Interlocking Angle Adjustable \$ 195 |
| | Seat Widths of 14"-21" only. |

| Footrest/Legrest Accessories | | |
|-----------------------------------------|----------------------------------------------------|--------|
| <input type="checkbox"/> | 1527015 Swing Away 70° Calf Pad | \$ 135 |
| | Not compatible with 70° Gel Pad. | |
| <input type="checkbox"/> | 1527020 Swing Away 70° Gel Pad | \$ 100 |
| | Not compatible with 70° Swing Away Calf Pad. | |
| <input type="checkbox"/> | 1527030 Swing Away ELR Gel Knee Pad | \$ 100 |
| <input type="checkbox"/> | 1527040 Adjustable Length Heel Loops E0951 | \$ 50 |
| <input type="checkbox"/> | 1527050 Bodypoint - Padded Calf Strap K0038 | \$ 105 |
| Miller's Dynamic Footrest Mounts | | |
| <input type="checkbox"/> | 1527060 For Plug In Front Rigging | \$ 190 |
| <input type="checkbox"/> | 1527065 For Swing Away 70° Front Rigging | \$ 380 |
| | Not available on 13" STF | |

| Extension Tubes | | |
|-----------------------------------------------|----------------------------------------------|-------|
| Regular Extension Tubes | | |
| | Not compatible with Swing Away 90°. | |
| <input type="checkbox"/> | 1528010 Standard (13"-17") | |
| <input type="checkbox"/> | 1528020 Long (15"-19") | \$ 75 |
| Extra Strength Extension Tubes (Solid) | | |
| | Available on Plug In Front Rigging only. | |
| <input type="checkbox"/> | 1528110 Standard (13"-17") | |
| <input type="checkbox"/> | 1528120 Long (15"-19") | |
| <input type="checkbox"/> | 1528200 90° Extension Tubes (11"-16") | \$ 75 |
| | Not compatible with 70° or ELR. | |

| Positioning Straps | | |
|-----------------------------------|----------------------------------------------|--------|
| Bodypoint - Hip Belt E0978 | | |
| <input type="checkbox"/> | 1529010 Center-Pull Auto Style | \$ 105 |
| <input type="checkbox"/> | 1529015 Center-Pull Auto Style, Padded | \$ 155 |
| <input type="checkbox"/> | 1529020 4-Point, Padded | \$ 175 |
| <input type="checkbox"/> | 1529030 Bodypoint - Chest Strap E0960 | \$ 85 |

| Miscellaneous Options | | |
|-------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|
| One-Arm Drive E0958 \$ 3,500 | | |
| | Requires 24" Composite Rear Wheel. Forward Wheel Position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks. | |
| <input type="checkbox"/> | 1530010 Installed on Right side | |
| <input type="checkbox"/> | 1530011 Installed on Left side | |
| <input type="checkbox"/> | 1530060 Stealth Flip Down Abductor E0957 | \$ 495 |
| <input type="checkbox"/> | 1530080 QR Bracket - Accessory Mount | \$ 200 |
| | Not compatible with TTDS. | |
| <input type="checkbox"/> | 1530080 Gel Hip Guides | \$ 325 |
| <input type="checkbox"/> | 1530110 IV Pole E1029 | \$ 295 |
| <input type="checkbox"/> | 1530120 O2 Holder E2208 | \$ 295 |
| | May add to overall width of the Chair. | |
| <input type="checkbox"/> | 1530150 Ventilator Tray E1029 | \$ 295 |
| | Minimum Seat Width is 16". Reduces ground clearance. May become unstable when no occupant in the chair. Not compatible with O2 holder. | |
| <input type="checkbox"/> | 1530152 Ventilator Bracket E1029 | \$ 295 |
| | Compatible with LTV1 (lap tops) ventilator. | |
| <input type="checkbox"/> | 1530154 Ventilator Battery Tray E1029 | \$ 95 |
| | Minimum Seat Width is 16". Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with seat height restrictions. | |
| <input type="checkbox"/> | 1530160 Clear Lap Tray (Slide On) | \$ 490 |
| | Not compatible with 8"- 11" Dual Post Arms, Desk Length Pads, or Gel Arm Pads. Seat Widths of 15"-20" only. | |
| <input type="checkbox"/> | 1530170 Adhesive Back Hook Strip, 2" x 30" | \$ 10 |
| | For Cushion mounting. | |
| <input type="checkbox"/> | 1530300 Fuze T50 to T20 Conversion Kit K0108 | \$ 995 |