

FUZE T50 JR

effortless

tilt



- Short wheelbase provides excellent maneuverability
- ISO tested
- Transit-ready option
- Width adjustable
- Depth adjustable

pdg
MOBILITY

WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

WWW.PDGMOBILITY.COM

FUZE T50 JR

JUNIOR TILT-IN-SPACE AT ITS BEST!



Nine standard frame colors to choose from and custom colors also available. Contact us for details.



Adjustable seat height



Adjustable axle position



Adjustable seat depth and back rest



Transit tie-downs



Adjustable seat width

Seat width, adjustable 12" – 18"
Seat depth, adjustable 13" – 18"
Seat height, adjustable 13" – 20"
Back angle, adjustable 90° – 120°
Back height 20" and 25"
Dynamic tilt 50°

Overall width
 12" & 16" rear wheels
 20", 22", 24" rear wheels

seat width + 7.5"
 seat width + 9"

Overall Length
Weight Capacity
Chair Weight

Transit Options Available
Warranty

24" (no front rigging)
 180 lbs
 31 lbs (without front rigging, wheels, anti-tippers, armrests)
 ANSI/RESNA WC-19
 PDG's exclusive TGIF warranty (Thank Goodness It Fits)
 Lifetime guarantee on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product.

USA
 PDG Product Design Group
 PMB # 173 - 1124 Fir Avenue
 Blaine, WA
 98230, USA

CANADA
 PDG Product Design Group
 Unit 103 - 318 East Kent Avenue South
 Vancouver, BC
 V5X 4N6

T: +1 604 323 9220
 F: +1 604 323 9097
 Toll Free (in N. America): 888-858-4422
 Toll Free Fax (in N. America): 877-323-9093
www.pdgmobility.com