

| | |
|--|-----|
| Purchase Order | Tag |
| Dealer/Organization Name | |
| Purchaser Name and Contact Information | |
| Mailing Address | |
| Shipping Address | |

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

| |
|--|
| |
| |
| |

Frame Type

1701000 **Fuze JR 50° Manual Tilt Wheelchair** \$ 3,595

Frame Options

- 1702025 **Adjustable/Growing Base and Seat Frame**
180 lbs Weight Capacity.
- 1702300 **Transit Tie Down System (TTDS)**
TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.

Seat Width

Adjustable Seat Widths

- 1703112 12"-15" Range - Set at 12"
- 1703113 12"-15" Range - Set at 13"
- 1703114 12"-15" Range - Set at 14"
- 1703115 12"-15" Range - Set at 15"
- 1703215 15"-18" Range - Set at 15"
- 1703216 15"-18" Range - Set at 16"
- 1703217 15"-18" Range - Set at 17"
- 1703218 15"-18" Range - Set at 18"

Seat Depth

Standard Seat Depths (Adjustable)

- 1704013 13" Seat Depth
- 1704014 14" Seat Depth
- 1704015 15" Seat Depth
- 1704016 16" Seat Depth
- 1704017 17" Seat Depth
- 1704018 18" Seat Depth

Frame Colour

PDG Primary Colours

- 1705000 Tracer Red
- 1705001 Gloss Black
- 1705002 Silver
- 1705003 Black Shimmer
- 1705004 Orange Fusion
- 1705005 Karl's Imagination
- 1705006 Blue Shimmer
- 1705007 Gloss White
- 1705008 Dark Titanium

PDG Premium Colours

- 1705009 Geyser Green
- 1705010 Lollipop Blue
- 1705011 Lollipop Red

N/C

1705200 **Custom Colour**

Contact PDG for specifications.

N/C

Arm Type

1706000 **No Armrest**

1706100 **Single Post Fixed Height (12")**

1706111 **Single Post Adjustable Height (8"-11")**

\$ 255

1706200 **Dual Post Fixed Height (12")**

Dual Post Adjustable Height

\$ 255

1706211 8"-11"

1706214 10"-14"

Locking Flip-up

\$ 255

Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.

1706308 Set at 8"

1706309 Set at 9"

1706310 Set at 10"

1706311 Set at 11"

1706312 Set at 12"

1706313 Set at 13"

1706314 Set at 14"

Armrests

- 1707095 **Curved Armpads (pair)**
Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests.
- Regular Armpads**
- 1707100 Desk Length - 2" x 10" (pair)
- 1707101 Full Length - 2" x 14" (pair)
- Wide Armpads** \$ 115
May interfere with wheels in some tilt configurations.
- 1707110 Desk Length - 3" x 10" (pair)
- 1707111 Full Length - 3" x 12" (pair)
- Gel Armpads** \$ 115
May interfere with wheels in some tilt configurations.
- 1707120 Full Length - 2" x 12" (pair)
- 1707121 Full Length - 3.5" x 12" (pair)
- Contour Arm Trough with Hand Pad** \$ 125
Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configurations.
- 1707130 Right (4.5" x 18")
- 1707131 Left (4.5" x 18")

Seat Pan

- 1708000 **No Seat Pan**
- 1708001 **Solid Seat Pan** \$ 120

Cushion

- 1709000 **No Cushion**
- next BioFit**
- 1709101 12"-18" width, 12"-18" depth \$ 595
- Cushion Size**
Available in even number sizes only.
- 1709400 Specify Size (Width" x Depth") _____

Seat Height

Circle desired Seat-to-Floor Height or enter information below:

STF Height is approximate and can vary by as much as 1" depending on rear wheel type and pre-tilt setting (Standard JR Chairs are shipped in 3° pre-tilt). To resolve, STF Height can be adjusted down by removing the pre-tilt collar/bumper.

| RW | 5" Casters | | | | | |
|-----|------------|-----|-----|-----|--|--|
| 12" | 13" | 14" | | | | |
| 16" | | 14" | 15" | 16" | | |
| 20" | | | 16" | 17" | | |
| 22" | | | | 17" | | |

| RW | 6" Casters & 8" Forks | | | | |
|-----|-----------------------|-----|-----|-----|-----|
| 12" | | | | 18" | 19" |
| 16" | | 15" | 16" | | |
| 20" | | | 16" | 17" | 18" |
| 22" | | | | 17" | 18" |
| 24" | | | | | 18" |

| RW | 6" Casters & 6" Forks | | | | |
|-----|-----------------------|-----|-----|-----|-----|
| 12" | | 14" | | 18" | |
| 16" | | 14" | 15" | 16" | |
| 20" | | | 16" | 17" | 18" |
| 22" | | | | 17" | 18" |
| 24" | | | | | 18" |

| RW | 8" Casters & 8" Forks | | | | |
|-----|-----------------------|-----|-----|-----|-----|
| 12" | | | | 18" | 19" |
| 16" | | 16" | | | 20" |
| 20" | | | 16" | 17" | 18" |
| 22" | | | | 17" | 18" |
| 24" | | | | | 18" |

Wheel size _____ Caster size _____ STF _____



Scan this code or visit www.pdgmobility.com to learn more about the Fuze JR Manual Tilt Wheelchair.

Castors

Urethane

- 1711051 5" x 1.25" Caster size
- 1711061 6" x 1.25" Caster size
- 1711063 6" x 1.50" Caster size \$ 85
- 1711081 8" x 1.25" Caster size
- 1711082 8" x 2.00" Caster size \$ 135

Pneumatic

- 1711161 6" x 1.25" Caster size \$ 110
- 1711181 8" x 1.25" Caster size \$ 120
- 1711182 8" x 2.00" Caster size \$ 150

Rear Wheels

Composite

- 1712012 12" Wheel size
- 1712016 16" Wheel size
- 1712020 20" Wheel size
- 1712022 22" Wheel size
- 1712024 24" Wheel size

Spoke

- 1712120 20" Wheel size \$ 75
- 1712122 22" Wheel size
- 1712124 24" Wheel size

Axle Type

- 1713000 **Quick Release Rear Axles**
- 1713010 **Fixed Rear Axles**

Rear Tires

- 1714000 **Urethane**
- 1714010 **Pneumatic** \$ 60
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.
- 1714020 **Pneumatic with Flat-Free Inserts** \$ 95
Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.

Hand Rims

- 1715000 **No Hand Rims**
- 1715010 **Standard** 20", 22", and 24" Rear Wheels only.
- 1715020 **Plastic Coated** 20", 22", and 24" Rear Wheels only. \$ 95

| Wheel Locks | | |
|--------------------------|---------|--|
| <input type="checkbox"/> | 1716005 | Foot Operated Wheel Lock \$ 200 |
| <input type="checkbox"/> | 1716010 | Push to Lock |
| <input type="checkbox"/> | 1716020 | Pull to Lock |
| <input type="checkbox"/> | 1716030 | Attendant Wheel Lock (Push to Lock) \$ 85 Not compatible with Wheel Lock Extensions. |
| <input type="checkbox"/> | 1716100 | Wheel Lock Extensions \$ 60 Not compatible with Attendant Wheel Lock. |

| Anti-Tippers | | |
|-------------------------------------|---------|---|
| <input checked="" type="checkbox"/> | 1717010 | Rear Anti-Tippers \$ 60 Required Safety Option. |

| Back Post Style | | |
|---|---------|---------------|
| Push Bar Back | | |
| <input type="checkbox"/> | 1718010 | Straight |
| <input type="checkbox"/> | 1718030 | 8° Bend \$ 75 |
| Cane Handle Back Not compatible with Headrest Mounting Plate. | | |
| <input type="checkbox"/> | 1718020 | Straight |
| <input type="checkbox"/> | 1718040 | 8° Bend \$ 75 |
| Back Post Height Seat Pan to Handle | | |
| <input type="checkbox"/> | 1718420 | 20" Height |
| <input type="checkbox"/> | 1718425 | 25" Height |

| Back Options | | |
|---|---------|--|
| Angle/Height Adjustable Push Bar Handle \$ 150 Not compatible with Headrest Mounting Plate or Cane Handle Back. | | |
| <input type="checkbox"/> | 1719010 | Button on Outside |
| <input type="checkbox"/> | 1719011 | Button on Inside |
| <input type="checkbox"/> | 1719050 | Removable Cane Handle Push Bar \$ 140 Not compatible with Push Bar Back. |
| Seat to Back Angle | | |
| <input type="checkbox"/> | 1719090 | Set at 90° |
| <input type="checkbox"/> | 1719096 | Set at 96° |
| <input type="checkbox"/> | 1719102 | Set at 102° |
| <input type="checkbox"/> | 1719108 | Set at 108° |
| <input type="checkbox"/> | 1719114 | Set at 114° |
| <input type="checkbox"/> | 1719120 | Set at 120° |
| <input type="checkbox"/> | 1719220 | Fold Down Back |
| <input type="checkbox"/> | 1719221 | Non Fold Down Back \$ 75 |

| Back Upholstery | | |
|-------------------------------------|---------|---------------------------|
| <input checked="" type="checkbox"/> | 1720000 | No Back Upholstery |

| Backrest | | |
|--|---------------------------------------|--------------------------------------|
| <input type="checkbox"/> | 1721000 | No Backrest |
| next Optima HA Active Contour | | |
| <input type="checkbox"/> | 1721100 | 16"-18" width, 16"-20" height \$ 690 |
| next Optima HA Plus Contour | | |
| <input type="checkbox"/> | 1721101 | 16"-18" width, 16"-20" height \$ 690 |
| Backrest Size Available in even number sizes only. | | |
| 1721400 | Specify Size (Width" x Height") _____ | |

| Tilt/Recline Options | | |
|----------------------------|---------|--|
| Attendant Tilt only | | |
| <input type="checkbox"/> | 1723011 | Installed on Right Back Post |
| <input type="checkbox"/> | 1723012 | Installed on Left Back Post |
| Tilt Limit | | |
| <input type="checkbox"/> | 1723220 | Set at 20° |
| <input type="checkbox"/> | 1723225 | Set at 25° |
| <input type="checkbox"/> | 1723230 | Set at 30° |
| <input type="checkbox"/> | 1723235 | Set at 35° |
| <input type="checkbox"/> | 1723240 | Set at 40° |
| <input type="checkbox"/> | 1723245 | Set at 45° |
| <input type="checkbox"/> | 1723250 | Set at 50° |
| <input type="checkbox"/> | 1723410 | Reclining Backrest (Manual) \$ 830 Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back or Tension Adjustable Upholstery. Not Available with 13" Seat Depth. |

| Headrest | | |
|--|---------|--|
| <input type="checkbox"/> | 1724000 | No Headrest |
| <input type="checkbox"/> | 1724010 | Headrest Mounting Plate \$ 75 Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back. |
| <input type="checkbox"/> | 1724100 | next MultiFit Head Support Hardware \$ 155 |
| next MultiFit Head Support Pad \$ 195 | | |
| <input type="checkbox"/> | 1724105 | Regular (9" x 5") |
| <input type="checkbox"/> | 1724106 | Large (10" x 5") |

| Front Rigging | | |
|--------------------------|---------|--|
| <input type="checkbox"/> | 1725000 | No Front Rigging |
| <input type="checkbox"/> | 1725001 | Sector Blocks only |
| <input type="checkbox"/> | 1725003 | External Front Rigging Mount \$ 75 |
| <input type="checkbox"/> | 1725140 | Swing Away 70° |
| <input type="checkbox"/> | 1725160 | Swing Away 90° Not compatible with 8" Casters |
| <input type="checkbox"/> | 1725200 | Swing Away Elevating Legrest \$ 355 Not available for Seat Widths of 12" or 13". |

| Foot Plates | | |
|--------------------------|--|--------|
| <input type="checkbox"/> | 1726010 Composite Not available on 12" Seat Width. | |
| <input type="checkbox"/> | 1726020 Angle Adjustable | \$ 140 |
| <input type="checkbox"/> | 1726040 One Piece Flip-up | \$ 395 |
| <input type="checkbox"/> | 1726050 Interlocking Angle Adjustable Seat Widths of 14"-18" only. | \$ 195 |

| Footrest/Legrest Accessories | | |
|------------------------------|---|--------|
| <input type="checkbox"/> | 1727015 Swing Away 70° Calf Pad Not compatible with 70° Gel Pad or Seat Widths of 12" or 13". | \$ 135 |
| <input type="checkbox"/> | 1727020 Swing Away 70° Gel Pad Not compatible with 70° Swing Away Calf Pad. | \$ 100 |
| <input type="checkbox"/> | 1727030 Swing Away ELR Gel Knee Pad | \$ 100 |
| <input type="checkbox"/> | 1727040 Adjustable Length Heel Loops | \$ 50 |
| <input type="checkbox"/> | 1727050 Bodypoint - Padded Calf Strap | \$ 105 |
| <input type="checkbox"/> | 1727065 Miller's Dynamic Footrest Mounts | \$ 380 |

| Extension Tubes | | |
|---|---|-------|
| Regular Extension Tubes Not compatible with Swing Away 90°. | | |
| <input type="checkbox"/> | 1728010 Standard (13"-17") | |
| <input type="checkbox"/> | 1728020 Long (15"-19") | \$ 75 |
| <input type="checkbox"/> | 1728200 90° Extension Tubes (11"-16") Not compatible with 70° or ELR. | \$ 75 |

| Positioning Straps | | |
|--|--|--------|
| Bodypoint - Pediatric Hip Belts | | |
| <input type="checkbox"/> | 1729017 Center-Pull Auto Style, Padded | \$ 125 |
| <input type="checkbox"/> | 1729022 4-Point, Padded | \$ 150 |
| <input type="checkbox"/> | 1729030 Bodypoint - Chest Strap | \$ 85 |

| Miscellaneous Options | | |
|--------------------------|--|--------|
| <input type="checkbox"/> | 1730160 Clear Lap Tray (Slide On) Not compatible with 8" - 11" Dual Post Arms, Desk Length Pads, or Gel Arm Pads. Seat Widths of 15"-18" only. | \$ 490 |
| <input type="checkbox"/> | 1730170 Adhesive Back Hook Strip, 2" x 30" For Cushion mounting. | \$ 10 |