

|  |     |
|--|-----|
| Purchase Order                         | Tag |
| Dealer/Organization Name               |     |
| Purchaser Name and Contact Information |     |
| Mailing Address                        |     |
| Shipping Address                       |     |

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877-323-9093.

### Order Notes

|  |
|--|
|  |
|  |
|  |

### Frame Type

- 1901000 **Stellar GL 45° Manual Tilt Wheelchair** \$ 3,995

### Frame Options

- 1902000 **Standard Base and Seat Frame**  
250 lbs Weight Capacity.
- 1902300 **Transit Tie Down System (TTDS)**  
TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

### Seat Width

#### Standard Seat Widths

- 1903014 14" Seat Width
- 1903015 15" Seat Width
- 1903016 16" Seat Width
- 1903017 17" Seat Width
- 1903018 18" Seat Width
- 1903019 19" Seat Width
- 1903020 20" Seat Width
- 1903021 21" Seat Width \$ 215
- 1903022 22" Seat Width \$ 215

### Seat Depth

#### Standard Seat Depths (Adjustable)

- 1904016 16" Seat Depth
- 1904017 17" Seat Depth
- 1904018 18" Seat Depth
- 1904019 19" Seat Depth
- 1904020 20" Seat Depth
- Custom Seat Depths \$ 325**
- 1904021 21" Seat Depth
- 1904022 22" Seat Depth

### Frame Colour

#### PDG Primary Colours

- 1905000 Tracer Red
- 1905001 Gloss Black
- 1905002 Silver
- 1905003 Black Shimmer
- 1905004 Orange Fusion
- 1905005 Karl's Imagination
- 1905006 Blue Shimmer
- 1905007 Gloss White
- 1905008 Dark Titanium

#### PDG Premium Colours

- 1905009 Geysir Green \$ 150
- 1905010 Lollipop Blue
- 1905011 Lollipop Red

- 1905200 **Custom Colour** \$ 300  
Contact PDG for specifications.

### Arm Type

- 1906000 **No Armrest**
- 1906100 **Single Post Fixed Height (12")**  
Full Length Armpad only.
- Single Post Adjustable Height \$ 255**
- 1906111 8"-11"
- 1906114 10"-14"
- 1906200 **Dual Post Fixed Height (12")**  
Full Length Armpad only.
- Dual Post Adjustable Height \$ 255**
- 1906211 8"-11"
- 1906214 10"-14"
- 1906216 12"-16"
- Locking Flip-up \$ 255**  
Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.
- 1906310 Set at 10"
- 1906311 Set at 11"
- 1906312 Set at 12"
- 1906313 Set at 13"
- 1906314 Set at 14"



Scan this code or visit [www.pdgmobility.com](http://www.pdgmobility.com) to learn more about the Stellar GL Manual Tilt Wheelchair.

| Armrests   |   |
|--|---|
| <b>Regular Armrests</b>  |   |
| <input type="checkbox"/> 1907100   | Desk Length - 2" x 10" (pair)                         |
| <input type="checkbox"/> 1907101   | Full Length - 2" x 14" (pair)                         |
| <b>Wide Armrests</b> \$ 115  |   |
| May interfere with wheels in some tilt configurations.                           |   |
| <input type="checkbox"/> 1907110   | Desk Length - 3" x 10" (pair)                         |
| <input type="checkbox"/> 1907111   | Full Length - 3" x 12" (pair)                         |
| <b>Gel Armrests</b> \$ 115   |   |
| May interfere with wheels in some tilt configurations.                           |   |
| <input type="checkbox"/> 1907120   | Full Length - 2" x 12" (pair)                         |
| <input type="checkbox"/> 1907121   | Full Length - 3.5" x 12" (pair)                       |
| <b>Contour Arm Trough</b> \$ 125   |   |
| May interfere with wheels in some tilt configurations.                           |   |
| <input type="checkbox"/> 1907130   | Right (4.5" x 18") with Hand Pad                      |
| <input type="checkbox"/> 1907131   | Left (4.5" x 18") with Hand Pad                       |
| <input type="checkbox"/> 1907220   | <b>Armrest Wideners for Dual Post Armrests</b> \$ 375 |
| Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels or Clothing Guards. |   |
| <b>Clothing Guards</b>   |   |
| Not compatible with Flip-up Armrests.  |   |
| <input type="checkbox"/> 1907310   | For Single Post Armrests \$ 110                       |
| 10"-14" Height range only  |   |
| <input type="checkbox"/> 1907320   | For Dual Post Armrests \$ 165                         |
| <b>Armrest Extensions</b> \$ 60  |   |
| Not compatible with Flip-up or Single Post Armrests.                             |   |
| <input type="checkbox"/> 1907400   | 3" Extension - Forward                                |
| <input type="checkbox"/> 1907401   | 3" Extension - Rearward                               |

| Seat Pan                         |                              |
|----------------------------------|------------------------------|
| <input type="checkbox"/> 1908000 | <b>No Seat Pan</b>           |
| <input type="checkbox"/> 1908001 | <b>Solid Seat Pan</b> \$ 120 |

| Cushion                              |                                      |
|--------------------------------------|--------------------------------------|
| <input type="checkbox"/> 1909000     | <b>No Cushion</b>                    |
| <b>nxt BioFit</b>                    |                                      |
| <input type="checkbox"/> 1909101     | 14"-20" width, 16"-20" depth \$ 595  |
| <input type="checkbox"/> 1909102     | 22"-24" width, 16"-20" depth \$ 625  |
| <b>Cushion Size</b>                  |                                      |
| Available in even number sizes only. |                                      |
| 1909400                              | Specify Size (Width" x Depth") _____ |

| Seat Height  |            |            |     |     |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
|--|------------|------------|-----|-----|-----|-----|--|-----|-----|-----|-----|-----|--|--|-----|--|--|-----|-----|-----|-----|-----|-----|-----|-----|-----|--|--|-----|--|-----|-----|-----|-----|--|-----|--|--|-----|-----|-----|-----|--|----|------------|--|--|--|-----|--|-----|-----|-----|-----|-----|--|-----|-----|-----|-----|
| Circle desired Seat-to-Floor Height or enter information below:  |            |            |     |     |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| <table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>20"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td></td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table> | RW         | 6" Casters |     |     |     |     |  | 12" | 13" | 14" | 15" | 16" |  |  | 16" |  |  | 15" | 16" | 17" | 18" | 20" | 13" | 14" | 15" | 16" |  |  | 22" |  | 14" | 15" | 16" | 17" |  | 24" |  |  | 15" | 16" | 17" | 18" | <table border="1"> <thead> <tr> <th>RW</th> <th colspan="4">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>24"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table> | RW | 8" Casters |  |  |  | 16" |  | 15" | 16" | 17" | 18" | 24" |  | 15" | 16" | 17" | 18" |
| RW   | 6" Casters |            |     |     |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| 12"  | 13"        | 14"        | 15" | 16" |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| 16"  |            |            | 15" | 16" | 17" | 18" |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| 20"  | 13"        | 14"        | 15" | 16" |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| 22"  |            | 14"        | 15" | 16" | 17" |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| 24"  |            |            | 15" | 16" | 17" | 18" |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| RW   | 8" Casters |            |     |     |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| 16"  |            | 15"        | 16" | 17" | 18" |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| 24"  |            | 15"        | 16" | 17" | 18" |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| Wheel size _____ Caster size _____ STF _____   |            |            |     |     |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |
| STF Height is adjustable within the range of the selected wheel size.  |            |            |     |     |     |     |  |     |     |     |     |     |  |  |     |  |  |     |     |     |     |     |     |     |     |     |  |  |     |  |     |     |     |     |  |     |  |  |     |     |     |     |  |    |            |  |  |  |     |  |     |     |     |     |     |  |     |     |     |     |

| Casters   |                                |
|---|--------------------------------|
| <b>Urethane</b>   |                                |
| 8" Casters are compatible only with 16" and 24" Rear Wheels.                    |                                |
| <input type="checkbox"/> 1911061  | 6" x 1.25" Caster size         |
| <input type="checkbox"/> 1911063  | 6" x 1.50" Caster size \$ 85   |
| <input type="checkbox"/> 1911081  | 8" x 1.25" Caster size         |
| <input type="checkbox"/> 1911084  | 8" x 2.00" Caster size \$ 135  |
| <b>Pneumatic</b>  |                                |
| 8" Casters are compatible only with 16" and 24" Rear Wheels.                    |                                |
| <input type="checkbox"/> 1911161  | 6" x 1.25" Caster size \$ 110  |
| <input type="checkbox"/> 1911181  | 8" x 1.25" Caster size \$ 120  |
| <input type="checkbox"/> 1911184  | 8" x 2.00" Caster size \$ 150  |
| <input type="checkbox"/> 1911260  | <b>6" Soft Roll</b> \$ 95      |
| <input type="checkbox"/> 1911400  | <b>Ultra Sport Fork</b> \$ 395 |
| Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. |                                |

| Rear Wheels                      |                                    |
|----------------------------------|------------------------------------|
| <b>Composite</b>                 |                                    |
| <input type="checkbox"/> 1912012 | 12" Wheel size                     |
| <input type="checkbox"/> 1912016 | 16" Wheel size                     |
| <input type="checkbox"/> 1912020 | 20" Wheel size                     |
| <input type="checkbox"/> 1912022 | 22" Wheel size                     |
| <input type="checkbox"/> 1912024 | 24" Wheel size                     |
| <b>Spoke</b> \$ 75               |                                    |
| <input type="checkbox"/> 1912120 | 20" Wheel size                     |
| <input type="checkbox"/> 1912122 | 22" Wheel size                     |
| <input type="checkbox"/> 1912124 | 24" Wheel size                     |
| <input type="checkbox"/> 1912224 | <b>24" Heavy Duty Spoke</b> \$ 350 |

| Axle Type                        |                                 |
|----------------------------------|---------------------------------|
| <input type="checkbox"/> 1913000 | <b>Quick Release Rear Axles</b> |
| <input type="checkbox"/> 1913010 | <b>Fixed Rear Axles</b>         |

| Rear Tires  |   |
|---|---|
| <input type="checkbox"/> 1914000  | <b>Urethane</b>                               |
| <input type="checkbox"/> 1914010  | <b>Pneumatic</b> \$ 60                        |
| Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.                                    |   |
| <input type="checkbox"/> 1914020  | <b>Pneumatic with Flat-Free Inserts</b> \$ 95 |
| Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels. |   |

| Hand Rims                        |   |
|----------------------------------|---|
| <input type="checkbox"/> 1915000 | <b>No Hand Rims</b>   |
| <input type="checkbox"/> 1915010 | <b>Standard</b> 20", 22", and 24" Rear Wheels only.             |
| <input type="checkbox"/> 1915020 | <b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only. \$ 95 |

| Wheel Locks                      |  |
|----------------------------------|--|
| <input type="checkbox"/> 1916010 | <b>Push to Lock</b>  |
| <input type="checkbox"/> 1916020 | <b>Pull to Lock</b>  |
| <input type="checkbox"/> 1916030 | <b>Attendant Wheel Lock (Push to Lock)</b> \$ 85<br>Not compatible with Wheel Lock Extensions.   |
| <input type="checkbox"/> 1916040 | <b>User &amp; Attendant Wheel Lock</b> \$ 125<br>Available only on 20", 22", or 24" Rear Wheels. |
| <input type="checkbox"/> 1916100 | <b>Wheel Lock Extensions</b> \$ 60<br>Not compatible with Attendant Wheel Lock.                  |
| <input type="checkbox"/> 1916200 | <b>Hill Holders</b> \$ 350<br>Requires Pneumatic Rear Tires with Flat-Free Inserts.              |

| Anti-Tippers                                |   |
|---|---|
| <input checked="" type="checkbox"/> 1917010 | <b>Rear Anti-Tippers</b> \$ 60<br>Required Safety Option.   |
| <input type="checkbox"/> 1917020            | <b>Front Anti-Tippers</b> \$ 300  |
| <input type="checkbox"/> 1917030            | <b>Heavy Duty Rear Anti-Tipper</b> \$ 600<br>Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder. |

| Back Post Style                               |   |
|---|---|
| <input type="checkbox"/> 1918010              | <b>Push Bar Back</b>  |
| <input type="checkbox"/> 1918020              | <b>Cane Handle Back</b><br>Not compatible with Headrest Mounting Plate. |
| <b>Back Post Height</b><br>Seat Pan to Handle |   |
| <input type="checkbox"/> 1918425              | 25" Height  |
| <input type="checkbox"/> 1918431              | 31" Height  |

| Back Options                     |   |
|----------------------------------|---|
| <input type="checkbox"/> 1919010 | <b>Angle/Height Adjustable Push Bar Handle</b> \$ 150<br>Not compatible with Headrest Mounting Plate or Cane Handle Back. |
| <input type="checkbox"/> 1919010 | Button on Outside   |
| <input type="checkbox"/> 1919011 | Button on Inside  |
| <b>Seat to Back Angle</b>        |   |
| <input type="checkbox"/> 1919090 | Set at 90°  |
| <input type="checkbox"/> 1919098 | Set at 98°  |
| <input type="checkbox"/> 1919106 | Set at 106°   |
| <input type="checkbox"/> 1919114 | Set at 114°   |
| <input type="checkbox"/> 1919122 | Set at 122°   |
| <input type="checkbox"/> 1919210 | <b>Backrest Extension 10" high</b> \$ 120<br>Not compatible with Angle/Height Adjustable Push Bar Handle.                 |
| <input type="checkbox"/> 1919220 | <b>Fold Down Back</b><br>Not compatible with Reclining Back.  |
| <input type="checkbox"/> 1919230 | <b>Spring Back Suspension</b> \$ 105<br>Not compatible with Reclining Back.   |

|                                  |  |
|----------------------------------|--|
| <input type="checkbox"/> 1919240 | <b>Dynamic Rocker Back Assembly</b> \$ 2,995<br>Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller. Requires longer Frame. |
|----------------------------------|--|

| Back Upholstery  |                           |
|--|---------------------------|
| <input type="checkbox"/> 1920000   | <b>No Back Upholstery</b> |
| <b>Standard</b> \$ 125<br>Back height measured from Seat Pan to top of Upholstery              |                           |
| <input type="checkbox"/> 1920021   | Set at 21"                |
| <input type="checkbox"/> 1920022   | Set at 22"                |
| <b>Short (25" Back only)</b> \$ 125<br>Back height measured from Seat Pan to top of Upholstery |                           |
| <input type="checkbox"/> 1920118   | Set at 18"                |
| <input type="checkbox"/> 1920119   | Set at 19"                |
| <input type="checkbox"/> 1920120   | Set at 20"                |
| <b>Tall (31" Back only)</b> \$ 125<br>Back height measured from Seat Pan to top of Upholstery  |                           |
| <input type="checkbox"/> 1920224   | Set at 24"                |
| <input type="checkbox"/> 1920225   | Set at 25"                |
| <input type="checkbox"/> 1920226   | Set at 26"                |
| <b>Adjustable Tension</b> \$ 250<br>Back height measured from Seat Pan to top of Upholstery    |                           |
| <input type="checkbox"/> 1920321   | Set at 21"                |
| <input type="checkbox"/> 1920322   | Set at 22"                |
| <input type="checkbox"/> 1920323   | Set at 23"                |

| Backrest   |                                       |
|--|---------------------------------------|
| <input type="checkbox"/> 1921000                             | <b>No Backrest</b>                    |
| <b>next Optima HA Active Contour</b>                         |                                       |
| <input type="checkbox"/> 1921100                             | 16"-20" width, 16"-20" height \$ 690  |
| <input type="checkbox"/> 1921105                             | 22"-24" width, 16"-20" height \$ 720  |
| <b>next Optima HA Plus Contour</b>                           |                                       |
| <input type="checkbox"/> 1921101                             | 16"-20" width, 16"-20" height \$ 690  |
| <input type="checkbox"/> 1921106                             | 22"-24" width, 16"-20" height \$ 720  |
| <b>Backrest Size</b><br>Available in even number sizes only. |                                       |
| 1921400  | Specify Size (Width" x Height") _____ |

| Tilt/Recline Options             |   |
|----------------------------------|---|
| <b>Attendant Tilt only</b>       |   |
| <input type="checkbox"/> 1923011 | Installed on Right Back Post  |
| <input type="checkbox"/> 1923012 | Installed on Left Back Post   |
| <b>Tilt Limit</b>                |   |
| <input type="checkbox"/> 1923220 | Set at 20°  |
| <input type="checkbox"/> 1923230 | Set at 30°  |
| <input type="checkbox"/> 1923240 | Set at 40°  |
| <input type="checkbox"/> 1923410 | <b>Reclining Backrest (Manual)</b> \$ 830<br>Maximum recline = 30°. Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, 20"-22" seat depth. Requires longer Frame. |

| Headrest                         |  |
|----------------------------------|--|
| <input type="checkbox"/> 1924000 | <b>No Headrest</b>   |
| <input type="checkbox"/> 1924010 | <b>Headrest Mounting Plate</b> \$ 75<br>Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back. |
| <input type="checkbox"/> 1924100 | <b>nxt MultiFit Head Support Hardware</b> \$ 155   |
|                                  | <b>nxt MultiFit Head Support Pad</b> \$ 195  |
| <input type="checkbox"/> 1924105 | Regular (9" x 5")  |
| <input type="checkbox"/> 1924106 | Large (10" x 5")   |
| <input type="checkbox"/> 1924107 | X-Large (14" x 5")   |

| Front Rigging                    |   |
|----------------------------------|---|
| <input type="checkbox"/> 1925000 | <b>No Front Rigging</b>                     |
| <input type="checkbox"/> 1925001 | <b>Sector Blocks only</b>                   |
| <input type="checkbox"/> 1925110 | <b>Plug In 70°</b>                          |
| <input type="checkbox"/> 1925120 | <b>Plug In Adjustable Knee Angle</b> \$ 975 |
| <input type="checkbox"/> 1925130 | <b>Swing Away 60°</b>                       |
| <input type="checkbox"/> 1925140 | <b>Swing Away 70°</b>                       |
| <input type="checkbox"/> 1925200 | <b>Swing Away Elevating Legrest</b> \$ 355  |

| Foot Plates                      |   |
|----------------------------------|---|
| <input type="checkbox"/> 1926010 | <b>Composite</b>  |
| <input type="checkbox"/> 1926020 | <b>Angle Adjustable</b> \$ 140  |
| <input type="checkbox"/> 1926030 | <b>One Piece Fixed</b> \$ 350<br>Available on Plug In Front Rigging only.   |
| <input type="checkbox"/> 1926040 | <b>One Piece Flip-up</b> \$ 395   |
| <input type="checkbox"/> 1926050 | <b>Interlocking Angle Adjustable</b> \$ 195<br>Seat Widths of 14"-21" only. |

| Footrest/Legrest Accessories     |  |
|----------------------------------|--|
| <input type="checkbox"/> 1927010 | <b>Swing Away 60° &amp; 70° Calf Pad</b> \$ 135<br>Not compatible with 60° & 70° Gel Pad.            |
| <input type="checkbox"/> 1927020 | <b>Swing Away 60° &amp; 70° Gel Pad</b> \$ 100<br>Not compatible with 60° & 70° Swing Away Calf Pad. |
| <input type="checkbox"/> 1927030 | <b>Swing Away ELR Gel Knee Pad</b> \$ 100  |
| <input type="checkbox"/> 1927040 | <b>Adjustable Length Heel Loops</b> \$ 50  |
| <input type="checkbox"/> 1927050 | <b>Bodypoint - Padded Calf Strap</b> \$ 105  |
|                                  | <b>Miller's Dynamic Footrest Mounts</b>  |
| <input type="checkbox"/> 1927060 | For Plug In Front Rigging \$ 190   |
| <input type="checkbox"/> 1927065 | For Swing Away 70° Front Rigging \$ 380  |

| Extension Tubes                  |  |
|----------------------------------|--|
|                                  | <b>Regular Extension Tubes</b>   |
| <input type="checkbox"/> 1928010 | Standard (13"-17")   |
| <input type="checkbox"/> 1928020 | Long (15"-19") \$ 75   |
|                                  | <b>Extra Strength Extension Tubes (Solid)</b> \$ 100<br>Available on Plug In Front Rigging only. |
| <input type="checkbox"/> 1928110 | Standard (13"-17")   |
| <input type="checkbox"/> 1928120 | Long (15"-19")   |

| Positioning Straps               |                                       |
|----------------------------------|---------------------------------------|
|                                  | <b>Bodypoint - Hip Belts</b>          |
| <input type="checkbox"/> 1929010 | Center-Pull Auto Style \$ 105         |
| <input type="checkbox"/> 1929015 | Center-Pull Auto Style, Padded \$ 155 |
| <input type="checkbox"/> 1929020 | 4-Point, Padded \$ 175                |
| <input type="checkbox"/> 1929030 | <b>Bodypoint - Chest Strap</b> \$ 85  |

| Miscellaneous Options            |   |
|----------------------------------|---|
| <input type="checkbox"/> 1930110 | <b>IV Pole</b> \$ 295<br>Not compatible with HD Anti Tipper, Attendant Wheel Locks.   |
| <input type="checkbox"/> 1930120 | <b>O2 Holder</b> \$ 295<br>Not compatible with HD Anti Tipper, Attendant Wheel Locks. |
| <input type="checkbox"/> 1930170 | <b>Adhesive Back Hook Strip, 2" x 30"</b> \$ 10<br>For Cushion mounting.              |