

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877-323-9093.

**Order Notes**

---



---



---

Frame Type			
<input checked="" type="checkbox"/>	1901000 <b>Stellar GL 45° Manual Tilt Wheelchair</b>	<i>E1161</i>	\$ 2,995

Frame Options	
<input checked="" type="checkbox"/>	1902000 <b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.
<input type="checkbox"/>	1902300 <b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width			
Standard Seat Widths			
<input type="checkbox"/>	1903014	14" Seat Width	
<input type="checkbox"/>	1903015	15" Seat Width	
<input type="checkbox"/>	1903016	16" Seat Width	
<input type="checkbox"/>	1903017	17" Seat Width	
<input type="checkbox"/>	1903018	18" Seat Width	
<input type="checkbox"/>	1903019	19" Seat Width	
<input type="checkbox"/>	1903020	20" Seat Width	<i>E2201</i> \$ 325
<input type="checkbox"/>	1903021	21" Seat Width	<i>E2201</i> \$ 325
<input type="checkbox"/>	1903022	22" Seat Width	<i>E2201</i> \$ 325

Seat Depth			
Standard Seat Depths (Adjustable)			
<input type="checkbox"/>	1904016	16" Seat Depth	
<input type="checkbox"/>	1904017	17" Seat Depth	
<input type="checkbox"/>	1904018	18" Seat Depth	
<input type="checkbox"/>	1904019	19" Seat Depth	
<input type="checkbox"/>	1904020	20" Seat Depth	<i>E2203</i> \$ 475
Custom Seat Depths			
<input type="checkbox"/>	1904021	21" Seat Depth	<i>E2203</i> \$ 475
<input type="checkbox"/>	1904022	22" Seat Depth	<i>E2204</i> \$ 750

Frame Color		
PDG Primary Colors		
<input type="checkbox"/>	1905000 Tracer Red	
<input type="checkbox"/>	1905001 Gloss Black	
<input type="checkbox"/>	1905002 Silver	
<input type="checkbox"/>	1905003 Black Shimmer	
<input type="checkbox"/>	1905004 Orange Fusion	
<input type="checkbox"/>	1905005 Karl's Imagination	
<input type="checkbox"/>	1905006 Blue Shimmer	
<input type="checkbox"/>	1905007 Gloss White	
<input type="checkbox"/>	1905008 Dark Titanium	
PDG Premium Colors		\$ 150
<input type="checkbox"/>	1905009 Geyser Green	
<input type="checkbox"/>	1905010 Lollipop Blue	
<input type="checkbox"/>	1905011 Lollipop Red	
Custom Color		\$ 300
<input type="checkbox"/> 1905200 <b>Custom Color</b>		
Contact PDG for specifications.		

Arm Type		
<input type="checkbox"/>	1906000 <b>No Armrest</b>	
<input type="checkbox"/>	1906100 <b>Single Post Fixed Height (12")</b> Full Length Armpad only.	
Single Post Adjustable Height		<i>E0973</i> \$ 255
<input type="checkbox"/>	1906111 8"-11"	
<input type="checkbox"/>	1906114 10"-14"	
<input type="checkbox"/>	1906200 <b>Dual Post Fixed Height (12")</b> Full Length Armpad only.	
Dual Post Adjustable Height		<i>E0973</i> \$ 255
<input type="checkbox"/>	1906211 8"-11"	
<input type="checkbox"/>	1906214 10"-14"	
<input type="checkbox"/>	1906216 12"-16"	
Locking Flip-up		<i>E0973</i> \$ 255
<input type="checkbox"/> 1906310 Set at 10"		
<input type="checkbox"/>	1906311 Set at 11"	
<input type="checkbox"/>	1906312 Set at 12"	
<input type="checkbox"/>	1906313 Set at 13"	
<input type="checkbox"/>	1906314 Set at 14"	

**PDG Product Design Group Inc.**  
318 East Kent Avenue South  
Vancouver BC V5X 4N6 Canada  
Phone 604-323-9220 | 888-858-4422

Armrests			
<b>Regular Armrests</b>			
<input type="checkbox"/>	1907100	Desk Length - 2" x 10" (pair)	
<input type="checkbox"/>	1907101	Full Length - 2" x 14" (pair)	
<b>Wide Armrests</b> \$ 115			
May interfere with wheels in some tilt configurations.			
<input type="checkbox"/>	1907110	Desk Length - 3" x 10" (pair)	
<input type="checkbox"/>	1907111	Full Length - 3" x 12" (pair)	
<b>Gel Armrests</b> \$ 115			
May interfere with wheels in some tilt configurations.			
<input type="checkbox"/>	1907120	Full Length - 2" x 12" (pair)	
<input type="checkbox"/>	1907121	Full Length - 3.5" x 12" (pair)	
<b>Contour Arm Trough</b> E2209 \$ 125			
May interfere with wheels in some tilt configurations.			
<input type="checkbox"/>	1907130	Right (4.5" x 18") with Hand Pad	
<input type="checkbox"/>	1907131	Left (4.5" x 18") with Hand Pad	
<input type="checkbox"/>	1907220	<b>Armrest Wideners for Dual Post Armrests</b> \$ 375	
Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels or Clothing Guards.			
<b>Clothing Guards</b> K0108			
Not compatible with Flip-up Armrests.			
<input type="checkbox"/>	1907310	For Single Post Armrests 10"-14" Height range only	\$ 110
<input type="checkbox"/>	1907320	For Dual Post Armrests	\$ 165
<b>Armrest Extensions</b> \$ 60			
Not compatible with Flip-up or Single Post Armrests.			
<input type="checkbox"/>	1907400	3" Extension - Forward	
<input type="checkbox"/>	1907401	3" Extension - Rearward	

Seat Pan			
<input type="checkbox"/>	1908000	<b>No Seat Pan</b>	
<input type="checkbox"/>	1908001	<b>Solid Seat Pan</b> E2231	\$ 120

Cushion			
<input type="checkbox"/>	1909000	<b>No Cushion</b>	
<b>nxt BioFit</b>			
<input type="checkbox"/>	1909101	14"-20" width, 16"-20" depth E2607	\$ 550
<input type="checkbox"/>	1909102	22"-24" width, 16"-20" depth E2608	\$ 600
<b>Matrx VI</b>			
<input type="checkbox"/>	1909111	14"-20" width, 16"-22" depth E2607	\$ 550
<input type="checkbox"/>	1909112	22" width, 18"-22" depth E2608	\$ 600
<input type="checkbox"/>	1909115	HD - 18"-30" width, 18"-22" depth E2608	\$ 1,150
<b>ROHO Quadtro Select HP</b>			
<input type="checkbox"/>	1909121	14"-20" width, 16"-22" depth E2624	\$ 650
<input type="checkbox"/>	1909122	22"-24" width, 16"-22" depth E2625	\$ 750
<b>ROHO High Profile</b>			
<input type="checkbox"/>	1909125	14"-20" width, 16"-22" depth E2622	\$ 650
<input type="checkbox"/>	1909126	22"-24" width, 16"-22" depth E2623	\$ 750
<b>Varilite Reflex</b>			
<input type="checkbox"/>	1909131	14"-20" width, 16"-20" depth E2603	\$ 375
<b>Cushion Size</b>			
	1909400	Specify Size (Width" x Depth") _____	



Scan this code or visit [www.pdgmobility.com](http://www.pdgmobility.com) to learn more about the Stellar GL Manual Tilt Wheelchair.

Seat Height																																																												
Circle desired Seat-to-Floor Height or enter the information below:																																																												
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>20"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td></td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table>	RW	6" Casters						12"	13"	14"	15"	16"			16"			15"	16"	17"	18"	20"	13"	14"	15"	16"			22"		14"	15"	16"	17"		24"			15"	16"	17"	18"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="4">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>24"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table>	RW	8" Casters				16"		15"	16"	17"	18"	24"		15"	16"	17"	18"
RW	6" Casters																																																											
12"	13"	14"	15"	16"																																																								
16"			15"	16"	17"	18"																																																						
20"	13"	14"	15"	16"																																																								
22"		14"	15"	16"	17"																																																							
24"			15"	16"	17"	18"																																																						
RW	8" Casters																																																											
16"		15"	16"	17"	18"																																																							
24"		15"	16"	17"	18"																																																							
Wheel size _____ Caster size _____ STF _____																																																												
STF Height is adjustable within the range of the selected wheel size.																																																												

Casters	
<b>Urethane</b>	
8" Casters are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/>	1911061 6" x 1.25" Caster size
<input type="checkbox"/>	1911063 6" x 1.50" Caster size \$ 85
<input type="checkbox"/>	1911081 8" x 1.25" Caster size
<input type="checkbox"/>	1911084 8" x 2.00" Caster size \$ 135
<b>Pneumatic</b>	
8" Casters are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/>	1911161 6" x 1.25" Caster size \$ 110
<input type="checkbox"/>	1911181 8" x 1.25" Caster size \$ 120
<input type="checkbox"/>	1911184 8" x 2.00" Caster size \$ 150
<input type="checkbox"/>	1911260 <b>6" Soft Roll</b> \$ 95
<input type="checkbox"/>	1911400 <b>Ultra Sport Fork</b> E1015 \$ 395
Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	

Rear Wheels	
<b>Composite</b>	
<input type="checkbox"/>	1912012 12" Wheel size
<input type="checkbox"/>	1912016 16" Wheel size
<input type="checkbox"/>	1912020 20" Wheel size
<input type="checkbox"/>	1912022 22" Wheel size
<input type="checkbox"/>	1912024 24" Wheel size
<b>Spoke</b> \$ 75	
<input type="checkbox"/>	1912120 20" Wheel size
<input type="checkbox"/>	1912122 22" Wheel size
<input type="checkbox"/>	1912124 24" Wheel size
<input type="checkbox"/>	1912224 <b>24" Heavy Duty Spoke</b> \$ 350

Axle Type	
<input type="checkbox"/>	1913000 <b>Quick Release Rear Axles</b>
<input type="checkbox"/>	1913010 <b>Fixed Rear Axles</b>

Rear Tires			
<input type="checkbox"/>	1914000	<b>Urethane</b>	
<input type="checkbox"/>	1914010	<b>Pneumatic</b>	\$ 60
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.			
<input type="checkbox"/>	1914020	<b>Pneumatic with Flat-Free Inserts</b> <i>E2213</i>	\$ 95
Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.			

Hand Rims			
<input type="checkbox"/>	1915000	<b>No Hand Rims</b>	
<input type="checkbox"/>	1915010	<b>Standard</b> 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1915020	<b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only.	\$ 95

Wheel Locks			
<input type="checkbox"/>	1916010	<b>Push to Lock</b>	
<input type="checkbox"/>	1916020	<b>Pull to Lock</b>	
<input type="checkbox"/>	1916030	<b>Attendant Wheel Lock (Push to Lock)</b>	\$ 85
Not compatible with Wheel Lock Extensions.			
<input type="checkbox"/>	1916040	<b>User &amp; Attendant Wheel Lock</b>	\$ 125
Available only on 20", 22", or 24" Rear Wheels.			
<input type="checkbox"/>	1916100	<b>Wheel Lock Extensions</b> <i>E0961</i>	\$ 60
Not compatible with Attendant Wheel Lock.			
<input type="checkbox"/>	1916200	<b>Hill Holders</b>	\$ 350
Requires Pneumatic Rear Tires with Flat-Free Inserts.			

Anti-Tippers			
<input checked="" type="checkbox"/>	1917010	<b>Rear Anti-Tippers</b> <i>E0971</i>	\$ 60
Required Safety Option.			
<input type="checkbox"/>	1917020	<b>Front Anti-Tippers</b>	\$ 300
<input type="checkbox"/>	1917030	<b>Heavy Duty Rear Anti-Tipper</b>	\$ 600
Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.			

Back Post Style			
<input type="checkbox"/>	1918010	<b>Push Bar Back</b>	
<input type="checkbox"/>	1918020	<b>Cane Handle Back</b>	
Not compatible with Headrest Mounting Plate.			
<b>Back Post Height</b>			
Seat Pan to Handle			
<input type="checkbox"/>	1918425	25" Height	
<input type="checkbox"/>	1918431	31" Height	

Back Options			
<input type="checkbox"/>		<b>Angle/Height Adjustable Push Bar Handle</b>	\$ 150
Not compatible with Headrest Mounting Plate or Cane Handle Back.			
<input type="checkbox"/>	1919010	Button on Outside	
<input type="checkbox"/>	1919011	Button on Inside	

Seat to Back Angle			
<input type="checkbox"/>	1919090	Set at 90°	
<input type="checkbox"/>	1919098	Set at 98°	
<input type="checkbox"/>	1919106	Set at 106°	
<input type="checkbox"/>	1919114	Set at 114°	
<input type="checkbox"/>	1919122	Set at 122°	
<input type="checkbox"/>	1919210	<b>Backrest Extension 10" high</b>	\$ 120
Not compatible with Angle/Height Adjustable Push Bar Handle.			
<input type="checkbox"/>	1919220	<b>Fold Down Back</b>	
Not compatible with Reclining Back.			
<input type="checkbox"/>	1919230	<b>Spring Back Suspension</b>	\$ 105
Not compatible with Reclining Back.			
<input type="checkbox"/>	1919240	<b>Dynamic Rocker Back Assembly</b> <i>K0108</i>	\$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller. Requires longer Frame.			

Back Upholstery			
<input type="checkbox"/>	1920000	<b>No Back Upholstery</b>	
<input type="checkbox"/>		<b>Standard</b>	\$ 125
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920021	Set at 21"	
<input type="checkbox"/>	1920022	Set at 22"	
<input type="checkbox"/>		<b>Short (25" Back only)</b>	\$ 125
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920118	Set at 18"	
<input type="checkbox"/>	1920119	Set at 19"	
<input type="checkbox"/>	1920120	Set at 20"	
<input type="checkbox"/>		<b>Tall (31" Back only)</b>	\$ 125
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920224	Set at 24"	
<input type="checkbox"/>	1920225	Set at 25"	
<input type="checkbox"/>	1920226	Set at 26"	
<input type="checkbox"/>		<b>Adjustable Tension</b>	\$ 250
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920321	Set at 21"	
<input type="checkbox"/>	1920322	Set at 22"	
<input type="checkbox"/>	1920323	Set at 23"	

Backrest			
<input type="checkbox"/>	1921000	<b>No Backrest</b>	
<b>next Optima HA</b>			
<input type="checkbox"/>	1921100	AC - 16"-20" width, 16"-20" height	<i>E2620</i> \$ 690
<input type="checkbox"/>	1921101	PC - 16"-20" width, 16"-20" height	<i>E2620</i> \$ 690
<input type="checkbox"/>	1921105	AC - 22"-24" width, 16"-20" height	<i>E2621</i> \$ 720
<input type="checkbox"/>	1921106	PC - 22"-24" width, 16"-20" height	<i>E2621</i> \$ 720
<b>Matrx Elite</b>			
<input type="checkbox"/>	1921110	14"-20" width, 10"-20" height	<i>E2620</i> \$ 825
<input type="checkbox"/>	1921112	HD - 22"-26" width, 16"-20" height	<i>E2621</i> \$ 1,095
<input type="checkbox"/>	1921116	Deep - 14"-20" width, 12"-20" height	<i>E2620</i> \$ 950

See next page for more Backrest Options →

ROHO Agility			
<input type="checkbox"/>	1921120	Min - 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1921122	Mid - 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1921124	Max - 14"-20" width, 10"-18" height	E2620 \$ 950
Varilite Icon			
<input type="checkbox"/>	1921130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1921132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
Backrest Size			
	1921400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1923011	Installed on Right Back Post	
<input type="checkbox"/>	1923012	Installed on Left Back Post	
Tilt Limit			
<input type="checkbox"/>	1923220	Set at 20°	
<input type="checkbox"/>	1923230	Set at 30°	
<input type="checkbox"/>	1923240	Set at 40°	
<input type="checkbox"/>	1923410	<b>Reclining Backrest (Manual)</b> E1225	\$ 830
Maximum recline = 30°. Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, 20"-22" seat depth. Requires longer Frame.			

Headrest			
<input type="checkbox"/>	1924000	<b>No Headrest</b>	
<input type="checkbox"/>	1924010	<b>Headrest Mounting Plate</b>	\$ 75
Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	1924100	<b>nxt MultiFit Head Support Hardware</b> E1028	\$ 155
<input type="checkbox"/>		<b>nxt MultiFit Head Support Pad</b> E0955	\$ 195
<input type="checkbox"/>	1924105	Regular (9" x 5")	
<input type="checkbox"/>	1924106	Large (10" x 5")	
<input type="checkbox"/>	1924107	X-Large (14" x 5")	
Stealth Headrest Hardware E1028			
<input type="checkbox"/>	1924110	Multi-Axis HMO475	\$ 150
<input type="checkbox"/>	1924111	Link TWB480	\$ 325
<input type="checkbox"/>		<b>Stealth Headrest Pad</b> E0955	\$ 200
<input type="checkbox"/>	1924115	10" width	
<input type="checkbox"/>	1924116	14" width	

Front Rigging			
<input type="checkbox"/>	1925000	<b>No Front Rigging</b>	
<input type="checkbox"/>	1925001	<b>Sector Blocks only</b>	
<input type="checkbox"/>	1925110	<b>Plug In 70°</b>	
<input type="checkbox"/>	1925120	<b>Plug In Adjustable Knee Angle</b>	\$ 975
<input type="checkbox"/>	1925130	<b>Swing Away 60°</b>	
<input type="checkbox"/>	1925140	<b>Swing Away 70°</b>	
<input type="checkbox"/>	1925200	<b>Swing Away Elevating Legrest</b> E0990	\$ 355

Foot Plates			
<input type="checkbox"/>	1926010	<b>Composite</b>	
<input type="checkbox"/>	1926020	<b>Angle Adjustable</b> K0040	\$ 140
<input type="checkbox"/>	1926030	<b>One Piece Fixed</b>	\$ 350
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1926040	<b>One Piece Flip-up</b>	\$ 395
<input type="checkbox"/>	1926050	<b>Interlocking Angle Adjustable</b>	\$ 195
Seat Widths of 14"-21" only.			

Footrest/Legrest Accessories			
<input type="checkbox"/>	1927010	<b>Swing Away 60° &amp; 70° Calf Pad</b>	\$ 135
Not compatible with 60° & 70° Gel Pad.			
<input type="checkbox"/>	1927020	<b>Swing Away 60° &amp; 70° Gel Pad</b>	\$ 100
Not compatible with 60° & 70° Swing Away Calf Pad.			
<input type="checkbox"/>	1927030	<b>Swing Away ELR Gel Knee Pad</b>	\$ 100
<input type="checkbox"/>	1927040	<b>Adjustable Length Heel Loops</b> E0951	\$ 50
<input type="checkbox"/>	1927050	<b>Bodypoint - Padded Calf Strap</b> K0038	\$ 105
Miller's Dynamic Footrest Mounts			
<input type="checkbox"/>	1927060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1927065	For Swing Away 70° Front Rigging	\$ 380

Extension Tubes			
Regular Extension Tubes			
<input type="checkbox"/>	1928010	Standard (13"-17")	
<input type="checkbox"/>	1928020	Long (15"-19")	\$ 75
Extra Strength Extension Tubes (Solid)			
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1928110	Standard (13"-17")	
<input type="checkbox"/>	1928120	Long (15"-19")	

Positioning Straps			
Bodypoint - Hip Belt E0978			
<input type="checkbox"/>	1929010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1929015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1929020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1929030	<b>Bodypoint - Chest Strap</b> E0960	\$ 85

Miscellaneous Options			
<input type="checkbox"/>	1930110	<b>IV Pole</b> E1029	\$ 295
Not compatible with HD Anti Tipper, Attendant Wheel Locks.			
<input type="checkbox"/>	1930120	<b>O2 Holder</b> E2208	\$ 295
Not compatible with HD Anti Tipper, Attendant Wheel Locks.			
<input type="checkbox"/>	1930170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	\$ 10
For Cushion mounting.			