

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type			
<input checked="" type="checkbox"/>	1901000 Stellar GL 45° Manual Tilt Wheelchair	E1161	\$ 2,995

Frame Options	
<input checked="" type="checkbox"/>	1902000 Standard Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/>	1902300 Transit Tie Down System (TTDS) TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width			
Standard Seat Widths			
<input type="checkbox"/>	1903014	14" Seat Width	
<input type="checkbox"/>	1903015	15" Seat Width	
<input type="checkbox"/>	1903016	16" Seat Width	
<input type="checkbox"/>	1903017	17" Seat Width	
<input type="checkbox"/>	1903018	18" Seat Width	
<input type="checkbox"/>	1903019	19" Seat Width	
<input type="checkbox"/>	1903020	20" Seat Width	E2201 \$ 325
<input type="checkbox"/>	1903021	21" Seat Width	E2201 \$ 325
<input type="checkbox"/>	1903022	22" Seat Width	E2201 \$ 325

Seat Depth			
Standard Seat Depths (Adjustable)			
<input type="checkbox"/>	1904016	16" Seat Depth	
<input type="checkbox"/>	1904017	17" Seat Depth	
<input type="checkbox"/>	1904018	18" Seat Depth	
<input type="checkbox"/>	1904019	19" Seat Depth	
<input type="checkbox"/>	1904020	20" Seat Depth	E2203 \$ 475
Custom Seat Depths			
<input type="checkbox"/>	1904021	21" Seat Depth	E2203 \$ 475
<input type="checkbox"/>	1904022	22" Seat Depth	E2204 \$ 750

Frame Color	
PDG Primary Colors	
<input type="checkbox"/>	1905000 Tracer Red
<input type="checkbox"/>	1905001 Gloss Black
<input type="checkbox"/>	1905002 Silver
<input type="checkbox"/>	1905003 Black Shimmer
<input type="checkbox"/>	1905004 Orange Fusion
<input type="checkbox"/>	1905005 Karl's Imagination
<input type="checkbox"/>	1905006 Blue Shimmer
<input type="checkbox"/>	1905007 Gloss White
<input type="checkbox"/>	1905008 Dark Titanium
PDG Premium Colors	
<input type="checkbox"/>	1905009 Geyser Green \$ 150
<input type="checkbox"/>	1905010 Lollipop Blue
<input type="checkbox"/>	1905011 Lollipop Red
Custom Color	
<input type="checkbox"/>	1905200 Custom Color \$ 300 Contact PDG for specifications.

Arm Type	
<input type="checkbox"/>	1906000 No Armrest
<input type="checkbox"/>	1906100 Single Post Fixed Height (12") Full Length Armpad only.
Single Post Adjustable Height	
	E0973 \$ 255
<input type="checkbox"/>	1906111 8"-11"
<input type="checkbox"/>	1906114 10"-14"
<input type="checkbox"/>	1906200 Dual Post Fixed Height (12") Full Length Armpad only.
Dual Post Adjustable Height	
	E0973 \$ 255
<input type="checkbox"/>	1906211 8"-11"
<input type="checkbox"/>	1906214 10"-14"
<input type="checkbox"/>	1906216 12"-16"
Locking Flip-up	
	E0973 \$ 255
Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.	
<input type="checkbox"/>	1906310 Set at 10"
<input type="checkbox"/>	1906311 Set at 11"
<input type="checkbox"/>	1906312 Set at 12"
<input type="checkbox"/>	1906313 Set at 13"
<input type="checkbox"/>	1906314 Set at 14"

PDG Product Design Group Inc.
318 East Kent Avenue South
Vancouver BC V5X 4N6 Canada
Phone 604-323-9220 | 888-858-4422

Armrests

Regular Armrests

- 1907100 Desk Length - 2" x 10" (pair)
- 1907101 Full Length - 2" x 14" (pair)

Wide Armrests

\$ 115

May interfere with wheels in some tilt configurations.

- 1907110 Desk Length - 3" x 10" (pair)
- 1907111 Full Length - 3" x 12" (pair)

Gel Armrests

\$ 115

May interfere with wheels in some tilt configurations.

- 1907120 Full Length - 2" x 12" (pair)
- 1907121 Full Length - 3.5" x 12" (pair)

Contour Arm Trough *E2209*

\$ 125

May interfere with wheels in some tilt configurations.

- 1907130 Right (4.5" x 18") with Hand Pad
- 1907131 Left (4.5" x 18") with Hand Pad

Armrest Wideners for Dual Post Armrests

\$ 375

Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels or Clothing Guards.

Clothing Guards *K0108*

Not compatible with Flip-up Armrests.

- 1907310 For Single Post Armrests
10"-14" Height range only \$ 110
- 1907320 For Dual Post Armrests \$ 165

Armrest Extensions

\$ 60

Not compatible with Flip-up or Single Post Armrests.

- 1907400 3" Extension - Forward
- 1907401 3" Extension - Rearward

Seat Pan

- 1908000 **No Seat Pan**
- 1908001 **Solid Seat Pan** *E2231* \$ 120

Cushion

- 1909000 **No Cushion**
- next BioFit**
 - 1909101 14"-20" width, 16"-20" depth *E2607* \$ 550
 - 1909102 22"-24" width, 16"-20" depth *E2608* \$ 600
- Matrx VI**
 - 1909111 14"-20" width, 16"-22" depth *E2607* \$ 550
 - 1909112 22" width, 18"-22" depth *E2608* \$ 600
 - 1909115 HD - 18"-30" width, 18"-22" depth *E2608* \$ 1,150
- ROHO Quadtro Select HP**
 - 1909121 14"-20" width, 16"-22" depth *E2624* \$ 650
 - 1909122 22"-24" width, 16"-22" depth *E2625* \$ 750
- ROHO High Profile**
 - 1909125 14"-20" width, 16"-22" depth *E2622* \$ 650
 - 1909126 22"-24" width, 16"-22" depth *E2623* \$ 750
- Varilite Reflex**
 - 1909131 14"-20" width, 16"-20" depth *E2603* \$ 375

Cushion Size

1909400 Specify Size (Width" x Depth") _____



Scan this code or visit www.pdgmobility.com to learn more about the Stellar GL Manual Tilt Wheelchair.

Seat Height

Circle desired Seat-to-Floor Height or enter the information below:

RW	6" Casters					
12"	13"	14"	15"	16"		
16"			15"	16"	17"	18"
20"	13"	14"	15"	16"		
22"		14"	15"	16"	17"	
24"			15"	16"	17"	18"

RW	8" Casters				
16"		15"	16"	17"	18"
24"		15"	16"	17"	18"

Wheel size _____ Caster size _____ STF _____

STF Height is adjustable within the range of the selected wheel size.

Casters

Urethane

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1911061 6" x 1.25" Caster size
- 1911063 6" x 1.50" Caster size \$ 85
- 1911081 8" x 1.25" Caster size
- 1911084 8" x 2.00" Caster size \$ 135

Pneumatic

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 1911161 6" x 1.25" Caster size \$ 110
- 1911181 8" x 1.25" Caster size \$ 120
- 1911184 8" x 2.00" Caster size \$ 150

- 1911260 **6" Soft Roll** \$ 95

- 1911400 **Ultra Sport Fork** *E1015* \$ 395

Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.

Rear Wheels

Composite

- 1912012 12" Wheel size
- 1912016 16" Wheel size
- 1912020 20" Wheel size
- 1912022 22" Wheel size
- 1912024 24" Wheel size

Spoke

\$ 75

- 1912120 20" Wheel size
- 1912122 22" Wheel size
- 1912124 24" Wheel size

- 1912224 **24" Heavy Duty Spoke** \$ 350

Axle Type

- 1913000 **Quick Release Rear Axles**
- 1913010 **Fixed Rear Axles**

Rear Tires			
<input type="checkbox"/>	1914000	Urethane	
<input type="checkbox"/>	1914010	Pneumatic	\$ 60
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.			
<input type="checkbox"/>	1914020	Pneumatic with Flat-Free Inserts <i>E2213</i>	\$ 95
Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.			

Hand Rims			
<input type="checkbox"/>	1915000	No Hand Rims	
<input type="checkbox"/>	1915010	Standard 20", 22", and 24" Rear Wheels only.	
<input type="checkbox"/>	1915020	Plastic Coated 20", 22", and 24" Rear Wheels only.	\$ 95

Wheel Locks			
<input type="checkbox"/>	1916010	Push to Lock	
<input type="checkbox"/>	1916020	Pull to Lock	
<input type="checkbox"/>	1916030	Attendant Wheel Lock (Push to Lock)	\$ 85
Not compatible with Wheel Lock Extensions.			
<input type="checkbox"/>	1916040	User & Attendant Wheel Lock	\$ 125
Available only on 20", 22", or 24" Rear Wheels.			
<input type="checkbox"/>	1916100	Wheel Lock Extensions <i>E0961</i>	\$ 60
Not compatible with Attendant Wheel Lock.			
<input type="checkbox"/>	1916200	Hill Holders	\$ 350
Requires Pneumatic Rear Tires with Flat-Free Inserts.			

Anti-Tippers			
<input checked="" type="checkbox"/>	1917010	Rear Anti-Tippers <i>E0971</i>	\$ 60
Required Safety Option.			
<input type="checkbox"/>	1917020	Front Anti-Tippers	\$ 300
<input type="checkbox"/>	1917030	Heavy Duty Rear Anti-Tipper	\$ 600
Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.			

Back Post Style			
<input type="checkbox"/>	1918010	Push Bar Back	
<input type="checkbox"/>	1918020	Cane Handle Back	
Not compatible with Headrest Mounting Plate.			
Back Post Height			
Seat Pan to Handle			
<input type="checkbox"/>	1918425	25" Height	
<input type="checkbox"/>	1918431	31" Height	

Back Options			
		Angle/Height Adjustable Push Bar Handle	\$ 150
Not compatible with Headrest Mounting Plate or Cane Handle Back.			
<input type="checkbox"/>	1919010	Button on Outside	
<input type="checkbox"/>	1919011	Button on Inside	

Seat to Back Angle			
<input type="checkbox"/>	1919090	Set at 90°	
<input type="checkbox"/>	1919098	Set at 98°	
<input type="checkbox"/>	1919106	Set at 106°	
<input type="checkbox"/>	1919114	Set at 114°	
<input type="checkbox"/>	1919122	Set at 122°	
<input type="checkbox"/>	1919210	Backrest Extension 10" high	\$ 120
Not compatible with Angle/Height Adjustable Push Bar Handle.			
<input type="checkbox"/>	1919220	Fold Down Back	
Not compatible with Reclining Back.			
<input type="checkbox"/>	1919230	Spring Back Suspension	\$ 105
Not compatible with Reclining Back.			
<input type="checkbox"/>	1919240	Dynamic Rocker Back Assembly <i>K0108</i>	\$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller. Requires longer Frame.			

Back Upholstery			
<input type="checkbox"/>	1920000	No Back Upholstery	
		Standard	\$ 125
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920021	Set at 21"	
<input type="checkbox"/>	1920022	Set at 22"	
		Short (25" Back only)	\$ 125
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920118	Set at 18"	
<input type="checkbox"/>	1920119	Set at 19"	
<input type="checkbox"/>	1920120	Set at 20"	
		Tall (31" Back only)	\$ 125
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920224	Set at 24"	
<input type="checkbox"/>	1920225	Set at 25"	
<input type="checkbox"/>	1920226	Set at 26"	
		Adjustable Tension	\$ 250
Back height measured from Seat Pan to top of Upholstery			
<input type="checkbox"/>	1920321	Set at 21"	
<input type="checkbox"/>	1920322	Set at 22"	
<input type="checkbox"/>	1920323	Set at 23"	

Backrest			
<input type="checkbox"/>	1921000	No Backrest	
next Optima HA			
<input type="checkbox"/>	1921100	AC - 16"-20" width, 16"-20" height	<i>E2620</i> \$ 690
<input type="checkbox"/>	1921101	PC - 16"-20" width, 16"-20" height	<i>E2620</i> \$ 690
<input type="checkbox"/>	1921105	AC - 22"-24" width, 16"-20" height	<i>E2621</i> \$ 720
<input type="checkbox"/>	1921106	PC - 22"-24" width, 16"-20" height	<i>E2621</i> \$ 720
Matrx Elite			
<input type="checkbox"/>	1921110	14"-20" width, 10"-20" height	<i>E2620</i> \$ 825
<input type="checkbox"/>	1921112	HD - 22"-26" width, 16"-20" height	<i>E2621</i> \$ 1,095
<input type="checkbox"/>	1921116	Deep - 14"-20" width, 12"-20" height	<i>E2620</i> \$ 950

See next page for more Backrest Options →

ROHO Agility			
<input type="checkbox"/>	1921120	Min - 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1921122	Mid - 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1921124	Max - 14"-20" width, 10"-18" height	E2620 \$ 950
Varilite Icon			
<input type="checkbox"/>	1921130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1921132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
Backrest Size			
	1921400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1923011	Installed on Right Back Post	
<input type="checkbox"/>	1923012	Installed on Left Back Post	
Tilt Limit			
<input type="checkbox"/>	1923220	Set at 20°	
<input type="checkbox"/>	1923230	Set at 30°	
<input type="checkbox"/>	1923240	Set at 40°	
<input type="checkbox"/>	1923410	Reclining Backrest (Manual) E1225	\$ 830
Maximum recline = 30°. Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, 20"-22" seat depth. Requires longer Frame.			

Headrest			
<input type="checkbox"/>	1924000	No Headrest	
<input type="checkbox"/>	1924010	Headrest Mounting Plate	\$ 75
Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	1924100	nxt MultiFit Head Support Hardware E1028	\$ 155
<input type="checkbox"/>		nxt MultiFit Head Support Pad E0955	\$ 195
<input type="checkbox"/>	1924105	Regular (9" x 5")	
<input type="checkbox"/>	1924106	Large (10" x 5")	
<input type="checkbox"/>	1924107	X-Large (14" x 5")	
Stealth Headrest Hardware E1028			
<input type="checkbox"/>	1924110	Multi-Axis HMO475	\$ 150
<input type="checkbox"/>	1924111	Link TWB480	\$ 325
<input type="checkbox"/>		Stealth Headrest Pad E0955	\$ 200
<input type="checkbox"/>	1924115	10" width	
<input type="checkbox"/>	1924116	14" width	

Front Rigging			
<input type="checkbox"/>	1925000	No Front Rigging	
<input type="checkbox"/>	1925001	Sector Blocks only	
<input type="checkbox"/>	1925110	Plug In 70°	
<input type="checkbox"/>	1925120	Plug In Adjustable Knee Angle	\$ 975
<input type="checkbox"/>	1925130	Swing Away 60°	
<input type="checkbox"/>	1925140	Swing Away 70°	
<input type="checkbox"/>	1925200	Swing Away Elevating Legrest E0990	\$ 355

Foot Plates			
<input type="checkbox"/>	1926010	Composite	
<input type="checkbox"/>	1926020	Angle Adjustable K0040	\$ 140
<input type="checkbox"/>	1926030	One Piece Fixed	\$ 350
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1926040	One Piece Flip-up	\$ 395
<input type="checkbox"/>	1926050	Interlocking Angle Adjustable K0040	\$ 195
Seat Widths of 14"-21" only.			

Footrest/Legrest Accessories			
<input type="checkbox"/>	1927010	Swing Away 60° & 70° Calf Pad	\$ 135
Not compatible with 60° & 70° Gel Pad.			
<input type="checkbox"/>	1927020	Swing Away 60° & 70° Gel Pad	\$ 100
Not compatible with 60° & 70° Swing Away Calf Pad.			
<input type="checkbox"/>	1927030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1927040	Adjustable Length Heel Loops E0951	\$ 50
<input type="checkbox"/>	1927050	Bodypoint - Padded Calf Strap K0038	\$ 105
Miller's Dynamic Footrest Mounts			
<input type="checkbox"/>	1927060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1927065	For Swing Away 70° Front Rigging	\$ 380

Extension Tubes			
Regular Extension Tubes			
<input type="checkbox"/>	1928010	Standard (13"-17")	
<input type="checkbox"/>	1928020	Long (15"-19")	\$ 75
Extra Strength Extension Tubes (Solid)			
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1928110	Standard (13"-17")	
<input type="checkbox"/>	1928120	Long (15"-19")	

Positioning Straps			
Bodypoint - Hip Belt E0978			
<input type="checkbox"/>	1929010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1929015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1929020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1929030	Bodypoint - Chest Strap E0960	\$ 85

Miscellaneous Options			
<input type="checkbox"/>	1930110	IV Pole E1029	\$ 295
Not compatible with HD Anti Tipper, Attendant Wheel Locks.			
<input type="checkbox"/>	1930120	O2 Holder E2208	\$ 295
Not compatible with HD Anti Tipper, Attendant Wheel Locks.			
<input type="checkbox"/>	1930170	Adhesive Back Hook Strip, 2" x 30"	\$ 10
For Cushion mounting.			